



# THE MEMPHIS DEPOT TENNESSEE

---

## ADMINISTRATIVE RECORD COVER SHEET

AR File Number 877

*Part II of III*

Controlled Copy

Copy No: \_\_\_\_\_

Implementation Date 10/17/03

SOP No: CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03Page 1 of 84**STL NORTH CANTON STANDARD OPERATING PROCEDURE**

**TITLE: GAS CHROMATOGRAPHIC ANALYSIS BASED ON METHOD 8000B,  
8021B, 8081A, 8082, 8151A, 8310, 8141A, 8015B, and 615**

(SUPERSEDES: Revision 5.6 Dated 05/25/01)

Prepared by:	<u><i>Amey Dayal</i></u>	<u>10/9/03</u>
		Date
Approved by:	<u><i>Marcel Mol</i></u>	<u>10/9/03</u>
	Technology Specialist	Date
Approved by:	<u><i>Beth Lambert</i></u>	<u>10/17/03</u>
	Quality Assurance Manager	Date
Approved by:	<u><i>Deborah Dunn</i></u>	<u>10/16/03</u>
	Environmental Health and Safety	Date
Approved by:	<u><i>Carol Ann Felt</i></u>	<u>10/16/03</u>
	Laboratory Director	Date
Approved by:	<u><i>Mark Burre</i></u>	<u>10/17/03</u>
	Technical Director	Date

## Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions. THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

**TABLE OF CONTENTS**

1.	Scope and Application	4
2.	Summary of Method	4
3.	Definitions	4
4.	Interferences	4
5.	Safety	4
6.	Equipment and Supplies	5
7.	Reagents and Standards	5
8.	Sample Preservation and Storage	5
9.	Quality Control	6
10.	Calibration and Standardization	9
11.	Procedure	14
12.	Data Analysis and Calculations	16
13.	Method Performance	19
14.	Pollution Prevention	20
15.	Waste Management	20
16.	References	20
17.	Miscellaneous	20

**List of Appendices:**

Appendix A	Analysis of Volatile Organics based on Method 8021B
Appendix B	Analysis of Organochlorine Pesticides based on Method 8081A
Appendix C	Analysis of Organochlorine Pesticides and PCBs based on Method 8082
Appendix D	Analysis of Phenoxy Acid Herbicides based on Method 8151A
Appendix E	Analysis of Organophosphorus Pesticides Based on Method 8141A
Appendix F	Total Petroleum Hydrocarbons Based on 8015B
Appendix G	Analysis of Non-Halogenated Organic Compounds by Method 8015B, Direct Injection
Appendix H	Analysis of Phillips 66 Compounds

**List of Tables**

Table A1	Standard analyte lists for 8021B
Table A2-A4	Recommended conditions for method 8021B
Table A5	Surrogate and Internal standard concentrations for aqueous and low level soil samples, method 8021B
Table A6	Concentrations for LCS and MS/MSD compounds, low level soil and aqueous, method 8021B.
Table B1	Standard analyte and reporting limits, method 8081A
Table B2	Recommended conditions, method 8081A
Table B3	Calibration levels, method 8081A
Table B4	Column degradation evaluation mix, method 8081A
Table B5	LCS/matrix spike and surrogate levels, method 8081A
Table B6	LCS/matrix spike and surrogate levels for TCLP, method 8081A
Table B7	Suggested analytical sequence, method 8081A
Table B8	Performance limits, method 8081A
Table C1	Standard analyte list, method 8082
Table C2	Instrumental conditions, method 8082
Table C3	Calibration standards, method 8082
Table C4	LCS/Matrix spike and surrogate levels, method 8082
Table C5	Suggested Analytical Sequence, Method 8082
Table D1	Standard Analyte List, method 8151A
Table D2	Instrumental conditions, method 8151A
Table D3	LCS/Matrix spike and surrogate levels, method 8151B
Table D4	Performance limits, method 8151B
Table E1	Analyte and Reporting Limit List, Method 8141A
Table E2	Instrument Conditions, Method 8141A
Table E3	Initial Calibration, Method 8141A
Table E4	LCS/Matrix spike and surrogate compounds, Method 8141A
Table F1	Suggested GC Temperature Program for TPH Analysis
Table F2	Reporting Limits for TPH Analysis
Table G1	Non-Halogenated Organic Compounds Reporting Limits
Table G2	Non-Halogenated Organic Compounds Working Standards
Table H1	Phillips 66 Compounds Reporting Limits



**1. SCOPE AND APPLICATION**

This SOP describes procedures for analysis of organic analytes by Gas Chromatography (GC). The procedures are based on SW-846 methodology and are applicable for measurements made to comply with the Resource Conservation and Recovery Act (RCRA). Individual analytes and methods are described in the appendices. Appendix G describes procedures for the analysis of petroleum hydrocarbons by SW-846 8015B methodology. Appendix H includes criteria for the analysis of non-halogenated organic compounds by Method 8015B, Direct Injection. Appendix I describes the analysis of Phillips 66 analytes by Method 8015B.

**2. SUMMARY OF METHOD**

In general, semivolatile analytes in aqueous samples are prepared for analysis using continuous or separatory funnel liquid / liquid extraction or solid phase extraction (SOP # CORP-OP-0001NCNC) Solid samples are prepared using sonication, Soxhlet or pressurized fluid extraction (SOP # CORP-OP-0001NC). Volatile analytes are prepared for analysis using purge and trap methodology (Appendix A).

After the initial preparation step, the sample is introduced to the GC and concentrations of target analytes are measured by the detector response within a defined retention time window, relative to the response to standard concentrations. Internal or external standardization procedures are used as specified in the method appendices.

**3. DEFINITIONS**

Definitions of terms used in this SOP may be found in the glossary of the STL North Canton Laboratory Quality Manual (LQM), current version.

**4. INTERFERENCES**

Contamination by carryover can occur when a low concentration sample is analyzed after a high concentration sample. In addition, some purge and trap autosamplers are susceptible to port specific contamination. Co-elution of target analytes with non-targets can occur, resulting in false positives or biased high results. In particular, this is a problem with non-selective detectors such as the Flame Ionization Detector (FID). See the appendices for interferences specific to individual tests and suggested corrective actions

**5. SAFETY**

- 5.1. Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual and this document.
- 5.2. Eye protection that prevents splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have become contaminated will be removed and discarded; other gloves will be cleaned immediately. Refer to the STL North Canton Chemical Hygiene plan for a complete description of personal protection equipment. Latex, Nitrile and vinyl gloves all provide adequate protection against the methanol used in this method.
- 5.3. The following is a list of the materials used in this method, which have a serious or significant hazard rating **NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table.** A complete list of materials used in the method can be found in the reagents and materials section.

Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

Material	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
Methanol	Flammable Poison Irritant	200 ppm-TWA	A slight irritant to the mucous membranes. Toxic effects exerted upon nervous system, particularly the optic nerve. Symptoms of overexposure may include headache, drowsiness and dizziness. Methyl alcohol is a defatting agent and may cause skin to become dry and cracked. Skin absorption can occur; symptoms may parallel inhalation exposure. Irritant to the eyes.
Acetone	Flammable	1000 ppm-TWA	Inhalation of vapors irritates the respiratory tract. May cause coughing, dizziness, dullness, and headache.
Hexane	Flammable Irritant	500 ppm-TWA	Inhalation of vapors irritates the respiratory tract. Overexposure may cause lightheadedness, nausea, headache, and blurred vision. Vapors may cause irritation to the skin and eyes.
Methylene Chloride	Carcinogen Irritant	25 ppm-TWA 125 ppm-STEL	Causes irritation to respiratory tract. Has a strong narcotic effect with symptoms of mental confusion, light-headedness, fatigue, nausea, vomiting and headache. Causes irritation, redness and pain to the skin and eyes. Prolonged contact can cause burns. Liquid degrades the skin. May be absorbed through skin.
Sulfuric Acid	Corrosive Oxidizer Dehydrator Poison Carcinogen	1 Mg/M3-TWA	Inhalation produces damaging effects on the mucous membranes and upper respiratory tract. Symptoms may include irritation of the nose and throat, and labored breathing. Symptoms of redness, pain, and severe burn can occur. Contact can cause blurred vision, redness, pain and severe tissue burns. Can cause blindness.

MTBE	Flammable Irritant	None established	Excessive exposure may cause irritation to the nose, throat, lungs, and respiratory tract. Central nervous system effects may include headache, dizziness, loss of balance and coordination.
THF	Flammable Irritant Possible Carcinogen	200 ppm TWA	Causes skin and eye irritation. Irritating to mucous membranes and upper respiratory tract.
Isooctane	Flammable Irritant	None Established	May cause eye, respiratory tract and skin irritation.
Hydrochloric Acid	Corrosive Poison	5 ppm-Ceiling	Inhalation of vapors can cause coughing, choking, inflammation of the nose, throat, and upper respiratory tract, and in severe cases, pulmonary edema, circulatory failure, and death. Can cause redness, pain, and severe skin burns. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
Sodium Bisulfate	Corrosive	None Listed	Contact may cause skin/eye burns. Inhalation can cause irritation of the respiratory tract with burning pain in the nose and throat, coughing, wheezing and shortness of breath. Causes chemical burns to the respiratory tract. May cause fatal spasms, inflammation or pulmonary/respiratory edema.
1 – Always add acid to water to prevent violent reactions.			
limit refers to the OSHA regulatory exposure limit.			

- 5.4. Opened containers of neat standards will be handled in a fume hood.
- 5.5. Sample extracts and standards, which are in a flammable solvent, shall be stored in an explosion-proof refrigerator
- 5.6. When using hydrogen gas as a carrier, all precautions listed in the CSM shall be observed.
- 5.7. Standard preparation and dilution shall be performed inside an operating fume hood.
- 5.8. The gas chromatograph contains zones that have elevated temperatures. The analyst needs to be aware of the locations of those zones, and must cool them to room temperature prior to working on them.
- 5.9. There are areas of high voltage in both the gas chromatograph and the mass spectrometer. Depending on the type of work involved, either turn the power to the instrument off, or disconnect it from its source of power.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. An analytical system complete with a gas chromatograph is required. A data system capable of measuring peak area and/or height is required. Recommended equipment and supplies for individual methods are listed in each method appendix.

**7. REAGENTS AND STANDARDS****7.1. Stock Standards**

Stock standards are purchased as certified solutions or prepared from pure solutions. Stock standards for method 8021B are stored at -10 to -20°C. Other stock standard solutions are stored as recommended by the manufacturer. All stock standards must be protected from light. Stock standard solutions should be brought to room temperature before using.

Semivolatile stock standard solutions must be replaced after one year. Stock standards of gases must be replaced at least every week, unless the acceptability of the standard is demonstrated (Less than 20% drift from the initial calibration is an acceptable demonstration). Other volatile stock standards must be replaced every 6 months or sooner if comparison with check standards prepared from an independent source indicates a problem.

- 7.1.1. Expiration times for all standards are measured from the time the standard is prepared or from the time that the standard ampoule is opened, if the standard is supplied in a sealed ampoule. If a vendor-supplied standard has an earlier expiration date then that date is used.

**7.2. Calibration Standards****7.2.1. Volatile Calibration Standards**

The procedure for preparation of volatile standards is given in Appendix A.

**7.2.2 Semivolatile Calibration Standards**

Semivolatile calibration standards are prepared as dilutions of the stock standards. Surrogates and internal standards are used as specified in the method appendices. Semivolatile calibration solutions must be refrigerated at  $\leq 6^{\circ}\text{C}$  and protected from light. The standards must be replaced at least every six months or sooner if comparison with check standards indicates a problem.

- 7.3. Gases for carrier and make-up: Hydrogen, Helium, Nitrogen, Argon/Methane.

**7.4. Quality control (QC) Standards**

QC standards (matrix spiking and LCS standards) are prepared and stored in the same way as calibration standards. They must be made from a stock independent from the calibration standards.

**8. SAMPLE PRESERVATION AND STORAGE**

Semivolatile extracts must be refrigerated at  $\leq 6^{\circ}\text{C}$  and analyzed within 40 days of the end of the extraction. Volatile sample storage conditions and holding times are given in Appendix A.

**9. QUALITY CONTROL****9.1. Initial Demonstration of Capability**

- 9.1.1. For the standard analyte list, the initial demonstration and method detection limit (MDL) studies described in Section 13 must be acceptable before analysis of samples may begin.

- 9.1.2. For non-standard analytes, a MDL study must be performed and calibration curve generated before analyzing any samples, unless lesser requirements are previously agreed to with the client. In any event the minimum initial demonstration required is analysis of an extracted standard at the reporting limit and a single point calibration.

9.2. Batch Definition

Batches are defined at the sample preparation stage. Batches should be kept together through the whole analytical process as far as possible, but it is not mandatory to analyze prepared extracts on the same instrument or in the same sequence. Refer to the STL North Canton QC Program document (QA-003) for further details of the batch definition.

9.2.1. Quality Control Batch

The batch is a set of up to 20 samples of the same matrix processed using the same procedures and reagents within the same time period. The Quality Control batch must contain a matrix spike / spike duplicate (MS/MSD), a Laboratory Control Sample (LCS), and a method blank. Laboratory generated QC samples (Blank, LCS, MS/MSD) do not count towards the maximum 20 samples in a batch. Field QC samples are included in the batch count. In some cases, at client request, the MS/MSD may be replaced with a matrix spike and sample duplicate. If insufficient sample is available for an MS/MSD a LCSD may be substituted.

9.3. Control Limits

In-house historical control limits must be determined for surrogates, matrix spikes, and laboratory control samples (LCS). These limits must be determined at least annually. The recovery limits are mean recovery  $\pm 3$  standard deviations, unless that limit is tighter than the calibration criteria, in which case limits may be widened. Refer to policy QA-003 for more details

- 9.3.1. These limits do not apply to dilutions (except for tests without a separate extraction), but surrogate and matrix spike recoveries will be reported unless the dilution is more than 5X

- 9.3.2. All surrogate, LCS, and MS recoveries (except for dilutions) must be entered into QuantIMS (when available) or other database so that accurate historical control limits can be generated. For tests without a separate extraction, surrogates and matrix spikes will be reported for all dilutions.

- 9.3.3. Refer to the QC Program document (QA-003) for further details of control limits.

9.4. Surrogates

All methods must use surrogates to the extent possible. Surrogate recoveries in samples and QC samples must be assessed to ensure that recoveries are within established limits. If any surrogates are outside limits, the following corrective actions must take place (except for dilutions):

- Check all calculations for error.
- Ensure that instrument performance is acceptable.
- Recalculate the data and/or reanalyze the extract if either of the above checks reveal a problem.
- Reprepare and reanalyze the sample or flag the data as "Estimated Concentration" if neither of the above resolves the problem. Repreparation is not necessary if there is obvious chromatographic interference.
- The decision to reanalyze or flag the data should be made in consultation with the client. It is only necessary to reprepare / reanalyze a sample once to demonstrate that poor surrogate recovery is due to matrix effect, unless the analyst believes that the repeated out of control results are not due to matrix effect

**GAS CHROMATOGRAPHIC ANALYSIS BASED ON  
METHOD 8000B, SW-846**

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03Page 9 of 86

---

- 9.4.1. If dual column analysis is used the choice of which result to report is made in the same way as for samples (Section 12.1.2) unless one column is out of control, in which case the in-control result is reported.
- 9.4.2. If the surrogates are out of control for the sample, matrix spike, and matrix spike duplicate, then matrix effect has been demonstrated for that sample and reparation is not necessary. If the sample is out of control and the MS and/or MSD is in control, then reparation or flagging of the data is required.
- 9.4.3. Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective actions.

**9.5. Method Blanks**

For each batch of samples, analyze a method blank. The method blank consists of reagent water for aqueous semivolatile samples, and sodium sulfate for semivolatile soils tests (Refer to SOP No. CORP-OP-0001NCNC for details). For low level volatiles, the method blank consists of reagent water. For medium level volatiles, the method blank consists of methanol as described in Appendix A. Surrogates are added and the method blank is carried through the entire analytical procedure. The method blank must not contain any analyte of interest at or above the reporting limit (except common laboratory contaminants, see below) or at or above 5% of the measured concentration of that analyte in the associated samples, whichever is higher.

If the analyte is a common laboratory contaminant (methylene chloride, acetone, 2-butanone, phthalate esters) the data may be reported with qualifiers if the concentration of the analyte is less than five times the reporting limit. Such action must be taken in consultation with the client.

Re-extraction and reanalysis of samples associated with an unacceptable method blank is required when reportable concentrations are determined in the samples.

If there is no target analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. Such action should be taken in consultation with the client.

- 9.5.1. Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective actions.

**9.6. Instrument Blanks**

- 9.6.1. An instrument blank must be analyzed during any 12-hour period of analysis that does not contain a method blank.
- 9.6.2. An instrument blank consists of the appropriate solvent with internal standards added. If internal standards are not used the surrogates should be added.
- 9.6.3. Control criteria are the same as for the method blank, except that only reanalysis of affected samples would be required, not re-extraction.

**9.7. Laboratory Control Samples (LCS)**

For each batch of samples, analyze a LCS. The LCS contains a representative subset of the analytes of interest, and must contain the same analytes as the matrix spike. The LCS may also contain the full set of analytes with a subset of control analytes. If any control analyte or surrogate is outside established control limits, the system is out of control and corrective action must occur. Corrective action will normally be reparation and reanalysis of the batch; however, if the matrix spike and matrix spike duplicate are within limits; the batch may be acceptable.

9.7.1. Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective action.

9.7.2. If dual column analysis is used the choice of which result to report is made in the same way as for samples (Section 12.1.2) unless one column is out of control, in which case the in control result is reported.

9.7.3. LCS compound lists are included in the appendices.

9.8. Matrix Spikes

For each QC batch, analyze a matrix spike and matrix spike duplicate. Spiking compounds and levels are given in the appendices. Compare the percent recovery and relative percent difference (RPD) to those in the laboratory specific historically generated limits.

- If any individual recovery or RPD falls outside the acceptable range, corrective action must occur. The initial corrective action will be to check the recovery of that analyte in the Laboratory Control Sample (LCS). Generally, if the recovery of the analyte in the LCS is within limits, then the laboratory operation is in control and analysis may proceed.
- If the recovery for any component is outside QC limits for both the Matrix spike / spike duplicate and the LCS, the laboratory is out of control and corrective action must be taken. Corrective action will normally include reparation and reanalysis of the batch.
- If a MS/MSD is not possible due to limited sample, then a LCS duplicate should be analyzed.
- The matrix spike / duplicate must be analyzed at the same dilution as the unspiked sample, unless the matrix spike components would then be above the calibration range.

9.8.1 If dual column analysis is used the choice of which result to report is made in the same way as for samples (Section 12.1 2) unless one column is out of control, in which case the in control result is reported.

9.9. Quality Assurance Summaries

Certain clients may require specific project or program QC that may supersede these method requirements. Quality Assurance Summaries should be developed to address these requirements.

9.10. STL North Canton QC Program

Further details of QC and corrective action guidelines are presented in the STL QC Program document (QA-003). Refer to this document if in doubt regarding corrective actions.

10. CALIBRATION AND STANDARDIZATION

Internal or external calibration may be used. Internal calibration is recommended unless the sample matrix is likely to interfere with the quantitation of the internal standard. In either event prepare standards containing each analyte of interest at a minimum of five concentration levels. The low level standard should be at or below the reporting limit. The other standards define the working range of the detector. Recommended calibration levels are given in the appendices.

- 10.1. A new calibration curve must be generated after major changes to the system or when the continuing calibration criteria cannot be met. Major changes include new columns, changing PID lamps or FID jets or replacing the ECD detector. A new calibration is not required after clipping the column, replacing the septum or syringe, or other minor maintenance.

- 10.2. With the exception of 10.3 below, it is NOT acceptable to remove points from a calibration curve for the purpose of meeting criteria, unless the points are the highest or lowest on the curve AND the reporting limit and/or linear range is adjusted accordingly. In any event, at least 5 points must be included in the calibration curve. Quadratic (second order) calibrations require at least six points. Third order calibrations require at least seven points.
- 10.3. A level may be removed from the calibration if the reason can be clearly documented, for example a broken vial or no purge run. A minimum of five levels must remain in the calibration. The documentation must be retained with the initial calibration. Alternatively, if the analyst believes that a point on the curve is inaccurate, the point may be reanalyzed and the reanalysis used for the calibration. All initial calibration points must be analyzed without any changes to instrument conditions, and all points must be analyzed within 24 hours.
- 10.4. External standard calibration
- Quantitation by the external standard method assumes a proportional relationship between the calibration run and the analyte in the sample. To use this approach, introduce each calibration standard into the GC using the technique that will be used for samples. The ratio of the peak height or area response to the mass or concentration injected may be used to prepare a calibration curve.

$$\text{Calibration Factor (CF)} = \frac{\text{Area or Height of Peak}}{\text{Mass Injected (ng)}}$$

Some data systems may use the inverse of this formula. This is acceptable so long as the same formula is used for standards and samples. It is also possible to use the concentration of the standard rather than the mass injected. (This would require changes in the equations used to calculate the sample concentrations). Use of peak area or height must be consistent. However, if matrix interferences would make quantitation using peak area inaccurate for a particular sample, then peak height may be used as a substitute.

- 10.5. Internal standard calibration
- 10.5.1. The internal standard approach assumes that variations in instrument sensitivity, amount injected etc. can be corrected by determining the ratio of the response of the analyte to the response of an internal standard that has been added to the extract. To use this approach, select one or more internal standard(s) that are similar in analytical behavior to the compounds of interest. Recommended internal standards are given in the appendices. The analyst must demonstrate that the measurement of the internal standard is not affected by method or matrix interferences. If the sample matrix interferes with quantitation of the internal standard, then the external standard approach must be used instead. In this event use the response factors from the previous continuing calibration to quantitate the analytes in the sample with the interference (applies only to the sample with the interference).
- 10.5.2. Introduce each calibration standard into the GC using the technique that will be used for samples. Response factors (RF) for each compound are calculated as follows:



$$RF = \frac{A_s \times C_{is}}{A_{is} \times C_s}$$

Where:

$A_s$  = Response for the analyte to be measured

$A_{is}$  = Response for the internal standard

$C_{is}$  = Concentration of internal standard

$C_s$  = Concentration of the analyte to be determined in the standard

#### 10.6. Calibration curve fits

Average response factor, linear regression, or quadratic curves may be used to fit the data. Average response factor may be used if the average % RSD of the response factors or calibration factors of all the analytes in the calibration standard taken together is  $\leq 20\%$ . The average %RSD is calculated by summing the RSD value for each analyte and dividing by the total number of analytes.

- 10.6.1. In general, for environmental analysis, average response factors are the most appropriate calibration model. Linear or curved regression fits should only be used if the analyst has reason to believe that the average RF model does not fit the normal concentration/response behavior of the detector.

#### 10.6.2. Average response factor

The average response factor may be used if the average percent relative standard deviation (%RSD) of all the response factors taken together is  $\leq 20\%$

The equation for average response factor is:

$$\text{Average response factor} = \overline{RF} = \frac{\sum_{i=1}^n RF_i}{n}$$

Where:  $n$  = Number of calibration levels

$$\sum_{i=1}^n RF_i = \text{Sum of response factors for each calibration level}$$

## 10.6.3. Linear regression

The linear fit uses the following functions:

## 10.6.3.1. External Standard

$$y = ax + b$$

or

$$x = \frac{(y - b)}{a}$$

Where:  $y$  = Instrument response

$x$  = Concentration

$a$  = Slope

$b$  = Intercept

## 10.6.3.2. Internal Standard

$$C_s = \frac{\left[ \frac{A_s C_{is}}{A_{is}} - b \right]}{a}$$

Where:  $C_s$  = Concentration in the sample

$A_s$  = Area of target peak in the sample

$A_{is}$  = Area of internal standard in the sample

$C_{is}$  = Concentration of the internal standard

## 10.6.4. Quadratic curve

The quadratic curve uses the following functions:

## 10.6.4.1. External standard

$$y = ax + cx^2 + b$$

Where  $c$  is the curvature

## 10.6.4.2. Internal Standard

$$y = a \left( \frac{A_s \times C_{is}}{A_{is}} \right) + c \left( \frac{A_s \times C_{is}}{A_{is}} \right)^2 + b$$

## 10.7. Evaluation of calibration curves

- 10.7.1. The percent relative standard error (%RSE) from the calibration curve is used to evaluate the initial calibration. This provides a measure of how much error is associated with using the calibration curve for quantitation.
- 10.7.2. The least squares regression line is calculated and used to calculate the predicted concentration for each level. The percent relative standard error is calculated as follows:

$$\% RSE = 100\% \times \sqrt{\frac{\sum_{i=1}^N \left[ \frac{C_i - PC_i}{C_i} \right]^2}{(N - P)}}$$

Where:

$N$  = Number of points in the curve

$P$  = Number of parameters in the curve (= 1 for average response factor, 2 for linear, 3 for quadratic)

$C_i$  = True concentration for level  $i$

$PC_i$  = Predicted concentration for level  $i$

Note that when average response factors are used, %RSE is equivalent to %RSD.

## 10.8. The following requirements must be met for any calibration to be used:

- Response must increase with increasing concentration.
- If a curve is used, the intercept of the curve at zero response must be less than  $\pm$  the reporting limit for the analyte.
- The average Relative Standard Error (RSD for average response factors) of the calibration points from the curve used must be  $\leq 20\%$ .
- Some data systems will not measure the %RSE from a linear or quadratic fit. For the linear case, the correlation coefficient may be used as an alternative to the %RSE, and must be greater than or equal to 0.990. For the quadratic case the Coefficient of Determination may be used, and must be greater or equal to 0.990.

**Note:** The Relative Standard Error (RSE) is superior to the Correlation Coefficient ( $r$ ) and Coefficient of Determination ( $r^2$ ) for testing the fit of a set of calibration points to a line. The lower points on a curve have little effect on  $r$ . As a result a curve may have a very good correlation coefficient ( $>0.995$ ), while also having  $> 100\%$  error at the low point.

## 10.9. Weighting of data points

10.9.1. In linear and quadratic calibration fits, the points at the lower end of the calibration curve have less absolute variance than points at the high concentration end of the curve. This can cause severe errors in quantitation at the low end of the calibration. However, in environmental analysis, accuracy at the low end of the curve is very important. For this reason it is preferable to increase the weighting of the lower concentration points.  $1/\text{Concentration}^2$  weighting (often called  $1/X^2$  weighting) will improve accuracy at the low end of the curve and should be used if the data system has this capability.

10.10. Non-standard analytes are sometimes requested. For these analytes, it may be acceptable to analyze a single standard at the reporting limit with each continuing calibration rather than a five point initial calibration. This action must be with client approval. If the analyte is detected in any of the samples, a five point initial calibration must be generated and the sample(s) reanalyzed for quantitation.

## 10.11. Calibration Verification

## 10.11.1. 12 hour Calibration

The working calibration curve or RF must be verified by the analysis of a mid point calibration standard at the beginning, after every 12 hours, and at the end of the analysis sequence. The center of each retention time window is updated with each 12-hour calibration or calibration verification.

## 10.11.2. Calibration Verification

It may be appropriate to analyze a mid point standard more frequently than every 12 hours. If these calibration verification standards are analyzed, requirements are the same as the 12-hour calibration with the exception that retention times are not updated

10.11.3. Any individual compounds with  $\%D < 15\%$  meet the calibration criteria. The calibration verification is also acceptable if the average of the  $\%D$  for all the analytes is  $< 15\%$ . This average is calculated by summing the entire absolute  $\%D$  results in the calibration (including surrogates) and dividing by the number of analytes. Any analyte that is reportable as found must have a  $\%$  difference of  $< 15\%$  in the calibration verification or 12 hour calibration, on the column used for quantitation. Refer to section 12.1.2 for which result to report.

10.11.4. It is not necessary to run a calibration verification standard at the beginning of the sequence if samples are analyzed immediately after the completion of the initial calibration.

10.11.5. Samples quantitated by external standard methods must be bracketed by calibration verification standards that meet the criteria listed above. The bracketing standards on the column used for calibration must meet the same criteria as the opening standards. Bracketing is not necessary for internal standard methods

10.11.6. If the analyst notes that a CCV has failed and can document the reason for failure (e.g. no purge, broken vial, carryover from the previous sample etc.) then a second CCV may be analyzed without any adjustments to the instrument. If this CCV meets criteria then the preceding samples have been successfully bracketed. If adjustments to the instrument are performed before the repeat CCV then the proceeding samples have not been successfully bracketed but analysis may continue.

10.11.7. In general, it is not advisable to analyze repeat CCVs on unattended runs. If repeat CCVs are analyzed then the first will serve as the bracketing standard for the preceding samples and the last will serve as the CCV for the following samples

10.11.8. If highly contaminated samples are expected it is acceptable to analyze blanks or primers at any point in the run.

10.11.9. % Difference calculation

% Difference for internal and external methods is calculated as follows

Internal Standard:

External standard:

$$\%D = \frac{RF_c - \overline{RF}}{\overline{RF}} \times 100$$

$$\%D = \frac{CF_c - \overline{CF}}{\overline{CF}} \times 100$$

Where  $RF_c$  and  $CF_c$  are the response and calibration factors from the continuing calibration

$\overline{RF}$  and  $\overline{CF}$  are the average response and calibration factors from the initial calibration

10.11.10. % Drift calculation

% Drift is used for comparing the continuing calibration to a linear or quadratic curve. The criteria for % drift are the same as for % difference

$$\% \text{ Drift} = \frac{\text{Calculated Conc.} - \text{Theoretical Conc.}}{\text{Theoretical Conc.}} \times 100\%$$

10.11.11. Corrective Actions for Continuing Calibration

If the overall average %D of all analytes is greater than  $\pm 15\%$  corrective action must be taken. This may include clipping the column, changing the liner or other minor instrument adjustments, followed by reanalyzing the standard. If the overall average %D still varies by more than  $\pm 15\%$ , a new calibration curve must be prepared.

10.11.12. Corrective Action for Samples

For internal standard methods, any samples injected after a standard not meeting the calibration criteria must be re-injected.

For external standard methods, any samples injected after the last good continuing calibration standard must be re-injected.

If the average %D for all the analytes in the calibration is over 15%, but all of the analytes requested for a particular sample have  $\%D \leq 15\%$ , then the analysis is acceptable for that sample.

## 11. PROCEDURE

### 11.1. Extraction

Extraction procedures are referenced in the appendices.

### 11.2. Cleanup

Cleanup procedures are referenced in the appendices

## 11.3. Gas Chromatography

Chromatographic conditions for individual methods are presented in the appendices.

## 11.4. Sample Introduction

In general, volatile analytes are introduced using purge and trap as described in Appendix A. Semivolatile analytes are introduced by direct injection of the extract. Samples, standards, and QC must be introduced using the same procedure.

## 11.5. Analytical Sequence

An analytical sequence starts with an initial calibration or a daily calibration. Refer to the individual method appendices for method specific details of daily calibrations and analytical sequences.

11.5.1. The daily calibration includes analysis of standards containing all single response analytes and updating the retention time windows.

11.5.2. If there is a break in the analytical sequence of greater than 12 hours, a new analytical sequence must be started with a daily calibration.

## 11.6. Retention Time Windows

11.6.1. Retention time windows must be determined for all analytes. Make an injection of all analytes of interest each day over a three-day period. Calculate the standard deviation of the three retention times for each analyte (relative retention times may also be used). For multi-response analytes (e.g., Aroclors) use the retention time of major peaks. Plus or minus three times the standard deviation of the retention times of each analyte defines the retention time window.

11.6.2. The center of the retention time window is the retention time from the last of the three standards. The centers of the windows are updated with the mid-point of the initial calibration and each 12-hour calibration. The widths of the windows will remain the same until new windows are generated following the installation of a new column.

11.6.3. If the retention time window as calculated above is less than  $\pm 0.05$  minutes, use  $\pm 0.05$  minutes as the retention time window. This allows for slight variations in retention times caused by sample matrix.

11.6.4. The laboratory must calculate new retention time windows each time a new column is installed. The new windows must be generated within one week of the installation of the new column. Until these standards have been run on the new column, the retention time windows from the old column may be used, updated with the retention times from the new initial calibration.

## 11.6.5. Corrective Action for Retention Times

The retention times of all compounds in the 12 hour calibration or calibration verification standard must be within the retention time window. If this condition is not met, all samples analyzed after the last compliant standard must be reanalyzed unless the following conditions are met for any compound that elutes outside the retention time window:

The retention time of that compound in the standard must be within a retention time range equal to twice the original window.

No peak that would be reportable may be present on the sample chromatogram within an elution time range equal to three times the original retention time window.

## 11.7. Daily Retention Time Windows

The center of the retention time windows determined in Section 11.6 are adjusted to the retention time of each analyte as determined in the 12 hour calibration standards or continuing calibration verification standards. (See the method 8081A and 8082 appendices for exceptions for multi-response components.) The retention time windows must be updated at the beginning of each analytical sequence and with each 12-hour calibration or continuing calibration verification.

#### 11.8 Percent Moisture

Analytical results may be reported as dry or wet weight, as required by the client. Percent moisture must be determined if results will be reported as dry weight. Refer to SOP CORP-OP-0001NC for determination of percent moisture.

#### 11.9. Procedural Variations

Procedural variations are allowed only if deemed necessary in the professional judgment of the supervisor to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and approved by a supervisor and QA/QC manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file. The nonconformance is also addressed in the case narrative. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

### 12. DATA ANALYSIS AND CALCULATIONS

#### 12.1. Qualitative Identification

- 12.1.1. Tentative identification occurs when a peak is found within the retention time window for an analyte, at a concentration above the reporting limit, or above the MDL if J flags are required. Normally confirmation is required on a second column, but if the detector is sufficiently specific or if the sample matrix is well enough defined, single column analysis may be adequate. In some cases GC/MS confirmation may be required. Client specific requirements may also define the need for second column confirmation and / or GC/MS confirmation. Refer to the appendices for test specific requirements for confirmation. Identification is confirmed if a peak is also present in the retention time window for that analyte on the confirmatory column, at a concentration greater than the reporting limit (MDL if J flag confirmation required)

#### 12.1.2. Dual column quantitation

For confirmed results, two approaches are available to the analyst.

A) The primary column approach

Or

B) The better result approach

Both are acceptable to avoid the reporting of erroneous or unconfirmed data.

##### 12.1.2.1. Primary column approach:

The result from the primary column is normally reported. The result from the secondary column is reported if any of the following three bulleted possibilities are true.

- There is obvious chromatographic interference on the primary column
- The result on the primary column is > 40% greater than the result on the secondary column
- Continuing or bracketing standard fails on the primary column but is acceptable on the secondary column (If the primary column result is > 40% higher than the secondary and the primary column calibration fails, then the sample must be evaluated for reanalysis.)

12.1.2.2. Better result approach

The lower of the two results is normally reported. The lower result is considered better because the higher result is generally higher because of chromatographic interference. The higher result is reported if any of the following two bulleted possibilities are true.

- There is obvious chromatographic interference on the column with the lower result
- The continuing or bracketing calibration on the column with the lower result fails. (If the higher result is > 40% higher and the calibration on the column with the lower result fails, then the sample must be evaluated for reanalysis.)

- 12.1.3. If the Relative percent difference (RPD) between the response on the two columns is greater than 40%, or if the opinion of an experienced analyst is that the complexity of the matrix is resulting in false positives, the confirmation is suspect and the results are qualified. RPD is calculated using the following formula:

$$RPD = \frac{R_1 - R_2}{\frac{1}{2}(R_1 + R_2)}$$

Where R=Result

12.1.4. Multi-response Analytes

For multi-response analytes, the analyst should use the retention time window, but should rely primarily on pattern recognition. The pattern of peaks will normally serve as confirmation.

- 12.1.5. The experience of the analyst should weigh heavily in the interpretation of the chromatogram. For example, sample matrix or laboratory temperature fluctuation may result in variation of retention times.

12.2. Calibration Range

If concentrations of any analytes exceed the working range as defined by the calibration standards, then the sample must be diluted and reanalyzed. Dilutions should target the most concentrated analyte in the upper half (over 50% of the high level standard) of the calibration range. It may be necessary to dilute samples due to matrix.

12.3. Dilutions

Samples may be screened to determine the appropriate dilution for the initial run. If the initial diluted run has no hits or hits below 20% of the calibration range and the matrix allows for analysis at a lesser dilution, then the sample must be reanalyzed at a dilution targeted to bring the largest hit above 50% of the calibration range.

12.3.1. Guidance for Dilutions Due to Matrix

If the sample is initially run at a dilution and only minor matrix peaks are, then the sample should be reanalyzed at a more concentrated dilution. Analyst judgement is required to determine the most concentrated dilution that will not result in instrument contamination

12.3.2. Reporting Dilutions

The most concentrated dilution with no target compounds above the calibration range will be reported. Other dilutions may be reported at client request, if the lower dilutions will not cause detector saturation, column overload, or carryover. Analyst judgement and client site history will factors in the reporting of dual dilutions



## 12.4. Interferences

If peak detection is prevented by interferences, further cleanup should be attempted. If no further cleanup is reasonable, then elevation of reporting levels and/or lack of positive identification must be addressed in the case narrative.

## 12.5. Internal Standard Criteria for Continuing Calibration

If internal standard calibration is used, then the internal standard response in a continuing calibration standard must be within 50 to 150% of the response in the mid level of the initial calibration.

## 12.6. Calculations

Capabilities of individual data systems may require the use of different formulas than those presented here. When this is the case, the calculations used must be shown to be equivalent and must be documented in an appendix attached to this document.

## 12.6.1. External Standard Calculations

## 12.6.1.1. Aqueous samples

$$\text{Concentration (mg / L)} = \frac{(A_s \times V_i \times D_f)}{(CF \times V_t \times V_s)}$$

Where:

$A_s$  = Response for the analyte in the sample

$V_i$  = Volume of extract injected,  $\mu\text{L}$

$D_f$  = Dilution factor

$V_t$  = Volume of total extract,  $\mu\text{L}$

$V_s$  = Volume of sample extracted or purged, mL

$CF$  = Calibration factor, area or height/ng, Section 10.1

## 12.6.1.2. Non-aqueous Samples

$$\text{Concentration (mg / kg)} = \frac{(A_s \times V_i \times D_f)}{(CF \times V_t \times W \times D)}$$

Where:

$W$  = Weight of sample extracted or purged, g

$$D = \frac{100 - \% \text{ Moisture}}{100} \quad (D = 1 \text{ if wet weight is required})$$

## 12.6.2. Internal Standard Calculations

## 12.6.2.1. Aqueous Samples

$$\text{Concentration (mg / L)} = \frac{(A_s \times C_R \times D_f)}{(A_R \times RF \times V_s)}$$

Where:

$C_{is}$  = Amount of internal standard added, ng

$A_{is}$  = Response of the internal standard

$RF$  = Response factor for analyte

#### 12.6.2.2. Non-aqueous Samples

$$\text{Concentration (mg / kg)} = \frac{(A_s \times C_{is} \times D_f)}{(A_{is} \times RF \times W \times D)}$$

#### 12.6.3. Surrogate Recovery

Concentrations of surrogate compounds are calculated using the same equations as for the target compounds. The response factor from the initial calibration is used. Surrogate recovery is calculated using the following equation:

$$\% \text{ Recovery} = \frac{\text{Concentration (or amount) found}}{\text{Concentration (or amount) spiked}} \times 100$$

### 13. METHOD PERFORMANCE

#### 13.1. Method Detection Limit

Each laboratory must generate a valid method detection limit for each analyte of interest. The MDL must be below the reporting limit for each analyte. The procedure for determination of the method detection limit is given in 40 CFR Part 136, Appendix B, and further defined in QA Policy #: QA-005.

#### 13.2. Initial Demonstration

Each laboratory must make a one time initial demonstration of capability for each individual method. Demonstration of capability for both soils and water matrices is required. This requires analysis of QC check samples containing all of the standard analytes for the method. For some tests it may be necessary to use more than one QC check mix to cover all analytes of interest.

13.2.1. Four aliquots of the QC check sample are analyzed using the same procedures used to analyze samples, including sample preparation. The concentration of the QC check sample should be equivalent to a mid-level calibration.

13.2.2. Calculate the average recovery and standard deviation of the recovery for each analyte of interest. Compare these results with the acceptance criteria given in each appendix.

13.2.3. If any analyte does not meet the acceptance criteria, the test must be repeated. Only those analytes that did not meet criteria in the first test need to be evaluated. Repeated failure for any analyte indicates the need for the laboratory to evaluate the analytical procedure and take corrective action.

## 13.3. Training Qualification

The group/team leader has the responsibility to ensure that an analyst who has been properly trained in its use and has the required experience performs this procedure.

## 14. POLLUTION PREVENTION

This method does not contain any specific modifications that serve to minimize or prevent pollution.

## 15. WASTE MANAGEMENT

15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."

## 15.2. Waste Streams Produced by the Method.

15.2.1. The following waste streams are produced when this method is carried out.

15.2.1.1. **Vials containing sample extracts:** These vials are placed in the vial waste located in the GC/MS laboratory.

15.2.1.2. **Tubes containing sample extracts for TPH, Pesticides, PCBs and Herbicides:** These capped tubes are placed in the PCB/flammable waste located in the GC prep laboratory.

15.2.1.3. **Samples, standards, and all extraction materials contaminated with high levels (>50ppm) of PCB's must be segregated into their own waste stream.** PCB wastes are collected in one of three waste streams, solid PCB, liquid PCB and PCB vial waste. PCB containing samples are located through a LIMS query and disposed of as PCB containing.

15.2.1.4. **Extracted solid samples contaminated with methylene chloride/acetone or acetone/hexane.** These materials are disposed of in the solid waste and debris in a red container located in the extractions lab.

15.2.1.5. **Discarded samples.** These samples are collected in the solid debris drum.

## 16. REFERENCES

Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW846, 3rd Edition, Final Update III, December 1996, and Section 8000B

## 17. MISCELLANEOUS

## 17.1. Modifications from Reference Method

17.1.1. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the Method Detection Limit. This SOP states that the Method Blank must not contain any analyte of interest at or above the reporting limit. Common lab contaminants are allowed to be up to 5 times the reporting limit in the blank following consultation with the client.

**GAS CHROMATOGRAPHIC ANALYSIS BASED ON  
METHOD 8000B, SW-846**

**SOP No. CORP-GC-0001NC**

Revision No. 5.7

Revision Date: 10/01/03

Page 23 of 86

---



**17.2. Modifications from Previous Revision**

The calibration criteria in section 10.11 have been rewritten to improve consistency with SW-846 and to improve clarity.


**17.3. Facility Specific SOPs**

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none.

## 17.4. Flow Diagrams

17.4.1. Initial demonstration and MDL<sup>1</sup>

<sup>1</sup> This flow diagram is for guidance and cannot cover all eventualities. Consult the SOP text and a supervisor if in doubt.

17.4.2. Sample Analysis<sup>1</sup>

<sup>1</sup> This flow diagram is for guidance and cannot cover all eventualities. Consult the SOP text and a supervisor if in doubt.

**ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B****1. SCOPE AND APPLICATION**

- 1.1. This method describes sample preparation and extraction for the analysis of volatile organics by a purge and trap procedure, following method 8021B. However, where required by a client QAPP this section may also be used to analyze aromatic volatiles by discontinued methods 8020A and 8010B. All requirements of the 8000B section of this SOP must be met except when superseded by this Appendix. Refer to Table A-1 for the individual analytes normally determined by these procedures.
- 1.2. Compounds within the scope of this method have boiling points below 200°C and are soluble or slightly soluble in water. Classes of compounds best suited to purge-and-trap analysis include low molecular weight halogenated hydrocarbons, aromatics, ketones, nitriles, acetates, acrylates, ethers, and sulfides.
- 1.3. Water samples and soils samples with low levels of contamination may be analyzed directly by purge-and-trap extraction and gas chromatography. Higher concentrations of these analytes in soil may be determined by the medium level methanol extraction procedure.
- 1.4. This method also describes the preparation of water-miscible liquids, non-water-miscible liquids, solids, wastes, and soils/sediments for analysis by the purge-and-trap procedure.
- 1.5. The associated LIMS method code is QR

**2. SUMMARY OF METHOD**

- 2.1. An inert gas is bubbled through the sample at ambient temperature or at 40°C (40°C required for low-level soils), and the volatile components are transferred from the aqueous phase to the vapor phase. The vapor is swept through a sorbent column where the volatile components are adsorbed. After purging is completed, the sorbent column is heated and back-flushed with inert gas to desorb the components onto a gas chromatographic column. Analytes are detected using a photoionization Detector, an electrolytic conductivity detector or a combination of both.
- 2.2. For soil samples, a portion of the sample is dispersed in methanol to dissolve the volatile organic constituents. A portion of the methanolic solution is combined with water. It is then analyzed by purge-and-trap GC following the normal water method. If very low detection limits are needed for soil samples then direct purge using sodium bisulfate preservation may be necessary.

**3. DEFINITIONS**

Refer to the STL North Canton Laboratory Quality Manual (LQM), current version, for definitions of terms used in this SOP.

**4. INTERFERENCES**

- 4.1. Refer to section 4 of the method 8000B part of this SOP for general information on chromatographic interferences
- 4.2. Impurities in the purge gas, and from organic compounds out-gassing from the plumbing ahead of the trap, account for the majority of contamination problems. The analytical system must be demonstrated to be free from contamination under the conditions of the analysis by running laboratory reagent blanks. The use of non-TFE plastic tubing, non-TFE thread sealants, or flow controllers with rubber components in the purging device should be avoided.
- 4.3. Samples can be contaminated by diffusion of volatile organics (particularly methylene chloride and fluorocarbons) through the septum seal of the sample vial during shipment and storage. A trip blank

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A2 of A17

**ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**

prepared from organic-free reagent water and carried through sampling and handling protocols serves as a check on such contamination.

- 4.4. Contamination by carryover can occur whenever high-concentration and low-concentration samples are analyzed sequentially. Whenever an unusually concentrated sample is analyzed, it should be followed by an analysis of organic-free reagent water to check for cross-contamination. The trap and other parts of the system are subject to contamination. Therefore, frequent bake-out and purging of the system may be required.
- 4.5. When utilizing an autosampler system, which has multiple ports for sample analysis, it is likely that only a single stage or port may be contaminated by a highly concentrated sample. If a port is suspect, a water blank should be analyzed to verify lack of contamination. If the water blank and subsequent blanks on that port show contamination consistent with the concentrated sample, further maintenance is required. This may include replacing or cleaning the multi-port valve, transfer lines, etc.
- 4.6. A holding blank is kept in the sample refrigerator. This is analyzed and replaced every 14 days. If the holding blank does not meet the method blank criteria, the source of contamination must be found and corrected. Evaluation of all samples analyzed in the 14-day period prior to the analysis of the contaminated holding blank is required.
- 4.7. Acidification of samples may result in hydrolysis of 2-chloroethyl-vinyl ether.

**5. SAFETY**

- 5.1. Refer to section 5 of the Method 8000B section of this SOP for general safety requirements.
- 5.2. Often, purge vessels on purge-and-trap instrumentation are pressurized by the time analysis is completed. Therefore, vent the pressure prior to removal of these vessels to prevent the contents from spraying out.
- 5.3. The toxicity or carcinogenicity of each chemical used in this procedure has not been fully defined. Additional health and safety information can be obtained from the MSDS files maintained in the laboratory. The following specific hazards are known.

Methanol -- Flammable and toxic
- 5.4. The following method analytes have been tentatively classified as known or suspected human or mammalian carcinogens: Benzene, Carbon Tetrachloride, 1,4-Dichlorobenzene, 1,2-Dichloroethane, Hexachlorobutadiene, 1,1,2,2-Tetrachloroethane, 1,1,2-Trichloroethane, Chloroform, 1,2-Dibromoethane, Tetrachloroethene, Trichloroethene, Vinyl Chloride. Pure standard materials and stock standard solutions of these compounds should be handled in a hood.
- 5.5. Methanol shall not be used in a CaptAir hood.
- 5.6. GC VOA instruments use an ultraviolet (UV) light source, which must be shielded from view. There should also be a warning label/sticker on each instrument that identifies it as a UV light source.
- 5.7. Sodium bisulfate creates Sulfuric Acid when mixed with water.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Microsyringes -- 10µL, 25µL, 100µL, 250µL, 500µL, and 1000µL. These should be equipped with a 20 gauge (0.006" ID) needle. These will be used to measure and dispense methanolic solutions and aqueous samples.



**ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**

- 6.2. Gas tight syringes -- 5 mL and 25 mL. Used for measuring sample volumes.
- 6.3. Purge and Trap Apparatus -- A device capable of extracting volatile compounds, trapping on a sorbent trap, and introducing onto a gas chromatograph.
- 6.4. Purge and Trap Autosampler -- In order to maintain high sample throughput, an autosampler is highly recommended.
- 6.5. Trap -- The trap used is dependent on the class of compound to be analyzed. Refer to Table A-2 for suggested traps for specific tests.
- 6.6. Purge Vessels -- These are dependent on the purge and trap unit/autosampler used. Both disposable culture tubes (needle sparge units) and specially designed vessels with fritted bottoms may be used. Follow the manufacturer's suggestions for configuration.
- 6.7. Columns - Refer to Table A-2 for details of columns.
- 6.8. Volumetric flasks, Class A: 5 mL to 250 mL
- 6.9. pH paper
- 6.10. Balance capable of weighing to 0.01g for samples

**7. REAGENTS AND SUPPLIES**

- 7.1. Refer to the method 8000B section of this SOP for general requirements for reagents and supplies.
- 7.2. Organic Free Water

Organic free water is defined as water in which an interferent is not observed at the reporting limit of the compounds of interest. Suggested methods for generating organic free water include

  - Filtration through a carbon bed.
  - Continuously sparging water with helium or nitrogen.
  - Use of commercial water purification systems.

Other methods may be used, so long as the requirement that the water show no interference is met. The procedure used should be documented in a lab specific attachment.
- 7.3. Sodium Bisulfate
- 7.4. Methanol -- Purge and Trap Grade
- 7.5. Standards

Refer to tables A-5 and A-6 for details of surrogate, matrix spiking and internal standards. Calibration standard levels are not specified, since they may depend on the sensitivity and linear range of specific detectors. However, the low level standard must be equivalent to the reporting limits specified in Table A-1.

  - 7.5.1. Volatile standards are prepared by injecting a measured volume of the stock standard into a syringe containing the appropriate volume of organic free water. The calibration standard is then loaded into the purge device.

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A4 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

## 8. SAMPLE COLLECTION, PRESERVATION, AND STORAGE

- 8.1. Holding times for all volatile analysis are 14 days from sample collection.
- 8.2. Water samples are normally preserved at pH < 2 with 1:1 hydrochloric acid. If residual chlorine is present, 2 drops of 10% sodium thiosulfate are added.
- 8.3. Solid samples are field preserved with sodium bisulfate solution for low level analysis, or with methanol for medium level analysis. Soil samples can also be taken using the EnCore™ sampler and preserved in the lab within 48 hours of sampling. At specific client request, unpreserved soil samples may be accepted.
- 8.4. There are several methods of sampling soil. The recommended method, which provides the minimum of field difficulties, is to take an EnCore sample. (The 5 g or 25 g sampler can be used, depending on client preference). Following shipment back to the lab the soil is preserved in methanol. This is the medium level procedure. If very low detection limits are needed (< 50 µg/kg for most analytes) then it will be necessary to use two additional 5 g EnCore samplers or to use field preservation.
- 8.5. Sample collection for medium level analysis using EnCore samplers
  - 8.5.1. Ship one 5 g (or 25 g) EnCore sampler per field sample position.
  - 8.5.2. An additional bottle must be shipped for percent moisture determination
  - 8.5.3. When the samples are returned to the lab, extrude the (nominal) 5g (or 25 g) sample into a tared VOA vial containing 5 mL methanol (25 mL methanol for the 25 g sampler). Obtain the weight of the soil added to the vial and note on the label.
  - 8.5.4. Add the correct amount of surrogate spiking mixture. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample )
  - 8.5.5. Add the correct amount of matrix spiking solution to the matrix spike and matrix spike duplicate samples. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.) The addition of spike introduces a slight error, (0.4%) which can be neglected, into the calculations.
  - 8.5.6. Prepare an LCS for each batch by adding the correct amount of matrix spiking solution to clean methanol. (100 µL of spike to 25 mL methanol or 20 µL spike to 5 mL methanol).
  - 8.5.7. Shake the samples for two minutes to distribute the methanol throughout the soil.
  - 8.5.8. Allow to settle, then remove a portion of methanol and store in a clean Teflon capped vial at 4±2°C until analysis
- 8.6. Sample collection for medium level analysis using field methanol preservation
  - 8.6.1. Prepare a VOA vial by adding 5 mL purge and trap grade methanol. (If a 25 g sample is to be used, add 25 mL methanol to the VOA vial).
  - 8.6.2. Seal the bottle and attach a label.
  - 8.6.3. Weigh the bottle to the nearest 0.01g and note the weight on the label.
  - 8.6.4. Ship with appropriate sampling instructions.

## APPENDIX A

SOP No CORP-GC-000INC

Revision No. 5.7Revision Date: 10/01/03

Page A 5 of A 17

## ANALYSIS OF VOLATILE ORGANICS BASED ON

## METHOD 8021B

- 
- 8.6.5. Each sample will require an additional bottle with no preservative for percent moisture determination.
- 8.6.6. At client request, the methanol addition and weighing may also be performed in the field.
- 8.6.7. When the samples are returned to the lab, obtain the weight of the soil added to the vial and note on the label.
- 8.6.8. Add the correct amount of surrogate spiking mixture. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.)
- 8.6.9. Add the correct amount of matrix spiking solution to the matrix spike and matrix spike duplicate samples. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.) The addition of spike introduces a slight error, (0.4%) which can be neglected, into the calculations.
- 8.6.10 Prepare an LCS for each batch by adding the correct amount of matrix spiking solution to clean methanol. (100 µL of spike to 25 mL methanol or 20 µL spike to 5 mL methanol).
- 8.6.11 Shake the samples for two minutes to distribute the methanol throughout the soil.
- 8.6.12. Allow to settle, then remove a portion of methanol and store in a clean Teflon capped vial at 4±2°C until analysis.
- 8.7. Low level procedure
- 8.7.1. If low detection limits are required (typically < 50 µg/kg) sodium bisulfate preservation must be used. However, it is also necessary to take a sample for the medium level (methanol preserved) procedure, in case the concentration of analytes in the soil is above the calibration range of the low-level procedure.
- 8.7.2 A purge and trap autosampler capable of sampling from a sealed vial is required for analysis of samples collected using this method (Varian Archon or O.I. 4552)
- 8.7.3. The soil sample is taken using a 5g EnCore sampling device and returned to the lab. It is recommended that two EnCore samplers be used for each field sample position, to allow for any reruns than may be necessary. A separate sample for % moisture determination is also necessary.
- 8.7.4. Prepare VOA vials by adding a magnetic stir bar, approximately 1 g of sodium bisulfate and 5 mL of reagent water
- 8.7.5. Seal the vial and attach a label. The label must not cover the neck of the vial or the autosampler will malfunction.
- 8.7.6. Weigh the vial to the nearest 0.01g and note the weight on the label
- 8.7.7. Extrude the soil sample from the EnCore sampler into the prepared VOA vial. Reweigh the vial to obtain the weight of soil and note on the label
- 8.7.8. **Note:** Soils containing carbonates may effervesce when added to the sodium bisulfate solution. If this is the case at a specific site, add 5 mL of water instead, and freeze at ≥ -10°C until analysis.
- 8.7.9. Alternatively the sodium bisulfate preservation may be performed in the field. Ship at least two vials per sample. The field samplers must determine the weight of soil sampled. Each sample will

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03



Page A 6 of A 17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

require an additional bottle with no preservative for percent moisture determination, and an additional bottle preserved with methanol for the medium level procedure.


- 8.8. Aqueous samples are stored in glass containers with Teflon lined septa at 4°C +/- 2°C, with minimum headspace.
- 8.9. Medium level solid extracts are aliquoted into 2 - 5 mL glass vials with Teflon lined caps and stored at 4°C +/- 2°C. The extracts are stored with minimum headspace.
- 8.10. The maximum holding time is 14 days from sampling until the sample is analyzed. (Samples that are found to be unpreserved still have a 14 day holding time. However they should be analyzed as soon as possible. The lack of preservation should be addressed in the case narrative). Maximum holding time for the EnCore sampler (before the sample is added to methanol or sodium bisulfate) is 48 hours.
- 8.11. A holding blank is stored with the samples. This is analyzed and replaced if any of the trip blanks show any contamination. Otherwise it is replaced every 14 days.

**Note:** Freezing is not allowed for Ohio VAP solids.

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**EnCore procedure when low level is not required (field steps  
in gray)****EnCore procedure when low level is required**

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

## Field methanol extraction procedure (field steps in gray)



## APPENDIX A

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03Page A 9 of A 17

---

**9. QUALITY CONTROL**

- 9.1. Refer to the method 8000B section of this SOP, section 9, for general quality control procedures, including batch definition, requirements for method blanks, LCS, matrix spikes, surrogates, and control limits.

**10. CALIBRATION AND STANDARDIZATION**

- 10.1. Refer to the method 8000B section of this SOP, section 10, for general calibration procedures.

10.2. Gas Chromatograph Operating Conditions

Various column configurations are possible. If dual column confirmation is necessary, the sample may be split using a Y splitter at the injector end to direct the sample to two columns and two detectors. For simultaneous determination of aromatic and halogenated volatiles, a single column is used and the PID and ELCD detectors are connected in series.

- 10.2.1. Refer to Table A-2, A-3 and A-4 for GC operating conditions.

10.3. Initial Calibration

- 10.3.1 Refer to Section 10 of the 8000B section of this SOP for details of initial calibration criteria.

- 10.3.2 Low level soil samples must be purged at 40°C; therefore the calibration curve must also be purged at 40°C. In addition, the low level soil calibration solutions should contain approximately the same amount of sodium bisulfate as the samples.

- 10.3.3. The low level calibration must be at the reporting limit or below. The remaining standards encompass the working range of the detector.

- 10.3.4 Calibrate the instrument using the same volume that will be used during sample analysis.

10.4 Calibration Verification

- 10.4.1. A mid level calibration standard is used for the calibration verification. The gases have 20 % D criteria rather than the 15% used for other analytes.

- 10.4.2. A calibration verification run is performed after every 10 samples for this method.

- 10.4.3 Bracketing of samples with calibration verification runs is only necessary for external standard analysis.

**11. PROCEDURE**

- 11.1 Refer to the method 8000B section of this SOP for general procedural requirements.

11.2. Analytical Sequence

The analytical sequence starts with an initial calibration of at least five points, or a 12 hour calibration that meets % difference criteria from an existing initial calibration

11.3 Confirmation

The PID and ELCD detectors are sufficiently selective that second column confirmation is not always necessary. Requirements for second column confirmation should be decided in consultation with the

client. If the PID and ELCD are used in series confirmatory information for many analytes can be gained by comparing the relative response from the two detectors.

11.4. Aqueous Sample Analysis (Purge and Trap units using sparge vessels)

- 11.4.1. Depending on the sensitivity of the instrument and capabilities of the purge and trap device, 5, 10, 20, or 25 mL sample volumes may be analyzed. A 5 mL sample volume is recommended.
- 11.4.2. Rinse a 5 mL (or 25 mL for larger sample volumes) gas-tight syringe with organic free water. Fill the syringe with the sample to be analyzed, and compress to volume
- 11.4.3. Check and document the pH of the sample remaining in the VOA vial after loading the syringe.
- 11.4.4. This procedure invalidates the contents of the VOA vial for further analysis, unless an aliquot is transferred to a smaller VOA vial with no headspace (e.g., 20 mL) at the same time the analysis aliquot is removed.
- 11.4.5. Spike with the appropriate volume of surrogate/internal standard solution and spike solution (if required) through the barrel of the syringe. The method blank is spiked with surrogates only, the LCS and matrix spikes with the surrogate and matrix spiking solutions. Refer to Tables A-5 and A-6 for volumes and concentrations of spiking solutions.
- 11.4.6. Load onto the purge and trap device and start the run.
- 11.4.7. If the initial analysis of a sample or a dilution of the sample has a concentration of analytes that exceeds the initial calibration range, the sample must be reanalyzed at a higher dilution. When a sample has a high response for a compound, an organic free water blank should follow analysis. It is recognized that during automated unattended analysis, this may not occur. If any potential carryover hits are present in samples following highly contaminated samples, the sample must be reanalyzed to determine if any of these hits are a result of carryover or are actually present in the sample.
- 11.4.8. Dilutions may be made in gas tight syringes unless the volume of sample used is less than 5  $\mu$ L, in which case dilution in volumetric flasks will be necessary.
  - 11.4.8.1. Spike with the same volume of surrogate/internal standard solution as used for undiluted samples prior to loading onto the purge and trap device
  - 11.4.8.2. For Matrix spike / matrix spike duplicates where the sample requires dilution, the sample is spiked after the dilution is performed.

11.5. Aqueous and Soil Sample Analysis (Purge and Trap units that sample directly from the VOA vial)

- 11.5.1. Units, which sample from the VOA vial, should be equipped with a module, which automatically adds surrogate and internal standard solution to the sample prior to purging the sample
- 11.5.2. If the autosampler uses automatic IS/SS injection, no further preparation of the VOA vial is needed. Otherwise the internal and surrogate standards must be added to the vial. *Note:* Aqueous samples with high amounts of sediment present in the vial may not be suitable for analysis on this instrumentation, or they may need to be analyzed as soils.
- 11.5.3. Sample remaining in the vial after sampling with one of these mechanisms is no longer valid for further analysis. A fresh VOA vial must be used for further sample analysis.
- 11.5.4. Check the pH of the sample remaining in the VOA vial after analysis is completed.



11.6. *Low-Level Solids Analysis using discrete autosamplers*

**Note:** This technique may seriously underestimate analyte concentration and must not be used except at specific client request for the purpose of comparability with previous data. It is no longer part of SW-846.

*This method is based on purging a heated sediment/soil sample mixed with reagent water containing the surrogate and, if applicable, internal and matrix spiking standards. Analyze all reagent blanks and standards under the same conditions as the samples (e.g., heated). The calibration curve is also heated during analysis. Purge temperature is 40°C.*

11.6.1. *Do not discard any supernatant liquids. Mix the contents of the container with a narrow metal spatula.*

11.6.2. *Weigh out 5 g (or other appropriate aliquot) of sample into a disposable culture tube or other purge vessel. Record the weight to the nearest 0.1 g. If method sensitivity is demonstrated, a smaller aliquot may be used. Do not use aliquots less than 1.0 g. If the sample is contaminated with analytes such that a purge amount less than 1.0 g is appropriate, use the medium level method described in section 11.7.*

11.6.3. *Connect the purge vessel to the purge and trap device.*

11.6.4. *Rinse a 5 mL gas-tight syringe with organic free water, and fill. Compress to 5 mL. Add surrogate/internal standard (and matrix spike solutions if required.) (See Tables A-5, A-6, A-7 and A-8.) Add directly to the sample from 11.6.2.*

11.6.5. *The above steps should be performed rapidly and without interruption to avoid loss of volatile organics.*

11.6.6. *Add the heater jacket or other heating device and start the purge and trap unit.*

11.6.7. *Soil samples that have low IS recovery when analyzed (<50%) should be reanalyzed once to confirm matrix effect. If external standard calibration is used, samples with surrogate recovery below the control limit should be reanalyzed once to confirm matrix effect.*

11.7. *Methanol Extract Soils*

11.7.1. *Rinse a gas-tight syringe with organic free water. Fill the syringe with the same volume of organic free water as used in the calibrations. Add no more than 2% (v/v) (100 µL for a 5 mL purge) methanolic extract (from Section 8.5 or 8.6) to the syringe. Add internal standard (if used). Load the sample onto the purge and trap device and analyze as for aqueous samples. If less than 5 µL of methanolic extract is to be added to the water, dilute the methanolic extract such that a volume greater than 5 µL will be added to the water in the syringe.*

12. **DATA ANALYSIS AND CALCULATIONS**

Refer to section 12 of the 8000B section of this SOP.

13. **METHOD PERFORMANCE**

13.1. Performance limits for the four replicate initial demonstration of capability required under Section 13.1 of the 8000B section of this SOP.

14. **POLLUTION PREVENTION**

This method does not contain any specific modifications that serve to minimize or prevent pollution.

**15. WASTE MANAGEMENT**

- 15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."
- 15.2. Waste streams produced by the method.
  - 15.2.1. The following waste streams are produced when this method is carried out.
    - 15.2.1.1. **Acidic material from the auto-sampler.** Waste stream must be collected and neutralized before discharge to a sewer system if the pH is less than 4.
    - 15.2.1.2. **Methanol waste from rinses and standards.** Methanol waste is discarded as a flammable liquid.
    - 15.2.1.3. **All samples including purged and extracted soils and waters:** Samples are collected in boxes and removed from the lab to storage. The waste coordinator handles crushing the vials and proper disposal.

**16. REFERENCES**

- 16.1. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW846, 3rd Edition, Final Update III, December 1996, Sections 5000, 5030B, 5035 and 8021B

**17. MISCELLANEOUS**

- 17.1. Modifications from Reference Method
- 17.2. Modifications from previous revision
  - 17.2.1. No revisions were made to this appendix.
- 17.3. Facility Specific SOPs

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none.

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

## 17.4. Tables

Table A-1					
Standard Analyte List					
Test	Compound	CAS number	Reporting Limit, µg/L or µg/kg		
			Aqueous	Low Soil	Medium Soil
Halogenated volatiles by 8021B	Bromodichloromethane	75-27-4	1.0	1.0	50
	Bromoform	75-25-2	1.0	1.0	50
	Bromomethane	74-83-9	1.0	1.0	50
	Carbon Tetrachloride	56-23-5	1.0	1.0	50
	Chlorobenzene	108-90-7	1.0	1.0	50
	Chloroethane	70-00-3	1.0	1.0	50
	2-Chloroethyl vinyl ether	110-75-8	5.0	5.0	250
	Chloroform	67-66-3	1.0	1.0	50
	Chloromethane	74-87-3	1.0	1.0	50
	Dibromochloromethane	124-48-1	1.0	1.0	50
	1,2-Dichlorobenzene	95-50-1	1.0	1.0	50
	1,3-Dichlorobenzene	541-73-1	1.0	1.0	50
	1,4-Dichlorobenzene	106-46-7	1.0	1.0	50
	Dichlorodifluoromethane	75-71-8	1.0	1.0	50
	1,1-Dichloroethane	75-34-3	1.0	1.0	50
	1,2-Dichloroethane	107-06-2	1.0	1.0	50
	1,1-Dichloroethene	75-45-4	1.0	1.0	50
	cis-1,2-Dichloroethene	156-59-4	1.0	1.0	50
	trans-1,2-Dichloroethene	156-60-5	1.0	1.0	50
	Dichloromethane(DCM)	75-09-2	5.0	5.0	250
	1,2-Dichloropropane	78-87-5	1.0	1.0	50
	cis-1,3-Dichloropropene	10061-01-5	1.0	1.0	50
	trans-1,3-Dichloropropene	10061-02-6	1.0	1.0	50
	1,1,2,2-Tetrachloroethane	79-34-5	1.0	1.0	50
	Tetrachloroethene	127-18-4	1.0	1.0	50
	1,1,1-Trichloroethane	71-55-6	1.0	1.0	50
	1,1,2-Trichloroethane	79-00-5	1.0	1.0	50
	Trichloroethene	79-01-6	1.0	1.0	50
	Trichlorofluoromethane	75-69-4	1.0	1.0	50
	Vinyl Chloride	75-01-4	1.0	1.0	50
Additional halogenated volatiles	Benzyl Chloride	100-44-7	5.0	5.0	250
	Bromobenzene	108-86-1	1.0	1.0	50
	Dibromomethane	74-95-3	1.0	1.0	50
	1,1,1,2-Tetrachloroethane	630-20-6	1.0	1.0	50
	Freon 113	76-13-1	1.0	1.0	50
	1,2,3-Trichloropropane	96-18-4	1.0	1.0	50
BTEX by	Benzene	71-43-2	1.0	1.0	50
8021B	Ethyl Benzene	100-41-4	1.0	1.0	50
	Toluene	108-88-3	1.0	1.0	50
	Xylenes (total)	1330-20-7	1.0	1.0	50
Aromatic	Benzene	71-43-2	1.0	1.0	50

## APPENDIX A

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

SOP No. CORP-GC-0001NC

Revision No 5.7

Revision Date: 10/01/03

Page A14 of A17

877 741

Table A-1					
Standard Analyte List					
Test	Compound	CAS number	Reporting Limit, µg/L or µg/kg		
			Aqueous	Low Soil	Medium Soil
volatiles by 8021B					
	Chlorobenzene	108-90-7	1.0	1.0	50
	1,2-Dichlorobenzene	75-34-3	1.0	1.0	50
	1,3-Dichlorobenzene	107-06-2	1.0	1.0	50
	1,4-Dichlorobenzene	75-45-4	1.0	1.0	50
	Ethyl Benzene	100-41-4	1.0	1.0	50
	Toluene	108-88-3	1.0	1.0	50
	Xylenes (total)	1330-20-7	1.0	1.0	50
Additional aromatic and unsaturated volatiles	1,2,4 Trimethylbenzene	95-63-6	1.0	1.0	50
	1,3,5 Trimethylbenzene	108-67-8	1.0	1.0	50
	Acetone	67-64-1	10	10	500
	MEK (2-butanone)	78-93-3	5.0	5.0	250
	MIBK (4-methyl-2-pentanone)	108-10-1	5.0	5.0	250
	Naphthalene	91-20-3	2.0	2.0	250
	Styrene	100-42-5	1.0	1.0	50
	Methyl tert-butyl ether (MTBE)	1634-04-4	1.0	1.0	50
Combined halogenated and aromatic volatiles by 8021B	Benzene	71-43-2	1.0	1.0	50
	Bromobenzene	108-86-1	1.0	1.0	50
	Bromochloromethane	74-97-5	1.0	1.0	50
	Bromodichloromethane	75-27-4	1.0	1.0	50
	Bromoform	75-25-2	1.0	1.0	50
	Bromomethane	74-83-9	1.0	1.0	50
	n-butylbenzene	104-51-8	1.0	1.0	50
	sec-Butylbenzene	135-98-8	1.0	1.0	50
	tert-Butylbenzene	98-06-6	1.0	1.0	50
	Carbon Tetrachloride	56-23-5	1.0	1.0	50
	Chlorobenzene	108-90-7	1.0	1.0	50
	Chlorodibromomethane	124-48-1	1.0	1.0	50
	Chloroethane	75-00-3	1.0	1.0	50
	Chloroform	67-66-3	1.0	1.0	50
	Chloromethane	74-87-3	1.0	1.0	50
	2-Chlorotoluene	95-49-8	1.0	1.0	50
	4-Chlorotoluene	106-43-4	1.0	1.0	50
	1,2-Dibromo-3-Chloropropane(DBCP)	96-12-8	1.0	1.0	50
	1,2-Dibromoethane(EDB)	106-93-4	1.0	1.0	50
	Dibromomethane	74-95-3	1.0	1.0	50
	1,2-Dichlorobenzene	95-50-1	1.0	1.0	50
	1,3-Dichlorobenzene	541-73-1	1.0	1.0	50
	1,4-Dichlorobenzene	106-46-7	1.0	1.0	50
	Dichlorodifluoromethane	75-71-8	1.0	1.0	50
	1,1-Dichloroethane	75-34-3	1.0	1.0	50
	1,2-Dichloroethane	107-06-2	1.0	1.0	50
	1,1-Dichloroethene	75-35-4	1.0	1.0	50

877 742

## APPENDIX A

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

SOP No. CORP-GC-0001NC

Revision No 5.7

Revision Date: 10/01/03

Page A15 of A17

Table A-1					
Standard Analyte List					
Test	Compound	CAS number	Reporting Limit, µg/L or µg/kg		
			Aqueous	Low Soil	Medium Soil
	cis-1,2-Dichloroethene	156-59-4	1.0	1.0	50
	trans-1,2-Dichloroethene	156-60-5	1.0	1.0	50
	1,2-Dichloropropane	78-87-5	1.0	1.0	50
	1,3-Dichloropropane	142-28-9	1.0	1.0	50
	2,2-Dichloropropane	590-20-7	1.0	1.0	50
	1,1-Dichloropropene	563-58-6	1.0	1.0	50
	cis-1,3-Dichloropropene	10061-01-5	1.0	1.0	50
	trans-1,3-Dichloropropene	10061-02-6	1.0	1.0	50
	Ethylbenzene	100-41-4	1.0	1.0	50
	Hexachlorobutadiene	87-68-3	1.0	1.0	50
	Isopropylbenzene	98-82-8	1.0	1.0	50
	p-Isopropyltoluene	99-87-6	1.0	1.0	50
	Methylene Chloride	75-09-2	5.0	5.0	250
	Naphthalene	91-20-3	2.0	2.0	250
	n-Propylbenzene	10306501	1.0	1.0	50
	Styrene	100-42-5	1.0	1.0	50
	1,1,1,2-Tetrachloroethane	630-20-6	1.0	1.0	50
	1,1,2,2-Tetrachloroethane	79-34-5	1.0	1.0	50
	Tetrachloroethene	127-18-4	1.0	1.0	50
	Toluene	108-88-3	1.0	1.0	50
	1,2,3-Trichlorobenzene	87-61-6	1.0	1.0	50
	1,2,4-Trichlorobenzene	120-82-1	1.0	1.0	50
	1,1,1-Trichloroethane	71-55-6	1.0	1.0	50
	1,1,2-Trichloroethane	79-00-5	1.0	1.0	50
	Trichloroethene	79-01-6	1.0	1.0	50
	Trichlorofluoromethane	75-69-4	1.0	1.0	50
	1,2,3-Trichloropropane	96-18-4	1.0	1.0	50
	1,2,4-Trimethylbenzene	95-63-6	1.0	1.0	50
	1,3,5-Trimethylbenzene	108-67-8	1.0	1.0	50
	Vinyl Chloride	75-01-4	1.0	1.0	50
	Xylenes (total)	1330-20-7	1.0	1.0	50

Table A-2	
Recommended Conditions for Aromatic Volatiles	
Parameter	Recommended Conditions
Temperature program	50°C, 1min, 10°C/min to 200°C, 1min
Column 1	Rtx-502.2 or DB-502.2 60m x 0.53mm 3.0µm
Column 2	Rtx-1 or DB-1 60m x 0.53mm 3.0 µm
Carrier gas	Helium or hydrogen
Purge Flow / time	40 mL/min, 11 minutes
Desorb Temp / time	180°C, 2 minutes (220°C for Vocab 3000)
Bake Time / temp	200°C, 12 minutes (230°C for Vocab 3000)
Transfer line / valve temp	115°C

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A 16 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

Table A-3 Recommended Conditions for Method Halogenated Volatiles	
Parameter	Recommended Conditions
Temperature program	35°C, 12 min, then 4°C/min to 200°C, hold for 5 min
Column 1	DB-VRX or RTX-502.2 105m x 0.53 mm id df = 3.0um
Column 2	DB-1 or RTX-1 105m x 0.53 mm ID df = 3.0um
Column 3	Rtx - Volatiles 120m x 0.53mm ID df=2.0um
Carrier gas	Helium
Purge Flow / time	40 mL/min, 11 minutes
Desorb Temp / time	180°C, 2 minutes (220°C for Vocarb 3000)
Bake Time / temp	200°C, 12 minutes (230°C for Vocarb 3000)
Transfer line / valve temp	115°C

Table A-4 Recommended Conditions for Method Combined Aromatic and Halogenated Volatiles	
Parameter	Recommended Conditions
Temperature program	35°C, 12 min, then 4°C/min to 200°C, hold for 5 min
Column 1	DB-VRX or RTX-502.2 105m x 0.53 mm id df = 3.0um
Column 2	DB-1 or RTX-1 105m x 0.53 mm ID df = 3.0um
Column 3	Rtx - Volatiles 120m x 0.53mm ID df=2.0um
Carrier gas	Helium
Purge Flow / time	40 mL/min, 11 minutes
Desorb Temp / time	180°C, 2 minutes (220°C for Vocarb 3000)
Bake Time / temp	200°C, 12 minutes (230°C for Vocarb 3000)
Transfer line / valve temp	115°C

## APPENDIX A

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

SOP No. CORP-GC-0001NC

Revision No. 5.7

Revision Date: 10/01/03

Page A17 of A17

Table A-5 Surrogate and Internal Standard Concentrations				
Standard	Components	Working Solution µg/mL	Spike amount µL (for 5 mL purge)	Final concentration µg/L (µg/kg)
Aromatic volatiles IS/SS	4-Chlorotoluene (SS)	20	5	20
	1-Chloro-4-fluorobenzene (IS)	40		40
Halogenated volatiles IS/SS	4-chlorotoluene (SS)	20	5	20
	1-Chloro-4-fluorobenzene (IS)	40		40
Combined Aromatic and halogenated volatiles IS/SS	Fluorobenzene (SS)	20	5	20
	1,4-Dichlorobutane (SS)	20		20
	1-Chloro-4-fluorobenzene (IS)	40		40

It may be necessary to select different surrogates in order to minimize sample interferences. 1-chloro-4-fluorobenzene and 4-chlorotoluene are fairly well resolved from analytes listed in this SOP. However 4-chlorotoluene may sometimes be requested as a target analyte. Other surrogates that may be considered, and issues associated with their use are.

Bromochloromethane: Elutes very close to chloroform and cis-1, 2-dichloroethene on the 502.2 column.  
May be a target analyte.

1,2-Bromochloroethane

1-Chloro-2-fluorobenzene: Elutes close to ethylbenzene on DB-1 or Rtx-1 and close to m,p-xylene on 502.2

a,a,a-Trifluorotoluene. Good for aromatic volatiles, coelutes or very close to trichloroethene

Bromofluorobenzene Close to 1,1,2,2-trichloroethane and 1,2,3-trichloropropane on the 502.2 column.

Good on DB-1 or Rtx-1.

2-Bromo-1-chloropropane: May coelute with 1,1,2-trichloroethane

Table A-6 Concentrations for LCS and MS/MSD compounds				
Standard	Components	Working Solution µg/mL	Spike amount µL (5 mL purge)	Final concentration µg/L (µg/kg)
Aromatic	Benzene	20	5	20
	Toluene	20		20
	Chlorobenzene	20		20
Halogenated	Chlorobenzene	20	5	20
	1,1-Dichloroethene	20		20
	Trichloroethene	20		20
Combination aromatic / halogenated	Benzene	20	5	20
	Toluene	20		20
	Chlorobenzene	20		20
	1,1-Dichloroethene	20		20
	Trichloroethene	20		20

## APPENDIX B

SOP No. CORP-GC-0001NC

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision No: 5.7  
Revision Date: 10/01/03  
Page B1 of B12**1. SCOPE AND APPLICATION**

This SOP Appendix describes procedures to be used when SW-846 Method 8081A is applied to the analysis of organochlorine pesticides by GC/ECD. This Appendix may also be applied when discontinued SW-846 Method 8080A is requested, and is applicable to extracts derived from any matrix which are prepared according to the appropriate STL North Canton sample extraction SOPs. (CORP-OP-0001NC)

Table B-1 lists compounds, which are routinely determined by this method, and gives the Reporting Limits (RL) for each matrix. RLs given are based on the low level standard and the sample preparation concentration factors. Matrix interferences may result in higher RLs than those listed.

At client request, this method may also be used for the analysis of PCBs (Arochlors) in combination with pesticides, although these are normally analyzed following method 8082, as described in Appendix C of this SOP. In any event, if samples for PCB analysis do not need the acid clean up procedure, then the same injection may be used for method 8081B and 8082, assuming all calibration and QC requirements for both methods are met. Extracts that have been acid cleaned may not be analyzed for pesticides, since several of the pesticides will be degraded.

- 1.1. The associated LIMS method code is QJ (8081A).

**2. SUMMARY OF METHOD**

This method presents conditions for the analysis of prepared extracts of organochlorine pesticides. The pesticides are injected onto the column and separated and detected by electron capture detection. Quantitation may be by internal or external standard methods

**3. DEFINITIONS**

Refer to the STL North Canton Laboratory Quality Manual (LQM), current version, for definitions of terms used in this document.

**4. INTERFERENCES**

- 4.1. Refer to the method 8000B section of this SOP for information regarding chromatographic interferences.
- 4.2. Interferences in the GC analysis arise from many compounds amenable to gas chromatography that give a measurable response on the electron capture detector. Phthalate esters, which are common plasticizers, can pose a major problem in the determinations. Avoiding contact with any plastic materials minimizes interferences from phthalates.
- 4.3. Sulfur will interfere and can be removed using procedures described in SOP CORP-OP-0001NC.
- 4.4. Interferences co-extracted from samples will vary considerably from source to source. The presence of interferences may raise quantitation limits for individual samples. Specific cleanups may be performed on the sample extracts, including florisil cleanup (Method 3620), Gel Permeation Chromatography (Method 3640), and Sulfur cleanup (Method 3660). These cleanup procedures are included in SOP # CORP-OP-0001NC. Using hexane / acetone as the extraction solvent (rather than hexane / methylene chloride) will reduce the amount of interferences extracted

**5. SAFETY**

- 5.1. Refer to section 5 of the Method 8000B SOP for general safety requirements.



## APPENDIX B

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03

Page B2 of B12

**ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081A**

- 
- 5.2. Aroclors have been classified as a potential carcinogen under OSHA. Concentrated solutions of Aroclors must be handled with extreme care to avoid excess exposure. Contaminated gloves and clothing must be removed immediately. Contaminated skin surfaces must be washed thoroughly.
- 5.3. The following parameters covered by this method have been tentatively classified as known or suspected, human or mammalian carcinogens: 4,4'-DDT, 4,4'-DDD, and the BHCs. Primary standards of these toxic compounds should be prepared in a hood.
- 5.4. All  $^{63}\text{Ni}$  sources shall be leak tested every six months, or in accordance with the manufacturer's general radioactive material license.
- 5.5. All  $^{63}\text{Ni}$  sources shall be inventoried every six months. If a detector is missing, the Director, EH&S shall be immediately notified and a letter sent to the NRC or local state agency
- 6. EQUIPMENT AND SUPPLIES**
- 6.1. Refer to Section 6 of the 8000B section of this SOP. A  $^{63}\text{Ni}$  electron capture detector is required
- 6.2. Refer to Table B-2 for analytical columns.
- 6.3. Microsyringes, various sizes, for standards preparation, sample injection, and extract dilution.
- 7. REAGENTS AND STANDARDS**
- 7.1. Refer to the method 8000B section of this SOP for general requirements for reagents and supplies.
- 7.2. Refer to Table B-3 for details of calibration standards.
- 7.3. Surrogate Standards  
Tetrachloro-m-xylene and decachlorobiphenyl are the surrogate standards. Refer to tables B-5 and B-6 for details of surrogate standards.
- 7.4. Column Degradation Evaluation Mix  
A mid-level standard containing 4,4'-DDT and Endrin and not containing any of their breakdown products must be prepared for evaluation of degradation of these compounds by the GC column and injection port. This mix must be replaced after one year, or whenever corrective action to columns fails to eliminate the breakdown of the compounds, whichever is shorter. This solution also contains the surrogates. Refer to Table B-4 for details of the column degradation evaluation mix.
- 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**  
Refer to Section 8 of the 8000B section of this SOP.
- 9. QUALITY CONTROL**  
Refer to Section 9 of the 8000B section of this SOP.
- 10. CALIBRATION AND STANDARDIZATION**
- 10.1. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements

## APPENDIX B

SOP No. CORP-GC-0001NC

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision No: 5.7Revision Date: 10/01/03

Page B3 of B12

- 10.2. Refer to Table B-2 for details of GC operating conditions. The conditions listed should result in resolution of all analytes listed in Table B-1 on both columns. Closely eluting pairs are DDE and Dieldrin on the Rtx-5 or DB-5 column and Endosulfan II and DDD on the 1701 column.
- 10.3. Column Degradation Evaluation
- Before any calibration runs, either initial or 12 hour, The column evaluation mix must be injected before each initial or daily calibration. The degradation of DDT and endrin must be calculated (see equations 9 and 10) and each shown to be less than 15% before calibration can proceed. This is only necessary if the target compound list includes DDT, Endrin, or any of their degradation products.
- If the breakdown of DDT and/or endrin exceeds the limits given above, corrective action must be taken. This action may include:
- Replacement of the injection port liner or the glass wool.
  - Cutting off a portion of the injection end of a capillary column.
  - Replacing the GC column.
- 10.4. Initial Calibration
- Refer to Section 10 of the 8000B section of this SOP for details of calibration procedures.
- 10.4.1. Refer to Table B-7 for the initial calibration analytical sequence.
- 10.4.2. The response for each single-peak analyte will be calculated by the procedures described in the general method for GC analysis.
- 10.4.3. The surrogate calibration curve is calculated from the Individual AB mix. Surrogates in the other calibration standards are used only as retention time markers. If there are resolution problems, then the A and B mixes may be analyzed separately
- 10.4.4. For multi-component pesticides:
- Single point calibration is used for multi-component pesticides (typically toxaphene and technical chlordane). Two options are possible; the same quantitation option must be used for standards and samples. Refer to section 12.3 for guidance on which option to use.
- 10.4.5. For multi-component analytes, the mid level standard must be analyzed as part of the initial calibration. This single point calibration is used to quantitate multi-component analytes.
- 10.4.6. The analyst may include a full 5 point calibration for any of the multi-component analytes with the initial calibration.
- 10.5. 12 hour Calibration Verification
- The 12 hour calibration verification sequence must be analyzed within 12 hours of the start of the initial calibration and at least once every 12 hours thereafter if samples are being analyzed. If more than 12 hours have elapsed since the injection of the last sample in the analytical sequence, a new analytical sequence must be started with a 12 hour calibration. A mid level calibration standard is used for the 12 hour calibration. Refer to the 8000B section of this SOP for acceptance criteria.
- 10.5.1. At a minimum, the 12 hour calibration includes analysis of the breakdown mix followed by mid level standards of any single and multi-component analytes.
- 10.5.2. The retention time windows for any analytes included in the 12 hour calibration are updated.

**ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081A**Revision No: 5.7Revision Date: 10/01/03

Page B4 of B12

**10.6. Continuing Calibration**

The AB calibration mix is analyzed as the continuing calibration standard. At a minimum, this is analyzed after every 20 samples, including matrix spikes, LCS, and method blanks. If 12 hours elapse analyze the 12 hour standard sequence instead. The continuing calibration standard need not include multi-component analytes. If instrument drift is expected due to sample matrix or other factors, it may be advisable to analyze the continuing calibration standard more frequently.

10.6.1. A mid level calibration standard is used for the continuing calibration.

**11. PROCEDURE**

11.1. Refer to the method 8000B section of this SOP for general procedural requirements.

**11.2. Extraction**

The extraction procedure is described in SOP No. CORP-OP-0001NC.

**11.3. Cleanup**

Cleanup procedures are described in SOP No. CORP-OP-0001NC.

11.4. Suggested gas chromatographic conditions are given in Table B-2.

11.5. Allow extracts to warm to ambient temperature before injection.

11.6. The suggested analytical sequence is given in Table B-7.

**12. DATA ANALYSIS AND CALCULATIONS**

12.1. Refer to the 8000B section of this SOP for identification and quantitation of single component analytes.

**12.2. Identification of Multi-component Analytes**

Retention time windows are also used for identification of multi-component analytes, but the "fingerprint" produced by major peaks of those compounds in the standard is used in tandem with the retention times to identify the compounds. The ratios of the areas of the major peaks are also taken into consideration. Identification of these compounds may be made even if the retention times of the peaks in the sample fall outside of the retention time windows of the standard, if in the analyst's judgment the fingerprint (retention time and peak ratios) resembles the standard chromatogram.

**12.3. Quantitation of Multi-component Analytes**

Use 3-10 major peaks or total area for quantitation as described in section 10.4.4, initial calibration of multi-component analytes.

12.3.1. If there are no interfering peaks within the envelope of the multi-component analyte, the total area of the standards and samples may be used for quantitation. Any surrogate or extraneous peaks within the envelope must be subtracted from the total area.

**12.3.1.1. Multiple peak option**

This option is particularly valuable if toxaphene is identified but interferences make quantitation based on total area difficult. Select 3-10 major peaks in the analyte pattern. Calculate the response using the total area or total height of these peaks. Alternatively, find the response of each of the 3-10 peaks per multi-peak pesticide, and use these responses independently, averaging the resultant concentrations found in samples for a final concentration result. When using this option, it is appropriate to remove peaks that appear to be

## APPENDIX B

SOP No. CORP-GC-0001NC

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision No: 5.7  
Revision Date: 10/01/03  
Page B5 of B12

coeluting with contaminant peaks from the quantitation. (i.e. peaks which are significantly larger than would be expected from the rest of the pattern.)

Chlordane may be quantitated either using the multiple peak option (12.3.1.1) total area option (12.3.1.2.) or by quantitation of the major components,  $\alpha$ -chlordane,  $\gamma$ -chlordane and heptachlor

## 12.3.1.2. Total area option

The total area of the standards and samples may be used for quantitation of multi-component analytes. Any surrogate or extraneous peaks within the envelope must be subtracted from the total area. This option should not be used if there are significant interference peaks within the multi-component pattern in the samples. The retention time window for total area measurement must contain at least 90% of the area of the analyte.

- 12.4. Second column confirmation multi-component analytes will only be performed when requested by the client, because the appearance of the multiple peaks in the sample usually serves as a confirmation of analyte presence.
- 12.5. Surrogate recovery results are calculated and reported for decachlorobiphenyl (DCB) and tetrachloro-m-xylene (TCMX). Corrective action is only necessary if DCB and TCMX are both outside of acceptance limits.
- 12.6. Calculation of Column Degradation/% Breakdown (%B)

Equation 9

$$DDT \%B = \frac{A_{DDD} + A_{DDE}}{A_{DDD} + A_{DDE} + A_{DDT}} \times 100$$

where:

$A_{DDD}$ ,  $A_{DDE}$ , and  $A_{DDT}$  = the response of the peaks for 4,4'-DDD, 4,4'-DDE, and 4,4'-DDT in the column degradation evaluation mix.

Equation 10

$$Endrin \%B = \frac{A_{EK} + A_{EA}}{A_{EK} + A_{EA} + A_E} \times 100$$

where:

$A_{EK}$ ,  $A_{EA}$ , and  $A_E$  = the response of endrin ketone, endrin aldehyde, and endrin in the column degradation evaluation mix.

## 13. METHOD PERFORMANCE

- 13.1. Performance limits for the four replicate initial demonstration of capability required under Section 13.1 of the main body of this SOP. Example performance limits are listed in Table B-8. The spiking level should be equivalent to a mid level calibration.

## 14. POLLUTION PREVENTION

Refer to section 14 of the 8000B section of this SOP.

## APPENDIX B

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03Page B6 of B12

---

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081A

## 15. WASTE MANAGEMENT

- 15.1. Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

## 16. REFERENCES

- 16.1. SW846, Update III, December 1996, Method 8081A

## 17. MISCELLANEOUS

- 17.1. Modifications from Reference Method  
None
- 17.2. Modifications from Previous Revisions
- 17.2.1. No revisions were made to this appendix

## APPENDIX B

SOP No. CORP-GC-0001NC

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision No: 5.7  
Revision Date: 10/01/03  
Page B7 of B12

## 17.3. Tables

Table B-1 Standard Analyte list and Reporting Limits			
Compound	Reporting Limit, µg/L or µg/kg		
	water	soil	waste
Aldrin	0.05	1.7	50
α-BHC	0.05	1.7	50
β-BHC	0.05	1.7	50
δ-BHC	0.05	1.7	50
γ-BHC (Lindane)	0.05	1.7	50
α-Chlordane	0.05	1.7	50
γ-Chlordane	0.05	1.7	50
Chlordane (technical)	0.5	17	500
4,4'-DDD	0.05	1.7	50
4,4'-DDE	0.05	1.7	50
4,4'-DDT	0.05	1.7	50
Dieldrin	0.05	1.7	50
Endosulfan I	0.05	1.7	50
Endosulfan II	0.05	1.7	50
Endosulfan Sulfate	0.05	1.7	50
Endrin	0.05	1.7	50
Endrin Aldehyde	0.05	1.7	50
Heptachlor	0.05	1.7	50
Heptachlor Epoxide	0.05	1.7	50
Methoxychlor	0.1	3.3	100
Toxaphene	2.0	67	2000
APPENDIX IX ADD ONs			
Diallate	1.0	33	1000
Isodrin	0.1	3.3	100
Chlorobenzilate	0.1	3.3	100
Kepone <sup>1</sup>	1.0	33	1000

<sup>1</sup> Kepone is sometimes requested for analysis by method 8081A. However kepone may produce peaks with broad tails that elute later than the standard by up to a minute (presumably due to hemi-acetal formation). As a result kepone analysis by 8081A is unreliable and not recommended. Analysis by method 8270C is a possible alternative. Note: alpha chlordane, gamma chlordane, and endrin ketone are not required for some projects

The following concentration factors are assumed in calculating the Reporting Limits.

	Extraction Vol.	Final Vol.
Ground water	1000 mL	10 mL
Low-level Soil	30 g	10 mL
High-level soil / waste	1 g	10 mL

877 752

APPENDIX B

SOP No. CORP-GC-0001NC

Revision No: 5.7

ANALYSIS OF ORGANOCHLORINE PESTICIDES

Revision Date: 10/01/03

BASED ON METHOD 8081A

Page B8 of B12

Table B-2	
Parameter	Recommended Conditions
Injection port temp	220°C
Detector temp	325°C
Temperature program	120°C for 1 min, 8.5°C/min to 285°C, , 6 min hold
Column 1	Rtx-CL Pesticides 30m x 0.32mm id, 0.5µm
Column 2	Rtx-35 30m x 0.32 mm id, 0.5µm
Column 3	DB-608, 30m X 0.32 mm, 0.25µm
Injection	2µL
Carrier gas	Helium or Hydrogen
Make up gas	Nitrogen
Y splitter	Restek or J&W or Supelco glass tee

## APPENDIX B

SOP No. CORP-GC-0001NC

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision No: 5.7  
Revision Date: 10/01/03  
Page B9 of B12

Table B-3						
Calibration Levels ng/mL						
	Level 1	Level 2	Level 3	Level 4	Level 5	Level 6 <sup>2</sup>
<b>Individual Mix AB<sup>1</sup></b>						
Aldrin	5	10	25	50	100	200
g-BHC (Lindane)	5	10	25	50	100	200
Heptachlor	5	10	25	50	100	200
Methoxychlor	10	20	50	100	200	400
Dieldrin	5	10	25	50	100	200
Endosulfan I	5	10	25	50	100	200
Endosulfan II	5	10	25	50	100	200
4,4'-DDT	5	10	25	50	100	200
Endrin Aldehyde	5	10	25	50	100	200
Endrin Ketone	5	10	25	50	100	200
β-BHC	5	10	25	50	100	200
δ-BHC	5	10	25	50	100	200
α-BHC	5	10	25	50	100	200
4,4'-DDD	5	10	25	50	100	200
4,4'-DDE	5	10	25	50	100	200
Endosulfan Sulfate	5	10	25	50	100	200
Endrin	5	10	25	50	100	200
α-Chlordane <sup>3</sup>	5	10	25	50	100	200
γ-Chlordane <sup>3</sup>	5	10	25	50	100	200
<b>Multi-component Standards</b>						
Chlordane (Technical)			250 <sup>4</sup>			
Toxaphene			1000 <sup>5</sup>			
Surrogates are included with all the calibration mixes at the following levels:						
Tetrachloro-m-xylene	5	10	25	50	100	200
Decachlorobiphenyl	5	10	25	50	100	200

<sup>1</sup> Standards may be split into an A and B mix if resolution of all compounds on both columns is not obtained.<sup>2</sup> Level 6 is optional and should only be used if linearity can be maintained on the instrument to this level.<sup>3</sup> Compounds may be used in lieu of running a daily technical Chlordane standard for samples that are non-detect for technical Chlordane.<sup>4</sup> This standard may be used for quantitation of technical chlordane between 50 and 1000 ng/mL. If the chlordane is more concentrated, the extract must be diluted and reanalyzed.<sup>5</sup> This standard may be used for quantitation of toxaphene between 200 and 4000 ng/mL. If the toxaphene is more concentrated, the extract must be diluted and reanalyzed.



## APPENDIX B

SOP No. CORP-GC-0001NC

Revision No: 5.7ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision Date: 10/01/03

Page B10 of B12

Table B-4	
Column Degradation Evaluation Mix ng/mL	
Component	Concentration
4,4'-DDT	25
Endrin	25
Tetrachloro-m-xylene (Surrogate)	20
Decachlorobiphenyl (Surrogate)	20

Table B-5			
LCS/Matrix Spike and Surrogate Spike levels µg/L or µg/kg			
	Aqueous	Soil	Waste
gamma BHC (Lindane)	0.20	33.3	200
Aldrin	0.20	33.3	200
Heptachlor	0.20	33.3	200
Dieldrin	0.50	33.3	500
Endrin	0.50	33.3	500
4,4'-DDT	0.50	33.3	500
Tetrachloro-m-xylene (Surrogate)	0.20	33.3	200
Decachlorobiphenyl (Surrogate)	0.20	33.3	200

Table B-6		
LCS/Matrix Spike and Surrogate Spike levels for TCLP µg/L or µg/kg		
	Aqueous	Waste
Heptachlor	5	500
Heptachlor epoxide	5	500
Lindane	5	500
Endrin	5	500
Methoxychlor	10	1000

## APPENDIX B

SOP No. CORP-GC-0001NC

ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision No: 5.7  
Revision Date: 10/01/03  
Page B11 of B12**Table B-7**  
**Suggested Analytical Sequence****Initial Calibration**

Solvent blank (optional)	
Breakdown Mix	
Individual mix AB	All levels
Technical Chlordane	Level 3 <sup>1</sup>
Toxaphene	Level 3 <sup>1</sup>
Up to 20 samples unless 12 hours comes first)	
Solvent blank (optional)	
Individual mix AB	Mid level (Continuing calibration)
Samples	
After 12 hours:	
Breakdown mix	
Individual mix AB	
Any other single component analytes	

Any multi-component analytes

<sup>1</sup> A five point curve for any of the multi-component analytes may be included. If Aroclors are included, a 5 point calibration for Arochlor 1016/1260 should be included with the initial calibration and a single point for the other Aroclors. The mid point 1016/1260 mix is included with the daily calibration (every 12 hours).

Note: A solvent blank or primer may be analyzed at any time during the sequence when highly contaminated samples are expected. A solvent blank or primer may not be analyzed as routine immediately prior to standards.

**12 hour Calibration**

At least every 12 hours, counting from the start of the initial calibration, or from the start of the last daily calibration, the retention time windows must be updated using the Individual mix AB, and the breakdown mix must be run before the continuing calibration.

## APPENDIX B

SOP No. CORP-GC-0001NC

Revision No: 5.7ANALYSIS OF ORGANOCHLORINE PESTICIDES  
BASED ON METHOD 8081ARevision Date: 10/01/03

Page B12 of B12

Table B-8		
Example Performance limits, four replicate initial demonstration of capability		
Compound	Initial demonstration, mean recovery limits	Initial demonstration, RSD limits
Aldrin	46-112	21
alpha-BHC	51-122	24
beta-BHC	61-120	32
delta-BHC	49.5-118.5	36
gamma-BHC	57-116	23
Chlordane	44.8-108.6	20
4,4'-DDD	52-126	28
4,4'-DDE	46-120	27.5
4,4'-DDT	54-137	36
Dieldrin	42.5-124.5	38
Endosulfan I	43-141	24.5
Endosulfan II	78-171	61
Endosulfan Sulfate	62-132	27
Endrin	49-126	37
Heptachlor	57-100	20
Heptachlor Epoxide	43.5-131.5	25.4
Toxaphene	44.4-111.2	20

## APPENDIX C

SOP No. CORP-GC-0001NC

## ANALYSIS OF PCBs BASED ON METHOD 8082

Revision No: 5.7Revision Date: 10/01/03

Page C1 of C8

**1. SCOPE AND APPLICATION**

1.1. This SOP Appendix describes procedures to be used when SW-846 Method 8000B is applied to the analysis of polychlorinated biphenyls (PCB) by GC/ECD. This Appendix is to be applied when SW-846 Method 8082 is requested, and is applicable to extracts derived from any matrix which are prepared according to the appropriate STL sample extraction SOPs. (CORP-OP-0001NC). The PCBs are determined and quantitated as Arochlor mixes.

1.2. Table C-1 lists compounds, which are routinely determined by this method, and gives the Reporting Limits (RL) for each matrix. RLs given are based on the low level standard and the sample preparation concentration factors. Matrix interferences may result in higher RLs than those listed.

**Note: SW-846 method 8082 provides incomplete guidance for determination of individual PCB congeners. This SOP does not include directions for congener specific analysis.**

1.3. The associated LIMS method code is QH (8082).

**2. SUMMARY OF METHOD**

This method presents conditions for the analysis of prepared extracts of PCBs. The PCBs are injected onto the column and separated and detected by electron capture detection. Quantitation is by the external standard method.

**3. DEFINITIONS**

Refer to the STL North Canton Laboratory Quality Manual (LQM), current version, for definitions of terms used in this document.

**4. INTERFERENCES**

4.1. Refer to the method 8000B section of this SOP for information regarding chromatographic interferences.

4.2. Interferences in the GC analysis arise from many compounds amenable to gas chromatography that give a measurable response on the electron capture detector. Phthalate esters, which are common plasticizers, can pose a major problem in the determinations. Avoiding contact with any plastic materials minimizes interferences from phthalates.

4.3. Sulfur will interfere and can be removed using procedures described in SOP CORP-OP-0001NC

4.4. Interferences co-extracted from samples will vary considerably from source to source. The presence of interferences may raise quantitation limits for individual samples. Specific cleanups may be performed on the sample extracts, including florisil cleanup (Method 3620), Gel Permeation Chromatography (Method 3640), and Sulfur cleanup (Method 3660). These cleanup procedures are included in SOP # CORP-OP-0001NC.

**5. SAFETY**

5.1. Refer to section 5 of the Method 8000B SOP for general safety requirements

5.2. Aroclors have been classified as a potential carcinogen under OSHA. Concentrated solutions of Aroclors must be handled with extreme care to avoid excess exposure. Contaminated gloves and clothing must be removed immediately. Contaminated skin surfaces must be washed thoroughly.

5.3. All  $^{63}\text{Ni}$  sources shall be leak tested every six months, or in accordance with the manufacturer's general radioactive material license.

- 5.4. All  $^{63}\text{Ni}$  sources shall be inventoried every six months. If a detector is missing, the Director, EH&S shall be immediately notified and a letter sent to the NRC or local state agency.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Refer to Section 6 of the 8000B section of this SOP. A  $^{63}\text{Ni}$  electron capture detector is required.
- 6.2. Refer to Table C-2 for analytical columns
- 6.3. Microsyringes, various sizes, for standards preparation, sample injection, and extract dilution.

## 7. REAGENTS AND STANDARDS

- 7.1. Refer to the method 8000B section of this SOP for general requirements for reagents and supplies. The standards must be replaced at least every six months or sooner if comparison with check standards indicates a problem
- 7.2. Refer to Table C-3 for details of calibration standards.
- 7.3. Surrogate Standards  
Tetrachloro-m-xylene and decachlorobiphenyl are the surrogate standards. Other surrogates may be used at client request. Refer to Table C-4 for details of surrogate standards

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

Refer to Section 8 of the 8000B section of this SOP.

## 9. QUALITY CONTROL

Refer to Section 9 of the 8000B section of this SOP.

## 10. CALIBRATION AND STANDARDIZATION

- 10.1. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements.
- 10.2. Initial Calibration
- 10.2.1. Refer to Table C-5 for the initial calibration analytical sequence
- 10.2.2. The response for each Arochlor will be calculated by the procedures described in the general method for GC analysis, with the following modifications.
- 10.2.3. A five point calibration of the Arochlor 1016/1260 mix is generated with at least mid level single points for the other Arochlor mixes. The average response factor is used to quantitate Arochlors 1260 and 1016, other Arochlors are quantitated from the mid level single point.
- 10.2.4. The analyst may include a full 5 point calibration for any of the Arochlors with the initial calibration.
- 10.2.5. The high and low standards for the initial 5 point calibration of 1016 / 1260 define the acceptable quantitation range for the other Arochlors. If any Arochlor is determined above this concentration the extract must be diluted and reanalyzed.
- 10.2.6. If the analyst knows that a specific Arochlor is of interest for a particular project, that Arochlor may be used for the five point calibration rather than the 1016 / 1260 mix

## APPENDIX C

SOP No. CORP-GC-0001NC

## ANALYSIS OF PCBs BASED ON METHOD 8082

Revision No: 5.7Revision Date: 10/01/03

Page C3 of C8

10.2.7. The surrogate calibration curve is calculated from the Aroclor 1016/1260 mix. Surrogates in the other calibration standards are used only as retention time markers.

10.2.8. Two options are possible for quantitation of Aroclors. The same quantitation option must be used for standards and samples.

10.2.8.1. Multiple peak option

Select 3-10 major peaks in the analyte pattern. Calculate the response using the total area or total height of these peaks. Alternatively, find the response of each of the 3-10 peaks per Aroclor, and use these responses independently, averaging the resultant concentrations found in samples for a final concentration result. When using this option, it is appropriate to remove peaks that appear to be coeluting with contaminant peaks from the quantitation. (i.e. peaks which are significantly larger than would be expected from the rest of the pattern.)

10.2.8.2. Total area option

The total area of the standards and samples may be used for quantitation of multi-component analytes. Any surrogate or extraneous peaks within the envelope must be subtracted from the total area. This option should not be used if there are significant interference peaks within the multi-component pattern in the samples. The retention time window for total area measurement must contain at least 90% of the area of the analyte.

10.3. 12 hour Calibration

The 12 hour calibration verification must be analyzed within 12 hours of the start of the initial calibration and at least once every 12 hours thereafter if samples are being analyzed. If there is a break in the analytical sequence of greater than 12 hours, then a new continuing calibration run must be analyzed before proceeding with the sequence. If more than 12 hours have elapsed since the injection of the last sample in the analytical sequence, a new analytical sequence must be started with a 12 hour calibration.

10.3.1. At a minimum, the 12 hour calibration includes analysis of the Aroclor 1260 / 1016 mix

10.3.2. It is adequate to verify calibration with a mixture of Aroclors 1016 and 1260. If a specific Aroclor is expected, it may be included in the daily calibration check

10.3.3. The retention time windows for any analytes included in the daily calibration and CCVs are updated.

10.3.4. For this method samples must be bracketed with successful calibration verification runs.

10.4. Calibration verification

The Aroclor 1260/1016 calibration mix is analyzed as the calibration verification standard. This is analyzed after every 20 samples, including matrix spikes, LCS, and method blanks. (Depending on the type of samples, it may be advisable to analyze verifications more frequently in order to minimize reruns.).

10.4.1. A mid level standard is used for the calibration verification

11. PROCEDURE

11.1. Refer to the method 8000B section of this SOP for general procedural requirements.

## APPENDIX C

SOP No. CORP-GC-0001NC

## ANALYSIS OF PCBs BASED ON METHOD 8082

Revision No: 5.7Revision Date: 10/01/03Page C4 of C8

---

## 11.2. Extraction

The extraction procedure is described in SOP No. CORP-OP-0001NC.

## 11.3. Cleanup

Cleanup procedures are described in SOP No. CORP-OP-0001NC.

## 11.4. Suggested gas chromatographic conditions are given in Table C-2.

## 11.5. Allow extracts to warm to ambient temperature before injection.

## 11.6. The suggested analytical sequence is given in Table C-5.

**12. DATA ANALYSIS AND CALCULATIONS**

## 12.1. Identification of Arochlors

Retention time windows are used for identification of Arochlors, but the "fingerprint" produced by major peaks of those analytes in the standard is used in tandem with the retention times for identification. The ratios of the areas of the major peaks are also taken into consideration. Identification may be made even if the retention times of the peaks in the sample fall outside of the retention time windows of the standard, if in the analyst's judgment the fingerprint (retention time and peak ratios) resembles the standard chromatogram.

A clearly identifiable Arochlor pattern serves as confirmation of single column GC analysis. Dual column confirmation may be used for specific program requirements or by client request.

## 12.2. Quantitation of Arochlors

Use 3-10 major peaks or total area for quantitation

If the analyst believes that a combination of Arochlor 1254 and 1260, or a combination of 1242, 1248 and 1232 is present, then only the predominant Arochlor is quantitated and reported, but the suspicion of multiple Arochlors is discussed in the narrative. If well separated Arochlor patterns are present, and then both Arochlors are quantitated and reported.

## 12.3. If there are no interfering peaks within the envelope of the Arochlor, the total area of the standards and samples may be used for quantitation. Any surrogate or extraneous peaks within the envelope must be subtracted from the total area.

## 12.4. Second column confirmation of Arochlors will only be performed when requested by the client. The appearance of the multiple peaks in the sample usually serves as a confirmation of Arochlor presence.

## 12.5. Surrogate recovery results are calculated and reported for decachlorobiphenyl (DCB) and tetrachloro-m-xylene (TCMX). Corrective action is only necessary if DCB and TCMX are both outside of acceptance limits

**13. METHOD PERFORMANCE**

## 13.1. Performance limits for the four replicate initial demonstration of capability are required as referenced under Section 13.1 of the main body of this SOP.

## 13.2. Method detection limits (MDL) are determined for all Arochlors.

**14. POLLUTION PREVENTION**

Refer to section 14 of the 8000B section of this SOP

## APPENDIX C

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03

Page C5 of C8

ANALYSIS OF PCBs BASED ON METHOD 8082

---

**15. WASTE MANAGEMENT**

- 15.1. Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**

- 16.1. SW846, Update III, December 1996, Method 8082

**17. MISCELLANEOUS**

- 17.1. Modifications from Reference Method

- 17.1.1. Method 8082 includes limited direction for congener specific quantitation. This is outside the scope of this SOP.

- 17.2. Modifications from Previous Revisions

No changes were made to this Appendix



## APPENDIX C

SOP No. CORP-GC-0001NC

## ANALYSIS OF PCBs BASED ON METHOD 8082

Revision No: 5.7

Revision Date: 10/01/03

Page C6 of C8

## 17.3. Tables

Table C-1 Standard Analyte list and Reporting Limits			
Compound	Reporting Limit, µg/L or µg/kg		
	water	soil	waste
Aroclor-1016	1.0	33	1000
Aroclor-1221	1.0	33	1000
Aroclor-1232	1.0	33	1000
Aroclor 1242	1.0	33	1000
Aroclor-1248	1.0	33	1000
Aroclor-1254	1.0	33	1000
Aroclor-1260	1.0	33	1000

The following concentration factors are assumed in calculating the Reporting Limits:

	Extraction Vol	Final Vol
Ground water	1000 mL	10 mL
Low-level Soil	30 g	10 mL
High-level soil / waste	1 g	10 mL

Table C-2	
Parameter	Recommended Conditions
Injection port temp	220°C
Detector temp	325°C
Temperature program	70°C for 0.5 min, 30°C/min to 190°C, 2.5°C/min to 225, 18°C/min to 280°C, 3 min hold
Column 1	DB-5 or Rtx-5 30m x 0.32mm id, 0.5µm
Column 2	DB-1701 or Rtx 1701 30m x 0.32 mm id, 0.25µm
Column 3	DB-608, 30m X 0.32 mm, 0.25µm
Injection	1-2µL
Carrier gas	Helium or Hydrogen
Make up gas	Nitrogen
Y splitter	Restek or J&W or Supelco glass tee

## APPENDIX C

SOP No. CORP-GC-0001NC

## ANALYSIS OF PCBs BASED ON METHOD 8082

Revision No: 5.7

Revision Date: 10/01/03

Page C7 of C8

Table C-3						
Calibration Levels ng/mL						
	Level 1	Level 2	Level 3	Level 4	Level 5	Level 6 <sup>1</sup>
Aroclor 1016/1260	100	200	500	1000	2000	4000
Aroclor 1242 <sup>2</sup>			500			
Aroclor 1221 +1254 <sup>2</sup>			500			
Aroclor 1232 <sup>2</sup>			500			
Aroclor 1248 <sup>2</sup>			500			
Surrogates are included with all the calibration mixes at the following levels:						
Tetrachloro-m-xylene	5	10	25	50	100	200
Decachlorobiphenyl	5	10	25	50	100	200

<sup>1</sup> Level 6 is optional and should only be used if linearity can be maintained on the instrument to this level.<sup>2</sup> Aroclors may be quantitated within the range 100 to 2000 ng/mL (4000ng/mL if the level 6 1016/1260 standard is included). If the Aroclor is more concentrated, it must be reanalyzed at a dilution.

Table C-4			
LCS/Matrix Spike and Surrogate Spike levels for Aroclor analysis with Acid Cleanup			
µg/L or µg/kg			
	Aqueous	Soil	Waste
Aroclor 1016/1260	10	333	10,000
Tetrachloro-m-xylene (Surrogate)	0.20	6.67	200
Decachlorobiphenyl (Surrogate)	0.20	6.67	200

Table C-5			
Michigan Analyte List and Reporting Limits <sup>1</sup>			
Compound	Reporting Limit		
	water (µg/L)	soil (µg/Kg)	
Aroclor-1016	0.2	330	
Aroclor-1221	0.2	330	
Aroclor-1232	0.4	330	
Aroclor 1242	0.2	330	
Aroclor-1248	0.2	330	
Aroclor-1254	0.2	330	
Aroclor-1260	0.2	330	

<sup>1</sup> Reporting Limits are only for samples performed under the Michigan program

## APPENDIX C

SOP No. CORP-GC-0001NC

## ANALYSIS OF PCBs BASED ON METHOD 8082

Revision No: 5.7Revision Date: 10/01/03Page C8 of C8

---

Table C-5

## Suggested Analytical Sequence

## Initial Calibration

## Injection #

1	Solvent blank (optional)	
2	Aroclor 1016/1260	Level 1
3	Aroclor 1016/1260	Level 2
4	Aroclor 1016/1260	Level 3
5	Aroclor 1016/1260	Level 4
6	Aroclor 1016/1260	Level 5
7	Aroclor 1232	Level 3
8	Aroclor 1242	Level 3
9	Aroclor 1248	Level 3
10	Aroclor 1221/1254	Level 3
11-30	Sample 1-20 (or as many samples as can be analyzed in 12 hours)	
	Solvent blank (optional)	
32	Aroclor 1016/1260	Level 3

etc

Note: A solvent blank or primer may be analyzed at any time during the sequence when highly contaminated samples are expected. A solvent blank or primer may not be analyzed as routine immediately prior to standards.

## 12 hour Calibration

At least every 12 hours, counting from the start of the initial calibration, or from the start of the last daily calibration, the retention time windows must be updated using the Aroclor 1260 / 1016 mix. Mid level standards of any other Aroclors expected to be present in the samples are also injected.

## APPENDIX D

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03

Page D1 of D5

**ANALYSIS OF PHENOXY ACID HERBICIDES BASED ON  
SW-846 METHOD 8151A and 615****1. SCOPE AND APPLICATION**

This method is applicable to the gas chromatographic determination of Chlorinated phenoxy acid herbicides in extracts prepared by SOP CORP-OP-0001NC. The herbicides listed in Table D1 are routinely analyzed. Other chlorinated acids may be analyzed by this method if the quality control criteria in Section 9 and the initial demonstration of method performance in Section 13 are met.

- 1.1. The associated LIMS method code is QS.

**2. SUMMARY OF METHOD**

This method presents conditions for the analysis of prepared extracts of phenoxy acid herbicides by gas chromatography. The herbicides, as their methyl esters, are injected onto the column, separated, and detected by electron capture detectors. Quantitation is by the external standard method.

**3. DEFINITIONS**

Refer to the STL North Canton Laboratory Quality Manual (LQM), current version, for definitions of terms used in this document.

**4. INTERFERENCES**

- 4.1. Refer to the method 8000B section of this SOP for general information regarding chromatographic interferences.
- 4.2. Chlorinated acids and phenols cause the most direct interference with this method.
- 4.3. Sulfur may interfere and may be removed by the procedure described in SOP#CORP-OP-0001NC.

**5. SAFETY**

- 5.1. Refer to section 5 of the Method 8000B SOP for general safety requirements.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Refer to Section 6 of the 8000B section of this SOP. A Ni<sub>63</sub> electron capture detector is required.
- 6.2. Refer to Table D2 for analytical columns.
- 6.3. Microsyringes, various sizes, for standards preparation, sample injection, and extract dilution.

**7. REAGENTS AND STANDARDS**

- 7.1. Refer to section 7 of the 8000B section of this SOP for general information on reagents and standards.

**8. SAMPLE PREPARATION, PRESERVATION AND STORAGE**

Refer to Section 8 of the 8000B section of this SOP.

**9. QUALITY CONTROL**

- 9.1. Refer to Section 9 of the 8000B section of this SOP for quality control requirements, including the initial demonstration of capability, definition of a batch, surrogate limits, method blanks, laboratory control samples (LCS), and matrix spikes (MS).

## APPENDIX D

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03

Page D2 of D5

**ANALYSIS OF PHENOXY ACID HERBICIDES BASED ON  
SW-846 METHOD 8151A and 615**

- 9.2. Refer to Table D-3 for the components and levels of the LCS and MS mixes.

**10. CALIBRATION AND STANDARDIZATION**

- 10.1. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements.
- 10.2. Calibration standards are prepared from purchased standards in the methyl ester form.
- 10.3. The low level standard must be at or below the laboratory reporting limit. Other standards are chosen to bracket the expected range of concentrations found in samples, without saturating the detector or leading to excessive carryover.
- 10.4. Refer to Table D-2, for details of GC operating conditions.

**11. PROCEDURE**

- 11.1. Refer to the method 8000B section of this SOP for procedural requirements.
- 11.2. Extraction  
The extraction procedure is described in SOP #CORP-OP-0001NC
- 11.3. Cleanup  
The alkaline hydrolysis and subsequent extraction of the basic solution described in the extraction procedure provides an effective cleanup.
- 11.4. Analytical Sequence  
The analytical sequence starts with an initial calibration of at least five points, or a daily calibration that meets % difference criteria from an existing initial calibration.
- 11.4.1. The daily calibration must be analyzed at least once every 24 hours when samples are being analyzed. If there is a break in the analytical sequence of greater than 12 hours, then a new continuing calibration run must be analyzed before proceeding with the sequence. If more than 24 hours have elapsed since the injection of the last sample in the analytical sequence, a new analytical sequence must be started with a daily calibration.
- 11.4.2. The daily calibration consists of mid level standards of all analytes of interest. Retention time windows must be updated with the daily calibration
- 11.4.3. After every 12 hours a continuing calibration is analyzed. The continuing calibration consists of mid level standards of all analytes of interest. Retention time windows are updated with continuing calibrations
- 11.5. Gas Chromatography  
Chromatographic conditions are listed in Table D-2.

**12. DATA ANALYSIS AND CALCULATIONS**

- 12.1. Refer to the 8000B section of this SOP for identification and quantitation of single component analytes.
- 12.2. The herbicides are analyzed as their methyl esters, but reported as the free acid. For this reason it is necessary to correct the results for the molecular weight of the ester versus the free acid. This is achieved through the concentrations of the calibration standards. For example the 20µg/L calibration standard for 2,4-D contains 21.3 µg/L of the methyl ester. No further correction is necessary.

**APPENDIX D**

SOP No CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03

Page D3 of D5

**ANALYSIS OF PHENOXY ACID HERBICIDES BASED ON  
SW-846 METHOD 8151A and 615**

---

- 12.3. A routine 10X dilution occurs on final extracts for all samples. Due to a QuantIMS limitation, the dilution factor field in QuantIMS cannot be used when a dilution is routine, because the dilution factor is automatically applied to all reference values creating reporting problems. For the herbicide analysis, the extract volume will be 10mL and an aliquot at 10X dilution will be analyzed. The final extract volume recorded on the laboratory bench sheet will be recorded as 100mL to avoid using the dilution factor field in QuantIMS.

**13. METHOD PERFORMANCE**

- 13.1. The EPA for this method has not published multiple laboratory performance data. Performance limits for the four replicate initial demonstration of capability are required as referenced under Section 13.1 of the main body of this SOP.

**14. POLLUTION PREVENTION**

This method does not contain any specific modifications that serve to minimize or prevent pollution

**15. WASTE MANAGEMENT**

Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**

Method 8151A, SW-846, Update III, December 1996

**17. MISCELLANEOUS**

- 17.1 Modifications from Reference Method

Refer to the method 8000B section of this SOP for modifications from the reference method

- 17.2 Modifications from Previous Revision

The calibration procedure has been changed to require esterification of the calibration standards

## APPENDIX D

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page D4 of D5

ANALYSIS OF PHENOXY ACID HERBICIDES BASED ON  
SW-846 METHOD 8151A and 615

## 17.3. Tables

Table D-1				
Standard Analyte list				
Compound	CAS Number	Reporting Limit, µg/L or µg/kg		
		Aqueous	Soil	Waste
2,4-D	94-75-7	4	80	4000
2,4-DB	94-82-6	4	80	4000
2,4,5-TP (Silvex)	93-72-1	1	20	1000
2,4,5-T	93-76-5	1	20	1000
Dalapon	75-99-0	2	40	2000
Dicamba	1918-00-9	2	40	2000
Dichloroprop	120-36-5	4	80	4000
Dinoseb	88-85-7	0.6	12	600
MCPA	94-74-6	400	8000	400,000
MCP	93-65-2	400	8000	400,000

The following concentration factors are assumed in calculating the Reporting Limits.

	Extraction Vol.	Final Vol.	Dilution Factor
Ground water	1000 mL / 10 mL	20	
Low-level Soil without GPC	50 g	10 mL	20
High-level soil / waste	1 g	10 mL	20

Specific reporting limits are highly matrix dependent. The reporting limits listed above are provided for guidance only and may not always be achievable. For special projects, the extracts may be analyzed without any dilution, resulting in reporting limits 20 times lower than those in Table D-1.

Table D-2	
Instrumental Conditions	
PARAMETER	Recommended conditions
Injection port temp	220°C
Detector temp	325°C
Temperature program	80,2/30/170,0/1/180,1
Column 1	DB-5MS or RTX 5 30x0.32, 0.5µm
Column 2	DB-1701 or Rtx-1701
Injection	1-2µL
Carrier gas	Helium / Hydrogen
Make up gas	Nitrogen

Recommended conditions should result in resolution of all analytes listed in Table D-1.

The reporting limits listed in Table D-1 will be achieved with these calibration levels and a 20 fold dilution of the sample extract. Lower reporting limits can be achieved with lesser dilutions of the sample extract.

## APPENDIX D

SOP No. CORP-GC-0001NC

ANALYSIS OF PHENOXY ACID HERBICIDES BASED ON  
SW-846 METHOD 8151A and 615Revision No: 5.7Revision Date: 10/01/03

Page D5 of D5

Table D-3			
LCS/Matrix Spike and Surrogate Spike levels $\mu\text{g/L}$ or $\mu\text{g/kg}$ <sup>1</sup>			
	Aqueous	Soil	Waste
2,4-D	16	800	16000
Silvex	4	200	4000
2,4,5-T	4	200	4000
2,4-DB	16	800	16000
Dalapon	8	400	8000
DCAA (surrogate)	16	800	16000

<sup>1</sup> LCS, MS and SS spikes are as the free acid.



APPENDIX E  
ORGANOPHOSPHORUS PESTICIDES BASED ON METHOD 8141A

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page E1 of E6

---

1. SCOPE AND APPLICATION

- 1.1. This method is applicable to the determination of the concentration of certain Organophosphorus Pesticides in waters, wastewaters, oils, soils, and sludges. It is based on SW846 Method 8141A. Table E1 shows reporting limits for compounds routinely analyzed by this method. The compounds include, but are not limited to, those shown in Table E1.
- 1.2. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary by the laboratory QA department.
- 1.3. The associated LIMS method code is P2.

2. SUMMARY OF METHOD

- 2.1. An aliquot of prepared sample is injected in a gas chromatograph (GC) and compounds in the effluent are detected by a flame photometric detector. Appropriate preparation techniques are described in SOP CORP-OP-0001NC. Ultrasonic Extraction (Method 3550) is **NOT** an appropriate sample preparation for Method 8141 and should not be used because of the potential for destruction of target analytes during the ultrasonic extraction process.

3. DEFINITIONS

Refer to the STL North Canton Laboratory Quality Manual (LQM), current version, for definitions of terms used in this document.

4. INTERFERENCES

- 4.1. Refer to the method 8000B section of this SOP for general information regarding chromatographic interferences.
- 4.2. Analytical difficulties encountered for target analysis include:
  - 4.2.1. The water solubility of Dichlorvos (DDVP) is 10 g/L at 20°C, and recovery is poor from aqueous solution.
  - 4.2.2. Naled is converted to Dichlorvos (DDVP) on column by debromination. This reaction may also occur during sample workup.
  - 4.2.3. Trichlorfon rearranges and is dehydrochlorinated in acidic, neutral, or basic media to form Dichlorvos (DDVP) and hydrochloric acid. If this method is to be used for the determination of organophosphates in the presence of Trichlorfon, the analyst should be aware of the possibility of rearrangement to Dichlorvos to prevent misidentification.
  - 4.2.4. Merphos is a single component pesticide that is readily oxidized to Merphos oxone. Chromatographic analysis of Merphos almost always results in two peaks.

5. SAFETY

- 5.1. Refer to Section 5 of the Method 8000B SOP for general safety requirements.

**APPENDIX E**  
**ORGANOPHOSPHORUS PESTICIDES BASED ON METHOD 8141A**

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page E2 of E6

---

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Refer to Section 6 of the 8000B section of this SOP.
- 6.2. Refer to Table E2 for Instrument settings.
- 6.3. Microsyringes, various sizes, for standards preparation, sample injection, and extract dilution.

**7. REAGENTS AND STANDARDS**

- 7.1. Refer to Section 7 of the 8000B section of this SOP for general information on reagents and standards.
- 7.2. Refer to Table E-3 for details of calibration and other standards.
- 7.3. Surrogate Standards  
Triphenyl phosphate is the surrogate standard. Refer to Table E-4 for details of the surrogate standard.

**8. SAMPLE PREPARATION, PRESERVATION AND STORAGE**

Refer to section 8 of the 8000B section of this SOP.

**9. QUALITY CONTROL**

- 9.1. Refer to Section 9 of the 8000B section of this SOP for quality control requirements, including the initial demonstration of capability, definition of a batch, surrogate limits, method blanks, laboratory control samples (LCS), and matrix spikes (MS).
- 9.2. Refer to Table E-4 for the components and levels of the LCS and MS mixes.

**10. CALIBRATION AND STANDARDIZATION**

- 10.1. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements.
- 10.2. Calibration standards are made up from purchased solutions. Table E-3 lists the calibration levels.
- 10.3. The low level standard must be at or below the laboratory reporting limit. Other standards are chosen to bracket the expected range of concentrations found in samples, without saturating the detector or leading to excessive carryover.

**11. PROCEDURE**

- 11.1. Refer to the method 8000B section of this SOP for procedural requirements.
- 11.2. Extraction  
The extraction procedure is described in SOP #CORP-OP-0001NC.
- 11.3. Cleanup
- 11.4. Analytical Sequence  
The analytical sequence starts with an initial calibration of at least five points, or a daily calibration that meets % difference criteria from an existing initial calibration.

**APPENDIX E**  
**ORGANOPHOSPHORUS PESTICIDES BASED ON METHOD 8141A**

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page E3 of E6

---

- 11.4.1 The daily calibration must be analyzed at least once every 24 hours when samples are being analyzed. If there is a break in the analytical sequence of greater than 12 hours, then a new continuing calibration run must be analyzed before proceeding with the sequence. If more than 24 hours have elapsed since the injection of the last sample in the analytical sequence, a new analytical sequence must be started with a daily calibration.
- 11.4.2 The daily calibration consists of mid level standards of all analytes of interest. Retention time windows must be updated with the daily calibration.
- 11.4.3 After every 12 hours a continuing calibration is analyzed. The continuing calibration consists of mid level standards of all analytes of interest. Retention time windows are updated with continuing calibrations.

**12. DATA ANALYSIS AND CALCULATIONS**

- 12.1 Refer to the 8000B section of this SOP for identification and quantitation of single component analytes.
- 12.2 Surrogate recovery results are calculated and reported for Triphenylphosphate unless it is determined that sample interference has adversely affected the quantitation of one of the surrogates. The surrogate must be within QC criteria. Corrective action is only necessary if Triphenylphosphate is outside of acceptance limits.

**13. METHOD PERFORMANCE**

- 13.1 Performance limits for the four replicate initial demonstration of capability are required as referenced under Section 13.1 of the main body of this SOP

**14. POLLUTION PREVENTION**

Refer to section 14 of the 8000B section of this SOP.

**15. WASTE MANAGEMENT**

- 15.1 Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required

**16. REFERENCES**

SW846, Update III, December 1996, Method 8141A.

**17. MISCELLANEOUS**

- 17.1 Reporting limits
  - 17.1.1 The lower standard reporting limits are listed in Table E-1
  - 17.1.2 If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.
    - 17.1.2.1 The nature of the FPD detector contributes to high dilutions for Method 8141A. There is a phenomenon known as quenching that occurs. This happens when light absorption occurs in the flame of the FPD due to hydrocarbons, sulfur, and certain light absorbing compounds. When this happens the analytes of interest do not reach the photomultiplier tube and are not detected even though they may be present.

**APPENDIX E**  
**ORGANOPHOSPHORUS PESTICIDES BASED ON METHOD 8141A**

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03

Page E4 of E6

## 17.2. Elution time

17.2.1. For any positive detection in a sample, the chromatogram is overlaid on screen with the nearest standard for elution time matching and pattern recognition using the Target system.

## 17.3. Troubleshooting guide

17.3.1. Consult the instrument manufacturer's operating manual for guidance.

Table E1: Organophosphorus Pesticides Routinely Analyzed and Reporting Limits

Compound	CAS Number	Reporting Limits	
		Water, µg/L	Solid, µg/kg
Azinphos methyl	86-50-0	1.0	33
Bolstar (Suprofos)	35400-43-2	1.0	33
Chlorpyrifos	2921-88-2	1.0	33
Coumaphos	56-72-4	1.0	33
Demeton, O and S	8065-48-3	1.0	33
Diazinon	333-41-5	1.0	33
Dichlorvos	62-73-7	1.0	33
Disulfoton	298-04-4	1.0	33
Ethoprop	13194-48-4	1.0	33
Fensulfothion	115-90-2	1.0	33
Fenthion	55-38-9	1.0	33
Malathion	121-75-5	1.0	33
Merphos	150-50-5	1.0	33
Methyl Parathion	298-00-0	1.0	33
Mevinphos	7786-34-7	1.0	33
Naled	300-76-5	1.0	33
Phorate	298-02-2	1.0	33
Ronnel	299-84-3	1.0	33
Stirophos	22248-79-9	1.0	33
Tokuthion	34643-46-4	1.0	33
Trichloronate	327-98-0	1.0	33
O,O,O-Tricntyl phosphorothioate	126-68-1	1.0	33
Thionazin	297-97-2	1.0	33
Sulfotepp	3689-24-5	1.0	33
Dimethoate	60-51-5	1.0	33
Parathion	56-38-2	1.0	33
Famphur	52-85-7	1.0	33

**APPENDIX E**  
**ORGANOPHOSPHORUS PESTICIDES BASED ON METHOD 8141A**

SOP No. CORP-GC-0001NC

Revision No: 5.7Revision Date: 10/01/03Page E5 of E6

---

Table E2: Instrumental Conditions	
PARAMETER	Recommended conditions
Injection port temp	175°C
Detector temp	230°C
Initial temp	135°C
Temperature program	(A) 5°C/minute (B) 20°C/minute
Final Temp	(A) 245°C (B) 295°C
Final Hold Time	(A) 0 minutes (B) 7 minutes
Column 1	RTX-OPP, 30 meter, 0.32 mm, 0.5 µm film
Column 2	RTX-1, 30 meter, 0.32 mm, 0.5 µm film
Injection	1-2µL
Carrier gas	Helium / Hydrogen
Make up gas	Nitrogen

APPENDIX E  
ORGANOPHOSPHORUS PESTICIDES BASED ON METHOD 8141A

SOP No. CORP-GC-0001NC

877 775

Revision No: 5.7  
Revision Date: 10/01/03  
Page E6 of E6

Table E3 Initial Calibration Concentrations (ng/uL)

COMPOUND	LEVEL1	LEVEL2	LEVEL3	LEVEL4	LEVEL5	LEVEL6	LEVEL7
o,o,o-Triethylphosphate	0.2	0.5	1	2	5	10	20
Dichlorvos	0.2	0.5	1	2	5	10	20
Mevinphos	0.2	0.5	1	2	5	10	20
Thionazin	0.2	0.5	1	2	5	10	20
Ethoprop	0.2	0.5	1	2	5	10	20
Naled	0.2	0.5	1	2	5	10	20
Sulfotepp	0.2	0.5	1	2	5	10	20
Phorate	0.2	0.5	1	2	5	10	20
Demeton	0.2	0.5	1	2	5	10	20
Dimethoate	0.2	0.5	1	2	5	10	20
Diazinon	0.2	0.5	1	2	5	10	20
Disulfoton	0.2	0.5	1	2	5	10	20
Methyl Parathion	0.2	0.5	1	2	5	10	20
Rommel	0.2	0.5	1	2	5	10	20
Fenthion	0.2	0.5	1	2	5	10	20
Chlorpyrifos	0.2	0.5	1	2	5	10	20
Parathion	0.2	0.5	1	2	5	10	20
Malathion	0.2	0.5	1	2	5	10	20
Trichloronate	0.2	0.5	1	2	5	10	20
Merphos	0.2	0.5	1	2	5	10	20
Stirophos	0.2	0.5	1	2	5	10	20
Tokuthion	0.2	0.5	1	2	5	10	20
Fensulfthion	0.2	0.5	1	2	5	10	20
Bolstar	0.2	0.5	1	2	5	10	20
Famphur	0.2	0.5	1	2	5	10	20
Azinphos methyl	0.2	0.5	1	2	5	10	20
Coumaphos	0.2	0.5	1	2	5	10	20

Table E4: LCS/Matrix Spike and Surrogate Spike Compounds – 20 ug/mL

Compound	Compound
Thiazin	Sulfotepp
Phorate	Disulfoton
Methyl Parathion	Parathion
Famphur	O,O,O-Triethylphosphate
Dimethoate	Triphenyl Phosphate - Surrogate

**APPENDIX F**  
**TOTAL PETROLEUM HYDROCARBONS BASED ON 8015B**

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page F1 of F4

---

**1. Scope and Application**

- 1.1. This method is applicable to the determination of the concentration and **tentative** identification of petroleum hydrocarbon mixes in waters, wastewaters, soils, and sludges.
- 1.2. This SOP is based on SW-846 Method 8015B, Modified, Revision 3, December 1996.
- 1.3. The associated LIMS method codes are HS (8015 MOD) and KI (8015B).

**2. Summary of Method**

- 2.1. This method provides gas chromatographic conditions for detection and identification of total petroleum hydrocarbons. Prior to the use of this method, appropriate sample preparation techniques are used.
- 2.2. An aliquot of the prepared sample is injected into a gas chromatograph (GC) and compounds in the effluent are detected by a flame ionization detector (FID)
- 2.3. The carbon range for Ohio VAP and BUSTR projects is C10-C20 and C20-C34.

**3. Definitions**

- 3.1. Refer to the glossary in the STL North Canton Laboratory Quality Manual (LQM), current version.

**4. INTERFERENCES**

- 4.1. Refer to the method 8000B section of this SOP for general information regarding chromatographic interferences

**5. SAFETY**

- 5.1. Refer to Section 5 of the Method 8000B SOP for general safety requirements.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Refer to Section 6 of the 8000B section of this SOP.
- 6.2. Microsyringes, various sizes, for standards preparation, sample injection, and extract dilution.

**7. REAGENTS AND STANDARDS**

- 7.1. Refer to Section 7 of the 8000B section of this SOP.
- 7.2. The petroleum hydrocarbons are purchased from a chemical supplier when available. When no chemical supplier is available, the fuels are purchased from public sources.
- 7.3. The OVAP and BUSTR standard is a commercially prepared standard containing alkanes from C10-C34.

**8. SAMPLE PREPARATION, PRESERVATION AND STORAGE**

- 8.1. Refer to Section 8 of the 8000B section of this SOP.

**APPENDIX F**  
**TOTAL PETROLEUM HYDROCARBONS BASED ON 8015B**

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page F2 of F4

---

**9. QUALITY CONTROL**

9.1. Refer to Section 9 of the 8000B section of this SOP for quality control requirements, including the initial demonstration of capability, definition of a batch, surrogate limits, method blanks, laboratory control samples (LCS), and matrix spikes (MS).

9.2. MS/MSD recoveries are calculated from a Diesel calibration.

9.3. Surrogates

9.3.1. Because of the nature of the TPH analysis, whereas certain petroleum mixtures can override the C9 surrogate, the C9 surrogate recoveries are advisory. Re-extraction due to surrogate recoveries is determined by analyst judgement.

NOTE: Ohio VAP rules require reanalysis when surrogate recoveries are outside of control limits.

**10. CALIBRATION AND STANDARDIZATION**

10.1. Recommended Instrument Conditions

10.1.1. Hydrogen carrier gas - flow rate 5 - 6 mL/min

10.1.2. Detector gas mixture - air hydrogen mixture in a 10:1 ratio, air 80 - 120 mL/min, hydrogen 8 - 12 mL/min

10.1.3. Temperature Program - refer to Appendix

10.1.4. Injection volume - 1 µL

10.2. Initial Calibration

10.2.1. Analyze a five point Diesel calibration standard referring to the recommended instrument conditions. The calibration concentrations are 100, 200, 500, 1000, and 2000 ng/uL. A 5000ng/uL standard may be analyzed if needed. The retention time window of C10-C32 shall be used for the Diesel calibration.

10.2.2. For Ohio VAP and BUSTR calibrations, analyze a five point calibration for the carbon range C10-C20. The concentrations are 60, 120, 240, 600 and 1200 ug/mL. Analyze a five point of the carbon range C20-C34. The concentration ranges are 80, 160, 320, 800, and 1600 ug/mL.

10.3. Continuing Calibration

10.3.1. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements

**11. PROCEDURE**

11.1. Refer to the method 8000B section of this SOP for procedural requirements.

11.2. Extraction

The extraction procedure is described in SOP #CORP-OP-0001NC.



**APPENDIX F**  
**TOTAL PETROLEUM HYDROCARBONS BASED ON 8015B**

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page F3 of F4

---

11.3. Cleanup

11.4. Analytical Sequence – Refer to Section 11 in the 8000B Section of this SOP.

11.5. Petroleum Hydrocarbon Identification and/or Fingerprinting

11.5.1 To identify the type of petroleum hydrocarbon, compare the chromatographic peak pattern to the patterns of known petroleum hydrocarbons analyzed under identical chromatographic conditions. Samples are quantified against diesel, but fingerprinting may be done when client requested.

11.5.2. Positive matching may not be possible, even using site-specific hydrocarbons. Degradation of the pattern can occur during environmental exposure of the fuel. See Table 2 for possible fingerprints

11.6 Sample Quantification

11.6.1. Samples are quantified against the initial calibration of diesel or DRO on a single column.

11.6.2. The total height or area of the hydrocarbon is determined in the same manner used for the hydrocarbon standard.

11.6.3. If the amount of sample injected into the GC exceeds the working range of the calibration curve, an appropriate dilution is performed before reanalysis

**12. DATA ANALYSIS AND CALCULATIONS**

12.1. Refer to the 8000B section of this SOP for identification and quantitation of single component analytes.

12.2. Surrogate recovery results are calculated and reported for Nonane (C-9) unless it is determined that sample interference has adversely affected the quantitation of one of the surrogates. The surrogate must be within QC criteria. Corrective action is only necessary if Nonane (C-9) is outside of acceptance limits

**13. METHOD PERFORMANCE**

13.1. Performance limits for the four replicate initial demonstration of capability are required as referenced under Section 13.1 of the main body of this SOP.

**14. POLLUTION PREVENTION**

Refer to section 14 of the 8000B section of this SOP.

**15. WASTE MANAGEMENT**

15.1 Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required

**16. REFERENCES**

16.1 SW846, Method 8015B, Nonhalogenated Organics Using GC/FID, Test Methods for Evaluating Solid Waste, Third Edition, USEPA

16.2. Related SOP

**APPENDIX F**  
**TOTAL PETROLEUM HYDROCARBONS BASED ON 8015B**

SOP No. CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page F4 of F4

16.2.1. CORP-OP-0001NC, Extraction and Cleanup of Organic Compounds from Waters and Soils, Based on SW846 3500 Series, 3600 Series, 8150, 8151, and 600 Series Methods

**Table F1: Suggested GC Temperature Program for TPH analysis**

Initial Temperature	40°C
Initial Hold Time	4 minutes
Temperature Program	10°C/minute
Final Temperature	280°C
Final Hold Time	10 minutes

**Table F2: Reporting Limits for TPH Analysis**

Analyte	Reporting Limits		
	Water (µg/L)	Solids (mg/kg)	Waste Dilution (mg/kg)
TPH (as Diesel) or DRO	100	3.3	200
C10-C20 (OVAP & BUSTR)	60	2.0	
C20-C34 (OVAP & BUSTR)	80	2.3	
Fingerprint Compounds <sup>1</sup>			
Mineral Spirits	Kerosene	Motor Oil	
Hydraulic Oil	Jet Fuel	Stoddard Solvent	
DRO Spiking Solution			
Decane	Dodecane	Tetradecane	
Hexadecane	Octobecane	Eicosane	
Docosane	Tetracosane	Hexacosane	
Octacosane			

<sup>1</sup> This list represents most of the common petroleum hydrocarbons. The list may be expanded to include other petroleum hydrocarbons.

## APPENDIX G

SOP No. CORP-GC-0001NC

**MODIFIED NON-HALOGENATED ORGANIC COMPOUNDS  
BY 8015B, DIRECT INJECTION**Revision No: 5.7Revision Date: 10/01/03Page G1 of G5

---

**1. Scope and Application**

- 1.1. This method is applicable to the determination of the concentration of various Non-halogenated Organic Compounds in waters, wastes, sludges, and solids. It is based on SW-846 Method 8015B
- 1.2. The applicable LIMs method codes are J7 (GC/FID 8015) and QU (Semivolatile Organics, 8015B). The preparation code is 88.

**2. Summary of Method**

- 2.1 This method provides gas chromatographic conditions for the detection of various nonhalogenated organic compounds. Samples are introduced to the GC by direct injection. Detection is achieved by a flame ionization detector (FID).

**3. Definitions**

- 3.1. Refer to the glossary in the STL North Canton Laboratory Quality Manual (LQM).

**4. INTERFERENCES**

- 4.1. Refer to the method 8000B section of this SOP for general information regarding chromatographic interferences.

**5. SAFETY**

- 5.1. Refer to section 5 of the Method 8000B SOP for general safety requirements.

**6. EQUIPMENT AND SUPPLIES**

- 6.1 Refer to Section 6 of the 8000B section of this SOP
- 6.2. Recommended Columns:
  - 6.2.1. RTX Stabilwax, fused silica, 60 m x 0.53, 0.5  $\mu$ m film thickness, or equivalent column.
  - 6.2.2. Stabilwax-DA, fused silica, 60 m x 0.32, 0.5  $\mu$ m film thickness, or equivalent column.
- 6.3. Detectors: Flame ionization (FID)
- 6.4. Microsyringes, various sizes, for standards preparation, sample injection, and extract dilution.

**7. REAGENTS AND STANDARDS**

- 7.1. Refer to Section 7 of the 8000B section of this SOP.
  - 7.1.1. Reagent water
- 7.2. Standards

## APPENDIX G

SOP No. CORP-GC-0001NC

**MODIFIED NON-HALOGENATED ORGANIC COMPOUNDS  
BY 8015B, DIRECT INJECTION**Revision No: 5.7Revision Date: 10/01/03Page G2 of G5

---

7.2.1. Refer to Table G2.

**8. SAMPLE PREPARATION, PRESERVATION AND STORAGE**

8.1. Refer to section 8 of the 8000B section of this SOP.

**9. QUALITY CONTROL**

9.1. Refer to Section 9 of the 8000B section of this SOP for quality control requirements, including the initial demonstration of capability, definition of a batch, surrogate limits, method blanks, laboratory control samples (LCS), and matrix spikes (MS).

**10. CALIBRATION AND STANDARDIZATION**

10.1. Recommended Instrument Conditions

10.1.1. The following conditions are recommended. The following table lists specific information for various compounds.

10.1.2. RTX- Stabilwax and Stabilwax-DA Columns:

Carrier Gas:	Hydrogen
Initial Temp:	45 °C
Initial Hold:	3 mins
Ramp Rate A:	5 °C/min
Ramp Rate B:	30 °C/min
Final Hold A:	0 min
Final Hold B:	3 mins
Analysis Time:	17.83 mins
Injector Temp:	275 °C
FID Temp:	300 °C
Injection Vol:	1 µL

10.2. Initial Calibration

10.2.1. For each non-halogenated organic compound and surrogate standard, analyze five or more calibration standards referring to the recommended GC conditions in Section 10.1. One of the standards analyzed should be at or near the concentration which corresponds to the calibration range.

10.2.2. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements.

**11. PROCEDURE**

11.1. Refer to the method 8000B section of this SOP for procedural requirements.

11.2. Sample Preparation Summary

11.2.1. Samples received fall into three general categories: waters, soils, or wastes.

## APPENDIX G

SOP No. CORP-GC-0001NC

**MODIFIED NON-HALOGENATED ORGANIC COMPOUNDS  
BY 8015B, DIRECT INJECTION**Revision No: 5.7Revision Date: 10/01/03Page G3 of G5

---

## 11.3. Sample Preparation Procedure

## 11.3.1. Waters

11.3.1.1. Add surrogate solution to the sample to achieve a concentration of 100 mg/L.

## 11.3.2. Sediments/soils and waste

11.3.2.1. Mix the contents of the sample container with a narrow metal or wood spatula. Weigh 5 g (wet weight) into a tarred culture tube. Use a top-loading balance. Record the weight to 0.01 gram.

11.3.2.2. Quickly add 5 mL of reagent water. Add surrogate standard. Cap the vial and vortex to mix for two minutes.

11.3.2.3. If extract is cloudy or has suspended sediment particles, refrigerate and allow sample to sit for a maximum of 24 hours. Filter sample if necessary.

## 11.4. Sample Analysis

## 11.4.1. Preliminary Evaluation

11.4.1.1. The sample or sample extract is introduced to the GC column by direct inject techniques. The concentration of the sample components is then calculated from the resulting chromatograms

11.4.2. Analytical Sequence – Refer to Section 11 in the 8000B Section of this SOP.

11.4.3. Inject 1 µL of the sample extract or diluted sample into the GC using the same operating conditions and techniques as those used in the calibration of the instrument.

## 11.5. Analytical Documentation

11.5.1. Record all analytical information in the analytical logbook/logsheet, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.5.2. All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. Logbooks are reviewed by the supervisor or designee.

11.5.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

**12. DATA ANALYSIS AND CALCULATIONS**

12.1. Refer to the 8000B section of this SOP for identification and quantitation of single component analytes.

12.2. Surrogate recovery results are calculated and reported unless it is determined that sample interference has adversely affected the quantitation of the surrogate. The surrogate must be within QC criteria. Corrective action is only necessary if the surrogate is outside of acceptance limits.

## APPENDIX G

SOP No. CORP-GC-0001NC

**MODIFIED NON-HALOGENATED ORGANIC COMPOUNDS**  
**BY 8015B, DIRECT INJECTION**

Revision No: 5.7Revision Date: 10/01/03

Page G4 of G5

**13. METHOD PERFORMANCE**

- 13.1. Performance limits for the four replicate initial demonstration of capability are required as referenced under Section 13.1 of the main body of this SOP.

**14. POLLUTION PREVENTION**

Refer to section 14 of the 8000B section of this SOP.

**15. WASTE MANAGEMENT**

- 15.1. Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**

- 16.1. SW846, Method 8015B, Nonhalogenated Organics Using GC/FID, Test Methods for Evaluating Solid Waste, Third Edition, USEPA

**17. Miscellaneous (Tables, Appendices, Etc...)****Table G1 Non-halogenated Organic Compounds, Reporting Limits<sup>1</sup>**

Compound	CAS Number	Reporting Limits	
		Water, mg/L	Solid, mg/kg
2-Methoxyethanol	109-86-4	1.0	1.0
Methanol	67-56-1	1.0	0.5
Isopropyl alcohol	67-63-0	1.0	0.5
n-Propyl alcohol	71-23-8	1.0	0.5
Ethanol	64-17-5	1.0	0.5
n-Butanol	71-36-3	1.0	0.5
1,4-Dioxane	123-91-1	1.0	0.5
Ethylene oxide	75-21-8	1.0	0.5
iso-Butanol	78-83-1	1.0	0.5

<sup>1</sup> If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.

Table G2 Non-halogenated Organic Compounds Working Standards

COMPOUND	LEVEL1	LEVEL2	LEVEL3	LEVEL4	LEVEL5	LEVEL 6
2-Methoxyethanol	1.0	2.0	5.0	10.0	20.0	50.0
Methanol	1.0	2.0	5.0	10.0	20.0	50.0
Isopropyl alcohol	1.0	2.0	5.0	10.0	20.0	50.0
n-Propyl alcohol	1.0	2.0	5.0	10.0	20.0	50.0
Ethanol	1.0	2.0	5.0	10.0	20.0	50.0
n-Butanol	1.0	2.0	5.0	10.0	20.0	50.0
1,4-Dioxane	1.0	2.0	5.0	10.0	20.0	50.0
Ethylene oxide	1.0	2.0	5.0	10.0	20.0	50.0
iso-Butanol	1.0	2.0	5.0	10.0	20.0	50.0

**1. Scope and Application**

- 1.1. This method is applicable to the determination of the concentration of Sulfolane and N-Methyl-2-pyrrolidone in water and solid samples. It is based on SW846 Method 8015B. The working linear range is 50 to 1000 µg/L. Table H1 lists the reporting limits associated with this method.
- 1.2. The applicable LIMS method code is KU.

**2. Summary of Method**

- 2.1. This method provides gas chromatographic conditions for the detection of mg/L levels of Phillips 66 compounds in water. Prior to use of this method, appropriate sample extraction techniques must be used.

**3. Definitions**

- 3.1. Refer to the glossary in the STL North Canton Laboratory Quality Manual (LQM), current version.

**4. INTERFERENCES**

- 4.1. Refer to the method 8000B section of this SOP for general information regarding chromatographic interferences

**5. SAFETY**

- 5.1. Refer to section 5 of the Method 8000B SOP for general safety requirements.
- 5.2. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:
  - 5.2.1. Chemicals that have been classified as **carcinogens**, or **potential carcinogens**, under OSHA include: **Methylene Chloride**.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Refer to Section 6 of the 8000B section of this SOP.
- 6.2. Gas Chromatograph
  - 6.2.1. Gas Chromatograph: Modified to accept capillary columns
  - 6.2.2. Data System: Capable of peak integration
  - 6.2.3. Gas Chromatographic Column: 30 m x 0.32 mm ID RTX-5 fused silica capillary column
  - 6.2.4. Autosampler: Capable of reproducible injections
  - 6.2.5. Carrier Gas: Hydrogen
  - 6.2.6. Detector: Flame ionization (FID)



## APPENDIX H

## Phillips 66 Compounds

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03Page H2 of H4

---

- 6.3. Volumetric Flasks: 10, 50, and 100 mL
- 6.4. Microsyringe: 10  $\mu$ L
- 6.5. Pipettes: Disposable  $\mu$ L, Pasteur
- 6.6. Autosampler Vials: 1 mL with 11 mm crimp cap, Teflon®/silicone septum liner.

**7. REAGENTS AND STANDARDS**

- 7.1. Refer to Section 7 of the 8000B section of this SOP.
- 7.2. Reagents
  - 7.2.1. Methylene Chloride: Pesticide grade or equivalent
- 7.3. Standards
  - 7.3.1. Refer to Section 7 of the 8000B section of this SOP.

**8. SAMPLE PREPARATION, PRESERVATION AND STORAGE**

- 8.1. Refer to Section 8 of the 8000B section of this SOP.

**9. QUALITY CONTROL**

- 9.1. Refer to Section 9 of the 8000B section of this SOP for quality control requirements, including the initial demonstration of capability, definition of a batch, surrogate limits, method blanks, laboratory control samples (LCS), and matrix spikes (MS).
- 9.2. Surrogates are not used for this analysis.

**10. CALIBRATION AND STANDARDIZATION**

- 10.1. Initial Calibration
  - 10.1.1. Refer to Section 10 of the 8000B section of this SOP for general calibration requirements.

**11. PROCEDURE**

- 11.1. Refer to the method 8000B section of this SOP for procedural requirements.
- 11.2. Sample Analysis
  - 11.2.1. Summary
    - 11.2.1.1. The sample extract is injected onto the GC column. The compounds are then identified and quantitated.
  - 11.2.2. Recommended Instrument Conditions

## 11.2.2.1. GC Conditions

Initial Temperature: 45°C  
Initial Hold Time: 4 minutes  
Temperature Program: 15°C/minute  
Final Temperature: 300°C, hold 2 minutes  
Final Time: 23 minutes  
Carrier Gas: Hydrogen  
Injection Volume: 1 µL

## 11.2.3. Sample Analysis Procedure

## 11.2.3.1. Preliminary Evaluation

The sample extracts may be screened to determine the level of analyte present. If the level of analyte exceeds the working range of the calibration curve, an appropriate dilution is performed to bring the level within the calibration range.

11.2.4. Inject 1 µL of the sample extract or diluted sample into the GC using the same conditions as those used in calibration.

## 11.2.5. Identification

11.2.5.1. Analytes of interest are identified by comparing retention times with known standards.

11.2.5.2. A single column is used for identification.

## 11.2.6. Sample Quantification

11.2.6.1. Refer to Section 11 in the 8000B Section of this SOP

## 11.3. Analytical Documentation

11.3.1. Record all analytical information in the analytical logbook/logsheets, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.3.2. All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. Logbooks are reviewed by the supervisor or designee.

11.3.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

## 12. DATA ANALYSIS AND CALCULATIONS

12.1. Refer to the 8000B section of this SOP for identification and quantitation of single component analytes.

**13. METHOD PERFORMANCE**

- 13.1. Performance limits for the four replicate initial demonstration of capability are required as referenced under Section 13.1 of the main body of this SOP.

**14. POLLUTION PREVENTION**

- 14.1. Refer to section 14 of the 8000B section of this SOP.

**15. WASTE MANAGEMENT**

- 15.1. Waste generated in this procedure will be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**

- 16.1. SW846, Method 8015B, Nonhalogenated Organics Using GC/FID, Test Methods for Evaluating Solid Waste, Third Edition, USEPA

**17. Miscellaneous (Tables, Appendices, Etc...)**

- 17.1. Reporting limits

17.1.1. The lower reporting limits are shown in Table H1

17.1.2. If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.

**TABLE H1**  
**PHILLIPS 66 REPORTING LIMITS**

Compound	Reporting Limits, µg/L
Tetramethylene sulfone (Sulfolane)	50
N-Methyl-2-pyrrolidone	50

Controlled Copy No. \_\_\_\_\_

Implementation Date: 10/13/03

SOP No: CORP-0P-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 1 of 82

877 789

## STL NORTH CANTON STANDARD OPERATING PROCEDURE

**TITLE:**      **EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS FROM  
WATERS AND SOILS, BASED ON SW-846 3500 SERIES, 3600 SERIES,  
8151A AND 600 SERIES METHODS.**

(SUPERSEDES: Revision 4.0 (Dated 02/04/03))

Prepared by:

Thomas R. Fawcett

10/9/03

Date

Reviewed by:

Jim S. Hill

10/8/03

Date

Technology Specialist

Approved by:

Beth Lambert

10/8/02

Date

Quality Assurance Manager

Approved by:

Dorothy Dunn

10/8/03

Date

Environmental Health and Safety Coordinator

Approved by:

Paul H. [Signature]

10/9/03

Date

Laboratory Director

Approved by:

Mark Benne

10/9/03

Date

Technical Director

### Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 2 of 82

---

**TABLE OF CONTENTS**

1. Scope and Application.....	4	
2. Summary of Method.....	4	
3. Definitions.....	5	
4. Interferences.....	5	
5. Safety.....	5	
6. Equipment and Supplies.....	6	
7. Reagents and Standards.....	8	
8. Sample Preservations and Storage.....	10	
9. Quality Control.....	11	
10. Calibrations.....	13	
11. Procedure.....	13	
11.1 Procedural Variations.....	13	
11.2 Separatory Funnel Liquid/Liquid Extraction of Water Samples.....		13
11.3 Short Term Continuous Liquid/Liquid Extraction.....	16	
11.4 Continuous Liquid/Liquid Extraction from Water Samples.....	16	
11.5 Sonication Extraction from Soil Samples.....	20	
11.6 Soxhlet.....	23	
11.7 Accelerated Soxhlet.....	24	
11.8 Waste Dilution.....	26	
11.9 Concentration.....	26	
11.10 Nitrogen Evaporation to Final Concentration.....	28	
11.11 Turbovap Method.....	30	
11.12 Cleanup Technique.....	31	
11.13 Gel Permeation Chromatography (GPC).....	32	
11.14 GPC Column Preparation.....	32	
11.15 Initial Calibration of the GPC Column.....	33	
11.16 GPC Calibration Check.....	35	
11.17 GPC Recovery for Pest/PCBs.....	35	
11.18 GPC Blank for Semivolatiles.....	36	
11.19 Sample Extract Cleanup.....	36	
11.20 Calculations.....	37	
11.21 Florisil Cartridge Cleanup.....	38	
11.22 Modifications For Other Analyte Classes.....	40	
11.23 Sulfur Removal.....	40	
11.24 Sulfuric Acid Cleanup.....	42	
11.25 Acid/Base Partition Cleanup.....	42	
11.26 Fluorocarbon Cleanup.....	43	
12. Data Analysis and Calculations.....	44	
13. Method Performance.....	44	
14. Pollution Prevention.....	44	
15. Waste Management.....	44	
16. References.....	45	

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 3 of 82

---

17. Miscellaneous.....	45
------------------------	----

**List of Appendices:**

Appendix A	Extraction Procedure for Chlorinated Acid Herbicides Based on Method 8151A .....	55
------------	---	----

**List of Tables**

Table 1	Liquid / liquid extraction conditions
Table 2	Exchange solvents and final volumes
Table 3	Surrogate spiking solutions
Table 4	Matrix spike and LCS solutions
Table 5	Surrogate spike components
Table 6	Matrix spike components
Table A1	Herbicide surrogate spiking solutions
Table A2	Herbicide matrix spike and LCS solutions
Table A3	Herbicide surrogate spike components
Table A4	Herbicide matrix spike and LCS components

---

## **1. SCOPE AND APPLICATION**

This SOP describes procedures for preparation (extraction and cleanup) of semivolatile organic analytes in aqueous, TCLP leachate, and soil matrices for analysis by Gas Chromatography (GC) and Gas Chromatography / Mass Spectrometry (GC/MS). The procedures are based on SW-846 and 600 series methodology and are applicable for measurements made to comply with the Resource Conservation and Recovery Act (RCRA) and for wastewater testing.

- 1.1. Extraction procedures for the following determinative methods are covered: 8081A, 8082, 8141A, 8151B, 8270C, 8310, 8015B, 8100, 608, 610, and 625
- 1.2. The extraction procedures here may be appropriate for other determinative methods when appropriate spiking mixtures are used.
- 1.3. Other extraction procedures (PFE, SPE, accelerated Soxhlet, etc.) may be used but are not currently covered in this SOP.
- 1.4. The applicable LIMS method codes are: P2 (8141A), QJ (8081A), QS (8151A), QL (8270C), SG (8310), DM (608), VT (610), DP (625), A1 (8100), HS (8015B)

## **2. SUMMARY OF METHOD**

### **2.1. Separatory Funnel Extraction**

A measured volume of sample, typically 1 liter, is adjusted, if necessary, to a specified pH and serially extracted with methylene chloride using a separatory funnel.

### **2.2. Continuous Liquid/Liquid Extraction**

A measured volume of sample, typically 1 liter, is placed into a continuous liquid/liquid extractor, adjusted, if necessary, to a specific pH and extracted with methylene chloride for 18-24 hours.

### **2.3. Sonication Extraction**

A measured weight of sample, typically 30 g, is mixed with anhydrous sodium sulfate to form a free flowing powder. This is solvent extracted three times using an ultrasonic horn.

### **2.4. Soxhlet Extraction (Accelerated and Traditional)**

A 30 g sample is mixed with anhydrous sodium sulfate to form a free flowing powder. This is extracted with refluxing solvent.

### **2.5. Cleanup and Concentration**

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 5 of 82

---

Procedures are presented for removing interferences from sample extracts, and for drying and concentration of the extract to final volume for analysis.

**2.6. Phenoxy acid herbicide extractions**

Procedures for the extraction and cleanup of phenoxy acid herbicides are presented in Appendix A.

**3. DEFINITIONS**

Definitions of terms used in this SOP may be found in the glossary of the STL North Canton Laboratory Quality Manual (LQM), current version .

**4. INTERFERENCES**

- 4.1. Method interferences may be caused by contaminants in solvents, reagents, glassware, and other processing apparatus. All these materials must be routinely demonstrated to be free from interferences under conditions of the analysis by running laboratory method blanks as described in the Quality Control section. Specific selection of reagents may be required to avoid introduction of contaminants.
- 4.2. Visual interferences or anomalies (such as foaming, emulsions, odor, etc.) must be documented.

**5. SAFETY**

- 5.1. Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual, Lab Specific Addendum to the CSM, and this document.
- 5.2. Eye protection that protects against splash, laboratory coat and appropriate gloves must be worn while samples, standards, solvents and reagents are being handled. Disposable gloves that have become contaminated will be removed and discarded, other gloves will be cleaned immediately. Viton gloves may be worn when halogenated solvents are used for extractions or sample preparation. Nitrile gloves may be used when other solvents are handled. [Note: VITON is readily degraded by acetone]. When good manual dexterity is needed, for example, when handling small quantities/containers, disposable gloves (such as latex or N-DEX®) shall be used. While these gloves protect against splashes, they give little or no protection against contact with large quantities of solvent, and no protection against spills or immersion.
- 5.3. The following analytes have been tentatively classified as known or suspected, human or mammalian carcinogens: benzo(a)anthracene, benzidine, 3,3'dichlorobenzidine, benzo(a)pyrene, alpha-BHC, beta-BHC, gamma-BHC,



**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 6 of 82

delta-BHC, dibenz(a,h)anthracene, N-nitrosodimethylamine, 4,4'-DDT, and polychlorinated biphenyl compounds. Primary standards of these toxic compounds should be prepared in hood.

- 5.4. The following is a list of the materials used in this method, which have a serious or significant hazard rating. NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table. A complete list of materials used in the method can be found in the reagents and materials section. Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

<b>Material (1)</b>	<b>Hazards</b>	<b>Exposure Limit (2)</b>	<b>Signs and symptoms of exposure</b>
Methylene Chloride	Carcinogen Irritant	25 ppm-TWA 125 ppm-STEL	Causes irritation to respiratory tract. Has a strong narcotic effect with symptoms of mental confusion, light-headedness, fatigue, nausea, vomiting and headache. Causes irritation, redness and pain to the skin and eyes. Prolonged contact can cause burns. Liquid decreases the skin. May be absorbed through skin.
Sodium Hydroxide	Corrosive Poison	2 ppm, 5 mg/m <sup>3</sup>  2 Mg/M3-Ceiling	This material will cause burns if comes into contact with the skin or eyes. Severe irritant. Effects from inhalation of dust or mist vary from mild irritation to serious damage of the upper respiratory tract, depending on severity of exposure. Symptoms may include sneezing, sore throat or runny nose. Contact with skin can cause irritation or severe burns and scarring with greater exposures. Causes irritation of eyes, and with greater exposures it can cause burns that may result in permanent impairment of vision, even blindness.
Sulfuric Acid (1)	Corrosive Oxidizer Dehydrator Poison Carcinogen	1 mg/m <sup>3</sup>	This material will cause burns if comes into contact with the skin or eyes. Inhalation of vapors will cause irritation of the nasal and respiratory system. Symptoms may include irritation of the nose and throat, and labored breathing. Symptoms of redness, pain, and severe burn can occur. Contact can cause blurred vision, redness, pain, and severe tissue burns. Can cause blindness.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 7 of 82

Acetone	Flammable	1000 ppm-TWA	Inhalation of vapors irritates the respiratory tract. May cause coughing, dizziness, dullness, and headache.
Acetonitrile	Flammable Poison	40 ppm-TWA	Early symptoms may include nose and throat irritation, flushing of the face, and chest tightness. Prolonged exposure to high levels of vapors may cause formation of cyanide anions in the body.
Hexane	Flammable Irritant	500 ppm-TWA	Inhalation of vapors irritates the respiratory tract. Overexposure may cause lightheadedness, nausea, headache, and blurred vision. Vapors may cause irritation to the skin and eyes.
1 – Always add acid to water to prevent violent reactions.			
2 – Exposure limit refers to the OSHA regulatory exposure limit.			

- 5.5. Exposure to hazardous chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples should be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.6. The preparation of standards and reagents and glassware cleaning procedures that involve solvents such as methylene chloride will be conducted in a fume hood with the sash closed as far as the operations will permit. Use of methylene chloride for glassware cleaning should be avoided as far as possible. If more than 500 mL of Methylene chloride is spilled, evacuate the area until the area has been cleaned by EH&S.
- 5.7. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.
- 5.8. During Kuderna-Danish (KD) concentration, do not allow the extract to boil to dryness. The solvent vapors remaining in the KD apparatus may superheat and create an explosion or fire hazard. The KD apparatus and glass separatory funnels have ground glass joints which can become stuck. Technicians must use Kevlar or other cut/puncture resistant gloves when separating stuck joints.

### 5.9. 3510 Separatory Funnel

5.9.1. The use of separatory funnels to extract aqueous samples with Methylene Chloride creates excessive pressure very rapidly. Initial venting should be done immediately after the sample container has been sealed and inverted, periodic venting may be necessary during the extraction. Vent the funnel into the hood away from people and other samples. This is considered a high-risk activity, the use of a face shield over safety glasses or goggles is recommended. Keep the sash on the fume hood as low as reasonably possible.

### 5.10. 3520 Extraction Cont Liq/Liq

5.10.1. All personnel are to ensure liquid-liquid area is clear of unnecessary items. Heating mantles used with liquid- liquid extractions generate temperatures that could ignite some materials that come in contact with the heating mantles.

5.10.2. Ensure all solvents are away from liquid-liquid extractor. Increased temperatures near solvents can cause the pressure in the containers to increase.

5.10.3. Ensure all boiling flasks have cooled to room temperature before disconnecting liquid-liquid bodies from boiling flasks to prevent any burns.

## 6. EQUIPMENT AND SUPPLIES

6.1. Glassware should be cleaned with soap and water, rinsed with water and dried in an oven at 400°C for at least 2 hours. Alternatively the glassware can be solvent rinsed with acetone or methanol followed by methylene chloride after the water rinse.

6.2. Equipment and supplies for extraction procedures

EQUIPMENT AND SUPPLIES	Sep fun.	CLLE	Soni	Sox	Conc
Separatory Funnel: 2 L	√				
Separatory Funnel Rack	√				
Balance: >1400 g capacity, accurate ±1 g	√	√			
pH indicator paper, wide-range covers extraction pH	√	√			
Graduated cylinder: 1 liter. (other sizes may be used)	√	√			
Erlenmeyer Flask or Fleaker: 125 & 300 mL (other sizes optional)	√		√		
Centrifuge	√				
Auto-Shaker	√				

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 9 of 82

EQUIPMENT AND SUPPLIES	Sep fun	CLLE	Soni	Sox	Conc
Methylene Chloride Collection Tank	√	√			
Solvent Dispenser Pump or 100 mL Graduated Cylinder	√		√		
Continuous Liquid/Liquid Extractor		√			
Round or flat Bottom: 250, 500 mL or 1 L		√			
Boiling Chips: Contaminant free, approximately 10/40 mesh (Teflon® PTFE, carbide or equivalent).		√		√	√
Cooling Condensers		√		√	
Heating Mantle: Rheostat controlled		√		√	
Auto-timer for heating mantle		√		√	
Beakers: 250 & 400 mL, graduated			√		
450mL wide-mouth glass jars			√		
Balance: >100 g capacity, accurate ±0.1 g			√	√	
Soxhlet Extractor				√	
Cellulose and Glass Thimbles				√	
Accelerated Soxhlet Extractor (Soxtherm-trade name)				√	
Sonicator (at least 300 watts)			√		
Sonicator horn, 3/4 inch			√		
Kuderna-Danish (K-D) Apparatus: 500 mL					√
Concentrator Tube: 10 mL, attached to K-D with clips					√
Snyder Column: Three-ball macro					√
Water Bath: Heated, with concentric ring cover, capable of temperature control (± 5°C) up to 95°C. The bath must be used in a hood or with a solvent recovery system					√
Turbo Vap II					√
Turbo Vap Tube: 200mL capacity					√
Vials: Glass, 2 mL, 4 mL, and 10 mL capacity with Teflon®-lined screw-cap					√
Nitrogen Blowdown Apparatus					√
Nitrogen: reagent grade.					√
Culture tubes: 10 mL, 16 mmx100 mm					√
Syringe: 1 mL	√	√	√	√	
Phase Separation Paper	√	√	√	√	
Glass Wool	√	√	√	√	
Glass Funnel: 75 X 75 mm	√	√	√	√	√
Disposable Pipets	√	√	√	√	√
Aluminum foil	√	√	√	√	√
Paper Towels	√	√	√	√	√
Wrist Shaker					

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date 10/07/03

Page 10 of 82

### 6.3. Equipment and Supplies for Cleanup Procedures

EQUIPMENT AND SUPPLIES	GPC	Florisil	Sulfur	Acid	A/B	FC
Gel permeation chromatography system (J2 Scientific Accuprep System)	✓					
Bio Beads: (S-X3) -200-400 mesh, 70 gm (Bio-Rad Laboratories, Richmond, CA, Catalog 152-2750 or equivalent).	✓					
Chromatographic column: 700 mm x 25 mm ID glass column. Flow is upward	✓					
Ultraviolet detector: Fixed wavelength (254 nm) and a semi-prep flow-through cell.	✓					
Strip chart recorder, recording integrator, or laboratory data system.	✓					
Syringe: 10 mL with Luerlok fitting.	✓					
Syringe filter assembly, with disposable 5 um filter discs, Millipore No. LSWP 01300 or equivalent.	✓					
Chromatographic column: 250 mm long x 10 mm ID; with Pyrex glass wool at the bottom and a Teflon stopcock (for silica gel cleanup)	✓					
Vacuum system for eluting multiple cleanup cartridges Supelco (or equivalent). The manifold design must ensure that there is no contact between plastics containing phthalates and sample extracts		✓				
Vacuum trap made from a 500 mL sidearm flask fitted with a one-hole stopper and glass tubing		✓				
Vacuum pressure gauge.		✓				
Rack for holding 10 mL volumetric flasks in the manifold.		✓				
Mechanical shaker or mixer: Vortex Genie or equivalent.			✓	✓		
Separatory Funnels with Ground-Glass Stoppers: 250 mL					✓	
Erlenmeyer Flasks: 125 mL					✓	
Disposable Pipets		✓	✓	✓	✓	✓
Culture tubes: 10 mL, 16 mmx100 mm	✓	✓	✓	✓	✓	✓

## 7. REAGENTS AND STANDARDS

### 7.1. Reagents for Extraction Procedures

All reagents must be ACS reagent grade or better unless otherwise specified.

REAGENTS	Sep fun.	CLLE	Soni	Sox	Conc
Sodium hydroxide (NaOH), Pellets: Reagent Grade	✓	✓			
Sodium hydroxide solution, 10 N: Dissolve 40 g of NaOH in reagent water and dilute to 100 mL.	✓	✓			
Sulfuric acid (H <sub>2</sub> SO <sub>4</sub> ), Concentrated. Reagent Grade	✓	✓			
Sulfuric acid (1:1): Carefully add 500 mL of H <sub>2</sub> SO <sub>4</sub> to 500 mL of reagent water. Mix well.	✓	✓			

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 11 of 82

REAGENTS	Sep fun.	CLLE	Soni	Sox	Conc
Hydrochloric Acid (HCl)					
Organic free reagent water.	√	√			
Sodium sulfate (Na <sub>2</sub> SO <sub>4</sub> ), Granular, Anhydrous: Purify by heating at 400°C a minimum of two hours.	√	√	√	√	√
Extraction/Exchange Solvents: Methylene chloride, hexane, acetonitrile, acetone, pesticide quality or equivalent	√	√	√	√	√
Acetone: Used for cleaning	√	√	√	√	√

## 7.2. Reagents for Cleanup Procedures

REAGENTS	GPC	Florisil	Sulfur	Acid	A/B	FC
Florisil: 500 mg or 1 g cartridges with stainless steel or Teflon frits (catalog 694-313, Analytichem, 24201 Frampton Ave., Harbor City, CA, or equivalent.)		√				
Mercury: triple distilled			√			
Tetrabutylammonium hydrogen sulfate			√			
Sodium sulfite			√			
Tetrabutylammonium (TBA) sulfite reagent: Prepare reagent by dissolving 3.39 g of Tetrabutylammonium hydrogen sulfate in 100 mL organic-free reagent water. Extract this solution 3 times with 20 mL portions of hexane. Discard the hexane extracts. Add 25 g sodium sulfite to the water solution.			√			
2-Propanol			√			
Nitric acid: 1N			√			
Copper powder: remove oxides (if powder is dark) by treating with 1N nitric acid, rinse with organic-free reagent water to remove all traces of acid, rinse with acetone, and dry under a stream of nitrogen.			√			
Sulfuric acid, Concentrated				√		
Sodium hydroxide, Pellets					√	
Sodium hydroxide, 10N: Dissolve 40 g of NaOH in 100 mL of reagent water					√	
Sulfuric acid (H <sub>2</sub> SO <sub>4</sub> ), Concentrated: Reagent Grade					√	
Sulfuric acid (1:1): Carefully add 500 mL of H <sub>2</sub> SO <sub>4</sub> to 500 mL of reagent water. Mix well.					√	
Fluorocarbon: PF-5080, 3M						√

## 7.3. Standards

### 7.3.1. Stock Standards

Stock standards are purchased as certified solutions or prepared from neats. Semivolatile stock standards are stored at ≤ 6°C. All stock standards must be protected from light. Stock standard solutions must be replaced after one year (from the time of preparation, if prepared in house, or from the

time the ampoule is opened if purchased.) Standards must be allowed to come to room temperature before use.

#### 7.3.2. Surrogate Spiking Standards

Prepare or purchase surrogate spiking standards at the concentrations listed in Table 5. Surrogate spiking standards are purchased or prepared as dilutions of the stock standards. Surrogate spiking solutions must be refrigerated and protected from light. The standards must be replaced at least every six months or sooner if there is reason to believe that the standard has degraded or concentrated.

#### 7.3.3. Matrix Spiking and Laboratory Control Spiking Standards.

The same spiking solution is used for the matrix spike and the Laboratory Control Sample. Prepare MS/LCS spiking standards at the concentrations listed in Table 6. Spiking standards are purchased or prepared as dilutions of the stock standards. Spiking solutions must be refrigerated and protected from light. The standards must be replaced at least every six months or sooner if there is reason to believe that the standard has degraded or concentrated.

#### 7.3.4. GPC calibration solution - prepare or purchase a solution in methylene chloride that contains the following analytes in the concentrations listed below:

Analyte	mg/mL
Corn Oil	25.0
Bis (2-ethylhexyl) phthalate	1.0
Methoxychlor	0.2
Perylene	0.02
Sulfur	0.08

NOTE: Sulfur is not very soluble in methylene chloride, however, it is soluble in warm corn oil. Therefore, one approach is to weigh out the corn oil, warm it, and transfer the weighed amount of sulfur into the warm corn oil. Mix it and then transfer into a volumetric flask with methylene chloride, along with the other calibration compounds. This standard has a lifetime of 6 months.

## 8. SAMPLE COLLECTION PRESERVATION AND STORAGE

8.1. Samples are not chemically preserved.

8.2. Samples are stored at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$  in glass containers with Teflon®-lined caps.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 13 of 82

---

**8.3. Holding Times**

8.3.1 Extraction is initiated within 7 days of the sampling date for aqueous samples, 14 days for solid and waste samples.

8.3.2 For TCLP leachates, extraction is initiated within seven days from when the TCLP Leach tumbling has been completed, excluding the filtration step. If the filtration step requires extended times, this time counts as part of the seven day holding time.

8.3.3. Analysis of the extracts is completed within forty days of extraction.

**9. QUALITY CONTROL**

**9.1. Quality Control Batch**

9.1.1. The batch is a set of up to 20 samples that are of the same matrix and are processed together using the same procedures and reagents. The batch must contain a method blank, an LCS and a matrix spike / matrix spike duplicate. (In some cases, at client request, it may be appropriate to process a matrix spike and sample duplicate in place of the MS / MSD). If clients specify specific samples for MS/MSD, the batch may contain multiple MS/MSD. See policy QA-003 for further definition of the batch.

**9.2. Definition of matrix**

9.2.1. The possible matrix types are aqueous, soil, waste and TCLP leachate.

**9.3. Insufficient Sample**

9.3.1. If insufficient sample is available to process a MS/MSD, then a second LCS must be processed. The LCS pair is then evaluated according to the MS/MSD criteria. Use of a LCS pair in place of a MS/MSD must be documented.

**9.4. Sample count**

9.4.1. Laboratory generated QC samples (method blanks, LCS, MS/MSD) are not included in the sample count. Field samples are included.

**9.5. Method Blank**



- 
- 9.5.1. A method blank consisting of all reagents added to the samples must be prepared and analyzed with each batch of samples. Surrogates are spiked into the method blank at the same level as the samples. The method blank is used to identify any background interference or contamination of the analytical system which may lead to the reporting of elevated concentration levels or false positive data.
- 9.5.2. Aqueous Method Blanks use 1000 mL of reagent water spiked with the surrogates. The method blank goes through the entire analytical procedure, including any cleanup steps.
- 9.5.3. Solid method blanks use 30 g of sodium sulfate spiked with the surrogates. The method blank goes through the entire analytical procedure, including any cleanup steps.
- 9.5.4. TCLP method blanks use 250 mL of leachate fluid (100 mL for herbicides) spiked with the surrogates. The leachate may optionally be diluted to 1000 mL with reagent water. The method blank goes through the entire analytical procedure, including any cleanup steps.
- 9.6. Laboratory Control Sample (LCS)
- 9.6.1. Laboratory Control Samples are well-characterized, laboratory generated samples used to monitor the laboratory's day to day performance of routine analytical methods. The LCS, spiked with a group of target compounds representative of the method analytes, is used to monitor the accuracy of the analytical process, independent of matrix effects. On-going monitoring of the LCS results provides evidence that the laboratory is performing the method within accepted QC guidelines for accuracy and precision. The LCS goes through the entire analytical procedure, including any cleanup steps.
- 9.6.2. The LCS is made up in the same way as the method blank (See sections 9.5.1 - 9.5.4) but spiked with the LCS standard and the surrogates.
- 9.7. Surrogates
- 9.7.1. Surrogates are organic compounds which are similar to the target analyte(s) in chemical composition and behavior in the analytical process, but which are not normally found in environmental samples.
- 9.7.2. Each applicable sample, blank, LCS and MS/MSD is spiked with surrogate standards. Surrogate spike recoveries must be evaluated by determining

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 15 of 82

---

whether the concentration (measured as percent recovery) falls within the required recovery limits.

**9.7. Matrix Spike/Matrix Spike Duplicate (MS/MSD)**

9.7.1. A matrix spike is an environmental sample to which known concentrations of target analytes have been added. A matrix spike duplicate is a second spiked aliquot of the same sample which is prepared and analyzed along with the sample and matrix spike.

**9.8. Initial Demonstration of Capability**

9.8.1. The initial demonstration and method detection limit studies described in section 13 must be acceptable before analysis of samples may begin.

**9.9. Quality Assurance Summaries**

9.9.1. Certain clients may require specific project or program QC which may supersede these method requirements. Quality Assurance Summaries (QAS) should be developed to address these requirements.

**9.10. STL North Canton QC Program**

9.10.1. Further details of QC and corrective action guidelines are presented in the STL North Canton QC Program document (QA-003). Refer to this document if in doubt regarding corrective actions.

**10. CALIBRATION AND STANDARDIZATION**

10.1. On a weekly basis, measure 1.0 mL of solvent into an autovial using a gastight syringe that is manufactured to a certified volume delivery tolerance of  $\pm 0.01$  mL (1 mL total volume). The "standard" autovial is sealed and the top and bottom of the meniscus are marked. The autovials containing the sample extracts are then compared against the "standard" vial to ensure that the final volume is consistently  $1.0 \pm 0.01$  mL. If a new box of autovials are used, then the steps are repeated to further ensure that variations due to vial size and shape are minimized. A log is kept of the lot number of the vials and the day the vials were prepared.

10.2. Refer to section 11.15. for calibration of the GPC. Otherwise this section is not applicable.

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 16 of 82

## 11. PROCEDURE

Procedures for separatory funnel liquid/liquid extraction (11.2), continuous liquid/liquid extraction (11.4), sonication extraction (11.5), soxhlet extraction (11.6), accelerated soxhlet (11.7), waste dilution (11.8), extract concentration (11.9), and extract cleanup (11.19) are presented in this section.

### 11.1. Procedural Variations

Procedural variations are allowed only if deemed necessary in the professional judgment of the supervisor to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance memo and approved by a supervisor and QA/QC manager. If contractually required, the client will be notified. The Nonconformance memo will be filed in the project file.

Any deviations from this procedure identified after the work has been completed must be documented as a nonconformance, with a cause and corrective action described. A Nonconformance memo shall be used for this documentation.

### 11.2. Separatory Funnel Liquid/Liquid Extraction of Water Samples.

A flow chart for this procedure is included in Section 17.

11.2.1. Measure the initial sample pH with wide-range pH paper and record on the extraction benchsheet. If sample is a leachate (c.g. TCLP), compare the current pH against leachate log, note on the benchsheet if there is any discrepancy.

11.2.2. The normal sample volume is 1 liter. Other sample volumes may be used to obtain specific reporting limits, and reduced sample volumes, diluted to 1 liter with reagent water, may be used for very dirty samples.

11.2.3. Weigh the sample container on a balance ( $\pm 1$  g), taring the sample and container. Transfer the sample to the separatory funnel. Add the appropriate volume of surrogate spiking solution (see Table 3). Also add appropriate volume of matrix spiking solution to any matrix spike / matrix spike duplicate samples (see Table 4) Mix well. Rinse the sample bottle with 60 mL methylene chloride and transfer to the separatory funnel. Reweigh the container. Assume a density of 1 g/mL and record the difference as the sample volume on the benchsheet to the nearest milliliter.

Note: Aqueous samples must be determined volumetrically for Ohio VAP samples.

**Warning:** Dichloromethane creates excessive pressure very rapidly! Therefore, initial venting should be done immediately after the sample container has been sealed and inverted. Vent into hood away from analysts and other samples.

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 17 of 82

11.2.3.1. If the entire sample bottle will not be used, transfer the aliquot to the separatory funnel, then add the spiking solutions to the sample in the separatory funnel.

**Note:** Alternative methods of measurement of sample volume include a) transferring the sample to a measuring cylinder and b) marking a meniscus on the sample bottle and then measuring the volume of water required to fill the bottle to the meniscus after the sample is transferred. The former method is not recommended because of the risk of cross contamination while the latter is not recommended because of poor accuracy. However, either method may be necessary for specific client programs.

11.2.4. Prepare a method blank, LCS and MS/MSD for each batch as specified in section 9 of this SOP. Use 1 L of reagent water for method blanks and LCS. The LCS is spiked with the surrogate and matrix spike solutions, the method blank only with the surrogates.

11.2.5. Use 250 mL of leachate for TCLP pesticides and TCLP semivolatiles, measured in a graduated cylinder. An alternative method would be to tare 400ml beaker on the balance and then add 250g to the beaker. Assume a density of 1g/mL and record the volume on the benchsheet. The leachate may be made up to 1 L in volume with reagent water.

11.2.6. For a TCLP method blank, LCS and LCS Dup measure 250 mL of the buffer solution used in the leaching procedure and transfer to the separatory funnel. Add 60 mL of methylene chloride to the separatory funnel. The TCLP leachate may be diluted to approximately 1 liter before extraction if desired.

11.2.7. Adjust sample pH as indicated in Table 1 for the initial extraction. Use the minimum amount of 1:1 H<sub>2</sub>SO<sub>4</sub> or 10 N NaOH necessary. Recheck the sample with pH paper by dipping a disposable pipette into the sample and wetting the pH paper. Record adjusted pH, spiking volumes and standard numbers on the benchsheet. Return spiking solutions to the refrigerator as soon as possible.

11.2.8. Seal and shake or rotate the separatory funnel vigorously for 2 minutes with periodic venting to release excess pressure. An autoshaker may be used to shake and rotate the separatory funnel.

**Warning:** Dichloromethane creates excessive pressure very rapidly! Therefore, initial venting should be done immediately after the separatory funnel has been sealed and inverted. Vent into hood away from analysts and other samples.

11.2.9. Allow the organic layer to separate from the water phase until complete visible separation has been achieved (approximately 10 minutes). If the

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 18 of 82

emulsion interface between layers is more than one-third the size of the solvent layer, the analyst must employ mechanical techniques to complete the phase separation. The optimum technique depends upon the sample and may include stirring, filtration of the emulsion through glass wool, centrifugation, or other physical methods. If the emulsion cannot be broken (recovery of <80% of the methylene chloride\*), transfer the sample, solvent, and emulsion into the extraction chamber of a continuous extractor and proceed as described in continuous liquid-liquid extraction (Section 11.4.). If this is done, the sample must be extracted as part of a valid CLLE batch.

\*Note: 15 - 20 mL of methylene chloride is expected to dissolve in 1 L of water. Thus, solvent recovery could be as low as 35 mL from the first shake and still be acceptable. Subsequent shakes should recover at least 50 mL of solvent.

- 11.2.10. Fill a funnel with 10-20 g of anhydrous sodium sulfate. The funnel can be plugged with glass wool or filter paper may be used to hold the sodium sulfate. Drain the solvent extract from the separatory funnel through the prepared filtration funnel into a clean glass container. The extract may be drained directly into the KD flask. Close the stopcock just before the water level begins draining out of the separatory funnel. If the sodium sulfate becomes saturated with water add more to the funnel or replace the existing sodium sulfate with fresh drying agent.
- 11.2.11. Repeat the extraction process two more times using fresh 60 mL portions of solvent, combining the three solvent extracts in the collection container.
- 11.2.12. If extraction at a secondary pH is required, adjust the pH of the sample in the separatory funnel to the pH indicated in Table 1 with a minimum amount of 10 N NaOH or 1:1 H<sub>2</sub>SO<sub>4</sub>. Measure with pH paper and record the adjusted pH on the benchsheet. Serially extract with three 60 mL portions of methylene chloride, as outlined in Steps 11.2.7 to 11.2.9. Collect these three extracts in the same container used for the previous fraction.
- 11.2.13. Rinse the extract residue from the sodium sulfate by pouring 20-30 mL of clean methylene chloride through the funnel and into the collection container.
- 11.2.14. Dispose of solvent and water remaining in the extractor into the appropriate waste container.
- 11.2.15. Cover with aluminum foil if the extract is not concentrated immediately. Refer to Section 11.9 for concentration and Section 11.12 for cleanup.

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 19 of 82

---

11.3. Short-Term Continuous Liquid-Liquid Extraction

- 11.3.1. Extraction personnel must first determine if one liter of sample exists in the sample bottle.
- 11.3.2. Mark the meniscus on the side of the sample bottle for later determination of sample volume.
- 11.3.3. Add 225 mL methylene chloride to the extractor, turn stopcock to the open position. Connect extractor to a 250 ml boiling flask with boiling chips.
- 11.3.4. Check the pH of the sample and adjust, if necessary, using 1:1 H<sub>2</sub>SO<sub>4</sub> or 10N sodium hydroxide. Quantitatively transfer the sample from the sample bottle to the extractor. \*Rinse the sample bottle with 60-ml of methylene chloride. Pour the rinsate from the sample bottle into the continuous extractor.
- 11.3.5. Adjust the volume of solvent in the boiling flask by adding reagent water to the extractor, if necessary.
- 11.3.6. Next, fill the sample bottle to the mark with reagent water. Pour the reagent water into a graduated cylinder. Measure the reagent water to  $\pm 10$  ml and document this sample volume on the extraction benchsheet.
- 11.3.7. Add the appropriate volume of surrogate to each sample and spiking solutions to the MS, MSD, LCS, and/or LCSD in the extractor and mix well.
- 11.3.8. Extract samples for 12 hours. 8270 extractions will still be performed in two pH extraction steps. Each step will run 6 hours, however, where timing is a factor the samples may be left to extract for longer than 12 hours.
- 11.3.9. Turn stopcocks to the closed position.
- 11.3.10. Allow the extractors to run another 5 minutes, then turn off the rheostats and mantles.
- 11.3.11. The samples will macro concentrate in the CLLE. After the flasks have cooled, continue to concentrate the extracts by KD or QES methodology.
- 11.3.12 This method of extraction **cannot** be used for method 625, 610 or 608.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 20 of 82

---

11.4 Continuous Liquid/Liquid Extraction from Water Samples.

A flow chart for this procedure is included in Section 17.

- 11.4.1. Assemble the apparatus. Add 250-400mL of methylene chloride to the extractor body. Add 3 to 5 boiling chips to the round-bottom distilling flask.
- 11.4.2. Measure the initial sample pH with wide-range pH paper and record on the extraction benchsheet. If sample is a leachate (e.g. TCLP), compare the current pH against leachate log, note on the benchsheet if there is any discrepancy.
- 11.4.3. Weigh the sample container on a balance ( $\pm 1$  g), taring the sample and container. Transfer the sample to the continuous liquid/liquid extractor. Add the appropriate volume of surrogate spiking solution (see Table 3). Also add appropriate volume of matrix spiking solution to any matrix spike / matrix spike duplicate samples (see Table 4) Mix well. Rinse the sample bottle with 60 mL methylene chloride and transfer to the continuous liquid/liquid extractor. Reweigh the container. Assume a density of 1 g/mL and record the difference as the sample volume on the benchsheet to the nearest milliliter.

Note: Aqueous samples must be determined volumetrically for Ohio VAP samples.

**Warning:** Dichloromethane creates excessive pressure very rapidly! Therefore, during the rinse for 600 series methodology venting should be done immediately after the sample container has been sealed and inverted. Vent into hood away from analysts and other samples.

- 11.4.3.1. If the entire sample bottle will not be used, transfer the aliquot to the extractor, then add the spiking solutions to the sample in the extractor.

**Note:** Alternative methods of measurement of sample volume include: a.) transferring the sample to a measuring cylinder and b.) marking a meniscus on the sample bottle and then measuring the volume of water required to fill the bottle to the meniscus after the sample is transferred. The former method is not recommended because of the risk of cross contamination while the latter is not recommended because of poor accuracy. However, either method may be necessary for specific client programs.

- 11.4.3.2. Prepare a method blank, LCS and MS/MSD for each batch as specified in Section 9 of this SOP. Use 1 L of reagent water for method blanks

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 21 of 82

---

and LCS. The method blank is spiked with the surrogates, the LCS and matrix spikes with the surrogates and matrix spiking solutions.

11.4.3.3. Use 250mL of leachate for TCLP for TCLP semivolatiles and TCLP pesticides measured in a graduated cylinder. An alternative method would be to tare 400ml beaker on the balance and then add 250g to the beaker. Assume a density of 1g/mL and record the volume on the benchsheet. Dilute to about 1 liter with reagent water.

11.4.3.4. For a TCLP method blank, LCS and LCS Dup measure 250 mL of the buffer solution used in the leaching procedure and transfer to the separatory funnel. Dilute to about 1 liter with reagent water.

Less than one liter of sample may be used, for highly contaminated samples, or if the reporting limit can be achieved with less than one liter of sample. In this event dilute the sample to about 1 liter with reagent water.



877 810


**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 22 of 82




EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 23 of 82



11.4.4. Adjust sample pH as indicated in Table 1 for the initial extraction. Use the minimum amount of 1:1  $\text{H}_2\text{SO}_4$  or 10 N NaOH necessary. Recheck the sample with pH paper. Record adjusted pH, spiking volumes and standard numbers on the benchsheet. Return spiking solutions to the refrigerator as soon as possible.

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 24 of 82

11.4.5. Add reagent water to the extractor body until approximately 250 mL of methylene chloride is pushed over into the round-bottomed flask to ensure proper operation and solvent cycling. Attach cold condenser (about 10°C). Turn on heating mantle. Inspect joints for leaks once solvent has begun cycling. Extract for 18-24 hours. (24 hours required for 600 series)

11.4.6. If extraction at a secondary pH is required, (see Table 1) turn off the heating mantle and allow the extractor to cool. Detach the condenser and adjust the pH of the sample in the extractor body to the pH indicated in Table 1 with a minimum amount of 10 N NaOH or 1:1 H<sub>2</sub>SO<sub>4</sub>. Measure with pH paper and record the adjusted pH on the benchsheet. Reattach the condenser and turn on heating mantle. Extract for 18-24 hours.

11.4.7. Turn off the heating mantle and allow the extractor to cool.

11.4.8. Place a funnel containing 10-20 g of anhydrous sodium sulfate on the Kuderna-Danish (K-D) apparatus or other glass container. The funnel can be plugged with glass wool enabling it to hold the granular anhydrous sodium sulfate or phase separation filter paper may be used.

11.4.9. Dry the extract in the round bottom flask by filtering it through the sodium sulfate filled funnel. Note that it is not necessary or advisable to attempt to add the solvent remaining in the continuous extractor body to the extract.

11.4.10 Collect the dried extract in a K-D or other glass container. Rinse the funnel with 20-30 mL of methylene chloride to complete the quantitative transfer. Dispose of solvent and water remaining in the extractor in the appropriate waste container.

**Note:** Some types of CLLE apparatus have built in drying columns. If this type of apparatus is used then a drying step subsequent to the extraction may not be necessary.

11.4.11. Cover with aluminum foil if the extract is not concentrated immediately. Refer to Section 11.9 for concentration and Section 11.12 for cleanup.

#### 11.5. Sonication

A flow chart for this procedure is included in Section 17.

11.5.1. Determination of percent moisture (Optional - if a different group performs this test, refer to the facility SOP.)

In some cases, sample results are needed on a dry weight basis. If this is the case, weigh 5-10 g of sample into a suitable tared container (typically an aluminum weigh pan. Determine the % moisture by drying overnight (at least 12 hours) at 105°C. Allow to cool in a desiccator before weighing.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 25 of 82

$$\% \text{ Moisture} = \frac{(\text{weight of wet sample} - \text{weight of dry sample})}{\text{weight of wet sample}} \times 100$$

- 11.5.2. Determination of pH (Optional - if a different group performs this test, refer to the facility SOP.)

If pH determination is required, transfer 10 g of soil to a beaker. Add 10 mL of water. Stir for five minutes then let stand for 1 hour. Determine the pH of the sample with a glass electrode and pH meter.

- 11.5.3. Decant and discard any water layer on a sediment/soil sample. Record and document if a water layer was discarded on the benchsheet. Homogenize the sample by mixing it thoroughly in the container. If this is not possible place the sample in clean beaker and homogenize. Upon completion of homogenization in beaker return sample to original container. Discard foreign objects such as sticks, leaves and rocks, unless extraction of this material is required by client. If the sample consists primarily of foreign materials consult with the client (via the Project Manager or Administrator).
- 11.5.4. Remove surrogate and matrix spiking solutions from refrigerator and allow to warm to room temperature.
- 11.5.5. Weigh 30 g of sample  $\pm$  0.2g into a 250 or 400 mL beaker. Record the weight to the nearest 0.01 g in the appropriate column on the benchsheet. Use 30 g of sodium sulfate for the method blank and 30 g of sodium sulfate with 30 g of reagent sand for the LCS.
- 11.5.6. Mix weighed sample with a spatula adding enough anhydrous sodium sulfate (approximately 30 g ) to be free flowing. (If the sample is not free flowing extraction efficiency may be reduced)
- 11.5.7. Prepare a method blank, LCS and MS/MSD for each batch as specified in Section 9 of this SOP.
- 11.5.8. Add the appropriate volume of surrogate spiking solution (see Table 3). Also add appropriate volume of matrix spiking solution to any matrix spike / matrix spike duplicate samples (See Table 4) Add 1 mL of the surrogate spiking solution to each sample, method blank, Laboratory Control Sample (LCS), and matrix spikes. Refer to Table 6 for details of the surrogate spiking solutions. Add 1 mL of the appropriate matrix spiking solution to each Matrix Spike/Matrix Spike Duplicate (MS/MSD) and LCS.

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 26 of 82

Refer to Tables 3 and 5 for details of the spiking solutions. Record spiking volumes and standard numbers on the benchsheet. Return spiking solutions promptly to refrigerator.

**Note:** The same volume of surrogate and matrix spiking solution is used if GPC is indicated since the final volume would be reduced to compensate for loss of extract during the GPC procedure.

- 11.5.9. Immediately add a minimum of 100 mL of solvent to the beaker.

Solvents:

Semivolatile GC/MS, TPH, Organochlorine pesticides and PCBs	1:1 Methylene Chloride / Acetone
--	----------------------------------

**Note:** Steps 11.5.5 - 11.5.9 should be performed rapidly to avoid loss of the more volatile extractables.

- 11.5.10. Place the bottom surface of the appropriate disrupter horn tip approximately ½ inch below the surface of the solvent, but above the sediment layer.

- 11.5.11. Sonicate for 3 minutes, making sure the entire sample is agitated. If the W-380 or W-385 sonicator is used the output should be set at 6 for the 3/4 inch high gain (Q) horn or 10 for the 3/4 inch standard horn. with mode switch on pulse, and percent-duty cycle knob set at 50%.

**Note:** Do *not* use *Microtip* probe.

- 11.5.12. Loosely plug the stem of a 75 mm x 75 mm glass funnel with glass wool and/or line the funnel with filter paper. Add 10-20 g of anhydrous sodium sulfate to the funnel cup.

- 11.5.13. Place the prepared funnel on a collection apparatus (beaker or K-D Apparatus).

- 11.5.14. Decant and filter extracts through the prepared funnel into a clean beaker or K-D Apparatus.

- 11.5.15. Repeat the extraction two more times with additional 100 mL minimum portions of solvent each time. Decant off extraction solvent after each sonication. On the final sonication pour the entire sample (sediment and solvent) into the funnel and rinse with an additional 10 mL-20 mL of the methylene chloride/acetone appropriate solvent (Refer to Table in 11.5.9).

**Note:** Alternatively, the three extracts may be collected together and then filtered through the sodium sulfate.

- 11.5.16. Cover with aluminum foil if the extract is not concentrated immediately. Refer to Section 11.9 for concentration and Section 11.12 for cleanup.

- 11.5.17. Sonicator Tuning: Tune the sonicator according to manufacturer's instructions. The sonicator must be tuned at least every time a new horn is installed.

#### 11.6. Soxhlet

##### 11.6.1. Determination of % moisture

In some cases, sample results are needed on a dry weight basis. If this is the case, weigh 5-10 g of sample into a suitable tared container (typically an aluminum weigh pan. Determine the % moisture by drying overnight at 105°C. Allow to cool in a desiccator before weighing.

$$\% \text{ Moisture} = \frac{(\text{weight of wet sample} - \text{weight of dry sample})}{\text{weight of wet sample}} \times 100$$

##### 11.6.2. Determination of pH

If pH determination is required, transfer 10 g of soil to a beaker. Add 10 mL of water. Stir for 1 hour. Determine the pH of the sample with a glass electrode and pH meter.

- 11.6.3. Decant and discard any water layer on a sediment/soil sample. Record and document if a water layer was discarded on the benchsheet. Homogenize the sample by mixing it thoroughly in the container. If this is not possible place the sample in clean beaker and homogenize. Upon completion of homogenization in beaker return sample to original container. Discard foreign objects such as sticks, leaves and rocks, unless extraction of this material is required by client. If the sample consists primarily of foreign materials consult with the client (via the Project Manager or Administrator).

- 11.6.4. Remove surrogate and matrix spiking solutions from refrigerator and allow to warm to room temperature.

- 11.6.5. Weigh 30 g of sample  $\pm$  0.2g into a beaker, recording the weight to the nearest 0.01 g on the benchsheet. Use 30 g of sodium sulfate for the method blank and LCS. Add 30 g of anhydrous sodium sulfate and mix well. The mixture should have a free flowing texture. If not, add more sodium sulfate. Add the sample/sodium sulfate mixture to a soxhlet thimble, but do not pack the thimble tightly. The extraction thimble must drain freely for the duration of the extraction period. A glass wool plug above and below the sample in the soxhlet extractor is an acceptable alternative for the thimble.

- 11.6.6. Add the appropriate amount of surrogate and matrix spiking solution as indicated in Tables 3, 4, 5, and 6.

Sample weights less than 30 g but over 5 g may be used if the appropriate reporting limits can be met.

- 11.6.7. Prepare a method blank, LCS and MS/MSD for each batch as specified in Section 9 of this SOP, using sodium sulfate as the matrix. The weight of sodium sulfate used should be approximately the weight of soil used in each sample.
- 11.6.8. Place approximately 200mL of solvent into a 250 mL flat bottom flask containing one or two clean boiling chips. Attach the flask to the extractor and extract the sample for 16-24 hours at 4-6 cycles per hour. Check the system for leaks at the ground glass joints after it has warmed up.

**NOTE:** If a reduced quantity of sample is extracted, it is usually necessary to increase the amount of sodium sulfate added or increase the solvent boiling rate to properly set the cycling rate.

Solvents:

Semivolatile GC/MS, OPP, PAH, TPH Organochlorine pesticides and PCBs	1:1 Methylene Chloride / Acetone
---	----------------------------------

- 11.6.9. Allow the extract to cool after the extraction is complete, then disassemble by gently twisting the soxhlet from the flask. Dry the extract in the flask by filtering it through a sodium sulfate filled funnel.
- 11.6.10. Collect the dried extract in a K-D or other glass container. Rinse the flask which contained the solvent extract with 20-30 mL of methylene chloride and add it to the funnel to complete the quantitative transfer.
- 11.6.11. Cover with aluminum foil if the extract is not concentrated immediately. Refer to Section 11.9 for concentration and Section 11.12 for cleanup.

11.7. Accelerated Soxhlet (Soxtherm Trade Name)

11.7.1. Demonstration of % moisture

In some cases, sample results are needed on a dry weight basis. If this is the case, weigh 5-10 g of sample into a suitable tared container (typically an aluminum weigh pan. Determine the % moisture by drying overnight at 105°C. Allow to cool in a desiccator before weighing.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No. 4.1

Revision Date: 10/07/03

Page 29 of 82

$$\% \text{Moisture} = \frac{(\text{Weight of wet sample} - \text{Weight of dry sample}) \times 100}{\text{Weight of wet sample}}$$

11.7.2. Determination of pH

If pH determination is required, transfer 10 g of soil to a beaker. Add 10 mL of water. Stir for 1 hour. Determine the pH of the sample with a glass electrode and Ph meter.

11.7.3. Decant and discard any water layer on a sediment/soil sample. Record and document if a water layer was discarded on the benchsheet. Homogenize the sample by mixing it thoroughly in the container. If this is not possible, place the sample in clean beaker and homogenize. Upon completion of homogenization in beaker, return sample to original container. Discard foreign objects such as sticks, leaves and rocks, unless extraction of this material is required by client. If the sample consists primarily of foreign materials, consult with the client (via the Project Manager or Administrator).

11.7.4. Remove surrogate and matrix spiking solutions from refrigerator and allow to return to room temperature.

11.7.5. Weigh 30 g of sample  $\pm$  0.2 g into a beaker, recording the weight to the nearest 0.01 g on the benchsheet. Use 30 g of sodium sulfate for the method blank and LCS. Add 30 g of anhydrous sodium sulfate and mix well. The mixture should have a free flowing texture. If not, add more sodium sulfate. Add the sample/sodium sulfate mixture to a soxhlet thimble, but do not pack the thimble tightly. The extraction thimble must drain freely for the duration of the extraction period. A glass wool plug above and below the sample in the thimble is required.

11.7.6. Add the appropriate amount of surrogate and matrix spiking solution as indicated in Tables 3, 4, 5, and 6.

11.7.7. Sample weights less than 30 g, but over 5 g may be used if the appropriate reporting limits can be met.

11.7.8. Prepare a method blank, LCS and MS/MSD for each batch as specified in Section 9 of this SOP, using sodium sulfate as the matrix. The weight of sodium sulfate used should be approximately the weight of soil used in each sample.

11.7.9. Place thimble in beaker containing clean boiling chips and add approximately 140 mL of solvent. Place beakers into positions on the accelerated soxhlet unit. Run appropriate program for the extraction solvent. Check the system for leaks at the joints periodically

Solvents:



**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 30 of 82

Semivolatile MS, PAH, TPH	1:1 Methylene Chloride / Acetone
Semivolatile GC PCB, PEST, OPP	1:1 Hexane/Acetone

11.7.10. Upon completion of the program, remove the beaker from the unit and dispose of the extracted sample. Dry the extract in the flask by filtering it through a sodium sulfate filled funnel.

11.7.11. Collect the dried extract in a K-D, TurboVap tube, or other glass container. Rinse the flask which contained the solvent extract with 10-15 mL of methylene chloride and pour it through the funnel. Rinse the funnel with 10-15 mL of methylene chloride to complete the quantitative transfer.

11.7.12. Cover with aluminum foil if the extract is not concentrated immediately. Refer to Section 11.9 for concentration and Section 11.12 for cleanup.

#### 11.8. Waste Dilution

11.8.1. This method is used for materials that are soluble in an organic solvent.

11.8.2. Transfer 10 mL of the solvent to be used for dilution into a Teflon capped vial. Mark the meniscus on the vial, then discard the solvent.

11.8.3. Tare the vial, then transfer approximately 1g of sample to the vial. Record the weight to the nearest 0.01 g.

11.8.4. Add appropriate volume of surrogate and spike solutions (Table 3).

11.8.5. Dilute to 10 mL with the appropriate solvent. (Methylene Chloride for GC/MS- and GCS-TPH. Add 10 mL of appropriate solvent (Hexane) for GCS-pesticide and/or PCB analysis.

11.8.6. Add 2 g + 0.1 g sodium sulfate to the sample. Cap and shake or vortex each extract for 2 minutes.

11.8.7. If H<sub>2</sub>O is still present, add 4-5 g sodium sulfate to a small pipette funnel. The funnel can be plugged with glass wool or phase separation filter paper may be used to hold the sodium sulfate.

11.8.8. Pour the sample through the funnel, collecting as much as possible in a clean vial. Do NOT rinse the funnel with additional solvent, and do NOT concentrate the sample. The final volume is defined as 10 mL.

11.8.9. Label the sample, which is now ready for cleanup or analysis.


#### 11.9. Concentration

According to the type of sample and any cleanup procedures needed, different final solvents and volumes will be required. Refer to Table 2 for the appropriate final volumes and concentrations.

11.9.1. Kuderna-Danish (KD) Method:

11.9.1.1. Assemble a Kuderna-Danish concentrator by attaching a 10 mL concentrator tube to the 500 mL KD flask. Transfer the sample to the K-D flask.

11.9.1.2. Add one or two clean boiling chips and the extract to be concentrated to the KD flask and attach a three ball Snyder Column. Add approximately 1 mL of clean methylene chloride to the top of the Snyder column (this is important to ensure that the balls are not stuck and that the column will work properly).



EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 32 of 82

- 11.9.1.3. Place the KD apparatus on a water bath (90-98°C) so that the tip of the concentrator tube is submerged. The water level should not reach the joint between the concentrator and the KD flask. At the proper rate of distillation, the balls will actively chatter but the chambers should not flood.
- 11.9.1.4. Concentrate to 5-15 mL. If the determinative method requires a solvent exchange add the appropriate exchange solvent (50 mL hexane, 20 mL hexane/acetone, or 4 mL acetonitrile, 4 mL toluene (Herbicide TCLP), or 10 mL toluene) to the top of the Snyder Column, and then continue the water bath concentration back down to 1-4 mL. Refer to Table 2 for details of exchange solvents and final volumes. The Snyder column may be insulated if necessary to maintain the correct rate of distillation.  
*Note:* Add an additional boiling chip with the addition of exchange solvent.  
An alternative technique for solvent exchange is to replace the macro Snyder column and KD flask with a micro Snyder column, concentrate to approximately 1 mL, add the appropriate solvent, and concentrate back down to 1 mL. The extract must be cool before the macro Snyder assembly is removed.  
*Note:* It is very important not to concentrate to dryness as analytes will be lost.
- 11.9.1.5. Remove the KD apparatus from the water bath and allow to cool for a minimum of 10 minutes. If the level of the extract is above the level of the concentrator tube joint, continue to distill the solvent as necessary. Again, allow the KD flask to cool for a minimum of 10 minutes.
- 11.9.1.6. If the final volume is 5 or 10 mL the extract may be made up to volume in the graduated KD tube or transferred to a 12 mL vial previously marked at the appropriate volume level. Document the final volume. Otherwise proceed to Section 11.10

11.10. Nitrogen Evaporation to Final Concentration

- 11.10.1. Transfer the entire extract to a calibrated evaporation tube. Rinse the concentrator tube with 1-2 mL of the appropriate solvent and transfer the solvent rinsate to the evaporation tube.
- 11.10.2. Place the tube in a warm water bath that is at least 5°C below the boiling temperature of the solvent being evaporated and evaporate the solvent using a gentle stream of nitrogen. The nitrogen flow will form a slight depression on the surface of the solvent, but should not create splattering of the extract.

Boiling points of commonly used solvents are:

---

Methylene chloride	40°C
Acetone	56°C
Hexane	69°C
Acetonitrile	82°C
Toluene	111°C

- 11.10.3. During the course of the evaporation rinse the sides of the evaporation tube twice with approximately 1 mL of clean solvent. The first rinse should be about half way through the process, with the second rinse when the solvent volume gets close to 1 mL. Concentrate the solvent accurately to the calibrated volume line and transfer the extract to the appropriate storage vial.

**Note:** It is very important not to concentrate to dryness as analytes will be lost.

An alternative technique is to follow the previous steps concentrating the solvent to slightly below the required final volume and then drawing the extract into a syringe. Rinse the evaporation tube with a small amount of solvent and draw additional solvent into the syringe to make up the accurate final volume.

**Note:** It is very important not to concentrate to dryness as analytes will be lost.

**Note:** The final concentration and volume measurement steps are critical. Use care when concentrating and make certain that the final volume measurement is accurate.

#### Alternative QES Concentration Method:

*This concentration method uses a hot water jacketed concentrator tube (CT) instead of the hot water bath and concentrator tube used in Section 11.9.1.1. The construction of the jacketed concentrator tube reduces the tendency of the extract to evaporate to dryness. Thus, low boiling analytes are retained in the extract better with less analyst monitoring of the concentration process.*

*Assemble the jacketed concentrator tube, KD body and hot water hoses. Add 1mL of exchange solvent (if needed) and one large, clean boiling chip.*

**NOTE:** *The boiling chips used in the jacketed concentrator tube must be large enough to prevent them falling down into the tip of the CT. If the boiling chip is not in the proper position, the*

---

*extract may superheat and bump vigorously if the extract solvent warms up slowly.*

*Pour the extract to be concentrated into the KD flask, and attach a three-ball Snyder column.*

*Turn on the hot (95°C) water flow to the jacketed concentrator tube.*

*Concentrate to 1 mL. If the determinative method requires a solvent exchange, add the appropriate exchange solvent (5-10 mL hexane, 1 mL acetonitrile, 4 mL acetone/hexane) to the top of the Snyder column. Continue to concentrate until the Snyder column balls stop chattering.*

*Cool the jacketed concentrator tube until it is cold to the touch.*

*Quantitatively transfer the extract and dilute to final volume, or continue concentration with nitrogen evaporation (Section 11.10).*

#### 11.11. Turbovap Method

- 11.11.1. Turn on the Turbovap and adjust the water temperature to 5-10°C less than the boiling point of the solvent to be evaporated.
- 11.11.2. Switch all endpoint sensors to the correct position.
- 11.11.3. Adjust the water bath level
- 11.11.4. Adjust the nitrogen gas pressure to approximately 14 psi.
- 11.11.5. Transfer the extract into the Turbovap tube and load into the Turbovap. Do not fill the Turbovap tubes over approximately 3/4 full.
- 11.11.6. Reset the sensor and close the lid.

**Note:** If the extract splashes when the nitrogen flow starts, reduce the nitrogen flow or transfer a portion of the extract back into the original extract container.

- 11.11.7 As the extract concentrates, transfer the remainder of the extract into the appropriate Turbovap tube. After all of the extract has been

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 35 of 82

transferred, rinse the flask with a few mL of methylene chloride and add to the Turbovap tube.

- 11.11.8. Concentrate the extract to slightly less than the required final volume.
- 11.11.9. If solvent exchange is required, concentrate to 1-4 mL and add appropriate volume of the exchange solvent. Concentrate back down to the appropriate volume. Refer to Table 2 for details of exchange solvents and final volumes.
- 11.11.10. Transfer the concentrated extract to volumetric glassware for adjustment of final volume, using a small amount of solvent to rinse the tube and complete the transfer.

*Note:* Water contamination from condensation during concentration is not acceptable. If water is present, remove the Turbovap tube and filter the extract through sodium sulfate. Transfer to a clean Turbovap tube and continue the concentration.

*Note 2:* Dark, opaque or turbid samples may not concentrate. If this occurs, supervise the entire concentration procedure.

#### 11.12. Cleanup Techniques

- 11.12.1. The following techniques may be used to remove interfering peaks, and /or to remove materials that may cause column deterioration and/ or loss of detector sensitivity.

- 11.12.2. Gel Permeation Chromatography (Section 11.13) is a generally applicable technique which can be used to prepare extracts for Semivolatiles (8270), and pesticides (8081) analysis. It is capable of separating high molecular weight material from the sample analytes, and so is particularly useful if tissue or vegetable matter is part of the sample, and for many soil samples. Note: GPC used only for CLP Projects

- 11.12.3. Florisil column cleanup (Section 11.21) is particularly useful for cleanup of pesticides for analysis by method 8081 and should normally be applied to these samples unless the matrix is clean. It separates compounds with a different polarity from the target analytes. Note: GPC used only for CLP Projects. Gel Permeation Chromatography and Florisil column cleanup may both be applied to samples for analysis by method 8081/8082. In this case the GPC should be performed first.

- 
- 11.12.4. Sulfur cleanup (Section 11.23) is generally applied to samples for analysis by method 8081/8082, since the Electron Capture Detector responds strongly to sulfur. It is performed after GPC and Florisil cleanup.
- 11.12.5. Sulfuric acid cleanup (Section 11.24) is applied to samples requiring analysis for Polychlorinated Biphenyls (PCBs) only. Most organic matter is destroyed by the sulfuric acid. WARNING: Sulfuric acid cleanup must not be performed on any matrix that may have water present as a violent reaction between the acid and water may result in acid exploding out of the vessel.
- 11.12.6. Acid Base Partition Cleanup (Section 11.25) is useful for separating organic acids and phenols from basic and neutral organics.
- 11.12.7. Fluorocarbon cleanup (Section 11.26) is used to remove hydrocarbons from water samples to be analyzed for water soluble alcohols.
- 11.13. Gel Permeation Chromatography (GPC)
- Note: GPC system used is J2 Scientific Accuprep.
- 11.14. GPC Column Preparation
- 11.14.1. Weigh out 70 g of Bio Beads (SX-3) into a 400-mL beaker.
- 11.14.2. Add approximately 300 mL of methylene chloride and stir gently.
- 11.14.3. Cover with aluminum foil and allow the beads to swell for a minimum of two hours. Maintain enough solvent to sufficiently cover the beads at all times.
- 11.14.4. Position and tighten the outlet bed support (top) plunger assembly in the tube by inserting the plunger and turning it clockwise until snug. Install the plunger near the column end but no closer than 5 cm (measured from the gel packing to the collar).
- 11.14.5. Turn the column upside down from its normal position with the open end up. Place the tubing from the top plunger assembly into a waste beaker below the column.
- 11.14.6. Swirl the bead/solvent slurry to get a homogeneous mixture and pour the mixture into the open end of the column. Transfer as much as possible with one continuous pour trying to minimize bubble formation. Pour enough to fill the column. Wait for the excess solvent to drain out before pouring in the rest. Add additional

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 37 of 82

methylene chloride to transfer the remaining beads and to rinse the beaker and the sides of the column. If the top of the gel begins to look dry, add more methylene chloride to rewet the beads.

- 11.14.7. Wipe any remaining beads and solvent from the inner walls of the column with a laboratory tissue. Loosen the seal slightly on the other plunger assembly (long plunger) and insert it into the column. Make the seal just tight enough so that any beads on the glass surface will be pushed forward, but loose enough so that the plunger can be pushed forward.

**CAUTION:** Do not tighten the seal if beads are between the seal and the glass surface because this can damage the seal and cause leakage.

- 11.14.8. Push the plunger until it meets the gel, then compress the column bed about 4 cm.
- 11.14.9. Connect the column inlet to the solvent reservoir and place the column outlet tube in a waste container. Pump methylene chloride through the column at a rate of 5 mL/min. for one hour.
- 11.14.10. After washing the column for at least one hour, connect the column outlet tube to the inlet side of the UV detector. Connect the system outlet to the outlet side of the UV detector. Placing a restrictor (made from a piece of capillary tubing of 1/16"OD x 10/1000"ID x 2") in the outlet tube from the UV detector will prevent bubble formation which causes a noisy UV baseline. The restrictor will not effect the flow rate. After pumping methylene chloride through the column for an additional 1-2 hours, adjust the inlet bed support plunger until approximately 6-10 psi back-pressure is achieved. Push the plunger in to increase pressure or slowly pull outward to reduce pressure. c
- 11.14.11. When the GPC column is not to be used for several days, connect the column inlet and outlet lines to prevent column drying and/or channeling. If channeling occurs, the gel must be removed from the column, re-swelled, and re-poured as described above. If drying occurs, pump methylene chloride through the column until the observed column pressure is constant and the column appears wet. Always recalibrate after column drying has occurred to verify that retention volumes have not changed.

11.15. Initial Calibration of the GPC Column



EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 38 of 82

- 
- 11.15.1. Before use, the GPC must be calibrated based on monitoring the elution of standards with a UV detector connected to the GPC column.
- 11.15.2. Pump solvent through the GPC column for 2 hours. Verify that the flow rate is 4.5-5.5 mL/min. Corrective action must be taken if the flow rate is outside this range. Record the column pressure (should be 6-10 psi) and room temperature (22°C is ideal).
- Note:** Changes in pressure, solvent flow rate, and temperature conditions can affect analyte retention times and must be monitored. If the flow rate and/or column pressure do not fall within the above ranges, a new column should be prepared.
- 11.15.3. Inject the calibration solution and retain a UV trace that meets the following requirements (See resolution calculation in Section 11.20.1.):
- Peaks must be observed and should be symmetrical for all compounds in the calibration solution.
  - Corn oil and phthalate peaks must exhibit >85% resolution.
  - Phthalate and methoxychlor peaks must exhibit >85% resolution.
  - Methoxychlor and perylene peaks must exhibit >85% resolution.
  - Perylene and sulfur peaks must not be saturated and must exhibit >90% baseline resolution.
- 11.15.4. A UV trace that does not meet the criteria in Section 11.15.1. indicates the need for system maintenance and/or the need for a new column.
- 11.15.5. Determine appropriate collect and dump cycles.
- 11.15.6. The calibrated GPC program for pesticides/PCB should dump >85% of the phthalate and should collect >95% of the methoxychlor and perylene. Use a wash time of 10 minutes.
- 11.15.7. For semivolatile extracts, initiate a column eluate collection just before the elution of bis (2-ethylhexyl) phthalate and after the elution of the corn oil. Stop eluate collection shortly after the elution of perylene. Stop collection before sulfur elutes. Use a wash time of 10 minutes after the elution of sulfur.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 39 of 82

- 
- 11.15.8. Reinject the calibration solution after appropriate dump and collect cycles have been set.
- 11.15.9. Measure and record the volume of collected GPC eluate in a graduated cylinder.
- 11.15.10. The retention times for both bis(2-ethylhexyl) phthalate and perylene must not vary more than +/- 5% between calibrations.
- 11.16. GPC calibration check
- 11.16.1. Check the calibration of the GPC immediately after the initial calibration and at least every 7 days thereafter, while the column is in use.
- 11.16.2. Inject the calibration solution, and obtain a UV trace. If the retention times of bis(2-ethylhexyl)phthalate or perylene have changed by more than  $\pm 5\%$  use this run as the start of a new initial calibration. Otherwise, proceed with the recovery check. Excessive retention time shifts may be caused by poor laboratory temperature control or system leaks, an unstabilized column, or high laboratory temperature causing outgassing of methylene chloride. Pump methylene chloride through the system and check the retention times each day until stabilized.
- 11.17. GPC Recovery Check for Pesticides/ PCBs
- 11.17.1. The recovery from the GPC must be verified immediately after the initial calibration and at least every 7 days, when the instrument is in use. Two recovery check solutions are used. The first mixture is prepared by diluting 1.0 2.0 mL of the pesticide matrix spiking solution (Table 6) to 10 mL in methylene chloride. The second mixture is prepared by diluting 12.0 mL of the PCB only matrix spiking solution (Table 6) to 10 mL with methylene chloride.
- 11.17.2. Load the pesticide matrix spike mixture, the PCB mixture, and a methylene chloride blank onto the GPC using the GC dump and collect values.
- Note:** If the analysis is for PCBs only, then the pesticide recovery check is not necessary.
- 11.17.3. After collecting the GPC calibration check fraction, concentrate, solvent exchanging to hexane. Adjust the final volume to 5.0 mL, and analyze by GC/EC. Refer to concentration, Section 11.9.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 40 of 82

- 
- 11.17.4. The methylene chloride blank may not exceed more than one half the reporting limit of any analyte. And if the recovery of each of the single component analytes is 80-10120% and if the Aroclor pattern is the same as previously run standards, then the analyst may use the column. If the above criteria are not met, there may be a need for system maintenance.

11.18. GPC Blank for Semivolatiles

- 11.18.1. The recovery from the GPC must be verified immediately after the initial calibration and at least every 7 days, when the instrument is in use.
- 11.18.2. Load a methylene chloride blank onto the GPC using the semivolatiles dump and collect values.
- 11.18.3. After collecting the GPC recovery check fraction, concentrate, and analyze by GC/MS. Refer to the concentration Section 11.9.

The blank should not contain any analytes at or above the reporting limit. If these conditions are met the column may be used for sample analysis. Otherwise correct the contamination problem, or extend the collect time to improve recovery of target analytes.

11.19. Sample Extract Cleanup

- 11.19.1. Reduce the sample extract volume to 1-2 mL, then adjust to 10 mL with methylene chloride prior to cleanup. This reduces the amount of acetone in the extract. Refer to Section 11.9.
- 11.19.2. Start the pump and let the flow stabilize for 2 hours. The solvent flow rate should be 4.5-5.5 mL/min. The ideal laboratory temperature to prevent outgassing of the methylene chloride is 22°C. The normal backpressure is 6-10 psi.
- 11.19.3. In order to prevent overloading of the GPC column, highly viscous sample extracts must be diluted prior to cleanup. Any sample extract with a viscosity greater than that of a 1:1 glycerol:water solution (by visual comparison) must be diluted and loaded into several loops.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC


Revision No: 4.1

Revision Date: 10/07/03

Page 41 of 82

- 
- 11.19.4. Samples being loaded onto the GPC should be filtered with a 5 micron (or less) filter disk. Attach a filter to a 10 mL Luerlok syringe and filter the 10 mL sample extract into the sample tube.
- 11.19.5. Load the filtered samples into the proper sample tubes and place on the GPC.
- 11.19.6. Set the collect, dump, and wash times determined by the calibration procedure.
- 11.19.7. Switch to the run mode and start the automated sequence. Process each sample using the collect and dump cycle times established by the calibration procedure.
- 11.19.8. Collect each sample in a suitable glass container. Monitor sample volumes collected.
- 11.19.9. Any samples that were loaded into 2 or more positions must be recombined.
- 11.19.10. Concentrate semivolatile sample extracts to 0.5 mL in methylene chloride. Refer to the concentration Section 11.9.
- 11.19.11. Solvent exchange pesticide/PCB sample extracts into hexane and concentrate to 5.0 mL. Refer to the concentration Section 11.9.
- 11.20. Calculations
- 11.20.1. Resolution
- To calculate the resolution between two peaks on a chromatograph, divide the depth of the valley between the peaks by the peak height of the smaller peak being resolved and multiply by 100.

## Resolution Calculation



$$\% \text{ Resolution} = \frac{A}{B} \times 100$$

Where: A = depth of valley to height of smaller peak

B = peak height of smaller peak

## 11.20.2. Dump Time

Mark on the chromatograph the point where collection is to begin. Measure the distance from the injection point. Divide the distance by the chart speed. Alternatively the collect and dump times may be measured by means of an integrator or data system.

$$\text{Dump time (min)} = \frac{\text{Distance (cm) from injection to collection start}}{\text{Chart speed (cm / min)}}$$

## 11.20.3. Collection Time

$$\text{Collection time (min)} = \frac{\text{Distance (cm) between collection start and stop}}{\text{Chart speed (cm / min)}}$$

## 11.21. Florisil Cartridge Cleanup

*Florisil cleanup is generally used for organochlorine pesticides, although it may be applied to other analytes. Sections 11.21.1 through 11.21.8 outline the procedure for organochlorine pesticides, while section 11.22 outlines modifications required for other analytes.*

**Note 1:** Systems for eluting multiple cleanup cartridges include the Supelco, Inc. Solid Phase Extraction (SPE) assembly or equivalent.

---

**Note 2:** Follow the lab specific procedure when using ABC Model 1002B.

- 11.21.1. Before Florisil cleanup sample volume must be reduced to 10 mL (5 mL if GPC cleanup was used) and the solvent must be hexane. Refer to Section 11.9 for details of concentration.
- 11.21.2. Attach a vacuum manifold to a vacuum pump or water aspirator with a trap installed between the manifold and the vacuum. Adjust the vacuum in the manifold to 5-10 psi.
- 11.21.3. Place one Florisil cartridge into the vacuum manifold for each sample extract. Prior to cleanup of samples, pre-elute each cartridge with 5 mL of hexane/acetone (9:1). Adjust the vacuum applied to each cartridge so that the flow through each cartridge is approximately 2 mL/min. Do not allow the cartridges to go dry.
- 11.21.4. Just before the cartridges go dry, release the vacuum to the manifold and remove the manifold top.
- 11.21.5. Place a rack of clean labeled 12 mL concentrator tubes into the manifold and replace the manifold top. Make sure that the solvent line from each cartridge is placed inside the appropriate tube.
- 11.21.6. After the clean tubes are in place, vacuum to the manifold is restored and 2.0 mL of the extract is added to the appropriate Florisil cartridge.
- 11.21.7. The pesticides/arocloris in the extract concentrates are then eluted through the column with 8 mL of hexane/acetone (90:10) and are collected into the 10 mL culture tube or concentrator tube held in the rack inside the vacuum manifold.
- 11.21.8. Transfer the extract to a graduated concentrator tube and concentrate the extract to 2 mL. Refer to the concentration Section. (11.9)

**Note 1:** A cartridge performance standard must be run with each lot of Florisil cartridges.

**Note 2:** Florisil cartridge performance check--every lot number of Florisil must be tested before use. Add 0.5 ug/mL of 2,4,5-trichlorophenol solution and 0.5 mL of GC Standard Mix A (midpoint concentration) to 4 mL hexane. Reduce volume to 0.5 mL. Add the concentrate to a pre-washed Florisil cartridge and elute with 9 mL hexane/acetone [(90:10)(v/v)]. Rinse cartridge with 1.0 mL hexane two additional times. Concentrate eluate to 1.0 mL final volume and transfer to vial. Analyze the solution by GC/EC. The test sample

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 44 of 82

must show 80 to 120% recovery of all pesticide analytes with <5% trichlorophenol recovery, and no peaks interfering with target compounds can be detected. This standard has a lifetime of six months. Alternatively, this standard may be purchased as a stock solution.

11.22. Modifications for other analyte classes

11.22.1. PCBs

Pre-elute the cartridge with 4 mL hexane. Add 2 mL of the sample extract and elute with 3 mL hexane. The eluant will contain the PCBs together with any heptachlor, aldrin, 4,4'DDE and part of any 4,4'DDT. Any BHC isomers, heptachlor epoxide, chlordane, endosulfan I and II, endrin aldehyde and endrin sulfate and methoxychlor will be retained on the column and can be eluted in a separate fraction with 8 mL 90:10 hexane:acetone if required.

11.23. Sulfur Removal

- 11.23.1. Sulfur can be removed by one of three methods: mercury, copper, or tetrabutylammonium sulfite (TBA) according to laboratory preference. If the sulfur concentration is such that crystallization occurs in the concentrated extract, centrifuge the extract to settle the crystals, and carefully draw off the sample extract with a disposable pipet, leaving the excess sulfur in the centrifuge tube. Transfer the extract to a clean concentrator tube before proceeding with further sulfur cleanup.

11.23.2. Sulfur Removal with Elemental Mercury

**Note:** Use Mercury sparingly in order to minimize exposure and disposal costs.

- 11.23.2.1. Transfer 2 mL of sample extract into a clean concentrator tube or Teflon sealed vial.

- 11.23.2.2. Add one to three drops of mercury to the extract vial and seal.

- 11.23.2.3. Shake well for 15-30 seconds. If prolonged shaking is required, use a mechanical shaker.

- 11.23.2.4. Remove the extract from the mercury using a disposable pipette and transfer to a clean vial.

- 11.23.2.5. If black precipitate forms, sulfur was present. Shake again, then centrifuge. After centrifugation, transfer the supernate to a clean test tube and repeat. Do this until relatively little

precipitate remains, or the screens indicate that cleanup is complete.

11.23.2.6. Properly dispose of the mercury waste.

11.23.3. Sulfur Removal with Copper

11.23.3.1. Transfer 1.0 mL of sample extract into a centrifuge or concentrator tube.

11.23.3.2. Add approximately 2 g of cleaned copper powder (see Section 7.2 for copper cleaning procedure) to the sample extract tube.

11.23.3.3. Mix for one minute on a mechanical shaker.

11.23.3.4. If the copper changes color, sulfur was present. Repeat the sulfur removal procedure until the copper remains shiny.

11.23.3.5. Transfer the supernate to a clean vial.

11.23.4. Sulfur Removal with Tetrabutylammonium (TBA) Sulfite Reagent

11.23.4.1. Transfer 1.0 mL of sample extract into a culture tube.

11.23.4.2. Add 1.0 mL TBA sulfite reagent and 2 mL 2-propanol to the sample extract. Cap and shake for 1 minute. If clear crystals (precipitated sodium sulfite) form, sufficient sodium sulfite is present.

11.23.4.3. If a precipitate does not form, add sodium sulfite in an approximately 0.1 g portions until a solid residue remains after repeated shaking.

11.23.4.4. Add 5 mL organic free reagent water and shake for 1 minute. Allow sample to stand for 5-10 minutes. (Centrifuge if necessary to separate the layers). Transfer the sample extract (top layer) to a vial. The final volume is defined as 1.0 mL in Section 11.23.4.1.

11.24. Sulfuric Acid Cleanup

11.24.1. Add approximately 2-5 mL of concentrated sulfuric acid to 2 mL of sample extract in a Teflon capped vial.



EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 46 of 82

---

**Caution:** There must be no water present in the extract or the reaction may shatter the sample container.

- 11.24.2. Shake or vortex for about thirty seconds and allow to settle. (Centrifuge if necessary)
- 11.24.3. Remove the sample extract (top layer) from the acid using a Pasteur pipette and transfer to a clean vial. **CAUTION:** It is not necessary to remove all the extract since the final volume is already determined. Transfer of small amounts of sulfuric acid along with the extract will result in extremely rapid degradation of the chromatographic column.
- 11.24.4. If the sulfuric acid layer becomes highly colored after shaking with the sample extract, transfer the hexane extract to a clean vial and repeat the cleanup procedure until color is no longer being removed by the acid, or a maximum of 5 acid cleanups.
- 11.24.5. Properly dispose of the acid waste.

11.25. Acid/Base Partition Cleanup

- 11.25.1. Place 10 mL of the solvent extract from a prior extraction procedure into a 125 mL separatory funnel.
- 11.25.2. Add 20 mL of methylene chloride to the separatory funnel.
- 11.25.3. Slowly add 20 mL of DI water which has been previously adjusted to a pH of 12 to 13 with 10 N sodium hydroxide.
- 11.25.4. Seal and shake the separatory funnel for at least two minutes with periodic venting to release excess pressure.

**CAUTION:** Methylene chloride creates excessive pressure very rapidly. Initial venting should be done immediately after the separatory funnel has been sealed.
- 11.25.5. Allow the organic layer to separate from the aqueous phase for a minimum of ten minutes.
- 11.25.6. If an emulsion interface between layers is more than one-third the size of the solvent layer, the analyst must employ mechanical techniques to complete the phase separation. The optimum technique depends upon the sample, and may include stirring, filtration of the emulsion through glass wool, centrifugation, or other physical methods.

- 
- 11.25.7. Separate the aqueous phase and transfer it to a 125 mL Erlenmeyer flask. Repeat the extraction two more times using fresh 20 mL aliquots of dilute sodium hydroxide (pH 12-13). Combine the aqueous extracts.
- 11.25.8. Water-soluble organic acids and phenols will be primarily in the aqueous phase. Base/neutral analytes will be in the methylene chloride. If the analytes of interest are only in the aqueous phase discard the methylene chloride and proceed to Section 11.25.9. If the analytes of interest are only in the methylene chloride, discard the aqueous phase and proceed to Section 11.25.11.
- 11.25.9. Externally cool the flask with ice while adjusting the aqueous phase to a pH of 1-2 with sulfuric acid (1:1). Transfer the cool aqueous phase to a clean 125 mL separatory funnel.
- 11.25.10. Add 20 mL of methylene chloride to the separatory funnel and shake for at least two minutes. Allow the methylene chloride to separate from the aqueous phase and collect the methylene chloride in an Erlenmeyer flask. Repeat the extraction two more times using fresh methylene chloride and extracting at pH 1-2. Combine the three extracts.
- 11.25.11. Dry the extract by passing through a funnel containing 10-20 g anhydrous sodium sulfate. Rinse the funnel with an additional 20-30 mL of clean methylene chloride.
- 11.25.12. Cover with aluminum foil if the extract is not concentrated immediately. Refer to Section 11.9 for concentration.
- 11.25.13. Dispose of solvent and water remaining in the separatory funnel into the appropriate waste container.
- 11.26. Fluorocarbon Cleanup
- This procedure is appropriate for the removal hydrocarbons from water samples prior to analysis for water soluble alcohols*
- 11.26.1. Transfer 1-2 mL of water sample to a 10 mL culture tube.
- 11.26.2. Add 1-2 mL of fluorocarbon solvent (PF-5080) to the culture and cap.
- 11.26.3. Shake for 30-60 seconds.

- 
- 11.26.4. After the two phases have separated, pipette about 1 mL of water sample (top layer) into an autosampler vial for analysis. If necessary, centrifuge to enhance the phase separation.

## **12. DATA ANALYSIS AND CALCULATIONS**

Not applicable

## **13. METHOD PERFORMANCE**

### **13.1. Method Detection Limit**

Each laboratory must generate a valid method detection limit for each analyte of interest. The procedure for the determination of the method detection limit is given in STL North Canton QA Policy #: QA-005

### **13.2. Initial Demonstration**

Each laboratory must make an initial demonstration of capability for each individual method. This requires analysis of four QC Check samples. The QC check sample is a well-characterized laboratory generated sample used to monitor method performance, which should contain all the analytes of interest. The spiking level should be equivalent to a mid-level calibration. (For certain tests more than one set of QC check samples may be necessary in order to demonstrate capability for the full analyte list.)

- 13.2.1. Four aliquots of the QC check sample are analyzed using the same procedures used to analyze samples, including sample preparation.

- 13.2.2. Calculations and acceptance criteria for the QC check samples are given in the determinative SOPs. (CORP-GC-0001, CORP-MS-0001, 0002)

### **13.3. Training Qualification**

The group/team leader has the responsibility to ensure that this procedure is performed by an analyst who has been properly trained in its use and has the required experience.

## **14. POLLUTION PREVENTION**

Within the constraints of following the methodology in this SOP, use of organic solvents should be minimized.

## **15. WASTE MANAGEMENT**

- 
- 15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."
- 15.2. The following waste streams are produced when this method is carried out.
- 15.2.1. **Extracted aqueous samples contaminated with methylene chloride.** Samples are drained into the liquid-liquid separation unit located in extractions. This tank is then periodically rolled to the tank room where the DCM is emptied a larger tank using tubing attached to the bottom of the tank. The remaining water is neutralized with sodium bicarbonate, the pH verified and the water discharged to the sanitary sewer. The waste is disposed of as Methylene chloride contaminated waste.
- 15.2.2. **Used sodium sulfate and glass wool or filter paper contaminated with methylene chloride/acetone or acetone/hexane from the extract drying step.** These materials are disposed of in the solid waste and debris in a red container located in the extractions lab.
- 15.2.3. **Assorted flammable solvent waste from various rinses.** These wastes are put into the halogenated/non-halogenated 25 gallon solvent waste container located under the fume hood in extractions.
- 15.2.4. **Methylene chloride waste from various rinses:** These wastes are disposed of in the liquid-liquid separation unit.
- 15.2.5. **Hexane- Hexane waste:** These samples are to be disposed in the flammable waste.
- 15.2.6. **Waste Hexane in vials.** These vials are placed in the vial waste located in the GC prep laboratory.
- 15.2.7. **Waste Methylene Chloride sample vials.** These vials are placed in the vial waste located in the GC prep laboratory.
- 15.2.8. **Extracted solid samples contaminated with methylene chloride/acetone or acetone/hexane.** These materials are disposed of in the solid waste and debris in a red container located in the extractions lab.
- 15.2.9. **Samples, standards, and all extraction materials contaminated with high levels (>50ppm) of PCB's must be segregated into their own waste stream.**

---

PCB wastes are collected in one of three waste streams, solid PCB, liquid PCB and PCB vial waste.

## 16. REFERENCES

### 16.1. References

16.1.1. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW846, 3rd Edition, Final Update III (December 1996). Sections 3500B, 3510C, 3520C, 3540C, 3550B, 3600C, 3610B, 3620B, 3640A, 3650B, 3660B, AND 3665A.

16.1.2. Corporate Quality Management Plan (QMP), current version.

16.1.3. STL Laboratory Quality Manual (LQM), current version.

### 16.2. Associated SOPs and Policies, latest version

16.2.1. QA Policy, QA-003

16.2.2. Glassware Washing, NC-QA-0014

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018

16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021

16.2.5. Navy/Army SOP, NC-QA-0016

16.2.6. Hazardous Waste Management SOP, NC-HS-0001.

## 17. MISCELLANEOUS

### 17.1. Modifications from Reference method

17.1.1. Some surrogate spiking concentrations are modified from those recommended in SW-846, in order to make the concentrations more consistent with the calibration levels in the determinative methods.

17.1.2. Aqueous sample volumes may be determined by weight.

17.1.3. Spiking levels for method 608 have been reduced by a factor of ten to bring the levels within the normal calibration range of the instrument.

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 51 of 82

- 17.1.4. 10 g of soil is used for pH determination, rather than the 50 g suggested in the reference method. The volume of water is also adjusted to maintain the sample / water ratio specified in the method.
- 17.2. Modifications from previous revision
- 17.2.1. SOP change forms are on record in the Quality Assurance Department.
- 17.3. Facility Specific SOPs
- Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none.
- 17.3. Tables

Table 1 Liquid /Liquid Extraction Conditions		
Determinative Method	Initial Ext. pH	Secondary Ext. pH
BNA: 8270C <sup>1</sup>	1-2	11-12
625	11-12	1-2
Pest/PCB: 8081A, 8082 & 608	5-9	None
Wisconsin DRO	2	
OPP: 8141A	as received	None
Hydrocarbons: 8015B	as received	None
PAH: 8310, & 610	as received	None

<sup>1</sup> If the laboratory has validated acid only 8270 extraction for the target compound list required then the base extraction step may be omitted. The required validation consists of a 4 replicate initial demonstration of capability and a method detection limit study. (See section 13). Additionally, either of the base or acid fractions of Method 8270 can be run first.

**Table 2  
Exchange Solvents and Final Volumes**

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 52 of 82

Final Volumes and Exchange Solvents if no cleanup is used		
Type	Exchange Solvent for Analysis	Final Volume for Analysis in mL
Semivolatiles	N/A	2.0 mL
PCB	50mL Hexane	10.0 for solids 5.0 for H <sub>2</sub> O 2.0 for H <sub>2</sub> O**
Pesticides	50mL Hexane	10.0 for solids 5.0 for H <sub>2</sub> O
Pesticides/TCLP	50mL Hexane	3.0 mL
PAH by HPLC	4mL Acetonitrile	1.0 mL
BNA – SIM	N/A	2.0 mL - Solids & H <sub>2</sub> O
TPH	N/A	1.0
OPP	20mL Hexane/Acetone	2.0

\*Requires high sensitivity mass spec tune (refer to NC-MS-0015 PAH by SIM)

\*\* Michigan work requires a final volume of 2 mL.

Final Volumes and solvents if GPC cleanup is used CLP ONLY – SOLIDS ONLY			
Type	Solvent for GPC	Final Volume for GPC	Final Volume and solvent for Analysis
Semivolatiles	CH <sub>2</sub> Cl <sub>2</sub>	10 mL <sup>1</sup>	0.5 mL CH <sub>2</sub> Cl <sub>2</sub> - OLM03.1
Semivolatiles	CH <sub>2</sub> Cl <sub>2</sub>	10 mL <sup>1</sup>	2.5 mL CH <sub>2</sub> Cl <sub>2</sub> - OLM04.2
Pesticides	CH <sub>2</sub> Cl <sub>2</sub>	10 mL <sup>1</sup>	5 mL, Hexane

<sup>1</sup> Final volume for GPC may be 4 mL if a 2 mL sample loop is used

Final volumes and solvents if Florisil cleanup is used CLP ONLY			
Type	Solvent for Florisil	Final Volume for Florisil	Final Volume and solvent for Analysis
Pesticides	Hexane	10 mL (2 mL aliquot used)	2 mL, hexane

Final volumes and solvents if both GPC and Florisil cleanup are used CLP ONLY					
Type	Solvent for GPC	Final Volume for GPC	Solvent for Florisil	Final Volume for Florisil	Final volume for analysis
Pesticides	Methylene Chloride	10 mL	Hexane	5 mL (2 mL aliquot used)	2 mL, hexane

Note: Different final volumes may be necessary to meet special client reporting limit requirements.

Table 3
Surrogate Spiking Solutions

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 53 of 82

Analyte Group	Surrogate Spike Solution ID	Volume (mL)
BNA	100/150 ppm BNA	0.2
BNA / SIM	100/150 ppm BNA	.2 / 0.02
PEST	0.2 ppm DCB/TCX	1.0
TPH	40ng C9	1.0
PCB	0.2 ppm DCB/TCX	1.0
PAH	1.0 ug/mL p-Terphenyl-d14 5.0 ug/mL Benzo(c)pyrene	1.0
OPP	10 ug/mL Triphenyl Phosphate	1.0

Table 4		
Matrix Spike and LCS Solutions		
Analyte Group	Matrix Spike Solution ID	Volume (mL)
BNA	100 ppm BNA All-Analyte Spike & Restek Spike	0.2
BNA / SIM	100 ppm BNA All-Analyte Spike & Restek Spike	0.2 / 0.02
PEST	Pest NPDES Spike	1.0
PEST TCLP	Pest TCLP Spike	1.0
PCB	10 ppm PCB Spike	1.0
PAH	See Spike List – Table 6	1.0
TPH	See Spike List – Table 6	1.0
OPP	See Spike List – Table 6	1.0



**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 54 of 82

<b>Table 5</b>			
<b>Surrogate Spike Components</b>			
Type	Compounds	Solvent	Conc. (µg/mL)
BNA	2-Fluorobiphenyl	Methanol	100
	Nitrobenzene-d5		100
	p-Terphenyl-d14		100
	2-Fluorophenol		150
	Phenol-d6		150
	2,4,6-Tribromophenol		150
	1,2-Dichlorobenzene-d4		100
	2-Chlorophenol-d4		150
Pest/PCB	Decachlorobiphenyl	Methanol Acetone	0.2
	Tetrachloro-m-xylenc		0.2
TPH	Nonane (C9)	Methanol	40.0
PAH	p-Terphenyl-d-14	CH3CN	1.0
	Benzo(e)pyrene		5.0
OPP	Triphenylphosphate	Acetone	10.0

<b>Table 6</b>			
<b>Matrix Spike Components</b>			
Type	Compounds	Solvent	Conc. (µg/mL)
TCL BNA	Acenaphthene	Methanol	100
	4-Chloro-3-Methylphenol		150
	2-Chlorophenol		150
	1,4-Dichlorobenzene		100
	2,4-Dinitrotoluene		100
	4-Nitrophenol		150

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 55 of 82

Table 6			
Matrix Spike Components			
Type	Compounds	Solvent	Conc. (µg/mL)
	N-Nitroso-Di-n-Propylamine		100
	Pentachlorophenol		150
	Phenol		150
	Pyrene		100
	1,2,4-Trichlorobenzene		100
BNA TCLP	1,4-Dichlorobenzene	Methanol	100
	2,4-Dinitrotoluene		100
	Hexachlorobenzene		100
	Hexachlorobutadiene		100
	Hexachloroethane		100
	2-Methylphenol		100
	3-Methylphenol		100
	4-Methylphenol		100
	Nitrobenzene		100
	Pentachlorophenol		100
	Pyridine		100
	2,4,5-Trichlorophenol		100
	2,4,6-Trichlorophenol		100
BNA NPDES		Methanol	
	Acenaphthene		100
	Acenaphthylene		100
	Anthracene		100
	Benzo(a)anthracene		100
	Benzo(b)fluoranthene	Methanol	100

877 844

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 56 of 82

Table 6			
Matrix Spike Components			
Type	Compounds	Solvent	Conc. (µg/mL)
	Benzo(k)fluoranthene		100
	Benzo(a)pyrene		100
	Benzo(ghi)perylene		100
	Benzyl butyl phthalate		100
	Bis(2-chloroethyl)ether		100
	Bis(2-chloroethoxy)methane		100
	Bis(2-ethylhexyl)phthalate		100
	Bis(2-chloroisopropyl)ether		100
	4-Bromophenyl phenyl ether		100
	2-Chloronaphthalene		100
	4-Chlorophenyl phenyl ether		100
	Chrysene		100
	Dibenzo(a,h)anthracene		100
	Di-n-butylphthalate		100
	1,3-Dichlorobenzene		100
	1,2-Dichlorobenzene		100
	1,4-Dichlorobenzene		100
	3,3'-Dichlorobenzidine		100
	Diethyl phthalate		100
	Dimethyl phthalate		100
	2,4-Dinitrotoluene		100
	2,6-Dinitrotoluene		100
	Di-n-octylphthalate		100
	Fluoranthene		100
	Fluorene		100

EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 57 of 82

Table 6			
Matrix Spike Components			
Type	Compounds	Solvent	Conc. (µg/mL)
	Hexachlorobenzene		100
	Hexachlorobutadiene		100
	Hexachloroethane		100
	Indeno(1,2,3-cd)pyrene		100
	Isophorone		100
	Naphthalene		100
	Nitrobenzene		100
	N-Nitrosodi-n-propylamine		100
	Phenanthrene		100
	Pyrene		100
	1,2,4-Trichlorobenzene		100
	4-Chloro-3-methylphenol		100
	2-Chlorophenol		100
	2,4-Dichlorophenol		100
	2,4-Dimethylphenol		100
	2,4-Dinitrophenol		100
	2-Methyl-4,6-dinitrophenol		100
	2-Nitrophenol		100
	4-Nitrophenol		100
	Pentachlorophenol		100
	Phenol		100
	2,4,6-Trichlorophenol		100
	Acetophenone		100
	Atrazine		100
	Caprolactum		100

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 58 of 82

<b>Table 6</b>			
<b>Matrix Spike Components</b>			
Type	Compounds	Solvent	Conc. (µg/mL)
	Benzaldehyde		100
	1,1'-Biphenyl		100
	Safrole		100
	1,4-Dioxane		100
	Pronamide		100
	p-Chlorobenzilate		100
	Phenacstin		100
	Ethyl methanesulfonate		100
	2-Picoline		100
	Phorate		100
	Quinoline		100
Pest TCLP	Heptachlor	Methanol Acetone	0.5
	Heptachlor epoxide		0.5
	Lindane		0.5
	Endrin		0.5
	Methoxychlor		1.0
Pest NPDES /Pest	Aldrin	Methanol Acetone	1.0
	alpha-BHC		1.0
	beta-BHC		1.0
	delta-BHC		1.0
	gamma-BHC (Lindane)		1.0
	4,4'-DDD		1.0
	4,4'-DDE		1.0

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 59 of 82

<b>Table 6</b>			
<b>Matrix Spike Components</b>			
Type	Compounds	Solvent	Conc. (µg/mL)
	4,4'-DDT		1.0
	Dieldrin		1.0
	alpha-Endosulfan		1.0
	beta-Endosulfan		1.0
	Endosulfan Sulfate		1.0
	Endrin		1.0
	Heptachlor		1.0
	Heptachlor Epoxide		1.0

<b>Diesel Range Organics (8015B) Spike</b>	
Compound	Final Concentration
n-decane	50 µg/ml
n-dodecane	50 µg/ml
n-tetradecane	50 µg/ml
n-hexadecane	50 µg/ml
n-octadecane	50 µg/ml
n-eicosane	50 µg/ml
n-docosane	50 µg/ml
n-tetracosane	50 µg/ml
n-hexacosane	50 µg/ml
n-octacosane	50 µg/ml

<b>Organophosphorous Pesticides (8141A)</b>	
Compound	Final Concentration
dimethoate	20 µg/mL
disulfoton	20 µg/mL
famphur	20 µg/mL
methyl parathion	20 µg/mL

877 848

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 60 of 82

parathion (ethyl)	20 µg/mL
phorate	20 µg/mL
sulfotepp	20 µg/mL
thionazin	20 µg/mL
o,o,o-triethyl phosphorothioate	20 µg/mL
triphenylphosphate (surrogate)	20 µg/mL

Polynuclear Aromatic Hydrocarbons (8310)	
Compound	Final Concentration
Acenaphthylene	10 µg/mL
Carbazole	10 µg/mL
Naphthalene	10 µg/mL
1-Methylnaphthalene	10 µg/mL
2-Methylnaphthalene	10 µg/mL
Acenaphthene	10 µg/mL
Fluorene	2 µg/mL
Phenanthrene	2 µg/mL
Anthracene	2 µg/mL
Fluoranthene	2 µg/mL
Pyrene	2 µg/mL
Benzo(a)anthracene	2 µg/mL
Chrysene	2 µg/mL
Benzo(a)pyrene	2 µg/mL
Benzo(k)fluoranthene	2 µg/mL
Benzo(a)pyrene(k)fluoranthene	2 µg/mL
Dibenzo(a,h)anthracene	2 µg/mL
Benzo(g,h,i)perylene	2 µg/mL
Indeno(1,2,3-cd)pyrene	2 µg/mL

**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

---

**877 849**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 61 of 82

17.5. Flow diagrams

17.5.1. Separatory funnel extraction



877 850


**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 62 of 82



17.5.2. Continuous liquid/liquid extraction


**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**

SOP No: CORP-OP-0001NC


Revision No: 4.1

Revision Date: 10/07/03

Page 63 of 82



## 17.5.3. Sonication Extraction



**EXTRACTION AND CLEANUP OF ORGANIC COMPOUNDS  
FROM WATERS AND SOILS, BASED ON SW-846 3500 SERIES,  
3600 SERIES, 8151A, AND 600 SERIES METHODS**


SOP No: CORP-OP-0001NC

Revision No: 4.1


Revision Date: 10/07/03

Page 65 of 82

17.5.4. Soxhlet extraction



## 17.5.5. Concentration and cleanup



## 1. SCOPE AND APPLICATION

- 1.1. This method is applicable to the extraction of chlorinated herbicides in waters, solids, oils, and TCLP extracts. Appropriate compounds for extraction by this method are listed in CORP-GC-0001, Appendix D, Gas Chromatography of Phenoxy Acid Herbicides based on Method 8151A.

## 2. SUMMARY OF METHOD

- 2.1. This method is based on SW846 method 8151A. Aqueous samples are hydrolyzed if esters and acids are to be determined, then washed with methylene chloride by a separatory funnel extraction. After acidifying the sample the free acids are extracted into diethyl ether. Solids are extracted into methylene chloride/ acetone by sonication. If esters and acids are to be determined, the extract is hydrolyzed and extracted into diethyl ether. For both soils and aqueous samples, the free acid herbicides in the ether extract are esterified. The final volume is adjusted to prepare the extract for gas chromatography.

## 3. DEFINITIONS

- 3.1. Refer to section 3 of the main body of this SOP.

## 4. INTERFERENCES

- 4.1. Refer to section 4 of the main body of this SOP.

## 5. SAFETY

- 5.1. Refer to section 5 of the main body this SOP for basic safety information.
- 5.2. DIAZOMETHANE is an extremely toxic gas with an explosion potential. Since the explosion potential is catalyzed by imperfections in glass, generation of diazomethane must be carried out in glassware free of scratches, cracks, chips and which does not have ground glass joints. Solutions of diazomethane will be kept at temperatures below 90°C. Diazomethane must be generated and handled in a fume hood.
- 5.3. Diethyl ether must be free of peroxides as demonstrated by EM (or equivalent) Quant test strips. This test can be done every time the ether is used or once per week if the bottle is marked with the test date(s).
- 5.4. Concentrated potassium hydroxide solution is highly caustic.

Material (1)	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
-----------------	---------	-----------------------	--------------------------------

APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A

SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 68 of 82

Nitric Acid	Oxidizer Poison	2 ppm-TWA 4 ppm-STEL	Nitric acid is extremely hazardous; it is corrosive, reactive, an oxidizer, and a poison. Inhalation of vapors can cause breathing difficulties and lead to pneumonia and pulmonary edema, which may be fatal. Other symptoms may include coughing, choking, and irritation of the nose, throat, and respiratory tract. Can cause redness, pain, and severe skin burns. Concentrated solutions cause deep ulcers and stain skin a yellow or yellow-brown color. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
Ethyl Ether	Flammable Irritant Peroxide Former	400 ppm-TWA	General anesthesia by inhalation can occur. Continued exposure may lead to respiratory failure or death. Early symptoms include irritation of nose and throat, vomiting, and irregular respiration, followed by dizziness, drowsiness, and unconsciousness. May cause irritation, redness and pain to the eyes. Irritating to the skin and mucous membranes by drying effect. Can cause dermatitis on prolonged exposure. May be absorbed through skin. <b>May form explosive peroxides on long standing or after exposure to air or light. This material must be disposed of within six months.</b>
1 – Always add acid to water to prevent violent reactions.			
2 – Exposure limit refers to the OSHA regulatory exposure limit.			

## 6. EQUIPMENT AND SUPPLIES

6.1. Refer to Section 6 of the main body of this SOP for basic extraction equipment and supplies. Additional equipment and supplies needed for this procedure are listed below.

6.2. Diazomethane generation apparatus

## 7. REAGENTS AND STANDARDS

7.1. Reagents are listed in Section 7 of the main body of this SOP. Additional reagents and standards needed for this procedure are listed below.

7.2. Reagents

APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A

877 857  
SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 69 of 82

- 
- 7.2.1. Potassium hydroxide solution, 37% aqueous solution, (w/v): Dissolve 37 g of potassium hydroxide pellets in reagent water and dilute to 100 mL. **Caution:** Considerable heat will be generated. Other volumes of solution may be made up as convenient.
- 7.2.2. Toluene, reagent grade
- 7.2.3. 2-(2-Ethoxyethoxy)ethanol, trade name Carbitol, 98%+purity
- 7.2.4. Diazald, 99% purity
- 7.2.5. Sodium sulfate,  $\text{Na}_2\text{SO}_4$ , Anhydrous, granular, acidified: Heat sodium sulfate in a shallow tray at  $400^\circ\text{C}$  for a minimum of 4 hours to remove phthalates and other interfering organic substances. In a large beaker, acidify by slurring 500 g sodium sulfate with just enough diethyl ether to cover. Add 20 mL of concentrated sulfuric acid and mix thoroughly. Place the mixture on a steam bath in a hood to evaporate the ether, or allow the ether to evaporate overnight. Larger or smaller batches of acidified sodium sulfate may be prepared using the reagents in the same proportions.
- 7.2.6. Sodium Chloride,  $\text{NaCl}$
- 7.2.7.  $\text{BF}_3$ -Methanol, Boron trifluoride-MeOH, lab use only
- 7.2.8. Diethyl ether, reagent grade.
- 7.2.9. Trimethylsilyldiaomethane
- 7.2.10. Methanol, reagent grade.
- 7.2.11. Silica gel
- 7.2.12. 2% methanolic KOH, semi-conductor grade
- 7.3. Standards
- 7.3.1. Surrogate Standard
- 7.3.1.1. See Table A3.
- 7.3.2. Matrix Spike and LCS standard
- 7.3.3. See Table A4.



APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A

SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 70 of 82

## 8. SAMPLE COLLECTION PRESERVATION AND STORAGE

8.1. Sample collection and storage is described in Section 8 of the main body of this SOP.

## 9. QUALITY CONTROL

9.1. Refer to Section 9 of the main body of this SOP for Quality control procedures.

## 10. CALIBRATION AND STANDARDIZATION

10.1. Not applicable

## 11. PROCEDURE

11.1. Preparation of Aqueous Samples

11.1.1. Weigh the sample bottle and pour approximately 500 ml (100 mL for TCLP leachates) into a 1 liter wide-mouth amber jar. Reweigh the bottle and record the sample volume on the benchsheet, assuming a density of 1.0. Alternatively, measure 500 ml in a graduated cylinder. If less than 500 mL was used, reagent water may be added to make the volume up to 500 mL. Note: Aqueous samples must be determined volumetrically for Ohio VAP samples.

11.1.2. Spike each sample blank, LCS and MS with 1.0 mL of DCAA surrogate solution. Spike matrix spikes and LCS with 1 mL of herbicide matrix spiking solution. (Refer to tables A1 and A2 )

11.1.3. Add 60-80 g of NaCl to the sample and shake to dissolve the salt.

11.1.4. Hydrolysis

11.1.4.1. Use this step only if herbicide esters in addition to herbicide acids are to be determined. This is normally the case. If the herbicide esters are not to be determined, omit this step and go to 11.4.1.8.

11.1.4.2. Add 3 mL of 10N NaOH to the sample, seal and shake. Check the pH of the sample with pH paper. If the pH of the sample is not  $\geq 12$  adjust to  $\geq 12$  by adding more NaOH. Let the sample sit at room temperature for 2 hours to complete the hydrolysis.

11.1.4.3. Add 300 mL of methylene chloride to the amber jar.

11.1.4.4. Prior to placing samples in tumbler, the samples should be shaken or rotated vigorously for 2 minutes, venting as necessary. Place the samples in tumbler and allow them to tumble for one (1) hour. Allow the organic layer to separate from the aqueous layer. If an emulsion layer greater than one third of the solvent layer forms, use mechanical techniques to complete the phase separation. Suggested techniques are stirring, filtration through glass wool and centrifugation.

11.1.4.5. Pour contents of amber jar into a pre-rinsed teflon sep funnel.

11.1.4.6. Discard the **methylene chloride** phase.

11.1.4.7. Add 6 mL of cold (4°C) 1:1 sulfuric acid to the sample. Seal, and shake to mix. Check the pH of the sample with pH paper. If the pH is not  $\leq 2$ , and more acid to adjust the pH to  $\leq 2$ .

**Caution:** *Addition of acid may cause heat and / or pressure build up.*

11.1.4.8. Add 100 mL diethyl ether to the sample and extract by shaking or rotating vigorously for 10 minutes, venting as necessary. Allow the organic layer to separate from the aqueous layer. If an emulsion layer greater than one third of the solvent layer forms, use mechanical techniques to complete the phase separation. Suggested techniques are stirring, filtration through glass wool and centrifugation.

11.1.4.9. Drain the aqueous layer into a clean flask or beaker. Collect the ether phase in a clean flask or bottle containing approximately 10g of acidified anhydrous sodium sulfate.

11.1.4.10. Return the aqueous phase to the separatory funnel, add 100 mL diethyl ether and repeat the extraction procedure a second time, combining the ether extracts. Repeat the extraction a third time with 100 mL diethyl ether. Discard the aqueous phase after the third extraction.

11.1.4.11. Allow the extract to remain in contact with the sodium sulfate for at least 2 hours, shaking periodically. (May be left overnight). The drying step is critical: if the sodium sulfate solidifies in a cake, add a few additional grams of acidified sodium sulfate. The amount of sodium sulfate is sufficient if some free flowing crystals are visible when the flask or bottle is swirled or shaken.

11.1.4.12. Proceed to Section 11.5, Concentration.

11.2. Extraction of soil and sediment samples

APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A

SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 72 of 82

- 
- 11.2.1. Decant and discard any water layer on a sediment/soil sample. Record and document if a water layer was discarded on the benchsheet. Homogenize the sample by mixing it thoroughly in the container. If this is not possible place the sample in clean beaker and homogenize. Upon completion of homogenization in beaker return sample to original container. Discard foreign objects such as sticks, leaves and rocks, unless extraction of this material is required by client. If the sample consists primarily of foreign materials consult with the client (via the Project Manager or Administrator).
- 11.2.2. Weigh 50.0 g of moist solid sample into a clean glass jar. Use 50 g of sodium sulfate for the Method Blank and the LCS. Acidify the sample with 5 mL of concentrated HCl.
- 11.2.3. There should be a small amount of liquid phase. If not, add reagent water until there is. Stir well with a spatula. (Note: This is not necessary for the method blank or LCS)
- 11.2.4. After 15 minutes, stir with a spatula and check the pH of the liquid phase. Add more acid if necessary to bring the pH to <2, repeating the stirring and standing time after each acid addition. (Note: The pH of the method blank and LCS are not determined.)
- 11.2.5. Add 60 g of acidified sodium sulfate and mix well. The sample should be free flowing. If not, add more sodium sulfate.
- 11.2.6. Spike each sample blank, LCS and MS with 1.0 mL of DCAA surrogate solution. Spike matrix spikes and LCS with 1 mL of herbicide matrix spiking solution. (Refer to tables A1 and A2)
- 11.2.7. Add a minimum of 100 mL of 1:1 methylene chloride:acetone to the beaker or 100 mL of methylene chloride for long list (dinoseb).
- 11.2.8. Place the bottom surface of the appropriate disrupter horn tip approximately ½ inch below the surface of the solvent, but above the sediment layer.
- 11.2.9. Sonicate for 3 minutes, making sure the entire sample is agitated.
- 11.2.10. Loosely plug the stem of a 75 mm x 75 mm glass funnel with glass wool and/or line the funnel with filter paper. Add 10-20 g of anhydrous sodium sulfate to the funnel cup.
- 11.2.11. Place the prepared funnel on a collection apparatus. If the herbicide esters are *not* to be determined, the collection apparatus is a bottle or flask containing approximately 10g of anhydrous acidified sodium sulfate. If the

herbicide esters *are* to be determined, (normally the case) the collection apparatus is glassware suitable for the hydrolysis step, typically a KD flask or Turbovap tube.

- 11.2.12. Decant and filter extracts through the prepared funnel into the collection apparatus.
- 11.2.13. Repeat the extraction two more times with additional 100 mL minimum portions of the appropriate solvent each time. Decant off extraction solvent after each sonication. On the final sonication pour the entire sample (sediment and solvent) into the funnel and rinse with an additional 10 mL-20 mL of the methylene chloride/acetone.  
**Note:** Alternatively, the three extracts may be collected together and then filtered through the sodium sulfate.
- 11.2.14. If the herbicide esters are not to be determined, dry the extract as described in Section 11.4 or go to cleanup, Section 11.3. If the herbicide esters are to be determined (normally the case) proceed to Section 28.2.15
- 11.2.15. Add 5 mL of 37% aqueous potassium hydroxide and 30 mL of water to the extract. Check the pH with pH paper. If the pH is not  $\geq 12$ , adjust with additional KOH.
- 11.2.16. Heat on a water bath at 60-65°C for 2 hours. Allow to cool. Higher temperatures, up to 90°C, may be used if needed to remove the ether layer within 2 hours.
- 11.2.17. Transfer the solution to a separatory funnel and extract three times with 100 mL portions of methylene chloride. **Discard the methylene chloride phase.** The aqueous solution contains the herbicides.
- 11.2.18. Adjust the pH of the solution to  $\leq 2$  with 1:1 sulfuric acid.
- 11.2.19. Extract three times with 60 mL diethyl ether.
- 11.2.20. Proceed to Section 11.3, Cleanup, if required, or Section 11.4, Extract drying.

### 11.3. Cleanup

- 11.3.1. This cleanup step may be necessary if the procedure for determining the herbicide acids only is being followed. (See Section 28.2.14) It is not normally required if the acids and esters

are being determined. (The usual case.) If cleanup is not required, proceed to Section 28.4, Extract drying.

- 11.3.2. Prepare 45 mL of basic extraction fluid by mixing 30 mL of reagent water with 15 mL of 37% KOH. Use three 15 mL portions of this fluid to partition the extract from section 28.2.12 or 28.2.20, using a small separatory funnel. **Discard the organic phase.**

- 11.3.3. Adjust the pH of the solution to  $\leq 2$  with cold (4°C) sulfuric acid. (1:1). Extract once with 40 mL diethyl ether and twice with 20 mL diethyl ether.

*Caution: Addition of acid may cause heat and / or pressure build up.*

11.4. Extract drying

- 11.4.1. Pour the extracts through a funnel containing acidified sodium sulfate into a flask or bottle containing approximately 10 g acidified sodium sulfate. Rinse the funnel with a little extra diethyl ether.
- 11.4.2. Allow the extract to remain in contact with the sodium sulfate for at least 2 hours, shaking periodically. (May be left overnight). The drying step is critical: if the sodium sulfate solidifies in a cake, add a few additional grams of acidified sodium sulfate. The amount of sodium sulfate is sufficient if some free flowing crystals are visible when the flask or bottle is swirled or shaken. Proceed to Section 28.5, concentration.

11.5. Concentration

- 11.5.1. Transfer the ether extract into a Turbovap concentrator tube or a 500 mL K-D flask **equipped with a 10 mL concentrator tube**. Use a stirring rod to crush the caked sodium sulfate during transfer. Rinse the flask or bottle with 20-30 mL ether to complete transfer.
- 11.5.2. Attach a three ball Snyder column to the K-D apparatus, prewet the column with a few mL of ether from the top, and place the apparatus on a water bath at approximately 60°C. At the proper rate of distillation, the balls of the column will chatter, but the chambers will not flood. When the apparent volume reaches 2 mL, remove from the water bath and allow to completely cool.
- 11.5.3. For TCLP extracts only, add 4 mL of toluene to the K-D (buffer reacts with diazald solution used in esterification).

- 
- 11.5.4. Remove TCLP extracts at approximately 2 mL when boiling slows and only toluene is remaining. Cool and cap in CT.
  - 11.5.5. Carefully disassemble the concentrator tube and rinse the lower glass joint with a small amount of diethyl ether.
  - 11.5.6. The extract is now ready for esterification by the Diazomethane Bubbler Method (Section 11.7) or the TCLP esterification by Boron Trifluoride Method (Section 11.8).
- 11.6. Esterification (Bubbler Method)
- 11.6.1. Assemble the diazomethane apparatus (see figure below) in a hood. Add 25-30 mL of diethyl ether to tube 1. Add 1.5 mL 37% KOH, 1 mL Carbitol, and 1 mL of ether and 3-4 g of diazald to tube 2.
  - 11.6.2. Place the tip of the disposable pipette into the vial containing the first sample extract. Apply nitrogen flow (approx. 10 mL/min) to bubble diazomethane through the sample extract for about 4-5 minutes. Replace the disposable pipette and place the tip into the vial containing the second extract. Replace Tube 2 solution after second extract.
  - 11.6.3. Allow the extracts to stand for 20 minutes, then add approximately 0.2 g of silica gel to each extract. Allow to stand for an additional 20 minutes.
  - 11.6.4. Adjust the volume to 10 mL with hexane. The sample is now ready for gas chromatography.
  - 11.6.5. A routine 10X dilution occurs on final extracts for all samples. Due to a QuantIMS limitation, the dilution factor field in QuantIMS cannot be used when a dilution is routine, because the dilution factor is automatically applied to all reference values creating reporting problems. For the herbicide analysis, the extract volume will be 10mL and an aliquot at 10X dilution will be analyzed. The final extract volume recorded on the laboratory bench sheet will be recorded as 100mL to avoid using the dilution factor field in QuantIMS.
- 11.7. Esterification by Boron trifluoride (TCLP extracts only).
- 11.7.1. To the concentrator tube with the extract, add 2 mL of Boron trifluoride. Place a two-ball micro-Snyder column on the concentrator tube and place in the Hot-Blok water bath adjusted to 35-40°C for 60 minutes. Remove and let cool for approximately ten minutes.

APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A


SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 76 of 82

- 
- 11.7.2. With a 10 mL graduated disposable 5-3/4" pipette 4.5 mL of 5% neutral sodium sulfate and place it in the concentrator tube. Seal with a tight fitting ground glass stopper. Vortex the mixture for one minute. Let stand for ten minutes to settle. With a 5-3/4" disposable pipette, withdraw the bottom aqueous layer into a 16 x 100 culture tube for proper disposal.
- 11.7.3. Prepare a clean up column in a 5-3/4" disposable pipette by placing a small amount of glass wool in the narrow end of the pipette and add about 1/2 inch – 1 inch of florisil and sodium sulfate each. Leave about a 1 inch gap at the top. Place the extract in the clean-up column and gently force it through by using a pipette bulb into a small test tube. **Care should be taken to avoid channeling.** Rinse the concentrator tube with 1-2 mL of toluene. Transfer the column rinsate into the test tube. Rinse the column with additional toluene. There should be approximately 4 mL collected in the test tube. Bring the final volume to 10 mL with toluene by visually comparing it to a calibrated collection tube.

**Note:** It is critical that all toluene is retained and no water should enter the column.

11.8 Esterification by Trimethylsilyldiazomethane

- 11.8.2 Exchange the extract into Hexane and concentrate to 2.0 mls. 200 uL of Methanol is added to the extract then 100 uL of the trimethylsilyldiazomethane solution. The extract then turns a yellow colour. If this does not occur, then an additional 100 uL aliquot is added until the yellow persists. The extract then sits for 1 hour at room temperature to allow the methylation reaction to proceed. After 1 hour, the reaction is halted and the diazomethane removed by the addition of silicic acid. The extract is then brought up to a final volume of 10.0 mls with Hexane and submitted for analysis.



## 12. DATA ANALYSIS AND CALCULATIONS

12.1. Not applicable

## 13. METHOD PERFORMANCE

13.1. Refer to CORP-GC-0001 for details of method performance.

## 14. POLLUTION PREVENTION

14.1. Refer to Section 14 of the main body of this SOP.

## 15. WASTE MANAGEMENT

15.1. Refer to Section 15 of the main body of this SOP.

15.2. The following waste streams are produced when this method is carried out.

15.2.1. **Aqueous acidic waste.** These wastes are disposed of in the liquid-liquid separation unit

15.2.2. **Non-hazardous sodium sulfate.** Non hazardous substances can be disposed of in the regular trash.



**APPENDIX A**  
**EXTRACTION PROCEDURE FOR CHLORINATED**  
**ACID HERBICIDES BASED ON METHOD 8151A**

**SOP No: CORP-OP-0001NC**

Revision No: 4.1

Revision Date: 10/07/03

Page 78 of 82

**16. REFERENCES**

- 16.1. SW846, Test Methods for Evaluating Solid Waste, Third Edition, Update III, December 1996, Chlorinated Herbicides, Method 8151A.

**17. MISCELLANEOUS**

- 17.1. Modifications from Reference Method

17.1.1. Directions to add sufficient reagent water to the soil sample so that the pH can be measured have been added (Section 11.5.2)

17.1.2. The bubbler esterification method uses methanolic KOH in place of the aqueous KOH / carbitol mixture recommended in method 8150B. This has been found to provide a more effective and reliable esterification.

- 17.2. Modifications from previous revisions

17.2.1. References have been updated

- 17.3. Tables

<b>Table A1</b>		
<b>Herbicide Surrogate Spiking Solutions</b>		
Analyte Group	Surrogate Spike Solution ID	Volume (mL)
Herbicides	Herbicides Water	1.0
Herbicides	Herbicides Soil	1.0

<b>Table A2</b>		
<b>Herbicide Matrix Spike and LCS Solutions</b>		
Analyte Group	Matrix Spike Solution ID	Volume (mL)
Herbicides	Herbicides MS-Soil	1.0
Herbicides	Herbicides MS-Water	1.0

APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A

877 867  
SOP No: CORP-OP-0001NC  
Revision No: 4.1  
Revision Date: 10/07/03  
Page 79 of 82

Table A3 Herbicide Surrogate Spike Components			
Type	Compounds <sup>1</sup>	Solvent	Conc. (ug/mL)
Herbicides WS	2,4-DCAA	Acetone	2
Herbicides SS	2,4-DCAA	Acetone Methanol	20

<sup>1</sup>The surrogate is spiked as the free acid

Table A4 Herbicide Matrix Spike Components				
Type	Compounds <sup>1</sup>	Compounds <sup>1</sup> Solvent	Water/Soil Conc. (ug/mL)	TCLP only Conc. (ug/mL)
Herbicides MS	2,4-D	Methanol	16	2
	2,4-DB		16	2
	2,4,5-TP (Silvex)		4	0.5
	Dalapon		8	1
	Dicamba		8	1
	Dichloroprop		16	1
	Dinoseb		2.5	0.3
	2,4,5-T		4	0.5
	MCPA		1600	
	MCPP		1600	
	Pentachlorophenol		2	

<sup>1</sup>The herbicide spiking solution contains the herbicides as the free acids.

APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A


SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03

Page 80 of 82

Extraction of Aqueous Samples



APPENDIX A  
EXTRACTION PROCEDURE FOR CHLORINATED  
ACID HERBICIDES BASED ON METHOD 8151A


SOP No: CORP-OP-0001NC

Revision No: 4.1

Revision Date: 10/07/03


Page 81 of 82

Extraction of Soils and Sediments



---

**Drying, Concentration and Esterification**



PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 1 of 41

Implementation Date: 2/5/04

## STL NORTH CANTON STANDARD OPERATING PROCEDURE

**TITLE: PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS SAMPLES BY  
COLD VAPOR ATOMIC ABSORPTION, SW846 7470A AND MCAWW 245.1**

(SUPERSEDES: REVISION 2.4, DATED 10/28/02)

Prepared by:	<u><i>Patricia O'Neil</i></u>	<u>1/21/04</u>	Date
Approved by:	<u><i>Larry R. Williams</i></u>	<u>1.21.04</u>	Date
	Technical Specialist		
Approved by:	<u><i>Beth Lambert</i></u>	<u>1/21/04</u>	Date
	Quality Assurance Manager		
Approved by:	<u><i>Debra L. Dunn</i></u>	<u>1/28/04</u>	Date
	Environmental, Health and Safety Coordinator		
Approved by:	<u><i>Paul M. Carter</i></u>	<u>1/29/04</u>	Date
	Laboratory Director		

### Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions. THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY. ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 2 of 41

---

**TABLE OF CONTENTS**

1. SCOPE AND APPLICATION.....	3
2. SUMMARY OF METHOD.....	3
3. DEFINITIONS.....	4
4. INTERFERENCES.....	4
5. SAFETY.....	5
6. EQUIPMENT AND SUPPLIES.....	7
7. REAGENTS AND STANDARDS.....	9
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE.....	10
9. QUALITY CONTROL.....	10
10. CALIBRATION AND STANDARDIZATION.....	14
11. PROCEDURE.....	15
12. DATA ANALYSIS AND CALCULATIONS.....	18
13. METHOD PERFORMANCE.....	20
14. POLLUTION PREVENTION.....	20
15. WASTE MANAGEMENT.....	20
16. REFERENCES.....	21
17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . ).....	22

**LIST OF APPENDICES:**

APPENDIX A - TABLES.....	26
APPENDIX B - MSA GUIDANCE.....	30
APPENDIX C - TROUBLESHOOTING GUIDE.....	33
APPENDIX D- CONTAMINATION CONTROL GUIDELINES.....	35
APPENDIX E- PREVENTATIVE MAINTENANCE.....	37
APPENDIX F - INSTRUMENT SET-UP.....	39

## 1. SCOPE AND APPLICATION

- 1.1. This procedure describes the preparation and analysis of mercury (Hg, CAS # 7439-97-6) by Cold Vapor Atomic Absorption Spectroscopy (CVAA) using SW-846 Method 7470A and MCAWW Method 245.1.
- 1.2. The associated LIMs method codes are BL (Method 245.1) and O8 (Method 7470A).
- 1.3. CVAA analysis provides for the determination of total mercury (organic and inorganic). The combination of the oxidants, potassium permanganate and potassium persulfate, has been found to give 100% recovery with both types of compounds. Detection limits, sensitivity and optimum concentration ranges for mercury analysis will vary with the matrices, instrumentation and volume of sample used.
- 1.4. Method 7470A is applicable to the preparation and analysis of mercury in ground water, aqueous samples, TCLP, and other leachates/extracts. Certain solid and sludge type wastes may also be analyzed, however Method 7471A (see CORP-MT-0007NC) is usually the method of choice. All matrices require sample preparation prior to analysis.
- 1.5. Method 245.1 is applicable to the determination of mercury in drinking, surface and saline waters and domestic and industrial wastes. All matrices require sample preparation prior to analysis.
- 1.6. The STL North Canton reporting limit for mercury in aqueous matrices is 0.0002 mg/L except for TCLP or SPLP leachates for which the reporting limit is 0.002 mg/L.

## 2. SUMMARY OF METHOD

- 2.1. This SOP describes a technique for the determination of mercury in solution. The procedure is a physical method based on the absorption of radiation at 253.7 nm by mercury vapor. A representative portion of the sample is digested in sulfuric and nitric acids. Organic mercury compounds are oxidized with potassium permanganate and potassium persulfate and the mercury reduced to its elemental state with stannous chloride and aerated from solution in a closed system. The mercury vapor passes through a cell positioned in the light path of an atomic absorption spectrophotometer. Absorbance is measured as a function of mercury concentration. Concentration of the analyte in the sample is determined by comparison of the sample absorbance to the calibration curve (absorbance vs. concentration).



### 3. DEFINITIONS

- 3.1. Dissolved Metals: Those elements which pass through a 0.45 um membrane. (Sample is acidified after filtration).
- 3.2. Suspended Metals: Those elements which are retained by a 0.45 um membrane.
- 3.3. Total Metals: The concentration determined on an unfiltered sample following digestion.

### 4. INTERFERENCES

Chemical and physical interferences may be encountered when analyzing samples using this method.

- 4.1. Potassium permanganate, which is used to breakdown organic mercury compounds also eliminates possible interferences from sulfide. Concentrations as high as 20 mg/L of sulfide as sodium sulfide do not interfere with the recovery of inorganic mercury from reagent water.
- 4.2. Copper has also been reported to interfere; however, copper concentrations as high as 10 mg/L had no effect on the recovery of mercury from spiked samples.
- 4.3. Chlorides can cause a positive interference. Sea waters, brines and industrial effluents high in chlorides require additional permanganate (as much as 25 mL) because, during the oxidation step, chlorides are converted to free chlorine, which also absorbs radiation at 253.7 nm. Care must be taken to ensure that free chlorine is absent before the mercury is reduced and swept into the cell. This is accomplished by adding excess hydroxylamine reagent (25 mL) and purging the sample headspace before stannous chloride is added. Both inorganic and organic mercury spikes have been quantitatively recovered from seawater using this technique.

**Note:** Sufficient addition of permanganate is apparent when the purple color persists at least 15 minutes. Some samples may require dilution prior to digestion due to extremely high concentrations of chloride.

- 4.4. Interference from certain volatile organic materials that absorb at this wavelength may also occur. While the possibility of absorption from certain organic substances present in the sample does exist, this problem is not routinely encountered. This is mentioned only to caution the analyst of the possibility.

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 5 of 41

- 4.5. Samples containing high concentrations of oxidizing organic materials, as evidenced by high COD levels, may not be completely oxidized by this procedure. When this occurs the recovery of mercury will be low.
- 4.6. The most common interference is laboratory contamination, which may arise from impure reagents, dirty glassware, improper sample transfers, dirty work areas, etc. Be aware of potential sources of contamination and take appropriate measures to minimize or avoid them.

## 5. SAFETY

- 5.1. Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual and this document.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The following is a list of the materials used in this method, which have a serious or significant hazard rating. NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table. A complete list of materials used in the method can be found in the reagents and materials section. Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

Material (1)	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
Sulfuric Acid	Corrosive  Oxidizer  Dehydrator  Poison	1 Mg/M3-TWA	Inhalation produces damaging effects on the mucous membranes and upper respiratory tract. Symptoms may include irritation of the nose and throat, and labored breathing. Symptoms of redness, pain, and severe burn can occur. Contact can cause blurred vision, redness, pain and severe tissue burns. Can cause blindness.
Nitric Acid	Corrosive  Oxidizer	2 ppm-TWA  4 ppm-STEL	Nitric acid is extremely hazardous; it is corrosive, reactive, an oxidizer, and a poison. Inhalation of vapors can cause breathing difficulties and lead to pneumonia and pulmonary edema, which may be

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 6 of 41

	Poison		fatal. Other symptoms may include coughing, choking, and irritation of the nose, throat, and respiratory tract. Can cause redness, pain, and severe skin burns. Concentrated solutions cause deep ulcers and stain skin a yellow or yellow-brown color. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
Hydrochloric Acid	Corrosive Poison	5 PPM-Ceiling	Inhalation of vapors can cause coughing, choking, inflammation of the nose, throat, and upper respiratory tract, and in severe cases, pulmonary edema, circulatory failure, and death. Can cause redness, pain, and severe skin burns. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
Potassium Permanganate	Oxidizer	5 Mg/M3 for Mn Compounds	Causes irritation to the respiratory tract. Symptoms may include coughing, shortness of breath. Dry crystals and concentrated solutions are caustic causing redness, pain, severe burns, brown stains in the contact area and possible hardening of outer skin layer. Diluted solutions are only mildly irritating to the skin. Eye contact with crystals (dusts) and concentrated solutions causes severe irritation, redness, and blurred vision and can cause severe damage, possibly permanent.
Mercury (1,000 ppm in Reagent)	Oxidizer Corrosive Poison	0.1 Mg/M3 Ceiling (Mercury Compounds)	Extremely toxic. Causes irritation to the respiratory tract. Causes irritation. Symptoms include redness and pain. May cause burns. May cause sensitization. Can be absorbed through the skin with symptoms to parallel ingestion. May affect the central nervous system. Causes irritation and burns to eyes. Symptoms include redness, pain, and blurred vision; may cause serious and permanent eye damage.
Potassium Persulfate	Oxidizer	None	Causes irritation to the respiratory tract. Symptoms may include coughing, shortness of breath. Causes irritation to skin and eyes. Symptoms include redness, itching, and pain. May cause dermatitis, burns, and moderate skin necrosis.
1 – Always add acid to water to prevent violent reactions.			
2 – Exposure limit refers to the OSHA regulatory exposure limit			

- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**. Therefore, unless they are known to be non-hazardous, all samples should be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.
- 5.6. Do not look directly into the beam of the Hg lamp. The UV light that these lamps radiate is harmful to the eyes.
- 5.7. Cylinders of compressed gas must be handled with caution, in accordance with local regulations. It is recommended that, wherever possible, cylinders are located outside the laboratory and the gas led to the instrument through approved lines.
- 5.8. The CVAA apparatus must be properly vented to remove potentially harmful fumes generated during sample analysis.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Temperature controlled water bath capable of maintaining a temperature of 90-95 °C.

6.2. Atomic Absorption Spectrophotometer equipped with:

- 6.2.1. Absorption Cell with quartz end windows perpendicular to the longitudinal axis. Dimensions of the cell must result in sufficient sensitivity to meet the SOP defined reporting limit. The quartz windows must be maintained to provide accurate measurements. Any scratches or fingerprints can alter the absorption of UV radiation.
- 6.2.2. Mercury specific hollow cathode lamp (HCL) or electrodeless discharge lamp (EDL).
- 6.2.3. Peristaltic pump, which can deliver 1 L/min, air.
- 6.2.4. Flowmeter capable of measuring an airflow of 1 L/min.
- 6.2.5. Recorder or Printer.
- 6.2.6. Aeration Tubing: A straight glass frit having a coarse porosity and Tygon tubing is used for the transfer of mercury vapor from the sample bottle to the absorption cell and return.
- 6.2.7. Drying device to prevent condensation in cell. The lamp is positioned to shine on the absorption cell maintaining the air temperature in the cell about 10 °C above room temperature. Other drying devices that achieve the same purpose are also acceptable (i.e., Gortex filter).

6.3. 8oz. HDPE Plastic bottles.

6.4. Nitrogen or argon gas supply, welding grade, or equivalent.

6.5. Calibrated automatic pipettes or Class A glass volumetric pipettes.

6.6. Class A volumetric flasks.

6.7. Thermometer (capable of accurate readings at 95 °C).

6.8. Disposable cups or tubes.

---

**7. REAGENTS AND STANDARDS**

- 7.1. Reagent water must be produced by a Millipore DI system or equivalent. Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.
- 7.2. Stock (10 ppm) mercury standards (in 10% HNO<sub>3</sub>) are purchased as custom STL North Canton solutions. All standards must be stored in FEP fluorocarbon or previously unused polyethylene or polypropylene bottles. Stock standard solutions must be replaced prior to the expiration date provided by the manufacturer. If no expiration date is provided, the stock solutions may be used for up to one year and must be replaced sooner if verification from an independent source indicates a problem.
- 7.3. Working mercury standard (0.1 ppm): Take 1 mL of the stock mercury standard (7.2) and dilute to 100 mL with reagent water. The working mercury standard must be made daily and must be prepared in a matrix of 0.15% HNO<sub>3</sub>. This acid (150 uL of concentrated HNO<sub>3</sub>) must be added to the flask/bottle before the addition of the stock standard aliquot.
- 7.4. The calibration standards listed in Table I must be prepared fresh daily from the working standard (7.3) by transferring 0, 0.2, 0.5, 1.0, 5.0 and 10.0 mL aliquots of the working mercury standard into 100 mL flasks and diluting to volume with reagent water.
- 7.5. The initial calibration verification standard must be made from a different stock solution than that of the calibration standards.
- 7.6. Refer to Table I (Appendix A) for details regarding the working standard concentrations for calibration, calibration verification and spiking solutions. All standards must be processed through the entire analytical procedure including sample preparation.
- 7.7. Nitric acid (HNO<sub>3</sub>), concentrated, trace metal grade or better.
- 7.8. Sulfuric acid (H<sub>2</sub>SO<sub>4</sub>), concentrated, trace metal grade or better.
- 7.9. Stannous chloride solution: Add 50 g of stannous chloride and 25 mL of concentrated hydrochloric acid to a 500mL volumetric flask and bring to volume with deionized water.
- 7.10. Sodium chloride-hydroxylamine hydrochloride solution: Add 240 g of sodium chloride and 240 g of hydroxylamine hydrochloride to every 2000 mL of reagent water.

7.11. Potassium permanganate, 5% solution (w/v): Dissolve 100 g of potassium permanganate for every 2000 mL of reagent water.

7.12. Potassium persulfate, 5% solution (w/v): Dissolve 100 g of potassium persulfate for every 2000 mL of reagent water.

## 8. **SAMPLE COLLECTION, PRESERVATION AND STORAGE**

8.1. Sample holding time for mercury is 28 days from time of collection to the time of analysis.

8.2. Aqueous samples are preserved with nitric acid to a pH of  $<2$  and may be stored in either plastic or glass. Refrigeration is not required. Preservation must be verified prior to analysis.

## 9. **QUALITY CONTROL**

Table II (Appendix A) provides a summary of quality control requirements including type, frequency, acceptance criteria and corrective action.

### 9.1. **Initial Demonstration of Capability**

Prior to the analysis of any analyte using 7470A or the 245.1, the following requirements must be met.

9.1.1. Method Detection Limit (MDL) - An MDL must be determined for each analyte/matrix prior to the analysis of any samples. The MDL is determined using seven replicates of reagent water, spiked with all the analytes of interest, that have been carried through the entire analytical procedure. MDLs must be redetermined on an annual basis in accordance with 40 CFR Part 136 Appendix B requirements. The spike level must be between the calculated MDL and 10X the MDL to be valid. The result of the MDL determination must be below the STL North Canton reporting limit.

9.1.2. Initial Demonstration Study - This requires the analysis of four QC check samples. The QC check sample is a well-characterized laboratory generated sample used to monitor method performance. The results of the initial demonstration study must be acceptable before analysis of samples may begin.

9.1.2.1. Four aliquots of the check sample (LCS) are prepared and analyzed using the procedures detailed in this SOP and the determinative SOPs.

9.2. Preparation Batch - A batch is a group of no greater than 20 samples excluding QC Samples (LCS, Method Blank, MS, MSD) which are processed similarly, with respect to the procedure. All sample setups must be initiated within a 24 hour period from the initial preparation or extraction and without interruption of the process. All samples within the batch must be treated with the same lots of reagents and the same processes. In some cases, at client request, it may be appropriate to process a matrix spike and sample duplicate in place of the MS/MSD. If clients specify specific samples for MS/MSD, the batch may contain multiple MS/MSD pairs.

9.3. Method Blank (MB) - One method blank must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting limit or at or above 10% of the measured concentration of that analyte in associated samples, whichever is higher (sample result must be a minimum of 20 times higher than the blank contamination level).

- Repreparation and reanalysis of all samples associated with an unacceptable method blank is required when reportable concentrations are determined in the samples (see exception noted above).
- If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. **Such action must be addressed in the project narrative.**
- If the above criteria are not met and reanalysis is not possible, then the sample data must be qualified. **This anomaly must be addressed in the project narrative.**

9.4. Laboratory Control Sample (LCS) - One aqueous LCS must be processed with each preparation batch. The LCS is used to monitor the accuracy of the analytical process. On going monitoring of the LCS results provides evidence that the laboratory is performing the



---

method within acceptable accuracy and precision guidelines. The LCS must be carried through the entire analytical procedure. If the LCS is outside established control limits the system is out of control and corrective action must occur. Until in-house control limits are established, a control limit of 80 - 120% recovery must be applied.

- In the instance where the LCS recovery is > 120% and the sample results are < RL, the data may be reported with qualifiers. Such action must be addressed in the project narrative.
- In the event that an MS/MSD analysis is not possible, a Laboratory Control Sample Duplicate (LCSD) must be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.
- Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.

9.5. Matrix Spike/Matrix Spike Duplicate (MS/MSD) - One MS/MSD pair must be processed for each preparation batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Samples identified as field blanks cannot be used for MS/MSD analysis. Spiking levels are provided in Table I (Appendix A).

- If analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. Until in-house control limits are established, a control limit of 75 - 125 % recovery and 20% RPD must be applied to the MS/MSD. If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include repreparation and reanalysis of the batch. MS/MSD results, which fall outside the control limits, must be addressed in the narrative.
- If the native analyte concentration in the MS/MSD exceeds 4 times the spike level for that analyte, the recovery data are reported as NC (i.e., not calculated). If the reporting software does not have the ability to report NC then the actual recovery must be reported and narrated as follows: "Results outside of limits do not

---

necessarily reflect poor method performance in the matrix due to high analyte concentrations in the sample relative to the spike level.”

- If an MS/MSD is not possible due to limited sample volume, then a laboratory control sample duplicate (LCSD) should be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.
- 9.6. Initial Calibration Verification (ICV/ICB) - Calibration accuracy is verified by analyzing a second source standard (ICV). The ICV result must fall within 10% of the true value for that solution. An ICB is analyzed immediately following the ICV to monitor low level accuracy and system cleanliness. The ICB result must fall within +/- the reporting limit (RL) from zero. If either the ICV or ICB fail to meet criteria, the analysis should be terminated, the problem corrected and the instrument recalibrated. (See Section 11.2.8 for required run sequence). If the cause of the ICV or ICB failure was not directly instrument related the corrective action will include reparation of the ICV, ICB, CRA, CCV, and CCB with the calibration curve.
- 9.7. Continuing Calibration Verification (CCV/CCB) - Calibration accuracy is monitored throughout the analytical run through the analysis of a known standard after every 10 samples. The CCV must be a mid-range standard at a concentration other than that of the ICV. The CCV result must fall within 20% of the true value for that solution. A CCB is analyzed immediately following each CCV. (See Section 11.2.8 for required run sequence.) The CCB result must fall within +/- RL from zero. Each CCV and CCB analyzed must reflect the conditions of analysis of all associated samples. Sample results may only be reported when bracketed by valid ICV/CCV and ICB/CCB pairs.
- 9.8. Detection Limit Standard (CRA)-To verify linearity at the reporting limit, a CRA standard is run at the beginning of each sample analysis run after the ICV/ICB. The CRA standard mercury concentration is 0.2 ug/L. Recovery must be  $\pm 50\%$  of the true value, or the standard is either rerun or the problem corrected and the instrument recalibrated. (See Section 11.2.8 for the required run sequence.)
- 9.9. Method of Standard Addition (MSA) -This technique involves adding known amounts of standard to one or more aliquots of the sample prior to preparation. This technique compensates for a sample interferent that may enhance or depress the analyte signal, thus producing a different slope from that of the calibration standards. It will not correct for additive interferences, which cause a baseline shift. Refer to Section 11.2.9 for additional information on when full 4 point MSA is required as well as Appendix B for specific MSA requirements.

---

## 10. CALIBRATION AND STANDARDIZATION

- 10.1. Calibration standards must be processed through the preparation procedure as described in Section 11.1.
- 10.2. Due to the differences in preparation protocols separate calibration and calibration verification standards must be prepared for aqueous and solid matrices.
- 10.3. Calibration must be performed daily (every 24 hours) and each time the instrument is set up. The instrument calibration date and time must be included in the raw data.
- 10.4. Set up the instrument with the operating parameters recommended by the manufacturer and listed in Appendix F. Allow the instrument to become thermally stable before beginning calibration (approximately 30 minutes of warm-up is required).
- 10.5. Calibrate the instrument according to instrument manufacturer's instructions, using a minimum of five standards and a blank. One standard must be at the STL North Canton reporting limit. Analyze standards in ascending order beginning with the blank. Refer to Section 7 and Table I for additional information on preparing calibration standards and calibration levels.
- 10.6. The calibration curve must have a correlation coefficient of  $\geq 0.995$  or the instrument shall be stopped and recalibrated prior to running samples. Sample results can not be reported from a curve with an unacceptable correlation coefficient.
- 10.7. Refer to Section 9.0 for calibration verification procedures, acceptance criteria and corrective actions.

## 11. PROCEDURE

### 11.1. Sample Preparation:

11.1.1. All calibration and calibration verification standards (ICV, ICB, CCV, CCB) are processed through the digestion procedure as well as the field samples. Transfer 0, 0.2, 0.5, 1.0, 5.0 and 10.0 mL aliquots of the working standard (7.3) into a series of 100 ml class A volumetrics, then dilute to volume. For the ICV, use a 2.5 ml aliquot of the working standard. The ICV working standard must be made from a source other than that used for the calibration standards.

11.1.2. Transfer 100 mL of well-mixed sample or standard to a clean sample digestion bottle.

**Note:** Reduced sample volumes can be used as long as a representative sample can be obtained and the reagent levels are adjusted to maintain the same sample to reagent ratio. All samples and standards must be processed similarly.

**Note:** Spiking is done before the addition of acids or reagents.

11.1.3. Add 5 mL of concentrated  $\text{H}_2\text{SO}_4$  and 2.5 mL of concentrated  $\text{HNO}_3$  mixing after each addition.

11.1.4. Add 15 mL of potassium permanganate solution. For samples high in organic materials or chlorides, additional permanganate may be added. Shake and add additional portions of permanganate solution until a purple color persists for at least 15 minutes. If after the addition of up to 25-mL additional permanganate the color does not persist, sample dilution prior to reanalysis may be required.

**Note:** When performing analyses using automated vs. manual techniques the sample dilution resultant from the addition of more than the original aliquot of permanganate solution must be compensated for by the addition of the same volume of permanganate to all other associated samples and standards in the run. In instances, where this is not feasible, the addition of excess reagent can be addressed through mathematical correction of the results to account for the resultant dilution effect.

11.1.5. Add 8 mL of potassium persulfate solution and heat for two hours in a water bath at 90 - 95 °C.

11.1.6. Cool samples.

11.2. Sample Analysis:

11.2.1. Refer to the Appendix F of this SOP for detailed setup and operation protocols.

11.2.2. When ready to begin analysis, add 6 mL of sodium chloride-hydroxylamine hydrochloride solution to the samples to reduce the excess permanganate (the permanganate has been reduced when no purple color remains). Add this solution in 6-mL increments until the permanganate is completely reduced.

11.2.3. Automated determination: Follow instructions provided by instrument manufacturer.

11.2.4. Perform a linear regression analysis of the calibration standards by plotting maximum response of the standards vs. concentration of mercury. Determine the mercury concentration in the samples from the linear regression fit of the calibration curve. Calibration using computer or calculation based regression curve fitting techniques on concentration/response data is acceptable.

11.2.5. All measurements must fall within the defined calibration range to be valid. Dilute and reanalyze all samples for analytes that exceed the highest calibration standard.

11.2.6. If the sample results are negative and the absolute value of the negative result is greater than the reporting limit, the sample must be diluted and reanalyzed.

11.2.7. The samples must be allowed to cool to room temperature prior to analysis or a decrease in the response signal can occur.

11.2.8. The following run sequence is consistent with 7470A, CLP and 245.1.

Instrument Calibration

ICV

ICB

CRA

CCV

CCB

10 samples

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 17 of 41

---

CCV

CCB

Repeat sequence of 10 samples between CCV/CCB pairs as required to  
complete run.

CCV

CCB

Refer to Quality Control Section 9.0 and Table II (Appendix A) for Quality Control  
criteria to apply to Methods 7470A and 245.1.

11.2.9. For TCLP samples, full four point MSA will be required if all of the following  
conditions are met:

- 1) recovery of the analyte in the matrix spike is not at least 50%,
- 2) the concentration of the analyte does not exceed the regulatory level, and,
- 3) the concentration of the analyte is within 20% of the regulatory level.

The reporting and matrix spike levels for TCLP analyses are detailed in Table I  
(Appendix A). Appendix B provides guidance on performing MSA analyses.  
For TCLP mercury determinations, MSA spikes must be added prior to sample  
preparation.

- 11.3. To facilitate the early identification of QC failures and samples requiring rerun it is strongly  
recommended that sample data are reviewed periodically throughout the run.
- 11.4. Guidelines are provided in the appendices on procedures to minimize contamination of  
samples and standards, preventive maintenance and troubleshooting.
- 11.5. Analytical Documentation

- 
- 11.5.1. Record all analytical information in the analytical logbook/logsheet which may be in an electronic format, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.
- 11.5.2. All standards are logged into a department standard logbook. All standards are assigned an unique number for identification. Logbooks are reviewed by the supervisor or designee.
- 11.5.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.
- 11.5.4. Sample results and associated QC are entered into the LIMs after final technical review.
- 11.6. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.
- 11.7. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

## 12. DATA ANALYSIS AND CALCULATIONS

- 12.1. ICV percent recoveries are calculated according to the equation:

$$\%R = 100 \left( \frac{\text{Found(ICV)}}{\text{True(ICV)}} \right)$$

- 12.2. CCV percent recoveries are calculated according to the equation:

$$\%R = 100 \left( \frac{\text{Found(CCV)}}{\text{True(CCV)}} \right)$$

- 12.3. Matrix spike recoveries are calculated according to the following equation:

$$\% R = 100 \left( \frac{SSR - SR}{SA} \right)$$

Where:

SSR = Spike Sample Result

SR = Sample Result

SA = Spike Added

- 12.4. The relative percent difference (RPD) of matrix spike/matrix spike duplicates or sample duplicates are calculated according to the following equations:

$$RPD = 100 \left[ \frac{|MSD - MS|}{\left( \frac{MSD + MS}{2} \right)} \right]$$

Where:

MS = determined spiked sample concentration

MSD = determined matrix spike duplicate concentration

$$RPD = 100 \left[ \frac{|DU1 - DU2|}{\left( \frac{DU1 + DU2}{2} \right)} \right]$$

Where:

DU1 = Sample result

DU2 = Sample duplicate result

- 12.5. The final concentration for an aqueous sample is calculated as follows:

$$mg/L = C \times D$$

Where:

C = Concentration (mg/L) from instrument readout

D = Instrument dilution factor

- 12.6. The LCS percent recovery is calculated according to the following equation:



$$\%R = 100 \left( \frac{Found(LCS)}{True(LCS)} \right)$$

12.7. Appropriate factors must be applied to sample values if dilutions are performed.

12.8. Sample results should be reported with up to three significant figures in accordance with the STL North Canton significant figure policy.

### 13. METHOD PERFORMANCE

13.1. Each laboratory must have initial demonstration of performance data on file for each analyte of interest as described in Section 9.0.

13.2. Method performance is determined by the analysis of method blanks, laboratory control samples, matrix spike and matrix spike duplicate samples. The matrix spike recovery should fall within +/- 25 % and the matrix spike duplicates should compare within 20% RPD. The method blanks must meet the criteria in Section 9.3. The laboratory control sample should recover within 20% of the true value until in house limits are established.

13.3. Training Qualification:

The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

### 14. POLLUTION PREVENTION

14.1. This method allows for the proportional reduction of sample and reagent volumes to decrease waste generation.

### 15. WASTE MANAGEMENT

15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."

15.2. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by an annual refresher training.

15.3. Waste Streams Produced by the Method

15.3.1. **Acid Waste- Aqueous waste generated by the analysis.** Samples are disposed of in the acid waste drum located in the metals lab. The contents of the drum are neutralized and released to the POTW.

## 16. REFERENCES

### 16.1. References

16.1.1. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW-846, 3rd Edition, Final Update II, Revision I, September 1994, Method 7470A (Mercury).

16.1.2. "Methods for the Chemical Analysis of Water and Wastes", EPA-600/4-79-020, U.S.EPA, August 1983, Method 245.1.

16.1.3. U.S.EPA Statement of Work for Inorganics Analysis, ILMO3.0 and ILMO4.0.

16.1.4. Corporate Quality Management Plan (QMP), current version.

16.1.5. STL Laboratory Quality Manual (LQM), current version.

16.2. Associated SOPs and Policies, latest version

16.2.1. QA Policy, QA-003.

16.2.2. Glassware Washing, NC-QA-0014.

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018.

16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021.

16.2.5. Navy/Army SOP, NC-QA-0016.

16.2.6. Preparation and Analysis of Mercury in Aqueous Samples by Cold Vapor Atomic Absorption, Method 245.1 CLP-M, SOW ILM03.0 and SOW ILM04.0 SOP, CORP-MT-0006NC, current version.

## 17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . )

17.1. Modifications/Interpretations from reference method.

17.1.1. Modifications from both 7470A and 245.1.

17.1.1.1. The 200 series methods and Chapter 1 of SW846 specify the use of reagent water with a purity equivalent to ASTM Type II water. This SOP specifies the use of a Millipore DI system or equivalent to produce reagent water. This SOP requires that reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.

17.1.1.2. This SOP allows for the use of reduced sample volumes to decrease waste generation. Reagent levels are adjusted to maintain the same ratios as stated in the source methods. According to a letter from Robert Booth of EPA EMSL-Cinn to David Payne of EPA Region V, "Reduction in sample size and appropriate corresponding reduction in sample volume is not considered a significant change in the methodology."

17.1.2. Modifications from Method 7470A

17.1.2.1. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the MDL. This SOP states that the

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 23 of 41

---

method blank must not contain any analyte of interest at or above the reporting limit.

17.1.3. Modifications from 245.1

- 17.1.3.1. Method 245.1 Section 9.3 states concentrations should be reported as follows: Between 1 and 10 ug/L, one decimal; above 10 ug/L, to the nearest whole number. STL North Canton reports all Hg results under this SOP to two significant figures.

**Figure 1.** Aqueous Sample Preparation - Mercury





Figure 2. CVAA Mercury Analysis



877 896

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1  
APPENDIX A -TABLES

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 26 of 41

---

## APPENDIX A

### TABLES

**TABLE I. MERCURY REPORTING LIMITS, CALIBRATION STANDARD\*, QC  
STANDARD AND SPIKING LEVELS (MG/L)**

Standard Aqueous RL	0.0002
TCLP RL	0.002
Std 0	0
Std 1/CRA	0.0002
Std 2	0.0005
Std 3	0.001
Std 4	0.005
Std 5	0.010
ICV	0.0025
LCS/CCV	0.005
Aqueous MS	0.001
TCLP MS	0.005



**TABLE II. Summary Of Quality Control Requirements**

QC PARAMETER	FREQUENCY	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
ICV	Beginning of every analytical run.	90-110 % recovery.	Terminate analysis; Correct the problem; Recalibrate or reprep with calibration curve. (See Section 9.6).
ICB	Beginning of every analytical run, immediately following the ICV.	The result must be within +/- RL from zero.	Terminate analysis; Correct the problem; Recalibrate or reprep with calibration curve. (See Section 9.6).
CRA	Beginning of every analytical run following the ICB and prior to sample analyses.	50-150% recovery	Rerun to verify; or correct problem and recalibrate or reprep with the calibration curve (See Sec. 9.8)
CCV	Every 10 samples and at the end of the run.	80 - 120 % recovery.	Terminate analysis; Correct the problem; Recalibrate and rerun all samples not bracketed by acceptable CCV or reprep with calibration curve. (See Section 9.7).
CCB	Immediately following each CCV.	The result must be within +/- RL from zero.	Terminate analysis; Correct the problem; Recalibrate and rerun all samples not bracketed by acceptable CCB or reprep with calibration curve. (See Section 9.7).

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1  
APPENDIX A -TABLES

877 899

SOP No. CORP-MT-0005NC  
Revision No. 2.5  
Revision Date: 01/08/04  
Page 29 of 41

Method Blank	One per sample preparation batch of up to 20 samples.	The result must be less than or equal to the RL. Sample results greater than 20x the blank concentration are acceptable. Samples for which the contaminant is < RL do not require redigestion (See Section 9.3).	Redigest and reanalyze samples.  Note exceptions under criteria section.  See Section 9.3 for additional requirements.
--------------	---	--	--

TABLE II. Summary of Quality Control Requirements (Continued)

QC PARAMETER	FREQUENCY	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
Laboratory Control Sample (LCS)	One per sample preparation batch of up to 20 samples.	Aqueous LCS must be within 80 - 120% recovery or in-house control limits.	Terminate analysis; Correct the problem; Redigest and reanalyze all samples associated with the LCS (see Section 9.4).
Matrix Spike	One per sample preparation batch of up to 20 samples.	75 - 125 % recovery or in-house control limits. If the MS/MSD is out for an analyte, it must be in control in the LCS.	In the absence of client specific requirements, flag the data; no flag required if the sample level is > 4x the spike added. (see Section 9.5)  For TCLP see Section 11.2.9
Matrix Spike Duplicate	See Matrix Spike	75 - 125 % recovery or in-house control limits; RPD $\leq$ 20%. (See MS)	See Corrective Action for Matrix Spike.

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1  
APPENDIX B - MSA GUIDANCE

877 90i  
SOP No. CORP-MT-0005NC  
Revision No. 2.5  
Revision Date: 01/08/04  
Page 31 of 41

---

**APPENDIX B**  
**MSA GUIDANCE**


## APPENDIX B. MSA GUIDANCE

### Method of Standard Addition

Four equal volume aliquots of sample are measured and known amounts of standards are added to three aliquots. The fourth aliquot is the unknown and no standard is added to it. The concentration of standard added to the first aliquot should be 50% of the expected concentration. The concentration of standard added to the second aliquot should be 100% of the expected concentration and the concentration of standard added to the third aliquot should be 150% of the expected concentration. The volume of the unspiked and spiked aliquots should be the same (i.e., the volume of the spike added should be negligible in relation to the volume of sample).

To determine the concentration of analyte in the sample, the absorbance (or response) of each solution is determined and a linear regression performed. On the vertical axis the absorbance (or response) is plotted versus the concentrations of the standards on the horizontal axis using 0 as the concentration of the unspiked aliquot. An example plot is shown in Figure 1. When the resulting line is extrapolated back to zero absorbance, the point of interception of the horizontal axis is the concentration of the unknown. Calculate the correlation coefficient ( $r$ ) and the x-intercept (where  $y=0$ ) of the curve. The concentration in the digestate is equal to the negative x-intercept.

Figure 1



- For the method of standard additions to be correctly applied, the following limitations must be taken into consideration.
- The plot of the sample and standards must be linear over the concentration range of concern. For best results, the slope of the curve should be similar to that of a plot of the aqueous standard curve.
- The effect of the interference should not vary as the ratio of the standard added to the sample matrix changes.

**APPENDIX C**  
**TROUBLESHOOTING GUIDE**

#### APPENDIX C. TROUBLESHOOTING GUIDE

Problem	Possible Cause
Poor or No Absorbance or Sensitivity Check failed	Incorrect wavelength Dirty windows Window loose Etched or dirty optics Wrong lamp Bad lamp Not enough or no sample introduced Empty sample cup Incorrectly made standards Gas leak EDL power supply set on "Continuous"
Erratic Readings	Source lamp not aligned properly Lamp not prewarmed Injection tip partially clogged Contaminated reagents Contaminated glassware Drying tube saturated Bad lamp Injection tip hitting outside of tube Injection tip coated or not set properly Leak in sample tubing Power fluctuations Air bubbles in tubing
EDL Won't Light	Lamp cable not plugged in Lamp power set at 0 Lamp is dead Power supply fuse is blown Short in cord
Standards reading twice or half normal absorbance or concentration	Incorrect standard used Incorrect dilution performed Dirty cell
Background Correction Light Blinking	Background screen or attenuator faulty



877 906

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846 METHOD 7470A AND MCAWW METHOD 245.1  
APPENDIX D - CONTAMINATION CONTROL GUIDELINES

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 36 of 41

## **APPENDIX D**

### **CONTAMINATION CONTROL GUIDELINES**

## APPENDIX D. CONTAMINATION CONTROL GUIDELINES

### The following procedures are strongly recommended to prevent contamination:

All work areas used to prepare standards and spikes should be cleaned before and after each use.

All glassware should be washed with detergent and tap water and rinsed with 1:1 nitric acid followed by deionized water.

Proper laboratory housekeeping is essential in the reduction of contamination in the metals laboratory. All work areas must be kept scrupulously clean.

Powdered gloves must not be used in the metals laboratory since the powder contains silica and zinc, as well as other metallic analytes. Alternatively, vinyl or nitrile gloves should be used in the metals laboratory.

Glassware should be periodically checked for cracks and etches and discarded if found. Etched glassware can cause cross contamination of any metallic analytes.

### The following are helpful hints in the identification of the source of contaminants:

Reagents or standards can contain contaminants or be contaminated with the improper use of a pipette.

Improper cleaning of glassware can cause contamination.

Separate glassware if an unusually high sample is analyzed and soak with sulfuric acid prior to routine cleaning.

**877 908**

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846 METHOD 7470A AND MCAWW METHOD 245.1  
APPENDIX E - PREVENTIVE MAINTENANCE

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 38 of 41

---

## **APPENDIX E**

### **PREVENTIVE MAINTENANCE**

## APPENDIX E. PREVENTIVE MAINTENANCE

A maintenance log is used to record when maintenance is performed on instruments. When an instrument problem occurs indicate the date, time and instrument number, then identify the problem and corrective action in the maintenance log.

The following procedures are required to ensure that that the instrument is fully operational.

### Cold Vapor Atomic Absorption (Leeman PS 200II)

Daily	As Needed
Check nitrogen flow.	Check Hg lamp intensity.
Check tubing.	Clean lens.
Check drain.	Check aperture.
	Replace drying tube.
	Change Hg lamp.
	Check liquid/gas separator.

877 910

PREPARATION AND ANALYSIS OF MERCURY IN AQUEOUS  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7470A AND MCAWW METHOD 245.1  
APPENDIX F – INSTRUMENT SET-UP

SOP No. CORP-MT-0005NC

Revision No. 2.5

Revision Date: 01/08/04

Page 40 of 41

---

## APPENDIX F

### INSTRUMENT SET UP

## Hg Analysis (Leeman PS200II)

### SYSTEM INITIALIZATION AND WARM UP

1. F1 Menu
2. Instrument
  - a. TASKMASTER
  - b. #4 Wake System Up Enter

The warming up period takes approximately 10 minutes.

### TO SET UP INSTRUMENT FOR ANALYSIS

1. F1 Menu
2. Autosampler
  - A. Rack Entry
  - B. Edit (ex. Rack 1), Enter
  - C. Cup ID - Enter (clears sample #'s)
  - D. Extended ID- type in matrix of sample (water or solid), Enter
  - E. Press Insert Key and move cursor with arrows to cup ID and begin typing labels.
  - F. F3 Print Screen
3. Press F2 Macro key and type in analyst's first name - Enter
  - A. Enter folder name - ex. HG0306, Enter If folder does not exist, type Y - Enter.
  - B. Type in - "Rack 1", "Rack 2" etc. , Enter.
  - C. Type in : FROM CUP TO CUP  
Ex. = 1 30

Do the same for position 2 if needed If not needed, you must press Enter 3 times to begin analysis.

Controlled Copy

Copy No. \_\_\_\_\_

Implementation Date: 12/3/02SOP No. CORP-MT-0007NCRevision No. 2.4Revision Date: 10/28/02Page 1 of 43**STL NORTH CANTON STANDARD OPERATING PROCEDURE**

**TITLE: PREPARATION AND ANALYSIS OF MERCURY IN SOLID SAMPLES BY  
COLD VAPOR ATOMIC ABSORPTION SPECTROSCOPY, SW846 7471A AND MCAWW  
245.5**

(SUPERSEDES: REVISION 2.3, DATED 05/15/01)

Prepared by:	<u>Karen R. Counts</u>	<u>11-25-02</u>
		Date
Approved by:	<u>Robert O'Meara</u>	<u>11/21/02</u>
	Technical Specialist	Date
Approved by:	<u>Beth Lambert</u>	<u>11/22/02</u>
	Quality Assurance Manager	Date
Approved by:	<u>Dorothy S Budd</u>	<u>11/22/02</u>
	Environmental Health and Safety Coordinator	Date
Approved by:	<u>Anthony J. Orlando</u>	<u>11/22/02</u>
	Laboratory Director	Date
Approved by:	<u>Phil Bue.</u>	<u>11/11/02</u>
	Corporate Technology	Date

**Proprietary Information Statement:**

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846

SOP No. CORP-MT-0007NC

Revision No. 2.4

Revision Date: 10/28/02

Page 2 of 42

**TABLE OF CONTENTS**

1. SCOPE AND APPLICATION .....	3
2. SUMMARY OF METHOD .....	3
3. DEFINITIONS .....	3
4. INTERFERENCES .....	3
5. SAFETY .....	5
6. EQUIPMENT AND SUPPLIES .....	6
7. REAGENTS AND STANDARDS .....	7
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE .....	9
9. QUALITY CONTROL .....	9
10. CALIBRATION AND STANDARDIZATION .....	13
11. PROCEDURE .....	13
12. DATA ANALYSIS AND CALCULATIONS .....	17
13. METHOD PERFORMANCE .....	19
14. POLLUTION PREVENTION .....	20
15. WASTE MANAGEMENT .....	20
16. REFERENCES .....	20
17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . ) .....	21

**LIST OF APPENDICES:**

APPENDIX A - TABLES .....	27
APPENDIX B - STL NORTH CANTON Hg DATA REVIEW CHECKLIST .....	31
APPENDIX C - MSA GUIDANCE .....	33
APPENDIX D - TROUBLESHOOTING GUIDE .....	36
APPENDIX E - CONTAMINATION CONTROL GUIDELINES .....	38
APPENDIX F - PREVENTATIVE MAINTENANCE .....	40



PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846

SOP No. CORP-MT-0007NC

Revision No. 2.4

Revision Date: 10/28/02

Page 3 of 42

---

APPENDIX G – INSTRUMENT SET-UP.....40

**1. SCOPE AND APPLICATION**

- 1.1. This procedure describes the preparation and analysis of mercury (Hg, CAS # 7439-97-6) by Cold Vapor Atomic Absorption Spectroscopy (CVAA) using SW-846 Method 7471A and MCAWW Method 245.5.
- 1.2. The associated LIMs method code is O9.
- 1.3. CVAA analysis provides for the determination of total mercury (organic and inorganic). The combination of the oxidants, potassium permanganate and potassium persulfate, has been found to give 100% recovery with both types of compounds. Detection limits, sensitivity and optimum concentration ranges for mercury analysis will vary with the matrices, instrumentation and volume of sample used.
- 1.4. Methods 7471A and 245.5 are applicable to the preparation and analysis of mercury in soils, sediments, bottom deposits, wastes, wipes and sludge-type materials. All matrices require sample preparation prior to analysis.
- 1.5. The STL North Canton reporting limit for mercury in solid matrices is 0.033 mg/kg based on a 0.6 g sample aliquot (wet weight).

**2. SUMMARY OF METHOD**

- 2.1. This SOP describes a technique for the determination of mercury in solution. The procedure is a physical method based on the absorption of radiation at 253.7 nm by mercury vapor. A representative portion of the sample is digested in hydrochloric and nitric acids. Organic mercury compounds are oxidized with potassium permanganate and potassium persulfate and the mercury reduced to its elemental state with stannous chloride and aerated from solution in a closed system. The mercury vapor passes through a cell positioned in the light path of an atomic absorption spectrophotometer. Absorbance is measured as a function of mercury concentration. Concentration of the analyte in the sample is determined by comparison of the sample absorbance to the calibration curve (absorbance vs. concentration).

### 3. DEFINITIONS

- 3.1. Total Metals: The concentration determined on an unfiltered sample following digestion.

### 4. INTERFERENCES

Chemical and physical interferences may be encountered when analyzing samples using this method.

- 4.1. Potassium permanganate which is used to breakdown organic mercury compounds also eliminates possible interferences from sulfide. Concentrations as high as 20 mg/L of sulfide as sodium sulfide do not interfere with the recovery of inorganic mercury from reagent water.
- 4.2. Copper has also been reported to interfere; however, copper concentrations as high as 10 mg/L had no effect on the recovery of mercury from spiked samples.
- 4.3. Chlorides can cause a positive interference. Samples high in chlorides require additional permanganate (as much as 25 mL) because, during the oxidation step, chlorides are converted to free chlorine, which also absorbs radiation at 253.7 nm. Care must be taken to ensure that free chlorine is absent before the mercury is reduced and swept into the cell. This is accomplished by adding excess hydroxylamine reagent (25 mL) and purging the sample headspace before stannous chloride is added. Both inorganic and organic mercury spikes have been quantitatively recovered from seawater using this technique.

**Note:** Sufficient addition of permanganate is apparent when the purple color persists at least 15 minutes. Some samples may require dilution prior to digestion due to extremely high concentrations of chloride.

- 4.4. Interference from certain volatile organic materials that absorb at this wavelength may also occur. If suspected, a preliminary run without stannous chloride can determine if this type of interference is present. While the possibility of absorption from certain organic substances present in the sample does exist, this problem is not routinely encountered. This is mentioned only to caution the analyst of the possibility. If this condition is found to exist, the mercury concentration in the sample can be determined by subtracting the result of the sample run without the reducing reagent (stannous chloride) from that obtained with the reducing reagent.
- 4.5. Samples containing high concentrations of oxidizable organic materials, as evidenced by high COD levels, may not be completely oxidized by this procedure. When this occurs the recovery of mercury will be low. The problem can be eliminated by reducing the volume of original sample used.

- 4.6. The most common interference is laboratory contamination which may arise from impure reagents, dirty glassware, improper sample transfers, dirty work areas, etc. Be aware of potential sources of contamination and take appropriate measures to minimize or avoid them.

## 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL North Canton associates.
- 5.2. Eye protection that satisfies ANSI Z87.1 (as per the Chemical Hygiene Plan), laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:
- 5.3.1. The following materials are known to be **corrosive**:
- hydrochloric acid, nitric acid and sulfuric acid.
- 5.3.2. The following materials are known to be **oxidizing agents**:
- nitric acid, potassium permanganate, potassium persulfate and magnesium perchlorate.
- 5.3.3. Mercury is a highly toxic element that must be handled with care. The analyst must be aware of the handling and clean-up techniques before working with mercury. Since mercury vapor is toxic, precaution must be taken to avoid its inhalation, ingestion or absorption through skin. All lines should be checked for leakage and the mercury vapor must be vented into a hood or passed through a mercury absorbing media such as:

---

5.3.3.1. Equal volumes of 0.1 M  $\text{KMnO}_4$  and 10%  $\text{H}_2\text{SO}_4$ , or

5.3.3.2. Iodine, 0.25%, in a 3% KI solution.

5.3.4. Magnesium sulfate is known to be a reproductive toxin (mutagen).

- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**. Therefore, unless they are known to be non-hazardous, all samples should be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.
- 5.6. Do not look directly into the beam of the Hg lamp. The UV light that these lamps radiate is harmful to the eyes.
- 5.7. Cylinders of compressed gas must be handled with caution, in accordance with local regulations. It is recommended that, wherever possible, cylinders be located outside the laboratory and the gas led to the instrument through approved lines.
- 5.8. The CVAA apparatus must be properly vented to remove potentially harmful fumes generated during sample analysis.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Temperature controlled water bath (capable of maintaining temperature of 90- 95 °C).
- 6.2. Atomic Absorption Spectrophotometer equipped with:
- 6.2.1. Absorption Cell with quartz end windows perpendicular to the longitudinal axis. Dimensions of the cell must result in sufficient sensitivity to meet the SOP defined reporting limit. The quartz windows must be maintained to provide accurate measurements. Any scratches or fingerprints can alter the absorption of UV radiation.
- 6.2.2. Mercury specific hollow cathode lamp (HCL) or electrodeless discharge lamp (EDL).

- 6.2.3. Peristaltic pump which can deliver 1 L/min air.
- 6.2.4. Flowmeter capable of measuring an airflow of 1 L/min.
- 6.2.5. Recorder or Printer.
- 6.2.6. Aeration Tubing: A straight glass frit having a coarse porosity and Tygon tubing is used for the transfer of mercury vapor from the sample bottle to the absorption cell and return.
- 6.2.7. Drying device (a drying tube containing magnesium perchlorate or magnesium sulfate and/or a lamp with a 60 W bulb) to prevent condensation in cell. The lamp is positioned to shine on the absorption cell maintaining the air temperature in the cell about 10 °C above room temperature. Other drying devices that achieve the same purpose are also acceptable (i.e., Gortex filter).

**Note:** Instruments designed specifically for the measurement of mercury using the cold vapor technique may be substituted for the atomic absorption spectrophotometer.

- 6.3. BOD bottles or equivalent.
- 6.4. Nitrogen or argon gas supply, welding grade or equivalent.
- 6.5. Calibrated automatic pipettes or Class A glass volumetric pipettes.
- 6.6. Class A volumetric flasks.
- 6.7. Top-loading balance, capable of reading up to two decimal places.
- 6.8. Thermometer (capable of accurate readings at 95 °C).
- 6.9. Disposable cups or tubes.

## 7. REAGENTS AND STANDARDS

- 7.1. Reagent water must be produced by a Millipore DI system or equivalent. Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.
- 7.2. Stock (10 ppm) mercury standards (in 10% HNO<sub>3</sub>) are purchased as custom STL North Canton solutions. All standards must be stored in FEP fluorocarbon or previously unused

---

polyethylene or polypropylene bottles. Stock standard solutions must be replaced prior to the expiration date provided by the manufacturer. If no expiration date is provided, the stock solutions may be used for up to one year and must be replaced sooner if verification from an independent source indicates a problem.

- 7.3. Working mercury standard (0.1 ppm): Take 1 mL of the stock mercury standard (7.2) and dilute to 100 mL with reagent water. The working mercury standard must be made daily and must be prepared in a matrix of 0.15% HNO<sub>3</sub>. This acid (150 uL of concentrated HNO<sub>3</sub>) must be added to the flask/bottle before the addition of the stock standard aliquot.

- 7.4. The calibration standards must be prepared fresh daily from the working standard (7.3) by transferring 0, 0.2, 0.5, 1.0, 5.0 and 10.0 mL aliquots of the working mercury standard into sample preparation bottles and proceeding as specified in Section 11.1

**Note:** Alternate approaches to standard preparation may be taken and alternate volumes of standard may be prepared as long as the accuracy and final standard concentrations as detailed in Table I are maintained. For example, automated mercury systems do not require 100 mL of standard and therefore smaller volumes may be generated to reduce waste generation.

- 7.5. The initial calibration verification standard must be made from a different stock solution than that of the calibration standards.
- 7.6. Refer to Table I (Appendix A) for details regarding the working standard concentrations for calibration, calibration verification and spiking solutions. All standards must be processed through the entire analytical procedure including sample preparation.

- 7.7. Nitric acid (HNO<sub>3</sub>), concentrated, trace metal grade or better.

**Note:** If a high reagent blank is obtained, it may be necessary to distill the nitric acid.

- 7.8. Sulfuric acid (H<sub>2</sub>SO<sub>4</sub>), concentrated, trace metal grade or better.

- 7.9. Hydrochloric acid (HCl), concentrated, trace metal grade or better.

- 7.10. Aqua Regia: Prepare immediately before use by carefully adding three volumes of concentrated HCl to one volume of concentrated HNO<sub>3</sub>.

- 7.11. Stannous sulfate solution: Add 25 g of stannous sulfate to 250 mL of 0.5 N sulfuric acid. This mixture is a suspension and should appear cloudy. This solution should be made daily and should be stirred continuously during use.

**Note:** Stannous chloride may be used in place of stannous sulfate. Prepare the stannous chloride solution according to the recommendations provided by the instrument manufacturer.

- 7.12. Sodium chloride-hydroxylamine hydrochloride solution: Add 12 g of sodium chloride and 12 g of hydroxylamine hydrochloride to every 100 mL of reagent water.

**Note:** Hydroxylamine sulfate may be used in place of hydroxylamine hydrochloride.

- 7.13. Potassium permanganate, 5% solution (w/v): Dissolve 5 g of potassium permanganate for every 100 mL of reagent water.

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

- 8.1. Sample holding time for mercury is 28 days from time of collection to the time of sample analysis.
- 8.2. Soil samples do not require preservation but must be stored at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$  until the time of analysis.

## 9. QUALITY CONTROL

Table II (Appendix A) provides a summary of quality control requirements including type, frequency, acceptance criteria and corrective action.

### 9.1. Initial Demonstration of Capability

Prior to the analysis of any analyte using 7471A or the 245.5, the following requirements must be met.

9.1.1. Method Detection Limit (MDL) - An MDL must be determined for each analyte/matrix prior to the analysis of any samples. The MDL is determined using seven replicates of reagent water, spiked with all the analytes of interest, that have been carried through the entire analytical procedure. MDLs must be redetermined on an annual basis in accordance with 40 CFR Part 136 Appendix B requirements. The spike level must be between the calculated MDL and 10X the MDL to be valid. The result of the MDL determination must be below the STL North Canton reporting limit.

9.1.2. Initial Demonstration Study - This requires the analysis of four QC check samples. The QC check sample is a well characterized laboratory generated sample used to

---

monitor method performance. The results of the initial demonstration study must be acceptable before analysis of samples may begin.

9.1.2.1. Four aliquots of the check sample (LCS) are prepared and analyzed using the procedures detailed in this SOP and the determinative SOPs.

- 9.2. Preparation Batch - A group of up to 20 samples that are of the same matrix and are processed together using the same procedures and reagents. The preparation batch must contain a method blank, a LCS and a matrix spike/matrix spike duplicate. In some cases, at client request, it may be appropriate to process a matrix spike and sample duplicate in place of the MS/MSD. If clients specify specific samples for MS/MSD, the batch may contain multiple MS/MSD pairs.
- 9.3. Sample Count - Laboratory generated QC samples (method blanks, LCS, MS/MSD) are not included in the sample count for determining the size of a preparation batch.
- 9.4. Method Blank (MB) - One method blank must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting or at or above 5% of the measured concentration of that analyte in associated samples, whichever is higher (sample result must be a minimum of 20 times higher than the blank contamination level).
- Repreparation and reanalysis of all samples associated with an unacceptable method blank is required when reportable concentrations are determined in the samples (see exception noted above).
  - If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. **Such action must be taken in consultation with the client and must be addressed in the project narrative.**
  - If the above criteria are not met and reanalysis is not possible, then the sample data must be qualified. **This anomaly must be addressed in the project narrative and the client must be notified.**
- 9.5. Laboratory Control Sample (LCS) - One aqueous LCS must be processed with each preparation batch. The LCS is used to monitor the accuracy of the analytical process. On-going monitoring of the LCS results provides evidence that the laboratory is performing the



---

method within acceptable accuracy and precision guidelines. The LCS must be carried through the entire analytical procedure. If the LCS is outside established control limits the system is out of control and corrective action must occur. Until in-house control limits are established, a control limit of 70-130% recovery must be applied.

- In the instance where the LCS recovery is > 130% and the sample results are < RL, the data may be reported with qualifiers. Such action must be taken in consultation with the client and must be addressed in the case narrative.”
- In the event that an MS/MSD analysis is not possible, a Laboratory Control Sample Duplicate (LCSD) must be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.
- Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.

9.6. Matrix Spike/Matrix Spike Duplicate (MS/MSD) - One MS/MSD pair must be processed for each preparation batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Samples identified as field blanks cannot be used for MS/MSD analysis. Spiking levels are provided in Table I (Appendix A).

- If analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. Until in-house control limits are established, a control limit of 70 - 130 % recovery and 20% RPD must be applied to the MS/MSD. If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include repreparation and reanalysis of the batch. MS/MSD results which fall outside the control limits must be addressed in the narrative.
- If the native analyte concentration in the MS/MSD exceeds 4 times the spike level for that analyte, the recovery data are reported as NC (i.e., not calculated). If the reporting software does not have the ability to report NC then the actual recovery must be reported and narrated as follows: “Results outside of limits do not

---

necessarily reflect poor method performance in the matrix due to high analyte concentrations in the sample relative to the spike level.”

- If an MS/MSD is not possible due to limited sample volume, then a laboratory control sample duplicate (LCSD) should be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.

- 9.7. Initial Calibration Verification (ICV/ICB) - Calibration accuracy is verified by analyzing a second source standard (ICV). The ICV result must fall within 10% of the true value for that solution. An ICB is analyzed immediately following the ICV to monitor low level accuracy and system cleanliness. The ICB result must fall within +/- the reporting limit (RL) from zero. If either the ICV or ICB fail to meet criteria, the analysis should be terminated, the problem corrected and the instrument recalibrated. (See Section 11.2.10 and Section 11.2.11 for required run sequence). If the cause of the ICV or ICB failure was not directly instrument related the corrective action will include re-preparation of the ICV, ICB, CRA, CCV and CCB with the calibration curve.
- 9.8. Continuing Calibration Verification (CCV/CCB) - Calibration accuracy is monitored throughout the analytical run through the analysis of a known standard after every 10 samples. The CCV must be a mid-range standard at a concentration other than that of the ICV. The CCV result must fall within 20% of the true value for that solution. A CCB is analyzed immediately following each CCV. (See Section 11.2.10 and 11.2.11 for required run sequence.) The CCB result must fall within +/- RL from zero. Each CCV and CCB analyzed must reflect the conditions of analysis of all associated samples. Sample results may only be reported when bracketed by valid ICV/CCV and ICB/CCB pairs.
- 9.9. Method of Standard Addition (MSA) - This technique involves adding known amounts of standard to one or more aliquots of the sample prior to preparation. This technique compensates for a sample interferent that may enhance or depress the analyte signal, thus producing a different slope from that of the calibration standards. It will not correct for additive interferences which cause a baseline shift. Refer to Section 11.2.12 for additional information on when full 4 point MSA is required as well as Appendix C for specific MSA requirements.

## 10. CALIBRATION AND STANDARDIZATION

- 10.1. Calibration standards must be processed through the preparation procedure as described in Section 11.1.
- 10.2. Due to the differences in preparation protocols separate calibration and calibration verification standards must be prepared for aqueous and solid matrices.

- 10.3. Calibration must be performed daily (every 24 hours) and each time the instrument is set up. The instrument calibration date and time must be included in the raw data.
- 10.4. Set up the instrument with the operating parameters recommended by the manufacturer. Allow the instrument to become thermally stable before beginning calibration (approximately 30 minutes of warm-up is required). Refer to the facility specific instrument SOP and CVAA instrument manual for detailed setup and operation protocols.
- 10.5. Calibrate the instrument according to instrument manufacturer's instructions, using a minimum of five standards and a blank. One standard must be at the STL North Canton reporting limit. Analyze standards in ascending order beginning with the blank. Refer to Section 7.5 and Table I for additional information on preparing calibration standards and calibration levels.
- 10.6. The calibration curve must have a correlation coefficient of  $\geq 0.995$  or the instrument shall be stopped and recalibrated prior to running samples. Sample results can not be reported from a curve with an unacceptable correlation coefficient.
- 10.7. Refer to Section 9.0 for calibration verification procedures, acceptance criteria and corrective actions.

## 11. PROCEDURE

### 11.1. Standard and Sample Preparation:

- 11.1.1. All calibration and calibration verification standards (ICV, ICB, CCV, CCB) are processed through the digestion procedure as well as the field samples.
- 11.1.2. Transfer 0, 0.2, 0.5, 1.0, 5.0 and 10.0 mL aliquots of the working standard (7.3) into a series of sample digestion bottles. For the ICV, transfer a 2.5 ml aliquot of the working standard. The ICV working standard must be made from a source other than that used for the calibration standards.  
  
**Note:** Alternate volumes of standard may be prepared as long as the accuracy and final standard concentrations as detailed in Table I are maintained.
- 11.1.3. Add reagent water to each standard bottle to make a total volume of 10 mL. Continue preparation as described under 11.1.5 below.
- 11.1.4. Transfer triplicate 0.2 g portions of a well mixed sample into a clean sample digestion bottle. Add 10mL of reagent water to each standard bottle. Continue

preparation as described under 11.1.5.

11.1.5. Water Bath protocol:

11.1.5.1. To each **standard** bottle: Add 5 mL of aqua regia.  
To each **sample** bottle: Add 10 mL of reagent water and 5 mL  
of aqua regia.

11.1.5.2. Heat for 2 minutes in a water bath at 90 - 95 ° C.

11.1.5.3. Add 40 mL of distilled water.

11.1.5.4. Add 15 mL of potassium permanganate solution.

11.1.5.5. Heat for 30 minutes in the water bath at 90 - 95 °C.

11.1.5.6. Cool.

11.1.5.7. Add 6 mL of sodium chloride-hydroxylamine sulfate solution to  
reduce the excess permanganate.

11.1.5.8. To each **standard** bottle: Add 50 mL of reagent water.  
To each **sample** bottle: Add 50 mL of reagent water.

11.1.5.9. Continue as described under Section 11.2.

11.2. Sample Analysis:

11.2.1. Because of differences between various makes and models of CVAA  
instrumentation, no detailed operating instructions can be provided. Refer to the  
facility specific instrument operating SOP and the CVAA instrument manual for  
detailed setup and operation protocols.

11.2.2. All labs are required to detail the conditions/programs utilized for each instrument  
within the facility specific instrument operation SOP.

11.2.3. Manual determination:

11.2.3.1. Treating each sample individually, purge the head space of the sample  
bottle for at least one minute.

- 11.2.3.2. Add 5 mL of stannous chloride solution and immediately attach the bottle to the aeration apparatus.
- 11.2.3.3. Allow the sample to stand quietly without manual agitation while the sample is aerated (1 L/min flow). Monitor the sample absorbance during aeration. When the absorbance reaches a maximum and the signal levels off, open the bypass valve and continue aeration until the absorbance returns to its baseline level. Close the bypass valve and remove the aeration device.
- 11.2.3.4. Place the aeration device into 100 mLs of 1% HNO<sub>3</sub> and allow to bubble rinse until the next sample is analyzed.
- 11.2.4. Automated determination: Refer to Appendix G for instrument setup and operation.
- 11.2.5. Perform a linear regression analysis of the calibration standards by plotting maximum response of the standards vs. ug of mercury. Determine the mercury concentration in the samples from the linear regression fit of the calibration curve. Calibration using computer or calculation based regression curve fitting techniques on concentration/response data is acceptable.
- 11.2.6. All measurements must fall within the defined calibration range to be valid. Dilute and reanalyze all samples for analytes that exceed the highest calibration standard.
- 11.2.7. If the sample results are negative and the absolute value of the negative result is greater than the reporting limit, the sample must be diluted and reanalyzed.
- 11.2.8. The samples must be allowed to cool to room temperature prior to analysis or a decrease in the response signal can occur.
- 11.2.9. Baseline correction is acceptable as long as it is performed after every sample or after the CCV and CCB; resloping is acceptable as long as it is immediately preceded and followed by a compliant CCV and CCB.
- 11.2.10. The following analytical sequence must be used with 7471A and 245.5:

Instrument Calibration

ICV

ICB

Maximum 10 samples

CCV

CCB

Repeat sequence of 10 samples between CCV/CCB pairs as required to  
complete run

CCV

CCB

Refer to Quality Control Section 9.0 and Table II (Appendix A) for quality  
control criteria to apply to Methods 7471A and 245.5.

**Note:** Samples include the method blank, LCS, MS, MSD, duplicate, field  
samples and sample dilutions.

- 11.2.11. The following run sequence is consistent with 7471A, CLP and 245.5 and may  
be used as an alternate to the sequence in 11.2.10. This run sequence is  
recommended if multiple method requirements must be accommodated in one  
analytical run:

Instrument Calibration

ICV

ICB

CRA\*

CCV

CCB

10 samples

CCV

CCB

Repeat sequence of 10 samples between CCV/CCB pairs as required to  
complete run.

CCV

CCB

Refer to the appropriate CLP SOP (CORP-MT-0008) for quality control  
requirements for QC samples.

\* Refer to the CLP SOP for information on the CRA.

11.2.12. For TCLP samples, full four point MSA will be required if all of the following conditions are met:

- 1) recovery of the analyte in the matrix spike is not at least 50%,
- 2) the concentration of the analyte does not exceed the regulatory level, and,
- 3) the concentration of the analyte is within 20% of the regulatory level.

Appendix E provides guidance on performing MSA analyses. For TCLP mercury determinations, MSA spikes must be added prior to sample preparation.

- 11.3. To facilitate the early identification of QC failures and samples requiring rerun it is strongly recommended that sample data be reviewed periodically throughout the run.
- 11.4. Guidelines are provided in the appendices on procedures to minimize contamination of samples and standards, preventive maintenance and troubleshooting.
- 11.5. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.
- 11.6. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

## 12. DATA ANALYSIS AND CALCULATIONS

12.1. ICV percent recoveries are calculated according to the equation:

$$\%R = 100 \left( \frac{Found(ICV)}{True(ICV)} \right)$$

12.2. CCV percent recoveries are calculated according to the equation:

$$\%R = 100 \left( \frac{Found(CCV)}{True(CCV)} \right)$$

12.3. Matrix spike recoveries are calculated according to the following equation:

---

$$\% R = 100 \left( \frac{SSR - SR}{SA} \right)$$

Where:

SSR = Spike Sample Result

SR = Sample Result

SA = Spike Added

- 12.4. The relative percent difference (RPD) of matrix spike/matrix spike duplicates or sample duplicates are calculated according to the following equations:

$$RPD = 100 \left[ \frac{|MSD - MS|}{\left( \frac{MSD + MS}{2} \right)} \right]$$

Where:

MS = determined spiked sample concentration

MSD = determined matrix spike duplicate concentration

$$RPD = 100 \left[ \frac{|DU1 - DU2|}{\left( \frac{DU1 + DU2}{2} \right)} \right]$$

Where:

DU1 = Sample result

DU2 = Sample duplicate result

- 12.5. For automated determinations, the final concentration determined in solid samples when reported on a dry weight basis is calculated as follows:

$$mg/kg, dry weight = (C \times V \times D) / (W \times S)$$

Where:

C = Concentration (ug/L) from instrument readout

V = Volume of digestate (L)

D = Instrument dilution factor

W = Weight in g of wet sample digested

S = Percent solids/100



---

**Note:** A Percent Solids determination must be performed on a separate aliquot when dry weight concentrations are to be reported. If the results are to be reported on a wet weight basis, the "S" factor should be omitted from the above equation.

- 12.6. For manual (total) determinations, the final concentration determined in solid samples when reported on a dry weight basis is calculated as follows:

$$mg/kg, dry weight = (C)/(W \times S)$$

Where:

C = Concentration (ug) from instrument readout

W = Weight in g of wet sample digested

S = Percent solids/100

**Note:** A Percent Solids determination must be performed on a separate aliquot when dry weight concentrations are to be reported. If the results are to be reported on a wet weight basis, the "S" factor should be omitted from the above equation.

- 12.7. The LCS percent recovery is calculated according to the following equation:

$$\%R = 100 \left( \frac{Found(LCS)}{True(LCS)} \right)$$

- 12.8. Sample results should be reported with up to three significant figures in accordance with the STL North Canton significant figure policy.

### 13. METHOD PERFORMANCE

- 13.1. Each laboratory must have initial demonstration of performance data on file for each analyte of interest as described in Section 9.1.
- 13.2. Method performance is determined by the analysis of method blank, laboratory control sample, matrix spike and matrix spike duplicate samples. The matrix spike recovery should fall within +/- 30 % and the matrix spike duplicates should compare within 20% RPD. The method blanks must meet the criteria in Section 9.4. The laboratory control sample should recover within 20% of the true value until in house limits are established.
- 13.3. Training Qualification:

The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

**14. POLLUTION PREVENTION**

- 14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

**15. WASTE MANAGEMENT**

- 15.1. Waste generated in the procedure must be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**

16.1. References

- 16.1.1. Test Methods for Evaluating Solid Waste , Physical/Chemical Methods, SW-846, 3rd Edition, Final Update II, Revision I, September 1994, Method 7471A (Mercury).
- 16.1.2. "Methods for the Chemical Analysis of Water and Wastes", EPA-600/4-79-020, U.S.EPA, August 1983, Method 245.5.
- 16.1.3. U.S.EPA Statement of Work for Inorganics Analysis, ILMO3.0 and IMLO4.0.
- 16.1.4. Corporate Quality Management Plan (QMP), current version.
- 16.2. STL Laboratory Quality Manual (LQM), current version. Associated SOPs and Policies, latest version

16.2.1. QA Policy, QA-003

16.2.2. Glassware Washing, NC-QA-0014

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018

16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021

16.2.5. Navy/Army SOP, NC-QA-0016

17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . )

17.1. Modifications/Interpretations from reference method.

17.1.1. Modifications from both 7471A and 245.5.

17.1.1.1. A potassium persulfate oxidation step has been included to facilitate the breakdown of organic mercurials which are not completely oxidized by potassium permanganate. Use of potassium persulfate in combination with the permanganate improves the recovery of mercury from organo-mercury compounds. The use of persulfate has been incorporated in several recent EPA mercury protocols.

17.1.1.2. The alternate run sequence presented in Section 11.2.11 is consistent with method requirements. An additional QC analysis (CRA) was added to accommodate the CLP protocol requirements.

17.1.2. Modifications from Method 7471A

17.1.2.1. Chapter 1 of SW846 specify the use of reagent water with a purity equivalent to ASTM Type II water. This SOP specifies the use of a Millipore DI system or equivalent to produce reagent water. This SOP requires that reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.

17.1.2.2. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the MDL. This SOP states that the method blank must not contain any analyte of interest at or above the reporting limit.

17.1.3. Modifications from 245.5

- 17.1.3.1. Method 245.5 Section 9.3 states concentrations should be reported as follows: Between 0.1 and 1 ug/g, to the nearest 0.01 ug; between 1 and 10 ug/g, to the nearest 0.1ug; above 10 ug/g, to the nearest ug. STL North Canton reports all Hg results under this SOP to two significant figures.

17.2. Modifications from previous SOP

- 17.2.1. Section 1.4 reporting limit changed from 0.1 mg/kg based on a 0.2 g to 0.3 mg/kg based on a 0.6 g sample aliquot
- 17.2.2. Section 9.3 added MS and MSDs as not counted in determination of preparation batches.

17.3. Facility Specific SOPs


Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none. Refer to the Appendices for any facility specific information required to support this SOP.

17.4. Documentation and Record Management


The following documentation comprises a complete CVAA raw data package:

- Raw data (direct instrument printout)
- Run log printout from instrument software where this option is available or manually generated run log. (A bench sheet may be substituted for the run log as long as it contains an accurate representation of the analytical sequence).
- Data review checklist - See Appendix B
- Standards Documentation (source, lot, date).
- Copy of digestion log.

Non-conformance summary (if applicable).

**Figure 1.** Solid Sample Preparation for Mercury - Water Bath Procedure

**Figure 2.** CVAA Mercury Analysis



877 936

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX A –TABLES

SOP No. CORP-MT-0007NC  
Revision No. 2.4  
Revision Date: 10/28/02  
Page 25 of 42

---

## APPENDIX A

### TABLES

**TABLE I. MERCURY REPORTING LIMITS, CALIBRATION STANDARD\*, QC  
STANDARD AND SPIKING LEVELS**

Soil RL (mg/kg)	0.1
Std 0 (mg/L)	0
Std 1 (mg/L)	0.0002
Std 2 (mg/L)	0.0005
Std 3 (mg/L)	0.001
Std 4 (mg/L)	0.005
Std 5 (mg/L) **	0.010
ICV (mg/L)	0.001 or 0.0025 ***
CCV/LCS (mg/L)	0.0025 or 0.005 ***
MS (mg/L)	0.001

- \* SOP specified calibration levels must be used unless prevented by the instrument configuration or client specific requirements. Deviations from specified calibration levels must be documented in the facility specific instrument operation SOP and must be approved by the facility technical manager and Quality Assurance Manager.
- \*\* Optional standard which may be used to extend the calibration range as allowed by the instrument configuration. If the instrument configuration prevents the use of 6 standards, the 2 ppb standard may be eliminated in favor of the 10 ppb standard.
- \*\*\* Concentration level dependent on high calibration standard used. CCV must be 50% of the high standard concentration and the ICV must be 20-25% of the high standard concentration.



**TABLE II. Summary Of Quality Control Requirements**

QC PARAMETER	FREQUENCY *	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
ICV	Beginning of every analytical run.	90-110 % recovery.	Terminate analysis; Correct the problem; Recalibrate or reprep with calibration curve (see Section 9.7).
ICB	Beginning of every analytical run, immediately following the ICV.	The result must be within +/- RL from zero.	Terminate analysis; Correct the problem; Recalibrate or reprep with calibration curve (see Section 9.7).
CCV	Every 10 samples and at the end of the run.	80 - 120 % recovery.	Terminate analysis; Correct the problem; Recalibrate and rerun all samples not bracketed by acceptable CCV or reprep with calibration curve (see Section 9.8).
CCB	Immediately following each CCV.	The result must be within +/- RL from zero.	Terminate analysis; Correct the problem; Recalibrate and rerun all samples not bracketed by acceptable CCB or reprep with calibration curve (see Section 9.8).
Method Blank	One per sample preparation batch of up to 20 samples.	<p>The result must be less than or equal to the RL.</p> <p>Sample results greater than 20x the blank concentration are acceptable.</p> <p>Samples for which the contaminant is &lt; RL do not require redigestion (See Section 9.4)</p>	<p>Redigest and reanalyze samples.</p> <p>Note exceptions under criteria section.</p> <p>See Section 9.4 for additional requirements.</p>

\*See Sections 11.2.10 and 11.2.11 for exact run sequence to be followed.

**TABLE II. Summary of Quality Control Requirements (Continued)**

QC PARAMETER	FREQUENCY	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
Laboratory Control Sample (LCS)	One per sample preparation batch of up to 20 samples.	Aqueous LCS must be within 80 - 120% recovery or in-house control limits.	Terminate analysis; Correct the problem; Redigest and reanalyze all samples associated with the LCS (see Section 9.5).
Matrix Spike	One per sample preparation batch of up to 20 samples.	75 - 125 % recovery or in-house control limits. If the MS/MSD is out for an analyte, it must be in control in the LCS.	In the absence of client specific requirements, flag the data; no flag required if the sample level is > 4x the spike added. (see Section 9.6)  For TCLP see Section 11.3.12
Matrix Spike Duplicate	See Matrix Spike	75 - 125 % recovery or in-house control limits; RPD ≤ 20%. (See MS)	See Corrective Action for Matrix Spike.

**877 940**

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX B - QUANTERRA HG DATA REVIEW CHECKLIST

SOP No. CORP-MT-0007NC

Revision No. 2.2

Revision Date: 02/06/01

Page 29 of 42

---

## **APPENDIX B**

### **EXAMPLE**

#### **STL NORTH CANTON Hg DATA REVIEW CHECKLIST**

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX B - QUANTERRA HG DATA REVIEW CHECKLIST

877 941

SOP No. CORP-MT-0007NC  
Revision No. 2.2  
Revision Date: 02/06/01  
Page 30 of 42

**Example**  
**STL North Canton Hg Data Review Checklist**

**Run/Project Information**

Run Date: \_\_\_\_\_ Analyst: \_\_\_\_\_ Instrument: \_\_\_\_\_  
Prep Batches Run: \_\_\_\_\_

Circle Methods used: 7470A / 245.1 : CORP-MT-0005 Rev 1      7471 / 245.5 : CORP-MT-0007 Rev 1  
CLP - AQ : CORP-MT-0006 Rev 0 CLP - SOL : CORP-MT-0008 Rev 0

**Review Items**

A. Calibration/Instrument Run QC	Yes	No	N/A	2ndLevel
1. Instrument calibrated per manufacturer's instructions and at SOP specified levels ?				
2. ICV/CCV analyzed at appropriate frequency and within control limits?				
3. ICB/CCB analyzed at appropriate frequency and within +/- RL or +/- CRDL (CLP)?				
4. CRA run (CLP only)?				
B. Sample Results				
1. Were samples with concentrations > the high calibration standard diluted and reanalyzed?				
2. All reported results bracketed by in control QC ?				
3. Sample analyses done within holding time?				
C. Preparation/Matrix QC				
1. LCS done per prep batch and within QC limits ?				
2. Method blank done per prep batch and < RL or CRDL (CLP) ?				
3. MS run at required frequency and within limits ?				
4. MSD or DU run at required frequency and RPD within SOP limits?				
D. Other				
1. Are all nonconformances documented appropriately ?				
2. Current IDL/MDL data on file?				
3. Calculations and Transcriptions checked for error ?				
4. All client/ project specific requirements met?				
5. Date of analysis verified as correct ?				

Analyst: \_\_\_\_\_ Date: \_\_\_\_\_  
Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

2nd Level Reviewer : \_\_\_\_\_ Date: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

877 942

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX C - MSA GUIDANCE

SOP No. CORP-MT-0007NC  
Revision No. 2.2  
Revision Date: 02/06/01  
Page 31 of 42

---

**APPENDIX C**  
**MSA GUIDANCE**


## APPENDIX C. MSA GUIDANCE

### Method of Standard Addition

Four equal volume aliquots of sample are measured and known amounts of standards are added to three aliquots. The fourth aliquot is the unknown and no standard is added to it. The concentration of standard added to the first aliquot should be 50% of the expected concentration. The concentration of standard added to the second aliquot should be 100% of the expected concentration and the concentration of standard added to the third aliquot should be 150% of the expected concentration. The volume of the unspiked and spiked aliquots should be the same (i.e., the volume of the spike added should be negligible in relation to the volume of sample).

To determine the concentration of analyte in the sample, the absorbance (or response) of each solution is determined and a linear regression performed. On the vertical axis the absorbance (or response) is plotted versus the concentrations of the standards on the horizontal axis using 0 as the concentration of the unspiked aliquot. An example plot is shown in Figure 1. When the resulting line is extrapolated back to zero absorbance, the point of interception of the horizontal axis is the concentration of the unknown. Calculate the correlation coefficient ( $r$ ) and the x-intercept (where  $y=0$ ) of the curve. The concentration in the digestate is equal to the negative x-intercept.

Figure 1



877 944

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX C - MSA GUIDANCE

SOP No. CORP-MT-0007NC  
Revision No. 2.2  
Revision Date: 02/06/01  
Page 33 of 42

---

- For the method of standard additions to be correctly applied, the following limitations must be taken into consideration.
- The plot of the sample and standards must be linear over the concentration range of concern. For best results, the slope of the curve should be similar to that of a plot of the aqueous standard curve.
- The effect of the interference should not vary as the ratio of the standard added to the sample matrix changes.

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
TROUBLESHOOTING GUIDE

SOP No. CORP-MT-0007NC

Revision No. 2.2

Revision Date: 02/06/01

Page 34 of 42

**APPENDIX D**  
**TROUBLESHOOTING GUIDE**



PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
TROUBLESHOOTING GUIDE

SOP No. CORP-MT-0007NC

Revision No. 2.2Revision Date: 02/06/01

Page 35 of 42

## APPENDIX D. TROUBLESHOOTING GUIDE

Problem	Possible Cause
Poor or No Absorbance or Sensitivity Check failed	Incorrect wavelength Dirty windows Window loose Etched or dirty optics Wrong lamp Bad lamp Not enough or no sample introduced Empty sample cup Incorrectly made standards Gas leak EDL power supply set on "Continuous"
Erratic Readings	Source lamp not aligned properly Lamp not prewarmed Injection tip partially clogged Contaminated reagents Contaminated glassware Drying tube saturated Bad lamp Injection tip hitting outside of tube Injection tip coated or not set properly Leak in sample tubing Power fluctuations Air bubbles in tubing
EDL Won't Light	Lamp cable not plugged in Lamp power set at 0 Lamp is dead Power supply fuse is blown Short in cord
Standards reading twice or half normal absorbance or concentration	Incorrect standard used Incorrect dilution performed Dirty cell
Background Correction Light Blinking	Background screen or attenuator faulty

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX E - CONTAMINATION CONTROL

---

SOP No. CORP-MT-0007NC

Revision No. 2.2

Revision Date: 02/06/01

Page 36 of 42

**APPENDIX E**  
**CONTAMINATION CONTROL GUIDELINES**

## APPENDIX E. CONTAMINATION CONTROL GUIDELINES

### **The following procedures are strongly recommended to prevent contamination:**

All work areas used to prepare standards and spikes should be cleaned before and after each use.

All glassware should be washed with detergent and tap water and rinsed with 1:1 nitric acid followed by deionized water.

Proper laboratory housekeeping is essential in the reduction of contamination in the metals laboratory. All work areas must be kept scrupulously clean.

Powdered or Latex Gloves must not be used in the metals laboratory since the powder contains silica and zinc, as well as other metallic analytes. Only vinyl or nitrile gloves should be used in the metals laboratory.

Glassware should be periodically checked for cracks and etches and discarded if found. Etched glassware can cause cross contamination of any metallic analytes.

Autosampler trays should be covered to reduce the possibility of contamination. Trace levels of elements being analyzed in the samples can be easily contaminated by dust particles in the laboratory.

### **The following are helpful hints in the identification of the source of contaminants:**

Reagents or standards can contain contaminants or be contaminated with the improper use of a pipette.

Improper cleaning of glassware can cause contamination.

Separate glassware if an unusually high sample is analyzed and soak with sulfuric acid prior to routine cleaning.

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX F - PREVENTIVE MAINTENANCE

---

877 949  
SOP No. CORP-MT-0007NC  
Revision No. 2.2  
Revision Date: 02/06/01  
Page 38 of 42

**APPENDIX F**  
**PREVENTIVE MAINTENANCE**

## APPENDIX F. PREVENTIVE MAINTENANCE

A maintenance log is used to record when maintenance is performed on instruments. When an instrument problem occurs indicate the date, time and instrument number, then identify the problem and corrective action in the maintenance log.

The following procedures are required to ensure that that the instrument is fully operational.

### Cold Vapor Atomic Absorption (Leeman PS 200) <sup>(1)</sup>

Daily	Semi-annually	Annually
Clean lens.	Check Hg lamp intensity.	Change Hg lamp.
Check aperture.		Check liquid/gas separator.
Check argon flow.		
Check tubing		
Check drain		
Replace drying tube.		

### Cold Vapor Atomic Absorption (PE 5000) <sup>(1)</sup>

Daily	Monthly
Clean aspirator by flushing with DI water.	Clean cell in aqua regia.
Check tubing and replace if needed.	Clean aspirator in aqua regia.
Clean windows with methanol.	
Change silica gel in drying tube	
Check argon gas supply.	
Adjust lamp	

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX G – INSTRUMENT SET-UP

---

SOP No. CORP-MT-0007NC

Revision No. 2.2

Revision Date: 02/06/01

Page 40 of 42

## APPENDIX G

### INSTRUMENT SET UP

## **Hg Analysis (Leeman PS200II)**

### **SYSTEM INITIALIZATION AND WARM UP**

1. F1 Menu
2. Instrument
  - a. Taskmaster
  - b. #4 Wake System Up   Enter

The warming up period takes approximately 10 minutes.

### **TO SET UP INSTRUMENT FOR ANALYSIS**

1. F1 Menu
2. Autosampler
  - A. Rack Entry
  - B. Edit (ex. Rack 1), Enter
  - C. Cup ID - Enter (clears sample #'s)
  - D. Extended ID- type in matrix of sample (water or solid), Enter.
  - E. Press Insert Key and move cursor with arrows to cup ID and begin typing labels.
  - F. F3 Print Screen
3. Press F2 Macro key and type in analyst's first name - Enter
  - A. Enter folder name - ex. HG0306, Enter. If folder does not exist, type Y - Enter.
  - B. Type in - "Rack 1", "Rack 2" etc. , Enter.

PREPARATION AND ANALYSIS OF MERCURY IN SOLID  
SAMPLES BY COLD VAPOR ATOMIC ABSORPTION, SW-846  
METHOD 7471A and MCAWW METHOD 245.5  
APPENDIX G – INSTRUMENT SET-UP

---

SOP No. CORP-MT-0007NC

Revision No. 2.2

Revision Date: 02/06/01

Page 42 of 42

D. Type in : FROM CUP    TO CUP

Ex.    =        1                    30

Do the same for position 2 if needed. If not needed, you must press Enter 3 times to begin analysis.



Controlled Copy

Copy No \_\_\_\_\_

Implementation Date: 1/22/04

SOP No. CORP-MT-0001NC

Revision No. 3.4Revision Date: 01/08/04Page 1 of 54**STL NORTH CANTON STANDARD OPERATING PROCEDURE**

**TITLE: INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION SPECTROSCOPY,  
SPECTROMETRIC METHOD FOR TRACE ELEMENT ANALYSES,  
SW-846 METHOD 6010B AND EPA METHOD 200.7**

(SUPERSEDES: REVISION 3.3, REVISION DATE (12/05/01))

Approved by:	<u><i>Patsy O'M</i></u>	<u>11/12/04</u>
		Date
Approved by:	<u><i>Karen J. Counts</i></u>	<u>1-12-2004</u>
	Technical Specialist	Date
Approved by:	<u><i>Beth Lambert</i></u>	<u>1-21-04</u>
	Quality Assurance Manager	Date
Approved by:	<u><i>Deborah Duran</i></u>	<u>1-19-2004</u>
	Environmental, Health and Safety Coordinator	Date
Approved by:	<u><i>Paul M. Cook</i></u>	<u>1/20/04</u>
	Laboratory Manager	Date
Approved by:	<u><i>Mark L. Bruce</i></u>	<u>1/22/04</u>
	Technical Director	Date

**Proprietary Information Statement:**

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

## TABLE OF CONTENTS

1. SCOPE AND APPLICATION .....	3
2. SUMMARY OF METHOD .....	3
3. DEFINITIONS .....	4
4. INTERFERENCES .....	4
5. SAFETY .....	5
6. EQUIPMENT AND SUPPLIES .....	7
7. REAGENTS AND STANDARDS .....	7
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE .....	8
9. QUALITY CONTROL .....	8
10. CALIBRATION AND STANDARDIZATION .....	16
11. PROCEDURE .....	17
12. DATA ANALYSIS AND CALCULATIONS .....	22
13. METHOD PERFORMANCE .....	23
14. POLLUTION PREVENTION .....	25
15. WASTE MANAGEMENT .....	25
16. REFERENCES .....	26
17. MISCELLANEOUS (TABLES, APPENDICES, ETC.) .....	27

## LIST OF APPENDICES:

APPENDIX A - TABLES .....	30
APPENDIX B - CROSS REFERENCE OF TERMS USED IN METHODS AND IN SOP .....	40
APPENDIX C - MSA GUIDANCE .....	42
APPENDIX D - TROUBLESHOOTING GUIDE .....	44
APPENDIX E - CONTAMINATION CONTROL GUIDELINES .....	46
APPENDIX F - PREVENTATIVE MAINTENANCE .....	48
APPENDIX G - ICP OPERATING INSTRUCTIONS .....	50

## 1. SCOPE AND APPLICATION

- 1.1. This procedure describes the analysis of trace elements including metals in solution by Inductively Coupled Plasma -Atomic Emission Spectroscopy (ICP-AES) using SW-846 Method 6010B and EPA Method 200.7. Table I of Appendix A lists the elements appropriate for analysis by Methods 6010B and 200.7. Additional elements may be analyzed under Methods 6010B and 200.7 provided that the method performance criteria presented in Section 13.0 are met.
- 1.2. ICP analysis provides for the determination of metal concentrations over several orders of magnitude. Detection limits, sensitivity and optimum concentration ranges of the metals will vary with the matrices and instrumentation used.
- 1.3. Method 6010B is applicable to the determination of dissolved, suspended, total recoverable and total elements in ground water, aqueous samples, soils, sludges, wastes, sediments, biological, and TCLP, EP and other leachates/extracts. All matrices require digestion prior to analysis with the exception of analyses for dissolved metals in filtered and acidified aqueous samples. Although digestion is not specifically required by the method, some clients and regulators may require digestion of **dissolved samples** and this must be clarified and documented before project initiation. Silver concentrations must be below 2.0 mg/L in aqueous samples and 100 mg/kg in solid matrix samples. Precipitation may occur in samples where silver concentrations exceed these levels and lead to the generation of erroneous data.
- 1.4. Method 200.7 is applicable to the determination of dissolved, suspended, total recoverable, and total elements in water, waste water, and solid wastes. All matrices require digestion prior to analysis with the exception of analyses for dissolved metals in filtered and acidified aqueous samples if the criteria in Section 11.1 are met. Silver concentrations must be below 0.1 mg/L in aqueous samples.
- 1.5. State-specific requirements may take precedence over this SOP for drinking water sample analyses.
- 1.6. The applicable LIMS method codes are QO (6010B), QM (6010B Trace), AS (200.7), JI (200.7 Trace). The Trivalent Chromium method code is GU.

## 2. SUMMARY OF METHOD

- 2.1. This method describes a technique for the determination of multi elements in solution using sequential or simultaneous optical systems and axial or radial viewing of the plasma. The basis of the method is the measurement of atomic emission by an optical spectroscopic technique. Samples are nebulized and the aerosol that is produced is transported to the plasma torch where excitation occurs. Characteristic atomic-line

emission spectra are produced by a radio frequency inductively coupled plasma (ICP). The spectra are dispersed by a grating spectrometer and the intensities of the emission lines are monitored by photomultiplier tubes. The photocurrents from the photomultiplier tubes are processed and controlled by a computer system. A background correction technique is required to compensate for variable background contribution to the determination of trace elements. Background must be measured adjacent to analyte lines during analysis. The position selected for the background intensity measurement, on either or both sides of the analytical line, will be determined by the complexity of the spectrum adjacent to the analyte line. The position used must be free of spectral interferences and reflect the same change in background intensity as occurs at the analyte wavelength measured. Background correction is not required in cases of line broadening where a background correction measurement would actually degrade the analytical result. The possibility of additional interferences should also be recognized and appropriate actions taken. Alternatively, multivariate calibration methods may be chosen for which point selection for background correction is superfluous since whole spectral regions are processed.

- 2.2. Refer to CORP-IP-0002NC, Acid Digestion of Soils, SW846 Method 3050B and CORP-IP-0003NC, Acid Digestion of Aqueous Samples by SW846 and MCAWW 200 Series Methods for details on sample preparation methods.

### 3. DEFINITIONS

- 3.1. Dissolved Metals: Those elements which pass through a 0.45 um membrane. (Sample is acidified after filtration).
- 3.2. Suspended Metals: Those elements which are retained by a 0.45 um membrane.
- 3.3. Total Metals: The concentration determined on an unfiltered sample following vigorous digestion.
- 3.4. Total Recoverable Metals: The concentration determined on an unfiltered sample following treatment with hot, dilute mineral acid.

### 4. INTERFERENCES

- 4.1. Spectral, physical and chemical interference effects may contribute to inaccuracies in the determinations of trace elements by ICP. Spectral interferences are caused by:
- Overlap of a spectral line from another element.
  - Unresolved overlap of molecular band spectra.

- Background contribution from continuous or recombination phenomena.
  - Stray light from the line emission of high concentration elements.
- 4.1.1. A background correction technique is required to compensate for variable background contribution to the determination of trace elements. Background correction is not required in cases where a background corrective measurement would actually degrade the analytical result.
- 4.1.2. Inter-element correction factors (IECs) are necessary to compensate for spectral overlap. Inter-element interferences occur when elements in the sample emit radiation at wavelengths so close to that of the analyte that they contribute significant intensity to the analyte channel. If such conditions exist, the intensity contributed by the matrix elements will cause an excessively high (or sometimes low) concentration to be reported for the analyte. Inter-element corrections IECs must be applied to the analyte to remove the effects of these unwanted emissions.
- 4.1.3. Physical interferences are generally considered to be effects associated with sample transport, nebulization and conversion within the plasma. These interferences may result in differences between instrument responses for the sample and the calibration standards. Physical interferences may occur in the transfer of solution to the nebulizer (e.g., viscosity effects), at the point of aerosol formation and transport to the plasma (e.g., surface tension) or during excitation and ionization processes within the plasma itself. Changes in viscosity and surface tension can cause significant inaccuracies, especially in samples containing high dissolved solids or high acid concentrations. If physical interferences are present, dilution of the sample, use of a peristaltic pump, mass flow controller, use of an internal standard and/or use of a high solids nebulizer can reduce the effect.
- 4.1.4. Chemical interferences are characterized by molecular compound formation, ionization effects and solute vaporization effects. Normally these effects are not significant with the ICP technique but if observed can be minimized by buffering the sample, matrix matching or standard addition procedures.

## 5. SAFETY

- 5.1. Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual and this document.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 6 of 54

that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.

- 5.3. The following is a list of the materials used in this method, which have a serious or significant hazard rating. **NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table.** A complete list of materials used in the method can be found in the reagents and materials section. Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

Material (1)	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
Nitric Acid	Corrosive Oxidizer Poison	2 ppm-TWA 4-ppm STEL	Nitric acid is extremely hazardous; it is corrosive, reactive, an oxidizer, and a poison. Inhalation of vapors can cause breathing difficulties and lead to pneumonia and pulmonary edema, which may be fatal. Other symptoms may include coughing, choking, and irritation of the nose, throat, and respiratory tract. Can cause redness, pain, and severe skin burns. Concentrated solutions cause deep ulcers and stain skin a yellow or yellow-brown color. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
Hydrochloric Acid	Corrosive Poison	5 ppm-Ceiling	Inhalation of vapors can cause coughing, choking, inflammation of the nose, throat, and upper respiratory tract, and in severe cases, pulmonary edema, circulatory failure, and death. Can cause redness, pain, and severe skin burns. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
1 – Always add acid to water to prevent violent reactions			
2 – Exposure limit refers to the OSHA regulatory exposure limit.			

- 5.3.1. The plasma emits strong UV light and is harmful to vision. **NOTE: AVOID looking directly at the plasma.**
- 5.3.2. The RF generator produces strong radio frequency waves, most of which are unshielded. People with pacemakers should not go near the instrument while in operation.
- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples should be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Metals digestates can be processed outside of a fume hood. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Inductively Coupled Plasma Atomic Emission Spectrometer equipped with autosampler and background correction.
- 6.2. Radio Frequency Generator.
- 6.3. Argon gas supply, welding grade or equivalent.
- 6.4. Coolflow or appropriate water cooling device.
- 6.5. Peristaltic Pump.
- 6.6. Calibrated automatic pipettes or Class A glass volumetric pipettes.
- 6.7. Class A volumetric flasks.
- 6.8. Autosampler tubes.

## 7. REAGENTS AND STANDARDS

- 7.1. Intermediate standards are purchased as custom STL North Canton multi-element mixes or as single-element solutions. All standards must be stored in FEP fluorocarbon or unused polyethylene or polypropylene bottles. Intermediate standard solutions must be replaced prior to the expiration date provided by the manufacturer. If no expiration date is provided,

the intermediate solutions may be used for up to one year and must be replaced sooner if verification from an independent source indicates a problem. Expiration dates can be extended provided that the acceptance criteria described in laboratory-specific SOPs are met.

- 7.2. Working calibration and calibration verification solutions may be used for up to 3 months and must be replaced sooner if verification from an independent source indicates a problem. Standards should be prepared in a matrix of 5% hydrochloric and 5% nitric acids. Refer to Tables III, IV, IVA, V and VI (Appendix A) for details regarding the working standard concentrations for calibration, calibration verification, interference correction and spiking solutions.
- 7.3. Concentrated nitric acid ( $\text{HNO}_3$ ), trace metal grade or better.
- 7.4. Concentrated hydrochloric acid ( $\text{HCl}$ ), trace metal grade or better.
- 7.5. Reagent water must be produced by a Millipore DI system or equivalent. Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

- 8.1. Sample holding times for metals are six months from time of collection to the time of analysis.
- 8.2. Aqueous samples are preserved with nitric acid to a pH of  $<2$  and may be stored in either plastic or glass. If boron or silica are to be determined, plastic containers are preferred. Refrigeration is not required. Preservation must be verified prior to analysis.
- 8.3. Soil samples do not require preservation but must be stored at  $4^\circ\text{C} \pm 2^\circ$  until the time of preparation.

## 9. QUALITY CONTROL

Table VII (Appendix A) provides a summary of quality control requirements including type, frequency, acceptance criteria and corrective action.

- 9.1. Initial Demonstration of Capability
  - 9.1.1. Prior to analysis of any analyte using either Method 200.7 or Method 6010B, the following requirements must be met.
  - 9.1.2. Instrument Detection Limit (IDL) - The IDL for each analyte must be determined for each analyte wavelength used for each instrument. The IDL must be determined



annually. If the instrument is adjusted in anyway that may affect the IDL, the IDL for that instrument must be redetermined. The IDL shall be determined by multiplying by 3, the standard deviation obtained from the analysis of a blank solution, with seven consecutive measurements. Each measurement must be performed as though it were a separate analytical sample (i.e., each measurement must be followed by a rinse and/or any other procedure performed between the analysis of separate samples). The result of the IDL determination must be below the STL North Canton reporting limit. The CLP IDL procedure can be used for this method.

- 9.1.3. Method Detection Limit (MDL) - An MDL must be determined for each analyte prior to the analysis of any client samples. Refer to STL North Canton SOP NC-QA-0021 for details on MDL analysis and criteria.
- 9.1.4. Linear Range Verification (LR) - The linear range must be verified every 6 months on at least an annual basis for each analyte wavelength used on each instrument. The linear range is the concentration above which results cannot be reported without dilution of the sample. The standards used to verify the linear range limit must be analyzed during a routine analytical run and must read within 5% of the expected value.

For the **initial** determination of the upper limit of the linear dynamic range (LDR) for each wavelength, determine the signal responses from a minimum of three to five different concentration standards across the estimated range. One standard should be near the upper limit of the estimated range. The concentration measured at the LDR must be no more than 10% less than the expected level extrapolated from lower standards. If the instrument is adjusted in any way that may affect the LRs, new dynamic ranges must be determined. The LR data must be documented and kept on file.

- 9.1.5. Background Correction Points - To determine the appropriate location for off-line background correction when establishing methods, the user must scan the area on either side adjacent to the wavelength and record the apparent emission intensity from all other method analytes. This spectral information must be documented and kept on file. The location selected for background correction must be either free of off-line interelement spectral interference or a computer routine must be used for automatic correction on all determinations. Tests to determine spectral interference must be done using analyte concentrations that will adequately describe the interference. Background correction points must be set prior to determining IECs. Refer to the ICP instrument manual for specific procedures to be used in setting background correction points.
- 9.1.6. Inter-element Corrections (IECs) - ICP interelement correction factors must be determined prior to the analysis of samples and every six months thereafter. If the instrument is adjusted in any way that may affect the IECs, the IECs must be redetermined. When initially determining IECs for an instrument, wavelength scans must be performed to ensure that solutions in use are free from contaminants. If an IEC varies significantly from the previously determined IEC then the possibility of contamination should be investigated. The purity of the IEC check solution can be verified by using a standard from a second source or an alternate method (i.e., GFAA or ICP-MS). Published wavelength tables (e.g. MIT tables, Inductively Coupled Plasma-Atomic Spectroscopy: Prominent Lines) can also be consulted to

evaluate the validity of the IECs. Refer to the instrument manufacturer's recommendations for specific procedures to be used in setting IECs. An IEC must be established to compensate for any interelement interference which results in a false analyte signal greater than  $\pm$  the RL as defined in Tables I, IA or II. For elements with a reporting limit 10 ug/L or less, the signal should be  $\pm$  2X the RL. To determine IECs, run a single element standard at the established linear range. To calculate an IEC, divide the observed concentration of the analyte by the actual concentration of the "interfering element."

Note: Trace ICP IECs are more sensitive to small changes in the plasma and instrument setup conditions. Adjustments in the IECs will be required on a more frequent basis for the Trace as reflected by the ICSA response. Additional spectral interference is present from easily ionizable elements such as potassium and sodium in axial viewing instruments.

9.1.7. Rinse Time Determination - Rinse times must be determined upon initial set-up of an ICP instrument. To determine the appropriate rinse time for a particular ICP system, the linear range verification standard (see 9.1.4) should be aspirated as a regular sample followed by the analysis of a series of rinse blanks. The length of time required to reduce the analyte signals to  $<$  RL will define the rinse time for a particular ICP system. For some analytes it may be impractical to set the rinse time based on the linear range standard result (i.e., analyte not typically detected in environmental samples at that level and an excessive rinse time would be required at the linear range level). Until the required rinse time is established, the method recommends a rinse period of at least 60 seconds between samples and standards. If a memory effect is suspected, the sample must be reanalyzed after a rinse period of sufficient length. Rinse time studies can be conducted at additional concentration levels. These additional studies must be documented and kept on file, if a concentration other than the linear range level is used to set the rinse time. The concentration levels used to establish the rinse time must be taken into consideration when reviewing the data.

9.2. Method Blank (MB) - One method blank must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting limit (exception: common laboratory contaminants, see below) or at or above 5% of the measured concentration of that analyte in associated samples, whichever is higher (sample result must be a minimum of 20x higher than the blank contamination level).

- If the analyte is a common laboratory contaminant (copper, iron, lead (Trace only) or zinc) the data may be reported with qualifiers if the concentration of the analyte in the method blank is less than two times the RL. Such action must be addressed in the project narrative.
  - Repreparation and reanalysis of all samples associated with an unacceptable method blank is required when reportable concentrations are determined in the samples (see exception noted above).
  - If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. Such action must be addressed in the project narrative.
  - If the above criteria are not met and reanalysis is not possible, then the sample data must be qualified. This anomaly must be addressed in the project narrative.
  - For dissolved metals samples which have not been digested, a CCB result is reported as the method blank. The CCB run immediately prior to the start of the dissolved sample analyses must be used for this purpose. No more than 20 samples can be associated with one CCB.
- 9.3. Laboratory Control Sample (LCS) - One aqueous LCS must be processed with each preparation batch. The LCS must contain all analytes of interest and must be carried through the entire analytical procedure. Aqueous LCS spike levels are provided in Table III (Appendix A). The LCS is used to monitor the accuracy of the analytical process. On-going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines.
- If any analyte is outside established control limits the system is out of control and corrective action must occur. Unless in-house control limits are established, a control limit of 80 - 120% recovery must be applied.
  - In the event that an MS/MSD analysis is not possible a Laboratory Control Sample Duplicate (LCSD) must be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.
  - In the instance where the LCS recovery is greater than 120% and the sample results are < RL, the data may be reported with qualifiers. Such action must be addressed in the report narrative.
  - Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.

- For dissolved metals samples which have not been digested, a CCV result is reported as the LCS. The CCV run immediately prior to the start of the dissolved sample analyses must be used for this purpose. No more than 20 samples can be associated with one CCV.

9.4. Matrix Spike/Matrix Spike Duplicate (MS/MSD) - One MS/MSD pair must be processed for each preparation batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSDs. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Samples identified as field blanks cannot be used for MS/MSD analysis. Spiking levels are provided in Tables III and VI (Appendix A).

- If any analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. For both methods 200.7 and 6010B, control limits of 75-125% recovery and 20% RPD or historical acceptance criteria must be applied to the MS/MSD. If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted. If the recovery of the LCS is outside limits corrective action must be taken. Corrective action will include repreparation and reanalysis of the batch. MS/MSD results which fall outside the control limits must be addressed in the narrative.
- If the native analyte concentration in the MS/MSD exceeds 4x the spike level for that analyte, the recovery data are reported as NC, MSB (i.e., not calculated). Two other narrative notes for metals analyses: Matrix spike/spike duplicate spike recovery/recoveries was/were outside the acceptance limits of some analytes. The acceptable LCS analysis data indicated that the analytical system was operating within control and this condition is most likely due to matrix interference. See the Matrix Spike Report for the affected analytes which will be flagged with N. Matrix spike/spike duplicate relative percent difference (RPD) exceeded the acceptance limits for some analytes. The imprecision may be attributed to sample heterogeneity. See the Matrix Spike Report for the affected analytes, which will be flagged with \*.
- If an MS/MSD is not possible due to limited sample volume then a laboratory control sample duplicate (LCSD) should be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.

- For dissolved metals samples which have not been digested, a MS/MSD must be performed per batch of up to 20 samples by spiking two aliquots of the sample at the levels specified in Table III (Appendix A).
- 9.5. Dilution test – A dilution test is performed to determine whether significant physical or chemical interferences exist due to the sample matrix. One sample per preparation batch must be processed as a dilution test. The test is performed by running a sample at a 5x (1:4) dilution. Samples identified as field blanks cannot be used for dilution tests. The results of the diluted sample, after correction for dilution, should agree within 10% of the original sample determination when the original sample concentration is greater than 50x the IDL. If the results are not within 10%, the possibility of chemical or physical interference exists and the data is flagged.
- 9.6. Initial Calibration Verification (ICV/ICB) - Calibration accuracy is verified by analyzing a second source standard (ICV). For analyses conducted under Method 200.7, the ICV result must fall within 5% of the true value for that solution with relative standard deviation <3% from replicate (minimum of two) exposures. For Method 6010B, the ICV must fall within 10% of the true value for that solution with relative standard deviation <5% from replicate (minimum of two) exposures. An ICB is analyzed immediately following the ICV to monitor low level accuracy and system cleanliness. The ICB result must fall within +/- the RL from zero. If either the ICV or ICB fail to meet criteria, the analysis should be terminated, the problem corrected, the instrument recalibrated and the calibration reverified. (See Section 11.8 or 11.11 for required run sequence).
- 9.7. Continuing Calibration Verification (CCV/CCB) - Calibration accuracy is monitored throughout the analytical run through the analysis of a known standard after every 10 samples and at the end of the sample run. The CCV is be a mid-range standard made from a dilution of the calibration standard. The CCV for both methods must fall within 10% of the true value for that solution with relative standard deviation <5% from replicate (minimum of two) exposures. A CCB is analyzed immediately following each CCV. (See Section 11.8 or 11.11 for required run sequence.) The CCB result must fall within +/- RL from zero. If the blank is less than 1/10 the concentration of the action level of interest, and no sample is within 10% of the action limit, reanalysis and recalibration are not required before continuation of the run. Sample results may only be reported when bracketed by valid CCV/CCB pairs. If a mid-run CCV or CCB fails, all affected samples will be re-analyzed with valid CCV/CCB pairs. (Refer to Section 11.13 for an illustration of the appropriate rerun sequence).
- 9.8. Interference Check Analysis (ICSA/ICSAB) - The validity of the interelement correction factors is demonstrated through the successful analysis of interference check solutions. The ICSA contains only interfering elements, the ICSAB contains analytes and interferents. Refer to Table V (Appendix A) for the details of ICSA and ICSAB composition. Custom

STL North Canton multielement ICS solutions must be used. All analytes should be spiked into the ICSAB solution, therefore, if a non-routine analyte is required then it should be manually spiked into the ICSAB using a certified ultra high purity single element solution or custom lab-specific mix. If the ICP will display overcorrection as a negative number then the non-routine elements can be controlled from the ICSA as described in section 9.8.3. Elements known to be interferents on a required analyte must be included in the ICP run when that analyte is determined. Aluminum, iron, calcium and magnesium must always be included in all ICP runs.

- 9.8.1. The ICSA and ICSAB solutions must be run at the beginning of the run. (See Section 11.10 or 11.13 for required run sequence.)
- 9.8.2. The ICSAB results for the interferents must fall within 80 - 120% of the true value. If any ICSAB interferent result fails criteria, the analysis should be terminated, the problem corrected, the instrument recalibrated and the samples rerun.
- 9.8.3. ICSA results for the non-interfering elements with reporting limits  $\leq 10$   $\mu\text{g/L}$  must fall within the STL North Canton guidelines of  $\pm 2x$  RL from zero. ICSA results for the non-interfering elements with  $\text{RLs} > 10$   $\mu\text{g/L}$  must fall within the STL North Canton guidelines of  $\pm 1x$  RL from zero. If the ICSA results for the non-interfering elements do not fall within  $\pm 2x$  RL ( $\text{RL} \leq 10$ ) or  $\pm 1x$  RL ( $\text{RL} > 10$ ) from zero the field sample data must be evaluated as follows:
  - If the non-interfering element concentration in the ICSA is the result of contamination versus a spectral interference, and this reason is documented, the field sample data can be accepted.
  - If the affected element was not required then the sample data can be accepted.
  - If the interfering elements are not present in the field sample at a concentration which would result in a false positive or negative result greater than  $\pm 2x$  RL from zero then the field sample data can be accepted.
  - If the interfering element is present in the field sample at a level which would result in a false analyte signal greater than  $\pm 2x$  RL from zero, the data can be accepted only if the concentration of the affected analyte in the field sample is more than  $10x$  the analyte signal in the ICSA.
  - If the data does not meet the above conditions then the IECs must be re-evaluated and corrected if necessary and the affected samples reanalyzed or the sample results manually corrected through application of the new IEC to the

raw results. If the results are recalculated manually the calculations must be clearly documented on the raw data.

- 9.9. CRI - To verify linearity near the RL for ICP analysis, a CRI standard is run at the beginning of each sample analysis run. Additionally, some projects may require CRI analysis at the end of the run. (See Section 11.10 or 11.13 for required run sequence.) Evaluate associated samples based upon advisory limits of  $\pm 50\%$  of true value.

Note: The custom STL North Canton CRI mix contains most analytes at a level near the standard lab reporting limit.

Note: For certified Ohio drinking water analysis, the CRI concentration must be at or below the reporting limit and must recover at  $\pm 30\%$  of true value.

- 9.10. Method of Standard Addition (MSA) -This technique involves adding known amounts of standard to one or more aliquots of the processed sample solution. This technique compensates for a sample interferent that may enhance or depress the analyte signal, thus producing a different slope from that of the calibration standards. It will not correct for additive interferences which cause a baseline shift. Refer to Section 11.17 for additional information on when MSA is required as well as Appendix C for specific MSA requirements.

## 10. CALIBRATION AND STANDARDIZATION

- 10.1. Set up the instrument with the operating parameters recommended by the manufacturer. Allow the instrument to become thermally stable before beginning calibration (approximately 30 minutes of warm-up is required). Refer to the instructions in Appendix G.
- 10.2. Profile and calibrate the instrument according to the instrument manufacturer's recommended procedures. Flush the system with the calibration blank between each standard or as the manufacturer recommends. The calibration curve must consist of a minimum of a blank and a standard. Refer to Appendix G for detailed set up and operation protocols.
- 10.3. Calibration must be performed daily and each time the instrument is set up. Instrument runs may be continued over periods exceeding 24 hours as long as all calibration verification (CCV) and interference check QC criteria are met. The instrument standardization date and time must be included in the raw data.
- 10.4. Refer to Section 9.0 for calibration verification procedures, acceptance criteria and corresponding corrective actions.



## 11. PROCEDURE

- 11.1. For 200.7 analyses, dissolved (preserved) samples must be digested unless it can be documented that the sample meets all of the following criteria:
- A. Visibly transparent with a turbidity measurement of 1 NTU or less.
  - B. Is of one liquid phase and free of particulate or suspended matter following acidification.
  - C. Is NOT being analyzed for silver.
- 11.2. A minimum of two exposures for each standard, field sample and QC sample is required. The average of the exposures is reported. For Trace ICP analyses, the results of the sum channel must be used for reporting.
- 11.3. Prior to calibration and between each sample/standard the system is rinsed with the calibration blank solution. The minimum rinse time between analytical samples is 60 seconds unless following the protocol outlined in 9.1.5 it can be demonstrated that a shorter rinse time may be used.
- 11.4. The use of automated QC checks through the instrument software is highly recommended for all calibration verification samples (ICV,CCV), blanks (ICB,CCB,PB), interference checks (ICSA,ICSAB) and field samples (linear range) to improve the data review process.
- 11.5. To facilitate the early identification of QC failures and samples requiring rerun it is strongly recommended that sample data be reviewed periodically throughout the run.
- 11.6. To facilitate the data review and reporting processes it is strongly recommended that all necessary dilutions be performed before closing out the instrument run.
- 11.7. The use of an internal standard is **required** on the Trace ICP unless the calibration and QC standards are matrix matched to each digestion procedure used as follows:

Preparation Method	% HNO <sub>3</sub>	% HCl
CLP Aqueous	1	5
CLP Soil	5	2.5
SW846 3050	10	10
SW846 3005	2	5

SW846 3010	6	5
------------	---	---

The following procedural guidelines must be followed when using an internal standard:

- 11.7.1. Typically used internal standards are: yttrium. (Note: Any element can be used that is not typically found in environmental samples at a high rate of occurrence.)
- 11.7.2. The internal standard (IS) must be added to every sample and standard at the same concentration. It is recommended that the IS be added to each analytical sample automatically through use of a third pump channel and mixing coil. Internal standards should be added to blanks, samples and standards in a like manner, so that dilution effects resulting from the addition may be disregarded.
- 11.7.3. The concentration of the internal standard should be sufficiently high to obtain good precision in the measurement of the IS analyte used for data correction and to minimize the possibility of correction errors if the IS analyte is naturally present in the sample.
- 11.7.4. The internal standard raw intensity counts must be printed on the raw data.
- 11.7.5. The analyst must monitor the response of the internal standard throughout the sample analysis run. This information is used to detect potential problems and identify possible background contributions from the sample (i.e., natural occurrence of IS analyte). The instrument automatically adjusts sample results based on comparison of the internal standard intensity in the sample to the internal standard intensity at calibration.
  - 11.7.5.1. If the internal standard counts fall within  $\pm 30\%$  of the counts observed in the ICB then the data is acceptable.
  - 11.7.5.2. If the internal standard counts in the field samples are more than  $\pm 30\%$  higher than the expected level, the field samples must then be diluted and reanalyzed.
- 11.8. The following analytical sequence must be used for Methods 6010B and 200.7:  
Instrument Calibration  
ICV  
ICB  
CRI

ICSA

ICSAB

7 samples

CCV

CCB

10 samples

CCV

CCB

Repeat sequence of up to 10 samples between CCV/CCB pairs as required to complete run

CRI (The CRI counts as a sample analysis.)

CCV

CCB

Refer to Quality Control Section 9.0 and Table VII (Appendix A) for Method 6010B and 200.7 quality control criteria.

- 11.9. Additional quality control analyses are necessary for analysis under the Contract Laboratory Program (CLP). If these are included then CLP, 6010 and 200.7 samples can be included in the same sequence. Refer to CORP-MT-0002NC for details.
- 11.10. Full method required QC must be available for each wavelength used in determining reported analyte results.
- 11.11. The following run sequence provides an illustration of a mid-run CCV or CCB failure and the appropriate corrective action run sequence as described in Section 9.7:

Original Run: Instrument Calibration

ICV

ICB

CRI

ICSA

ICSAB

7 samples

CCV1

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7

877 973  
SOP No. CORP-MT-0001NC  
Revision No. 3.4  
Revision Date: 01/08/04  
Page 20 of 54

---

CCB1

10 samples

CCV2

CCB2

10 samples \*\*

CCV3 \*                      \* Failure occurs at CCV3/CCB3

CCB3 \*                      \*\*Samples requiring rerun for affected analytes

10 samples \*\*

CCV4

CCB4

10 samples

CCV5

CCB5

Reanalysis:    Recalibrate

ICV

ICB

CRI

ICSA

ICSAB

CCV2

CCB2

10 samples

CCV3

CCB3

10 samples

CCV4

CCB4

Notes: Samples between CCV4 and CCV5 do not require reanalysis as they were bracketed by compliant QC samples.

See CORP-MT-0002NC for the appropriate reanalysis sequence if CLP requirements must also be met.

- 11.14 The instrument may be reprofiled between CCV/CCB pairs to correct for environment induced drift.
- 11.15 Guidelines are provided in the appendices on procedures to minimize contamination of samples and standards, preventive maintenance and troubleshooting.
- 11.16 All measurements must fall within the defined linear range where spectral interference correction factors are valid. Dilute and reanalyze all samples for required analytes that exceed the linear range or use an alternate wavelength for which QC data are established. If an interelement correction exists for an analyte which exceeds the linear range, the IEC may be inaccurately applied. Therefore, even if an overrange analyte may not be required to be reported for a sample, if that analyte is a interferent for any requested analyte in that sample, the sample must be diluted. Acid strength must be maintained in the dilution of samples.
- 11.17 For TCLP samples, full four-point MSA will be required if all of the following conditions are met:
- 1) recovery of the analyte in the matrix spike is not at least 50%,
  - 2) the concentration of the analyte does not exceed the regulatory level, and,
  - 3) the concentration of the analyte is within 20% of the regulatory level.
- The reporting and regulatory limits for TCLP analyses as well as matrix spike levels are detailed in Table VI (Appendix A). Appendix C provides guidance on performing MSA analyses.
- 11.18 Any variation in procedure shall be completely documented using instrument run logs, maintenance logs, report narratives, a Nonconformance Memo, or an anomaly report and

is approved by a Supervisor/Group Leader and QA Manager. If contractually required, the client shall be notified by the Project Manager.

11.19 Nonconformance documentation shall be filed in the project file.

11.20 Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

11.21 Analytical Documentation

11.21.1 Record all analytical information in the analytical logbook/logsheets which may be in an electronic format, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.21.2 All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. Logbooks are reviewed by the supervisor or designee.

11.21.3 Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

11.21.4 Sample results and associated QC are entered into the LIMs after final technical review.

## 12 DATA ANALYSIS AND CALCULATIONS

12.1. ICV percent recoveries are calculated according to the equation:

$$\%R = 100 \left( \frac{\text{Found(ICV)}}{\text{True(ICV)}} \right)$$

12.2. CCV percent recoveries are calculated according to the equation:

$$\%R = 100 \left( \frac{\text{Found(CCV)}}{\text{True(CCV)}} \right)$$

12.3. Matrix Spike Recoveries are calculated according to the following equation:

$$\%R = 100 \left( \frac{SSR - SR}{SA} \right)$$

Where:

SSR = Spike Sample Result

SR = Sample Result

SA = Spike Added

- 12.4. The relative percent difference (RPD) of matrix spike/matrix spike duplicates are calculated according to the following equations:

$$RPD = 100 \left[ \frac{|MSD - MS|}{\left( \frac{MSD + MS}{2} \right)} \right]$$

Where:

MS = determined spiked sample concentration

MSD = determined matrix spike duplicate concentration

- 12.5. The final concentration for a digested aqueous sample is calculated as follows:

$$mg / L = \frac{C \times V1 \times D}{V2}$$

Where:

C = Concentration (mg/L) from instrument readout

D = Instrument dilution factor

V1 = Final volume in liters after sample preparation

V2 = Initial volume of sample digested in liters

- 12.6. The final concentration determined in digested solid samples when reported on a dry weight basis is calculated as follows:

$$mg / Kg, dry weight = \frac{C \times V \times D}{W \times S}$$

Where:

C = Concentration (mg/L) from instrument readout

D = Instrument dilution factor

V = Final volume in liters after sample preparation

W = Weight in Kg of wet sample digested

S = Percent solids/100

Note: A Percent Solids determination must be performed on a separate aliquot when dry weight concentrations are to be reported. If the results are to be reported on wet weight basis the "S" factor should be omitted from the above equation.

- 12.7. The LCS percent recovery is calculated according to the following equation:

$$\%R = 100 \left( \frac{\text{Found}(LCS)}{\text{True}(LCS)} \right)$$

- 12.8. The dilution test percent difference for each component is calculated as follows:

$$\%Difference = \frac{|I - S|}{I} \times 100$$

Where:

I = Sample result (Instrument reading)

S = Dilution test result (Instrument reading  $\times$  5)

- 12.9. Appropriate factors must be applied to sample values if dilutions are performed.
- 12.10. Sample results should be reported with up to three significant figures in accordance with the STL North Canton significant figure policy.

#### 12.11. Trivalent Chromium

- 12.11.1. Trivalent chromium ( $\text{Cr}^{+3}$ ) is the result obtained by subtracting the hexavalent chromium ( $\text{Cr}^{+6}$ ) results for a sample from the total chromium result from that sample. The total chromium result is determined using the procedures in this SOP. The hexavalent chromium result is determined using the procedures in STL North Canton SOP NC-WC-0024.

#### 12.11.2. Reporting Limits

- 12.11.2.1. The STL North Canton water reporting limit for trivalent chromium is 0.02 mg/l.

- 12.11.2.2. The STL North Canton solid reporting limit for trivalent chromium is 2.0 mg/kg, wet weight.

- 12.11.3. Calculations:  $\text{Cr}^{+3} = \text{Cr, total} - \text{Cr}^{+6}$

### 13. METHOD PERFORMANCE



- 13.1. Each laboratory must have initial demonstration of performance data on file for each analyte of interest as described in Section 9.0.
- 13.2. Refer to Tables I, IA & II in Appendix A for the list of Method 6010B and 200.7 analytes as well as additional analytes that may be analyzed using this SOP.
- 13.3. Method performance is determined by the analysis of MS and MSD samples as well as method blanks and laboratory control samples. The MS or MSD recovery should fall within +/- 25 % and the MS/MSD should compare within 20% RPD or within the laboratory's historical acceptance limits. These criteria apply to analyte concentrations greater than or equal to 10xIDL. Method blanks must meet the criteria specified in Section 9.2. The laboratory control samples should recover within 20% of the true value or within the laboratory's historical acceptance limits.
- 13.4. Training Qualification:
  - 13.4.1. The group/team leader or the supervisor has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

#### 14. POLLUTION PREVENTION

- 14.2. This method does not contain any specific modifications that serve to minimize or prevent pollution.

#### 15. WASTE MANAGEMENT

- 15.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.
- 15.2. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."
- 15.3. Waste Streams Produced by this Method
  - 15.3.1 The following waste streams are produced when this method is carried out:
  - 15.3.2 Acid waste consisting of sample and rinse solution - Any sample waste generated must be collected and disposed of in the acid waste drum located in the metals lab.

15.3.3 Standards should be purchased and prepared in volumes consistent with laboratory use to minimize the volume of expired standards to be disposed.

## 16. REFERENCES

### 16.1. References

- 16.1.1. 40 CFR Part 136, Appendix B, 7-5-95, Determination of Method Detection Limits.
- 16.1.2. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW-846, 3rd Edition, Final Update III, Revision 2, December 1996. Method 6010B.
- 16.1.3. Determination of Metals and Trace Elements in Water and Wastes by Inductively Coupled Plasma-Atomic Emission Spectrometry, Revision 4.4, May 1994. Method 200.7.
- 16.1.4. Inductively Coupled Plasma – Atomic Emission Spectrometric Method for Trace Element Analysis of water and wastes Method 200.7, 40 CFR – Chapter I – Part 136 – Appendix C. Electronic version dated September 30, 2002.
- 16.1.5. CORP-MT-0002NC, Inductively Coupled Plasma-Atomic Emission Spectroscopy, Method 200.7 & CLP-M, SOW ILMO3.0 and ILMO4.0).
- 16.1.6. Corporate Quality Management Plan (QMP), current version
- 16.1.7. STL Laboratory Quality Manual (LQM), current version.

### 16.2. Associated SOPs and Policies, latest version

- 16.2.1. QA-003, STL North Canton QC Program.
- 16.2.2. Glassware Washing, NC-QA-0014
- 16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018
- 16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021
- 16.2.5. Navy/Army SOP, NC-QA-0016
- 16.2.6. QA-004, Rounding and Significant Figures.

- 16.2.7. NC-WC-0024, Hexavalent Chromium (Colorimetric)
- 16.2.8. CORP-IP-0002NC, Acid Digestion of Soils, SW846 Method 3050B
- 16.2.9. CORP-IP-0003NC, Acid Digestion of Aqueous Samples by SW846 and MCAWW 200 Series Methods

## 17. MISCELLANEOUS (TABLES, APPENDICES, ETC.)

### 17.1. Modifications/Interpretations from reference method

#### 17.1.1. Modifications/interpretations from both Methods 6010B and 200.7.

- 17.1.1.1. STL North Canton laboratories use mixed calibration standard solutions purchased from approved vendors instead of using individual mixes prepared in house as recommended by the subject methods.
- 17.1.1.2. Methods 200.7 and 6010B state that if the correction routine is operating properly, the determined apparent analyte(s) concentration from analysis of each interference solution should fall within a specific concentration range around the calibration blank. In determining IECs, because of lack of definition clarification for "concentration range around the calibration blank," STL North Canton has adopted the procedure in EPA CLP ILMO4.0.
- 17.1.1.3. Section 8.5 of Method 6010B and Section 9.5 of Method 200.7 recommend that whenever a new or unusual matrix is encountered, a series of tests be performed prior to reporting concentration data for that analyte. The dilution test helps determine if a chemical or physical interference exists. Because STL North Canton laboratories receive no prior information from clients regarding when to expect a new or unusual matrix, STL North Canton may select to perform a dilution test on one sample in each prep batch. According to the method, the post digestion spike (PDS) determines any potential matrix interferences. At STL North Canton, matrix interference is determined by evaluating data for the LCS and MS/MSD. STL North Canton REQUIRES documented, clear guidance when a new or unusual matrix will be received for a project and a request to perform the dilution test or PDS on a client-identified sample.

17.1.2. Modifications from Method 200.7.

- 17.1.2.1. Method 200.7 defines the IDL as the concentration equivalent to a signal, due to the analyte, which is equal to three times the standard deviation of a series of ten replicate measurements of the calibration blank signal at the same wavelength. STL North Canton labs utilize the IDL definition as defined in Section 9.1 of this SOP.
- 17.1.2.2. The calibration blank is prepared in an acid matrix of 5% HNO<sub>3</sub>/5% HCl instead of the specified 2% HNO<sub>3</sub>/10% HCl matrix as the former matrix provides for improved performance relative to the wide variety of digestate acid matrices which result from the various EPA preparation protocols applied.
- 17.1.2.3. Method section 9.3.4 specifies that "Analysis of the IPC (ICSA/AB) solution immediately following calibration must verify that the instrument is within  $\pm 5\%$  of calibration with a relative standard deviation  $<3\%$  from replicate integrations  $\geq 4$ ." STL North Canton uses a minimum of two exposures.
- 17.1.2.4. The 40 CFR version of Method 200.7 requires the instrument check standard to agree within  $\pm 5\%$  of expected values. Also, the 40 CFR version requires the interference check sample to be analyzed at the beginning, end, and at periodic intervals throughout the sample run. At STL North Canton, the instrument check standard equals the CCV, which must agree within  $\pm 10\%$  of expected values, and the ICSA standards are analyzed only at the beginning of a sample run. STL's procedures are in line with the Rev. 4.4, May 1994 version of Method 200.7.
- 17.1.2.5. Section 7.12 of 200.7 indicates that the QCS (ICV) should be prepared at a concentration near 1 ppm. The ICV specified in this SOP accommodates the 1 ppm criteria for the majority of analytes. For the remaining analytes, this SOP specifies ICV concentrations which are appropriate to the range of calibration. The intent of the ICV, verification of calibration standard accuracy, is independent of the ICV concentration used.
- 17.1.2.6. The ICS criteria applied by this SOP differ from those stated in the method. Method 200.7 section 10.4 states that results should fall within the established control limits of 3 times the standard deviation of the

calibration blank for that analyte. The control limits listed in this SOP are those applicable to the EPA designed solution.

- 17.1.2.7. Method 200.7 section 9.3.4 states the CCB should be less than the IDL, but > the lower 3-sigma control limit of the calibration blank. The intent of this requirement is to ensure that the calibration is not drifting at the low end. STL North Canton has adopted an absolute control limit of +/- RL from zero for calibration blank criteria. SOP section 9.7 provides the detailed corrective action criteria that must be followed.

17.1.3. Modifications from Method 6010B.

- 17.1.3.1. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the MDL. This SOP states that the method blank must not contain any analyte of interest at or above the reporting limit. Common lab contaminants are allowed up to two times the reporting limit in the blank following consultation with the client.
- 17.1.3.2. Method 6010B section 8.6.1.3 states that the results of the calibration blank are to agree within 3x the IDL. If not, repeat the analysis two or more times and average the results. If the average is not within three standard deviations of the background mean, terminate the analysis, correct the problem, recalibrate, and reanalyze the previous 10 samples. The intent of this requirement is to ensure that the calibration is not drifting at the low end. STL North Canton has adopted an absolute control limit of +/- RL from zero for calibration blank criteria. See SOP Section 9.7 for a detailed description of the required corrective action procedures.

**APPENDIX A**

**TABLES**

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP No. CORP-MT-0001NC  
Revision No. 3.4  
Revision Date: 01/08/04  
Page 31 of 54

**TABLE I. Method 200.7 and 6010B Target Analyte List**

ELEMENT	Symbol	CAS #	6010B analyte	200.7 analyte	Reporting Limit (ug/L) Water	Reporting Limit (mg/kg) Soil
Aluminum	Al	7429-90-5	X	X	200	20
Antimony	Sb	7440-36-0	X	X	60	6
Arsenic	As	7440-38-2	X	X	300	30
Barium	Ba	7440-39-3	X	X	200	20
Beryllium	Be	7440-41-7	X	X	5.0	0.5
Boron	B	7440-42-8	X	X	200	20
Cadmium	Cd	7440-43-9	X	X	5.0	0.5
Calcium	Ca	7440-70-2	X	X	5000	500
Chromium	Cr	7440-47-3	X	X	10	1
Cobalt	Co	7440-48-4	X	X	50	5
Copper	Cu	7440-50-8	X	X	25	2.5
Iron	Fe	7439-89-6	X	X	100	10
Lead	Pb	7439-92-1	X	X	100	10
Magnesium	Mg	7439-95-4	X	X	5000	500
Manganese	Mn	7439-96-5	X	X	15	1.5
Molybdenum	Mo	7439-98-7	X	X	40	4
Nickel	Ni	7440-02-0	X	X	40	4
Potassium	K	7440-09-7	X	X	5000	500
Selenium	Se	7782-49-2	X	X	250	25
Silver	Ag	7440-22-4	X	X	10	1
Sodium	Na	7440-23-5	X	X	5000	500
Thallium	Tl	7440-28-0	X	X	2000	200
Vanadium	V	7440-62-2	X	X	50	5
Zinc	Zn	7440-66-6	X	X	20	2

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

877 985  
SOP No. CORP-MT-0001NC  
Revision No. 3.4  
Revision Date: 01/08/04  
Page 32 of 54

**TABLE IA. Method 200.7 and 6010B Trace ICP Target Analyte List**

ELEMENT	Symbol	CAS #	Reporting Limit (ug/L) Water	Reporting Limit (mg/kg) Soil
Arsenic	As	7440-38-2	10	1.0
Lead	Pb	7439-92-1	3.0	0.3
Selenium	Se	7782-49-2	5.0	0.5
Thallium	Tl	7440-28-0	10	1.0
Antimony	Sb	7440-36-0	10	1.0
Cadmium	Cd	7440-43-9	2.0	0.2
Silver	Ag	7440-22-4	5.0	0.5
Chromium	Cr	7440-47-3	5.0	0.5

**TABLE II. Non-Routine Analyte List**

ELEMENT	Symbol	CAS #	Reporting Limit (ug/L) Water	Reporting Limit (mg/kg) Soil
Tin	Sn	7440-31-5	100	10
Titanium	Ti	7440-32-6	50	5



INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 33 of 54

**TABLE III. Matrix Spike and Aqueous Laboratory Control Sample Levels**

ELEMENT	LCS Level (ug/L)	Matrix Spike Level (ug/L)
Aluminum	2000	2000
Antimony	500	500
Arsenic	2000	2000
Barium	2000	2000
Beryllium	50	50
Cadmium	50	50
Calcium	50000	50000
Chromium	200	200
Cobalt	500	500
Copper	250	250
Iron	1000	1000
Lead	500	500
Magnesium	50000	50000
Manganese	500	500
Molybdenum	1000	1000
Nickel	500	500
Potassium	50000	50000
Selenium	2000	2000
Silver	50	50
Sodium	50000	50000
Thallium	2000	2000
Vanadium	500	500
Zinc	500	500
Boron	1000	1000
Tin	2000	2000
Titanium	1000	1000

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 34 of 54

**TABLE IV. ICP Calibration and Calibration Verification Standards**

Element	Calibration Level	RL (ug/L)	ICV (ug/L)	CCV (ug/L)
Aluminum	100000	200	25000	50000
Antimony	10000	60	1000	5000
Arsenic	10000	300	1000	5000
Barium	10000	200	1000	5000
Beryllium	10000	5	1000	5000
Cadmium	10000	5	1000	5000
Calcium	100000	5000	25000	50000
Chromium	10000	10	1000	5000
Cobalt	10000	50	1000	5000
Copper	10000	25	1000	5000
Iron	100000	100	25000	50000
Lead	10000	100	1000	5000
Magnesium	100000	5000	25000	50000
Manganese	10000	15	1000	5000
Molybdenum	10000	40	1000	5000
Nickel	10000	40	1000	5000
Potassium	100000	5000	25000	50000
Selenium	10000	250	1000	5000
Silver	2000	10	500	1000
Sodium	100000	5000	25000	50000
Thallium	20000	2000	5000	10000
Vanadium	10000	50	1000	5000
Zinc	10000	20	1000	5000
Boron	10000	200	1000	5000
Tin	10000	100	1000	5000
Titanium	10000	50	1000	5000

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 35 of 54

**TABLE IVA. Trace ICP Calibration and Calibration Verification Standards**

Element	Calibration Level	RL (ug/L)	ICV (ug/L)	CCV (ug/L)
Aluminum	50000	200	12500	25000
Antimony	1000	10	250	500
Arsenic	1000	10	250	500
Barium	4000	10	1000	2000
Beryllium	4000	5	1000	2000
Cadmium	1000	2	250	500
Calcium	100000	5000	25000	50000
Chromium	4000	5	1000	2000
Cobalt	4000	50	1000	2000
Copper	4000	25	1000	2000
Iron	50000	100	12500	25000
Lead	1000	3	250	500
Magnesium	100000	5000	25000	50000
Manganese	4000	15	1000	2000
Molybdenum	4000	40	1000	2000
Nickel	4000	40	1000	2000
Potassium	100000	5000	25000	50000
Selenium	1000	5	250	500
Silver	2000	5	500	1000
Sodium	100000	5000	25000	50000
Thallium	2000	10	500	1000
Vanadium	4000	50	1000	2000
Zinc	4000	20	1000	2000

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP<sup>1</sup>No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 36 of 54

**TABLE V. Interference Check Sample Concentrations \***

Element	ICSA (ug/L)	ICSAB (ug/L)
Aluminum	500000	500000
Antimony	-	1000
Arsenic	-	1000
Barium	-	500
Beryllium	-	500
Cadmium	-	1000
Calcium	500000	500000
Chromium	-	500
Cobalt	-	500
Copper	-	500
Iron	200000	200000
Lead	-	1000
Magnesium	500000	500000
Manganese	-	500
Molybdenum	-	1000
Nickel	-	1000
Potassium	-	10000
Selenium	-	1000
Silver	-	1000
Sodium	-	10000
Thallium	-	10000**
Vanadium	-	500
Zinc	-	1000
Tin	-	1000
Boron		1000
Titanium		1000

\* Custom STL North Canton solutions contain analytes common to all STL North Canton facilities. Non-routine elements not listed above should be spiked into the ICSAB at 1000 ug/L.

\*\* Thallium level for Trace ICP should be at 1000 ug/L.

877 990

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 37 of 54

---

**TABLE VI. TCLP Reporting Limits, Regulatory Limits and Matrix Spike Levels**

ELEMENT	Reporting Level (ug/L)	Regulatory Limit (ug/L)	Spike Level (ug/L)
Arsenic	500	5000	5000
Barium	10000	100000	50000
Cadmium	100	1000	1000
Chromium	500	5000	5000
Lead	500	5000	5000
Selenium	250	1000	1000
Silver	500	5000	1000

**TABLE VII. Summary Of Quality Control Requirements**

QC PARAMETER	FREQUENCY *	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
Two-point Initial Calibration	Beginning of every analytical run, every 24 hours, whenever instrument is modified, or CCV criterion is not met		Terminate analysis; Correct the problem; Prepare new standards; Recalibrate following system performance.
ICV	Beginning of every analytical run.	Method 200.7: 95 - 105 % recovery. RSD dupl. exp < 3%  Method 6010B: 90 - 110 % recovery. RSD dupl. exp < 5%	Terminate analysis; Correct the problem; Recalibrate.
ICB	Beginning of every analytical run, immediately following the ICV.	The result must be within +/- RL from zero.	Terminate analysis; Correct the problem; Recalibrate.
CCV	Every 10 samples and at the end of the run.	Method 200.7 & 6010B:  90 - 110 % recovery. RSD dupl. exp < 5%	Terminate analysis; Correct the problem; Recalibrate and rerun all samples not bracketed by acceptable CCV.
CCB	Immediately following each CCV.	The result must be within +/- RL from zero.	Terminate analysis; Correct the problem; Recalibrate and rerun all samples not bracketed by acceptable CCB.
CRI	Beginning of every run	50-150% recovery (advisory)	Evaluate associated samples.
ICSA	Beginning of every run	See Section 9.8.3	See Section 9.8.3
ICSAB	Immediately following each ICSA.	Results must be within 80 - 120% recovery.	See Section 9.8.2.

\* See Sections 11.10 and 11.13 for exact run sequence to be followed.

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX A - TABLES**

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 39 of 54

**TABLE VII. Summary of Quality Control Requirements (Continued)**

QC PARAMETER	FREQUENCY	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
Dilution Test	One per prep batch.	For samples > 50x IDL, dilutions must agree within 10%.	Narrate the possibility of physical or chemical interference per client request.
Method Blank	One per sample preparation batch of up to 20 samples.	<p>The result must be less than or equal to the RL.</p> <p>Common lab contaminants may be accepted up to 2x the RL (See 9.2).</p> <p>Sample results greater than 20x the blank concentration are acceptable.</p> <p>Samples for which the contaminant is &lt; RL may not require redigestion or reanalysis (see Section 9.2).</p>	<p>Redigest and reanalyze samples.</p> <p>Note exceptions under criteria section.</p> <p>See Section 9.2 for additional requirements.</p>
Laboratory Control Sample (LCS)	One per sample preparation batch of up to 20 samples.	<p>Aqueous LCS must be within 80 - 120% recovery or in-house control limits.</p> <p>Samples for which the contaminant is &lt; RL and the LCS results are &gt; 120% may not require redigestion or reanalysis (see Section 9.3)</p>	<p>Terminate analysis;</p> <p>Correct the problem;</p> <p>Redigest and reanalyze all samples associated with the LCS.</p>

**TABLE VII. Summary of Quality Control Requirements (Continued)**

<b>QC PARAMETER</b>	<b>FREQUENCY</b>	<b>ACCEPTANCE CRITERIA</b>	<b>CORRECTIVE ACTION</b>
Matrix Spike	One per sample preparation batch of up to 20 samples.	75 - 125 % recovery or in- house control limits. For TCLP See Section 11.17.	In the absence of client specific requirements, flag the data; no flag required if the sample level is > 4x the spike added. For TCLP see Section 11.17.
Matrix Spike Duplicate	See Matrix Spike	75 - 125 % recovery; RPD ≤ 20% .	See Corrective Action for Matrix Spike.



**877 994**

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX B- CROSS REFERENCE OF TERMS USED**

---

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 41 of 54

## **APPENDIX B**

**CROSS REFERENCE OF TERMS USED IN METHODS 6010B, 200.7, AND BY STL NORTH  
CANTON**

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX B- CROSS REFERENCE OF TERMS USED**

877 995  
SOP No. CORP-MT-0001NC  
Revision No. 3.4  
Revision Date: 01/08/04  
Page 42 of 54

**CROSS REFERENCE OF TERMS COMMONLY USED IN  
METHODS EPA 200.7, SW6010B, AND STL NORTH CANTON INC. SOP**

EPA 200.7	SW6010B	STL North Canton Inc. SOP
Calibration blank (CB)	Calibration blank	Initial and continuing calibration blanks (ICB/CCB)
Dilution test	Dilution test	Dilution Test
Instrument detection limit (IDL)	Instrument detection limit (IDL)	Instrument detection limit (IDL)
Instrument performance check (IPC)	Continuing calibration verification (CCV)	Continuing calibration verification (CCV)
Internal standard	Internal standard	Internal standard (IS)
Laboratory duplicates	n/a	n/a
Laboratory fortified blank (LFB)	n/a	Laboratory control sample (LCS)
Laboratory fortified sample matrix (LFM)	Matrix spike and matrix spike duplicate (MS/MSD)	Matrix spike and matrix spike duplicate (MS/MSD)
Laboratory reagent blank (LRB)	Method blank	Method or Prep blank (MB)
Linear dynamic range (LDR)	Linear dynamic range (LDR)	Linear dynamic range (LDR)
Method detection limit (MDL)	Method detection limit (MDL)	Method detection limit (MDL)
Quality control sample (QCS)	Check standard or Initial calibration verification (ICV)	Initial calibration verification (ICV)
Spectral interference check solution (SIC)	Interference check solution (ICS)	Interference check solution (ICSA/ICSAB)

**877 996**

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX C- MSA GUIDANCE**

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 43 of 54

---

## **APPENDIX C**

### **MSA GUIDANCE**

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX C- MSA GUIDANCE**

SOP No. CORP-MT-0001NC

Revision No. 3.4Revision Date: 01/08/04Page 44 of 54


## Appendix C. MSA Guidance

### Method of Standard Addition

Four equal volume aliquots of sample are measured and known amounts of standards are added to three aliquots. The fourth aliquot is the unknown and no standard is added to it. The concentration of standard added to the first aliquot should be 50% of the expected concentration. The concentration of standard added to the second aliquot should be 100% of the expected concentration and the concentration of standard added to the third aliquot should be 150% of the expected concentration. The volume of the unspiked and spiked standard should be the same.

In order to determine the concentration of analyte in the sample, the analytical value of each solution is determined and a plot or linear regression performed. On the vertical axis the analytical value is plotted versus the concentrations of the standards on the horizontal axis. An example plot is shown in Figure 1. When the resulting line is extrapolated back to zero absorbance, the point of interception of the horizontal

axis is the  
concentration of the  
unknown.



Conc. of Sample	Addn 0 No Addn	Addn 1 Addn of 50% of Expected Amount	Addn 2 Addn of 100% of Expected Amount	Addn 3 Addn of 150% of Expected Amount
--------------------	-------------------	--	---	---

For the method of standard additions to be correctly applied, the following limitations must be taken into consideration:

- The plot of the sample and standards must be linear over the concentration range of concern. For best results, the slope of the curve should be similar to that of a plot of the aqueous standard curve.

877 998

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX C- MSA GUIDANCE**

---

SOP No. CORP-MT-0001NC

Revision No. 3.4

Revision Date: 01/08/04

Page 45 of 54

- The effect of the interference should not vary as the ratio of the standard added to the sample matrix changes.

**APPENDIX D**  
**TROUBLESHOOTING GUIDE**

8771000

INDUCTIVELY COUPLED PLASMA-ATOMIC EMISSION  
SPECTROSCOPY, SPECTROMETRIC METHOD FOR TRACE  
ELEMENT ANALYSIS, METHOD 6010B AND METHOD 200.7  
**APPENDIX D - TROUBLESHOOTING GUIDE**

SOP No. CORP-MT-0001NC  
Revision No. 3.4  
Revision Date: 01/08/04  
Page 47 of 54

**APPENDIX D. TROUBLESHOOTING GUIDE**

Problem	Possible Cause/ Solution
High Blanks	Increase rinse time Clean or replace tip Clean or replace torch Clean or replace sample tubing Clean or replace nebulizer Clean or replace mixing chamber Lower Torch
Instrument Drift	RF not cooling properly Vacuum level is too low Replace torch (Crack) Clean or replace nebulizer (blockage) Check room temperature (changing) Replace pump tubing Room humidity too high Clean torch tip (salt buildup) Check for argon leaks Adjust sample carrier gas Reprofile Horizontal Mirror Replace PA tube
Erratic Readings, Flickering Torch or High RSD	Check for argon leaks Adjust sample carrier gas Replace tubing (clogged) Check drainage(back pressure changing) Increase uptake time (too short) Increase flush time (too short) Clean nebulizer, torch or spray chamber Increase sample volume introduced Check that autosampler tubes are full Sample or dilution of sample not mixed Increase integration time (too short) Realign torch Reduce amount of tubing connectors
Cu/Mn Ratio Outside Limits or Low Sensitivity	Plasma conditions changed Clean nebulizer, torch or spray chamber Replace tubing (clogged) Realign torch Check IECs
Standards reading twice normal absorbance or concentration	Incorrect standard used Incorrect dilution performed

**APPENDIX E**  
**CONTAMINATION CONTROL GUIDELINES**



## **APPENDIX E. CONTAMINATION CONTROL GUIDELINES**

### **The following procedures are strongly recommended to prevent contamination:**

All work areas used to prepare standards and spikes should be cleaned before and after each use.

All glassware should be washed with detergent and tap water and rinsed with 1:1 nitric acid followed by deionized water.

Proper laboratory housekeeping is essential in the reduction of contamination in the metals laboratory. All work areas must be kept scrupulously clean.

Powdered or Latex Gloves must not be used in the metals laboratory since the powder contains silica and zinc as well as other metallic analytes. Only vinyl or nitrile gloves should be used in the metals laboratory.

Glassware should be periodically checked for cracks and etches and discarded if found. Etched glassware can cause cross contamination of any metallic analytes.

Autosampler trays should be covered to reduce the possibility of contamination. Trace levels of elements being analyzed in the samples can be easily contaminated by dust particles in the laboratory.

### **The following are helpful hints in the identification of the source of contaminants:**

Yellow pipet tips and volumetric caps can sometimes contain cadmium.

Some sample cups have been found to contain lead.

The markings on glass beakers have been found to contain lead. If acid baths are in use for glassware cleaning, they should be periodically checked for contaminants since contaminant concentrations will increase over time.

New glassware especially beakers can be a source of silica and boron.

Reagents or standards can contain contaminants or be contaminated with the improper use of a pipette.

Improper cleaning of glassware can cause contamination.

Latex gloves contain over 500 ppb of zinc.

**APPENDIX F**  
**PREVENTIVE MAINTENANCE**

## **APPENDIX F. PREVENTIVE MAINTENANCE**

A maintenance log is used to record when maintenance is performed on instruments. When an instrument problem occurs indicate the date, time and instrument number, then identify the problem and corrective action in the maintenance log.

**The following procedures are required to ensure that that the instrument is fully operational.**

<b>Daily</b>	Change sample pump tubing and pump windings Check argon gas supply level Check rinse solution and fill if needed Check waste containers and empty if needed Check sample capillary tubing is clean and in good condition Check droplet size to verify nebulizer is not clogged. Check sample flow for cross flow nebulizer Check Cu/Mn ratio-should be 30% of value at date that IECs were performed Check pressure for vacuum systems
<b>As Needed</b>	Clean plasma torch assembly to remove accumulated deposits Clean nebulizer and drain chamber; keep free-flowing to maintain optimum performance Replace peristaltic pump tubing, sample capillary tubing and autosampler sipper probe
<b>Weekly</b>	Apply silicon spray on autosampler tracks Check water level in coolflow
<b>Monthly</b>	Clean air filters on back of power unit to remove dust Check D mirror for air instruments
<b>Bi-yearly</b>	Change oil for vacuum systems Replace coolant water filter (may require more or less frequently depending on quality of cooling water)

**APPENDIX G**  
**ICP OPERATING INSTRUCTIONS**

## ICP Analysis (TJA 61E)

(Example)

### 1. SETUP

- a. Plasma Control Panel (enter)
- b. (F1)-Startup
- c. (F9)-Continue
- d. (F2)-Levels
  1. Change auxiliary gas to low – use space bar to toggle
  2. Change nebulizer gas flow to 0.5 L/min.
  3. Change pump rate to 130
  4. Esc
  5. Allow instrument to warm up approximately 30 minutes.

### 2. DEVELOPMENT

- a. Methods (enter)
- b. Enter method name
- c. (F3)-Method Info.
- d. Change file name
- e. (F9)-Done
- f. (F9)-Done/Keep

## 18. OPERATION

- a. Analysis (enter)
- b. (F5)-Profile
  1. (F3)-Automatic
  2. (F1)-Run
  3. If peak position is greater than +/- .05 units from the center peak position, you must adjust the profile. If it is within +/- .05 units, press (F9)-done.
  4. To adjust select (F1)-CalcSS and enter current vernier position. (enter)
  5. Adjust to new vernier position(F9)-done
  6. Rerun profile until peak position is +/- .05 units.
  7. (F9)-Done
- c. Autosampler (F9)
  1. Enter method name (enter)
  2. Enter autosampler table name (enter)
  3. (F1)-Run

8771008

Controlled Copy

Copy No. \_\_\_\_\_

Implementation Date: 4/25/03

SOP No. CORP-IP-0002NC

Revision No. 2.4

Revision Date: 02/19/03

Page 1 of 26

**STL NORTH CANTON STANDARD OPERATING PROCEDURE**

**TITLE: ACID DIGESTION OF SOILS, SW846 METHOD 3050B**

**(SUPERSEDES: REVISION 2.3 DATE 01/18/02)**

Prepared by:	<u>Trish McNair</u>	<u>4.25.03</u>	Date
Reviewed by:	<u>Patrick Smith</u> Technology Specialist	<u>3/26/03</u>	Date
Reviewed by:	<u>Betty Lambert</u> Quality Assurance Manager	<u>3/4/03</u>	Date
Reviewed by:	<u>Deborah Buddel</u> Environmental Health and Safety	<u>3/31/03</u>	Date
Reviewed by:	<u>Paul M. Ford</u> Laboratory Director	<u>4/1/03</u>	Date
Reviewed by:	<u>Mark Bruce</u> <del>Corporate Technology</del> Technical Director	<u>4/11/03</u>	Date

**Proprietary Information Statement:**

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED

**TABLE OF CONTENTS**

1. SCOPE AND APPLICATION.....	3
2. SUMMARY OF METHOD.....	3
3. DEFINITIONS .....	3
4. INTERFERENCES .....	3
5. SAFETY .....	4
6. EQUIPMENT AND SUPPLIES.....	5
7. REAGENTS AND STANDARDS.....	6
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE.....	8
9. QUALITY CONTROL.....	8
10. CALIBRATION AND STANDARDIZATION .....	10
11. PROCEDURE.....	10
12. DATA ANALYSIS AND CALCULATIONS.....	13
13. METHOD PERFORMANCE.. .....	13
14. POLLUTION PREVENTION.....	14
15. WASTE MANAGEMENT. ....	14
16. REFERENCES.....	14
17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . ).....	15

**LIST OF APPENDICES:**

APPENDIX A - TABLES.....	19
APPENDIX B - CONTAMINATION CONTROL GUIDELINES .....	24



## 1. SCOPE AND APPLICATION

- 1.1. This procedure describes the preparation of soil samples for the analysis of certain metals by Graphite Furnace Atomic Absorption (GFAA), Flame Atomic Absorption (FLAA), Inductively Coupled Plasma Atomic Emission Spectroscopy (ICP), and Inductively Coupled Plasma-Mass Spectrometry (ICP/MS) as specified in SW846 Method 3050B.
- 1.2. Samples prepared by the protocols detailed in this SOP may be analyzed by ICP, ICP/MS, FLAA or GFAA for the elements listed in Table I (Appendix A). Other elements and matrices may be analyzed following digestion by these protocols provided that the method performance criteria specified in Section 13.0 of this SOP are met.
- 1.3. This method is not a total digestion, but will dissolve almost all metals that could become "environmentally available". By design, metals bound in silicate structures are not dissolved by this procedure, as they are not usually mobile in the environment. This SOP can be applied to metals in solids, sludges, wastes and sediments.
- 1.4. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary.

## 2. SUMMARY OF METHOD

A representative 1 gram (wet weight) portion of sample is digested in nitric acid and hydrogen peroxide. The digestate is refluxed with hydrochloric acid for ICP or FLAA analysis. The digestates are then filtered and diluted to 100 mL/100 g.

## 3. DEFINITIONS

- 3.1. Refer to the glossary in the Laboratory Quality Manual (LQM), latest version.
- 3.2. Total Metals: The concentration determined on an unfiltered sample following digestion. Note that this method is designed to determine the total *environmentally available* metals.

## 4. INTERFERENCES

- 4.1. There are numerous routes by which samples may become contaminated. Potential sources of trace metals contamination include: metallic or metal-containing labware (e.g., talc gloves which contain high levels of zinc), containers, impure reagents, dirty glassware, improper sample transfers, dirty work areas, atmospheric inputs such as dirt and dust, etc. Be aware of potential sources of contamination and take appropriate measures to minimize or avoid them.

- 4.2. The entire work area, including the bench top and fume hood, should be thoroughly cleaned on a routine schedule in order to minimize the potential for environmental contamination.
- 4.3. Boron and silica from the glassware will grow into the sample solution during and following sample processing. For critical low level determinations of boron and silica, only quartz and/or plastic labware should be used.
- 4.4. Physical interference effects may contribute to inaccuracies in the determinations of trace elements. Oils, solvents and other matrices may not be digested using these methods if they are not soluble with acids. If physical interferences are present, they should be documented.
- 4.5. Visual interferences or anomalies (such as foaming, emulsions, precipitates, etc.) must be documented.
- 4.6. Allowing samples to boil or go dry during digestion may result in the loss of volatile metals. If this occurs the sample must be reprepared. Antimony is easily lost by volatilization from hydrochloric media.
- 4.7. Specific analytical interferences are discussed in each of the determinative methods.

## 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL North Canton associates.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:

5.3.1. The following materials are known to be **corrosive**:

Hydrochloric acid and nitric acid.

5.3.2. The following materials are known to be **oxidizing agents**:

Nitric acid and hydrogen peroxide.

5.3.3. All heating of samples must be carried out in a fume hood.

5.4. The acidification of samples containing reactive materials may result in the release of toxic gases, such as cyanides or sulfides. Acidification of samples should be done in a fume hood. The analyst should also be aware of the potential for a vigorous reaction.

5.5. Exposure to chemicals must be maintained **as low as reasonably achievable**. Therefore, unless they are known to be non-hazardous, all samples should be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.

5.6. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit or under other means of mechanical ventilation.

5.7. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.

5.8. Always carry bulk concentrated acid bottles in appropriate impact proof containers.

5.9. Acid/peroxide spills must be neutralized immediately, flushed with water and cleaned up using appropriate spill kits.

5.10. Discard chipped or broken beakers to prevent injury. Chipped glassware may be fire polished as an alternative to disposal.

5.11. Any and all accidents and spills must be reported to the lab supervisor or EH&S coordinator.

## 6. EQUIPMENT AND SUPPLIES

6.1. Hot plate, digestion block, steam bath or other heating source capable of maintaining a temperature of 90-99°C.

- 6.2. Calibrated thermometer that covers a temperature range of 0-200°C.
- 6.3. Griffin beakers of assorted sizes or equivalent.
- 6.4. Vapor recovery device (Watch glasses, ribbed or other device).
- 6.5. Whatman No. 41 filter paper or equivalent.
- 6.6. Funnels or equivalent filtration apparatus.
- 6.7. Centrifugation equipment (if desired method of removing particulates is centrifugation).
- 6.8. Graduated cylinder or equivalent capable of measuring 100 mL within 3% accuracy.
- 6.9. Analytical balance capable of accurately weighing to the nearest 0.01 grams.
- 6.10. Repipetors or suitable reagent dispensers.
- 6.11. Calibrated automatic pipettes with corresponding pipette tips or Class A glass volumetric pipettes.
- 6.12. Class A volumetric flasks.
- 6.13. pH indicator strips (pH range 0 - 6).
- 6.14. Plastic bottles.

## **7. REAGENTS AND STANDARDS**

- 7.1. Reagent water must be produced by a Millipore DI system or equivalent. Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks as defined in the determinative SOPs.
- 7.2. Laboratory Control Sample (LCS) and matrix spike (MS) solutions are purchased as custom STL North Canton solutions. All standards must be stored in FEP fluorocarbon or previously unused polyethylene or polypropylene bottles. Stock standard solutions must be replaced prior to the expiration date provided by the manufacturer. If no expiration date is provided, the stock solutions may be used for up to one year and must be replaced sooner if verification from an independent source indicates a problem.

- 7.3. Working ICP LCS/MS spike solution: Prepare the ICP LCS/MS working spike solutions from custom stock standards to the final concentration listed in Table II. The working spike must be prepared in a matrix of 5% HNO<sub>3</sub>. This acid (5 mL of concentrated HNO<sub>3</sub> per 100 mL) must be added to the volumetric flask before the addition of the stock standard aliquot. The working ICP LCS solution must be made fresh every three months.
- 7.4. Working GFAA LCS/MS spike solution: Prepare the GFAA working LCS/MS spike solutions by diluting the custom stock solution (7.2) 200x. The working spike solution must be prepared in a matrix of 5% HNO<sub>3</sub>. This acid (5 mL of concentrated HNO<sub>3</sub> per 100 mL) must be added to the volumetric flask before the addition of the stock standard aliquot. The working GFAA LCS/MS solution must be made fresh every three months.
- 7.5. The LCS and MS samples must contain all the elements designated for analysis in each batch of samples. If a non-routine element is required that is not contained in the custom STL North Canton solution, the individual facility must purchase a solution from the designated vendor that will cover the additional analyte(s) of interest and provide for a final spike concentration that is appropriate to the determinative method.
- 7.6. Aqueous laboratory control samples (LCSW) and matrix spike samples are prepared as described in Sections 9.5 and 9.6. Refer to Tables II through IV (Appendix A) for details regarding the stock, working standard and final digestate spike concentrations for ICP and GFAA LCS and matrix spike preparations.
- 7.7. Nitric acid (HNO<sub>3</sub>), concentrated, trace metal grade or better.
- 7.8. Nitric acid, 1:1 - dilute concentrated HNO<sub>3</sub> with an equal volume of reagent water.
- Note:** When preparing diluted acids always add acid to water. If the water is added to the acid a violent reaction may occur.
- 7.9. Hydrochloric acid (HCl), concentrated, trace metal grade or better.
- 7.10. Hydrochloric acid, 1:1 - dilute concentrated HCl with an equal volume of reagent water.
- Note:** When preparing diluted acids always add acid to water. If the water is added to the acid a violent reaction may occur.
- 7.11. 30% Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), reagent grade.

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

- 8.1. Sample holding time for metals included under the scope of this SOP is 180 days from the date of collection to the date of analysis.
- 8.2. Soil samples do not require preservation but must be stored at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$  until the time of analysis.

## 9. QUALITY CONTROL

Table V (Appendix A) provides a summary of quality control requirements including type, frequency, acceptance criteria and corrective action.

### 9.1. Initial Demonstration of Capability

Prior to analysis of any analyte using Method 3050B the following requirements must be met.

- 9.1.1. Method Detection Limit (MDL) - An MDL must be determined for each analyte/matrix prior to the analysis of any samples. The MDL is determined using seven replicates of reagent water, spiked with all the analytes of interest, that have been carried through the entire analytical procedure. MDLs must be redetermined on an annual basis in accordance with 40 CFR Part 136 Appendix B requirements or verified as detailed in STL North Canton QA Policy QA-005. The spike level must be between the calculated MDL and 10X the MDL to be valid. The result of the MDL determination must be below the STL North Canton reporting limit.
- 9.1.2. Initial Demonstration Study- this requires the analysis of four QC check samples. The QC check sample is a well characterized laboratory generated sample used to monitor method performance, which should contain all the analytes of interest. The results of the initial demonstration study must be acceptable before analysis of samples may begin. The results of the initial demonstration study may be used to extend a method for the analysis of other elements provided all acceptance criteria are met.
  - 9.1.2.1. Four aliquots of the check sample (LCS) are prepared and analyzed using the procedures detailed in this SOP and the determinative SOPs.
  - 9.1.2.2. Calculations and acceptance criteria for QC check samples are given in the determinative SOPs (CORP-MT-0001, CORP-MT-0003).

- 9.2. Preparation Batch - A group of up to 20 samples that are of the same matrix and are processed together using the same procedures and reagents. The preparation batch must contain a method blank, a LCS and a matrix spike/matrix spike duplicate. In some cases, at client request, it may be appropriate to process a matrix spike and sample duplicate in place of the MS/MSD. If clients specify specific samples for MS/MSD, the batch may contain multiple MS/MSD pairs.
- 9.3. Sample Count - Laboratory generated QC samples (method blanks, LCS, MS/MSD) are not counted towards the maximum 20 samples in a batch. Field QC samples are included in the batch count.
- 9.4. Method Blank (MB) - One method blank must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. Criteria for the acceptance of blanks are contained within the individual analytical method SOP's. If the method blank does not meet the criteria contained within the analytical method SOPs, the blank and all associated samples in the batch must be redigested.
- 9.5. Laboratory Control Sample (LCS) - One aqueous LCS must be processed with each preparation batch. The LCS must contain all analytes of interest and must be carried through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines. Criteria for the acceptance of LCS results are contained within the individual analytical method SOP's. Corrective action when LCS results fail to meet control limits will be repreparation and reanalysis of the batch. Tables II and III provide the details regarding the stock, working standards and final spike concentrations for ICP and GFAA. Refer to Section 7.3 or 7.4 for instructions on preparation of the aqueous LCS.
- 9.5.1. The LCS is prepared by spiking a 1mL or 1 g aliquot of reagent water with 2 mL of the working LCS/MS spike solution (7.3 or 7.4). The LCS is then processed as described in either Section 11.10.
- 9.6. Matrix Spike/Matrix Spike Duplicate (MS/MSD) - One MS/MSD pair must be processed for each preparation batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives

(DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Samples identified as field blanks cannot be used for MS/MSD analysis. If any analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include repreparation and reanalysis of the batch. Corrective action when MS results fail to meet control limits does not include repreparation of samples unless the results indicate that a spiking error may have occurred. Tables II through IV provide the details regarding the stock, working standards and final matrix spike concentrations for ICP and GFAA. Refer to Sections 7.3 and 7.4 for instructions on preparation of the working matrix spike solutions.

9.6.1. The soil matrix spike sample is prepared by spiking a 1 g aliquot of a sample with 2mL of the working LCS/MS spike solution (7.3 or 7.4). The matrix spike sample is then processed as described in either Section 11.10.

9.7. Quality Assurance Summaries - Certain clients may require specific project or program QC which may supersede the SOP requirements. Quality Assurance Summaries (QAS) should be developed to address these requirements.

## 10. CALIBRATION AND STANDARDIZATION

10.1. Hotplate or block temperature must be verified daily for each unit used and must be recorded on either the metals preparation log or in a hotplate temperature logbook. The hotplate temperature should be verified by measuring the temperature of a beaker of reagent water placed on each hotplate. For block digestors, use a tube containing water.

## 11. PROCEDURE

11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

11.3. The heating procedures are carried out in a properly functioning hood.



- 11.4. All samples are to be checked out of sample control with the chain of custody documentation filled out completely.
- 11.5. Proper sample identification is extremely important in any preparation procedure. Labeling of beakers and bottles must be done in a manner to ensure connection with the proper sample. The use of automatic label printing programs is recommended to reduce transcription errors (QuantIMS option).
- 11.6. Samples are typically logged in as either waters or soils. Wastes such as organic liquids or sludges and tissues (animal/vegetable) are usually logged in with solid test codes. When initiating prep, examine the sample to see if the sample matches the matrix designation. If the sample is logged in as aqueous but it appears more like a waste (biphasic, sludge like, organic liquid, lots of sediment etc.) contact the lab supervisor or project administrator for further instructions. In some cases it may be more appropriate to process these samples as solids.
- 11.7. If possible prepare all the samples of a project at the same time to minimize the QC required and streamline the flow of the project through the lab and reporting group.
- 11.8. In some cases, both AA and ICP digests are required on each sample. It is recommended that both aliquots be weighed out and processed at the same time.
- 11.9. Guidelines are provided in the appendices on procedures to minimize contamination of samples and standards.
- 11.10. Preparation of Soils, Sediments and Sludges for Analysis by GFAA, ICP, ICP/MS and FLAA.
  - 11.10.1. Mix sample thoroughly by stirring with a clean plastic or wooden spoon or spatula. Some project plans may require, when possible, transfer of the entire sample from the original container to butcher paper or a clean beaker, mixing thoroughly, then returning the sample to its original container.
  - 11.10.2. For each digestion procedure required (i.e., ICP or GFAA), weigh a 1.0 portion of solid and record the exact weight to the nearest 0.01 g. A 2 g sample size may also be used if needed to meet the reporting limits.
  - 11.10.3. Measure additional aliquots of the designated samples for the MS and MSD analyses.

- 11.10.4. Spike each of the MS and MSD aliquots with 2 mL of the working LCS/MS spiking solution (7.3 or 7.4).
- 11.10.5. Prepare a beaker for the method blank.
- 11.10.6. Prepare a beaker for the LCS. Add 2 mL of the working LCS/MS spiking solution (7.3 or 7.4).
- 11.10.7. Add 10 mL of 1:1 HNO<sub>3</sub> and mix the sample.
- 11.10.8. Heat sample to 95° ± 4° C and reflux for 10 minutes without boiling, using a vapor recovery device.  
  
**Note: DO NOT ALLOW SAMPLE TO BOIL OR GO DRY during any part of the digestion.** Doing so will result in the loss of analyte and the sample must be reprepared.
- 11.10.9. Allow sample to cool, if necessary.
- 11.10.10. Add 5 mL of concentrated HNO<sub>3</sub> and replace vapor recovery device.
- 11.10.11. Reflux at 95° ± 4° C for 30 minutes. (Add reagent water as needed to ensure that the volume of solution is not reduced to less than 5 mL.)
- 11.10.12. If brown fumes are observed, repeat step 11.10.10 until no more fumes are evolved.
- 11.10.13. Allow the samples to cool, if necessary.
- 11.10.14. Add approximately 2 mL of reagent water and 1 mL of 30 % H<sub>2</sub>O<sub>2</sub>. Care must be taken to ensure that losses do not occur due to excessively vigorous effervescence.
- 11.10.15. Replace the vapor recovery device and heat sample until effervescence subsides.
- 11.10.16. Allow the sample to cool, if necessary.
- 11.10.17. Continue adding 30% H<sub>2</sub>O<sub>2</sub> in 1 mL aliquots with warming until effervescence is minimal or sample appearance is unchanged.

**Note:** Do not add more than a total of 10 mL of 30 % H<sub>2</sub>O<sub>2</sub>.

- 11.10.18. Continue heating at  $95^{\circ} \pm 4^{\circ}$  C for 75 minutes until the volume is reduced to approximately 5-10 mL. Alternatively the sample may be heated for 2 hours.
- 11.10.19. If the sample is being prepared for ICP or FLAA analyses add 10 mL of concentrated HCl and reflux for an additional 10 minutes without boiling. This step is omitted for analysis by ICP-MS and GFAA

**Note:** Antimony and silver have poor solubility in dilute nitric acid solution. Therefore it is strongly recommended that these elements are determined by the ICP procedure that includes HCl as the final digestion acid.

- 11.10.20. Allow the sample to cool.
- 11.10.21. Wash down beaker walls and vapor recovery device with reagent water.
- 11.10.22. Filter sample through Whatman 41 filter paper or equivalent into a graduated cylinder or equivalent or a pre-weighed bottle. Other measuring bottles (for example, Corning Snap Seals™) may be used if their accuracy is documented and is better than  $\pm 2\%$ . Rinse beaker and filter paper with reagent water to ensure complete sample transfer.

**Note:** In place of filtering, the samples, after dilution and mixing, may be centrifuged or allowed to settle by gravity overnight to remove insoluble material

- 11.10.23. Dilute sample to 100 mL or 100g with reagent water. The sample is now ready for analysis.

**Note:** This SOP allows for samples to be weighed instead of measured volumetrically. This assumes the density of the diluted sample is close to 1.0 g/mL (See Section 17.1.2).

## 12. DATA ANALYSIS AND CALCULATIONS

Not Applicable.

## 13. METHOD PERFORMANCE

- 13.1. Method performance is determined by the analysis of matrix spike and matrix spike duplicate samples as well as method blanks and laboratory control samples. In general, the matrix spike recovery should fall within  $\pm 20\%$  and the matrix spike duplicates should compare

within 20% RPD. Method blanks must meet the criteria specified in the determinative SOPs. The laboratory control samples should recover within 20% of the true value until in house control limits are established. Acceptance criteria are given in the determinative SOPs.

- 13.2. The initial demonstration study as detailed in Section 9.1.2 must be acceptable before the analysis of field samples under this SOP may begin. The results of the initial demonstration study may be used to extend a method for the analysis of other elements provided all acceptance criteria are met.

- 13.3. Training Qualification:

The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

#### 14. **POLLUTION PREVENTION**

- 14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

#### 15. **WASTE MANAGEMENT**

- 15.1. Waste generated in the procedure must be segregated and disposed according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.
- 15.2. Standards should be purchased and prepared in volumes consistent with laboratory use to minimize the volume of expired standards to be disposed.

#### 16. **REFERENCES**

- 16.1. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW-846, 3rd Edition, Final Update III, December 1996. Method 3050B.
- 16.2. CORP-MT-0001, Inductively Coupled Plasma-Atomic Emission Spectroscopy, Spectrometric Method for Trace Element Analysis of Water and Wastes, Method 6010B and Method 200.7.
- 16.3. CORP-MT-0003, Graphite Furnace Atomic Absorption Spectroscopy, SW846 Method 7000A and MCAWW 200 Series Methods.

- 16.4. NC-MT-0002, Inductively Coupled Plasma-Mass Spectrometry, EPA Methods 6020 and 200.8.
- 16.5. QA-003, STL North Canton QC Program.
- 16.6. QA-004, Rounding and Significant Figures.
- 16.7. QA-005, Method Detection Limits.
- 17. **MISCELLANEOUS (TABLES, APPENDICES, ETC. . . )**
  - 17.1. Modifications/Interpretations from reference method.
    - 17.1.1. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the MDL. This SOP states that the method blank must not contain any analyte of interest at or above the reporting limit. Common lab contaminants, as defined in the determinative SOPs, are allowed up to two times the reporting limit in the blank following consultation with the client.
    - 17.1.2. This SOP allows for aqueous samples to be weighed instead of measured volumetrically. This assumes the density of the sample is close to 1.0 g/mL. Samples with large amounts of sediment or suspended solids, sludges, non-aqueous liquids must be processed volumetrically. Weighing samples directly into the digestion vessel minimizes the potential for cross contamination, offers improved accuracy over the use of graduated cylinders (comparable to volumetric flask accuracy), uses less glassware and is more efficient.
  - 17.2. Modifications from previous SOP

17.2.1. ICP/MS has been added as an appropriate determinative technique.

17.2.2. The approved analyte list may be modified provided that additional elements meet the requirements in section 13.

17.2.3. Directions for digestion for set time periods rather than reduction to set volumes have been added.

17.2.4. The order of two steps in the digestion has been changed. (See section 11.10.20)

17.2.5. Definition of the method as determining total environmentally available metals has been added.

17.3. Facility Specific SOPs

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none. Refer to the Appendices for any facility specific information required supporting this SOP.

17.4. Documentation and Record Management

The preparation benchsheet should, at a minimum, include the following information:

- Preparation date, analyst initials, matrix, prep type (ICP or GFAA).
- Sample ID; initial weight/volume and final weight/volume.
- Standards Documentation (source, lot, prep date, volume added).
- Analyst initials.

Figure 1. Soil Sample Preparation (Section 11.10)





Figure 2. Soil Sample Preparation (continued)





8771026

ACID DIGESTION OF SOILS, SW846 METHOD 3050B

SOP No. CORP-IP-0002NC

Revision No. 2.4

Revision Date: 02/19/03

Page 19 of 26

---

APPENDIX A - TABLES

## APPENDIX A

### TABLES

## APPENDIX A - TABLES

TABLE I. Method 3050A Approved Analyte List

ELEMENT	Symbol	CAS Number
Aluminum	Al	7429-90-5
Antimony	Sb	7440-36-0
Arsenic	As	7440-38-2
Barium	Ba	7440-39-3
Beryllium	Be	7440-41-7
Cadmium	Cd	7440-43-9
Calcium	Ca	7440-70-2
Chromium	Cr	7440-47-3
Cobalt	Co	7440-48-4
Copper	Cu	7440-50-8
Iron	Fe	7439-89-6
Lead	Pb	7439-92-1
Magnesium	Mg	7439-95-4
Manganese	Mn	7439-96-5
Molybdenum	Mo	7439-98-7
Nickel	Ni	7440-02-0
Potassium	K	7440-09-7
Selenium	Se	7782-49-2
Silver	Ag	7440-22-4
Sodium	Na	7440-23-5
Thallium	Tl	7440-28-0
Vanadium	V	7440-62-2
Zinc	Zn	7440-66-6

TABLE II. ICP and FLAA Soil Matrix Spike and Aqueous LCS Levels

ELEMENT	Working LCS/MS Standard (mg/L)	Aqueous LCS/MS Level* (ug/L)	Soil MS Level ** (mg/Kg)
Aluminum	100	2000	200
Antimony	25	500	50
Arsenic	100	2000	200
Barium	100	2000	200
Beryllium	2.5	50	5
Cadmium	2.5	50	5
Calcium	2500	50000	5000
Chromium	10	200	20
Cobalt	25	500	50
Copper	12.5	250	25
Iron	50	1000	100
Lead	25	500	50
Lithium	50	1000	100
Magnesium	2500	50000	5000
Manganese	25	500	50
Molybdenum	50	1000	100
Nickel	25	500	50
Potassium	2500	50000	5000
Selenium	100	2000	200
Silver	2.5	50	5
Sodium	2500	50000	5000
Strontium	50	1000	100
Thallium	100	2000	200
Vanadium	25	500	50
Zinc	25	500	50
Boron	50	1000	100
Tin	100	2000	200
Titanium	50	1000	100

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 2.0 mL working spike (7.3) to 100 mL of sample.

\*\* Final soil spike concentration based on the addition of 2.0 mL working spike (7.3) to 1.0 g of sample/100 mL final volume (assumes 100% solids).

## APPENDIX A - TABLES

TABLE III. ICPMS Soil Matrix Spike and Aqueous LCS Levels

ELEMENT	Working LCS/MS Standard (mg/L)	Aqueous LCS/MS Level* (ug/L)	Soil MS Level ** (mg/Kg)
Aluminum	100	1000	100
Antimony	10	100	10
Arsenic	10	100	10
Barium	10	100	10
Beryllium	10	100	10
Cadmium	10	100	10
Calcium	100	1000	100
Chromium	10	100	10
Cobalt	10	100	10
Copper	10	100	10
Iron	100	1000	100
Lead	10	100	10
Magnesium	100	1000	100
Manganese	10	100	10
Molybdenum	10	100	10
Nickel	10	100	10
Potassium	100	1000	100
Selenium	10	100	10
Silver	10	100	10
Sodium	100	1000	100
Strontium	10	100	10
Thallium	10	100	10
Vanadium	10	100	10
Zinc	10	100	10
Boron	10	100	10
Tin	10	100	10
Titanium	10	100	10
Zirconium	10	100	10

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 1.0 mL working spike (7.4) to 100 mL of sample.

\*\* Final soil spike concentration based on the addition of 1.0 mL working spike (7.4) to 1.0 g of sample/100 mL final volume (assumes 100% solids).

TABLE IV. GFAA Soil Matrix Spike and Aqueous LCS Spike Levels

ELEMENT	Stock LCS/MS Standard (mg/L)	Working LCS/MS Standard (ug/L)	Aqueous LCS/ MS Level * (ug/L)	Soil MS Level** (mg/Kg)
Arsenic	400	2000	40	4
Selenium	400	2000	40	4
Lead	400	2000	40	4
Thallium	400	2000	40	4
Antimony	400	2000	40	4
Cadmium	40	200	4	0.4
Chromium	100	500	10	1
Silver	50	250	5	0.5

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 1.0 mL working spike (7.5) to 50 mL of sample.

\*\* Final soil spike concentration based on the addition of 2.0 mL working spike (7.5) to 1.0 g of sample/100 mL final volume (assumes 100% solids).

TABLE V. Summary of Quality Control Requirements

QC PARAMETER	FREQUENCY	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
Method Blank	One per sample preparation batch of up to 20 samples.	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	Redigest and reanalyze samples.
Laboratory Control Sample (LCS)	One per sample preparation batch of up to 20 samples.	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	Redigest and reanalyze all samples associated with the LCS.
Matrix Spike	One per sample preparation batch of up to 20 samples.	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	Reprep not required unless preparation error suspected.
Matrix Spike Duplicate	See Matrix Spike	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	See Corrective Action for Matrix Spike.

8771032

ACID DIGESTION OF SOILS, SW846 METHOD 3050B

SOP No. CORP-IP-0002NC

Revision No. 2.4

Revision Date: 02/19/03

Page 25 of 26

APPENDIX B – CONTAMINATION CONTROL GUIDELINES

---

## **APPENDIX B**

### **CONTAMINATION CONTROL GUIDELINES**

**APPENDIX B – CONTAMINATION CONTROL GUIDELINES**

---

**APPENDIX B. CONTAMINATION CONTROL GUIDELINES**

**The following procedures are strongly recommended to prevent contamination:**

All work areas used to prepare standards and spikes should be cleaned before and after each use.

All glassware should be washed with detergent and tap water and rinsed with 1:1 nitric acid followed by deionized water.

Proper laboratory housekeeping is essential in the reduction of contamination in the metals laboratory. All work areas must be kept scrupulously clean.

Powdered or latex gloves must not be used in the metals laboratory since the powder contains silica and zinc, as well as other metallic analytes. Only vinyl or nitrile gloves should be used in the metals laboratory.

Glassware should be periodically checked for cracks and etches and discarded if found. Etched glassware can cause cross contamination of any metallic analytes.

Autosampler trays should be covered to reduce the possibility of contamination. Trace levels of elements being analyzed in the samples can be easily contaminated by dust particles in the laboratory.

**The following are helpful hints in the identification of the source of contaminants:**

Reagents or standards can contain contaminants or be contaminated with the improper use of a pipette.

Improper cleaning of glassware can cause contamination.

Separate glassware if an unusually high sample is analyzed and soak with sulfuric acid prior to routine cleaning.



8771034

Controlled Copy

Copy No. \_\_\_\_\_

Implementation Date: 4/25/03

SOP No. CORP-IP-0003NC

Revision No. 1.4

Revision Date: 02/19/03

Page 1 of 37

## STL STANDARD OPERATING PROCEDURE

### TITLE: ACID DIGESTION OF AQUEOUS SAMPLES BY SW846 AND MCAWW 200 SERIES METHODS

(Supersedes: Revision 1.3 Dated 09/25/01)

Prepared by: Lisa McCall 4.25.03  
Date

Reviewed by: Patricia J. M. 3/26/03  
Technology Specialist Date

Reviewed by: Beth Lambert 3/4/03  
Quality Assurance Manager Date

Reviewed by: Deborah L. Bricker 3/31/03  
Environmental Health and Safety Date

Reviewed by: Chad M. Cole 4/1/03  
Laboratory Director Date

Reviewed by: Mark Bruma 4/1/03  
Corporate Technical Director Technology

Date

#### Proprietary Information Statement

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: © 2000 SEVERN TRENT LABORATORIES, INC ALL RIGHTS RESERVED

APPENDIX A - TABLES

---

**TABLE OF CONTENTS**

1. SCOPE AND APPLICATION.....	3
2. SUMMARY OF METHOD.....	4
3. DEFINITIONS.....	5
4. INTERFERENCES.....	5
5. SAFETY.....	6
6. EQUIPMENT AND SUPPLIES.....	7
7. REAGENTS AND STANDARDS.....	8
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE.....	10
9. QUALITY CONTROL.....	10
10. CALIBRATION AND STANDARDIZATION.....	13
11. PROCEDURE.....	13
12. DATA ANALYSIS AND CALCULATIONS.....	19
13. METHOD PERFORMANCE.....	19
14. POLLUTION PREVENTION.....	19
15. WASTE MANAGEMENT.....	19
16. REFERENCES.....	20
17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . )..	20

**LIST OF APPENDICES:**

APPENDIX A - TABLES.....	32
APPENDIX B - CONTAMINATION CONTROL GUIDELINES .....	38

**1. SCOPE AND APPLICATION**

- 1.1. This procedure describes the preparation of aqueous samples for the analysis of certain metals by Graphite Furnace Atomic Absorption (GFAA), Flame Atomic Absorption (FLAA), Inductively Coupled Plasma Atomic Emission Spectroscopy (ICP), and Inductively Coupled Plasma-Mass Spectrometry (ICP/MS) using the MCAWW 200 series methods (NPDES) and SW846 Methods 3005A, 3010A, 3020A and 7060A/7740 (RCRA).
- 1.2. The applicability of each of these preparation protocols to specific analytes is detailed in Tables I and II (Appendix A) and the applicable determinative methods are illustrated by Figures 6 and 7 (Section 17). Additional elements may be analyzed following digestion by these protocols provided that the method performance criteria specified in Section 13.0 of this SOP are met.
- 1.3. This SOP provides procedures applicable to the preparation of dissolved, suspended, total recoverable and total elements in ground water, aqueous samples, certain aqueous sludges, wastes, and biological tissues, and leachates/extracts.
- 1.4. SW-846 Method 3005A is used to prepare surface and groundwater samples for total recoverable and dissolved metals determination by FLAA, ICP and GFAA (antimony only).
- 1.5. MCAWW Method 200.7 Section 9.4 is used to prepare surface water, domestic and industrial waste samples for total recoverable and dissolved metals determination by ICP.
- 1.6. SW-846 Method 3010A is used to prepare aqueous samples, EP and mobility-procedure extracts, and wastes that contain suspended solids for total metals analysis by FLAA or ICP.
- 1.7. MCAWW Method 200.7 Section 9.3 is used to prepare surface water and wastes that contain suspended solids for total metals analysis by ICP.
- 1.8. SW-846 Method 3020A is used to prepare aqueous samples, EP and mobility-procedure extracts, and wastes that contain suspended solids for total metals by GFAA, or ICP/MS.
- 1.9. MCAWW Method 200.0 Section 4.1.3 is used to prepare surface water and wastes that contain suspended solids for total metals analysis by GFAA.

APPENDIX A - TABLES

- 1.10. MCAWW Method 200.0 Section 4.1.4 is used to surface water, domestic and industrial waste samples for total recoverable and dissolved metals determination by GFAA.
- 1.11. SW-846 Methods 7060A and 7740, respectively, contain the procedure for the preparation of aqueous samples for arsenic and selenium.
- 1.12. MCAWW Methods 206.2 and 270.2, respectively, contain the procedure for the preparation of aqueous samples for arsenic and selenium.
- 1.13. All matrices require digestion prior to analysis with the exception of analyses for dissolved metals in filtered and acidified aqueous samples. Although digestion is not specifically required by the method, some clients and regulators do require digestion of dissolved samples and this must be clarified before project initiation.

2. **SUMMARY OF METHOD**

- 2.1. Method 3005A / Method 200.7 Section 9.4 - Preparation for Total Recoverable or Dissolved Metals Analysis by FLAA or ICP Spectroscopy
  - 2.1.1. A representative aliquot of sample is heated with nitric and hydrochloric acids and substantially reduced in volume. The digestate is filtered (if necessary) and diluted to volume.
- 2.2. Method 3010A / Method 200.7 Section 9.3 - Preparation for Total Metals Analysis by FLAA or ICP Spectroscopy
  - 2.2.1. A representative aliquot of sample is refluxed with nitric acid. After the digestate has been reduced to a low volume, it is refluxed with hydrochloric acid, filtered (if necessary) and brought up to volume.
- 2.3. Method 3020A / Method 200.0 Section 4.1.3 - Preparation for Total Metals for Analysis by GFAA Spectroscopy and ICP/MS.

APPENDIX A - TABLES

---

2.3.1. A representative aliquot of sample is refluxed with nitric acid. After the digestate has been reduced to a low volume, it is cooled, filtered (if necessary) and brought up to volume.

2.4. Methods 7060A/206.2 and Methods 7740/270.2 - Preparation for Arsenic/Selenium Analysis by GFAA

2.4.1. A representative aliquot of sample is heated with nitric acid and peroxide until the digestate has been reduced to a low volume. The sample is cooled, filtered (if necessary) and brought up to volume.

2.5. Method 200.0 Section 4.1.4 - Total Recoverable GFAA Preparation (NPDES)

A representative aliquot of sample is heated with nitric acid and until the digestate has been reduced to a low volume. The sample is cooled, filtered (if necessary) and brought up to volume.

### 3. DEFINITIONS

Additional definitions of terms used in this SOP may be found in the glossary of the LQM.

3.1. Dissolved Metals: Those elements which pass through a 0.45 um membrane. (Sample is acidified after filtration).

3.2. Suspended Metals: Those elements which are retained by a 0.45 um membrane.

3.3. Total Metals: The concentration determined on an unfiltered sample following digestion.

3.4. Total Recoverable Metals: The concentration determined on an unfiltered sample following treatment with hot, dilute mineral acid.

### 4. INTERFERENCES

4.1. There are numerous routes by which samples may become contaminated. Potential sources of trace metals contamination include: metallic or metal-containing labware (e.g., talc gloves which contain high levels of zinc), containers, impure reagents, dirty glassware, improper sample transfers, dirty work areas, atmospheric inputs such as dirt and dust, etc. Be aware of potential sources of contamination and take appropriate measures to minimize or avoid them.

APPENDIX A - TABLES

---

- 4.2. The entire work area, including the bench top and fume hood, should be thoroughly cleaned on a routine schedule in order to minimize the potential for environmental contamination. Refer to Appendix B for additional contamination control guidelines.
- 4.3. Boron and silica from the glassware will migrate into the sample solution during and following sample processing. For critical low level determinations of boron and silica, only quartz and/or plastic labware should be used.
- 4.4. Physical interference effects may contribute to inaccuracies in the determinations of trace elements. Oils, solvents and other matrices may not be digested using these methods if they are not soluble with acids. If physical interferences are present, they should be documented.
- 4.5. Visual interferences or anomalies (such as foaming, emulsions, precipitates, etc.) must be documented.
- 4.6. Allowing samples to boil or go dry during digestion may result in the loss of volatile metals. If this occurs the sample must be reprepared. Antimony is easily lost by volatilization from hydrochloric acid media.
- 4.7. Precipitation of silver chloride ( $\text{AgCl}$ ) may occur when chloride ions and high concentrations of silver (i.e., greater than 1 mg/L) are present in the sample.
- 4.8. Specific analytical interferences are discussed in each of the determinative methods.

5. **SAFETY**

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL associates.
- 5.2. Eye protection that satisfies ANSI Z87.1 (as per the Chemical Hygiene Plan), laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:

APPENDIX A - TABLES

---

5.3.1. The following materials are known to be **corrosive**:

hydrochloric acid and nitric acid.

5.3.2. The following materials are known to be **oxidizing agents**:

nitric acid and hydrogen peroxide.

5.3.3. All sample digestions, including cooling of digestates, must be carried out in a fume hood.

5.4. The acidification of samples containing reactive materials may result in the release of toxic gases, such as cyanides or sulfides. Acidification of samples should be done in a fume hood. The analyst should also be aware of the potential for a vigorous reaction.

5.5. Exposure to chemicals must be maintained **as low as reasonably achievable**. Therefore, unless they are known to be non-hazardous, all samples should be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.

5.6. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL associate. The situation must be reported **immediately** to a laboratory supervisor.

5.7. Always carry bulk concentrated acid bottles in appropriate impact proof containers.

5.8. Acid/peroxide spills must be neutralized immediately, flushed with water and cleaned up using appropriate spill kits.

5.9. Discard chipped or broken beakers to prevent injury. Chipped glassware may be fire polished as an alternative to disposal.

5.10. Any and all accidents and spills must be reported to the lab supervisor or EH&S coordinator.

## 6. EQUIPMENT AND SUPPLIES

6.1. Hot plate, digestion block or other adjustable heating source capable of maintaining a temperature of 95°C ( $\pm 4$ ).

6.2. Calibrated thermometer that covers a temperature range of 0-200°C.

APPENDIX A - TABLES

- 6.3. Griffin beakers of assorted sizes or equivalent.
- 6.4. Watch glasses, ribbed or equivalent.
- 6.5. Whatman No. 4 filter paper or equivalent.
- 6.6. Funnels or equivalent filtration apparatus.
- 6.7. Centrifugation equipment (if desired method of removing particulates is centrifugation).
- 6.8. Graduated cylinder or equivalent capable of measuring 50 mL within 3% accuracy.
- 6.9. Analytical balance capable of accurately weighing to the nearest 0.01 grams.
- 6.10. Repipetors or suitable reagent dispensers.
- 6.11. Calibrated automatic pipettes with corresponding pipette tips or Class A glass volumetric pipettes.
- 6.12. Class A volumetric flasks.
- 6.13. pH indicator strips (pH range 0 - 6).
- 6.14. Plastic digestate storage bottles.

**7. REAGENTS AND STANDARDS**

- 7.1. Reagent water must be produced by a Millipore DI system or equivalent. Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks as defined in the determinative SOPs.
- 7.2. Laboratory Control Sample (LCS) and matrix spike (MS) solutions are purchased as custom STL solutions. All standards must be stored in FEP fluorocarbon or previously unused polyethylene or polypropylene bottles. Stock standard solutions must be replaced prior to the expiration date provided by the manufacturer. If no expiration date is provided, the stock solutions may be used for up to one year and must be replaced sooner if verification from an independent source indicates a problem.
- 7.3. Working ICP LCS/MS spike solution: Prepare the ICP LCS/MS working spike solution from custom stock standards to the final concentration listed in Table III. The working



APPENDIX A - TABLES

---

spike must be prepared in a matrix of 5%  $\text{HNO}_3$ . This acid (5 mL of concentrated  $\text{HNO}_3$  per 100 mL) must be added to the volumetric flask before the addition of the stock standard aliquot. The working ICP LCS solution must be made fresh every three months.

- 7.4. Working GFAA LCS/MS spike solution: Prepare the GFAA working LCS spike solution by diluting the custom stock solution (7.2) 200x. The working spike solution must be prepared in a matrix of 5%  $\text{HNO}_3$ . This acid (5 mL of concentrated  $\text{HNO}_3$  per 100 mL) must be added to the volumetric flask before the addition of the stock standard aliquot. The working GFAA LCS solution must be made fresh every three months.
- 7.5. The TCLP MS working spike solution is provided directly by the vendor, no further standard preparation is necessary.
- 7.6. The LCS and MS samples must contain all the elements designated for analysis in each batch of samples. If a non-routine element is required that is not contained in the custom STL solution, the individual facility must purchase a solution from the designated vendor that will cover the additional analyte(s) of interest and provide for a final spike concentration that is appropriate to the determinative method.
- 7.7. Aqueous laboratory control samples (LCSW) and matrix spike samples are prepared as described in Sections 9.5 and 9.6. Refer to Tables III and IV(Appendix A) for details regarding the stock, working standard and final digestate spike concentrations for ICP and GFAA LCS and matrix spike preparations.
- 7.8. Nitric acid ( $\text{HNO}_3$ ), concentrated, trace metal grade or better.
- 7.9. Nitric acid, 1:1 - dilute concentrated  $\text{HNO}_3$  with an equal volume of reagent water.
- Note:** When preparing diluted acids always add acid to water. If the water is added to the acid a violent reaction may occur.
- 7.10. Hydrochloric acid ( $\text{HCl}$ ), concentrated, trace metal grade or better.
- 7.11. Hydrochloric acid, 1:1 - dilute concentrated  $\text{HCl}$  with an equal volume of reagent water.
- Note:** When preparing diluted acids always add acid to water. If the water is added to the acid a violent reaction may occur.
- 7.12. 30% Hydrogen peroxide ( $\text{H}_2\text{O}_2$ ), reagent grade.

APPENDIX A - TABLES

8. **SAMPLE COLLECTION, PRESERVATION AND STORAGE**

- 8.1. Sample holding time for metals included under the scope of this SOP is 180 days from the date of collection to the date of analysis.
- 8.2. Aqueous samples are preserved with nitric acid to a pH of <2 and may be stored in either plastic or glass. If boron or silica are to be determined, plastic containers are preferred. Refrigeration is not required. Preservation must be verified prior to analysis.
- 8.3. For dissolved metals analysis, the samples should be filtered through a 0.45 um filter prior to preservation. Filtration must be done in the field or within 24 hours of collection.

**Note:** If a sample being analyzed for dissolved metals is found to contain sediment the analyst should contact their supervisor or group leader. The client should be notified of the problem to decide how to treat the sample.

9. **QUALITY CONTROL**

Table VI (Appendix A) provides a summary of quality control requirements including type, frequency, acceptance criteria and corrective action.

9.1. Initial Demonstration of Capability

Prior to analysis of any analyte using any method contained within this SOP the following requirements must be met:

- 9.1.1. Method Detection Limit (MDL) - An MDL must be determined for each analyte/matrix prior to the analysis of any samples. The MDL is determined using seven replicates of reagent water, spiked with all the analytes of interest, that have been carried through the entire analytical procedure. MDL's must be redetermined on an annual basis in accordance with 40 CFR Part 136 Appendix B requirements as detailed in STL QA Policy QA-005. The spike level must be between the calculated MDL and 10X the MDL to be valid. The result of the MDL determination must be below the STL reporting limit.
- 9.1.2. Initial Demonstration Study - This requires the analysis of four QC check samples. The QC check sample is a well-characterized laboratory generated sample used to monitor method performance, which should contain all the analytes of interest. The results of the initial demonstration study must be acceptable before analysis of samples may begin. The results of the initial demonstration study may be used to

## APPENDIX A - TABLES

extend a method for the analysis of other elements provided all acceptance criteria are met.

9.1.2.1. Four aliquots of the check sample (LCS) are prepared and analyzed using the procedures detailed in this SOP and the determinative SOPs.

9.1.2.2. Calculations and acceptance criteria for QC check samples are given in the determinative SOPs (CORP-MT-0001, CORP-MT-0003).

- 9.2. Preparation Batch - A group of up to 20 samples that are of the same matrix and are processed together using the same procedures and reagents. The preparation batch must contain a method blank, a LCS and a matrix spike/matrix spike duplicate. In some cases, at client request, it may be appropriate to process a matrix spike and sample duplicate in place of the MS/MSD. If clients specify specific samples for MS/MSD, the batch may contain multiple MS/MSD pairs.
- 9.3. Sample Count - Laboratory generated QC samples (method blanks, LCS) are not included in the sample count for determining the size of a preparation batch. MS/MSD are not included in the sample count unless there are multiple sets of MS/MSD per batch. In other words, the first MS/MSD are not counted; all additional MS and MSDs are counted as samples.
- 9.4. Method Blank (MB) - One method blank must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. Criteria for the acceptance of blanks are contained within the individual analytical method SOP's. If the method blank does not meet the criteria contained within the analytical method SOPs; the blank and all associated samples in the batch must be redigested.
- 9.4.1. Aqueous method blanks are prepared by taking 50 mL or 50 g of reagent water through the appropriate procedure as described in Section 11.
- 9.4.2. TCLP method blanks are prepared by taking 50 mL or 50 g of leachate fluid through the appropriate procedure as described in Section 11.
- 9.5. Laboratory Control Sample (LCS) - One aqueous LCS must be processed with each preparation batch. The LCS must contain all analytes of interest and must be carried

APPENDIX A - TABLES

---

through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines. Criteria for the acceptance of LCS results are contained within the individual analytical method SOP's. Corrective action when LCS results fail to meet control limits will be reparation and reanalysis of the batch. Refer to Section 7.3 and 7.4 for instructions on preparation of the aqueous LCS spike solution.

9.5.1. The aqueous LCS is prepared by spiking a 50 mL aliquot of reagent water with 1.0 mL of the working LCS/MS spike solution (7.3 or 7.4). The LCS is then processed through the appropriate procedure as described in Section 11.

9.6. Matrix Spike/Matrix Spike Duplicate (MS/MSD) - One MS/MSD pair must be processed for each preparation batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Samples identified as field blanks cannot be used for MS/MSD analysis. If any analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include reparation and reanalysis of the batch. Corrective action when MS results fail to meet control limits does not include reparation of samples unless the results indicate that a spiking error may have occurred.

9.6.1. The aqueous matrix spike sample is prepared by spiking a 50 mL aliquot of a sample with 1.0 mL of the working LCS/MS spike solution (7.3 or 7.4). The matrix spike sample is then processed as described in Section 11.

9.6.2. The TCLP matrix spike sample is prepared by spiking a 50 mL aliquot of a leachate with 0.5 mL of the working TCLP spike solution (7.5). The matrix spike sample is then processed as described in Section 11.

NOTE: The TCLP matrix spike must be added prior to preservation of the leachate.

9.6.3. If insufficient sample is available to process a MS/MSD, then a second LCS must be processed. The LCS pair is then evaluated according to the MS/MSD criteria.

APPENDIX A - TABLES

---

- 9.7. Quality Assurance Summaries - Certain clients may require specific project or program QC that may supersede the SOP requirements. Quality Assurance Summaries (QAS) should be developed to address these requirements.

**10. CALIBRATION AND STANDARDIZATION**

- 10.1 Hotplate temperature must be verified daily for each hotplate used and must be recorded on either the metals preparation log or in a hotplate temperature logbook. The hotplate temperature should be verified by measuring the temperature of a beaker of reagent water placed on each hotplate.

**11. PROCEDURE**

- 11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.
- 11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.
- 11.3. All digestion procedures must be carried out in a properly functioning hood.
- 11.4. All samples are to be checked out of sample control with the chain of custody documentation filled out completely.
- 11.5. Proper sample identification is extremely important in any preparation procedure. Labeling of beakers and bottles must be done in a manner to ensure connection with the proper sample.
- 11.6. Samples are typically logged in as either waters or soils. Wastes such as organic liquids or sludges and tissues (animal/vegetable) are usually logged in with solid test codes. When initiating prep, examine the sample to see if the sample matches the matrix designation. If the sample is logged in as aqueous but it appears more like a waste (biphasic, sludge like, organic liquid, lots of sediment etc.) contact the lab supervisor or project manager for further instructions. In some cases it may be more appropriate to process these samples as solids.

APPENDIX A - TABLES

---

- 11.7. If possible prepare all the samples of a project at the same time to minimize the QC required and streamline the flow of the project through the lab and reporting group.
- 11.8. In most cases, both AA and ICP digests are required on each sample. It is recommended that both aliquots be measured out and processed at the same time.
- 11.9. Guidelines are provided in the appendices on procedures to minimize contamination of samples and standards.
- 11.10. The following procedure must be followed for all aqueous sample preparations:
- 11.10.1. Measure sample pH with pH paper on a separate aliquot of sample.
- Note:** If the sample pH is  $> 2$  pH units, the client must be notified of the anomaly.
- Note:** If sample pH has already been verified and documented in sample receipt this step may be omitted.
- 11.10.2. Mix sample by shaking the container.
- 11.10.3. Measure and transfer 50 mL or 50 g of the sample into a beaker.
- Note:** This SOP allows for samples to be weighed instead of measured volumetrically.
- 11.10.4. Measure two extra aliquots of sample selected for the MS/MSD analysis. Spike each aliquot with the appropriate spiking solutions (7.3-7.5,9.6).
- 11.10.5. Measure and transfer 50 mL of reagent water into a beaker for the method blank.
- 11.10.6. Measure and transfer 50 mL of reagent water into a beaker for the LCS and add the appropriate spiking solutions (7.3-7.5,9.6).

## APPENDIX A - TABLES

11.11. Proceed to the appropriate Section for the desired method as follows:

Method 3005A or Method 200.7 Section 9.4	11.12
Method 3010A or Method 200.7 Section 9.3	11.13
Method 3020A or Method 200.0 Section 4.1.3	11.14
Method 7060A/7740 or Method 206.2/270.2	11.15
Method 200.0 Section 4.1.4	11.16

11.12. **Method 3005A / Method 200.7 Section 9.4 - Preparation for Total Recoverable or Dissolved Metals Analysis by FLAA or ICP (See Figures 1, 6 and 7)**

11.12.1. To the sample beaker, add 1 mL of concentrated  $\text{HNO}_3$  and 2.5 mL of concentrated  $\text{HCl}$ .

11.12.2. Cover with ribbed watch glass.

11.12.3. Heat at  $95^\circ\text{C}$  ( $\pm 4$ ) until volume is reduced to between 15 and 20 mL.

**NOTE: DO NOT ALLOW SAMPLE TO BOIL OR GO DRY.** Doing so will result in the loss of analyte and the sample must be reprepared.

11.12.4. Cool the beaker in a fume hood.

11.12.5. Wash down beaker walls and watch glass with reagent water.

11.12.6. Filter sample, if insoluble materials are present, through Whatman 4 filter paper that has been pre-rinsed with dilute nitric acid.

**Note:** If any samples in a preparation batch are filtered, the method blank and LCS associated with that batch must also be filtered.

**Note:** In place of filtering, the samples, after dilution and mixing, may be centrifuged or allowed to settle by gravity overnight to remove insoluble material.

APPENDIX A - TABLES

11.12.7. Rinse beaker and filter paper with reagent water to ensure complete sample transfer.

11.12.8. Adjust the final volume/mass to 50 mL or 50 g with reagent water. The sample is now ready for analysis

**11.13. Method 3010A / Method 200.7 Section 9.3 - Preparation for Total Metals Analysis by FLAA or ICP Spectroscopy (See Figures 2, 6 and 7)**

11.13.1. To the sample beaker, add 3.0mL of concentrated HNO<sub>3</sub>

11.13.2. Cover with ribbed watch glass.

11.13.3. Place beaker on hotplate 95°C (± 4) and evaporate for 4-5 hours or to low volume of 15-20 mL while ensuring that no portion of the bottom of the beaker is allowed to go dry.

**NOTE: DO NOT ALLOW SAMPLE TO BOIL OR GO DRY.** Doing so will result in the loss of analyte and the sample must be reprepared.

11.13.4. If necessary, add another 1.5 ml portion of concentrated HNO<sub>3</sub> and re-cover the beaker. Reflux 15 minutes.

11.13.5. Add 5 mL of 1:1 HCl.

11.13.6. Cover and reflux for an additional 15 minutes to dissolve precipitate or residue. Cool in a fume hood.

11.13.7. Wash down beaker walls and watch glass with reagent water.

11.13.8. Filter sample, if insoluble materials are present, through Whatman 4 filter paper.

**Note:** If any samples in the QC batch are filtered the method blank and LCS associated with that batch must also be filtered.

**Note:** In place of filtering, the samples, after dilution and mixing, may be centrifuged or allowed to settle by gravity overnight to remove insoluble material.



## APPENDIX A - TABLES

11.13.9. Rinse beaker and filter paper with reagent water to ensure complete sample transfer.

11.13.10. Adjust final volume/mass to 50 mL or 50 g with reagent water. The sample is now ready for analysis.

**11.14. Method 3020A / Method 200.0 Section 4.1.3 - Preparation for Total Metals Analysis by GFAA or Total Recoverable Metals by ICPMS (See Figures 3, 6 and 7)**

11.14.1. To the sample beaker, add 1.5 mL of concentrated  $\text{HNO}_3$ .

11.14.2. Cover with ribbed watch glass.

11.14.3. Place beaker on hotplate  $95^\circ\text{C}$  ( $\pm 4$ ) and evaporate to low volume of 15-20 mL while ensuring that no portion of the bottom of the beaker is allowed to go dry.

**NOTE: DO NOT ALLOW SAMPLE TO BOIL OR GO DRY.** Doing so will result in the loss of analyte and the sample must be reprepared.

11.14.4. If necessary, add another 1.5 mL portion of concentrated  $\text{HNO}_3$ . Recover, and reflux 15 minutes. Cool the beaker in a fume hood.

11.14.5. Filter sample, if insoluble materials are present, through Whatman 4 filter paper.

**Note:** If any samples in the QC batch are filtered the method blank and LCS associated with that batch must also be filtered.

**Note:** In place of filtering, the samples, after dilution and mixing, may be centrifuged or allowed to settle by gravity overnight to remove insoluble material.

11.14.6. Rinse beaker and filter paper with reagent water to ensure complete sample transfer.

11.14.7. Adjust final volume to 50 mL with reagent water. The sample is now ready for analysis.

**11.15. Method 7060A/7740 and Method 206.2/270.2 - Preparation for Arsenic and Selenium Analysis by GFAA (See Figures 4, 6 and 7)**

APPENDIX A - TABLES

- 11.15.1. To the sample beaker, add 1 mL of 30 %  $\text{H}_2\text{O}_2$  and 1.0 mL of 1:1  $\text{HNO}_3$ .
- 11.15.2. Heat, until the digestion is complete, at  $95^\circ\text{C} (\pm 4)$  or until the volume has been reduced to 15-20 mL.
- 11.15.3. Cool beaker.
- 11.15.4. Filter sample, if insoluble materials are present, through Whatman 4 filter paper that has been pre-rinsed with dilute nitric acid.

**Note:** If any samples in the QC batch are filtered the method blank and LCS associated with that batch must also be filtered.

**Note:** In place of filtering, the samples, after dilution and mixing, may be centrifuged or allowed to settle by gravity overnight to remove insoluble material.

- 11.15.5. Rinse beaker and filter paper with reagent water to ensure complete sample transfer.
- 11.15.6. Adjust final volume to 50 mL with reagent water. The sample is now ready for analysis.

**11.16. Method 200.0 Section 4.1.4 - Preparation for Total Recoverable GFAA Analyses.  
(See Figures 5 and 7)**

- 11.16.1. To the sample beaker, add 1.0 mL of 1:1  $\text{HNO}_3$ .
- 11.16.2. Heat, until the digestion is complete, at  $95^\circ\text{C} (\pm 4)$  or until the volume has been reduced to 15 - 20 mL.
- 11.16.3. Cool beaker.
- 11.16.4. Filter sample, if insoluble materials are present, though Whatman 4 filter paper.

**Note:** If any samples in the QC batch are filtered the method blank and LCS associated with that batch must also be filtered.

**Note:** In place of filtering, the samples, after dilution and mixing, may be centrifuged or allowed to settle by gravity overnight to remove insoluble material.

## APPENDIX A - TABLES

11.16.5. Rinse beaker and filter paper with reagent water to ensure complete sample transfer.

11.16.6. Adjust final volume to 50 mL with reagent water. The sample is now ready for analysis.

**12. DATA ANALYSIS AND CALCULATIONS**

Not Applicable.

**13. METHOD PERFORMANCE**

13.1. Method performance is determined by the analysis of matrix spike and matrix spike duplicate samples as well as method blanks and laboratory control samples. In general, the matrix spike recovery should fall within +/- 20 % and the matrix spike duplicates should compare within 20% RPD. Method blanks must meet the criteria specified in determinative SOPs. The laboratory control samples should recover within 20% of the true value until in house control limits are established. Acceptance criteria are given in the determinative SOPs.

13.2. The initial demonstration study as detailed in Section 9.1.2 must be acceptable before the analysis of field samples under this SOP may begin. The results of the initial demonstration study may be used to extend a method for the analysis of other elements provided all acceptance criteria are met.

13.3. Training Qualification:

The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

**14. POLLUTION PREVENTION**

14.1. This method allows for the proportional reduction of sample and reagent volumes to decrease waste generation.

**15. WASTE MANAGEMENT**

15.1. Waste generated in the procedure must be segregated and disposed according to the facility hazardous waste procedures. The facility EH & S coordinator should be contacted if additional information is required.

---

**APPENDIX A - TABLES**

---

- 15.2. Standards should be purchased and prepared in volumes consistent with laboratory use to minimize the volume of expired standards to be disposed.

**16. REFERENCES**

- 16.1. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW-846, 3rd Edition, Final Update I, Revision 1, July 1992. Methods 3005A, 3010A, 3020A, 7060A and 7740A.
- 16.2. Methods for the Chemical Analysis of Water and Waste (MCAWW), 1983.
- 16.3. CORP-MT-0001, Inductively Coupled Plasma-Atomic Emission Spectroscopy, Spectrometric Method for Trace Element Analysis of Water and Wastes, Method 6010A and Method 200.7.
- 16.4. CORP-MT-0003, Graphite Furnace Atomic Absorption Spectroscopy, SW846 Method 7000A and MCAWW 200 Series Methods.
- 16.5. QA-003, STL QC Program.
- 16.6. QA-004, Rounding and Significant Figures.
- 16.7. QA-005, Method Detection Limits.

**17. MISCELLANEOUS (TABLES, APPENDICES, ETC. . . )**

- 17.1. Modifications/Interpretations from reference methods.
- 17.1.1. Modifications applicable to SW-846 reference methods.
- 17.1.1.1. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the MDL. This SOP states that the method blank must not contain any analyte of interest at or above the reporting limit. Common lab contaminants are allowed up to two times the reporting limit in the blank following consultation with the client.
- 17.1.1.2. The referenced methods as well as Table 3-1 of SW-846 refer to the use of a 100 mL aliquot for digestion. This SOP requires the use of a 50 mL sample size to reduce waste generation. The use of reduced sample volumes are supported in EPA's document "Response to Public

APPENDIX A - TABLES

---

Comments Background Document, Promulgation of the Second Update to SW-846, Third Edition" dated November 3, 1994. This document stated "flexibility to alter digestion volumes is addressed and "allowed" by the table (3-1) and is also inherently allowed by specific digestion methods. Table 3-1 is only to be used as guidance when collecting samples..." EMSL-Ci has also taken the stance that "reduction in sample size and appropriate corresponding reduction in sample volume is not considered a significant change in the methodology." Additionally, in written correspondence from the Office of Solid Waste, Olliver Fordham stated "As a "representative sample" can be assured, scaling causes no loss of precision and accuracy in the analysis."

## 17.1.2. Modifications Specific to Method 3010A

17.1.2.1. Section 11.13.7 of this SOP requires the sample be reduced to a volume of 15 - 20 mL. Section 7.2 of Method 3010A states the volume should be reduced to 3 mL but also states that no portion of the bottom of the beaker should go dry. The SOP required volume is a closer approximation of the volume required to provide an adequate covering of the beaker so as to prevent the loss of critical analytes through volatilization.

17.1.2.2. The scope of 3010A has been expanded to include silver based on comparison studies with 7760A. Method 3010A consistently demonstrated improved accuracy and precision over Method 7760A in the matrices tested (reagent water, surface water and TCLP leachate) up to a concentration of 1 ppm silver.

## 17.1.3. Modifications Specific to Method 3020A

17.1.3.1. Section 11.14.3 of this SOP requires the sample be reduced to a volume of 15 - 20 mL. Section 7.2 of Method 3010A states the volume should be reduced to 3 mL but also states that no portion of the bottom of the beaker should go dry. The SOP required volume is a closer approximation of the volume required to provide an adequate covering of the beaker so as to prevent the loss of critical analytes through volatilization.

APPENDIX A - TABLES

17.1.4. Modifications Specific to Method 7060A/7740

17.1.4.1. Methods 7060A and 7740A incorporate the use of a two step dilution to accommodate the addition of a nickel nitrate modifier. This SOP performs the dilution directly in one step and omits the addition of the modifier. The modifier is added automatically at the instrument by direct injection into the furnace.

17.1.5. Modifications Specific to MCAWW Methods

It was determined by technical review that several of the MCAWW methods were equivalent to the SW-846 methods and therefore were combined under the scope of this SOP as described in Section 11.0. The nature of the differences were deemed insignificant in regards to the amount of acid added and the evaporative volume based on the flexibility allowed by the methods (i.e., add additional acid as required) and the subjective wording of the methods (i.e., evaporate to near dryness vs. an exact volume).

17.2. Modifications from previous SOP

17.2.1. Added ICP/MS to the digestion procedures.

17.3. Facility Specific SOPs

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none. Refer to the Appendices for any facility specific information required to support this SOP.


17.4. Documentation and Record Management

The preparation benchsheet should, at a minimum, include the following information:

- Preparation date, analyst name, matrix, prep type (ICP or GFAA), SOP reference.
- Sample ID; initial weight/volume and final weight/volume.
- Standards Documentation (source, lot, prep date, volume added).
- Analyst Signature.
- Reviewer's Signature and date.


## APPENDIX A - TABLES

Figure 1. Method 3005A / Method 200.7 Section 9.4 (Section 11.12)



APPENDIX A - TABLES


Figure 2. Method 3010A / Method 200.7 Section 9.3 (Section 11.13)






## APPENDIX A - TABLES

Figure 3. Method 3020A / Method 200.0 Section 4.1.3 (Section 11.14)



APPENDIX A - TABLES

Figure 4. Method 7060A/7740A and Method 206.2/270.2 (Section 11.15)



## APPENDIX A - TABLES

Figure 5. Method 200.0 Section 4.1.4 (Section 11.16)






Figure 6. Overview of SW846 Aqueous Preparation Methods by Determinative Method



TAL list: Al, Sb, As, Ba, Be, Cd, Ca, Cr, Co, Cu, Fe, Pb, Mg, Mn, Ni, K, Se, Ag, Na, Tl, V, Zn

**Figure 7. Overview of MCAWW Aqueous Preparation Methods by Determinative Technique**

TAL list: Al, Sb, As, Ba, Be, Cd, Ca, Cr, Co, Cu, Fe, Pb, Mg, Mn, Ni, K, Se, Ag, Na, Tl, V, Zn

APPENDIX A - TABLES

APPENDIX A

TABLES

## APPENDIX A - TABLES

**TABLE I. Approved Preparation Method Analytes - SW846**

ELEMENT	Symbol	CAS Number	3005A	3010A	3020A	7060A 7740
Aluminum	Al	7429-90-5	X	X		
Antimony	Sb	7440-36-0	X			
Arsenic	As	7440-38-2	X	X		X
Barium	Ba	7440-39-3	X	X		
Beryllium	Be	7440-41-7	X	X	X	
Cadmium	Cd	7440-43-9	X	X	X	
Calcium	Ca	7440-70-2	X	X		
Chromium	Cr	7440-47-3	X	X	X	
Cobalt	Co	7440-48-4	X	X	X	
Copper	Cu	7440-50-8	X	X		
Iron	Fe	7439-89-6	X	X		
Lead	Pb	7439-92-1	X	X	X	
Magnesium	Mg	7439-95-4	X	X		
Manganese	Mn	7439-96-5	X	X		
Molybdenum	Mo	7439-98-7	X	X	X	
Nickel	Ni	7440-02-0	X	X		
Potassium	K	7440-09-7	X	X		
Selenium	Se	7782-49-2	X	X		X
Silver	Ag	7440-22-4	X	X		
Sodium	Na	7440-23-5	X	X		
Thallium	Tl	7440-28-0	X	X	X	
Vanadium	V	7440-62-2	X	X	X	
Zinc	Zn	7440-66-6	X	X		

X - Designates that the preparation method is approved for an element

**Note:** Additional elements may be analyzed following digestion by these protocols provided the method performance criteria specified in Section 13.0 of the SOP are met.

APPENDIX A - TABLES

**TABLE II. Approved Preparation Method Analytes - NPDES**

ELEMENT	Symbol	CAS Number	200.7 (9.4)	200.7 (9.3)	200.0 (4.1.4)	200.0 (4.1.3)	206.2 270.2
Aluminum	Al	7429-90-5	X	X			
Antimony	Sb	7440-36-0	X	X	X	X	
Arsenic	As	7440-38-2	X	X			X
Boron	B	7440-42-8	X	X			
Barium	Ba	7440-39-3	X	X	X	X	
Beryllium	Be	7440-41-7	X	X	X	X	
Cadmium	Cd	7440-43-9	X	X	X	X	
Calcium	Ca	7440-70-2	X	X			
Chromium	Cr	7440-47-3	X	X	X	X	
Cobalt	Co	7440-48-4	X	X	X	X	
Copper	Cu	7440-50-8	X	X	X	X	
Iron	Fe	7439-89-6	X	X	X	X	
Lead	Pb	7439-92-1	X	X	X	X	
Magnesium	Mg	7439-95-4	X	X			
Manganese	Mn	7439-96-5	X	X	X	X	
Molybdenum	Mo	7439-98-7	X	X	X	X	
Nickel	Ni	7440-02-0	X	X	X	X	
Potassium	K	7440-09-7	X	X			
Selenium	Se	7782-49-2	X	X			X
Silicon	Si	7631-86-9	X	X			
Silver	Ag	7440-22-4	X	X	X	X	
Sodium	Na	7440-23-5	X	X			
Thallium	Tl	7440-28-0	X	X	X	X	
Vanadium	V	7440-62-2	X	X	X	X	
Zinc	Zn	7440-66-6	X	X			

X - Designates that the preparation method is approved for an element

**Note:** Additional elements may be analyzed following digestion by these protocols provided the method performance criteria specified in Section 13.0 of the SOP are met.



**8771066****ACID DIGESTION OF AQUEOUS SAMPLES BY SW846 AND  
MCAWW 200 SERIES METHODS**

SOP No. CORP-IP-0003NC

Revision No. 1.4

Revision Date: 02/19/03

Page 33 of 37

**APPENDIX A - TABLES****TABLE III. ICP and FLAA Matrix Spike and Aqueous Laboratory Control Sample Levels**

ELEMENT	Working LCS/MS Standard (mg/L)	Aqueous LCS/ MS Level * (ug/l)
Aluminum	100	2000
Antimony	25	500
Arsenic	100	2000
Barium	100	2000
Beryllium	2.5	50
Cadmium	2.5	50
Calcium	2500	50000
Chromium	10	200
Cobalt	25	500
Copper	12.5	250
Iron	50	1000
Lead	50	500
Magnesium	2500	50000
Manganese	25	500
Molybdenum	50	1000
Nickel	25	500
Phosphorous	500	10000
Potassium	2500	50000
Selenium	100	2000
Silver	2.5	50
Sodium	2500	50000
Thallium	100	2000
Vanadium	25	500
Zinc	25	500
Boron	50	1000
Tin	100	2000
Titanium	50	1000

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 1.0 mL working spike (7.3) to 50 mL of sample.

APPENDIX A - TABLES

**TABLE IV. GFAA Matrix Spike and Aqueous LCS Spike Levels**

ELEMENT	Stock LCS/MS Standard (mg/L)	Working LCS/MS Standard (ug/L)	Aqueous LCS/ MS Level * (ug/l)
Arsenic	400	2000	40
Selenium	400	2000	40
Lead	400	2000	40
Thallium	400	2000	40
Antimony	400	2000	40
Cadmium	40	200	4
Chromium	100	500	10
Silver	50	250	5

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 1.0 mL working spike (7.4) to 50 mL of sample.

**TABLE V. TCLP Reporting Limits, Regulatory Limits and Matrix Spike Levels**

ELEMENT	RL (ug/L)	Regulatory Limit (ug/L)	Spike Level (ug/L)*
Arsenic	500	5000	5000
Barium	10000	100000	50000
Cadmium	100	1000	1000
Chromium	500	5000	5000
Lead	500	5000	5000
Selenium	250	1000	1000
Silver	500	5000	1000

\* Levels shown indicate the spike concentration in the final digestate of the aqueous LCS or matrix spike based on the addition of 0.5 mL working spike (7.4) to 50 mL of sample

8771068

ACID DIGESTION OF AQUEOUS SAMPLES BY SW846 AND  
MCAWW 200 SERIES METHODS

SOP No. CORP-IP-0003NC

Revision No. 1.4

Revision Date: 02/19/03

Page 35 of 37

## APPENDIX A - TABLES

TABLE VI. Summary of Quality Control Requirements

QC PARAMETER	FREQUENCY	ACCEPTANCE CRITERIA	CORRECTIVE ACTION
Method Blank	One per sample preparation batch of up to 20 samples.	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	Redigest and reanalyze samples associated with the method blank.
Laboratory Control Sample (LCS)	One per sample preparation batch of up to 20 samples.	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	Redigest and reanalyze all samples associated with the LCS.
Matrix Spike	One per sample preparation batch of up to 20 samples.	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	Reprep not required unless preparation error suspected.
Matrix Spike Duplicate	See Matrix Spike	Refer to determinative SOPs: NC-MT-0002 CORP-MT-0001 CORP-MT-0003	See Corrective Action for Matrix Spike.

**APPENDIX B**  
**CONTAMINATION CONTROL GUIDELINES**

**APPENDIX B – CONTAMINATION CONTROL GUIDELINES****APPENDIX B. CONTAMINATION CONTROL GUIDELINES**

**The following procedures are strongly recommended to prevent contamination:**

All work areas used to prepare standards and spikes should be cleaned before and after each use.

All glassware should be washed with detergent and tap water and rinsed with 1:1 nitric acid followed by deionized water.

Proper laboratory housekeeping is essential in the reduction of contamination in the metals laboratory. All work areas must be kept scrupulously clean.

Powdered or Latex Gloves must not be used in the metals laboratory since the powder contains silica and zinc, as well as other metallic analytes. Only vinyl or nitrile gloves should be used in the metals laboratory.

Glassware should be periodically checked for cracks and etches and discarded if found. Etched glassware can cause cross contamination of any metallic analytes.

Autosampler trays should be covered to reduce the possibility of contamination. Trace levels of elements being analyzed in the samples can be easily contaminated by dust particles in the laboratory.

**The following are helpful hints in the identification of the source of contaminants:**

Reagents or standards can contain contaminants or be contaminated with the improper use of a pipette.

Improper cleaning of glassware can cause contamination.

Separate glassware if an unusually high sample is analyzed and soak with sulfuric acid prior to routine cleaning.



# STL

8771071

SOP No. COI-GC-005  
Revision No. 4  
Revision Date: 05/21/04  
Effective Date: 06/04/04  
Page: 1 of 19

## STANDARD OPERATING PROCEDURE

### **TITLE: SAMPLE PREPARATION AND THE DETERMINATION OF DISSOLVED GASES IN WATER BY USING GC HEADSPACE EQUILIBRATION TECHNIQUE (EPA RSKSOP-175 MODIFIED)**

(SUPERSEDES: REV. 3)

Prepared by: David Kammerer

Reviewed by: Lerylina C Yabut  
Technical Specialist, Lerylina Yabut

Approved by: Linda Scharpenberg  
Quality Assurance Manager, Linda Scharpenberg

Approved by: Linda Scharpenberg  
Environmental Health and Safety Coordinator, Linda Scharpenberg

Approved by: Elizabeth A. Winger  
Laboratory Director, Elizabeth A. Winger

#### Proprietary Information Statement

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use, and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the valuation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

**THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY:**

**COPYRIGHT 2004 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED**

Severn Trent Laboratories, Inc.

STL Los Angeles • 1721 South Grand Avenue, Santa Ana, CA 92705-4808

Tel 714 258 8610 Fax 714 258 0921 • [www.stl-inc.com](http://www.stl-inc.com)

## TABLE OF CONTENTS

1. SCOPE AND APPLICATION	3
2. SUMMARY OF METHOD	3
3. DEFINITIONS	3
4. INTERFERENCES	4
5. SAFETY	5
6. EQUIPMENT AND SUPPLIES	6
7. REAGENTS AND STANDARDS	7
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE	7
9. QUALITY CONTROL	8
10. CALIBRATION AND STANDARDIZATION	10
11. PROCEDURE	11
12. DATA ANALYSIS AND CALCULATIONS	12
13. REPORTING	15
14. METHOD PERFORMANCE	16
15. POLLUTION PREVENTION & WASTE MANAGEMENT	16
16. REFERENCES	16
17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)	17

## 1. SCOPE AND APPLICATION

- 1.1. This method is applicable to the preparation of water samples and analysis of the prepared headspace to quantify part-per-million to percent levels of dissolved gases in the water sample. Although this method is specifically for the determination of Methane, Ethane, and Ethene, it may also be modified for Nitrogen, Oxygen, Carbon Dioxide, and Carbon Monoxide. This procedure is based on a non-promulgated method, EPA RSKSOP 175.
- 1.2. This method is restricted to use by or under the supervision of analysts experienced in sample preparation and in the use of gas chromatography and the interpretation of chromatograms

## 2. SUMMARY OF METHOD

- 2.1. A water sample is collected in a 40ml Voa vial, free of headspace, and capped using a Teflon-faced septum and crimp or screw-on cap, or equivalent, of the appropriate size to fit the bottle.
- 2.2. A headspace is generated in the laboratory by replacing 10 percent of the water sample with high purity Helium.
- 2.3. The sample bottle is agitated for 5 minutes and a sample is taken of the headspace and injected onto a gas chromatographic system where the gaseous components of interest are separated and detected by flame ionization detector or thermal conductivity detector.
- 2.4. By using Henry's Law, the headspace concentration of the gas, the bottle volume and the temperature of the sample, the concentration of the dissolved gas in the original water sample can be determined.

## 3. DEFINITIONS

- 3.1. Batch - An Analytical Batch is defined as a set of up to 20 client samples of the same matrix processed using the same procedures and reagents within the same time period. A batch must contain a Laboratory Control Sample (LCS), a Laboratory Control Sample Duplicate (LCSD), and a Method Bank. Using this method, the first CCV analysis will normally start a new analytical batch.
- 3.2. The Quality Control Batch must contain a Laboratory Control Sample (LCS), a Laboratory Control Sample Duplicate (LCSD) and a Method Bank.



- 3.3. Method Blank - A Method Blank consisting of all reagents added to the samples must be analyzed with each batch of samples. The Method Blank is used to identify any background interference or contamination of the analytical system that may lead to the reporting of elevated concentration levels or false positive data.
- 3.4. Laboratory Control Samples (LCS/LCSD) - Laboratory Control Samples are laboratory-generated samples used to monitor the laboratory's day-to-day performance. The LCS, spiked with a group of target compounds representative of the method analytes, is used to monitor the accuracy of the analytical process, independent of matrix effects. Ongoing monitoring of the LCS results provides evidence that the laboratory is performing the method within accepted QC guidelines for accuracy and precision.
- 3.5. Method Detection Limits - The Method Detection Limit (MDL) is defined as the minimum concentration of a substance that can be measured and reported with 99% confidence that the analyte concentration is greater than zero and is determined from seven replicate analysis of low level standards in a typical representative matrix.
- 3.6. Ambient Air (Room Air) - "Room air" is defined as containing 78.08% v/v Nitrogen, 20.946% v/v Oxygen, and 0.934% v/v Argon.

#### 4. INTERFERENCES

- 4.1. This system is relatively free of interferences due to specificity of the multiple columns used and the back-flush capabilities of the system. Compounds of interest are well separated and there are no major baseline upsets at expected retention times to complicate correct peak integration.
- 4.2. Ambient air (in the laboratory) contains approximately 400ppmv Carbon Dioxide and several ppmv of Methane. These concentrations should be considered variable and be taken into consideration whenever "room air" is used as a standard.
- 4.3. Argon naturally occurs in ambient air at 0.934 percent, and co-elutes with Oxygen.
- 4.4. The referenced methods note that water vapor can interfere with the chromatographic baseline. STL Los Angeles' choice of columns and back-flush technique after elution of target compounds minimizes this problem.
- 4.5. Hydrochloric Acid used for preservation of the samples has been shown to contribute to the Carbon Dioxide result, either through chemical interaction with the sample or from Carbon Dioxide contained in the HCL. Ideally, samples analyzed for Carbon Dioxide should be unpreserved.

- 4.6. This method is dependent on sample volumes, so final results are affected by samples containing large amounts of sediment or solid material.

## 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL associates.
- 5.2. Eye protection that satisfies ANSI Z87.1 (as per the Corporate Safety Manual), laboratory coat and appropriate gloves must be worn while samples, standards, solvents and reagents are being handled. Employees must abide by the policies and procedures in the Corporate Safety Manual, Lab Specific Addendum to the CSM, and this document.
- 5.3. Specific Safety Concerns and Requirements
- 5.3.1. MAPA Blue Grip gloves should be worn when handling the VOA vials.
- 5.3.2. Hydrogen is flammable
- 5.3.3. All safety issues involving compressed gas cylinders should be followed. All compressed gas cylinders must be securely fastened to a bench or wall.
- 5.3.4. Tedlar bags may be used for standard preparation. These must be handled with care and must not be over pressurized to prevent them from rupturing. Make sure the valves are closed tightly when not in use.
- 5.3.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit
- 5.4. Primary Materials Used
- 5.4.1. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDSs) maintained in the laboratory. Analysts and Technicians need to be familiar with the information contained in the applicable MSDS's before beginning work on this method. The following specific hazards are known:

Material (1)	Hazards	OSHA Exposure Limit (2)	Signs and symptoms of exposure/Unusual Hazards
Hydrochloric Acid	Corrosive Poison	5 ppm- Ceiling	Inhalation of vapors can cause coughing, choking, inflammation of the nose, throat, and upper respiratory tract, and in severe cases, pulmonary edema, circulatory failure, and death. Can cause redness, pain, and severe skin burns. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.

- 5.5. Exposure to chemicals must be maintained as low as reasonably achievable, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.6. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL associate. The situation must be reported immediately to a Supervisor

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Gas Chromatograph - capable of temperature programming for the oven and automatic multiple valve control, and equipped with thermal conductivity and flame ionization detectors (Varian 3400)
- 6.2. Custom valving and reduction catalyst housing (Valco Instruments and Lotus Consulting)
- 6.3. Tube oven for oxidation catalyst capable of 600 degrees C (Lindberg or equivalent)
- 6.4. Chromatographic grade stainless steel tubing and stainless steel plumbing fittings
- 6.5. Stainless steel chromatographic columns - Tenax 1/8" x 4', Chromasorb 106 1/8" x 6', Mol Sieve 5A 1/8 x 7' (Supelco, Inc.)
- 6.6. An assortment of gas tight syringes from 0.1 ml to 1.0 liter volume - for standard preparation and sample injection (Hamilton Syringe or equivalent)

- 6.7. Pressure regulators for carrier gas, flame ionization detectors and standards - 2 stage, stainless steel diaphragm
- 6.8. Tedlar Bags in 1 or 3 liter sizes (SKC or equivalent).
- 6.9. Automated data systems capable of archiving instruments runs (Varian d-654, PE Nelson Turbochrom, or equivalent).

## **7. REAGENTS AND STANDARDS**

- 7.1. High purity Helium for carrier gas, standard preparation, making dilutions, and purging blank water
- 7.2. Compressed Air for valve actuation
- 7.3. High purity Hydrogen for the Flame Ionization Detector
- 7.4. High purity Air for the Flame Ionization Detector
- 7.5. Prepared calibration standards. Standards are available at various concentrations commercially and are analytically certified by the supplier (Scott Specialty or equivalent).
- 7.6. Pure cylinders of Nitrogen, Carbon Dioxide, and Methane Room air can be used for the source of Oxygen/Argon, and also for Nitrogen at or below 78%.
- 7.7. Working standards. Standards are prepared at appropriate levels by making dilutions of the pure and prepared standards in Nitrogen using syringes and Tedlar Bags and/or pressure gauges in Summa or Silco canisters.
- 7.8. Reagent grade (distilled) water for blank preparation.
- 7.9. Reagent grade Hydrochloric Acid (1:1) for sample preservation (if required) and preserved blank preparation.

## **8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

- 8.1. Water samples should be collected in the field or prepared in the lab by placing the water in a 40ml VOA vial. There should be no headspace present in the vial.

- 8.2. Add the water down the side of the bottle so as not to agitate or contaminate the sample. Fill to the top and cap using a Teflon-faced rubber septum and appropriate size crimp or screw-on cap. Care should be taken to eliminate or reduce the formation of bubbles.
- 8.3. Field samples should be fixed with 1:1 Hydrochloric Acid to a pH less than 2 before they are capped. Pre-preserved Voa vials may be used.

**NOTE: DO NOT ADD ACID IF CARBON DIOXIDE IS TO BE DETERMINED,** since it may convert inorganic carbon to carbon dioxide. The use of unpreserved Voa vials is recommended.

- 8.4. STL Los Angeles observes a holding-time of 14 days from the date of collection whether there is preservative added or not.
- 8.5. Samples must be stored at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$

## 9. QUALITY CONTROL

### 9.1. Initial Demonstration of Method Capability

- 9.1.1. Method Detection Limit (MDL) - A MDL must be determined prior to the analysis of any samples. The MDL is determined using seven replicates that have been spiked with the analyte of interest. These spiked replicates must be carried through the entire analytical procedure. MDL's must be determined and verified on an annual basis. The spike level must be between the calculated MDL and 10x the MDL to be valid. The result of the MDL determination must be below the reporting limit.
- 9.1.2. Demonstration of Capability (DOC) Study - Replicate LCS analysis. At the initiation of this method, four replicates of the LCS must be analyzed and evaluated to determine method accuracy and precision and the middle portion of the linear dynamic range. Results of the initial demonstration study must be acceptable before the analysis of samples may proceed. Annual DOC's are done to assure continued proficiency in the method for an analyst.
- 9.1.3. The department manager is responsible for ensuring that this procedure is performed by an analyst who has been properly trained and has the required experience.

## 9.2. Control Limits

- 9.2.1. When available, in-house historical control limits must be used for Laboratory Control Samples (LCS). These limits must be determined annually at a minimum. The recovery limits are determined as the mean recovery,  $\pm 3$  standard deviations for the LCS. Default limits are, at a minimum, 50% to 150% for recovery and 20% for the RPD between the LCS and the LCSD. In-house control limits may be tighter than the default limits.

## 9.3. Quality Control Batch

- 9.3.1. The QC batch is a set of up to 20 client samples of the same matrix processed using the same procedures and reagents within the same time period. The batch also must contain an LCS/LCSD pair and a method blank. Laboratory generated QC samples or instrument QC (MB, LCS/LCSD pair, calibration standards, ICV's and CCV's) do not count towards the maximum 20 samples in a batch. Field QC samples are included in the batch count.

## 9.4. Method Blanks

- 9.4.1. For each batch of up to 20 client samples, analyze a method blank consisting of Helium-purged reagent grade water with headspace prepared identically to the samples. The method blank is analyzed after the calibration standards, normally before any samples.
- Target analytes in the method blank must not exceed the reporting limit.
- 9.4.2. If there is no target analyte greater than the reporting limit (RL) in the samples associated with an unacceptable MB, the data may be reported with qualifiers. Such action should be done in consultation with the client and a NCM filed.

## 9.5. Laboratory Control Samples (LCS)

- 9.5.1. For each batch of samples, analyze an LCS/LCSD pair. The pair is analyzed after the calibration standard, and normally before any samples. The LCS contains all the analytes of interest in the mid-range of the calibration. If any analyte is outside established control limits, the system is out of control and corrective action must occur. Corrective action will normally include either the re-prep and/or reanalysis of the LCS, and sample reanalysis if there are samples that apply to the failed QC batch.

9.5.2. If the batch is not re-analyzed, the reasons for accepting the batch must be clearly presented in the project records, NCM and the report narrative.

- If reanalysis of the batch is not possible the LCS is reported, all associated samples are flagged, and appropriate comments are made in a narrative to provide further documentation.
- If any analyte in the LCS exceeds the upper control limit and that analyte is not detected in the sample, then no corrective action is required. Document in an NCM.

## **10. CALIBRATION AND STANDARDIZATION**

### **10.1. Standard Preparation**

10.1.1. Standards are prepared in Tedlar Bags and/or Summa or Silco canisters using certified standards and static dilution techniques. Nitrogen is the diluent gas of choice.

### **10.2. Calibration and Verification**

10.2.1. An initial calibration curve (ICAL) consisting of 5 points is run to determine the linear working range of the system for each compound. The low point should be at or below the reporting limit. A relative standard deviation (RSD) is calculated for each target analyte response factor using the calculation in Section 12. Each ICAL must have a %RSD less than 25% for each compound prior to analysis of samples.

10.2.2. Continuing Calibration Verification (CCV) standards for all target analytes are analyzed at a frequency of every 24 hours or 20 reportable continuously injected client samples, whichever is more frequent. Standards are prepared in Tedlar Bags using static dilution techniques and appropriate levels of gas standards certified by an outside supplier. Room air may be used for Nitrogen and Oxygen/Argon calibration. All target analyte RF's must be within 25%D of the initial calibration average RF's. Corrective action includes re-preparation and re-analysis of the standard, instrument maintenance and generation of a new initial calibration.

10.2.3. Once the above criteria have been met, sample analysis may begin. A closing CCV shall be injected if the sample batch cannot be finalized within the continuous injection period as defined above.

10.2.4. If sample analysis must be halted for more than an eight (8) hour period then an opening CCV must be analyzed to ensure that instrumentation conditions have remained stable.

10.2.4.1. If the opening CCV passes then sample analysis may proceed.

10.2.4.2. If the initial CCV fails, re-inject a second CCV. If the second CCV passes, analysis may proceed. If it fails, conduct maintenance and re-calibrate the instrument. Document the failure in the run log.

10.2.4.3. If the CCV's indicate valid instrument conditions, then the batch QC can be used for the data generated on two separate days.

## **11. PROCEDURE**

11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and/or QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

### **11.1. Sample Analysis Procedure**

11.1.1. GC parameters, nominal sample and QC volume, and sample temperature conditions are listed in Section 17.2.

11.1.2. All samples must be analyzed as part of a batch.

11.1.3. It is not necessary to reanalyze batch QC with reanalysis of samples. However, any reanalysis of samples must be associated with valid instrument QC.

11.1.4. Remove Voa's from refrigerator and allow to come to room temperature. Voa's should be at ambient temperature before they are extracted. Elevated sample temperatures are to be avoided, as they will affect sample results.

11.1.5. Fill a clean gas-tight syringe with a volume of contaminant-free Helium equal to 10% of the volume of the sample container. A standard 40ml VOA vial has been found to contain 43.5ml. Typically 4.4ml of sample is removed and replaced with Helium. Insert the needle of this syringe and another, empty gas-tight syringe through the septum of the sample container.



- 11.1.6. Slowly inject the Helium into the sample container while removing an equivalent volume of water with the other syringe.
- 11.1.7. Remove both syringes and agitate the bottle on a vortex for 5 minutes. **Caution: excessive handling of the sample should be avoided, as this will raise the temperature of the sample.**
- 11.1.7.1. If a vortex or automated shaker is not available, manual shaking can be used. The analyst must limit hand contact with the Voa, such as, only holding the ends of the vial with the fingertips.
- 11.1.8. After agitating the sample, use a 500 $\mu$ L gas-tight syringe to remove a 300 $\mu$ L aliquot of headspace sample. Inject the syringe's contents into the GC.

## 12. DATA ANALYSIS AND CALCULATIONS

### 12.1. Qualitative Identification

- 12.1.1. Retention time windows are set at +/- 0.5 minutes of the expected retention time as referenced to the CCV values. The reference method for RSK-175 does not address this issue.

### 12.2. Quantitative Analysis

- 12.2.1. The target analyte quantity is calculated by dividing the peak area of the target analyte by the average response factor of the ICAL, and applying any dilution factor generated, taking into account static dilutions and/or injection volumes.
- 12.2.2. The headspace concentration is calibrated and calculated as volume/volume (either %v/v or ppmv/v) and then used in the Henry's Law equation in section 12.4 to determine the concentration in water.

### 12.3. Calculations

- 12.3.1. The Turbochrom data system automatically quantitates the sample results based on a predetermined sample size. The results are in ppmv/v from the FID and %v/v from the TCD. The sample size is nominally 0.30 ml.

12.3.2. Calculation for Percent Relative Standard Deviation (%RSD):

$$\%RSD = \text{Std Dev of RF's} / \text{Mean of RF's} \times 100$$

12.3.3. Calculation for Percent Difference (%D):

$$\% D = [(RF \text{ cpd CCV} - \text{ave. RF cpd ICAL}) / \text{ave. RF cpd ICAL}] \times 100$$

12.3.4. Calculation for Determining Concentration of Compounds in Headspace

$$\text{Conc. Cpd} = \text{Area Cpd} / \text{ave. RF Cpd} \times \text{DF}$$

12.3.5. Calculation for Percent Recovery (%Rec):

$$\% \text{ Rec} = (\text{Amount Cpd recovered} / \text{Amount Cpd Spiked}) \times 100$$

12.3.6. Calculation for Relative Percent Difference (RPD):

$$\text{RPD} = |\text{Value A} - \text{Value B}| / \text{Average of Values} \times 100$$

12.3.7. Calculation of Decimal Fraction  $p_g$

$$p_g = \text{ppmv/v} \times 10^{-6} = \%v/v \times 10^{-2}$$

12.4. The following parameters are needed to perform Henry's Law calculations:

$p_g$  = headspace concentration of analyte in decimal fraction, eg. 10ppm = 0.00001

H = Henry's Law Constant

ST = temperature of the sample (assumed to be standard temp. 298 °K)

bv = volume of sample bottle

hv = headspace volume

1)  $p_g$  = result in ppmv divided by 1000000 or result in %v/v divided by 100

2) equilibrium mole fraction of the dissolved gas,  $x_g = p_g / H$

where H = Henry's Law Constant for the gas

3) Let  $n_g$  = moles analyte and  $n_w$  = moles water

$$\text{Then } x_g = n_g / (n_g + n_w) \text{ and } n_g = x_g (n_g + n_w)$$

If,  $n_g * x_g \ll n_g$ ,

Then,  $n_g = x_g * n_w$  or  $n_g = n_w (p_g/H)$

Therefore;  $n_g / V = n_w / V (p_g/H)$

4) One liter of water is 55.5 g-moles,

$$n_g / V = 55.5 / L (p_g/H)$$

5) Saturation concentration of the gas,

$$C = (n_g/V) (MW) (1000 \text{ mg/g})$$

where MW = molecular weight of the analyte

6) Density calculation

$$p = (MW) / (22.4 \text{ l/mole}) (ST \text{ in } ^\circ K / 273 ^\circ K)$$

where: p = density

ST = sample temperature

7)  $v = (bv_{mls} - hv_{mls}) (1L/1000mls)$

where: bv = bottle volume

hv = headspace volume

**Then:**

$$8) A_h = hv_{mls} * p_g$$

where  $A_h$  = ml of analyte in headspace

then liquid phase analyte ( $A_l$ ) is

$$A_l = (A_h/(v))(p_g)(1000 \text{ mg/g})(1L/1000ml)$$

**Then: TC =  $A_l$  + C**

Where: TC = Total Concentration of analyte in the original sample

$A_l$  = liquid phase analyte

C = saturation concentration

The result will be in units of milligrams of gas per liter of water

### 12.5. Example Calculation:

$p_g = 0.0018$  (1800 ppmv/v) methane

Henry's Law Constant,  $H = 4.13 \times 10^4$  for Methane

Sample Temperature 25 C (298 K)

Bottle Volume 60 ml

Headspace volume 6 ml

Using equation 2,  $x_g = 0.0018 / 4.13 \times 10^4$  or  $4.358 \times 10^{-8}$  mole  $\text{CH}_4$

Using equation 4 and the value above,

$n_g/V = (55.5)(4.358 \times 10^{-8})$  or  $2.42 \times 10^{-6}$  moles  $\text{CH}_4$  per liter of water.

$p = (16 \text{ g/mole}) / ((22.4 \text{ l/mole})(298/273)) = 0.654 \text{ g } \text{CH}_4 / \text{liter } \text{H}_2\text{O}$

$b_v = 60 \text{ ml}$  and  $h_v = 6 \text{ ml}$ ,  $v = (60 \text{ ml} - 6.0 \text{ ml})(1 \text{ L} / 1000 \text{ ml}) = 0.054 \text{ L}$

$A_h = 6 \text{ ml}$ ,  $0.0018 = 0.0108 \text{ ml } \text{CH}_4$

$A_l = (0.0108 \text{ ml} / 0.054 \text{ L})(0.654 \text{ g/L})(1 \text{ L} / 1000 \text{ ml})(1000 \text{ mg/g})$

$A_l = 0.1308 \text{ mg } \text{CH}_4 / \text{liter } \text{H}_2\text{O}$

Then  $\text{TC} = A_l + C = 0.038 \text{ mg/L} + 0.131 \text{ mg/L} = 0.169 \text{ mg} / \text{liter } \text{H}_2\text{O}$

## 13. REPORTING

13.1. Reporting Limits can be found in Section 18.5

13.2. Reporting Results

13.2.1. Estimates of uncertainty are based on historical control limits for the LCS.  
 These limits can be provided when requested.

13.1. No conversion of the analytical results to the standard conditions is made.

#### 14. METHOD PERFORMANCE

- 14.1. Method performance is controlled through the measurement of accuracy and precision. Analysis of a LCS and a LCSD are performed to measure both accuracy and precision. The control limits established for the LCS and LCSD are used to maintain method performance within a well-defined set of criteria.

#### 15. WASTE MANAGEMENT AND POLLUTION PREVENTION

- 15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by, and this method is set up in accordance with, section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention".
- 15.2. Waste Streams produced by this method.
- 15.2.1. Expired Standards – these are identified as expired, stored under manufacturer's recommended conditions and then lab packed for disposal.
- 15.2.2. Expired gas cylinders are emptied and returned to the manufacturer.
- 15.2.3. Sample Voa vials – maintained under refrigeration and then transferred to 55 gallon drums in the 90 day area.

#### 16. REFERENCES

- 16.1. Newell, B.S., RSKSOP-175, Revision Number 0, August 1994, R.S. Kerr Environmental Research Laboratory of the United States Environmental Protection Agency.
- 16.2. Kampbell, D.H., J.T. Wilson, S.A. Vandegrift, Dissolved Oxygen and Methane in Water by GC Headspace Equilibration Technique, International Journal of Environmental Analytical Chemistry, Volume 36, pp. 249-257, 1991.
- 16.3. Vandegrift, S.A., RSKSOP-114, Revision Number 0, January 1991.
- 16.4. Newell, B.S., RSKSOP-147, Revision Number 0, August 1993.
- 16.5. Perry, J.H., *Chemical Engineer's Handbook*, (McGraw-Hill, NY, 1978), 5th ed.

16.5.1. Varian manual 3400 GC

16.5.2. PE Nelson Turbochrom manual

**17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)**

17.1. Table 1 - Variables used in normal calculations

17.2. Table 2 - GC Conditions: Varian 3400

17.3. Table 3 - CCV and LCS Nominal Concentrations

17.4. Table 4 - Nominal ICAL Concentrations

17.5. Table 5 - Reporting Limits and Default Control Criteria

**Table 1 - Variables used in normal calculations**

Variable	Value	Units
Temperature	298	°K
Volume 40 ml VOA vial (bv)	43.5	ml
Volume headspace generated (hv)	4.4	ml
MW Methane	16	g/mole
MW Ethane	30	g/mole
MW Ethylene	28	g/mole
MW Carbon dioxide	44	g/mole
MW Oxygen	32	g/mole
MW Nitrogen	28	g/mole
MW Carbon Monoxide (25 deg C)	28	g/mole
Henry's Const. Methane (25 deg C)	4.13e+4	none
Henry's Const. Ethane (25 deg C)	3.02e+4	none
Henry's Const. Ethylene (25 deg C)	1.14e+4	none
Henry's Const. Carbon dioxide (25 deg C)	1.64e+3	none
Henry's Const. Oxygen (25 deg C)	4.38e+4	none
Henry's Const Nitrogen (25 deg C)	8.65e+4	none
Henry's Const. Carbon Monoxide (25 deg C)	5.80e+4	none

8771088

SOP No. COI-GC-005  
Revision No. 4  
Revision Date: 05/21/04  
Effective Date: 06/04/04  
Page: 18 of 19

**Table 2 - GC Conditions: Varian 3400**

Initial Column Temp	65	Det B range (init)	12
Hold min	5.10	Det B autozero	yes
Prog 1 Temp	100	Det B range at time 0.01 min	11 no A/Z
Prog 1 rate deg/min	35	Det B range at time 5.0 min	12 no A/Z
Prog 1 hold min	11.9	A/S vial mode	off
Inj temp	100	Wait for ready	yes
Aux. Temp. (FID and reduction catalyst)	380	Initial Relays	-1-2-3-4
Detector temp (TCD)	150	Relays 0.01 min	+2
Det A range	0.50	Relays 2.85 min (may vary)	+3
Det A autozero	yes	Relays 8.00 min	-2
Det A filament temp	200	Relays 12.00 min	-3
Det A polarity positive	no	External Ox catalyst oven	550
Det A time program	no		

**Table 3 - CCV and LCS Nominal Concentrations**

Compound	CCV Concentration	LCS Concentration
Carbon Dioxide (TCD)	2.0 %v/v	1.0% v/v
Oxygen (TCD)	5.0% v/v	2.188 %v/v (10 % v/v air)
Nitrogen (TCD)	20.0% v/v	7.808 %v/v (10 % v/v air)
Methane (TCD)	2.0 %v/v	--
Carbon Dioxide (FID)	1,000 ppm v/v	10,000 ppm v/v
Ethylene (FID)	100 ppm v/v	500 ppm v/v
Ethane (FID)	100 ppm v/v	500 ppm v/v
Methane (FID)	100 ppm v/v	500 ppm v/v
Carbon Monoxide (FID)	100 ppm v/v	500 ppm v/v

**Table 4 - Nominal ICAL Concentrations**

Compound	Level 1	Level 2	Level 3	Level 4	Level 5	Level 6
Carbon Dioxide (TCD) %v/v	0.10	0.50	1.0	2.0	10.0	50
Oxygen (TCD) %v/v	0.219	1.09	2.19	10.9	21.88	--
Nitrogen (TCD) %v/v	0.780	3.90	7.81	39.0	78.08	--
Methane (TCD) %v/v	0.10	0.50	1.0	5.0	10.0	50
Carbon Dioxide (FID) ppmv/v	100	500	1000	5000	10000	--
Ethylene (FID) ppm v/v	5	10	100	500	1000	--
Ethane (FID) ppm v/v	5	10	100	500	1000	--
Methane (FID) ppm v/v	10	100	500	1000	5000	10000
Carbon Monoxide (FID) ppm v/v	10	100	500	1000	5000	20000

**Table 5 - Reporting Limits and Default Control Criteria**

Compound	Reporting Limits	Control Limits % **	RPD %
Carbon Dioxide (TCD)	1.7 mg/l	50 - 150	20
Oxygen (TCD)	4.0 mg/l	50 - 150	20
Nitrogen (TCD)	10 mg/l	50 - 150	20
Methane (TCD)	0.10 mg/l	50 - 150	20
Carbon Dioxide (FID)	0.17 mg/l	50 - 150	20
Ethylene (FID)	0.0010 mg/l	50 - 150	20
Ethane (FID)	0.0020 mg/l	50 - 150	20
Methane (FID)	0.0010 mg/l	50 - 150	20
Carbon Monoxide (FID)	0.0010 mg/l	50 - 150	20

\*\* Default limits, actual control limits may be tighter



8771090

Controlled Copy

Copy No. \_\_\_\_\_

Implementation Date 1/24/03SOP No. NC-WC-0017Revision No. 2.2Revision Date: 01/24/03

Page 1 of 15

**STL North Canton  
STANDARD OPERATING PROCEDURE**

**TITLE: TOTAL ORGANIC CARBON (TOC)**

**(SUPERSEDES: REVISION 2.1, DATED 08/17/00)**

Prepared by:	<u>Alquita Colon</u>	<u>1/24/03</u> Date
Reviewed by:	<u>Angela Serra</u> Technology Specialist	<u>1/24/03</u> Date
Approved by:	<u>Beth Lambert</u> Quality Assurance Manager	<u>1/24/03</u> Date
Approved by:	<u>Debrah S. Buckle</u> Environmental Health and Safety Coordinator	<u>1/24/03</u> Date
Approved by:	<u>Golden Lee</u> Laboratory Director	<u>1/24/03</u> Date

**Proprietary Information Statement:**

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

---

## TABLE OF CONTENTS

<i>1. SCOPE AND APPLICATION .....</i>	<i>3</i>
<i>2. SUMMARY OF METHOD .....</i>	<i>3</i>
<i>3. DEFINITIONS.....</i>	<i>3</i>
<i>4. INTERFERENCES.....</i>	<i>3</i>
<i>5. SAFETY.....</i>	<i>3</i>
<i>6. EQUIPMENT AND SUPPLIES .....</i>	<i>4</i>
<i>7. REAGENTS AND STANDARDS.....</i>	<i>5</i>
<i>8. SAMPLE COLLECTION, PRESENTATION AND STORAGE.....</i>	<i>6</i>
<i>9. QUALITY CONTROL.....</i>	<i>6</i>
<i>10. CALIBRATION AND STANDARDIZATION.....</i>	<i>9</i>
<i>11. PROCEDURE.....</i>	<i>11</i>
<i>12. DATA ANALYSIS AND CALCULATIONS .....</i>	<i>12</i>
<i>13. METHOD PERFORMANCE.....</i>	<i>13</i>
<i>14. POLLUTION PREVENTION .....</i>	<i>14</i>
<i>15. WASTE MANAGEMENT.....</i>	<i>14</i>
<i>16. REFERENCES.....</i>	<i>14</i>
<i>17. MISCELLANEOUS (TABLES, APPENDICES, ETC...).....</i>	<i>15</i>

**1. SCOPE AND APPLICATION**

- 1.1. This method is applicable to the determination of Total Organic Carbon (TOC) in waters and similar matrices. It is based on SW846 Method 9060 and EPA Method 415.1. The working linear range is instrument dependent at 1 mg/L to 50 mg/L with a reporting limit of 1 mg/L.
- 1.2. QuantIMS method codes are DA (415.1) and FM (9060).
- 1.3. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary.

**2. SUMMARY OF METHOD**

- 2.1. Organic Carbon is converted to carbon dioxide (CO<sub>2</sub>) using chemical oxidation. The CO<sub>2</sub> is then measured by an infrared detector.

**3. DEFINITIONS**

- 3.1. Refer to the glossary in the Laboratory Quality Manual (LQM), latest version.

**4. INTERFERENCES**

- 4.1. Contaminants in solvents, reagents, glassware, and other processing apparatus that lead to discrete artifacts may cause Method interferences. All of these materials must be routinely demonstrated to be free from interferences under conditions of the analysis by running laboratory method blanks as described in the Quality Control section. Specific selection of reagents may be required to avoid introduction of contaminants.
- 4.2. Carbonate and bicarbonate interfere but are eliminated by the acidification and purging step of the instrument.

**5. SAFETY**

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL associates.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the

Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:

- 5.3.1. The following materials are known to be **corrosive: Phosphoric Acid, Sulfuric Acid.**
- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.6. It is recommended that neat standards be purchased only as a last resort. The preparation of standards from neat materials and reagents {as well as glassware cleaning procedures that involved solvents such as methylene chloride} should be conducted in a fume hood with the sash closed as far as the operations will permit.
- 5.7. Standards in solution may be diluted in the open laboratory when syringes and the like are utilized.
- 5.8. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL associate. The situation must be reported **immediately** to a laboratory supervisor.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. O-I Corporation Model 1010 TOC Analyzer with 1051 vial multisampler and Printer
- 6.2. Nitrogen Gas and Regulator
- 6.3. Volumetric flasks: Various sizes
- 6.4. Volumetric pipettes: Various sizes
- 6.5. Vials: 40 mL glass
- 6.6. Graduated cylinders: Various sizes
- 6.7. pH Strips
- 6.8. Whatman filter #4

6.9. Top loading balance: capable of accurately weighing  $\pm 0.01$  g

## 7. REAGENTS AND STANDARDS

### 7.1. Reagents

7.1.1. Sodium Persulfate: Reagent Grade

7.1.2. Sodium Persulfate Solution: Add 200 g sodium persulfate ( $\text{Na}_2\text{S}_2\text{O}_8$ ) to a 1 liter volumetric flask and dilute to volume with reagent water.

7.1.3. Phosphoric Acid, concentrated: Reagent Grade

7.1.4. Phosphoric Acid Solution: Carefully add 59 mL concentrated phosphoric acid ( $\text{H}_2\text{PO}_4$ ) to 900 mL of reagent water in a 1 liter volumetric flask. Dilute to volume with reagent water.

7.1.5. Sulfuric Acid, concentrated: Reagent Grade

### 7.2 Standards

7.2.1. All standards should be prepared in volumetric flasks, using volumetric pipettes, and diluted to volume with reagent water.

7.2.2. TOC Stock Standard

7.2.3. Primary and secondary sources are needed.

7.2.3.1. TOC 1000 mg/L

7.2.3.1.1. Dilute 1.06 g KHP (potassium acid phthalate) to volume in a 500 mL volumetric flask. A commercially prepared solution may also be used.

7.2.3.2. Prepare every six months.

7.2.4. TOC Calibration Standards

## TOTAL ORGANIC CARBON (TOC)

SOP No. NC-WC-0017Revision No. 2.2Revision Date: 01/24/03Page 6 of 15

---

- 7.2.5. Prepare the following standards from the primary stock standard described in Section 7.2.2.2.1.

Concentration (mg/L)	Volume (mL)	Stock Concentration (mg/L)	Final Volume (mL)
50	5	1000	100
25 (MS/MSD) (LCS)	12.5	1000	500
10	1	1000	100
1	0.1	1000	100

- 7.2.6. TOC Verification Standard (LCS)

7.2.6.1. A commercially prepared solution is used.

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

8.1 Samples are preserved to a pH <2 with sulfuric acid (H<sub>2</sub>SO<sub>4</sub>) or hydrochloric acid (HCl) and stored in plastic or glass containers at 4°C±2°C.

8.2 The holding time is twenty-eight days from sampling to analysis.

## 9. QUALITY CONTROL

### 9.1. Batch Definition

9.1.1. A batch is a group of no greater than 20 samples excluding QC samples (LCS, Method Blank, MS, MSD) which are processed similarly, with respect to the procedure. All sample setups must be initiated within a 24 hour period from the initial preparation or extraction and without interruption of the process. All samples within the batch must be treated with the same lots of reagents and the same processes.

### 9.2. Method Blank

9.2.1. One method blank (MB) must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including

preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting limit.

9.2.2. A method blank consisting of 40 mL of reagent water and all reagents added to the samples must be prepared and analyzed with each batch of samples. The method blank is used to identify any background interference or contamination of the analytical system, which may lead to the reporting of elevated concentration levels or false positive data.

9.2.3. Corrective Action for Blanks

9.2.3.1. If the analyte level in the method blank exceeds the reporting limit for the analytes of interest in the sample, all associated samples are reprepared and reanalyzed. If this is not possible due to limited sample quantity or other considerations, the corresponding sample data **must be taken in consultation with the client and must be addressed in the project narrative.**

9.2.3.2. If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. **Such action must be taken in consultation with the client and must be addressed in the project narrative.**

9.3. Laboratory Control Sample (LCS)

9.3.1. Laboratory Control Samples are well characterized; laboratory generated samples used to monitor the laboratory's day to day performance of routine analytical methods. The LCS is used to monitor the accuracy of the analytical process, independent of matrix effects. Ongoing monitoring of the LCS results provides evidence that the laboratory is performing the method within accepted QC guidelines for accuracy and precision.

9.3.2. A purchased LCS must be analyzed with each batch of samples.

9.3.3. Corrective Action for LCS

9.3.3.1. If any analyte is outside established control limits the system is out of control and corrective action must occur.

9.3.3.2. The only exception is that if the LCS recoveries are biased high and the associated sample is ND for the parameter(s) of interest, the batch is acceptable. **This must be addressed in the project narrative.**

9.3.3.3. Corrective action will be reparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.

#### 9.4. Matrix Spike/Matrix Spike Duplicate (MS/MSD)

9.4.1. One MS/MSD pair must be processed for each batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Samples identified as field blanks cannot be used for MS/MSD analysis.

9.4.2. An MS/MSD consisting of 20 mL of sample and 20 mL of the 25 mg/L standard will be analyzed with each analytical batch of samples.

##### 9.4.3. Corrective action for MS/MSDs

9.4.3.1. If the analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include reparation and reanalysis of the batch.

9.4.3.2. If the native analyte concentration in the MS/MSD exceeds 4x the spike level for that analyte, the recovery data is reported as "amount" MSB. The Exception Code is changed to NC. The following two footnotes will appear on the report page "NC The recovery and/or RPD were not calculated." "MSB The recovery and RPD were not calculated because the sample amount was greater than four times the spike amount."

9.4.3.3. If an MS/MSD is not possible due to limited sample volume then a laboratory control sample duplicate (LCSD) should be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.



## 9.5. Control Limits

9.5.1. Control limits are established by the laboratory as described in SOP, NC-QA-0018.

9.5.2. Laboratory control limits are internally generated and updated periodically unless method specified. Control limits are listed in the Laboratory Quality Manual (LQM) and the latest version is easily accessible via the LIMs (QC Browser program).

## 9.6. Method Detection Limits (MDLs) and MDL Checks

9.6.1. MDLs and MDL Checks are established by the laboratory as described in SOP, NC-QA-0021.

9.6.2. MDLs are listed in the Laboratory Quality Manual (LQM) and the latest version is easily accessible via the LIMs (QC Browser program).

## 9.7. Nonconformance and Corrective Action

9.7.1. Any deviations from QC procedures must be documented as a nonconformance, with applicable cause and corrective action approved by the facility QA Manager.

# 10. CALIBRATION AND STANDARDIZATION

## 10.1. Recommended Initial Setup

Contant Settings		
STD Mass	=	6.76 ug C
Sample Vol	=	2.0 mL
Acid Vol	=	4 x 100 uL
Oxidant Vol	=	10 x 100 uL

10.1.1 Adjust the nitrogen to 30 psi using the flow valve on the tank. The gauge should always be set at 30 psi when the instrument is not in use.

- 10.1.2 Remove the reagent bottles and fill with appropriate reagents (phosphoric acid solution and sodium persulfate). Do not fill bottles completely full; leave a small amount of air space. Loosely reconnect caps (and tubing), and replace the bottles into the instrument.
- 10.1.3 Blank and calibrate the instrument when CCVs and/or CCBs fail to meet acceptance criteria or when other problems are encountered.
  - 10.1.3.1. Choose "Calibration". Start a new file with the current date.
  - 10.1.3.2. Choose "Sequences" from the "databases" menu option, open up the calibration template.
  - 10.1.3.3. Confirm all information. If blanking is required, it is best if done before calibration. Enter the desired number of blanks (no less than five) in the "reagent blanks before" field.
  - 10.1.3.4. Analyze an ICV/ICB
  - 10.1.3.5. Save the file using the current date as the filename.
  - 10.1.3.6. Update data file information.
    - 10.1.3.6.1. Choose the setup menu option, go into win TOC output, change the log fill name and prefix counter.
  - 10.1.3.7. When analysis is complete, print the run from "/utilities/view run log"
  - 10.1.3.8. Evaluate the data. The correlation coefficient of the original curve must be  $\geq 0.995$  or recalibration is required.

## 10.2. Continuing Calibration

- 10.6.1. The run is checked at the beginning, after every ten samples, and at the end of the run of the same species using a midrange CCV made from a primary source (Section 7.2.4. or 7.2.7.) to verify continued linearity. A CCV cannot vary from the original curve by more than  $\pm 10\%$ , or recalibration is required.
- 10.6.2. System cleanliness is checked every ten samples and at the end of the run using Continuing Calibration Blank (CCB). A CCB cannot contain the analyte of interest above the reporting limit, or recalibration is required.

## 11. PROCEDURE

11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

11.2. Sample Preparation Procedure

11.2.1. If excess particulate matter exists, filter an aliquot of sample through a Whatman #4 filter into a TOC vial or decant.

11.3. Preparation Documentation

11.3.1. Record any sample preparation on the analytical logsheet.

11.4. Analytical Documentation

11.4.1 Record all analytical information in the analytical logbook/logsheet, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.4.2 All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. The supervisor or designee reviews logbooks.

11.4.3 Documentation, such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs, is available for each data file.

11.4.4 Sample results and associated QC are entered into the Laboratory Information Management System (LIMS) after final technical review.

11.4.5 Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

## 12 DATA ANALYSIS AND CALCULATIONS

12.1. Sample Analysis Procedure

12.1.1. Type a run protocol sequence into the computer using the run template, if desired. Update the data file information in the setup/win TOC output.

## TOTAL ORGANIC CARBON (TOC)

SOP No. NC-WC-0017Revision No. 2.2Revision Date: 01/24/03

Page 12 of 15

12.1.2. For method 9060, quadruplicate analysis is required. If quadruplicate reporting is requested, each of four results are reported. Replicate analysis should be taken from separate vials, if available. If only one reportable result is requested per sample, the four results should be taken from one vial, and the average of the four results are reported.

12.1.3. For method 415.1 only one analysis is required. The single analysis is reported directly from the instrument printout.

12.1.4. All samples and standards should be poured into 40 mL vials. Samples received in vials can be run in those containers, provided there is not an excess of solids.

12.1.5. Be sure the samples are loaded on the sampler, the first one positioned under the needle.

12.1.6. Click the "Start" button.

12.1.7. Samples that fall outside the linear range (>50 mg/L) of the instrument must be diluted and reanalyzed.

12.1.7.1. Samples following a high sample should be re-analyzed if carryover is a concern.

12.1.8. Print the run from "utilities/view run log".

12.1.9. When analysis is complete, properly dispose of or put away samples and standards.

## 12.2. Calculations for 9060 only

12.2.1. Total Organic Carbon, mg/L = Average Instrument Value x Dilution

*Where:*

*TOC, mg/L = average of the 4 instrument values\* x dilution, calculated (without dilution) by the instrument.*

$$12.2.2. \quad LCS \% Recovery = \frac{\text{Instrument values}^*}{\text{True Value}} \times 100$$

12.2.3. MS/MSD % recovery

$$\left( \frac{(\text{Instrument values} * \text{MS or MSD}) - (\text{Avg sample instrument value} * 2)}{12.5} \right) \times 100$$

*\*One of the values may be judged erroneous and disregarded if three of the four are consistent. If no consistency can be found in the four values, the sample must be rerun.*

### 13 METHOD PERFORMANCE

13.1. Each laboratory must have initial demonstration of performance data on file and corresponding method detection limit files.

13.2. Training Qualifications:

13.2.1. The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

13.2.2. Method validation information (where applicable) in the form of laboratory demonstrations of capabilities is maintained for this method in the laboratory QA files.

### 14. POLLUTION PREVENTION

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

### 15. WASTE MANAGEMENT

15.1. Acid waste must be collected in clearly labeled acid waste containers.

15.2. Solid materials (gloves, soiled paper products, etc.) are placed in the solid debris container. Do not put liquids in the solid waste container.

15.3. Refer to the Laboratory Sample and Waste Disposal plan.

15.4. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by an annual refresher training.

15.5. Solvent waste must be disposed of in clearly labeled waste cans.

## 16. REFERENCES

### 16.1. References

16.1.1. SW846, Test Methods for Evaluating Solid Waste, Third Edition, Total Organic Carbon, Method 9060.

16.1.2. EPA 600, Methods for Chemical Analysis of Water and Wastes, Organic Carbon, Method 415.1.

16.1.3. Corporate Quality Management Plan (QMP), current version.

16.1.4. STL Laboratory Quality Manual Plan (LQMP), current version.

### 16.2. Associated SOPs and Policies, latest version

16.2.1. QA Policy, QA-003

16.2.2. Glassware Washing, NC-QA-0014

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018

16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021

16.2.5. Navy/Army SOP, NC-QA-0016

## 17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)

### 17.1. Reporting limits

17.1.1. The lower reporting limit is 1 mg/L

17.1.2. If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.

### 17.2. Troubleshooting guide

8771104

TOTAL ORGANIC CARBON (TOC)

SOP No. NC-WC-0017

Revision No. 2.2

Revision Date: 01/24/03

Page 15 of 15

---

17.2.1. See the manufacturer's instructions for an instrument troubleshooting guide and maintenance requirements.

17.3.Method deviations

17.3.1. A blender is not used to homogenize samples.

17.3.2. For Method 9060, the calibration must be verified with an independently prepared check standard every 15 samples. The laboratory is verifying the calibration every 10 samples.

Controlled Copy  
Copy No. \_\_\_\_\_

Implementation Date: 4/10/01

8771105  
SOP No. NC-WC-0018  
Revision No. 2.1  
Revision Date: 03/30/01  
Page 1 of 12

**STL NORTH CANTON STANDARD OPERATING PROCEDURE**

**TITLE: TOTAL ORGANIC CARBON (TOC) ANALYSIS FOR NON-WATERS**

**(SUPERSEDES: REVISION 2.0 DATED 02/27/97)**

Prepared by:	<u>Deborah R Marcum</u>	<u>4-4-01</u>	Date
Approved by:	<u>Kurt C. All</u> Technology Specialist	<u>4/02/01</u>	Date
Approved by:	<u>Carol Ann Fink</u> Quality Assurance Manager	<u>4/6/01</u>	Date
Approved by:	<u>[Signature]</u> Environmental Health and Safety Coordinator	<u>4-9-01</u>	Date
Approved by:	<u>Christy K. Goolsby</u> Laboratory Director	<u>4-9-01</u>	Date

**Proprietary Information Statement:**

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.



**TABLE OF CONTENTS**

<b>1. SCOPE AND APPLICATION .....</b>	<b>3</b>
<b>2. SUMMARY OF METHOD .....</b>	<b>3</b>
<b>3. DEFINITIONS.....</b>	<b>3</b>
<b>4. INTERFERENCES .....</b>	<b>3</b>
<b>5. SAFETY .....</b>	<b>3</b>
<b>6. EQUIPMENT AND SUPPLIES .....</b>	<b>4</b>
<b>7. REAGENTS AND STANDARDS.....</b>	<b>5</b>
<b>8. SAMPLE COLLECTION, PRESENTATION AND STORAGE .....</b>	<b>5</b>
<b>9. QUALITY CONTROL .....</b>	<b>6</b>
<b>10. CALIBRATION AND STANDARDIZATION.....</b>	<b>8</b>
<b>11. PROCEDURE.....</b>	<b>8</b>
<b>12. DATA ANALYSIS AND CALCULATIONS.....</b>	<b>10</b>
<b>13. METHOD PERFORMANCE .....</b>	<b>11</b>
<b>14. POLLUTION PREVENTION.....</b>	<b>11</b>
<b>15. WASTE MANAGEMENT.....</b>	<b>11</b>
<b>16. REFERENCES.....</b>	<b>12</b>
<b>17. MISCELLANEOUS (TABLES, APPENDICES, ETC...).....</b>	<b>12</b>

## **1. SCOPE AND APPLICATION**

- 1.1. This method is applicable to the determination of Total Organic Carbon in liquid, oils, sludge, soil, and sediment samples. It is based on Methods of Soil Analysis, Walkley-Black. The working linear range is 100 to 15,000 mg/kg.
- 1.2. The associated QuantIMs method code is VR.
- 1.3. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary by the laboratory QA department.

## **2. SUMMARY OF METHOD**

- 2.1. An aliquot of a solid sample is treated with excess potassium dichromate and concentrated sulfuric acid. After treatment, the solution is backtitrated with ferrous sulfate to determine the amount of dichromate reduced during digestion.

## **3. DEFINITIONS**

- 3.1. Refer to the glossary in the STL Laboratory Quality Manual (LQM), current version.

## **4. INTERFERENCES**

- 4.1. Method interferences may be caused by contaminants in solvents, reagents, glassware, and other processing apparatus that lead to discrete artifacts. All of these materials must be routinely demonstrated to be free from interferences under conditions of the analysis by running laboratory method blanks as described in the Quality Control section. Specific selection of reagents may be required to avoid introduction of contaminants.
- 4.2. Chloride and iron give a positive interference. Chloride may be totally or partially eliminated by the addition of mercuric sulfate.

## **5. SAFETY**

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL North Canton associates.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.

- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazard is known:
- 5.3.1. The following material is known to be **corrosive: Sulfuric acid**.
- 5.4. Exposure to chemicals must be maintained as **low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.6. It is recommended that neat standards be purchased only as a last resort. The preparation of standards from neat materials and reagents {as well as glassware cleaning procedures that involved solvents such as methylene chloride} should be conducted in a fume hood with the sash closed as far as the operations will permit.
- 5.7. Standards in solution may be diluted in the open laboratory when syringes and the like are utilized.
- 5.8. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Buret: 25 mL Class A
- 6.2. Analytical balance: capable of weighing to  $\pm 0.0001$  g
- 6.3. Top loading balance: capable of weighing to  $\pm 0.01$  g
- 6.4. Amber bottles
- 6.5. Beakers: various
- 6.6. Graduated cylinders: various

TOTAL ORGANIC CARBON (TOC)  
ANALYSIS FOR NON-WATERS

SOP No. NC-WC-0018

Revision No. 2.1

Revision Date: 03/30/01

Page 5 of 12

---

6.7. Volumetric pipettes: various , Class A

6.8. Erlenmeyer flasks: various

6.9. Whatman #4 filter paper

**7. REAGENTS AND STANDARDS**

7.1. Reagents

7.1.1. Sulfuric Acid ( $H_2SO_4$ ): concentrated, Tracepur grade

7.1.2. Ferroin indicator, purchased

7.1.3. Potassium Dichromate ( $K_2Cr_2O_7$ ): primary standard grade

7.1.4. 1N Potassium Dichromate Solution: Accurately weigh 49.04 g of potassium dichromate (dried overnight at  $105^\circ C$ ) in a liter volumetric flask and dilute to volume. Store in amber bottle and refrigerate. Replace after six months.

7.1.5. Ferrous Sulfate ( $FeSO_4 \cdot 7H_2O$ ): reagent grade

7.1.6. 0.5 N Ferrous Sulfate Titrant: Accurately weigh 140 g of  $FeSO_4 \cdot 7H_2O$  into a 1 liter volumetric flask and dissolve with 500 mL reagent water. Carefully add 15 mL of concentrated sulfuric acid and allow to cool. Dilute to volume with reagent water. Store in amber bottle and refrigerate.

7.1.7. Mercuric Sulfate ( $HgSO_4$ ): reagent grade

7.2. Standards

7.2.1. Laboratory Control Sample

7.2.1.1. Potassium Hydrogen Phthalate ( $KHC_8H_4O_4$ ), purchased

**8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

8.1. Samples are stored in a glass container at  $4^\circ C \pm 2^\circ C$ .

8.2. Samples are not chemically preserved. In lieu of no guidance, holding time is based on water requirements.

8.3. The holding time is twenty-eight days from sampling to analysis.

## 9. QUALITY CONTROL

### 9.1. Batch Definition

9.1.1. A batch is a group of no greater than 20 samples excluding QC samples (LCS, Method Blank, MS, MSD) which are processed similarly, with respect to the procedure. All sample setups must be initiated within a 24-hour period from the initial preparation or extraction and without interruption of the process. All samples within the batch must be treated with the same lots of reagents and the same processes.

### 9.2. Method Blank

9.2.1. One method blank (MB) must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting limit.

9.2.2. A reagent water blank consisting of 200-mL reagent water is being prepared and analyzed with each analytical batch of samples.

#### 9.2.3. Corrective Action for Blanks

9.2.3.1. If the analyte level in the method blank exceeds the reporting limit for the analytes of interest in the sample, all associated samples are reprepared and reanalyzed. If this is not possible due to limited sample quantity or other considerations, the corresponding sample data **must be addressed in the project narrative.**

9.2.3.2. If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. **Such action must be taken in consultation with the client and must be addressed in the project narrative.**

### 9.3. Laboratory Control Sample (LCS)

9.3.1. One LCS must be processed with each preparation batch. The LCS must be carried through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines.

9.3.2. A midrange LCS using 0.02 g of potassium hydrogen phthalate is prepared and analyzed with each batch of samples.

9.3.3. Corrective Action for LCS

9.3.3.1. If any analyte is outside established control limits the system is out of control and corrective action must occur.

9.3.3.2. The only exception is that if the LCS recoveries are biased high and the associated sample is ND for the parameter(s) of interest, the batch is acceptable. **This must be addressed in the project narrative.**

9.3.3.3. Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.

9.4. Sample Duplicate

9.4.1. A sample duplicate (DU) is a second aliquot of an environmental sample, taken from the same sample container when possible, that is processed with the first aliquot of that sample. That is, sample duplicates are processed as independent samples within the same QC batch. The sample and DU results are compared to determine the effect of the sample matrix on the precision of the analytical process. As with the MS/MSD results, the sample/DU precision results are not necessarily representative of the precision for other samples in the batch.

9.4.2. Sample duplicates are performed at a frequency of 10% or one per batch which ever is more frequent and must meet laboratory-specific limits for precision.

9.5. Control Limits

9.5.1. Control limits are established by the laboratory as described in SOP, NC-QA-0018.

9.5.2. Laboratory control limits are internally generated and updated periodically unless method specified. Control limits are listed in the Laboratory Quality Manual

---

(LQM) and the latest is version easily accessible via the LIMs (QC Browser program).

9.6. Method Detection Limits (MDLs) and MDL Checks

9.6.1. MDLs and MDL Checks are established by the laboratory as described in SOP, NC-QA-0021.

9.6.2. MDLs are listed in the Laboratory Quality Manual (LQM) and the latest version is easily accessible via the LIMs (QC Browser program).

9.7. Nonconformance and Corrective Action

9.7.1. Any deviations from QC procedures must be documented as a nonconformance, with applicable cause and corrective action approved by the facility QA Manager.

**10. CALIBRATION AND STANDARDIZATION**

10.1. The ferrous sulfate titrant is standardized daily as follows.

10.1.1. Pipette 10.0 mL of 1.00 N potassium dichromate solution into a 250 mL Erlenmeyer flask and add 90-mL reagent water.

10.1.2. Carefully add 30 mL of concentrated sulfuric acid and allow cooling completely.

10.1.3. Add 2-3 drops of ferroin indicator.

10.1.4. Titrate with 0.5 N ferrous sulfate titrant to a reddish-brown endpoint or to the first color change after reaching an emerald-green color.

10.1.5. Calculate the normality using the following equation.

$$N = \frac{10}{\text{mL ferrous sulfate}}$$

10.1.6. Repeat steps 10.1.1 through 10.1.5 two more times.

10.1.7. The average of the triplicate standardization is used.

**11. PROCEDURE**

TOTAL ORGANIC CARBON (TOC)  
ANALYSIS FOR NON-WATERS

8771113

SOP No. NC-WC-0018  
Revision No. 2.1  
Revision Date: 03/30/01  
Page 9 of 12

---

- 11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.
- 11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.
- 11.3. Sample Preparation
  - 11.3.1. Physical Preparation
    - 11.3.1.1. Mix the sample thoroughly before selecting a portion for analysis.
    - 11.3.1.2. Discard any foreign objects such as sticks, leaves, and rocks.
  - 11.3.2. Analytical Preparation
    - 11.3.2.1. Weigh an aliquot of soil of 2.50 g to the nearest 0.01 g (use less sample if TOC is known to be high). Record the weight on the analytical logsheet.
    - 11.3.2.2. Place sample in a 500 mL Erlenmeyer flask and add 10.0 mL of 1 N potassium dichromate.
    - 11.3.2.3. Under a hood, carefully add 20 mL of concentrated sulfuric acid and gently swirl for one minute.
    - 11.3.2.4. Allow sample to cool for about 30 minutes.
    - 11.3.2.5. Add 200 mL of reagent water and swirl to mix. If necessary, filter sample through Whatman #4 filter.
- 11.4. Sample Analysis Procedure
  - 11.4.1. Add 2-3 drops ferroin indicator.
  - 11.4.2. Titrate with 0.5 N ferrous Sulfate Solution to a reddish-brown endpoint or first color change after reaching an emerald-green color.



8771114

TOTAL ORGANIC CARBON (TOC)  
ANALYSIS FOR NON-WATERS

SOP No. NC-WC-0018  
Revision No. 2.1  
Revision Date: 03/30/01  
Page 10 of 12

---

11.4.2.1. If the digestate of the sample is already green or reddish-brown after the addition of the ferroin indicator, the sample needs to be re-extracted with a smaller sample amount or if less than 5 mLs of titrant is used.

11.4.3. Document the amount of titrant on the analytical logsheet.

11.5. Analytical Documentation

11.5.1. Record all analytical information in the analytical logbook/logsheets, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.5.2. All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. Logbooks are reviewed by the supervisor or designee.

11.5.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

11.5.4. Sample results and associated QC are entered into the LIMs after final technical review.

12. DATA ANALYSIS AND CALCULATIONS

12.1. Total Organic Carbon, mg/kg =

$$\left[ \frac{[(10)(N \text{ Potassium Dichromate}) - (mL \text{ Ferrous Sulfate}) (N \text{ Ferrous Sulfate})]}{\text{Weight of Soil (g)}} \times 3000 \right] \times 1.3$$

Where: 1.3 = Correction factor recommended in method

---

$$12.2. \quad \text{TOC, \%} = \frac{\text{mg / kg}}{10,000}$$

$$12.3. \quad \text{LCS, \%} = \frac{\text{TOC, \%}}{61.152 \text{ (true)}} \times 100$$

### 13. METHOD PERFORMANCE

13.1. Each laboratory must have initial demonstration of performance data on file and corresponding method detection limit files.

13.2. Training Qualifications:

13.2.1. The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

13.2.2. Method validation information (where applicable) in the form of laboratory demonstrations of capabilities is maintained for this method in the laboratory QA files.

### 14. POLLUTION PREVENTION

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

### 15. WASTE MANAGEMENT

15.1. Solvent waste must be disposed of in clearly labeled waste cans.

15.2. Acid waste must be collected in clearly labeled acid waste containers.

15.3. Solid materials (gloves, soiled paper products, etc.) are placed in the solid debris container. Do not put liquids in the solid waste container.

15.4. Refer to the Laboratory Sample and Waste Disposal plan.

15.5. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL North Canton. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by annual refresher training.

**16. REFERENCES**

## 16.1. References

16.1.1. Methods of Soil Analysis, 1982 Second Edition Method 29-3.5.2 Walkley-Black Procedure.

16.1.2. STL North Canton Laboratory Quality Manual (LQM), current version.

## 16.2. Associated SOPs and Policies, latest version

16.2.1. QA Policy, QA-003

16.2.2. Glassware Washing, NC-QA-0014

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018

16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021

16.2.5. Navy/Army SOP, NC-QA-0016

**17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)**

## 17.1. Reporting limits

17.1.1. The lower reporting limit (RL) for undiluted samples is 100 mg/kg.

17.1.2. If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.

## 17.2. Troubleshooting guide

17.2.1. When interferences as described in Section 4 are encountered or suspected, treat the sample as specified in that section.

17.2.2. If a high level of TOC is suspect (black sample), a smaller amount will be required.

Controlled Copy  
Copy No: \_\_\_\_\_

SOP No: NC-WC-0084

Revision No: 3

Revision Date: 10/03/00

Page 1 of 22

Implementation Date: 10/16/00

## STL STANDARD OPERATING PROCEDURE

TITLE: DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

(SUPERSEDES: REVISION (2))

Prepared by: Deborah N. Marcum 10-5-00  
Date

Reviewed by: Kurt C. Hill 10/5/00  
Technology Specialist Date

Approved by: Paul H. Feh 10/5/00  
Quality Assurance Manager Date

Approved by: [Signature] 10-10-00  
Environmental Health and Safety Coordinator Date

Approved by: Chris M. Prandi 10/11/00  
Laboratory Director Date

## Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00Page 2 of 22

---

1. SCOPE AND APPLICATION.....	3
2. SUMMARY OF METHOD .....	3
3. DEFINITIONS .....	3
4. INTERFERENCES .....	3
5. SAFETY.....	4
6. EQUIPMENT AND SUPPLIES .....	5
7. REAGENTS AND STANDARDS .....	5
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE.....	12
9. QUALITY CONTROL.....	12
10. CALIBRATION AND STANDARDIZATION.....	14
11. PROCEDURE.....	15
12. DATA ANALYSIS AND CALCULATIONS .....	16
13. METHOD PERFORMANCE.....	17
14. POLLUTION PREVENTION.....	17
15. WASTE MANAGEMENT .....	17
16. REFERENCES .....	18
17. MISCELLANEOUS (TABLES, APPENDICES, ETC...).....	18

**DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY**SOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00Page 3 of 22

---

**1. SCOPE AND APPLICATION**

- 1.1. This method covers the determination of fluoride, chloride, nitrite, bromide, nitrate, ortho-phosphate and sulfate in drinking water, surface water, mixed domestic and industrial wastewaters, groundwater, reagent waters, solids (after extraction 11.7) and leachates (when no acetic acid is used).
- 1.2. A listing of associated LIMs method codes is located in Section 8.2.
- 1.3. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary.

**2. SUMMARY OF METHOD**

- 2.1. A 25 uL volume of sample is introduced into the ion chromatograph. The sample is pumped through three different ion exchange columns and into a conductivity detector. The first two columns, a precolumn or guard column and a separator column, are packed with low-capacity, strongly basic anion exchange resin. Ions are separated into discrete bands based on their affinity for the exchange sites of the resin. The last column is a suppresser column that reduces the background conductivity of the eluent to a low or negligible level and converts the anions in the sample to their corresponding acids. The separated anions in their acid form are measured using an electrical conductivity cell. Anions are identified based on their retention times compared to known standards. Quantitation is accomplished by measuring the peak height or area and comparing it to a calibration curve generated from known standards.

**3. DEFINITIONS**

- 3.1. Refer to the glossary in the Laboratory Quality Manual (LQM).

**4. INTERFERENCES**

- 4.1. Interferences can be caused by substances with retention times that are similar to and overlap those of the anion of interest. Large amounts of an anion can interfere with the peak resolution of an adjacent anion. Sample dilution and/or fortification can be used to solve most interference problems associated with retention times.
- 4.2. The water dip or negative peak that elutes near, and can interfere with, the fluoride peak can usually be eliminated by the addition of concentrated eluent to each standard and sample.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

SOP No. NC-WC-0084

Revision No: 3

Revision Date: 10/03/00

Page 4 of 22

- 
- 4.3. Method interferences may be caused by contaminants in the reagent water, reagents, glassware and other sample processing apparatus that lead to discrete artifacts or an elevated baseline in the ion chromatograms.
  - 4.4. Any anion that is not retained by the column or only slightly retained will elute in the area of fluoride and interfere. Known coelution is caused by carbonate and other small organic anions. At concentrations of fluoride above 1.5 mg/L, this interference may not be significant; however, it is the responsibility of the user to generate precision and accuracy information in each sample matrix.
  - 4.5. The acetate anion elutes early during the chromatographic run. The retention times of the anions also seem to differ when large amounts of acetate are present. Therefore, this method is not recommended for leachates of solid samples when acetic acid is used for pH adjustment.

## 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL associates.
- 5.2. Eye protection that satisfied ANSI Z87.1 (as per the Chemical Hygiene Plan), laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory.
- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore; unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation when possible. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.6. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL associate. The situation must be reported **immediately** to a laboratory supervisor.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00

Page 5 of 22

---

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Balance -- Analytical, capable of accurately weighing to the nearest 0.0001 g.
- 6.2. Ion Chromatograph -- Analytical system complete with ion chromatograph and all required accessories including analytical columns, compressed gases and detectors.
  - 6.2.1. Anion guard column: A protector of the separator column. If omitted from the system the retention times will be shorter. Usually packed with same substrate as the separator column. 4 x 50 mm, Dionex IonPac AG14 P/N 46134, or equivalent.
  - 6.2.2. Anion separator column: The separation shown in Figure 1 was generated using a Dionex IonPac AS14 column (P/N 46134). Equivalent column may be used if comparable resolution is obtained, and the requirements of Sect. 9.2 can be met.
  - 6.2.3. Anion suppresser device: Dionex anion micro membrane suppresser (P/N 37106) or ASRS-Ultra Self-Regenerating Suppressor (4mm) P/N 53946 or equivalent.
  - 6.2.4. Detector -- Conductivity cell: approximately 1.25 uL internal volume, Dionex, or equivalent.
  - 6.2.5. Dionex --PeakNet 5.1 Data Chromatography Software or equivalent.
- 6.3. Assorted laboratory glassware (pipettes, volumetric flasks, etc.).

**7. REAGENTS AND STANDARDS**

- 7.1. Reagent grade chemicals shall be used in all tests. Unless otherwise indicated, it is intended that all reagents shall conform to the specifications of the Committee on Analytical Reagents of the American Chemical Society, where such specifications are available. Other grades may be used, provided it is first ascertained that the reagent is of sufficiently high purity to permit its use without lessening the accuracy of the determination.
- 7.2. Reagent water: Distilled or deionized water, free of the anions of interest. Water should contain particles no larger than 0.20 microns.
- 7.3. Eluent solution: sodium bicarbonate (CASRN 144-55-8) 1.0 mM, sodium carbonate (CASRN 497-19-8) 3.5 mM. Dissolve 1.680 g sodium bicarbonate ( $\text{NaHCO}_3$ ) and 7.417 g



of sodium carbonate ( $\text{Na}_2\text{CO}_3$ ) in reagent water (7.2) and dilute to 100 mL in a volumetric flask. Take 10 mL of this concentrated eluent solution and dilute to 2 L for use as the working eluent solution or dissolve the entire bicarbonate/carbonate amount in 20 L of reagent water.

7.4. Stock solutions (1,000 mg/L): All stocks may be prepared as described below or purchased from commercial sources. Primary and secondary sources are required for each target analyte.

7.4.1. Fluoride stock solution (1.00 mL = 1.00 mg  $\text{F}^-$ ): In a 1 liter volumetric flask, dissolve 2.2100 g of sodium fluoride ( $\text{NaF}$ ) in reagent water, and dilute to volume with reagent water. Store in chemical-resistant glass or polyethylene.

7.4.2. Chloride stock solution (1.00 mL = 1.00 mg  $\text{Cl}^-$ ): Dry sodium chloride ( $\text{NaCl}$ ) for 12 hours at  $105^\circ\text{C}$ , and cool in a desiccator. In a 1 liter volumetric flask, dissolve 1.6485 g of the dry salt in reagent water and dilute to volume with reagent water.

7.4.3. Nitrite stock solution (1.00 mL = 1.00 mg  $\text{NO}_2^- - \text{N}$ ): Place approximately 10.0 g of sodium nitrite ( $\text{KNO}_2$ ) in a 125 mL beaker and dry to constant weight (about 24 hours) in a desiccator. In a 1 liter volumetric flask, dissolve 6.0790 g of the dried salt in reagent water and dilute to volume with reagent water. Store in a sterilized glass bottle. Refrigerate and prepare monthly.

- Nitrite is easily oxidized, especially in the presence of moisture, and only fresh reagents are to be used.
- Prepare sterile bottles for storing nitrite solutions by heating for 1 hour at  $170^\circ\text{C}$  in an air oven.

7.4.4. Bromide stock solution (1.00 mL = 1.00 mg  $\text{Br}^-$ ): Dry approximately 5.0 g of sodium bromide ( $\text{NaBr}$ ) for 12 hours at  $105^\circ\text{C}$ , and cool in a desiccator. In a 1 liter volumetric flask, dissolve 1.2876 g of the dried salt in reagent water and dilute to volume with reagent water.

7.4.5. Nitrate stock solution (1.00 mL = 1.00 mg  $\text{NO}_3^- - \text{N}$ ): Dry approximately 10.00 g of sodium nitrate ( $\text{KNO}_3$ ) at  $105^\circ\text{C}$  for 24 hours. In a 1 liter volumetric flask, dissolve 7.2200 g of the dried salt in reagent water and dilute to volume with reagent water.

7.4.6. Phosphate stock solution (1.00 mL = 1.00 mg  $\text{PO}_4 - \text{P}$ ): Dry approximately 10.00 g of potassium dihydrogen phosphate ( $\text{KH}_2\text{PO}_4$ ) for 1 hour at  $105^\circ\text{C}$  and cool in a

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

8771123

SOP No. NC-WC-0084  
Revision No: 3  
Revision Date: 10/03/00  
Page 7 of 22

---

desiccator. In a 1 liter volumetric flask, dissolve 4.3937 g of the dry salt in reagent water and dilute to volume with reagent water.

- 7.4.7. Sulfate stock solution ( 1.00 mL = 1.00 mg  $\text{SO}_4^{--2}$ ): Dry approximately 5.00 g of potassium sulfate ( $\text{K}_2\text{SO}_4$ ) at 105°C for 1 hour and cool in a desiccator. In a 1 liter volumetric flask, dissolve 1.8141 g of the dried salt in reagent water and dilute to volume with reagent water.
- 7.4.8. Commercial stock solution A:  $\text{F}^-$  - 25 mg/L,  $\text{Cl}^-$  - 500 mg/L,  $\text{Br}^-$  - 100 mg/L,  $\text{NO}_3^-$  - N- 25 mg/L,  $\text{PO}_4 - \text{P}$  - 25 mg/L,  $\text{SO}_4^{--2}$  - 500 mg/L
- 7.4.9. Commercial stock solution B:  $\text{NO}_2^-$  - N- 25 mg/L
- 7.4.10. Commercial IC Spike solution A: :  $\text{F}^-$  - 125 mg/L,  $\text{Cl}^-$  - 2500 mg/L,  $\text{Br}^-$  - 500 mg/L,  $\text{NO}_3^-$  - N- 125 mg/L,  $\text{PO}_4 - \text{P}$  - 125 mg/L,  $\text{SO}_4^{--2}$  - 2500 mg/L
- 7.4.11. Commercial IC Spike solution B:  $\text{NO}_2^-$  - N- 125 mg/L
- 7.5. Working standards: Prepare calibration standard #5 in a 10 mL volumetric flask and transfer to a vial. Adjust the amount of stock solution used to prepare the working standards if the stock concentration differs from 1000 mg/L as assumed Alternatively prepare Cal standard #5 by mixing 4.0 mL commercial stock A, 4.0 mL commercial stock B and 2.0 mL of reagent water.

---

**Calibration Standard #5**

Analyte	mL of Stock	Final Conc.
Fluoride	0.10mL	10.0 mg/L
Chloride	2.0 mL	200. mg/L
Nitrite	0.10 mL	10.0 mg/L
Bromide	0.40 mL	40.0 mg/L
Nitrate	0.10 mL	10.0 mg/L
Ortho-Phosphate	0.10 mL	10.0 mg/L
Sulfate	2.0 mL	200. mg/L

7.5.1. In 5 mL PolyVials prepare the following calibration standards in reagent grade water.  
Final concentrations of working standards are shown below.

Calibration Standard #4: take 2.50 mL of calibration standard #5 and add 2.50 mL of reagent water.

Calibration Standard #2: take 250  $\mu$ L of calibration standard #5 and add 4.75 mL of reagent water.

Calibration Standard #1: take 25.0  $\mu$ L of calibration standard #5 and add 4.95 mL of reagent water.

Calibration Standard #3: take 1.25 mL of calibration standard #5 and add 3.75 mL of reagent water.

### Calibration Standard #1

Analyte	25.0 $\mu$ L of Cal Std #5	Final Conc
Fluoride		0.05 mg/L
Chloride		1.0 mg/L
Nitrite		0.05 mg/L
Bromide		0.20 mg/L
Nitrate		0.05 mg/L
Ortho-Phosphate		0.05 mg/L
Sulfate		1.0 mg/L

### Calibration Standard #2

Analyte	250 $\mu$ L of Cal Std #5	Final Conc.
Fluoride		0.5 mg/L
Chloride		10. mg/L
Nitrite		0.5 mg/L
Bromide		2.0 mg/L
Nitrate		0.5 mg/L
Ortho-Phosphate		0.5 mg/L
Sulfate		10. mg/L

**Calibration Standard #3**

Analyte	1.25 mL of Cal Std #5	Final Conc.
Fluoride		2.5 mg/L
Chloride		50. mg/L
Nitrite		2.5 mg/L
Bromide		10. mg/L
Nitrate		2.5 mg/L
Ortho-Phosphate		2.5 mg/L
Sulfate		50. mg/L

**Calibration Standard #4**

Analyte	2.5 mL of Cal Std #5	Final Conc.
Fluoride		5.0 mg/L
Chloride		100 mg/L
Nitrite		5.0 mg/L
Bromide		40. mg/L
Nitrate		5.0 mg/L
Ortho-Phosphate		5.0 mg/L
Sulfate		100 mg/L

- 7.5.2. Prepare or purchase a secondary stock standard(s) using a standards source other than that used for the primary standards as described in Section 7.5. Dilute these stock standards to as indicated in the table below to prepare the mixture to be used for the LCS and CCV solution. Alternatively prepare this solution by mixing 0.50 mL commercial stock A, 0.50 mL commercial stock B and 4.0 mL of reagent water.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

SOP No. NC-WC-0084

Revision No: 3

Revision Date: 10/03/00

Page 11 of 22

**LCS & Continuing Calibration Verification Solution**

Analyte	Final Conc. (V <sub>f</sub> =5ml)
Fluoride	2.5 mg/L
Chloride	50. mg/L
Nitrite	2.5 mg/L
Bromide	10. mg/L
Nitrate	2.5 mg/L
Ortho Phosphate	2.5 mg/L
Sulfate	50. mg/L

- 7.5.3. Prepare or purchase a secondary stock standard(s) using a standards source other than that used for the primary standards as described in Section 7.5. Dilute these stock standards to prepare the mixture to be used for the Matrix Spike solution. Alternatively purchase these mixes (ready to use) from a commercial source. Add 100 uL of each IC Spike solution to 5 mL of sample when preparing the MS/MSD. Dilute as needed after spiking the sample.

**Matrix Spike "True" Values**

Analyte	Final Conc.
Fluoride	2.5 mg/L
Chloride	50. mg/L
Nitrite	2.5 mg/L
Bromide	10. mg/L
Nitrate	2.5 mg/L
Ortho Phosphate	2.5 mg/L
Sulfate	50. mg/L

**NOTE:** Stock standards, calibration standard #5 and LCS standard should be stored in the dark at 4° ± 2°C. Replace these standards when instrument response indicates target analyte

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00

Page 12 of 22

degradation may have occurred or after the standard has expired (12 months commercial mix or 6 months in house mix), which ever occurs first. Nitrite and ortho-phosphate are particularly light and oxygen sensitive.

**8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

- 8.1. Samples should be collected in plastic or glass bottles. All bottles must be thoroughly cleaned and rinsed with reagent water. The volume collected should be sufficient to ensure a representative sample, allow for replicate analysis, if required, and minimize waste disposal.
- 8.2. Sample preservation and holding times for the anions that can be determined by this method for water samples are as follows:

QuanTIMs Method Code		Analyte	Preservation	Holding Time
EPA 300.0A	SW846 9056A	Fluoride	4° ± 2°C	28 days
C8	3C			
CX	3D	Chloride	4° ± 2°C	28 days
GO	3E	Nitrite	4° ± 2°C	48 hours
GM	3F	Bromide	4° ± 2°C	28 days
C9	3G	Nitrate	4° ± 2°C	48 hours
DO	3H	Ortho Phosphate	4° ± 2°C	48 hours
CY	3I	Sulfate	4° ± 2°C	28 days

**Note:** Soil leachates will follow the same preservation and holding times as the water samples; starting from the time of extraction.

**9. QUALITY CONTROL**

- 9.1. The STL QC Program document provides further details of the QC and corrective action guidelines presented in this SOP. Refer to this document if additional guidance is required.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

8771129

SOP No. NC-WC-0084  
Revision No: 3  
Revision Date: 10/03/00  
Page 13 of 22

---

- 9.2. Table I provides a summary of quality control requirements including type, frequency, acceptance criteria and corrective action.
- 9.3. Initial Demonstration of Capability
- 9.3.1. Prior to the analysis of any samples by Ion Chromatography, the following requirements must be met:
- 9.3.1.1. Method Detection Limit (MDL): An MDL must be determined prior to analysis of any samples. The MDL is determined using seven replicates of reagent water spiked with the anions of interest that has been carried through the entire analytical procedure. MDLs must be redetermined on an annual basis. The spike level must be greater than the calculated MDL but less than or equal to 10x the MDL. The result of the MDL determination must be below the STL reporting limit.
- 9.4. Batch definition: Preparation and QC batch definitions are provided in the STL QC Policy.
- 9.5. Method Blank (MB): One method blank must be processed with each preparation batch. The method blank consists of reagent grade water that has been taken through the entire preparation and analytical process. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest above the reporting limit.
- 9.6. Laboratory Control Sample (LCS): One LCS must be processed with each preparation batch and must be carried through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On-going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines. If the result is outside established control limits the system is out of control and corrective action must occur. Until in-house limits are established, a control limit of 90 - 110% recovery must be applied. Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable. The LCS consists of reagent grade water containing a known amount of target analytes that has been injected into the ion chromatography system. The LCS is prepared from a separate stock standard, or neat material, of a different manufacturer than the stock, or neat material, used to prepare the calibration standard.
- 9.7. Matrix Spike/Matrix Spike Duplicate (MS/MSD): One MS/MSD pair must be processed for each QC batch. A matrix spike (MS) is a field sample to which a known concentration of target analyte has been added. A matrix spike duplicate (MSD) is a second aliquot of the



same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific DQO's may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Spiking levels will be the same as the LCS values.

- If the MS/MSD recovery or RPD falls outside the acceptance range, the recovery of the analyte must be in control for the LCS. Until in-house control limits are established, a control limit of 90-110% recovery and 20% RPD must be applied to the MS/MSD.
- If the native analyte concentration in the MS/MSD exceeds 4x the spike level for that analyte, the recovery data are reported as NC (i.e. not calculated).
- If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted.
- If the recovery of the LCS is outside the limits, corrective action must be taken. Corrective action will include reparation and reanalysis of the batch.
- If a MS/MSD is not possible due to limited sample volume then a LCS duplicate must be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike limits.

9.8. Continuing Calibration Verification (CCV/CCB): Continuing calibration is verified by analyzing the calibration standard after every ten (10) samples. The CCV must fall within  $\pm 10\%$  of the true value for each target analyte. A CCB is analyzed immediately following the CCV to monitor low level accuracy and system cleanliness. The CCB result must be below the reporting limit for that analyte. If either the CCV or CCB fail to meet criteria, the analysis must be terminated, the problem corrected and reparation and analysis of all samples following the last CCV and CCB which were in control.

## 10. CALIBRATION AND STANDARDIZATION

- 10.1. Establish ion chromatographic operating parameters equivalent to those indicated in table 2. Refer to Table 3 for typical standard run retention times. Other than the presence of the analytical column the instrument conditions are the same.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00

Page 15 of 22

- 
- 10.2. For each analyte of interest, prepare a **minimum** of 3 calibration standards and a blank by adding accurately measured volumes of one or more stock standards to a volumetric flask and dilution to volume with reagent water. If a sample analyte concentration exceeds the calibration range the sample may be diluted to fall within the range. If this is not possible then three new calibration concentrations must be chosen, two of which must bracket the concentration of the sample analyte of interest. Each attenuation range of the instrument used to analyze a sample must be calibrated individually.
- 10.3. Using an injection volume of 25 uL of each calibration standard, tabulate peak height or area responses against the concentration. The results are used to prepare a calibration curve for each analyte. During this procedure, retention times must be recorded. All analytes will be calibrated using a quadratic regression forced through the origin. Correlation coefficients ( $R^2$ ) must be 0.995 or better.

**11. PROCEDURE**

- 11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.
- 11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.
- 11.3. Table 2 summarizes the recommended operating conditions for the ion chromatograph. Included in this table are estimated retention times that can be achieved by this method. Other columns, chromatographic conditions, or detectors may be used if the requirements of Sect. 9.2 are met.
- 11.4. Check system calibration daily as outlined in Table 1 and, if required, recalibrate as described in Sect 10.
- 11.5. Load and inject a fixed amount (25 uL) of settled & filtered sample. If the sample is cloudy then it should be filtered prior to loading into the autosampler polyvial. Flush injection loop thoroughly, using each new sample. Use the same size loop for standards and samples. Record the resulting peak size in area or peak height units. An automated constant volume injection system may also be used.
- 11.6. The width of the retention time window used to make identifications should be based upon measurements of actual retention time variations of various concentration. Three times the standard deviation of a retention time can be used to calculate a suggested window size for each analyte.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00Page 16 of 22

---

However, the experience of the analyst should weigh heavily in the interpretation of chromatograms since retention time is concentration dependent for most analytes..

- 11.7. If the response for the peak exceeds the working range of the system, dilute the sample with an appropriate amount of reagent water and reanalyze.
- 11.8. If the resulting chromatogram fails to produce adequate resolution, or if identification of specific anions is questionable, fortify the sample with an appropriate amount of standard and reanalyze.
- NOTE: Retention time is affected by concentration. Nitrate and sulfate exhibit the greatest amount of change, although all anions are affected to some degree. See Table 3. In some cases this peak migration may produce poor resolution or identification.
- 11.9. The following extraction should be used for solid materials: Add an amount of reagent water equal to ten times the weight of dry solid material taken as a sample. This slurry is mixed for one hour using a magnetic stirring device or tumbler. Filter the resulting slurry before injecting using a 0.45 um membrane type filter. This can be the type that attaches directly to the end of the syringe
- 11.10. Should more complete resolution be needed between peaks the eluent (7.3) can be diluted. This will spread out the run but will also cause the later eluting anions to be retained longer. The analyst must determine to what extent the eluent is diluted. This dilution should not be considered a deviation from the method.

**12. DATA ANALYSIS AND CALCULATIONS**

- 12.1. Prepare a calibration curve for each analyte by plotting instrument response against standard concentration. Compute sample concentration by comparing sample response with the standard curve. Multiply answer by appropriate dilution factor.
- 12.2. Report only those values that fall between the lowest and the highest calibration standards. Samples exceeding the highest standard should be diluted and reanalyzed.
- 12.3. Report results in mg/L for aqueous samples and mg/Kg for 1 hour leachates and mg/L for 18 hour leachates of solid samples.

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00Page 17 of 22

---

12.4. Report  $\text{NO}_2^-$  as N $\text{NO}_3^-$  as N $\text{HPO}_4^{=}$  as P**13. METHOD PERFORMANCE**

13.1. The reporting limits for the following analytes are based on a 25 uL injection volume:

Analyte	Water RL	Soil RL
Fluoride	1.0 mg/L	10 mg/kg
Chloride	1.0 mg/L	10 mg/kg
Nitrite	0.5 mg/L	5 mg/kg
Bromide	0.5 mg/L	5 mg/kg
Nitrate	0.05 mg/L	0.5 mg/kg
O-Phosphate	0.5 mg/L	5 mg/kg
Sulfate	1.0 mg/L	10 mg/kg

13.2. The group/team leader has the responsibility to ensure that this procedure is performed by an analyst who has been properly trained in its use and has the required experience. The analyst must be given two blind performance samples to analyze or process for analysis. Upon successful completion of the performance evaluation (PE) samples, these analyses will be documented as initial qualification,. Requalification must be performed annually thereafter for this procedure. The group/team leader must document the training and PE performance and submit the results to the QA Manager for inclusion in the associate's training files.

**14. POLLUTION PREVENTION**

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

**15. WASTE MANAGEMENT**

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00Page 18 of 22

---

- 15.1. Waste generated in this procedure must be segregated, and disposed of according to the facility hazardous waste procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**


- 16.1. Method 300.0, "Determination of Inorganic Anions by Ion Chromatography", Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. Environmental Protection Agency, Cincinnati, Ohio, Revision 2.1, August 1993.
- 16.2. Method 9056A, "Determination of Inorganic Anions by Ion Chromatography", SW846, Test Methods for Evaluating Solid Waste, Third Edition, Draft Revision 1, September 1999.
- 16.3. STL North Canton Laboratory Quality Manual (LQM), current version.

**17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)**

- 17.1. Attachment #1, method Flow Chart
- 17.2. Table 1, Quality Control Samples
- 17.3. Table 2, Standard Instrument Operating Parameters
- 17.4. Table 3, Retention Time Matrix
- 17.5. Figure 1, Example Chromatogram

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHYSOP No. NC-WC-0084Revision No: 3Revision Date: 10/03/00

Page 19 of 22

**Determination Of Inorganic Anions By Ion Chromatography**

ATTACHMENT #1

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

SOP No. NC-WC-0084

Revision No: 3

Revision Date: 10/03/00

Page 20 of 22

**TABLE 1**  
**QUALITY CONTROL SAMPLES**

<b>QC Samples</b>	<b>Frequency</b>	<b>Acceptance Criteria</b>	<b>Corrective Action</b>
Initial Calibration Verification (ICV)	At the start of each day following calibrating prior to sample analysis	+/- 10% of true value	Recalibrate and reanalyze
Initial Calibration Blank (ICB)	After Initial Calibration Verification and prior to sample analysis	< the Reporting Limit	Reprepare and reanalyze
Laboratory Control Sample (LCS)	1 per batch of 20 samples	Meets laboratory historical limits	Reanalyze all samples associated with unacceptable LCS
Matrix Spike Sample (MS/MSD)	1 MS/MSD pair per batch or 20 samples	Meets laboratory historical limits	Supervisor's technical judgment
Continuing Calibration Verification (CCV)	Between each group of 10 injections and at the end of the analytical sequence	+/- 10% of true value	Recalibrate and reanalyze all samples since the last acceptable CCV
Continuing Calibration Blank (CCB)	Between each group of 10 injections and at the end of the analytical sequence	< the Reporting Limit	Recalibrate and reanalyze all samples since the last acceptable CCB

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

SOP No. NC-WC-0084

Revision No: 3

Revision Date: 10/03/00

Page 21 of 22

**TABLE 2****Standard Instrument Operating Parameters**

Standard Conditions:

Eluent Pump Rate:	1.20 mL/min (DX-120 and DX-320)
Sample Loop:	25 uL
Eluent:	1.0mM sodium bicarbonate, 3.5mM sodium carbonate
Detector output	Baseline conductivity should be 15 - 20 uS prior to sample analysis.

**TABLE 3****Standard Run Retention Time Matrix (minutes)\***

Analyte	Concentration (mg/L)													RT window
	0.05	0.2	0.5	1	2	2.5	5	10	20	40	50	100	200	
F <sup>-</sup>	2.75	2.75				2.75	2.75	2.75						
Cl <sup>-</sup>				3.97				3.98			4.03	4.08	4.17	
NO <sub>2</sub> <sup>-</sup>	4.80	4.80				4.78	4.78	4.80						
Br <sup>-</sup>		6.15			6.13			6.10	6.08	6.07				
NO <sub>3</sub> <sup>-</sup>	7.33	7.27				7.17	7.13	7.07						
o-PO <sub>4</sub> <sup>2-</sup>	9.53	9.53				9.52	9.50	9.48						
SO <sub>4</sub> <sup>2-</sup>				11.50				11.48			11.43	11.38	11.27	

\* Analyte retention time is concentration dependent for most anions. Retention time increases with increasing concentration for chloride. Retention time decreases with increasing concentration for bromide, nitrate, ortho-phosphate and sulfate.



8771138

DETERMINATION OF INORGANIC ANIONS BY ION  
CHROMATOGRAPHY

SOP No. NC-WC-0084


Revision No: 3

Revision Date: 10/03/00

Page 22 of 22

**EXAMPLE ION CHROMATOGRAM**

*cal std 4 IC stds 9001/9002*



**Figure #1**

Controlled Copy

Copy No. \_\_\_\_\_

SOP No. NC-WC-0006

Revision No. 4

Revision Date: 02/06/01

Page 1 of 13

Implementation Date: 03/02/01

## STL NORTH CANTON STANDARD OPERATING PROCEDURE

TITLE: ALKALINITY (TOTAL)

(SUPERSEDES: REVISION 3, REVISION DATE 04/04/00)

Prepared by: Angela Silva 02/08/01  
Date

Reviewed by: K. A. Celli 02/08/01  
Technology Specialist Date

Approved by: Boulton 2/15/01  
Quality Assurance Manager Date

Approved by: [Signature] 2-16-01  
Environmental Health and Safety Coordinator Date

Approved by: Charles McQuinn 2-20-01  
Laboratory Director Date

## Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

---

**TABLE OF CONTENTS**

<i>1. Scope and Application.....</i>	<i>3</i>
<i>2. Summary of Method .....</i>	<i>3</i>
<i>3. Definitions.....</i>	<i>3</i>
<i>4. Interferences.....</i>	<i>3</i>
<i>5. SAFETY.....</i>	<i>4</i>
<i>6. Equipment and Supplies.....</i>	<i>4</i>
<i>7. Reagents and Standards.....</i>	<i>5</i>
<i>8. Sample Collection, PRESERVATION and Storage.....</i>	<i>6</i>
<i>9. Quality Control .....</i>	<i>6</i>
<i>10. Calibration and Standardization.....</i>	<i>8</i>
<i>11. Procedure .....</i>	<i>9</i>
<i>12. Data Analysis and Calculations.....</i>	<i>11</i>
<i>13. Method Performance .....</i>	<i>12</i>
<i>14. Pollution Prevention .....</i>	<i>12</i>
<i>15. Waste Management.....</i>	<i>12</i>
<i>16. References .....</i>	<i>12</i>
<i>17. Miscellaneous (Tables, Appendices, Etc...) .....</i>	<i>13</i>

## ALKALINITY (TOTAL)

SOP No. NC-WC-0006

Revision No. 4

Revision Date: 02/06/01Page 3 of 13

---

**1. SCOPE AND APPLICATION**

- 1.1. This method is applicable for the determination of total alkalinity in drinking, surface, saline, domestic, and industrial waters and wastewaters. It is also applicable to the determination of water-soluble alkalinity in solid samples if they have been prepared according to NC-WC-0073. It is based on EPA Method 310.1 and Standard Methods 2320B. The working linear range is 5 to 1000 mg/L.
- 1.2. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary by the laboratory QA department.
- 1.3. QuantIMs reference for total alkalinity is VC (310.1) and LV (2320B).

**2. SUMMARY OF METHOD**

- 2.1. An unaltered sample is titrated to an electrometrical endpoint of pH 4.5. **The sample must not be filtered, concentrated, or altered in any way.**

**3. DEFINITIONS**

- 3.1. Refer to the glossary in the STL North Canton Laboratory Quality Manual (LQM), current version.

**4. INTERFERENCES**

- 4.1. Method interference's may be caused by contaminants in solvents, reagents, glassware, and other processing apparatus that lead to discrete artifacts. All of these materials must be routinely demonstrated to be free from interferences under conditions of the analysis by running laboratory method blanks as described in the Quality Control section. Specific selection of reagents may be required to avoid introduction of contaminants.
- 4.2. Samples with salts of weak organic and inorganic acids and greases or oils will interfere with pH measurements.
- 4.3. The method is suitable for all concentration ranges of alkalinity; however, appropriate aliquots should be used to avoid a titration volume greater than 50 mL.

## ALKALINITY (TOTAL)

SOP No. NC-WC-0006Revision No. 4Revision Date: 02/06/01Page 4 of 13

---

**5. SAFETY**

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL North Canton associates.
- 5.2. Eye protection that satisfies ANSI Z87.1 (as per the Chemical Hygiene Plan), laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:
  - 5.3.1. The following materials are known to be **corrosive: Sulfuric Acid**
- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.6. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL NORTH CANTON associate. The situation must be reported **immediately** to a laboratory supervisor.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. Alkalinity – Manual
  - 6.1.1. Stir plate and stir bars
  - 6.1.2. Graduated cylinders: various
  - 6.1.3. Beakers: various
  - 6.1.4. Buret: Class A 25 mL or 50 mL (preferred)
- 6.2. Alkalinity - Automated

**ALKALINITY (TOTAL)**SOP No. NC-WC-0006Revision No. 4Revision Date: 02/06/01Page 5 of 13

---

6.2.1. Autotitrator

6.2.2. 50 mL centrifuge tubes

## 6.3. Alkalinity – Manual and Automated

6.3.1. pH meter and electrode(s) with temperature compensation

6.3.2. Volumetric pipettes: various

6.3.3. Autopipettor and disposable tips

6.3.4. Top loading balance: Capable of accurately weighing  $\pm 0.01$  g

6.3.5. Volumetric flasks: various

6.3.6. Oven

6.3.7. Desiccator

**7. REAGENTS AND STANDARDS**

## 7.1. Reagents

7.1.1. 0.02 N Sulfuric Acid: reagent grade, purchased, standardized monthly.

7.1.2. Sodium Carbonate ( $\text{Na}_2\text{CO}_3$ ): standard grade, dry overnight in  $180^\circ\text{C}$  oven and cool in a desiccator or purchased primary standard grade.7.1.3. Sodium Carbonate Solution: Add 2.50g of  $\text{Na}_2\text{CO}_3$  (record exact weight of  $\text{Na}_2\text{CO}_3$  used) to a 1000mL volumetric flask and dilute to volume with reagent water. Mix well.

## 7.2. Standards

## 7.2.1. Target Calibration Standard

7.2.1.1. pH Buffers: 4, 7, and 10 (Manufactured)

## 7.2.2. Laboratory Control Sample

7.2.2.1. Alkalinity Standard, 25,000 mg/L  $\text{CaCO}_3$  , purchased or other commercially available reference solutions.

## 7.2.3. Matrix Spike Standard

7.2.3.1. Alkalinity Standard, 25,000 mg/L  $\text{CaCO}_3$  , purchased

**8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

8.1. Samples are not chemically preserved.

8.2. Samples are stored in plastic or glass containers at  $4^\circ\text{C} \pm 2^\circ\text{C}$ .

8.3. The holding time is fourteen days from sampling to analysis.

8.4. The bottle must be filled with no headspace and provided in a separate container.

8.5. Do not open sample bottle before analysis. If other tests are to be performed from the same bottle, Alkalinity must be determined first. This is dependent on the client actually sending a separate bottle for alkalinity.

**9. QUALITY CONTROL**

## 9.1. Batch Definition

9.1.1. A batch is a group of no greater than 20 samples excluding QC samples (LCS, Method Blank, MS, MSD) which are processed similarly, with respect to the procedure. All sample setups must be initiated within a 24-hour period from the initial preparation or extraction and without interruption of the process. All samples within the batch must be treated with the same lots of reagents and the same processes.

## 9.2. Method Blank

9.2.1. One method blank (MB) must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including

preparation and analysis. The method blank is used to identify any system and process interference's or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting limit.

9.2.2. A reagent water blank consisting of 50 mL of reagent water and all other reagents added to samples within the analytical batch is analyzed with each analytical batch of samples.

9.2.3. Corrective Action for Blanks

9.2.3.1. If the analyte level in the method blank exceeds the reporting limit for the analytes of interest in the sample, all associated samples are reprepared and reanalyzed. If this is not possible due to limited sample quantity or other considerations, the corresponding sample data **must be addressed in the project narrative**.

9.2.3.2. If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. **Such action must be taken in consultation with the client and must be addressed in the project narrative.**

9.3. Laboratory Control Sample (LCS)

9.3.1. One aqueous LCS must be processed with each preparation batch. The LCS must be carried through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines.

9.3.2. An LCS consisting of 1mL of the 25,000 mg/L alkalinity standard and 50 mL reagent water or other commercially available reference solution is analyzed with each analytical batch of samples.

9.3.3. Corrective Action for LCS

9.3.3.1. If any analyte is outside established control limits the system is out of control and corrective action must occur.

9.3.3.2. Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.



#### 9.4. Matrix Spike/Matrix Spike Duplicate (MS/MSD)

9.4.1. One MS/MSD pair must be processed for each batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Samples identified as field blanks cannot be used for MS/MSD analysis.

9.4.2. An MS/MSD consisting of 1 mL of the 25,000 mg/L alkalinity standard and 50 mL of the sample will be analyzed.

##### 9.4.3. Corrective action for MS/MSDs

- 9.4.3.1. If the analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include reparation and reanalysis of the batch.
- 9.4.3.2. If the native analyte concentration in the MS/MSD exceeds 4x the spike level for that analyte, the recovery data is reported as DIL (diluted out).
- 9.4.3.3. If an MS/MSD is not possible due to limited sample volume then a laboratory control sample duplicate (LCSD) should be analyzed. The RPD of the LCS and LCSD must be compared to the matrix spike RPD limits.

#### 9.5. QC Acceptance Criteria

- 9.5.1.1. Control limits are established by the laboratory as described in NC-QA-0018.

### 10. CALIBRATION AND STANDARDIZATION

#### 10.1. Instrument Directions

10.1.1. Calibrate the pH meter according to the manufacturer's specifications. See pH Electrode Method SOP # NC-WC-0010.

10.2. Initial Calibration

10.2.1. The pH meter is calibrated everyday with the 4 and the 7 calibration buffers and is verified at the beginning of the run by using the 10 buffer. The pH buffers should bracket the sample concentration.

10.3. Continuing Calibration

10.3.1. The pH meter is checked every ten readings with a midrange (pH 7) buffer to ensure the calibration remain linear. The acceptance range for the calibration check is  $7 \pm 0.05$  pH units or recalibration is necessary.

**11. PROCEDURE**

11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

11.3. Sample Preparation

11.3.1. For solids preparation, see SOP NC-WC-0073.

11.3.2. No preparation is necessary for water samples.

11.4. Standardization

11.4.1. To standardize 0.02 N sulfuric acid, titrate 50 mL reagent water and 0.125 g sodium carbonate (weighed accurately and recorded) with 0.02 N  $H_2SO_4$  to a pH of 4.5. This should be performed monthly or on a new lot of acid (whichever is more frequent). Calculate as follows:

8771148

ALKALINITY (TOTAL)

SOP No. NC-WC-0006

Revision No. 4

Revision Date: 02/06/01

Page 10 of 13

---

$$N = \frac{A \times 1000}{53.00 \times B} \text{ (manual)}$$

$$N = \frac{A \times 1000}{53.00 \times B} \times \frac{20}{50} \text{ (Autotitrator)}$$

A=gNa<sub>2</sub>CO<sub>3</sub>

B=mL .02 N H<sub>2</sub>SO<sub>4</sub> titrant

11.5. Repeat standardization two or three more times. Record standardization in the calibration logbook.

11.6. Sample Analysis – Manual

11.6.1. Do not shake sample.

11.6.2. Record the initial pH prior to sample analysis.

11.6.3. Record the initial pH prior to sample analysis. Use a sufficiently large volume of titrant (>20 mL in a 50 mL buret) to obtain good precision while keeping the volume low enough to permit a sharp end point

11.6.4. Place 50 mL of sample, or an aliquot diluted to 50 mL with reagent water, in a beaker. Begin mixing; measure and record initial pH of the sample. Titrate the sample to an endpoint of pH 4.5 with 0.02 N H<sub>2</sub>SO<sub>4</sub>. Record the volume of the titrant on the analytical logsheet. Samples requiring >50 mL titrant (> 75 mL using autotitrator) should be re-analyzed using less sample volume. Record dilution.

11.6.5. Use a sufficiently large volume of titrant (>20 mL in a 50 mL buret) to obtain good precision while keeping the volume low enough to permit a sharp end point.

11.7. Sample Analysis – Automated – Summary

11.7.1. The samples are analyzed on the autotitrator for Total Alkalinity.

11.7.2. Do not shake sample.

11.7.3. Place 50 mL of sample, or an aliquot diluted to 50 mL with reagent water, in a 50 mL centrifuge tube. See Manufacturer's information for operating instructions.

11.7.4. If a dilution of the sample was done change the volume on the schedule to reflect the dilution.(based on a 20 mL sample inject) For alkalinity the dilution factor will be taken into account in the final calculation. Do not manually multiply the dilution unless it was not typed into the schedule.

11.7.5. After the results have been gathered from the instrument make sure to check the pH of all the samples. If a sample has an initial pH of >4.5 and the Total Alkalinity is zero, the sample must be diluted and reanalyzed.

#### 11.8. Analytical Documentation

11.8.1. Record all analytical information in the analytical logbook/logsheets, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.8.2. All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. The supervisor or designee reviews logbooks.

11.8.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

11.8.4. Sample results and associated QC are entered into the LIMs after final technical review.

## 12. DATA ANALYSIS AND CALCULATIONS

$$12.1. \text{ Alkalinity, mg/L CaCO}_3 \text{ to pH 4.5} = \frac{A \times N \times 50,000}{\text{mL of sample}}$$

$$12.2. \text{ LCS \%} = \frac{\text{mg / L}}{500 \text{ (true)}} \times 100$$

$$\text{MS/MSD \%} = \frac{B - C}{500 \text{ (true)}} \times 100$$

Where:

A = mL of Titrant

N = Normality of Titrant

B = MS/MSD, mg/L

C = Sample, mg/L

### **13. METHOD PERFORMANCE**

13.1. Each laboratory must have initial demonstration of performance data on file and corresponding method detection limit files.

13.2. Training Qualifications:

13.2.1. The group/team leader has the responsibility to ensure that an associate who has been properly trained in its use and has the required experience performs this procedure.

### **14. POLLUTION PREVENTION**

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

### **15. WASTE MANAGEMENT**

15.1. Acid waste must be collected in clearly labeled acid waste containers.

15.2. Solid materials (gloves, soiled paper products, etc.) are placed in the solid debris container. Do not put liquids in the solid waste container.

15.3. Refer to the Laboratory Sample and Waste Disposal plan.

15.4. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by an annual refresher training.

### **16. REFERENCES**

16.1. References

## ALKALINITY (TOTAL)

SOP No. NC-WC-0006

Revision No. 4

Revision Date: 02/06/01

Page 13 of 13

8771151

16.1.1. EPA-600/4-79-020, Methods for Chemical Analysis of Water and Wastes, Revised March 1983, Alkalinity, Method 310.1.

16.1.2. Standard Methods for the Examination of Water and Wastewater, 18th Edition, 1992, Alkalinity Methods, 2320B.

16.1.3. EPA 600, Methods for Chemical Analysis of Water and Wastes, pH, Method 150.1.

### 16.2. Associated SOPs

16.2.1. Solid Extraction for Wet Chemistry Parameters, NC-WC-0073.

16.2.2. pH Electrode Method for Wet Chemistry Parameters, NC-WC-0010.

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018.

## 17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)

### 17.1. Reporting limits

17.1.1. The lower reporting limit (RL) for undiluted samples is 5 mg/L  $\text{CaCO}_3$ .

17.1.2. If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.

### 17.2. Method Deviation

17.2.1. A fixed endpoint of 4.5 is used for all samples since the sample concentration is often unknown.

17.2.2. The Sodium Carbonate ( $\text{Na}_2\text{CO}_3$ ) is dried at 180°C overnight instead of at 250° C for 4 hours.

17.2.3. The standard acid solution is not boiled gently for 3-5 minutes under a watch glass cover.

8771152

Controlled Copy

Copy No. \_\_\_\_\_

Implementation Date 7-18-03SOP No. NC-WC-0060Revision No. 3Revision Date: 06/17/03Page 1 of 19**STL NORTH CANTON STANDARD OPERATING PROCEDURE****TITLE: SULFIDE****(SUPERSEDES: REVISION 2; DATED 09/17/02)**

Prepared by: *Stephanie A. Fuller-Gustaf* 7/14/03  
Date

Reviewed by: *Angela Suria* 7/11/03  
Technology Specialist Date

Approved by: *Beth Gumbert* 7/14/03  
Quality Assurance Manager Date

Approved by: *Deborah Budd* 7/17/03  
Environmental Health and Safety Coordinator Date

Approved by: *Pauline Corbett* 7/18/03  
Laboratory Director Date

**Proprietary Information Statement:**

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions. THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

## TABLE OF CONTENTS

1. <i>Scope and Application</i> .....	3
2. <i>Summary of Method</i> .....	3
3. <i>Definitions</i> .....	3
4. <i>Interferences</i> .....	3
5. <i>SAFETY</i> .....	4
6. <i>Equipment and Supplies</i> .....	5
7. <i>Reagents and Standards</i> .....	5
8. <i>Sample Collection, PRESERVATION and Storage</i> .....	8
9. <i>Quality Control</i> .....	8
10. <i>Calibration and Standardization</i> .....	11
11. <i>Procedure</i> .....	11
12. <i>Data Analysis and Calculations</i> .....	16
13. <i>Method Performance</i> .....	17
14. <i>Pollution Prevention</i> .....	17
15. <i>Waste Management</i> .....	17
16. <i>References</i> .....	17
17. <i>Miscellaneous (Tables, Appendices, Etc...)</i> .....	18



## 1. SCOPE AND APPLICATION

- 1.1. This method is applicable to the determination of the concentration of Sulfide in waters, liquids, solids, and sludges. It is based on SW846 Method 9030B and Methods for Chemical Analysis of Water and Wastes (MCAWW) 376.1. The working range is 1 to 30 mg/L for waters and 10-650 mg/kg for solids and sludges.
- 1.2. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary by the laboratory QA department.
- 1.3. The associated QuantIMS method codes are TV (9030B) and CT (376.1).

## 2. SUMMARY OF METHOD

- 2.1. For acid soluble sulfide samples, separation of sulfide from the sample matrix is accomplished by the addition of sulfuric acid to the sample. The sample is heated to 70<sup>0</sup>C and the hydrogen sulfide (H<sub>2</sub>S) which is formed, is distilled under acidic conditions and carried by a nitrogen stream into zinc acetate scrubbing bottles where it is precipitated as zinc sulfide.
- 2.2. For acid-insoluble sulfide samples, separation of sulfide from the sample matrix is accomplished by suspending the sample in concentrated hydrochloric acid by vigorous agitation. Tin (II) chloride is present to prevent oxidation of sulfide to sulfur by the metal ion (as in copper (II)), by the matrix, or by dissolved oxygen in the reagents. The prepared sample is distilled under acidic conditions at 100<sup>0</sup>C under a stream of nitrogen. Hydrogen sulfide gas is released from the sample and collected in gas scrubbing bottles containing zinc(II) and a strong acetate buffer. Zinc sulfide precipitates.
- 2.3. An excess of iodine is added to a sample which oxidizes the Sulfide to sulfur under acidic conditions. The excess iodine is back titrated with sodium thiosulfate.

## 3. DEFINITIONS

- 3.1. Refer to the glossary in the STL North Canton Laboratory Quality Manual (LQM).

## 4. INTERFERENCES

- 4.1. Method interferences may be caused by contaminants in solvents, reagents, glassware, and other processing apparatus that lead to discrete artifacts. All of these materials must be routinely demonstrated to be free from interferences under conditions of the analysis

---

by running laboratory method blanks as described in the Quality Control section. Specific selection of reagents may be required to avoid introduction of contaminants.

- 4.2. Reducing substances such as thiosulfite, sulfites, and various organic compounds cause interferences, but treatment with zinc acetate solution will eliminate some of these interferences. (Use approximately 15 drops of 2 N zinc acetate per 500 mL of sample if not already preserved with it.)
- 4.3. Samples that contain strong oxidizers or reducers will interfere with this method.

## 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL North Canton associates.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:
  - 5.3.1. The following materials are known to be **corrosive: Hydrochloric Acid, Sodium Hydroxide, and Sulfuric Acid.**
- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.6. It is recommended that neat standards be purchased only as a last resort. The preparation of standards from neat materials and reagents should be conducted in a fume hood with the sash closed as far as the operations will permit.
- 5.7. Standards in solution may be diluted in the open laboratory when syringes and the like are utilized.

- 5.8. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Volumetric pipettes: various
- 6.2. Buret: 25 mL Class A
- 6.3. Erlenmeyer flasks: 500 mL
- 6.4. Graduated cylinder: 250 mL
- 6.5. Top loading balance: capable of accurately weighing  $\pm 0.01$  g
- 6.6. Volumetric flasks: various
- 6.7. Vacuum pump, filter and flask
- 6.8. Whatman 934-AH filters
- 6.9. Distillation apparatus containing: 250 mL addition funnel, 500 mL – 3-neck reaction flask, sparging tube and 1 – 500 mL Erlenmeyer Flask.
- 6.10. Stirring / Hot Plates
- 6.11. Crystallizing dishes

## 7. REAGENTS AND STANDARDS

### 7.1. Reagents

- 7.1.1. (1:1) Hydrochloric Acid: Add 250 mL concentrated hydrochloric acid (HCl) to 250 mL of reagent water.
- 7.1.2. Starch Indicator: Add 10 mL of reagent water to 5 g starch (potato) and mix. Add starch mixture to 500 mL of boiling reagent water. Mix, cool, and store in a well-labeled squirt bottle. Alternately, use purchased starch solution.

7.1.3. 0.025 N Sodium Thiosulfate (stored in dessicator): Add 0.4 g NaOH and 6.205 g of sodium thiosulfate ( $\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$ ) to 500 mL of reagent water in a 1 liter volumetric. Dilute to volume with reagent water. Store in a dark container. Also available commercially.

7.1.3.1. Standardization of 0.025 N Sodium Thiosulfate Solution: \*To make 0.025N Biodate Solution, dissolve 0.462 g  $\text{KH}(\text{IO}_3)_2$  in 500 mL with reagent water. Weigh 2 g KI in a 500 mL Erlenmeyer flask. Add 100 to 150 mL reagent water, 5 drops  $\text{H}_2\text{SO}_4$  and 20 mL biodate solution using a volumetric pipet. Dilute to 200 mL with reagent water. Titrate with Sodium Thiosulfate. When a pale straw yellow color is reached, add 1-2 mL starch. Continue titrating from a blue to a clear end point.

\*Note: Biodate Solution may be purchased.

Calculation

$$\text{Na}_2\text{S}_2\text{O}_3 \text{ Normality} = \frac{(a)(b)}{c}$$

*a* = mLs Biodate (20 mL)

*b* = Normality Biodate (0.025N)

*c* = mLs of  $\text{Na}_2\text{S}_2\text{O}_3$  used to titrate Repeat two more times

7.1.4. 0.0282 N Iodine Solution: Add 20 g KI (potassium iodide) and 3.2 g iodine to a 1 liter volumetric flask. Add 500 - 700 mL of reagent water and dissolve. Dilute to volume with reagent water. Store in a dark container. Also available commercially.

7.1.4.1. Standardization 0.025 N Iodine Solution: Perform three method blanks daily. Refer to method blank section in SOP.

#### Calculation

$$\text{Normality Iodine} = \frac{(\text{Normality Na}_2\text{S}_2\text{O}_3)(\text{mL of titrant Na}_2\text{S}_2\text{O}_3)}{20 \text{ mL Iodine}}$$

7.1.5. 2N Zinc Acetate: Dissolve 220 g of zinc acetate in 870 mL of reagent water and dilute to 1 liter with reagent water.

7.1.6. Formaldehyde (37% solution), CH<sub>2</sub>O. This solution is commercially available.

7.1.7. Zinc Acetate for the Erlenmeyer flasks:

7.1.7.1. For acid-soluble Zinc: Zinc acetate solution (approximately 0.5M). Dissolve 110g Zinc acetate, dihydrate NaC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>, in 800 mL of reagent water. Add 1 mL concentrated hydrochloric acid and dilute to 1 liter.

7.1.7.2. For acid-insoluble sulfides: Zinc acetate/sodium acetate buffer. Dissolve 100 g sodium acetate, NaC<sub>2</sub>H<sub>3</sub>O<sub>2</sub>, and 11 g zinc acetate dihydrate in 800 mL of reagent water. Add 1 mL concentrated hydrochloric acid and dilute to 1 liter. The resulting pH should be 6.8

7.1.8. Sulfuric acid – 50%. Place 450 mL of reagent water in a volumetric flask. **Slowly** add 500 mL concentrated Sulfuric Acid (H<sub>2</sub>SO<sub>4</sub>). **Use extreme caution – this is an exothermic reaction and will create excess heat.** This solution is commercially available.

7.1.9. Hydrochloric Acid, 9.8N, for **acid-insoluble sulfides**: Place 200 mL of reagent water in a 1-liter beaker. Slowly add concentrated HCl to bring the total volume to 1 liter.

7.1.10. Tin (II) chloride, SnCl<sub>2</sub>, granular

---

## 7.2. Standards

### 7.2.1. Laboratory Control Sample

7.2.1.1. 2000 ppm Sulfide: Add 3.5 g of sodium sulfide to 100 mL of reagent water in a 250mL volumetric flask. Dilute to volume with reagent water. The sulfide standard must be verified each working day. If the resulting value is <75% of the original standard, a new solution must be prepared.

### 7.2.2. Matrix Spike Standard

7.2.2.1. Prepare a midrange matrix spike standard as described in 7.2.1 for use as a MS/MSD.

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

8.1. Waters are preserved to a pH > 9 with NaOH and zinc acetate. Non-water samples are unpreserved. All matrices are stored at 4°C ± 2°C in plastic or glass containers.

8.2. The holding time is seven days from sampling to analysis.

## 9. QUALITY CONTROL

### 9.1. Batch Definition

9.1.1. A batch is a group of no greater than 20 samples excluding QC samples (LCS, Method Blank, MS, MSD) which are processed similarly, with respect to the procedure. All sample setups must be initiated within a 24 hour period from the initial preparation or extraction and without interruption of the process. All samples within the batch must be treated with the same lots of reagents and the same processes.

### 9.2. Method Blank

9.2.1. One method blank (MB) must be processed with each preparation batch. The method blank consists of reagent water containing all reagents specific to the method that is carried through the entire analytical procedure, including preparation and analysis. The method blank is used to identify any system and process interferences or contamination of the analytical system that may lead to the reporting of elevated analyte concentrations or false positive data. The method blank should not contain any analyte of interest at or above the reporting limit.

9.2.2. A reagent water blank consisting of 250 mL of reagent water must be analyzed with each analytical batch of samples.

9.2.3. Corrective Action for Blanks

9.2.3.1. If the analyte level in the method blank exceeds the reporting limit for the analytes of interest in the sample, all associated samples are reprepared and reanalyzed. If this is not possible due to limited sample quantity or other considerations, the corresponding sample data **must be addressed in the project narrative.**

9.2.3.2. If there is no analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. **Such action must be taken in consultation with the client and must be addressed in the project narrative.**

9.3. Laboratory Control Sample (LCS)

9.3.1. One aqueous LCS must be processed with each preparation batch. The LCS must be carried through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On-going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines.

9.3.2. A midrange LCS is prepared by adding 2.5 mL for Method 376.1 and 1.0 mL for Method 9030B (water samples) or 1.0 mL (solid samples) of 2000ppm sulfide standard to a flask. This standard must be analyzed with each analytical batch of samples.

9.3.3. Corrective Action for LCS

9.3.3.1. If any analyte is outside established control limits the system is out of control and corrective action must occur.

9.3.3.2. The only exception is that if the LCS recoveries are biased high and the associated sample is ND for the parameter(s) of interest, the batch is acceptable. **This must be addressed in the project narrative.**

9.3.3.3. Corrective action will be repreparation and reanalysis of the batch unless the client agrees that other corrective action is acceptable.

#### 9.4. Matrix Spike/Matrix Spike Duplicate (MS/MSD)

9.4.1. One MS/MSD pair must be processed for each batch. A matrix spike (MS) is a field sample to which known concentrations of target analytes have been added. A matrix spike duplicate (MSD) is a second aliquot of the same sample (spiked identically as the MS) prepared and analyzed along with the sample and matrix spike. Some client specific data quality objectives (DQO's) may require the use of sample duplicates in place of or in addition to MS/MSD's. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, these results may have immediate bearing only on the specific sample spiked. Samples identified as field blanks cannot be used for MS/MSD analysis.

9.4.2. A MS/MSD consisting of 30 g or 250 mL of the sample and 2.5 mL for Method 376.1 or 1.0 mL for Method 9030B of 2000 ppm sulfide standard should be analyzed with every 20 samples.

##### 9.4.3. Corrective action for MS/MSDs

9.4.3.1. If the analyte recovery or RPD falls outside the acceptance range, the recovery of that analyte must be in control for the LCS. If the LCS recovery is within limits, then the laboratory operation is in control and the results may be accepted. If the recovery of the LCS is outside limits, corrective action must be taken. Corrective action will include reparation and reanalysis of the batch.

9.4.3.2. If the native analyte concentration in the MS/MSD exceeds 4x the spike level for that analyte, the recovery data is reported as "amount" MSB. The Exception Code is changed to NC. The following two footnotes will appear on the report page "NC The recovery and/or RPD were not calculated." "MSB The recovery and RPD were not calculated because the sample amount was greater than four times the spike amount."

9.4.3.3. If an MS/MSD is not possible due to limited sample volume then a laboratory control sample duplicate (LCSD) should be analyzed. The RPD of the LCS and LCSD must be compared to the laboratory limits.

#### 9.5. Control Limits

9.5.1. Control limits are established by the laboratory as described in SOP, NC-QA-0018.



9.5.2. Laboratory control limits are internally generated and updated periodically unless method specified. Control limits are listed in the Laboratory Quality Manual (LQM) and the latest is version easily accessible via the LIMs (QC Browser program).

9.6. Method Detection Limits (MDLs) and MDL Checks

9.6.1. MDLs and MDL Checks are established by the laboratory as described in SOP, NC-QA-0021.

9.6.2. MDLs are listed in the Laboratory Quality Manual (LQM) and the latest version is easily accessible via the LIMs (QC Browser program).

9.7. Nonconformance and Corrective Action

9.7.1. Any deviations from QC procedures must be documented as a nonconformance, with applicable cause and corrective action approved by the facility QA Manager.

**10. CALIBRATION AND STANDARDIZATION**

10.1. Not Applicable

**11. PROCEDURE**

11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

11.3. Sample Preparation

11.3.1. Water Samples Method 376.1

11.3.1.1. Using the vacuum apparatus, filter 250 mL (or a lesser aliquot) of well homogenized sample. Discard the filtrate and analyze the filter pad and solids following section 11.4. Record volume filtered on the analytical logsheet.

### 11.3.2. Aqueous Sample Method 9030B

11.3.2.1. For an efficient distillation, the mixture in the distillation flask must be of such a consistency that the motion of the stirring bar is sufficient to keep the solids from settling. The mixture must be free of solid objects that could disrupt the stirring bar. Prepare the sample using one of the procedures in this section then proceed with the distillation step.

11.3.2.2. If the sample is aqueous, shake the sample container to suspend any solids, then quickly decant the appropriate volume (250 mL) of the sample to a graduated cylinder. Transfer the contents of the graduated cylinder into the reaction flask.

11.3.2.3. If the sample is aqueous, but contains a large proportion of solids, the sample may be roughly separated by phase and the amount of each phase measured and weighed to the nearest milligram into the distillation flask in proportion to their abundance in the sample. Reagent water may be added up to a total volume of 250 mL.

### 11.3.3. Solids and Waste Samples, Method 9030B

11.3.3.1. Weigh out 30 g +/- 0.1 g of homogenized sample and put into the reaction flask.

11.3.3.2. Samples that are not water miscible (oils, various solvents) cannot be analyzed using this method.

### 11.3.4. 9030B Distillation

#### 11.3.4.1. Acid Soluble Sulfides:

11.3.4.1.1. Add 30g +/- 0.1g of solid sample plus 200 mL of reagent water or 250 mL of aqueous sample and a stir bar to the reaction flask

11.3.4.1.2. Attach the reaction flask to the distillation apparatus such that the bottom of the reaction flask does not touch the bottom of the crystallizing dish (submerge approximately 1/3 of the flask in the warm water).

11.3.4.1.3. Add 100 mL of 50% H<sub>2</sub>SO<sub>4</sub> (7.1.8) to the addition flask and place in the center neck of the reaction flask. Attach the nitrogen flow line to the top of the addition flask.

11.3.4.1.4. Prepare one (1) gas trap bottles for each distillation setup by adding 200 mL of reagent water, 40 mL of zinc acetate buffer (7.1.7.1) and 10 mL of formaldehyde (7.1.6) to each trap.

11.3.4.1.5. Connect the trap and insert the trap arm into the right neck of the reaction flask. Turn on the nitrogen to approximately 5 psi and adjust the flow to 3 – 5 bubbles per second and purge for 15 minutes.

11.3.4.1.6. Add spike and LCS solution by removing the trap arm and pipetting the spike solution, below the surface of the water. Replace the trap arm and re-establish the flow of nitrogen, purge for 5 minutes.

11.3.4.1.7. Open the addition funnel to add the sulfuric acid (7.1.8) drop by drop. Do not open the stop cock fully. Distill the sample for 90 minutes, maintaining 70<sup>0</sup>C in the water bath.

11.3.4.1.8. Fill the crystallizing dish with 650 – 700 mL of reagent water and place on the hotplate/stirrer. Turn on the hotplate to a setting of “6” for approximately 15 – 20 minutes. Turn down to a setting of “4” to maintain a temperature of 70<sup>0</sup>C, +/- 5<sup>0</sup>C.

#### 11.3.4.2. Acid-Insoluble Sulfide:

11.3.4.2.1. As the concentration of HCl during distillation must be within a narrow range for successful distillation of H<sub>2</sub>S, the water content must be controlled. It is imperative that the final concentration of HCl in the distillation flask be about 6.5N and that the sample is mostly suspended in the fluid by the action of the stirring bar. This is achieved by adding 50 mL of reagent water, including water in the sample, 100 mL of 9.8N HCl, and the sample to the distillation flask. Solids which absorb water and swell will restrict fluid motion and, therefore, lower recovery will be obtained. Such samples should be limited to 25 g dry weight.

11.3.4.2.2. If the matrix is aqueous, then a maximum of 50 g of the sample may be used. No additional water may be added.

11.3.4.2.3. If the matrix is dry solid, use 30 g +/- 0.1 g of the sample and add 50 mL of reagent water.

11.3.4.2.4. Add 5 g SnCl to each distillation flask

11.3.4.2.5. Assemble the distillation apparatus. Place 200 +/- 4.0 mL of zinc acetate/sodium acetate buffer solution and 10.0 mL +/- 2.0 mL of 37% formaldehyde in each Erlenmeyer flask. Add 40 mL DI water.

11.3.4.2.6. Add 100 +/- 1.0 mL of 9.8N HCL to the addition funnel. Connect the nitrogen line to the top of the funnel and turn the nitrogen on to pressurize the dropping funnel headspace.

11.3.4.2.7. Set the nitrogen flow at 25 mL/min. The nitrogen in the Erlenmeyer flask should bubble at about five (5) bubbles per second. Purge the oxygen from the system for about 15 minutes.

11.3.4.2.8. Turn on the magnetic stirrer. Set the stirring bar to spin as fast as possible. The fluid should form a vortex. If not, the distillation will exhibit poor recovery. Add all the HCl from the dropping funnel to the flask.

11.3.4.2.9. Heat the water bath to the boiling point (100°C). the sample may or may not be boiling. Allow the purged distillation to proceed for 90 minutes at 100°C.

\*Note: Watch water dishes so they don't go dry!

#### 11.4. Sample Analysis

##### 11.4.1. Summary

11.4.1.1. Excess iodine is added to a sample and under acidic conditions is back titrated with sodium thiosulfate from a blue to a clear color.

##### 11.4.2. Sample Analysis Procedure

11.4.2.1. Method 376.1: Place the filter pad from the filtered water sample in a 500 mL Erlenmeyer flask. At this time, add any spiking solutions if necessary for the filtered water samples. Add 20.0 mL .028 N Iodine solution and 1-2 mL 1:1 HCl solution (watch for fumes). Check the pH prior to titration to make sure it is less than 2. If the pH is not <2, add additional acid. Add 250 mL reagent water only for the water samples and mix. Add 1 squirt (1-2 mL) of starch indicator and mix. Titrate from blue to clear

with .025 N sodium thiosulfate titrant. Record the amount of titrant used on the analytical logsheet.

Note: Some matrices may be turbid or colored and the color change from blue to clear may not be easily seen. In this case, look for a shade change.

11.4.2.2. 9030B/9034 Titrations: Titrate the scrubber solution in the Erlenmeyer. Add 20 mL of 0.028 N iodine solution under the liquid level in the Erlenmeyer, and 1-2 mL of 1:1 HCl solution. (More HCl is needed for insoluble 9030B.) Check the pH prior to titration to make sure it is less than 2. Add 1 squirt of starch indicator and mix. Titrate from blue to clear with .025 N sodium thiosulfate titrant.

\*Note: Because of the buffer and the formaldehyde, the titration will take between 15 – 30 minutes. Add titrant slowly to allow the reaction to take place. Overtitrating could be an issue if the titrant is added too quickly.

11.4.2.3. After adding 20 mL iodine, the color should be orange/red. If the color remains yellow, add additional 10 mL aliquots until the orange/red color persists (adjust the calculation accordingly). If the sample requires more than 60 mL of iodine, the sample must be re-prepped at a smaller dilution. The iodine should turn a yellow-orange color when added to the sample after the addition of the reagent water. If it does not, the sample may be high in sulfide and less sample should be used, or more iodine for 9030B.

## 11.5. Analytical Documentation

11.5.1. Record all analytical information in the analytical logbook/logsheet, which may be in electronic format, including the analytical data from standards, blanks, LCSs, MS/MSDs, and any corrective actions or modifications to the method.

11.5.2. All standards are logged into a department standard logbook. All standards are assigned a unique number for identification. Logbooks are reviewed by the supervisor or designee.

11.5.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

11.5.4. Sample results and associated QC are entered into the LIMs after final technical review.

---

## 12. DATA ANALYSIS AND CALCULATIONS

### 12.1. Calculations

$$\text{Sulfide, mg / L or mg / kg} = \frac{[(A \times B) - (C \times D)] \times 16,000}{\text{mL or g of sample used}}$$

*Where:*

*A* = *mL of iodine solution*

*B* = *Normality of iodine solution*

*C* = *mL of sodium thiosulfate titrant*

*D* = *Normality of sodium thiosulfate titrant*

$$\text{Sulfide, mg / L or mg / kg} = \frac{(20 - \text{mL titrant}) \times 400}{\text{mL or g of sample used}}$$

#### 12.1.1.

$$\text{LCS \% Recovery} = \frac{\text{mg/L (from 12.1.1)}}{20 (\text{true})} \times 100$$

**Note:** The true value of the standard is determined daily.

#### 12.1.2.

$$\text{MS/MSD \% Recovery} = \frac{A - B}{20 (\text{waters}) \text{ or } 1000 (\text{solids})} \times 100$$

*Where:*

*A* = *(20 - mL titrant for MS/MSD) x 400*

*B* = *Concentration from 12.1.1 x mL or g of sample used*

**13. METHOD PERFORMANCE**

13.1. Each laboratory must have initial demonstration of performance data on file and corresponding method detection limit files.

13.2. Training Qualifications:

13.2.1. The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

13.2.2. Method validation information (where applicable) in the form of laboratory demonstrations of capabilities is maintained for this method in the laboratory QA files.

**14. POLLUTION PREVENTION**

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

**15. WASTE MANAGEMENT**

15.1. All aqueous sample preparations can be rinsed down the drain with copious amounts of water.

15.2. Solvent waste must be disposed of in clearly labeled waste cans.

15.3. Acid waste must be collected in clearly labeled acid waste containers.

15.4. Solid materials (gloves, soiled paper products, etc.) are placed in the solid debris container. Do not put liquids in the solid waste container.

15.5. Refer to the Laboratory Sample and Waste Disposal plan.

15.6. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL North Canton. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by an annual refresher training.

**16. REFERENCES**

---

## 16.1. References

- 16.1.1. SW846, Test Methods of Evaluating Solid Waste, Third Edition, Sulfide, Method 9030B.
- 16.1.2. EPA 600, Methods for Chemical Analysis of Waters and Wastes, Sulfide (Titrimetric, Iodine), Method 376.1
- 16.1.3. Standard Methods Eighteenth Edition, Sulfide Method 4500-5<sup>2</sup>-E
- 16.1.4. STL North Canton Laboratory Quality Manual (LQM), current version.
- 16.1.5. STL Quality Management Plan (QMP), current version.
- 16.1.6. SW846, Test Methods of Evaluating Solid Waste, Third Edition, Titrimetric Procedure for Acid-soluble and Acid-insoluble Sulfides, Method 9034.

## 16.2. Associated SOPs and Policies, latest version

- 16.2.1. QA Policy, QA-003
- 16.2.2. Glassware Washing, NC-QA-0014
- 16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018
- 16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021
- 16.2.5. Navy/Army SOP, NC-QA-0016

## 17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)

### 17.1. Reporting limits

- 17.1.1. The lower reporting limits (RL) are 1 mg/L for waters and 50 mg/kg for solids.
- 17.1.2. If samples require dilution or smaller volumes than specified in this method, the RL will be elevated.

### 17.2. Method Deviations



8771170

SULFIDE

SOP No. NC-WC-0060

Revision No. 3

Revision Date: 06/17/03

Page 19 of 19

---

17.2.1. The lab does not perform the distillation procedure described in Method SW846 9030A.

17.2.2. The laboratory uses one collection flask instead of two.

17.2.3. The laboratory does not titrate the sample in the original container as specified in Method 376.1

## LABORATORY STANDARD OPERATING PROCEDURES

Once printed, this is considered an uncontrolled document

SOP No.	Revision No.	Date	Page
AWC-VFA-20	1	February 4, 2004	1 of 13

TITLE: Volatile Fatty Acids

Supersedes: Revision 0

REVIEWED & APPROVED BY:	Signature	Date
Christopher Oprandi, Laboratory Director	<i>Christopher Oprandi</i>	2/5/04
Verl D. Preston, Quality Manager	<i>Verl D. Preston</i>	2/5/04
Peggy Gray-Erdmann, Supervisor	<i>Peggy Gray-Erdmann</i>	2/10/04

## Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories, Inc. (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it upon STL's request and not to reproduce, copy, lend or otherwise disclose its contents, directly or indirectly, and not to use it for any purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to those parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY SEVERN TRENT LABORATORIES IS PROTECTED BY STATE AND FEDERAL LAWS OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR, THE FOLLOWING NOTICE SHALL APPLY:

©COPYRIGHT 2004 SEVERN TRENT LABORATORIES, INC., INC. ALL RIGHTS RESERVED.

## 1.0 IDENTIFICATION OF TEST METHODS

- 1.1. Volatile Fatty Acids are determined using Dionex proprietary method. This method is a modification to Standards Methods for the Examination of Water and Waste Water 20<sup>th</sup> ed., Method 5560 Organic and Volatile Acids.

## 2.0 APPLICABLE MATRIX

- 2.1. This method is applicable to surface water, groundwater, wastewater, drinking waters and soils.

## 3.0 REPORTING LIMIT

- 3.1. The reporting limit for each acid is listed below:

3.1.1. Acetic Acid – 1 mg/L

3.1.2. Propionic Acid – 1 mg/L

3.1.3. Butyric Acid – 1 mg/L

3.1.4. Lactic Acid – 1 mg/L

3.1.5. Formic Acid – 1 mg/L

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	Page 2 of 13
AWC-VFA-20	1	February 4, 2004	

**TITLE:** Volatile Fatty Acids

**Supercedes:** Revision 0

3.1.6. Pyruvic Acid – 1 mg/L

- 3.2. MDLs are calculated every year in accordance with method specification and kept on file with the QA department.

#### 4.0 SCOPE AND APPLICATION

- 4.1. Ion Chromatography provides a single instrumental technique that may be used for the measurement in environmental samples of the common Volatile Fatty acids or Organic acids. Fatty acids are low-molecular weight carboxylic acids. Identification and quantification is performed by Ion chromatography. Separation is accomplished by an ion-exclusion column. Peak sensitivity and lower reporting limits are obtained by using an eluent suppressor. The acids are then identified by conductivity detection.

#### 5.0 SUMMARY OF TEST METHOD

- 5.1. A filtered aqueous sample is injected into an ion chromatograph with the use of an automated sampler. The sample merges with an eluent stream and is pumped through the system. The ion exchanger separates the acids of interest. Ions are separated based on their affinity for the exchange sites of the resin. The separated anions in their acid form are measured using an electrical conductivity cell. Anions are identified based on their retention times compared to known standards. Quantitation is accomplished by measuring the peak area and comparing it to a calibration curve generated from known standards.

#### 6.0 DEFINITIONS

- 6.1. Standard definitions can be found in section 3.0 of the STL Buffalo Laboratory Quality Manual.
- 6.2. VFA: Volatile Fatty Acid

#### 7.0 INTERFERENCES

- 7.1. Interferences can be caused by substances with retention times that are similar to and overlap those of the acid of interest. Acids of high concentrations can interfere with the peak resolution of an adjacent acid. Diluting the sample can minimize overlap.
- 7.2. Method interferences may be caused by contaminants in the reagent water, reagents, glassware and other sample processing apparatus that lead to discrete artifacts or an elevated baseline in the ion chromatograms.

<b>SOP No.</b> AWC-VFA-20	<b>Revision No.</b> 1	<b>Date</b> February 4, 2004	<b>Page 3 of 13</b>
------------------------------	--------------------------	---------------------------------	---------------------

**TITLE: Volatile Fatty Acids**

**Supersedes: Revision 0**

- 7.3.1 All samples must be filtered through a 20um filter before injection. If particles contaminate the guard or analytical columns, follow the manufacturer's suggestions for cleaning, or simply replace the column.

## 8.0 SAFETY

- 8.1. Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual and this document.

### 8.2. SPECIFIC SAFETY CONCERNS OR REQUIREMENTS

- 8.2.1. This method uses weak carboxylic acids and heptafluorobutyric acid. All acids will be poured into water.

- 8.2.2. Exercise caution when using syringes with attached filter assemblies. Application of excessive force has, upon occasion, caused a filter disc to burst during the process.

### 8.3. PRIMARY MATERIALS USED

- 8.3.1. The following is a list of the materials used in this method, which have a serious or significant hazard rating. **NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table.** A complete list of materials used in the method can be found in the reagents and materials section. Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

Material (1)	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
Heptafluor-Butyric acid	Corrosive	NA	Causes irritation to the respiratory tract, skin and eyes. Symptoms may include coughing, shortness of breath. Symptoms include redness, itching, and pain.
Fatty Acid Custom Mix 1000ppm in water	Corrosive	NA	May irritate eyes and/or skin. Irritates Respiratory tract. Low blood pressure. Dermatitis. Pulmonary edema. Lung damage.
1 – Always add acid to water to prevent violent reactions.			
2 – Exposure limit refers to the OSHA regulatory exposure limit.			

## 9.0 EQUIPMENT AND SUPPLIES

- 9.1. Ion chromatograph complete with all required accessories:

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	Page
AWC-VFA-20	1	February 4, 2004	4 of 13

**TITLE:** Volatile Fatty Acids

**Supersedes:** Revision 0

9.1.1. ICE-AS6 separator column capable of resolving weak organic acids.

9.1.2. AMMS-ICE suppressor.

9.1.3. Conductivity detector with temperature control and separate working and reference electrodes.

9.1.4. Pump able to deliver 1.2 ml/min of constant flow rate.

9.1.5. Data collection and analysis system.

9.1.6. Automated sampler.

9.2. Various laboratory glassware such as Class A graduated cylinders, syringes, volumetric flasks and pipettes.

9.3. 10 ml syringes and 0.2 um syringe filters for colored samples

9.4. Analytical balance, capable of weighing to the nearest 0.0001g.

9.5. Filter caps for clean samples purchased from Dionex

9.6. 5 ml sample vials purchased from Dionex

#### 10.0 REAGENTS AND STANDARDS

10.1. Sample bottles: Glass or polyethylene bottles of sufficient volume to allow replicate analyses of anions of interest.

10.2. Reagent water: Distilled or deionized water free of the anions of interest. Water should contain particles no larger than 0.20 microns.

10.3. Regenerant Concentrate (TetraButylAmmonium Hydroxide TBAOH) from VWR catalog# (JTV365-7) 0.4M water solution.

10.3.1 Regenerant solution (10mM): Add 50ml of 0.4M TBAOH into 2L reagent water.

10.4. Eluent Concentrate (Heptafluorobutyric acid): 99% Heptafluorobutyric acid is purchased from Aldrich. Catalog # 16,419-4

10.4.1. Eluent Solution (1.0mM.): Weigh 0.21g of the Heptafluorobutyric acid (10.4) to 1 liter with reagent water or 250ul into 2 liter

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No. AWC-VFA-20	Revision No. 1	Date February 4, 2004	Page 5 of 13
-----------------------	-------------------	--------------------------	--------------

**TITLE: Volatile Fatty Acids**

**Supersedes: Revision 0**

- 10.5. Multi analyte Ion Chromatography Custom Standards purchased from Restek & Supelco Scientific. Custom Mix contains 1000mg/l of the following acids: Acetic, Propionic, Butyric, Lactic, Formic, and Pyruvic. The standard source from Supelco is used for the calibration curve (10.6) and the standard source from Restek is used for the ICV/CCV/LCS Solutions and the matrix spikes (10.7, 10.8, 10.9, 10.10).

- 10.6. Calibration standards: all are made from dilutions of the Multi Element IC Standards in reagent water.

- 10.6.1. Prepare the calibration standards for a 5-point curve by measuring the following volumes into a 25 ml Class A volumetric. Bring to the final volume of 25 ml with reagent water.

	Level 1	Level 2	Level 3	Level 4	level 5
Stock solution (1000 mg/L)	25ul	50ul	125ul	250ul	1.25ml
Final Volume	25ml	25ml	25ml	25ml	25ml

- 10.6.2 The final concentrations of each anion in the 5 calibration points are summarized below.

	Level 1 (mg/L)	Level 2 (mg/L)	Level 3 (mg/L)	Level 4 (mg/L)	Level 5 (mg/L)
Formic	1	2	5	10	50
Propionic	1	2	5	10	50
Butyric	1	2	5	10	50
Lactic	1	2	5	10	50
Formic	1	2	5	10	50
Pyruvic	1	2	5	10	50

- 10.7. ICV/CCV/LCS, MS and SD (MSD) solution: 500ul of the Multi Element IC Standard diluted to 100ml with reagent water. The final concentration for each acid in the solution is as follows:

## ICV &amp; CCV/LCS

Acetic	5 mg/L
Propionic	5 mg/L
Butyric	5 mg/L
Lactic	5 mg/L
Formic	5 mg/L
Pyruvic	5 mg/L

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	Page
AWC-VFA-20	1	February 4, 2004	6 of 13

**TITLE:** Volatile Fatty Acids

**Supercedes:** Revision 0

## 11.0 SAMPLE COLLECTION, PRESERVATION, SHIPMENT AND STORAGE

- 11.1. Samples should be shipped and stored in 40 mL glass amber vials at  $4 \pm 2$  degrees C. Samples should be analyzed for VFA within 28 days of collection.
- 11.2. Soil will follow the same preservation and holding times as the water samples, starting from the time of extraction. Soils should be collected in 4 oz. or 8 oz. wide mouth amber jars.

## 12.0 QUALITY CONTROL

- 12.1. Before analyzing samples, the laboratory must establish a method detection limit (MDL). The MDL is repeated every year.
- 12.2. Each group of sample analyses must be bracketed by an acceptable calibration verification sample and calibration blank. All quality control data should be maintained and available for easy reference or inspection.
- 12.3. Initial and Continuing Calibration Blank (ICB, CCB): To determine freedom from contamination, prepare one calibration blank (ICB) at the beginning of the analytical procedure and another (CCB) after every ten samples and at the end of the analytical procedure. The blank consists of 5 ml reagent water that gets the same treatment as the samples and standards. The blanks must be free of the analytes of concern at levels less than the STL Buffalo quantitation limit.
  - 12.3.1. All blanks associated with USACE samples should be less than half the STL Buffalo quantitation limit for each anion.
- 12.4. Initial and Continuing Calibration Verification/Laboratory Control Sample (ICV/CCV/LCS): Prepare an ICV at the beginning of the analytical procedure and additional CCV/LCS after every ten samples and again at the end of the procedure. The recovery of the ICV/CCV/LCSs must be within 80-120% of the true value.
- 12.5. Sample Duplicate: Analyze either a Matrix Duplicate (MD) or a Matrix Spike Duplicate (MSD) with every batch of twenty or fewer samples. Acceptable RPD between replicate analyses should be less than 20%.
- 12.6. Matrix Spikes (MS) are to be run with every batch of 20 or fewer samples. Deviations may occur due to specific client, state, or protocol requirements. Spike 10ml of sample with 50ul of Multi-analyte Stock solution (1000mg/L) section 10.5.

## 13.0 CALIBRATION AND STANDARDIZATION

- 13.1. Prepare standard 5-point curve by plotting instrument response against concentration values.

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	
AWC-VFA-20	1	February 4, 2004	Page 7 of 13

**TITLE: Volatile Fatty Acids**

**Supersedes: Revision 0**

- 13.1.1. Generate a linear regression curve. Do not force through zero and do not average in the origin.
- 13.1.2. A calibration curve may be fitted to the calibration solution concentration/response data using the manufacturer's software.
- 13.1.3. Acceptance criteria for the calibration curve is a correlation coefficient (R value)  $\geq 0.995$ . If the R-value is less than 0.995, the calibration standards must be remade and a new curve analyzed.
- 13.1.4. New calibration curves must be run every three months or if the instrument falls out of calibration whichever is sooner
- 13.2. Initial Calibration Verification Solution prepared from a different standard source is analyzed immediately after the calibration curve to verify the accuracy of the curve. The recovery of the ICV must be within 80-120%.

#### 14.0 PROCEDURE

##### 14.1. System Equilibrium:

- 14.1.1. Set up the ion chromatograph as specified in the manufacturer's instructions.
- 14.1.2 Turn on and prime the pump.
- 14.1.3 Adjust the eluent flow rate to  $1.2 \pm 0.1$  ml/min. Adjust Regenerant pressure to 5PSI. Adjust restrictor at the end of the regenerant waste line to allow for 3-5ml/min. flow.
- 14.1.4 Allow the system to come to equilibrium (15-20 minutes). A stable baseline indicates system equilibrium.

##### 14.2. Sample analysis:

- 14.2.1. For dirty samples filter sample through a pre-washed 0.2um pore diameter membrane filter. If sample is clean use 20um filter autosampler caps.
- 14.2.2. Fill autosampler vials with the sample to the fill line marked on the vial body (approximately 5 ml). Place vial cap into vial.
- 14.2.3. Place the filled vial into the sampler cassette and fully insert the cap using the insertion tool.
- 14.2.4. Place the filled cassettes into the automated sampler and start the run. Set sample run time to 64 minutes.



## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	Page
AWC-VFA-20	1	February 4, 2004	8 of 13

**TITLE:** Volatile Fatty Acids

**Supersedes:** Revision 0

14.2.5. Check data for any needed dilutions and calculate percent recovery of check standards and sample spikes. Any data from samples that were diluted will have to be multiplied by the dilution factor before reporting.

**14.3. Column cleanup procedure**

14.3.1. Disconnect the suppressor from the analytical column. Reverse the order of the guard and analytical column in the eluent flow path. Double check that the eluent flows in the direction designated on each of the column labels.

14.3.1.1. **CAUTION:** When cleaning an analytical column and a guard column in series, ensure that the guard column is placed after the analytical column in the eluent flow path. Contaminants that have accumulated on the guard column can be eluted onto the analytical column and irreversibly damage it. If in doubt, clean each column separately.

14.3.2. Set the pump flow rate to 1.0 ml/min for the ICE-AS6 analytical.

14.3.3. Rinse the column for 15 minutes with reagent water before pumping the 0.1M oxalic acid over the column. Acetone can also be used to remove organics. Start with a 5% solution and increase up to 15%, but do not exceed 15%.

14.3.4. Pump the cleanup solution through the column for at least 60 minutes.

14.3.5. Rinse the column for 30 minutes with reagent water before pumping eluent over the column.

14.3.6. Equilibrate the columns with eluent before resuming normal operations for at least 30 minutes.

14.3.7. Reconnect the ASRS-Ultra to the analytical column.

14.4 Retention time (migration time) is the expected time retention time or migration time in minutes for the component. If the retention time is unknown, enter any number greater than zero. The correct retention time can be determined later from the first calibration run, and the component table then updated. In subsequent calibrations, PeakNet will automatically update the retention time. The Update Retention Time setting must be selected in the calibration Parameters dialog box.

**15.0 CALCULATIONS**

15.1. Using the computer and software packages, prepare a linear regression calibration curve for each analyte by plotting instrument response against standard concentration. Compute sample concentration by comparing sample response with the standard curve.

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	Page
AWC-VFA-20	1	February 4, 2004	9 of 13

**TITLE:** Volatile Fatty Acids

**Supersedes:** Revision 0

The response factor produced from a linear equation best fits the detector's response. The equation used is shown below.

$$Y = K0 + K1 \times X$$

At least four points are needed to fit the equation: thus, the calibration must have at least four levels for all components.

The following values, used to calculate component amount, are determined automatically by the Method Editor and cannot be edited.

**X** = area

**K0** indicates the Y intercept of the calibration curve.

**K1** is the coefficient for the first-degree variable. When the fit type is linear, K1 indicates the slope of the calibration curve for the selected calibration level.

The equation for the calibration curve fit used to calculate the component amount is displayed at the bottom of the replicate page. The  $r^2$  value (Coefficient of Determination) for the component is shown at the bottom of the replicate page.

- 15.2. The analyst corrects the results for and dilution factors:

$$X_f = X_j \times \text{Dilution Factor}$$

Where:

$X_f$  = Final sample concentration

$X_j$  = calculated concentration of sample at instrument

- 15.3. Report only those values that are less than the highest calibration standards. Samples exceeding the highest standard should be diluted and reanalyzed.

#### 16.0 METHOD PERFORMANCE

- 16.1. The method detection limit (MDL) is to be performed every year in accordance with the specifications in 40 CFR 136, appendix B, and must demonstrate the ability to quantitate at or below the reporting limit for each acid. The current MDL is on file with the department supervisor and the QA Department.

#### 17.0 DATA ASSESSMENT AND ACCEPTANCE CRITERIA FOR QUALITY CONTROL MEASURES

- 17.1. Obtained ICV and CCV/LCS values must be within 80-120% of the true value.
- 17.2. Acceptance limits for sample spike recovery are based on the historical data and are statistically derived annually. They are maintained in the laboratory LIMS system. If the

## LABORATORY STANDARD OPERATING PROCEDURES

<b>SOP No.</b> AWC-VFA-20	<b>Revision No.</b> 1	<b>Date</b> February 4, 2004	<b>Page</b> 10 of 13
------------------------------	--------------------------	---------------------------------	----------------------

**TITLE:** Volatile Fatty Acids

**Supercedes:** Revision 0

lab calculated limits are wider than the method limits, the method limits of 80-120% are used for evaluation of sample spike acceptance.

17.3. Sample duplicates are required to have a calculated RPD  $\leq 20$ .

17.4. ICB and CCB values must be less than the STL quantitation limit.

17.4.1. All blanks associated with USACE samples should be less than  $\frac{1}{2}$  the STL Buffalo quantitation limit for each anion.

## 18.0 CORRECTIVE ACTIONS FOR OUT-OF-CONTROL DATA

18.1. If acceptance criteria are exceeded for any QC element, all related samples and check standards must be repeated.

## 19.0 CONTINGENCIES FOR HANDLING OUT-OF-CONTROL OR UNACCEPTABLE DATA

19.1 If acceptable data can not be obtained, a Job Exception Form is to be filled out and turned in to the appropriate project manager in order to notify of the client.

19.2 Historical data review may be used to evaluate sample results.

## 20.0 WASTE MANAGEMENT/ POLLUTION PREVENTION

20.1 All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."

20.2 Waste Streams Produced by the Method -The following waste streams are produced when this method is carried out.

20.2.1 Alkaline and/or acidic waste generated by the analysis. Dispose of this waste in the "A" waste container.

20.2.2 Contaminated plastic materials such as IC syringes, filters, caps and vials utilized for sample preparation. All plastic materials should be disposed of in the recycling containers located throughout the lab.

## 21.0 REFERENCE

21.1. Method 300.0, "Determination of Inorganic Anions by Ion Chromatography", Environmental Monitoring Systems Laboratory, Office of Research and Development, U.S. EPA, Cincinnati, Ohio, Revision 2.1, August 1993

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	
AWC-VFA-20	1	February 4, 2004	Page 11 of 13

**TITLE: Volatile Fatty Acids**

**Supersedes: Revision 0**

21.2. Method 5560, Standard Methods for the Examination of Water and Wastewater, 20th Edition,

21.3. Dionex application note #45, "Fatty Acid Analysis"

**22.0 TABLES, DIAGRAMS, FLOWCHARTS AND VALIDATION DATA**

22.1. Analytical Run sequence

22.2. Wet Chemistry Batch Summary & Data Review Checklist

**23.0 CHANGES FROM PREVIOUS REVISION**

23.1 Changed the reporting limit of Butyric acid from 2 mg/L to 1 mg/L

23.2 Changed the concentration of the regenerant solution from 4 mM to 10 mM

23.3 Changed the concentration of the eluent solution from 0.6 mM to 1.0 mM

23.4 Removed MSD reference and added specific recipe for MS in section 12.6

23.5 Added 64 minute run time to section 14.2.4

## LABORATORY STANDARD OPERATING PROCEDURES

SOP No.	Revision No.	Date	
AWC-VFA-20	1	February 4, 2004	Page 12 of 13

**TITLE:** Volatile Fatty Acids

**Supersedes:** Revision 0

## 22.1 Analytical Run Sequence

LCS/CCV

ICB

Sample

Sample

Sample

Sample

Sample

Sample

Sample

Sample

Sample

Sample

CCV

CCB

Sample

Sample

Sample

Sample

Sample

Sample

Sample

Sample

Sample dup

Sample Spike

CCV

CCB

## LABORATORY STANDARD OPERATING PROCEDURES

<b>SOP No.</b> AWC-VFA-20	<b>Revision No.</b> 1	<b>Date</b> February 4, 2004	<b>Page 13 of 13</b>
------------------------------	--------------------------	---------------------------------	----------------------

TITLE: Volatile Fatty Acids

Supercedes: Revision 0

## 22.2 Wet Chemistry Batch Summary &amp; Data Review Summary

## WET CHEMISTRY BATCH SUMMARY

Parameter \_\_\_\_\_ Method \_\_\_\_\_ Batch # \_\_\_\_\_

Comment #	Comment
1	NA
2	Sample(s) was diluted for matrix interference.
3	Sample(s) was diluted for excessive foaming/
4	Sample(s) was diluted for turbidity.
5	NA
6	NA
7	NA
8	Sample(s) was diluted for high concentration of target analyte
9	Sample(s) was diluted for turbidity.
10	Sample(s) was diluted for color.
11	There was insufficient volume for a lower dilution.
12	Sample(s) was diluted for viscosity.
13	Sample(s) was diluted for other reason (detail required)
14	Sample(s) required re-run to verify result.
15	Sample(s) requires re-run to verify deviation from historical result.
16	Sample(s) requires re-run for CCB failure.
17	Sample(s) affected by elevated CCB are greater than 10x detection limit.
18	Sample was colored.
19	Sample(s) was received outside of Holding Times.
20	Sample(s) contained a high amount of settleable material.
21	Sample(s) contained a high amount of suspended material.
22	Sample(s) were centrifuged for turbidity.
23	There was insufficient volume for analysis of sample at method required volume.
24	There was insufficient volume for re-analysis of the sample(s).
25	There was insufficient volume for dilution of the sample(s).
26	There was insufficient volume for Dup/Spk.
27	Sample(s) was cloudy
28	See accompanying Job Exception Report.

## Comments and Corrective actions

# \_\_\_\_\_ Sample(s) \_\_\_\_\_

# \_\_\_\_\_ Sample(s) \_\_\_\_\_

# \_\_\_\_\_ Sample(s) \_\_\_\_\_

# \_\_\_\_\_ Sample(s) \_\_\_\_\_

CCV/CCB Compliant? NA \_\_\_\_\_ YES \_\_\_\_\_ NO \_\_\_\_\_ (see reason below)

Other \_\_\_\_\_

\_\_\_\_\_

Technician \_\_\_\_\_ Date \_\_\_\_\_

2<sup>nd</sup> Review \_\_\_\_\_ Date \_\_\_\_\_

Number of Reanalysis for this batch: \_\_\_\_\_

Review \_\_\_\_\_ Date \_\_\_\_\_

WC Summary Rev 1 7-2003

8771184

Controlled Copy

Copy No. \_\_\_\_\_

Implementation Date: 10/18/00

SOP No. CORP-IP-0004NC

Revision No. 1.1Revision Date: 10/10/00

Page: 1 of 48

**STL STANDARD OPERATING PROCEDURE****TITLE: TOXICITY CHARACTERISTIC LEACHING PROCEDURE****AND SYNTHETIC PRECIPITATION LEACHING PROCEDURE****(SUPERSEDES: REVISION 1, REVISION DATE 02/01/00)**

Prepared by:

Mark H. Bure10/17/00  
Date

Reviewed by:

Phil Bure

Technology Specialist

10/17/00  
Date

Approved by:

Galina Cahn

Quality Assurance Manager

10/17/00  
Date

Approved by:

Heidi

Environmental Health and Safety Coordinator

10-17-00  
Date

Approved by:

Christy McQuade

Laboratory Director

10-17-00  
Date

Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

## **Table of Contents**

<b>1. SCOPE AND APPLICATION .....</b>	<b>4</b>
<b>2. SUMMARY OF METHOD.....</b>	<b>5</b>
<b>3. DEFINITIONS.....</b>	<b>6</b>
<b>4. INTERFERENCES .....</b>	<b>6</b>
<b>5. SAFETY.....</b>	<b>7</b>
<b>6. EQUIPMENT AND SUPPLIES .....</b>	<b>8</b>
<b>7. REAGENTS AND STANDARDS.....</b>	<b>9</b>
<b>8. SAMPLE COLLECTION, PRESERVATION AND STORAGE .....</b>	<b>11</b>
<b>TABLE 1 - HOLDING TIMES (DAYS) .....</b>	<b>12</b>
<b>9. QUALITY CONTROL .....</b>	<b>13</b>
<b>10. CALIBRATION AND STANDARDIZATION .....</b>	<b>14</b>
<b>11. PROCEDURE .....</b>	<b>14</b>
<b>11.1. GENERAL COMMENTS.....</b>	<b>14</b>
<b>11.2. PRELIMINARY EVALUATIONS.....</b>	<b>14</b>
<b>TABLE 2 - MINIMUM REQUIRED LEACHATE VOLUMES .....</b>	<b>14</b>
<b>11.3. PROCEDURE: NON-VOLATILE CONSTITUENTS: SEMI-VOLATILES, PESTICIDES, HERBICIDES, METALS (BOTTLE LEACH).....</b>	<b>20</b>
<b>11.4. PROCEDURE: VOLATILE CONSTITUENTS (ZHE).....</b>	<b>22</b>
<b>12. DATA ANALYSIS AND CALCULATIONS.....</b>	<b>26</b>
<b>13. METHOD PERFORMANCE.....</b>	<b>28</b>
<b>14. POLLUTION PREVENTION.....</b>	<b>28</b>
<b>15. WASTE MANAGEMENT.....</b>	<b>29</b>



16. REFERENCES .....	29
17. MISCELLANEOUS (METHOD INTERPRETATIONS, RECORD MANAGEMENT, ETC...).....	32

**LIST OF APPENDICES**

APPENDIX A - TABLES .....	35
APPENDIX B - FIGURES.....	37
APPENDIX C - WORKSHEETS.....	41
APPENDIX D - FLOW CHARTS.....	48

## 1. SCOPE AND APPLICATION

1.1. This SOP describes the application of the Toxicity Characteristic Leaching Procedure (TCLP), SW846 Method 1311. The Toxicity Characteristic (TC) of a waste material is established by determining the levels of 8 metals and 31 organic chemicals in the aqueous leachate of a waste. The TC is one of four criteria in 40 CFR Part 261 to determine whether a solid waste is classified as a hazardous waste. The other three are corrosivity, reactivity and ignitability. The TC Rule utilizes the TCLP method to generate the leachate under controlled conditions that were designed to simulate leaching through a landfill. EPA's "worst case" waste disposal model assumes mismanaged wastes will be exposed to leaching by the acidic fluids generated in municipal landfills. The EPA's model also assumes the acid/base characteristics of the waste will be dominated by the landfill fluids. The TCLP procedure directs the testing laboratory to use a more acidic leaching fluid if the sample is an alkaline waste, again in keeping with the model's assumption that the acid fluids will dominate leaching chemistry over time.

1.2. The specific list of TC analytes and regulatory limits may be found in Appendix A.

**Note:** The list in Appendix A does not include the December 1994 EPA rule for Universal Treatment Standards for Land Disposal Restrictions. Those requirements include 216 specific metallic and organic compounds and, in some cases, lower detection limit requirements (see 40 CFR 268.40). TCLP leachates are part of the new Universal Treatment Standards, but the conventional analytical methods will not necessarily meet the new regulatory limits. Consult with the client and with STL Technical Specialists before establishing the instrumental methods for these regulations.

1.3. This SOP also describes the application of the Synthetic Precipitation Leaching Procedure (SPLP) which was designed to simulate the leaching that would occur if a waste was disposed in a landfill and exposed only to percolating rain water. The procedure is based on SW846 Method 1312. The list of analytes for SPLP may extend beyond the toxicity characteristic compounds shown in Appendix A. With the exception of the use of a modified extraction fluid, the SPLP and TCLP protocols are essentially equivalent. Where slight differences may exist between the SPLP and TCLP they are distinguished within this SOP.

1.4. The procedure is applicable to liquid, solid, and multiphase wastes.

1.5. The results obtained are highly dependent on the pH of the extracting solution, the length of time that the sample is exposed to the extracting solution, the temperature during extraction,

TCLP and SPLP Leaching Procedure

---

and the particle size/surface area of the sample. These parameters must be carefully controlled.

- 1.6. The reporting limits are based on the individual samples as well as the individual analysis techniques. However, the sample is determined to be hazardous if it contains any analyte at levels greater than or equal to the regulatory limits.
- 1.7. If a total analysis of the waste demonstrates that individual analytes are not present in the waste, or that they are present but at such low concentrations that the appropriate regulatory levels could not possibly be exceeded, the procedure need not be run. If the total analysis results indicate that TCLP is not required, the decision to cease TCLP analysis should be remanded to the client.
- 1.8. If an analysis of any one of the liquid fractions of the procedure leachate indicates that a regulated compound is present at such a high concentration that, even after accounting for dilution from the other fractions of the leachate, the concentration would be equal to or above the regulatory level for that compound, then the waste is hazardous and it may not be necessary to analyze the remaining fractions of the leachate. However, the remaining analyses should not be terminated without the approval of the client.
- 1.9. Volatile organic analysis of the leachate obtained using a bottle extraction, normally used for extractable organics and metals, can be used to demonstrate that a waste is hazardous, but only the ZHE option can be used to demonstrate that the concentration of volatile organic compounds is below regulatory limits due to potential analyte loss into the headspace during the bottle extraction.

## 2. SUMMARY OF METHOD

- 2.1. For liquid wastes that contain less than 0.5% dry solid material, the waste, after filtration through 0.6 to 0.8  $\mu\text{m}$  glass fiber filter, is defined as the TCLP leachate.
- 2.2. For wastes containing greater than or equal to 0.5% solids, the liquid, if any, is separated from the solids and stored for later analysis or recombination with the leachate. The particle size of the remaining solid phase is reduced, if necessary. The solid phase is extracted with an amount of extraction fluid equal to 20 times the weight of the solid phase. For TCLP, the extraction fluid employed for extraction of non-volatile analytes is a function of the alkalinity of the solid phase of the waste. For SPLP, the extraction fluid employed is a function of the region of the country where the sample site is located if the sample is a soil. If the sample is a waste or wastewater the extraction fluid employed is a pH 4.2 solution. Two leachates may be generated: a) one for analysis of non-volatile constituents (semi-volatile organics, pesticides, herbicides and metals and/or b) one from a Zero Headspace Extractor (ZHE) for

---

TCLP and SPLP Leaching Procedure

---

analysis of volatile organic constituents. Following extraction, the liquid leachate is separated from the solid phase by filtration through a 0.6 to 0.8  $\mu\text{m}$  fiber filter.

- 2.3. If compatible (i.e., multiple phases will not form on combination), the initial liquid filtrate of the waste is added to the liquid leachate and these are prepared and analyzed together. If incompatible, the liquids are analyzed separately and the results are mathematically combined to yield a volume-weighted average concentration.

### 3. DEFINITIONS

- 3.1. "Leachate" is used to refer to the solutions generated from these procedures (TCLP, SPLP, deionized water leach).
- 3.2. "Wet Solids" is that fraction of a waste sample from which no liquid may be forced out by pressure filtration.

### 4. INTERFERENCES

- 4.1. Oily wastes may present unusual filtration and drying problems. If requested by the client and as recommended by EPA (see Figure 3), oily wastes can be assumed to be 100% liquid and analysis for total concentrations of contaminants will be performed. This applies specifically to samples containing viscous non-aqueous liquids that would be difficult to filter. Alternately, the oil may be subjected to pressure filtration. The portion that passes through the filter will be prepared and analyzed separately as an organic waste. The "wet solid" portion that remains behind on the filter will be subjected to leaching, prepared and analyzed. The results will then be mathematically combined.
- 4.2. Wastes containing free organic liquids (e.g., oil, paint thinner, fuel) usually require dilution prior to analysis to address the matrix interferences. In most instances this results in reporting limits elevated above the TCLP regulatory limits.
- 4.3. Solvents, reagents, glassware and other sample processing hardware may yield artifacts and/or interferences to sample analysis. All these materials must be demonstrated to be free from interferences under the conditions of the analysis by analyzing method blanks as described in Section 9 and the individual determinative SOPs.
- 4.4. Glassware and equipment contamination may result in analyte degradation. Soap residue on glassware and equipment may contribute to this. All glassware and equipment should be rinsed very carefully to avoid this problem.

TCLP and SPLP Leaching Procedure

---

- 4.5. Phthalates may be eliminated by proper glassware cleanup and by avoiding plastics. Only glass, Teflon or Type 316 stainless steel tumblers may be used for leachates to be analyzed for organics. Plastic tumblers may be used for leachates to be analyzed for the metals.
- 4.6. Overexposure of the sample to the environment will result in the loss of volatile components.
- 4.7. Potential interferences that may be encountered during analysis are discussed in the individual analytical methods.

**5. SAFETY**

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL associates.
- 5.2. Eye protection that satisfies ANSI Z87.1 (as per the Chemical Hygiene Plan), laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory. The following specific hazards are known:
  - 5.3.1. Chemicals that have been classified as carcinogens, or potential carcinogens, under OSHA include:  
Methylene chloride
  - 5.3.2. Chemicals known to be flammable are:  
Methanol
  - 5.3.3. The following materials are known to be corrosive:  
Hydrochloric acid, nitric acid, sulfuric acid, acetic acid, sodium hydroxide
  - 5.3.4. The following materials are known to be oxidizing agents:  
Nitric Acid.
- 5.4. Gas pressurized equipment is employed in this procedure. Be sure all valves and gauges are operating properly and that none of the equipment, especially tubing, is over-pressurized. CAUTION: Do not open equipment that has been pressurized until it has returned to ambient pressure.

TCLP and SPLP Leaching Procedure

---

- 5.5. A rotary agitation apparatus is used in this procedure. Certain samples may break the glass jars used in the procedure. For these samples, extra caution, including plastic or polyethylene overwraps of the glass jar, may be necessary. Turning the jar or bottle sideways rather than tumbling end over end may also reduce the chance of breakage. If sideways tumbling is used, note this change in the logbook comment section.
- 5.6. Secure tumbler and extraction apparatus before starting rotary agitation apparatus.
- 5.7. During sample rotation, pressure may build up inside the bottle. Periodic venting of the bottle will relieve pressure.
- 5.8. Exposure to hazardous chemicals must be maintained as low as reasonably achievable, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.9. The preparation of standards and reagents and glassware cleaning procedures that involve solvents such as methylene chloride will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.10. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL associate. The situation must be reported immediately to a laboratory supervisor.
- 5.11. Due to the potential for ignition, flammability or production of noxious fumes, do not attempt to dry non-aqueous liquid samples in an oven. Use extended drying in a ventilation hood.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Extraction vessels
  - 6.1.1. For volatile analytes - zero-headspace extraction (ZHE) vessel, gas-pressure actuated, Millipore YT3009OHW or equivalent (see Figure 2).
  - 6.1.2. For metals - either borosilicate glass jars (2.5 L, with Teflon lid inserts) or 2.5 L HDPE (Nalgene or equivalent) bottles may be used.
  - 6.1.3. For non-volatile organics - only borosilicate glass may be used.
- 6.2. Vacuum filtration apparatus and stainless steel pressure filtration apparatus (142 mm diameter), capable of 0 - 50 psi.

TCLP and SPLP Leaching Procedure

---

- 6.3. Borosilicate glass fiber filters, 0.6 - 0.8  $\mu\text{m}$  (Whatman GF/F 14.2 cm, 9.0 cm, 4.7 cm, 0.7  $\mu\text{m}$  or equivalent). When analyzing for metals, wash the filters with 1 N nitric acid and de-ionized water prior to use, or purchase pre-washed filters. Glass fiber filters are fragile and should be handled with care.
- 6.4. Rotary agitation apparatus, multiple-vessel, Associated Design and Manufacturing Company 3740-6 or equivalent (see Figure 1). The apparatus must be capable of rotating the extraction vessel in an end-over-end fashion at  $30 \pm 2$  rpm.
- 6.5. ZHE Extract Collection Device: Gas-tight syringes, 50 or 100 mL capacity, Hamilton 0158330 or equivalent.
- 6.6. Top loading balance, capable of  $0 - 4000 \pm 0.01$  g (all measurements are to be within  $\pm 0.1$  grams).
- 6.7. pH meter and probe capable of reading to the nearest 0.01 unit, and with automatic temperature compensation.
- 6.8. pH probes.
- 6.9. Magnetic stirrer/hotplate and stirring bars.
- 6.10. VOA vials, 40 mL, with caps and septa.
- 6.11. Glass bottles, 1 liter, with Teflon lid-inserts.
- 6.12. Nalgene plastic bottles or equivalent, 1 liter.
- 6.13. Miscellaneous laboratory glassware and equipment.

## 7. REAGENTS AND STANDARDS

- 7.1. Reagent water for non-volatile constituents must be produced by a Millipore DI system or equivalent. For volatile constituents, water must be passed through an activated carbon filter bed (Milli-Q or tap water passed through activated carbon). Reagent water must be free of the analytes of interest as demonstrated through the analysis of method blanks.
- 7.2. Hydrochloric acid, 1 N: Carefully add 83 mL concentrated reagent grade HCl to 800 mL reagent water, cool and dilute to 1 liter with reagent water. Cap and shake to mix well.

TCLP and SPLP Leaching Procedure

---

- 7.3. Sodium hydroxide, 1 N: Carefully add 40 g reagent grade NaOH pellets to 800 mL reagent water, stir until the pellets are completely dissolved, cool and dilute to 1 liter with reagent water.

CAUTION: Heat is generated during this process.

- 7.4. Acetic acid, glacial: concentrated, reagent grade liquid (HOAc).
- 7.6. pH calibration solutions: buffered to a pH of 4, 7, and 10. Commercially available.
- 7.5. TCLP Leaching Fluids

7.5.1. General Comments

7.5.1.1. The pH of both solutions listed below should be monitored daily and the pH probes are to be calibrated prior to use.

7.5.1.2. The leaching fluids MUST be prepared correctly. If the desired pH range is not achieved and maintained, the TCLP may yield erroneous results due to improper leaching. If the pH is not within the specifications, the fluid must be discarded and fresh extraction fluid prepared.

7.5.1.3. Additional volumes of extraction fluids listed below may be prepared by multiplying the amounts of acetic acid and NaOH by the number of liters of extraction fluid required.

7.5.2. TCLP Fluid #1: Carefully add 5.7 mL glacial acetic acid and 64.3 mL of 1 N NaOH to 500 mL reagent water in a 1 liter volumetric flask. Dilute to a final volume of 1 L with reagent water, cap and shake to mix well. For 8 L of fluid use 45.6 mL glacial acetic acid and 514 mL 1N NaOH, dilute to 8 L with reagent water. When correctly prepared, the pH of this solution is  $4.93 \pm 0.05$ . The density of TCLP fluid #1 is 0.997 g/mL.

7.5.3. TCLP Fluid #2: Carefully add 5.7 mL glacial acetic acid to 500 mL reagent water in a 1 liter volumetric flask. Dilute to a final volume of 1 L with reagent water, cap and shake to mix well. For 8 L of fluid use 45.6 mL glacial acetic acid, dilute to 8 L with reagent water. When correctly prepared, the pH of this solution is  $2.88 \pm 0.05$ . The density of TCLP fluid #2 is 0.997 g/mL.

- 7.6. Nitric acid, 50% solution: Slowly and carefully add 500 mL concentrated  $\text{HNO}_3$  to 500 mL reagent water. Cap and shake to mix well.



TCLP and SPLP Leaching Procedure

---

- 7.7. Sulfuric acid / nitric acid (60/40 weight percent mixture)  $\text{H}_2\text{SO}_4/\text{HNO}_3$ . Cautiously mix 60 g of concentrated sulfuric acid with 40 g of concentrated nitric acid.
- 7.8. SPLP Leaching fluids
- 7.8.1. SPLP solutions are unbuffered and exact pH may not be attained. The pH of TCLP and SPLP fluids should be checked prior to use. If not within specifications, the fluid should be discarded and fresh fluid prepared.
- 7.8.2. SPLP fluid #1: Add 60/40 weight percent mixture of sulfuric and nitric acids to reagent water until the pH is  $4.20 \pm 0.05$ . This fluid is used for soils from a site that is east of the Mississippi River and for wastes and wastewaters.
- 7.8.3. SPLP fluid #2: Add 60/40 weight percent mixture of sulfuric and nitric acids to reagent water until the pH is  $5.00 \pm 0.05$ . This fluid is used for soils from a site that is west of the Mississippi River.
- 7.8.4. SPLP fluid #3: This fluid is reagent water and is used for leaching of volatiles. Additionally, any cyanide-containing waste or soil is leached with fluid #3 because leaching of cyanide containing samples under acidic conditions may result in the formation of hydrogen cyanide gas.
- 7.9. Methanol and methylene chloride - used to aid in cleaning oil contaminated equipment.

**8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

- 8.1. Samples being analyzed for non-volatile organic compounds should be collected and stored in glass containers with Teflon lid liners. Chemical preservatives shall NOT be added UNTIL AFTER leachate generation.
- 8.2. Samples being analyzed for metals only can be collected in either glass or polyethylene containers.
- 8.3. When the waste is to be evaluated for volatile analytes, care should be taken to minimize the loss of volatiles. Samples shall be collected and stored in a manner intended to prevent the loss of volatile analytes (e.g., samples should be collected in Teflon lined septum capped vials with minimal headspace and stored at  $4 \pm 2^\circ\text{C}$ ). Samples should be opened only immediately prior to extraction.

## TCLP and SPLP Leaching Procedure

- 8.4. Samples should be refrigerated to  $4 \pm 2^{\circ}\text{C}$  unless refrigeration results in irreversible physical changes to the waste. If precipitation occurs, the entire sample (including precipitate) should be extracted.
- 8.5. The minimum TCLP sample collection size is determined by the physical state or states of the waste and the analytes of concern. The amount of waste required varies with the percent solids. The lower the percent solids, the more waste will be required for preliminary and final testing. For aqueous samples containing between 0.5 and 10% solids, several kilograms of sample are required to complete the analyses. The general minimal requirements when the samples are 100% solids include: 1 - 32 oz jar for semi-volatile organic analysis and metals, and 1 - 4 oz jar for volatile organic analysis. Low density sample materials, such as rags or vegetation, will require larger volumes of sample. For liquid samples (less than 50% solids), minimum requirements are 2 - 32 oz jars for semi-volatile organic analysis and metals, and 2 - 8 oz jars for volatile organic analysis. If volatile organic analysis is the only requested parameter, 2 separate jars are required. If matrix spike or duplicate control samples are requested, additional sample volume is required. If sufficient sample volumes were not received, analyses cannot be started and the client should be notified as soon as possible.
- 8.6. TCLP leachates should be prepared for analysis and analyzed as soon as possible following extraction. Leachates or portions of leachates for metallic analyte determinations must be acidified with nitric acid to a pH less than 2, unless precipitation occurs. If precipitation occurs upon addition of nitric acid to a small aliquot of the leachate, then the remaining portion of the leachate shall not be acidified and the leachate shall be analyzed as soon as possible. All other leachates should be stored under refrigeration ( $4 \pm 2^{\circ}\text{C}$ ) until analyzed. ZHE leachates must be stored in VOA vials filled to eliminate all headspace.
- 8.7. Samples are subject to appropriate treatment within the following time periods:

**Table 1 – Holding Times (days)**

Parameter	Collection to Start of TCLP Leach	End of TCLP Tumble to Preparation	Start of TCLP Leach or Semi-volatile Prep Extraction to Analysis	Total Elapsed Time
Volatiles:	14	N/A	14	28
Semi-volatiles:	14	7	40	61
Mercury:	28	N/A	28	56
Other Metals:	180	N/A	180	360

**NOTE:** The initial holding time is measured from date of collection to date TCLP extraction started. (This should be the TCLP extraction date in QuanTims.) Semi-volatile method prep holding time is measured from the day tumbling is complete

TCLP and SPLP Leaching Procedure

---

to the start of method extraction. Subsequent analysis holding times are measured from the date extraction (TCLP or method prep) starts. If sample holding times are exceeded, the values obtained will be considered minimal concentrations. Exceeding holding times is not acceptable in establishing that a waste does not exceed the regulatory level. Exceeding the holding time will not invalidate characterization if the waste exceeds the regulatory limit. The Total Elapsed Time is to be used as guidance. If preps are initiated at the last possible moment of a holding time, the elapsed times may be exceeded.

## 9. QUALITY CONTROL

- 9.1. Quality Control Batch (QC Batch) - QA-003 defines a QC Batch as a set of up to 20 field samples of similar matrix that behave similarly and are processed using the same procedures, reagents and standards within the same time period. The same lot of reagents must be used within a batch. A minimum of one TCLP extraction blank (Method Blank), one Laboratory Control Sample (LCS), one Matrix Spike (MS), and one Matrix Spike Duplicate (MSD) will be prepared with each TCLP leachate batch.
- 9.2. TCLP Extraction Blanks - A minimum of one blank (using the same extraction fluid as used for the samples) must be prepared and analyzed for every batch of samples extracted in a particular vessel type. The blanks are generated in the same way as the samples (i.e., blanks will be tumbled and filtered with the samples). ZHE Extraction vessels will be uniquely numbered. Consult the STL QC Program and the individual analysis SOPs for blank acceptance criteria.
- 9.3. Laboratory Control Sample (LCS) - A LCS is required with each batch of 20 or fewer samples. The LCS shall be generated after a batch of TCLP leachates have been generated (i.e., at the time of the preparative digestion or extraction) by spiking an aliquot of the appropriate extraction fluid used for that batch or reagent water. Consult the individual analysis SOPs for additional LCS guidance (i.e., spike amounts, spike levels, recovery criteria, etc.).
- 9.4. Matrix Spike (MS/MSD) - Matrix spikes are used to monitor the performance of the analytical methods on the matrix and to assess the presence of interferences. A MS/MSD pair are required with each batch of 20 or fewer samples.
  - 9.4.1. Matrix spikes are to be added after filtration of the TCLP leachate. Spikes are not to be added prior to the TCLP leaching. For metals, matrix spikes are to be added before preservation with nitric acid.
  - 9.4.2. Consult the individual analysis SOPs for additional guidance on spike compounds and levels.

TCLP and SPLP Leaching Procedure

---

## 9.5. Corrective Actions

9.5.1. Consult the STL QC Program and individual analysis SOPs for corrective action for blanks and LCS

9.5.2. Method of Standard Additions (MSA) shall be used for metals if all of the following conditions are met:

- Recovery of the analyte in matrix spike is not at least 50%,
- the concentration of the analyte does not exceed the regulatory level, and
- the concentration of the analyte measured in the sample is within 20% of the appropriate regulatory level.

If the matrix spike recovery is 5% or less due to dilution or matrix interference, contact the project manager and client for guidance. The client should also be contacted prior to initiation of any MSA steps. Refer to the individual analysis SOPs for details on how to perform MSA analysis.

**10. CALIBRATION AND STANDARDIZATION**

10.1. Refer to appropriate analysis SOPs.

**11. PROCEDURE**

## 11.1. GENERAL COMMENTS

11.1.1. One-time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented on a Nonconformance Memo kept in the project file and described in the final report. The variation must be approved by a project manager, Technical Specialist and QA Manager. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

11.1.2. All masses should be recorded to the nearest 0.1 g.

11.2. PRELIMINARY SAMPLE EVALUATIONS (Refer to Flow Chart #1, Appendix D)

## TCLP and SPLP Leaching Procedure

- 11.2.1. Determine the total volume of TCLP leachate (solid phase leachate + liquid filtrate) that needs to be generated for analysis according to the following:

**Table 2. Recommended TCLP Leachate Volume**

Analysis	TCLP Required Volume (mL)	SPLP Required Volume (mL)
Volatiles	3 x 40	3 x 40
Semi-volatiles	500	1000
Pesticides	500	1000
Herbicides	500	1000
Metals	300	300

- 11.2.1.1. For TCLP and SPLP samples used for matrix spike and matrix spike duplicate analysis, two to three times the listed volumes are required.

11.2.2. Sample Description (determine sample matrix)

- 11.2.2.1. Solid - If the waste will obviously yield no free liquid when subjected to pressure filtration, then proceed to Section 11.2.5 or 11.4 (Bottle Extraction Procedure or ZHE Procedure).

- 11.2.2.2. Liquid - If the sample is a monophasic liquid, proceed to Section 11.2.3 (Percent Solid Determination).

- 11.2.2.3. Multiphasic – The sample has discernible layers (liquid/liquid or liquid/solid). If more than one container of multi-phasic materials is received from the field, each container might show different amounts of each phase. Consult client to determine sample selection alternatives (composite all sample containers, select one, resample, etc.) if this occurs.

11.2.3. Solids Determination

- 11.2.3.1. Determine Type of Filtration Apparatus and Process

- 11.2.3.1.1. Percent Solids and ZHE Extractions - The ZHE filtration apparatus cannot accurately determine percent solids less than 5%. If an extraction is to be performed solely for volatile organic

TCLP and SPLP Leaching Procedure

---

compounds and the percent solids concentration is apparently greater than 5%, proceed to Section 11.4 (ZHE Extraction Procedure, Volatile Constituents). Otherwise, continue with Section 11.2.3.2. The aliquot of sample used here cannot be used again for the ZHE extraction.

11.2.3.1.2.If the sample is mostly a non-viscous liquid (water or non-viscous organic liquid) of low solids content (expected to be < 0.5%) , vacuum filtration should be used initially. Proceed to determination of percent dry solids (Section 11.2.3.2)

11.2.3.1.3.If the sample is viscous (sludge, oil, or is expected to have solids content > 0.5%), use pressure filtration. Proceed to determination of wet solids (Section 11.2.3.3).

11.2.3.2.Determination of percent dry solids

11.2.3.2.1.Measure and record the weight of the filter. Load the filter into the filter holder and assemble vacuum filter apparatus.

11.2.3.2.2.Homogenize the waste, then transfer 100 g subsample to a glass beaker. Record the sample weight in the percent dry solids section of the logbook.

11.2.3.2.3.Turn on vacuum source. Transfer the sample to the vacuum filtration device attempting to spread the waste sample evenly over the surface of the filter. Be sure to transfer all particulates from the beaker to the filter. Use a reagent water rinse if necessary.

11.2.3.2.4.Once all liquid has been pulled though the filter, remove the filter with the wet solids from the vacuum filtration apparatus.

11.2.3.2.5.Dry the filter and solid phase at  $100 \pm 20^{\circ} \text{C}$  for approximately 15 minutes.

11.2.3.2.6.Remove the filter from the oven and allow to cool in a desiccator.

11.2.3.2.7.Weigh and record the dry weight of filter + particulates.

11.2.3.2.8.Calculate and record the percent dry solids.

TCLP and SPLP Leaching Procedure

---

11.2.3.2.9.If the percent dry solids is  $\geq 0.5\%$ , repeat the drying step. Weigh and record the second filter + particulates dry weight. If the two weighings do not agree within 1%, perform additional drying and weighing until successive weighings agree within 1%.

11.2.3.2.10.If the dry solids result is  $\geq 0.5\%$ , proceed to Section 11.2.3.3 using a fresh wet portion of the multiphase waste.

11.2.3.2.11.If the percent solids result is less than 0.5%, discard the solid phase. No leaching will be necessary. Filter sufficient sample with either the pressure filtration system or ZHE system as described in Sections 11.3 and 11.4. The filtrate is the TCLP leachate.

#### 11.2.3.3.Determination of wet solids

11.2.3.3.1.Assemble the pressure filtration apparatus (use blunt forceps to handle the 0.6 to 0.8  $\mu\text{m}$  filter membrane).

11.2.3.3.2.Homogenize the waste, transfer a minimum of a 100 mL subsample to the glass beaker. Measure and record the gross weight (logbook column A).

11.2.3.3.3.Measure and record the tare weight of the filtrate collection bottle (logbook column D).

11.2.3.3.4.Transfer the sample to the filtration device attempting to spread the waste sample evenly over the surface of the filter. Measure and record the tare weight of the empty glass beaker and any residual sample (logbook column B).

11.2.3.3.5.Calculate and record the net weight of sample used for testing (logbook column C).

11.2.3.3.6.Slowly apply gentle pressure of 10 psi to the filtration apparatus. Allow the sample to filter until no SIGNIFICANT additional liquid has passed through the filter during a 2 minute period.

11.2.3.3.7.If necessary, repeat previous step by increasing the pressure in 10 psi increments until a maximum of 50 psi is reached. Stop the filtration when no additional filtrate is generated within a 2 minute period.

TCLP and SPLP Leaching Procedure

---

**Note:** Some samples will contain liquid material that does not filter (e.g., oil). Do not attempt to filter the sample again by exchanging filters. Viscous oils or any wastes which do not pass through the filter are classified by the method as a solid.

11.2.3.3.8. Remove the filtrate collection bottle, weigh and record the gross weight (logbook column E).

11.2.3.3.9. Calculate and record the net weight of filtrate (logbook column F). This result will be used in the percent solids calculation.

11.2.3.3.10. Pour the filtrate into an appropriately sized graduated cylinder. Measure and record the volume of the filtrate in the logbook.

11.2.3.3.11. Retain the filtrate for possible recombination with the leachate in Section 11.3.6. Retain the filter and wet solids for the leaching in Section 11.3.

11.2.3.3.12. For multiphase sample preparations, calculate the total weight of wet solids and record the result in logbook column G.

11.2.4. Particle-size Reduction for Fluid Selection

11.2.4.1. The subsample used for fluid selection must consist of particles less than approximately 1 mm in diameter (versus the less than 1 cm requirement for the material used for the actual extraction). The method requires a smaller particle size to partially compensate for the shorter duration of contact time with the leachate solution as compared to the full extraction. Inappropriate use of coarser materials could result in the selection of the wrong fluid type.

11.2.4.2. Surface area exclusion - size reduction is not required if the sample surface area is greater than or equal to 3.1 cm<sup>2</sup> per gram.

11.2.4.3. If the sample contains particles greater than approximately 1 mm in diameter, crush, cut, or grind the solids to the required size.

11.2.4.4. Consult a supervisor or manager when dealing with unusual sample matrices (e.g., wood, cloth, metal, brick).

11.2.5. Determination of Appropriate Extraction Fluid



TCLP and SPLP Leaching Procedure

---

11.2.5.1. If the solid content is greater than or equal to 0.5%, and if the sample is being analyzed for metals or nonvolatile organic compounds, the type of leaching solution must be determined.

11.2.5.2. Follow times, temperature, and particle size specified in this section as closely as possible. If reaction time between the acid solution and solid waste is too short or too long, the procedure may produce false pH readings.

11.2.5.3. For SPLP, refer to Section 7.8 for fluid selection. Record the fluid type in the logbook.

11.2.5.4. The TCLP leaching fluid for all volatiles is TCLP Fluid #1.

11.2.5.5. TCLP leach fluid determination for non-volatile analytes

11.2.5.5.1. Calibrate the pH meter with fresh buffer solution in accordance with the pH SOP.

11.2.5.5.2. Weigh out a  $5.0 \pm 0.1$  g subsample (less than 1 mm particle size) of the solid phase into a glass container, and record in the logbook.  
Note: If sample quantity is limited, consult supervisor or manager.

Note: Many multiphase samples have limited solids quantity. In these instances use a 5 g aliquot of the whole sample. Document this difference in the logbook comment section.

11.2.5.5.3. Add  $96.5 \pm 1.0$  mL of reagent water, cover with a watchglass, and stir for 5 minutes.

11.2.5.5.4. Measure and record the pre-test sample pH in the logbook.

**Note:** To avoid damaging a glass pH probe when organic liquid is present, use narrow range pH indicator paper or an ISFET pH meter.

11.2.5.5.5. If the pH is less than or equal to 5.0, use TCLP Fluid #1.

11.2.5.5.6. If the fluid pH is greater than 5.0, add 3.5 mL 1 N HCl. Slurry the sample briefly. Clip thermometer to the inside edge of one sample in each pre-test sample group to monitor the temperature.

TCLP and SPLP Leaching Procedure

---

All samples in the group must go on the hot plate at the same time in order for the temperature of the one monitored sample to represent the others. Heat to 50°C and maintain for 10 minutes.

**Note:** The heating cycle is a critical step. If the solid waste does not remain in contact with the acidic solution under specified time and temperature conditions, an erroneous pH may be measured.

11.2.5.5.7. Cool to room temperature.

11.2.5.5.8. Measure and record the pH immediately after the sample has reached room temperature.

11.2.5.5.8.1. If the pH is less than or equal to 5.0, use TCLP Fluid #1. Record the buffer in the logbook.

11.2.5.5.8.2. If the pH is greater than 5.0, use TCLP Fluid #2. Record the buffer in the logbook.

11.2.6. For samples requiring analysis for semi-volatile organics, pesticides, herbicides or metals proceed to Section 11.3.

11.2.7. For samples requiring analysis for volatile organics (ZHE), proceed to Section 11.4.

11.3. BOTTLE EXTRACTION PROCEDURE: NON-VOLATILE CONSTITUENTS: SEMI-VOLATILES, PESTICIDES, HERBICIDES, METALS (Refer to Flow Chart #2, Appendix D)

11.3.1. Evaluate the solid portion of the waste for particle size. If it contains particles greater than 1 cm in size, prepare the solid portion of the waste for leaching by crushing, cutting, or grinding such that all particles are less than 1 cm in size (i.e., capable of passing through a 9.5 mm, 0.375 inch, standard sieve). Size reduction is not required if the sample surface area is greater than or equal to 3.1 cm<sup>2</sup> per gram. If particle size reduction was required, record this in comments column in logbook.

11.3.1.1. Consult your supervisor or manager when dealing with unusual sample matrices (e.g., wood, cloth, metal, brick). Scissors or shears may be used to cut cloth, plastic or sheet metal. Saws may be used for wood or solid metal. Determination of particle size reduction tools should take into account the requested analytes (e.g. avoid chromium steel tools when TCLP metals have been requested). Bricks, rocks, or other solids amenable to grinding may be

subcontracted out for particle size reduction. (Contact PA or PM.) Note that size reduction to fine powder is not appropriate, and could invalidate results. If necessary, consult client for guidance.

11.3.2. Determine the minimum total volume of solid phase leachate that needs to be generated. Refer to Section 11.2.1.

11.3.3. Use 100 g of solid unless sample quantity is limited. If limited sample, divide the total volume of solid phase leachate required by 20 to determine the minimum mass of solid phase required for leaching. Round this mass UP to the nearest 5g. Client must be notified if less than 100 g of solid material is used.

**Note:** Solid phase material is often in limited quantity from multiphase samples. Generally all the *solid* phase material and the filter from Section 11.2.3.3.11 are transferred to the leaching bottle

11.3.4. Weigh the required mass of solid phase into an appropriate bottle (plastic for metals only, glass for all others) and **slowly** add 20 times its mass of appropriate leaching fluid (e.g., 100 g of sample would require 2000 mL of leaching fluid). Record the weight of the sample aliquoted for the extraction. Record the volume of extraction fluid added in the logbook if other than 2000 mL.

11.3.5. Ensure any effervescence has stopped before capping the bottle tightly. Secure in a rotary agitator and rotate end-over-end at 28-32 rpm for 16-20 hours. The temperature of the room should be  $23 \pm 2^{\circ}\text{C}$ . Record the rotary agitator I.D. and the date and time extraction is started and completed in the logbook.

**NOTE:** As agitation continues, pressure may build up within the bottle for some types of wastes. To relieve excessive pressure, the bottle may be removed and opened periodically in a properly vented hood to relieve any built-up pressure

11.3.6. After tumbling in the rotary agitator is completed, remove the bottle and allow the solids to settle. Record the date and time the extraction is completed in the logbook. If sample was multiphase with an initial filtrate, drop a few drops of the filtrate (with a disposable glass pipette) into the extraction bottle and observe whether the filtrate is insoluble or forms a precipitate with the leachate. If so then the filtrate is not compatible with the leachate and must be bottled and analyzed separately. The results are normally mathematically recombined (Section 12.1.2). If the filtrate is compatible with the leachate (ie completely soluble) then pour the entire filtrate into the leachate bottle, recap and mix. Proceed with the leachate filtration step in the next section.

TCLP and SPLP Leaching Procedure

---

- 11.3.7. Filter the sample using pressure filtration by filtering through a new glass fiber filter. For final filtration of the TCLP leachate, the glass fiber filter may be changed, if necessary, to facilitate filtration. Filters must be acid washed if metals are to be determined (see Section 6.3). The entire sample need not be filtered; however, sufficient volume should be generated to support the required analyses.
- 11.3.8. Measure the pH of the TCLP leachate and record in the logbook. (Use narrow range pH paper or ISFET pH meter to measure the pH of oily samples as a glass pH probe may be damaged.)
- 11.3.9. Prepare subsamples for metals for MS/MSD quality control testing using the appropriate TCLP spiking solution (do not spike for organics). Refer to the appropriate determinative SOPs for further guidance on the spike components, levels and action criteria.
- 11.3.10. Immediately preserve the leachate as follows:

Metals	pH < 2 w/ HNO <sub>3</sub> for aqueous filtrates and leachates (do not acidify oils and other non-aqueous liquids)
All others	Refrigerate to 4 ± 2 °C

**Note:** Refer to Section 8.6 if precipitation occurs upon preservation.

- 11.3.11. Label each sample with the appropriate information and submit to the appropriate analytical groups for prep and analysis. For multiphase samples requiring mathematical recombination provide copies of the TCLP preparation logbook sheets to the sample preparation and analysis groups. Most mathematically recombined samples will require data entry for the filtrate and leachate portions as well as for the mathematically recombined results. Contact the project manager to ensure the proper sample login is completed.
- 11.4. **ZHE EXTRACTION PROCEDURE: VOLATILE CONSTITUENTS** (Refer to Flow Chart #3, Appendix D)
- 11.4.1. Use the ZHE device to obtain a TCLP leachate for analysis of volatile compounds only. Leachate resulting from the use of the ZHE shall NOT be used to evaluate the mobility of non-volatile analytes (e.g., metals, pesticides, herbicides and semi-volatile organics).
- 11.4.2. Due to some shortcomings of the method, losses of volatile compounds may occur. Extra care should be observed during the ZHE procedure to ensure that such losses

TCLP and SPLP Leaching Procedure

---

are minimized. Charge the ZHE with sample only once and do not open the device until the final extract has been collected. Do not allow the waste, the initial liquid phase, or the extract to be exposed to the atmosphere any longer than necessary.

11.4.3. Install new O-rings and adjust the ZHE piston in the ZHE body to the appropriate height (slightly moisten the O-rings with leaching fluid if necessary).

11.4.4. If the preliminary evaluations indicated the need for particle size reduction, homogenize the waste, weigh out a sufficient size subsample and prepare for leaching by crushing, cutting, or grinding such that all particles are less than 1 cm in size as measured with a ruler (Do NOT sieve the sample). Size reduction is not required if the sample surface area is greater than or equal to 3.1 cm<sup>2</sup> per gram. If particle size reduction was required, record this in the comments column of the logbook.

**Note:** To minimize loss of volatiles, samples for volatiles that require particle size reduction should be kept in sample storage (at 4°C) until immediately before size reduction. Aggressive reduction which would generate heat should be avoided and exposure of the waste to the atmosphere should be avoided to the extent possible. Size reduction to a fine powder is not appropriate. Also see Section 11.3.1.

11.4.4.1. Consult your supervisor or manager when dealing with unusual sample matrices (e.g., wood, cloth, metal, brick). Scissors or shears may be used to cut cloth, plastic or sheet metal. Saws may be used for wood or solid metal. Bricks, rocks, or other solids amenable to grinding may be subcontracted out for particle size reduction (Contact PM).

11.4.5. Homogenize and transfer an appropriate size subsample of the waste into the ZHE and record the mass in the logbook.

11.4.5.1. For wastes that are solid, a 15 g sample is used.

11.4.5.2. For wastes containing < 0.5% solids, the liquid portion of the waste, after filtration, is defined as the TCLP leachate. Filter enough of the sample to support all of the volatile analyses required.

11.4.5.3. If the sample has  $\geq$  0.5% solids and has non-volatile TCLP/SPLP requested, the appropriate sample size should be estimated based on the wet solids content determined in Section 11.2.3.3. If ZHE only, use visual wet solids estimate to sample subaliquot.

TCLP and SPLP Leaching Procedure

---

**Note:** For wastes containing greater than 0.5% wet or dry solids, the "solids" value from the ZHE filtration process may be used to determine the volume of fluid to load into the ZHE. This approach is recommended since the solids value from Section 11.2.3.3 may differ from the ZHE filtration solids due to sample variability or differences in the filtration apparatus.

- 11.4.6. Carefully place the glass fiber filter between the support screens and secure to the ZHE. Tighten all the fittings.
- 11.4.7. Place the ZHE in a vertical position; open both the gas AND liquid inlet/outlet valves. Attach a gas line to the gas inlet/outlet valve.
- 11.4.8. If the waste is solid, slowly increase the pressure to a maximum of 50 psi to force out as much headspace as possible and proceed to Section 11.4.13.
- 11.4.9. If this is a multiphase sample, carefully apply gentle pressure of 10 psi (or more, if necessary) to force all headspace slowly out of the ZHE. At the FIRST appearance of liquid from the liquid inlet/outlet valve, quickly close the valve and discontinue gas pressure.
- 11.4.10. Assemble a syringe and place the plunger in all the way. Attach the pre-weighed syringe to the liquid inlet/outlet valve and open the valve. Record the tare weight of the collection syringe in column D of the logbook.
- 11.4.11. Carefully apply gas pressure of no more than 10 psi to force out the liquid phase. Allow the sample to filter until no SIGNIFICANT additional filtrate has passed in a 2 minute period.

**Note:** If the capacity of the syringe is reached, close the liquid inlet/outlet valve, discontinue gas pressure, remove the syringe, weigh, record weight in column E and filtrate volume in the logbook. Return to Section 11.4.10.

- 11.4.12. Repeat previous step increasing the pressure in 10 PSI increments until 50 psi is reached and no significant liquid has passed in a 2 minute period. Close the valve and discontinue gas pressure. Remove the collection device and record the total weight of the collection device with filtrate in column E and filtrate volume in the logbook. Transfer the filtrate to VOA vials and label appropriately. Calculate the weight of filtrate collected and record in column F in the logbook.

**Notes:** If the original waste contained less than 0.5% solids (Section 11.2.3.2), this

TCLP and SPLP Leaching Procedure

---

filtrate is defined as the TCLP leachate and you may proceed to Section 11.4.22. Otherwise, save the vials by storing at 4°C under minimal headspace conditions, for recombination as in Section 11.4.21.

The material remaining in the ZHE is defined to be the "solid" phase. Calculate the weight of the solid phase and record in column G of the logbook by subtracting the weight of the filtrate from the weight of the sample.

- 11.4.13. Determine the amount of buffer to use. Solid samples use 300 mL of leach fluid (20 X 15 g). For multiphase samples use the wet solids (column G) amount and multiply by 20. Record the leach fluid volume in column H of the logbook.

**Note:** The TCLP ZHE prep uses only TCLP fluid #1; the SPLP ZHE prep uses only SPLP fluid #3.

- 11.4.14. Load the fluid transfer reservoir with an excess of Fluid #1 and preflush the transfer line to eliminate air pockets. Be sure the required volume remains.

- 11.4.15. Attach the transfer line to the liquid inlet/outlet valve and open the valve. Carefully pump the required volume into the ZHE and close the valve. Disconnect the transfer line.

- 11.4.16. Check the ZHE to make sure all the valves are closed and manually rotate the ZHE (end-over-end) 2 or 3 times. Reposition the ZHE in the vertical position.

- 11.4.17. Pressurize the ZHE to 5-10 psi. If the ZHE appears to be leaking, follow the corrective action protocols recommended by the manufacturer and repeat the analysis.

- 11.4.18. Slowly open the liquid inlet/outlet valve to bleed out any headspace that may have been introduced during the introduction of the Fluid. Upon the first sign of liquid from the valve, close the valve.

- 11.4.19. Repressurize the ZHE to 5-10 psi and place in the rotary agitator. Rotate at 28-32 rpm for 16-20 hours. Room temperature should be  $23 \pm 2$  °C. The room temperature is recorded using a continuous temperature monitor.

- 11.4.20. Confirm that the pressure of 5-10 psi was maintained throughout the leaching. If it was NOT maintained, return to Section 11.4 and repeat the leachate with a new aliquot of sample.

TCLP and SPLP Leaching Procedure

---

11.4.21. If there is an initial liquid filtrate (Sec 11.4.12) determine if it is compatible with the leachate if the filtrate has not been previously tested (Sec. 11.3.6).

11.4.21.1. Remove the plunger from the syringe and attach the barrel to the ZHE vessel. Open the outlet valve and pressurize as necessary to transfer about 1 mL of leachate into the syringe. Close the outlet valve.

11.4.21.2. With a glass pipette transfer a few drops of initial filtrate into the open syringe barrel. Formation of separate layers or a precipitate indicates the filtrate and leachate are not compatible. Bottle the filtrate for separate preparation and analysis. The results are normally mathematically recombined.

11.4.21.3. If the filtrate is compatible gently pour the remainder of the filtrate into the syringe barrel. Install the plunger. Bleed any pressure in the ZHE piston. Open the inlet/outlet valve and depress the syringe plunger to inject the filtrate into the ZHE vessel. Do not inject the air bubble (if present) from the syringe.

11.4.21.4. Close the valve and rotate a few times to mix. Proceed with leachate filtration as described in the next section.

11.4.22. Attach an empty syringe to the outlet valve. Open the valve and pressurize the piston to expel the leachate from the ZHE vessel. Following collection, store the TCLP leachate in 2 or 3 40-mL VOA vials with minimal headspace at  $4 \pm 2$  °C and prepare for analysis as soon as possible using the appropriate organic analysis procedure (see Section 16.3).

11.4.23. If the individual phases are analyzed separately, combine the results mathematically by using the recombination calculation in Section 12.1.2. . Provide copies of the TCLP preparation logbook sheets to the sample preparation and analysis groups. Most mathematically recombined samples will require data entry for the filtrate and leachate portions as well as for the mathematically recombined results. Contact the project manager to ensure the proper sample login is completed.

## 12. DATA ANALYSIS AND CALCULATIONS

### 12.1. Calculations



TCLP and SPLP Leaching Procedure

---

## 12.1.1. Calculation of weight of extraction fluid to use:

Volume of extraction fluid = 20 X weight of wet solids to be extracted

## 12.1.2. Mathematical recombination of analytical results:

$$\text{Final Analyte Concentration} = \frac{(V_1 \times C_1) + (V_2 \times C_2)}{V_1 + V_2}$$

$V_1$  = total volume of the initial filtrate phase (L).

$C_1$  = analyte concentration in initial filtrate phase (mg/L).

$V_2$  = volume of the theoretical solid phase leachate (L).

$C_2$  = analyte concentration in solid phase leachate (mg/L).,

## 12.2. REPORTING REQUIREMENTS

## 12.2.1. Follow these reporting conventions for multi-phase samples:

12.2.1.1. If both phases have positive results, use the values from each phase to calculate the recombined result. Use the reporting limit for each phase to calculate the recombined reporting limit.

12.2.1.2. If both phases are "ND" (not detected) the recombined result is "ND," and the reporting limit is calculated from the reporting limit for each phase.

12.2.1.3. If one phase is "ND" and the other phase has a positive result, use the zero for the "ND" phase and the positive value for the other phase to calculate the combined result. This will produce a minimum known concentration. Alternatively, at client request, the maximum possible concentration can be calculated by using the reporting limit for the "ND" phase rather than zero. The combined reporting limit is based on the reporting limit for both phases

12.2.2. Units - regardless of the nature of the sample, all TCLP and SPLP results are reported in units of mg/L.

12.2.3. For limits and significant figures, consult the appropriate analytical methods (Section 16.3).

---

TCLP and SPLP Leaching Procedure

---

12.2.4. Anomalies - all anomalies observed during the leach procedure must be noted on the worksheet or an NCM form. Some examples of such anomalies are:

12.2.4.1. Sample was monolithic - particle size reduction not possible due to nature of matrix.

12.2.4.2. Insufficient sample - less than the required 100 g minimum was available.

### 12.3. REVIEW REQUIREMENTS

12.3.1. Review all applicable holding times. If a holding time was exceeded, confirm that a holding time violation was properly documented in an NCM.

12.3.2. If Total analysis results are available, those results may be compared with the TCLP analysis results according to the following:

$$Total \geq 20 \times TCLP$$

**NOTE:** Assumes the sample is 100% Solids.

12.3.3. Total constituent analysis results can be used to demonstrate the TCLP protocol is unnecessary. In performing a TCLP analysis, there is a 20:1 dilution of the original sample with the leaching solution. Thus, if the "total constituent" result is less than 20 times the TC level, it is impossible for the leachate to "fail" and the TCLP does not need to be performed. For example, the TC level for lead is 5.0 mg/L (ppm). Therefore, if a sample of lead-contaminated soil contains less than 100 ppm total lead, a TCLP test need not be run to demonstrate that lead is less than the TCLP limit.

## 13. METHOD PERFORMANCE

13.1. Refer to individual analysis SOPs.

13.2. Training Qualification:

The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

## 14. POLLUTION PREVENTION

TCLP and SPLP Leaching Procedure

---

- 14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

**15. WASTE MANAGEMENT**

- 15.1. Waste generated in this procedure must be segregated and disposed according to the facility's hazardous wastes procedures. The Environmental Health and Safety Director should be contacted if additional information is required.

**16. REFERENCES**

- 16.1. Method 1311, Toxicity Characteristic Leaching Procedure, Revision 0, July 1992, SW-846 Final Update I.
- 16.2. Method 1312, Synthetic Precipitation Leaching Procedure, Revision 0, November 1994, SW-846 Update II.
- 16.3. Related Documents
- 16.3.1. Toxicity Characteristic: Corrections to Final Rule. Method 1311, Federal Register, Vol. 55, No. 126, Friday, June 29, 1990.
- 16.3.2. Toxicity Characteristic: Final Rule. Method 1311, Federal Register, Vol. 55, No. 61, Thursday, March 29, 1990.
- 16.3.3. Technical Background Document and Response To Comments, Method 1311, Toxicity Characteristic Leaching Procedure, USEPA/OSW, April, 1989.
- 16.3.4. QA-003, Quality Control Program
- 16.3.5. CORP-IP-0003NC: Acid Digestion of Aqueous Samples by SW846 and MCAWW 200 Series Methods.
- 16.3.6. CORP-MT-0001NC: Inductively Coupled Plasma-Atomic Emission Spectroscopy, Spectrometric Method for Trace Element Analysis, Method 6010B and Method 200.7.
- 16.3.7. CORP-MT-0003NC: Graphite Furnace Atomic Absorption Spectroscopy, (Thallium only)

TCLP and SPLP Leaching Procedure

---

- 16.3.8. CORP-MT-0005NC: Preparation and Analysis of Mercury in Aqueous Samples by Cold Vapor Atomic Absorption, SW-846 7470A and MCAWW 245.1.
- 16.3.9. CORP-MS-0002NC: Determination of Volatile Organics by GC/MS based on Methods 8260B, 624, and 524.2.
- 16.3.10. CORP-MS-0001NC : GC/MS Analysis Based on Method 8270C and 625.
- 16.3.11. CORP-GC-0001NC: Gas Chromatographic Analysis Based on Methods 8000B, 8021B, 8081A, 8082, 8151A, 8310, 610, and 8141A.
- 16.3.12. CORP-OP-0001NC: Extraction and Cleanup of Organic Compounds from Waters and Soils, Based on SW846 3500 Series, 3600 Series, 8151A and 600 Series Methods.

**17. MISCELLANEOUS****17.1. Modifications/Interpretations from Reference Methods**

- 17.1.1. Section 11.2: Preliminary Evaluations. Section 7.1 of the source method 1311 states that the sample aliquot used for the preliminary evaluation "...may not actually undergo TCLP extraction." Section 7.1.5 of the source method indicates that the portion used for the preliminary evaluation may be used for either the ZHE or non-volatile extraction if the sample was 100% solid. Section 7.1.5 further indicates that if the sample was subjected to filtration (i.e., < 100% solid) that this aliquot may be used for the non-volatile extraction procedure only as long as sufficient sample is available (minimum 100 g). Samples which have been subjected to the oven drying step may not be used for TCLP extraction because solid phase degradation may result upon heating.
- 17.1.2. Sections 11.3.6 and 11.4.21: Determination of Filtrate/Extraction Fluid Compatibility. Section 7.2.13 of the source method provides no guidance as to how to make this determination. As a result, the procedure herein was developed.
- 17.1.3. Section 9.2: TCLP Extraction Blanks. Section 8.1 of the source method states that a minimum of one blank for every 20 extractions "...that have been conducted in an extraction vessel." STL has interpreted this to mean one blank per twenty samples leached per TYPE of leaching vessel (i.e., Bottle or ZHE) per leach fluid used.

17.1.4. Section 11.2.5.5.8.1: Determination of Appropriate Extraction Fluid. Method 1311 does not address the appropriate approach to take if the pH equals 5.0. This SOP requires that Fluid #1 must be used if the pH is less than or equal to 5.0.

17.1.5. Section 9.4: QA/QC - Matrix Spikes. Section 8.2 of the source method states "A matrix spike shall be performed for each waste type..." and "A minimum of one matrix spike must be analyzed for each analytical batch." Further, Section 8.2.3 of the source method also states "The purpose of the matrix spike is to monitor the performance of the analytical methods used, and to determine whether matrix interferences exist." The standard STL LQM is designed to address the performance monitoring of analytical methodology through the LCS program. A minimum of one MS and MSD will be prepared for each TCLP leachate batch. The MS/MSD results are used to determine the effect of a matrix on the precision and accuracy of the analytical process. Due to the potential variability of the matrix of each sample, the MS/MSD results have immediate bearing only on the specific sample spiked and not all samples in the batch.

17.1.6. Section 8.2.2 of the source method states that "In most cases, matrix spikes should be added at a concentration equivalent to the corresponding regulatory level." The method also states "If the analyte concentration is less than one half the regulatory level, the spike concentration may be as low as one half of the analyte concentration but may not be less than five times the method detection limit". For several analytes, spiking at the regulatory level is inappropriate to the range of analysis afforded by the determinative methods. Due to the wide range in these levels, STL spikes at the levels specified in the determinative SOPs.

## 17.2. Modifications from Previous SOP

## 17.3. Facility Specific SOPs

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none. Refer to the Appendices for any facility specific information required to support this SOP.

8771215

SOP No. CORP-IP-0004NC

Revision No. 1.1

Revision Date: 10/10/00

Page: 32 of 51

---

## APPENDIX A - TABLES

### APPENDIX A

### TABLES

8771216

SOP No. CORP-IP-0004NC

Revision No. 1.1

Revision Date: 10/10/00

Page: 33 of 51

APPENDIX A - TABLES

---

**Table 3 - Toxicity Characteristic Analytes and Regulatory Levels (Final Rule)**

## APPENDIX A - TABLES

Contaminant	mg/L
Arsenic	5.0
Barium	100.0
Benzene	0.5
Cadmium	1.0
Carbon tetrachloride	0.5
Chlordane	0.03
Chlorobenzene	100.0
Chloroform	6.0
Chromium	5.0
o-Cresols	200.0
m-Cresols	200.0
p-Cresols	200.0
Total Cresols (used if isomers not resolved)	200.0
2,4-D	10.0
1,4-Dichlorobenzene	7.5
1,2-Dichloroethane	0.5
2,4-Dinitrotoluene	0.13
1,1-Dichloroethylene	0.7
Endrin	0.02
Heptachlor (& epoxide)	0.008
Hexachlorobenzene	0.13
Hexachlorobutadiene	0.5
Hexachloroethane	3.0
Lead	5.0
Lindane	0.4
Mercury	0.2
Methoxychlor	10.0
Methyl ethyl ketone	200.0
Nitrobenzene	2.0
Pentachlorophenol	100.0
Pyridine	5.0
Selenium	1.0
Silver	5.0
Tetrachloroethylene	0.7
Toxaphene	0.5
Trichloroethylene	0.5
2,4,5-Trichlorophenol	400.0
2,4,6-Trichlorophenol	2.0
2,4,5-TP (Silvex)	1.0
Vinyl chloride	0.2



8771218

SOP No. CORP-IP-0004NC

Revision No. 1.1

Revision Date: 10/10/00

Page: 35 of 51

---

APPENDIX B - FIGURES

**APPENDIX B**

**FIGURES**

## APPENDIX B - FIGURES

Figure 1 &amp; 2 - Rotary Agitation Apparatus and Zero Headspace Extraction Vessel (ZHE)





Figure 1. Rotary Agitation Apparatus



## APPENDIX B - FIGURES

Figure 3 - US Environmental Protection Agency Memorandum #35, Page 1

UNITED STATES ENVIRONMENTAL PROTECTION AGENCY  
WASHINGTON, D.C. 20460

OFFICE OF  
SOLID WASTE AND EMERGENCY RESPONSE

MEMORANDUM # 35

DATE: June 12, 1992  
SUBJECT: Notes on RCRA Methods and QA Activities  
From: Gail Hansen, Chief *Gail Hansen*  
Methods Section (OS-331)

This memo addresses the following topics:

- o 1992 Symposium on Waste Testing and Quality Assurance
- o SW-846 Update
  - Final Rule for January 23, 1989 Proposed Rule
  - Notice, Proposed Rulemaking for the Second Update to the Third Edition
- o Chlorofluorocarbon 113 (CFC-113) Solvent Replacement Update
- o Environmental Monitoring Methods Index (EMMI)
- o Sampling Work Group Formation
- o MICE Update
- o Oily Waste Analysis
- o Electronic SW-846 Availability.

## APPENDIX B - FIGURES

Figure 3 - US Environmental Protection Agency Memorandum #35, Page 10

Oil Waste Analysis

One of the most frequently asked questions on the MICE Service concerns the application of the TCLP, Method 1311, to oily wastes. Many callers request technical guidance on the extraction of oily wastes due to the difficulty in the filtration on these types of waste. In many cases, an oily waste does not filter completely due to premature clogging of the glass fiber filter. This can result in the retention of standing liquid on the glass fiber filter. Material that do not pass through the glass fiber filter at the conclusion of the filtration step is defined by the method as the solid phase of the waste. The solid phase is then subjected to the leaching procedure of the TCLP. For oily wastes, clogging of the glass fiber filter can result in an overestimation of the amount of solid material available for leaching.

To solve this problem, the Agency recommends a conservative approach, one that probably will overestimate the amount of leaching. Rather than performing the TCLP extraction on the unfiltered portion of the oily waste, assume the waste is 100% liquid (e.g., will pass through the glass fiber filter) and perform a totals analysis on the oily waste to determine if the oil exceeds the appropriate regulatory level.

Filterable waste oil generated during the TCLP must be analyzed for a variety of organic and inorganic analytes. The OSW recognizes the difficulty in achieving acceptable performance for the analysis of waste oil using methods currently provided in SW-846. As a result, the Agency will provide several new methods for the preparation and analysis of oil samples to the Organic Methods Workgroup in July. In addition, a microwave assisted digestion procedure should improve the analysis of metals and will be proposed as part of the Second Update of the Third Edition of SW-846. Brief descriptions of these techniques are provided below, for additional information on the organic procedures contact Barry Lesnik at (202) 260-7459. For additional information on microwave digestion contact Ollie Fordham (202) 260-4778.

The use of purge-and-trap (Method 5030) for volatiles in oil generally results in severe contamination of analytical instrumentation. Traps, transfer lines and chromatography columns may become contaminated with oil. This leads to elevated baselines, hydrocarbon background in subsequent analyses, and cross-contamination. Headspace (Method 3810) is currently allowed only as a screening procedure in SW-846. The Agency is evaluating the use of headspace in conjunction with isotope dilution mass spectrometry for the quantitative analysis of volatiles in oil. Headspace reduces interference problems encountered with purge-and-trap. However, headspace quantitation can be questionable because the distribution of analytes is not

SOP No. CORP-IP-0004NC  
Revision No. 1.1  
Revision Date: 10/10/00  
Page: 39 of 51

APPENDIX C – Logbook Pages

---

8771222

## APPENDIX C

### Logbook Pages

**TCLP/SPLP ZHE Logbook**

## STL NORTH CANTON ZHE LOGBOOK

[illegible]

8771223

8771224

SOP No. CORP-IP-0004NC

Revision No. 1.1

Revision Date: 10/10/00

Page: 41 of 51

APPENDIX C – Logbook Pages

N:\QAQC\SOPS\Inorgprep\CPIP0004NC TCLP\_SPLP.doc 09/28/00





SOP No. CORP-IP-0004NC  
Revision No. 1.1  
Revision Date: 10/10/00  
Page: 43 of 51

APPENDIX C – Logbook Pages

\* Filtrate volume only needed when filtrate is incompatible with leachate

N:\QAQC\SOPS\Inorgprep\CPIP0004NC TCLP\_SPLP.doc 09/28/00

8771226

SOP No.: NC-IP-0005

STL NORTH CANTON TCLP LOGBOOK

[illegible]

SOP No. CORP-IP-0004NC  
Revision No. 1.1  
Revision Date: 10/10/00  
Page: 45 of 51

APPENDIX C – Logbook Pages

N:\QAQC\SOPS\Inorgprep\CPIP0004NC TCLP\_SPLP.doc 09/28/00

## APPENDIX C – Logbook Pages

[illegible]

\* Filtrate volume only needed when filtrate is incompatible with leachate

N:\QAQC\SOPS\Unorgprep\CPIP0004NC TCLP SPLP.doc 09/28/00

8771230

SOP No. CORP-IP-0004NC

Revision No. 1.1

Revision Date: 10/10/00

Page: 47 of 51

---


APPENDIX D - FLOW CHARTS

**APPENDIX D**

**FLOW CHARTS**


APPENDIX D - FLOW CHARTS

**Flow Chart 1. Preliminary Sample Evaluation  
 (Section 11.2)**




## APPENDIX D - FLOW CHARTS

**Flow Chart 2. Bottle Extraction, Non-Volatile Constituents  
(Section 11.3)**




## APPENDIX D - FLOW CHARTS

Flow Chart 3. ZHE Extraction, Volatile Constituents  
(Section 11.4)



## APPENDIX D - FLOW CHARTS

**Flow Chart 3. ZHE Extraction  
(Continued)**



Controlled Copy  
Copy No. \_\_\_\_\_

SOP No. NC-WC-0010  
Revision No. 5  
Revision Date: 02/05/03  
Page 1 of 11

Implementation Date \_\_\_\_\_

## STL NORTH CANTON STANDARD OPERATING PROCEDURE

### TITLE: PH ELECTROMETRIC METHOD

(SUPERSEDES: REVISION 4.1, DATED 11/28/00)

Prepared by: *William P. Hill - L. Hill* 2/7/03  
Date

Reviewed by: *Angela Serra* 2/7/03  
Technology Specialist Date

Approved by: *Beth Hamlet* 2/10/03  
Quality Assurance Manager Date

Approved by: *Deborah Budd* 2/7/03  
Environmental Health and Safety Coordinator Date

Approved by: *Paul Hill* 2/8/03  
Laboratory Director Date

#### Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions. THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

**TABLE OF CONTENTS**

<b>1. Scope and Application.....</b>	<b>3</b>
<b>2. Summary of Method.....</b>	<b>3</b>
<b>3. Definitions .....</b>	<b>3</b>
<b>4. Interferences.....</b>	<b>3</b>
<b>5. SAFETY.....</b>	<b>3</b>
<b>6. Equipment and Supplies .....</b>	<b>4</b>
<b>7. Reagents and Standards .....</b>	<b>4</b>
<b>8. Sample Collection, PRESERVATION and Storage .....</b>	<b>5</b>
<b>9. Quality Control.....</b>	<b>5</b>
<b>10. Calibration and Standardization.....</b>	<b>6</b>
<b>11. Procedure .....</b>	<b>7</b>
<b>12. Data Analysis and Calculations .....</b>	<b>9</b>
<b>13. Method Performance.....</b>	<b>9</b>
<b>14. Pollution Prevention .....</b>	<b>10</b>
<b>15. Waste Management.....</b>	<b>10</b>
<b>16. References .....</b>	<b>10</b>
<b>17. Miscellaneous (Tables, Appendices, Etc...) .....</b>	<b>11</b>

**PH ELECTROMETRIC METHOD**SOP No. NC-WC-0010Revision No. 5Revision Date: 02/05/03Page 3 of 11

---

**1. SCOPE AND APPLICATION**

- 1.1. This method is applicable to the determination of pH in waters, wastewaters, and solids. It is based on SW846 Methods 9040B and 9045C and EPA Method 150.1. The approximate working range is 1 - 14 pH units. Samples with a pH of < 1 are reported as < 1.
- 1.2. The associated method codes are PU (9040B), OZ (9045C), and AJ (150.1). The preparation codes are 88 and 1C.
- 1.3. This document accurately reflects current laboratory standard operating procedures (SOP) as of the date above. All facility SOPs are maintained and updated as necessary.

**2. SUMMARY OF METHOD**

- 2.1. The pH is determined electrometrically by using an electrode. The pH meter is calibrated with a series of known pH buffers.

**3. DEFINITIONS**

- 3.1. Refer to the glossary in the Laboratory Quality Manual (LQM), latest version.

**4. INTERFERENCES**

- 4.1. Samples with very low or very high pH may give incorrect readings on the meter. For samples with a true pH of >10, the measured pH may be incorrectly low. This error can be minimized by using a low-sodium-error electrode.

**5. SAFETY**

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL North Canton associates.
- 5.2. Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined. Additional health and safety information can be obtained from the Material Safety Data Sheets (MSDS) maintained in the laboratory.

## PH ELECTROMETRIC METHOD

SOP No. NC-WC-0010Revision No. 5Revision Date: 02/05/03Page 4 of 11

---

- 5.4. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples must be opened, transferred and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
- 5.5. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.6. It is recommended that neat standards be purchased only as a last resort. The preparation of standards from neat materials and reagents {as well as glassware cleaning procedures that involved solvents such as methylene chloride} should be conducted in a fume hood with the sash closed as far as the operations will permit.
- 5.7. Standards in solution may be diluted in the open laboratory when syringes and the like are utilized.
- 5.8. All work must be stopped in the event of a known or potential compromise to the health and safety of an STL North Canton associate. The situation must be reported **immediately** to a laboratory supervisor.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. pH meter with electrode(s) and temperature compensation
- 6.2. Beakers: various
- 6.3. Top loading balance: Capable of accurately weighing  $\pm 0.01$  g
- 6.4. Stir plate and stir bars
- 6.5. Shaker or mechanical tumbler
- 6.6. Autotitrator
- 6.7. Centrifuges tubes

**7. REAGENTS AND STANDARDS**

- 7.1. Standards

**PH ELECTROMETRIC METHOD**SOP No. NC-WC-0010Revision No. 5Revision Date: 02/05/03Page 5 of 11

---

**7.1.1. Target Calibration Standards**

7.1.1.1. pH 4, 7, and pH 10 buffers, purchased

7.1.1.2. Fresh buffers are poured and used each working day.

**8. SAMPLE COLLECTION, PRESERVATION AND STORAGE**

8.1. Samples are not chemically preserved.

8.2. Samples are stored in plastic or glass containers at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$ .

8.3. Samples should be analyzed as soon as possible after sampling, but not to exceed twenty-four hours.

**9. QUALITY CONTROL****9.1. Batch Definition**

9.1.1. A batch is a group of no greater than 20 samples excluding QC samples (LCS and Sample Duplicate) which are processed similarly, with respect to the procedure. All sample setups must be initiated within a 24 hour period from the initial preparation or extraction and without interruption of the process. All samples within the batch must be treated with the same lots of reagents and the same processes.

**9.2. Sample Duplicate**

9.2.1. A sample duplicate (DU) is a second aliquot of an environmental sample, taken from the same sample container when possible, that is processed with the first aliquot of that sample. That is, sample duplicates are processed as independent samples within the same QC batch. The sample and DU results are compared to determine the effect of the sample matrix on the precision of the analytical process. As with the MS/MSD results, the sample/DU precision results are not necessarily representative of the precision for other samples in the batch.

9.2.2. Sample duplicates are performed at a frequency of 10% and must meet laboratory-specific limits for precision.

**PH ELECTROMETRIC METHOD**SOP No. NC-WC-0010Revision No. 5Revision Date: 02/05/03Page 6 of 11

---

**9.3. Laboratory Control Sample (LCS)**

9.3.1. One aqueous LCS must be processed with each analytical batch. The LCS must be carried through the entire analytical procedure. The LCS is used to monitor the accuracy of the analytical process. On-going monitoring of the LCS results provides evidence that the laboratory is performing the method within acceptable accuracy and precision guidelines.

9.3.2. A commercially available (Environmental Resource Associates or equivalent) control standard will be analyzed. Recovery must be within +/- 2% of true value.

**9.3.3. Corrective Action for LCS**

9.3.3.1. If the pH is outside the established control limits the system is out of control and corrective action must occur.

9.3.3.2. Corrective action consists of identification and correction of the cause for the out of control situation and reanalysis of all effected samples.

**10. CALIBRATION AND STANDARDIZATION****10.1. Initial Calibration**

10.1.1. Refer to the manufacturer's manual for instrumental calibration.

10.1.2. The following procedure is applicable for use with the Orion 250 pH meter.

10.1.2.1. Rinse the electrodes with reagent water and place in the pH 4.0 buffer. Press "Cal". Allow the value to stabilize and then, using the arrow keys, adjust the value up or down to read 4.00. Press Enter.

10.1.2.2. Rinse the electrodes and place in the pH 7.0 buffer. Allow the value to stabilize and then, using the arrow keys, adjust the value up or down to read 7.00. Press Enter.

10.1.2.3. Calibration Check: Rinse the electrodes and place in the pH 10.0 buffer. Allow value to stabilize. The pH should be between 9.95 and 10.05 or recalibration is necessary.

**NOTE:** When analyzing drinking water samples, calibrate as described in section 10.1.2, using the pH 7.0 and pH 10.0 buffers for calibration and the pH 4.0 buffer for the calibration check.

10.1.3 The pH meter should be calibrated daily. The calibration is recorded on the analytical logsheet.

10.1.4 If the pH meter has been turned off, it must be calibrated prior to use.

10.2. Continuing Calibration

10.2.1. A pH 7 buffer is analyzed before analysis and every ten samples to ensure the calibration remains linear.

10.2.2. The pH meter must be recalibrated if the buffer deviates by more than  $\pm 2\%$ . If this range is exceeded, reanalyze all samples analyzed since the last pH buffer that met criteria.

11. PROCEDURE

11.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and is approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

11.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

11.3. Sample Preparation

11.3.1. Waters

11.3.1.1. No preparation necessary for waters and wastewaters.

11.3.2. Solids and Soils

11.3.2.1. Place 20 g ( $\pm 0.5$  g) of sample in a beaker or other suitable container.



11.3.2.2. Add 20 mL of reagent water and mix for five minutes.

11.3.2.3. Allow sample to stand for one hour to allow the solids to settle out.

#### 11.4. Sample Analysis

##### 11.4.1. Manual Procedure

###### 11.4.1.1. Waters

11.4.1.1.1. Place the sample in a clean beaker using a sufficient volume to cover the sensing elements of the electrode(s). Allow the pH to stabilize (swirling or stirring may quicken stabilization). Record the pH on the analytical logsheet. Remove the electrodes from the sample. Rinse and gently dab off the electrodes between each measurement. Store the electrodes in pH 7 buffer when not in use.

###### 11.4.1.2. Solids

11.4.1.2.1. Immerse the pH electrodes in the supernatant layer of the sample - be careful not to stir up solids. Allow pH to stabilize and record it on the analytical logsheet. Remove and rinse the electrodes between each measurement. Store electrodes in the pH 7.0 buffer.

**NOTE:** If the sample contains oil or other substances that will coat or damage the electrodes, the pH should be analyzed following SOP# NC-WC-0009, pH - Paper Method.

11.4.2. Automated Procedure

11.4.2.1. Load the appropriate schedule on the autotitrator, starting with the pH calibration.

11.4.2.2. Pour a homogenized sample into the centrifuge tubes and place the tubes in the appropriate position on the autosampler. Remember to include a pH 7 buffer check after every ten positions.

11.4.2.3. Start the autotitrator.

11.5. Analytical Documentation

11.5.1. Record all analytical information in the analytical logbook/logsheet, including the analytical data from standards and any corrective actions or modifications to the method.

11.5.2. All standards are logged into a department standard logbook. All standards are assigned an unique number for identification. Logbooks are reviewed by the supervisor or designee.

11.5.3. Documentation such as all associated instrument printouts (final runs, screens, reruns, QC samples, etc.) and daily calibration data corresponding to all final runs is available for each data file.

11.5.4. Sample results and associated QC are entered into the LIMs after final technical review.

**12. DATA ANALYSIS AND CALCULATIONS**

12.1. Not Applicable

**13. METHOD PERFORMANCE**

13.1. Each laboratory must have initial demonstration of performance data on file and corresponding method detection limit files.

13.2. Training Qualifications:

**PH ELECTROMETRIC METHOD**SOP No. NC-WC-0010Revision No. 5Revision Date: 02/05/03Page 10 of 11

---

13.2.1. The group/team leader has the responsibility to ensure that this procedure is performed by an associate who has been properly trained in its use and has the required experience.

13.2.2. Method validation information (where applicable) in the form of laboratory demonstrations of capabilities is maintained for this method in the laboratory QA files.

**14. POLLUTION PREVENTION**

14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

**15. WASTE MANAGEMENT**

15.1. Solvent waste must be disposed of in clearly labeled waste cans.

15.2. Acid waste must be collected in clearly labeled acid waste containers.

15.3. Solid materials (gloves, soiled paper products, etc.) are placed in the solid debris container. Do not put liquids in the solid waste container.

15.4. Refer to the Laboratory Sample and Waste Disposal plan.

15.5. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL North Canton. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by an annual refresher training.

**16. REFERENCES****16.1. References**

16.1.1. SW846, Test Methods for Evaluating Solid Waste, 3rd Edition, pH Electrometric Measurement, Method 9040B

16.1.2. EPA 600, Methods for Chemical Analysis of Water and Wastes, pH (Electrometric), Method 150.1

16.1.3. SW846, Test Methods for Evaluating Solid Waste, 3rd Edition, Soil pH, method 9045C.

**PH ELECTROMETRIC METHOD**SOP No. NC-WC-0010Revision No. 5Revision Date: 02/05/03Page 11 of 11

---

16.1.4. STL North Canton Laboratory Quality Manual (LQM), current version

16.1.5. Corporate Quality Management Plan (QMP), current version.

16.2. Associated SOPs and Policies, latest version

16.2.1. QA Policy, QA-003

16.2.2. Glassware Washing, NC-QA-0014

16.2.3. Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018

16.2.4. Method Detection Limits and Instrument Detection Limits, NC-QA-0021

16.2.5. Navy/Army SOP, NC-QA-0016

**17. MISCELLANEOUS (TABLES, APPENDICES, ETC...)**

17.1. Reporting limits

17.1.1. A minimum reporting limit is not listed in LIMS. Units are reported as No Units.

8771246

Controlled Copy

Copy No: \_\_\_\_\_

Implementation Date: 06/05/01

SOP No. CORP-MS-0002NC

Revision No. 2.3

Revision Date: 05/23/01

Page 1 of 67

### STL STANDARD OPERATING PROCEDURE

#### TITLE: DETERMINATION OF VOLATILE ORGANICS BY GC/MS BASED ON METHOD 8260B, 8260A, AND 624

(SUPERSEDES: REVISION 2.2, DATED 11/28/00)

Prepared by: Thomas E. Stullen 5/25/01  
Date

Approved by: Ramona Evans 5/25/01  
Technology Specialist Date

Approved by: Paul A. Cash 5/25/01  
Quality Assurance Manager Date

Approved by: [Signature] 5-29-01  
Environmental Health and Safety Coordinator Date

Approved by: Christopher O'Connell 5-25-01  
Lab Manager Date

Approved by: [Signature] 6/4/01  
Corporate Technology and/or Corporate Quality Assurance Date

#### Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

## TABLE OF CONTENTS

1. SCOPE AND APPLICATION .....	4
2. SUMMARY OF METHOD .....	4
3. DEFINITIONS .....	5
4. INTERFERENCES .....	7
5. SAFETY.....	7
6. EQUIPMENT AND SUPPLIES .....	9
7. REAGENTS AND STANDARDS.....	11
8. SAMPLE COLLECTION, PRESERVATION AND STORAGE .....	12
9. QUALITY CONTROL .....	19
10. CALIBRATION AND STANDARDIZATION .....	22
11. PROCEDURE.....	27
12. DATA ANALYSIS AND CALCULATIONS .....	32
13. METHOD PERFORMANCE .....	39
14. POLLUTION PREVENTION .....	40
15. WASTE MANAGEMENT.....	40
16. REFERENCES .....	40
17. MISCELLANEOUS .....	40
APPENDIX A.....	
ANALYSIS OF VOLATILE ORGANICS BY METHOD 624 .....	55

**LIST OF TABLES**

<b>Table 1</b>	<b>STL Primary Standard and Reporting Limits</b>
<b>Table 2</b>	<b>STL Primary Standard Calibration Levels</b>
<b>Table 2A</b>	<b>STL Primary Standard Calibration Levels- Low Level</b>
<b>Table 3</b>	<b>STL Appendix IX Standard and Reporting Limits</b>
<b>Table 4</b>	<b>STL Appendix IX Standard Calibration Levels</b>
<b>Table 5</b>	<b>Reportable Analytes for STL Standard Tests, Primary Standard</b>
<b>Table 6</b>	<b>Reportable Analytes for STL Standard Tests, Appendix IX Standard</b>
<b>Table 7</b>	<b>Internal Standards</b>
<b>Table 8</b>	<b>Surrogate Standards</b>
<b>Table 9</b>	<b>Matrix Spike / LCS Control Standards</b>
<b>Table 10</b>	<b>BFB Tune Criteria</b>
<b>Table 11</b>	<b>SPCC Compounds and Minimum Response Factors</b>
<b>Table 12</b>	<b>CCC Compounds</b>
<b>Table 13</b>	<b>Characteristic Ions</b>
<b>Table A-1</b>	<b>Method 624 Analytes and Reporting Limits</b>
<b>Table A-2</b>	<b>Method 624 QC Acceptance Criteria</b>

## 1.0 SCOPE AND APPLICATION

- 1.1. This method is applicable to the determination of Volatile Organic Compounds in waters, wastewater, soils, sludges and other solid matrices. Standard analytes are listed in Tables 5 and 6.
- 1.2. This SOP is applicable to method 8260B. It may also be used for analysis following method 8260A. Appendix A presents modifications to the procedures in the main SOP that are necessary for analysis of wastewater by method 624. The associated LIMS method codes are QK (8260B), DN (624), and MZ (8260A). Ohio VAP projects are distinguished by Program Code 2J. The following Prep Codes are used: 15 (5 mL purge), 25 (25 mL purge), 4B (Methanol preservation, EnCore™), 4D (Sodium Bisulfate preservation, EnCore™), 4P (Frozen, EnCore™), and 73 (5030A Methanol Prep).
- 1.3. This method can be used to quantify most volatile organic compounds that have boiling points below 200°C and are insoluble or slightly soluble in water. Volatile water soluble compounds can be included in this analytical technique; however, for more soluble compounds, quantitation limits are approximately ten times higher because of poor purging efficiency.
- 1.4. The method is based upon a purge and trap, gas chromatograph/mass spectrometric (GC/MS) procedure. The approximate working range is 5 to 200 µg/L for 5 mL waters, 1 to 40 µg/L for 25 mL purge waters, 5 to 200 µg/kg for low-level soils, and 250 to 25,000 µg/kg for medium-level soils. Reporting limits are listed in Tables 1 and 3.
- 1.5. Method performance is monitored through the use of surrogate compounds, matrix spike/matrix spike duplicates, and laboratory control spike samples.

## 2. SUMMARY OF METHOD

- 2.1. Volatile compounds are introduced into the gas chromatograph by the purge and trap method. The components are separated via the chromatograph and detected using a mass spectrometer, which is used to provide both qualitative and quantitative information.
- 2.2. Aqueous samples are purged directly. Generally, soils are preserved by extracting the volatile analytes into methanol. If especially low detection limits are required, soil samples may be preserved with sodium bisulfate and purged directly.
- 2.3. In the purge and trap process, an inert gas is bubbled through the solution at ambient temperature or at 40°C (40°C required for low level soils) and the volatile components are efficiently transferred from the aqueous phase to the vapor phase. The vapor is swept through a sorbant



column where the volatile components are trapped. After purging is completed, the sorbant column (trap) is heated and backflushed with inert gas to desorb the components onto a gas chromatographic column. The gas chromatographic column is then heated to elute the components which are detected with a mass spectrometer.

- 2.4. Qualitative identifications are confirmed by analyzing standards under the same conditions used for samples and comparing the resultant mass spectra and GC retention times. Each identified component is quantified by relating the MS response for an appropriate selected ion produced by that compound to the MS response for another ion produced by an internal standard.

### 3. DEFINITIONS

#### 3.1. Batch

The batch is a set of up to 20 samples of the same matrix processed using the same procedures and reagents within the same time period. Using this method, each BFB analysis will normally start a new batch. Batches for medium level soils are defined at the sample preparation stage and may be analyzed on multiple instruments over multiple days, although reasonable effort should be made to keep the samples together.

- 3.1.1. The Quality Control batch must contain a matrix spike/spike duplicate (MS/MSD), a Laboratory Control Sample (LCS), and a method blank. In some cases, at client request, the MS/MSD may be replaced with a matrix spike and sample duplicate. Refer to the STL QC Program document (QA-003) for further details of the batch definition.

3.2. Method Blank

- 3.2.1. A method blank consisting of all reagents added to the samples must be analyzed with each batch of samples. The method blank is used to identify any background interference or contamination of the analytical system which may lead to the reporting of elevated concentration levels or false positive data.

3.3. Laboratory Control Sample (LCS)

- 3.3.1. Laboratory Control Samples are well characterized, laboratory generated samples used to monitor the laboratory's day-to-day performance of routine analytical methods. The LCS, spiked with a group of target compounds representative of the method analytes, is used to monitor the accuracy of the analytical process, independent of matrix effects. Ongoing monitoring of the LCS results provides evidence that the laboratory is performing the method within accepted QC guidelines for accuracy and precision.

3.4. Surrogates

- 3.4.1. Surrogates are organic compounds which are similar to the target analyte(s) in chemical composition and behavior in the analytical process, but which are not normally found in environmental samples. Each sample, blank, LCS, and MS/MSD is spiked with surrogate standards. Surrogate spike recoveries must be evaluated by determining whether the concentration (measured as percent recovery) falls within the required recovery limits.

3.5. Matrix Spike/Matrix Spike Duplicate (MS/MSD)

- 3.5.1. A matrix spike is an environmental sample to which known concentrations of target analytes have been added. A matrix spike duplicate is a second aliquot of the same sample which is prepared and analyzed along with the sample and matrix spike. Matrix spikes and duplicates are used to evaluate accuracy and precision in the actual sample matrix.

3.6. Calibration Check Compound (CCC)

- 3.6.1. CCCs are a representative group of compounds which are used to evaluate initial calibrations and continuing calibrations. Relative percent difference for the initial calibration and % drift for the continuing calibration response factors are calculated and compared to the specified method criteria.

3.7. System Performance Check Compounds (SPCC)

SPCCs are compounds which are sensitive to system performance problems and are used to evaluate system performance and sensitivity. A response factor from the continuing calibration is calculated for the SPCC compounds and compared to the specified method criteria.

#### 4. INTERFERENCES

- 4.1. Method interferences may be caused by contaminants in solvents, reagents, glassware, and other processing apparatus that lead to discrete artifacts. All of these materials must be routinely demonstrated to be free from interferences under conditions of the analysis by running laboratory method blanks as described in the Quality Control section. The use of ultra high purity gases, pre-purged purified reagent water, and approved lots of purge and trap grade methanol will greatly reduce introduction of contaminants. In extreme cases the purging vessels may be pre-purged to isolate the instrument from laboratory air contaminated by solvents used in other parts of the laboratory.
- 4.2. Samples can be contaminated by diffusion of volatile organics (particularly methylene chloride and fluorocarbons) into the sample through the septum seal during shipment and storage. A field blank prepared from reagent water and carried through the sampling and handling protocol can serve as a check on such contamination.
- 4.3. Matrix interferences may be caused by non-target contaminants that are coextracted from the sample. The extent of matrix interferences will vary considerably from source to source depending upon the nature and diversity of the site being sampled.
- 4.4. Cross-contamination can occur whenever high-level and low-level samples are analyzed sequentially or in the same purge position on an autosampler. Whenever an unusually concentrated sample is analyzed, it should be followed by one or more blanks to check for cross-contamination. The purge and trap system may require extensive bake-out and cleaning after a high-level sample.
- 4.5. Some samples may foam when purged due to surfactants present in the sample. When this kind of sample is encountered an antifoaming agent (e.g., J.T. Baker's Antifoam B silicone emulsion) can be used. A blank spiked with this agent must be analyzed with the sample because of the non-target interferences associated with the agent.

#### 5. SAFETY

- 5.1. Procedures shall be carried out in a manner that protects the health and safety of all STL associates.

- 
- 5.2. The Chemical Hygiene Plan (CHP) gives details about the specific health and safety practices which are to be followed in the laboratory area. Personnel must receive training in the CHP, including the written Hazard Communication plan, prior to working in the laboratory. Consult the CHP, the STL Health and Safety Policies and Procedures Manual, and available Material Safety Data Sheets (MSDS) prior to using the chemicals in the method.
  - 5.3. Consult the STL Health and Safety Policies and Procedures Manual for information on Personal Protective Equipment. Eye protection that protects against splash and a laboratory coat must be worn in the lab. Appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have been contaminated will be removed and discarded; other gloves will be cleaned immediately. Disposable gloves shall not be reused.
  - 5.4. The health and safety hazards of many of the chemicals used in this procedure have not been fully defined, therefore each chemical compound should be treated as a potential health hazard. Additional health and safety information can be obtained from the MSDS files maintained in the laboratory. The following specific hazards are known:
    - 5.4.1. Chemicals that have been classified as carcinogens, or potential carcinogens, under OSHA include: Acrylonitrile, benzene, carbon tetrachloride, chloroform, 1,2-dibromo-3-chloropropane, 1,4-dichlorobenzene, and vinyl chloride.
    - 5.4.2. Chemicals known to be flammable are: **Methanol**.
  - 5.5. Exposure to chemicals must be maintained **as low as reasonably achievable**, therefore, unless they are known to be non-hazardous, all samples should be opened, transferred, and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers will be kept closed unless transfers are being made.
  - 5.6. The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operations will permit.
  - 5.7. All work must be stopped in the event of a known or potential compromise to the health and safety of a STL associate. The situation must be reported **immediately** to a laboratory supervisor.
  - 5.8. Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices outlined in the STL Health and Safety Manual. These employees must have training on the hazardous waste disposal practices initially upon assignment of these tasks, followed by an annual refresher training.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Microsyringes: 10  $\mu$ L and larger, 0.006 inch ID needle.
- 6.2. Syringe: 5 or 25 mL glass with luerlok tip, if applicable to the purging device.
- 6.3. Balance: Analytical, capable of accurately weighing 0.0001 g, and a top-loading balance capable of weighing 0.1 g
- 6.4. Glassware:
  - 6.4.1. Vials: 20 mL with screw caps and Teflon liners.
  - 6.4.2. Volumetric flasks: 10 mL and 100 mL, class A with ground-glass stoppers.
- 6.5. Spatula: Stainless steel.
- 6.6. Disposable pipets: Pasteur.
- 6.7. pH paper: Wide range.
- 6.8. Gases:
  - 6.8.1. Helium: Ultra high purity, gr. 5, 99.999%.
  - 6.8.2. Nitrogen: Ultra high purity, from cylinders or gas generators, may be used as an alternative to helium for purge gas.
  - 6.8.3. Compressed air: Used for instrument pneumatics.
  - 6.8.4. Liquid nitrogen: Used for cryogenic cooling if necessary.
- 6.9. Purge and Trap Device: The purge and trap device consists of the sample purger, the trap, and the desorber.
  - 6.9.1. Sample Purger: The recommended purging chamber is designed to accept 5 mL samples with a water column at least 3 cm deep. The purge gas must pass through the water column as finely divided bubbles, each with a diameter of less than 3 mm at the origin. The purge gas must be introduced no more than 5 mm from the base of the water column. Alternative sample purge devices may be used provided equivalent performance is demonstrated. Low level soils are purged directly from a VOA vial.

- 
- 6.9.2. Trap: A variety of traps may be used, depending on the target analytes required. For most purposes the Vocab 3000 trap is suitable. Other traps, such as Vocab 4000, or Tenax / Silica gel / Charcoal may be used if the Quality Control criteria are met.
- 6.9.3. Desorber: The desorber should be capable of rapidly heating the trap to 180°C. Many such devices are commercially available.
- 6.9.4. Sample Heater: A heater capable of maintaining the purge device at 40°C is necessary for low level soil analysis.
- 6.10. Gas Chromatograph/Mass Spectrometer System:
- 6.10.1. Gas Chromatograph: The gas chromatograph (GC) system must be capable of temperature programming.
- 6.10.2. Gas Chromatographic Columns: Capillary columns are used. Some typical columns are listed below:
- 6.10.2.1. Column 1: 105m x 0.53 ID Rtx-624 with 3 µm film thickness.
- 6.10.2.2. Column 2: 75 m x 0.53 ID DB-624 widebore with 3 µm film thickness.
- 6.10.2.3. Mass Spectrometer: The mass spectrometer must be capable of scanning 35-300 AMU every two seconds or less, using 70 volts electron energy in the electron impact mode and capable of producing a mass spectrum that meets the required criteria when 50 ng of 4-Bromofluorobenzene (BFB) are injected onto the gas chromatograph column inlet.
- 6.10.3. GC/MS interface: In general glass jet separators are used but any interface (including direct introduction to the mass spectrometer) that achieves all acceptance criteria may be used.
- 6.10.4. Data System: A computer system that allows the continuous acquisition and storage on machine readable media of all mass spectra obtained throughout the duration of the chromatographic program. The computer must have software that allows searching any GC/MS data file for ions of a specified mass and plotting such ion abundances versus time or scan number. This type of plot is defined as an Extracted Ion Current Profile (EICP). Software must also be available that allows integrating the abundances in any EICP between the specified time or scan-number limits. Also, for the non-target compounds, software must be available that allows for the comparison of sample spectra against reference library spectra. The most recent release of the NIST/EPA mass spectral library should be used as

the reference library. The computer system must also be capable of backing up data for long-term off-line storage.

6.10.5. Cryogenic Cooling: Some columns require the use of liquid nitrogen to achieve the subambient temperature required for the proper separation of the gases.

## 7. REAGENTS AND STANDARDS

### 7.1. Reagents

7.1.1. Methanol: Purge and Trap Grade, High Purity

7.1.2. Reagent Water: High purity water that meets the requirements for a method blank when analyzed. (See section 9.4) Reagent water may be purchased as commercial distilled water and prepared by purging with an inert gas overnight. Other methods of preparing reagent water are acceptable.

### 7.2. Standards

#### 7.2.1. Calibration Standard

7.2.1.1. Stock Solutions: Stock solutions may be purchased as certified solutions from commercial sources or prepared from pure standard materials as appropriate. These standards are prepared in methanol and stored in Teflon-sealed screw-cap bottles with minimal headspace at -10° to -20°C.

7.2.1.2. Working standards: A working solution containing the compounds of interest prepared from the stock solution(s) in methanol. These standards are stored in the freezer or as recommended by the manufacturer. Working standards are monitored by comparison to the initial calibration curve. If any of the calibration check compounds drift in response from the initial calibration by more than 20% then corrective action is necessary. This may include steps such as instrument maintenance, preparing a new calibration verification standard or tuning the instrument. If the corrective actions do not correct the problem then a new initial calibration must be performed.

7.2.1.3. Aqueous Calibration Standards are prepared in reagent water using the secondary dilution standards. These aqueous standards must be prepared daily.

7.2.1.4. If stock or secondary dilution standards are purchased in sealed ampoules they may be used up to the manufacturers expiration date.

- 
- 7.2.2. Internal Standards: Internal standards are added to all samples, standards, and blank analyses. Refer to Table 7 for internal standard components.
  - 7.2.3. Surrogate Standards: Refer to Table 8 for surrogate standard components and spiking levels.
  - 7.2.4. Laboratory Control Sample Spiking Solutions: Refer to Table 9 for LCS components and spiking levels.
  - 7.2.5. Matrix Spiking Solutions: The matrix spike contains the same components as the LCS. Refer to Table 9.
  - 7.2.6. Tuning Standard: A standard is made up that will deliver 50 ng on column upon injection. A recommended concentration of 25 ng/ $\mu$ L of 4-Bromofluorobenzene in methanol is prepared as described in Sections 7.2.1.1 and 7.2.1.2.

## 8. SAMPLE COLLECTION, PRESERVATION AND STORAGE

- 8.1. Holding times for all volatile analysis are 14 days from sample collection.
- 8.2. Water samples are normally preserved at pH  $\leq$  2 with 1:1 hydrochloric acid. If residual chlorine is present, 2 drops of 10% sodium thiosulfate are added.
- 8.3. Solid samples are field preserved with sodium bisulfate solution for low level analysis, or with methanol for medium level analysis. Soil samples can also be taken using the EnCore™ sampler and preserved in the lab within 48 hours of sampling. At specific client request, unpreserved soil samples may be accepted.
- 8.4. There are several methods of sampling soil. The recommended method, which provides the minimum of field difficulties, is to take an EnCore™ sample. (The 5 g or 25 g sampler can be used, depending on client preference). Following shipment back to the lab the soil is preserved in methanol. This is the medium level procedure. If very low detection limits are needed (< 50  $\mu$ g/kg for most analytes) then it will be necessary to use two additional 5 g EnCore™ samplers or to use field preservation.
- 8.5. Sample collection for medium level analysis using EnCore™ samplers.
  - 8.5.1. Ship one 5 g (or 25 g) EnCore™ sampler per field sample position.
  - 8.5.2. An additional bottle must be shipped for percent moisture determination.



- 8.5.3. When the samples are returned to the lab, extrude the (nominal) 5g (or 25 g) sample into a tared VOA vial containing 5 mL methanol (25 mL methanol for the 25 g sampler). Obtain the weight of the soil added to the vial and note on the label.
- 8.5.4. Add the correct amount of surrogate spiking mixture. (Add 25  $\mu$ L of 2500  $\mu$ g/mL solution for a nominal 25 g sample, 5  $\mu$ L for a nominal 5 g sample.) Refer to Section 17.8 for Michigan project criteria.
- 8.5.5. Add the correct amount of matrix spiking solution to the matrix spike and matrix spike duplicate samples. (Add 500  $\mu$ L of 50  $\mu$ g/mL solution for a nominal 25 g sample, 100  $\mu$ L for a nominal 5 g sample.) Reduce the volume of methanol added to ensure the final volume is 25 mL for nominal 25 g sample or 5 mL methanol for a nominal 5 g sample. Refer to Section 17.8 for Michigan project criteria.
- 8.5.6. Prepare an LCS for each batch by adding the correct amount of matrix spiking solution to clean methanol. (50  $\mu$ L of spike to 25 mL methanol or 10  $\mu$ L spike to 5 mL methanol). Refer to Section 17.8 for Michigan project criteria.
- 8.5.7. Shake the samples for two minutes to distribute the methanol throughout the soil.
- 8.5.8. Allow to settle, then remove a portion of methanol and store in a clean Teflon capped vial at  $4 \pm 2^\circ\text{C}$  until analysis.
- 8.6. Sample collection for medium level analysis using field methanol preservation
  - 8.6.1. Prepare a 2 oz sample container by adding 25 mL purge and trap grade methanol. (If a 5 g sample is to be used, add 5 mL methanol to a 2 oz container or VOA vial).
  - 8.6.2. Seal the bottle and attach a label.
  - 8.6.3. Weigh the bottle to the nearest 0.01g and note the weight on the label.
  - 8.6.4. Ship with appropriate sampling instructions.
  - 8.6.5. Each sample will require an additional bottle with no preservative for percent moisture determination.
  - 8.6.6. At client request, the methanol addition and weighing may also be performed in the field.
  - 8.6.7. When the samples are returned to the lab, obtain the weight of the soil added to the vial and note on the label.

- 
- 8.6.8. Add the correct amount of surrogate spiking mixture. (Add 25  $\mu\text{L}$  of 2500  $\mu\text{g/mL}$  solution for a nominal 25 g sample, 5  $\mu\text{L}$  for a nominal 5 g sample.) Refer to Section 17.8 for Michigan project criteria.
- 8.6.9. Add the correct amount of matrix spiking solution to the matrix spike and matrix spike duplicate samples. (Add 25  $\mu\text{L}$  of 50  $\mu\text{g/mL}$  solution for a nominal 25 g sample, 100  $\mu\text{L}$  for a nominal 5 g sample.) Reduce the volume of methanol added to ensure the final volume is 25 mL for nominal 25 g sample or 5 mL methanol for a nominal 5 g sample. Refer to Section 17.8 for Michigan project criteria.
- 8.6.10. Prepare an LCS for each batch by adding the correct amount of matrix spiking solution to clean methanol. (500  $\mu\text{L}$  of spike to 25 mL methanol or 100  $\mu\text{L}$  spike to 5 mL methanol). Refer to Section 17.8 for Michigan project criteria.
- 8.6.11. Shake the samples for two minutes to distribute the methanol throughout the soil.
- 8.6.12. Allow to settle, then remove a portion of methanol and store in a clean Teflon capped vial at  $4\pm 2^\circ\text{C}$  until analysis.
- 8.7. Low level procedure
- 8.7.1. If low detection limits are required (typically  $< 50 \mu\text{g/kg}$ ) sodium bisulfate preservation must be used. However, it is also necessary to take a sample for the medium level (field methanol preserved or using the EnCore™ sampler) procedure, in case the concentration of analytes in the soil is above the calibration range of the low level procedure.
- 8.7.2. A purge and trap autosampler capable of sampling from a sealed vial is required for analysis of samples collected using this method. (Varian Archon or O.I. 4552).
- 8.7.3. The soil sample is taken using a 5g EnCore™ sampling device and returned to the lab. It is recommended that two EnCore™ samplers be used for each field sample position, to allow for any reruns that may be necessary. A separate sample for % moisture determination is also necessary.
- 8.7.4. Prepare VOA vials by adding a magnetic stir bar, approximately 1 g of sodium bisulfate and 5 mL of reagent water.
- 8.7.5. Seal and label the vial. It is strongly recommended that the vial is labeled with an indelible marker rather than a paper label, since paper labels may cause the autosampler to bind and malfunction. The label absolutely must not cover the neck of the vial or the autosampler will malfunction.

8.7.6. Weigh the vial to the nearest 0.1g and note the weight on the label.

8.7.7. Extrude the soil sample from the EnCore™ sampler into the prepared VOA vial. Reweigh the vial to obtain the weight of soil and note on the label.

**Note:** Soils containing carbonates may effervesce when added to the sodium bisulfate solution. If this is the case at a specific site, add 5 mL of water instead, and freeze at  $<-10^{\circ}\text{C}$  within 48 hours, analyzed within 12 days after preserving with water, and stored at a 45 degree angle in the freezer.

**Note:** Freezing is not allowed for Ohio VAP soil samples.

8.7.9. Alternatively the sodium bisulfate preservation may be performed in the field. This is not recommended because of the many problems that can occur in the field setting. Ship at least two vials per sample. The field samplers must determine the weight of soil sampled. Each sample will require an additional bottle with no preservative for percent moisture determination, and an additional bottle preserved with methanol for the medium level procedure. Depending on the type of soil it may also be necessary to ship vials with no or extra preservative.

#### 8.8. *Unpreserved soils*

8.8.1. *At specific client request unpreserved soils packed into glass jars or brass tubes may be accepted and subsampled in the lab. This is the old procedure based on method 5030A and method 8260A. It is no longer included in SW846 and is likely to generate results that are biased low, possibly be more than an order of magnitude.*


8.9. Aqueous samples are stored in glass containers with Teflon lined septa at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$ , with minimum headspace.

8.10. Medium level solid extracts are aliquoted into 2 - 5 mL glass vials with Teflon lined caps and stored at  $4^{\circ}\text{C} \pm 2^{\circ}\text{C}$ . The extracts are stored with minimum headspace.


8.11. The maximum holding time is 14 days from sampling until the sample is analyzed. (Samples that are found to be unpreserved still have a 14 day holding time. However they should be analyzed as soon as possible. The lack of preservation should be addressed in the case narrative). Maximum holding time for the EnCore™ sampler (before the sample is added to methanol or sodium bisulfate) is 48 hours.


8.12. A holding blank is stored with the samples. This is analyzed and replaced if any of the trip blanks show any contamination. Otherwise it is replaced every 14 days.

**EnCore procedure when low level is not required (field steps in gray)**




**EnCore procedure when low level is required**



**Field methanol extraction procedure (field steps in gray)**

### Field bisulfate preservation procedure (field steps in gray)



## 9. QUALITY CONTROL

### 9.1. Initial Demonstration of Capability

- 9.1.1. For the standard analyte list, the initial demonstration described in Section 13 and method detection limit (MDL) studies must be acceptable before analysis of samples may begin. MDLs should be analyzed for low and medium soils and aqueous samples.
- 9.1.2. For non-standard analytes, a MDL study must be performed and calibration curve generated before analyzing any samples, unless lesser requirements are previously agreed to with the client. In any event, the minimum initial demonstration required is analysis of a standard at the reporting limit and a single point calibration.

### 9.2. Control Limits

In-house historical control limits must be determined for surrogates, matrix spikes, and laboratory control samples (LCS). These limits must be determined at least annually. The recovery limits are mean recovery  $\pm$  3 standard deviations for surrogates, matrix spikes and LCS. Precision limits for matrix spikes / matrix spike duplicates are 0 to mean relative percent difference + 3 standard deviations.

- 9.2.1. All surrogate, LCS, and MS recoveries (except for dilutions) must be entered into QuantIMS (when available) or other database so that accurate historical control limits can be generated. For tests without a separate extraction, surrogates and matrix spikes will be reported for all dilutions.

- 9.2.2. Refer to the QC Program document (QA-003) for further details of control limits.

### 9.3. Surrogates

Every sample, blank, and QC sample is spiked with surrogates. Surrogate recoveries in samples, blanks, and QC samples must be assessed to ensure that recoveries are within established limits. The compounds included in the surrogate spiking solutions are listed in Table 8. If any surrogates are outside limits, the following corrective actions must take place (except for dilutions):

- Check all calculations for error.
- Ensure that instrument performance is acceptable.

- Recalculate the data and/or reanalyze if either of the above checks reveal a problem.
- Reprepare and reanalyze the sample or flag the data as "Estimated Concentration" if neither of the above resolves the problem.

The decision to reanalyze or flag the data should be made in consultation with the client. It is only necessary to reprepare/reanalyze a sample once to demonstrate that poor surrogate recovery is due to matrix effect, unless the analyst believes that the repeated out of control results are not due to matrix effect.

- 9.3.1. If the surrogates are out of control for the sample, matrix spike, and matrix spike duplicate, then matrix effect has been demonstrated for that sample and reparation is not necessary. If the sample is out of control and the MS and/or MSD is in control, then reanalysis or flagging of the data is required.
- 9.3.2. Refer to the STL QC Program document (QA-003) for further details of the corrective actions.

#### 9.4. Method Blanks

- 9.4.1. For each batch of samples, analyze a method blank. The method blank is analyzed after the calibration standards, normally before any samples. For low-level volatiles, the method blank consists of reagent water. For medium-level volatiles, the method blank consists of 25.0 mL of methanol. Surrogates are added and the method blank is carried through the entire analytical procedure. The method blank must not contain any analyte of interest at or above the reporting limit (except common laboratory contaminants, see below) or at or above 5% of the measured concentration of that analyte in the associated samples, whichever is higher.
  - If the analyte is a common laboratory contaminant (methylene chloride, acetone, 2-butanone) the data may be reported with qualifiers if the concentration of the analyte is less than five times the reporting limit. Such action must be taken in consultation with the client.
  - Reanalysis of samples associated with an unacceptable method blank is required when reportable concentrations are determined in the samples.
  - If there is no target analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. Such action should be done in consultation with the client.
- 9.4.2. The method blank must have acceptable surrogate recoveries. If surrogate recoveries are not acceptable, the data must be evaluated to determine if the method blank has served the



purpose of demonstrating that the analysis is free of contamination. If surrogate recoveries are low and there are reportable analytes in the associated samples re-extraction of the blank and affected samples will normally be required. Consultation with the client should take place.

9.4.3. If reanalysis of the batch is not possible due to limited sample volume or other constraints, the method blank is reported, all associated samples are flagged with a "B," and appropriate comments may be made in a narrative to provide further documentation.

9.4.4. Refer to the STL QC Program document (QA-003) for further details of the corrective actions.

#### 9.5. Laboratory Control Samples (LCS)

9.5.1. For each batch of samples, analyze a LCS. The LCS is analyzed after the calibration standard, and normally before any samples. The LCS contains a representative subset of the analytes of interest (See Table 9), and must contain the same analytes as the matrix spike. If any analyte or surrogate is outside established control limits, the system is out of control and corrective action must occur. Corrective action will normally be re-preparation and reanalysis of the batch.

- If the batch is not re-extracted and reanalyzed, the reasons for accepting the batch must be clearly presented in the project records and the report. (Examples of acceptable reasons for not reanalyzing might be that the matrix spike and matrix spike duplicate are acceptable, and sample surrogate recoveries are good, demonstrating that the problem was confined to the LCS.)
- If re-extraction and reanalysis of the batch is not possible due to limited sample volume or other constraints, the LCS is reported, all associated samples are flagged, and appropriate comments are made in a narrative to provide further documentation.

9.5.2. Refer to the STL QC Program document (QA-003) for further details of the corrective action.

9.5.3. If full analyte spike lists are used at client request, it will be necessary to allow a percentage of the components to be outside control limits as this would be expected statistically. These requirements should be negotiated with the client. Refer to Section 17.5 for Ohio VAP specific analytes.

#### 9.6. Matrix Spikes

9.6.1. For each QC batch, analyze a matrix spike and matrix spike duplicate. Spiking compounds and levels are given in Table 9. Compare the percent recovery and relative percent difference

(RPD) to that in the laboratory specific historically generated limits. See Section 17.5 for Ohio VAP specific analytes.

- If any individual recovery or RPD falls outside the acceptable range, corrective action must occur. The initial corrective action will be to check the recovery of that analyte in the Laboratory Control Sample (LCS). Generally, if the recovery of the analyte in the LCS is within limits, then the laboratory operation is in control and analysis may proceed. The reasons for accepting the batch must be documented.
- If the recovery for any component is outside QC limits for both the matrix spike/ spike duplicate and the LCS, the laboratory is out of control and corrective action must be taken. Corrective action will normally include reanalysis of the batch.
- If a MS/MSD is not possible due to limited sample, then a LCS duplicate should be analyzed. RPD of the LCS and LCSD are compared to the matrix spike limits.
- The matrix spike/duplicate must be analyzed at the same dilution as the unspiked sample, even if the matrix spike compounds will be diluted out.

#### 9.7. Nonconformance and Corrective Action

- 9.7.1. Any deviations from QC procedures must be documented as a nonconformance, with applicable cause and corrective action approved by the facility QA Manager.

#### 9.8. Quality Assurance Summaries

Certain clients may require specific project or program QC which may supersede these method requirements. Quality Assurance Summaries should be developed to address these requirements.

#### 9.9. STL QC Program

Further details of QC and corrective action guidelines are presented in the STL QC Program document (QA-003). Refer to this document if in doubt regarding corrective actions.

### 10. CALIBRATION AND STANDARDIZATION

#### 10.1. Summary

- 10.1.1. Prior to the analysis of samples and blanks, each GC/MS system must be tuned and calibrated. Hardware tuning is checked through the analysis of the 4-Bromofluorobenzene (BFB) to establish that a given GC/MS system meets the standard mass spectral abundance criteria. The

GC/MS system must be calibrated initially at a minimum of five concentrations (analyzed under the same BFB tune), to determine the linearity of the response utilizing target calibration standards. Once the system has been calibrated, the calibration must be verified each twelve hour time period for each GC/MS system.

#### 10.2.1. General

Electron Energy:	70 volts (nominal)
Mass Range:	35–300 AMU
Scan Time:	to give at least 5 scans/peak, but not to exceed 2 second/scan
Injector Temperature:	200–250°C
Source Temperature:	According to manufacturer's specifications
Transfer Line	Temperature: 250–300°C
Purge Flow:	40 mL/minute
Carrier Gas	Flow: 15 mL/minute
Make-up Gas Flow:	25–30 mL/minute

#### 10.2.2. Gas chromatograph suggested temperature program

##### 10.2.2.1. BFB Analysis

Isothermal:	170°C
-------------	-------

##### 10.2.2.2. Sample Analysis

Initial Temperature:	40°C
Initial Hold Time:	4 minutes
Temperature Program:	8°C/minute
Final Temperature:	184°C
Second Temperature	Program: 40°C/minute
Final Temperature:	240°C
Final Hold Time:	2.6 minutes

#### 10.3. Instrument Tuning

10.3.1. Each GC/MS system must be hardware-tuned to meet the abundance criteria listed in Table 10 for a maximum of a 50 ng injection or purging of BFB. Analysis must not begin until these criteria are met. These criteria must be met for each twelve-hour time period. The twelve-hour time period begins at the moment of injection of BFB.

#### 10.4. Initial Calibration

- 10.4.1. A series of five initial calibration standards is prepared and analyzed for the target compounds and each surrogate compound. Six standards must be used for a quadratic least squares calibration. Suggested calibration levels for a 5 mL purge are: 5, 20, 50, 100, and 200 µg/L. Certain analytes are prepared at higher concentrations due to poor purge performance. Suggested calibration levels for a 25 mL purge are 1, 5, 10, 20, and 40 µg/L. Again, some analytes are prepared at higher levels. Tables 2 and 4 list the calibration levels for each analyte. Other calibration levels and purge volumes may be used depending on the capabilities of the specific instrument. (For example, adequate sensitivity can be obtained on the Agilent 5973 instruments to use a 5 mL purge volume to reach the same reporting limits that once required a 25 mL purge. The calibration levels will still be the same 1, 5, 10, 20, 40 µg/L.) However, the same purge volume must be used for calibration and sample analysis, and the low level standard must be at or below the reporting limit.
- 10.4.2. It may be necessary to analyze more than one set of calibration standards to encompass all of the analytes required for same tests. For example, the Appendix IX list requires the Primary standard (Table 5) and the Appendix IX standard (Table 6). If acceptable analytical performance can be obtained the primary and appendix IX standards may be analyzed together.
- 10.4.3. Internal standard calibration is used. The internal standards are listed in Table 7. Target compounds should reference the nearest internal standard. Each calibration standard is analyzed and the response factor (RF) for each compound is calculated using the area response of the characteristic ions against the concentration for each compound and internal standard. See equation 1, Section 12, for calculation of response factor.
- 10.4.4. The % RSD of the calibration check compounds (CCC) must be less than 30%. Refer to Table 12 for the CCCs.
- 10.4.4.1. If none of the CCCs are required analytes, project specific calibration specifications must be agreed with the client.
- 10.4.5. The average RF must be calculated for each compound. A system performance check is made prior to using the calibration curve. The five system performance check compounds (SPCC) are checked for a minimum average response factor. Refer to Table 11 for the SPCC compounds and required minimum response factors.
- 10.4.6. If the average of all the %RSDs in the calibration is  $\leq 15\%$ , then all analytes may use average response factor for calibration.

10.4.6.1. If the software in use is capable of routinely reporting curve coefficients for data validation purposes, and the necessary calibration reports can be generated, then the analyst should evaluate analytes with %RSD > 15% for calibration on a curve. If it appears that substantially better accuracy would be obtained using quantitation from a curve then the appropriate curve should be used for quantitation. If Relative Standard Error (RSE) is used to evaluate the curve it must be better than 15%. Otherwise the correlation coefficient (coefficient of determination for non-linear curves) must be  $\geq 0.990$ .

10.4.6.2. If the average of all the %RSDs in the calibration is > 15% then calibration on a curve must be used for all analytes with %RSD > 15%. Linear or quadratic curve fits may be used. The analyst should consider instrument maintenance to improve the linearity of response. If Relative Standard Error (RSE) is used to evaluate the curve it must be better than 15%. If the % RSD is > 15%, the analyst may drop the low or high in the ICAL, as long as a minimum of 5 points are maintained and the quantitation range is adjusted accordingly. Otherwise the correlation coefficient, (coefficient of determination,  $r^2$  for non-linear curves) must be  $\geq 0.990$ . If the correlation coefficient is < 0.990, then any hit for these compounds must be flagged as estimated.

10.4.6.3. Refer to Section 17.5 for specific Ohio VAP criteria.

#### 10.4.7. Weighting of data points

In a linear or quadratic calibration fit, the points at the lower end of the calibration curve have less weight in determining the curve generated than points at the high concentration end of the curve. However, in environmental analysis, accuracy at the low end of the curve is very important. For this reason it is preferable to increase the weighting of the lower concentration points.  $1/\text{Concentration}^2$  weighting (often called  $1/X^2$  weighting) will improve accuracy at the low end of the curve and should be used if the data system has this capability.

10.4.8. If time remains in the 12-hour period initiated by the BFB injection before the initial calibration, samples may be analyzed. Otherwise, proceed to continuing calibration.

10.4.9. The calibration standards for the initial 5-point calibration for low level soils that are not preserved in sodium bisulfate (i.e. are preserved by freezing, or not preserved) must be heated to 40°C for purging. Using this calibration curve for water samples is acceptable as long as all calibration, QC, and samples are also heated to 40°C. A separate five point calibration must be prepared for analysis of low level soils that are preserved with sodium bisulfate. Low level soils analysis requires the use of a closed vial autosampler such as the Varian Archon, O.I. 4552 or Tekmar Precept. Each standard for analysis of sodium

bisulfate preserved samples is prepared by spiking the methanolic standard solution through the septum of a VOA vial containing 5 mL of water and 1 g sodium bisulfate. The standards are heated to 40°C for purging. All low-level soil samples, standards, and blanks must also be heated to 40°C for purging. Medium soil extracts should be analyzed using the water (unheated or optionally heated) calibration curve as long as all calibration standards, samples, and QC samples are purged at the same temperature.

10.4.10. Non-standard analytes are sometimes requested. For these analytes, it is acceptable to analyze a single standard at the reporting limit with each continuing calibration rather than a five point initial calibration. If the analyte is detected in any of the samples, a five point initial calibration must be generated and the sample(s) reanalyzed for quantitation. However, if the analyte is not detected, the non-detect may be reported and no further action is necessary.

**Note:** This procedure is may not be used for Ohio VAP samples.

10.5. Continuing Calibration: The initial calibration must be verified every twelve hours.

10.5.1. Continuing calibration begins with analysis of BFB as described in Section 10.3. If the system tune is acceptable, the continuing calibration standard(s) are analyzed. The level 3 calibration standard is used as the continuing calibration.

10.5.2. The RF data from the standards are compared with the average RF from the initial five-point calibration to determine the percent drift of the CCC compounds. The calculation is given in equation 4, Section 12.3.4.

10.5.3. The % drift of the CCCs must be  $\leq 20\%$  for the continuing calibration to be valid. The SPCCs are also monitored. The SPCCs must meet the criteria described in Table 11. In addition, the % drift of all analytes must be  $\leq 50\%$  with allowance for up to six target analytes to have % drift  $> 50\%$ .

10.5.3.1. If none of the CCCs are required analytes, project specific calibration specifications must be agreed with the client.

10.5.3.2. Cyclohexanone, one of the components of the Appendix IX standard, is unstable in the calibration solution, forming 1,1-dimethoxycyclohexane. No calibration criteria are applied to cyclohexanone and quantitation is tentative. Cyclohexanone is included on the Universal Treatment Standard and FO-39 regulatory lists (but not on Appendix IX).

10.5.3.3. Refer to Table 12 for specific Ohio VAP analytes.

10.5.4. If the CCCs and or the SPCCs do not meet the criteria in Sections 10.5.3 and 10.5.4, the system must be evaluated and corrective action must be taken. The BFB tune and continuing calibration must be acceptable before analysis begins. Extensive corrective action such as a different type of column will require a new initial calibration.

10.5.5. Once the above criteria have been met, sample analysis may begin. **Initial calibration average RFs (or the calibration curve) will be used for sample quantitation, not the continuing calibration RFs.** Analysis may proceed until 12 hours from the injection of the BFB have passed. (A sample *desorbed* less than or equal to 12 hours after the BFB is acceptable.)

## 11. PROCEDURE

### 11.1. Procedural Variations

11.1.1. One time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation shall be completely documented using a Nonconformance Memo and approved by a Supervisor or group leader and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file.

11.1.2. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

### 11.2. Preliminary Evaluation

11.2.1. Where possible, samples are screened by headspace or GC/MS off-tune analysis to determine the correct aliquot for analysis. Alternatively, an appropriate aliquot can be determined from sample histories.

11.2.2. Dilutions should be done just prior to the GC/MS analysis of the sample. Dilutions are made in volumetric flasks or in a Lucrlok syringe. Calculate the volume of reagent water required for the dilution. Fill the syringe with reagent water, compress the water to vent any residual air and adjust the water volume to the desired amount. Adjust the plunger to the mark and inject the proper aliquot of sample into the syringe. If the dilution required would use less than 1  $\mu$ L of sample then serial dilutions must be made in volumetric flasks.

11.2.2.1. The diluted concentration is to be estimated to be in the upper half of the calibration range.

### 11.3. Sample Analysis Procedure

- 11.3.1. All analysis conditions for samples must be the same as for the continuing calibration standards (including purge time and flow, desorb time and temperature, column temperatures, multiplier setting etc.).
- 11.3.2. All samples must be analyzed as part of a batch. The batch is a set of up to 20 samples of the same matrix processed using the same procedures and reagents within the same time period. The batch also must contain a MS/MSD, a LCS, and a method blank.
  - 11.3.2.1. If there is insufficient time in the 12-hour tune period to analyze 20 samples, the batch may be continued into the next tune period. However, if any re-tuning of the instrument is necessary, or if a period of greater than 24 hours from the preceding BFB tune has passed, a new batch must be started. For medium level soils the batch is defined at the sample preparation stage.
  - 11.3.2.2. Laboratory generated QC samples (Blank, LCS, MS/MSD) do not count towards the maximum 20 samples in a batch. Field QC samples are included in the batch count.
  - 11.3.2.3. It is not necessary to reanalyze batch QC with reanalyses of samples. However, any reruns must be as part of a valid batch.

### 11.4. Water Samples

- 11.4.1. All samples and standard solutions must be at ambient temperature before analysis.
- 11.4.2. Fill a syringe with the sample. If a dilution is necessary it may be made in the syringe if the sample aliquot is  $\geq 5 \mu\text{L}$ . Check and document the pH of the remaining sample.
- 11.4.3. Add 250 ng of each internal and surrogate standard (10  $\mu\text{L}$  of a 25  $\mu\text{g/mL}$  solution, refer to Tables 7, 8 and 16). The internal standards and the surrogate standards may be mixed and added as one spiking solution (this results in a 50  $\mu\text{g/L}$  solution for a 5 mL sample, and a 10  $\mu\text{g/L}$  solution for a 25 mL sample). Inject the sample into the purging chamber.
  - 11.4.3.1. For TCLP samples use 0.5 mL of TCLP leachate with 4.5 mL reagent water and spike with 10  $\mu\text{L}$  of the 25  $\mu\text{g/mL}$  TCLP spiking solution. (Note that TCLP reporting limits will be 10 times higher than the corresponding aqueous limits).
- 11.4.4. Purge the sample for eleven minutes (the trap must be below 35°C).



11.4.5. After purging is complete, desorb the sample, start the GC temperature program, and begin data acquisition. After desorption, bake the trap for approximately 3-10 minutes to condition it for the next analysis. When the trap is cool, it is ready for the next sample.

11.4.6. Desorb and bake time and temperature are optimized for the type of trap in use. The same conditions must be used for samples and standards.

#### 11.5. Methanol Extract Soils

11.5.1. Rinse a gas-tight syringe with organic free water. Fill the syringe with the same volume of organic free water as used in the calibrations. Add no more than 2% (v/v) (100  $\mu$ L for a 5 mL purge) methanolic extract (from Section 8.5 or 8.6) to the syringe. Add internal standard (if used). Load the sample onto the purge and trap device and analyze as for aqueous samples. If less than 5  $\mu$ L of methanolic extract is to be added to the water, dilute the methanolic extract such that a volume greater than 5  $\mu$ L will be added to the water in the syringe. Refer to Section 17.8 for Michigan project requirements.

#### 11.6. Liquid wastes that are soluble in methanol and insoluble in water.

11.6.1. Pipet 2 mL of the sample into a tared vial. Use a top-loading balance. Record the weight to the nearest 0.1 gram.

11.6.2. Quickly add 7 mL of methanol, then add 1 mL of surrogate spiking solution to bring the final volume to 10 mL. Cap the vial and shake for 2 minutes to mix thoroughly. For a MS/MSD or LCS, 6 mL of methanol, 1 mL of surrogate solution, and 1 mL of matrix spike solution is used.

11.6.3. Rinse a gas-tight syringe with organic free water. Fill the syringe with the same volume of organic free water as used in the calibrations. Add no more than 2% (v/v) (100  $\mu$ L for a 5 mL purge) methanolic extract (from Section 8.5 or 8.6) to the syringe. Add internal standard (if used). Load the sample onto the purge and trap device and analyze as for aqueous samples. If less than 5  $\mu$ L of methanolic extract is to be added to the water, dilute the methanolic extract such that a volume greater than 5  $\mu$ L will be added to the water in the syringe.

#### 11.7. Aqueous and Low level Soil Sample Analysis (Purge and Trap units that sample directly from the VOA vial)

11.7.1. Units which sample from the VOA vial should be equipped with a module which automatically adds surrogate and internal standard solution to the sample prior to purging the sample.

- 11.7.2. If the autosampler uses automatic IS/SS injection, no further preparation of the VOA vial is needed. Otherwise the internal and surrogate standards must be added to the vial. *Note:* Aqueous samples with high amounts of sediment present in the vial may not be suitable for analysis on this instrumentation, or they may need to be analyzed as soils.
- 11.7.3. Soil samples must be quantitated against a curve prepared with standards containing about the same amount of sodium bisulfate as the samples (1 g in 5 mL).
- 11.7.4. Sample remaining in the vial after sampling with one of these mechanisms is no longer valid for further analysis. A fresh VOA vial must be used for further sample analysis.
- 11.7.5. For aqueous samples, check the pH of the sample remaining in the VOA vial after analysis is completed.

11.8. *Low-Level Solids Analysis using discrete autosamplers, Method 8260A, 5030A.*

**Note: This technique may seriously underestimate analyte concentration and must not be used except at specific client request for the purpose of comparability with previous data. It is no longer part of SW-846.**

*This method is based on purging a heated soil/sediment sample mixed with reagent water containing the surrogates and internal standards. Analyze all reagent blanks and standards under the same conditions as the samples (e.g., heated). The calibration curve is also heated during analysis. Purge temperature is 40°C.*

- 11.8.1. *Do not discard any supernatant liquids. Mix the contents of the container with a narrow metal spatula.*
- 11.8.2. *Weigh out 5 g (or other appropriate aliquot) of sample into a disposable culture tube or other purge vessel. Record the weight to the nearest 0.1 g. If method sensitivity is demonstrated, a smaller aliquot may be used. Do not use aliquots less than 1.0 g. If the sample is contaminated with analytes such that a purge amount less than 1.0 g is appropriate, use the medium level method. For the medium level method, add 4g soil to 10 mL methanol containing the surrogates, mix for two minutes, allow to settle then remove a portion of the methanol and store in a clean Teflon capped vial at 4°C until analysis. Analyze as described in section 11.5.*
- 11.8.3. *Connect the purge vessel to the purge and trap device.*
- 11.8.4. *Rinse a 5 mL gas-tight syringe with organic free water, and fill. Compress to 5 mL. Add surrogate/internal standard (and matrix spike solutions if required.). Add directly to the sample from 11.5.2.*

*11.8.5. The above steps should be performed rapidly and without interruption to avoid loss of volatile organics.*

*11.8.6. Add the heater jacket or other heating device and start the purge and trap unit.*

*11.8.7. Soil samples that have low IS recovery when analyzed (<50%) should be reanalyzed once to confirm matrix effect.*

#### 11.9. Medium-Level Soil/Sediment and Waste Samples

11.9.5. Sediments/soils and waste that are insoluble in methanol.

11.9.5.1. Sediments/soils and waste that are insoluble in methanol.

11.9.5.1.1. Gently mix the contents of the sample container with a narrow metal or wood spatula. Weigh 4 g (wet weight) into a tared vial. Use a top-loading balance. Record the weight to 0.1 gram. Do not discard any supernatant liquids.

11.9.5.1.2. Quickly add 9 mL of methanol, and 1 mL of surrogate spiking solution to bring the final volume of methanol to 10 mL. For an LCS or MS/MSD sample add 8 mL of methanol, 1 mL of surrogate spike solution, and 1 mL of matrix spike solution. Cap the vial and vortex to mix thoroughly.

NOTE: Sections 11.9.5.1.1 and 11.9.5.1.2 must be performed rapidly and without interruption to avoid the loss of volatile organics.

#### 11.10. Initial review and corrective actions

11.10.1. If the retention time for any internal standard in the continuing calibration changes by more than 0.5 minutes from the mid-level initial calibration standard, the chromatographic system must be inspected for malfunctions and corrected. Reanalysis of samples analyzed while the system was malfunctioning is required.

11.10.2. If the internal standard response in the continuing calibration is more than 200% or less than 50% of the response in the mid-level of the initial calibration standard, the chromatographic system must be inspected for malfunctions and corrected. Reanalysis of samples analyzed while the system was malfunctioning is required.

11.10.2.1. Any samples that do not meet the internal standard criteria for the continuing calibration must be evaluated for validity. If the change in sensitivity is a matrix effect confined to an individual sample reanalysis is not necessary. If the change in sensitivity is

due to instrumental problems all affected samples must be reanalyzed after the problem is corrected.

- 11.10.3. The surrogate standard recoveries are evaluated to ensure that they are within limits. Corrective action for surrogates out of control will normally be to reanalyze the affected samples. However, if the surrogate standard response is out high and there are no target analytes or tentatively identified compounds, reanalysis may not be necessary. Out of control surrogate standard response may be a matrix effect. It is only necessary to reanalyze a sample once to demonstrate matrix effect, but reanalysis at a dilution should be considered.

#### 11.11. Dilutions

If the response for any compound exceeds the working range of the GC/MS system, a dilution of the extract is prepared and analyzed. An appropriate dilution should be in the upper half of the calibration range. Samples may be screened to determine the appropriate dilution for the initial run. If the initial diluted run has no hits or hits below 20% of the calibration range and the matrix allows for analysis at a lesser dilution, then the sample must be reanalyzed at a dilution targeted to bring the largest hit above 50% of the calibration range.

##### 11.10.1. Guidance for Dilutions Due to Matrix

If the sample is initially run at a dilution and the baseline rise is less than half the height of the internal standards, or if individual non target peaks are less than twice the height of the internal standards, then the sample should be reanalyzed at a more concentrated dilution. This requirement is approximate and subject to analyst judgement.

##### 11.10.2. Reporting Dilutions

The most concentrated dilution with no target compounds above the calibration range will be reported. Other dilutions will only be reported at client request.

## 12. DATA ANALYSIS AND CALCULATIONS

### 12.1. Qualitative identification

An analyte is identified by retention time and by comparison of the sample mass spectrum with the mass spectrum of a standard of the suspected compound (standard reference spectrum). Mass spectra for standard reference may be obtained on the user's GC/MS by analysis of the calibration standards or from the NIST Library. Two criteria must be satisfied to verify identification: (1) elution of sample component at the same GC retention time as the standard component; and (2) correspondence of the sample component and the standard component

---

characteristic ions. (Note: Care must be taken to ensure that spectral distortion due to co-elution is evaluated.)

- The sample component retention time must compare to within  $\pm 0.2$  min. of the retention time of the standard component. For reference, the standard must be run within the same twelve hours as the sample.
- All ions present in the standard mass spectra at a relative intensity greater than 10% (most abundant ion in the spectrum equals 100%) should be present in the sample spectrum.
- The relative intensities of ions should agree to within  $\pm 30\%$  between the standard and sample spectra. (Example: For an ion with an abundance of 50% in the standard spectra, the corresponding sample abundance must be between 20 and 80 percent.)

12.1.1. If a compound cannot be verified by all the above criteria, but in the technical judgment of the analyst, the identification is correct, then the analyst shall report that identification and proceed with quantitation.

## 12.2. Tentatively Identified Compounds (TICs)

12.2.1. If the client requests components not associated with the calibration standards, a search of the NIST library may be made for the purpose of tentative identification. Guidelines are:

- 12.2.1.1. Relative intensities of major ions in the reference spectrum (ions  $> 10\%$  of the most abundant ion) should be present in the sample spectrum.
- 12.2.1.2. The relative intensities of the major ions should agree to within 20%. (Example: If an ion shows an abundance of 50% in the standard spectrum, the corresponding sample ion abundance must be between 30% and 70%).
- 12.2.1.3. Molecular ions present in the reference spectrum should be present in the sample spectrum.
- 12.2.1.4. Ions present in the sample spectrum but not in the reference spectrum should be reviewed for possible background contamination or presence of coeluting compounds.
- 12.2.1.5. Ions present in the reference spectrum but not in the sample spectrum should be reviewed for possible subtraction from the spectrum because of background contamination or coeluting peaks. (Data system reduction programs can sometimes create these discrepancies.)

12.2.1.6. Computer-generated library search routines should not use normalization routines that would misrepresent the library or unknown spectra when compared to each other. Only after visual inspection of the sample with the nearest library searches should the analyst assign a tentative identification.

### 12.3. Calculations.

#### 12.3.1. Response factor (RF):

##### **Equation 1**

$$RF = \frac{A_x C_{is}}{A_{is} C_x}$$

Where:

$A_x$  = Area of the characteristic ion for the compound to be measured

$A_{is}$  = Area of the characteristic ion for the specific internal standard

$C_{is}$  = Concentration of the specific internal standard, ng

$C_x$  = Concentration of the compound being measured, ng

#### 12.3.2. Standard deviation (SD):

##### **Equation 2**

$$SD = \sqrt{\sum_{i=1}^N \frac{(X_i - X)^2}{N - 1}}$$

$X_i$  = Value of X at i through N

$N$  = Number of points

$X$  = Average value of  $X_i$

#### 12.3.3. Percent relative standard deviation (%RSD):

##### **Equation 3**

$$\%RSD = \frac{\text{Standard Deviation}}{RF_i} \times 100$$

$RF_i$  = Mean of RF values in the curve

12.3.4. Percent drift between the initial calibration and the continuing calibration:

**Equation 4**

$$\% \text{ Drift} = \frac{C_{\text{expected}} - C_{\text{found}}}{C_{\text{expected}}} \times 100$$

Where

$C_{\text{expected}}$  = Known concentration in standard

$C_{\text{found}}$  = Measured concentration using selected quantitation method

12.3.5. Target compound and surrogate concentrations:

Concentrations in the sample may be determined from linear or second order (quadratic) curve fitted to the initial calibration points, or from the average response factor of the initial calibration points. Average response factor may only be used when the % RSD of the response factors in the initial calibration is  $\leq 15\%$ .

12.3.5.1. Calculation of concentration using Average Response Factors

**Equation 5**

$$\text{Concentration } \mu\text{g} / \text{L} = \frac{x}{RF}$$

12.3.5.2. Calculation of concentration using Linear fit

**Equation 6**

$$\text{Concentration } \mu\text{g} / \text{L} = A + Bx$$

12.3.5.3. Calculation of concentration using Quadratic fit

**Equation 7**

$$\text{Concentration } \mu\text{g} / \text{L} = A + Bx + Cx^2$$

$x$  is defined in equations 8, 9 and 10

$A$  is a constant defined by the intercept

$B$  is the slope of the curve

$C$  is the curvature

12.3.5.4. Calculation of  $x$  for Water and water-miscible waste:

**Equation 8**

$$x = \frac{(A_x)(I_s)(D_f)}{(A_{is})(V_o)}$$

Where:

$A_x$  = Area of characteristic ion for the compound being measured (secondary ion quantitation is allowed only when there are sample interferences with the primary ion)

$A_{is}$  = Area of the characteristic ion for the internal standard

$I_s$  = Amount of internal standard added in ng

$$\text{Dilution Factor} = D_f = \frac{\text{Total volume purged (mL)}}{\text{Volume of original sample used (mL)}}$$

$V_o$  = Volume of water purged, mL

12.3.5.5. Calculation of  $x$  for Medium level soils:

**Equation 9**

$$x = \frac{(A_x)(I_s)(V_t)(1000)(D_f)}{(A_{is})(V_a)(W_s)(D)}$$

Where:

$A_x$ ,  $I_s$ ,  $D_f$ ,  $A_{is}$ , same as for water.

$V_t$  = Volume of total extract, mL (Typically 25 mL)

$V_a$  = Volume of extract added for purging,  $\mu\text{L}$

$W_s$  = Weight of sample extracted, g

$$D = \frac{100 - \% \text{moisture}}{100}$$



12.3.5.6. Calculation of  $x$  for Low level soils:

**Equation 10**

$$x = \frac{(A_x)(I_s)}{(A_{is})(W_s)(D)}$$

Where:

$A_x$ ,  $I_s$ ,  $A_{is}$ , same as for water.

$D$  is as for medium level soils

$W_s$  = Weight of sample added to the purge vessel, g

12.3.5.7. Calculation of TICs: The calculation of TICs (tentatively identified compounds) is identical to the above calculations with the following exceptions:

$A_x$  = Area in the total ion chromatogram for the compound being measured

$A_{is}$  = Area of the total ion chromatogram for the nearest internal standard without interference

$RF = 1$

In other words, the concentration is equal to  $x$  as defined in equations 8, 9 and 10.

12.3.6. MS/MSD Recovery

**Equation 11**

$$\text{Matrix Spike Recovery, \%} = \frac{SSR - SR}{SA} \times 100$$

$SSR$  = Spike sample result

$SR$  = Sample result

$SA$  = Spike added

12.3.7. Relative % Difference calculation for the MS/MSD

**Equation 12**

$$RPD = \frac{|MSR - MSDR|}{\frac{1}{2}(MSR + MSDR)} \times 100$$

Where:

**RPD** = Relative percent difference

**MSR** = Matrix spike result

**MSDR** = Matrix spike duplicate result

### 13. METHOD PERFORMANCE

#### 13.1. Method Detection Limit

- 13.2. Generally, each laboratory must generate a valid method detection limit for each analyte of interest. The MDL must be below the reporting limit for each analyte. The procedure for determination of the method detection limit is given in 40 CFR Part 136, Appendix B, and further defined in QA Policy #: QA-005. When non-standard compounds are analyzed at client request, lesser requirements are possible with client agreement. At a minimum, a standard at the reporting limit must be analyzed to demonstrate the capability of the method.

#### 13.3. Initial Demonstration

- 13.4. Each laboratory must make a one time initial demonstration of capability for each individual method. Demonstration of capability for both soil and water matrices is required. This requires analysis of QC check samples containing all of the standard analytes for the method. For some tests it may be necessary to use more than one QC check mix to cover all analytes of interest. The QC check sample is made up at 20 µg/L. (Some compounds will be at higher levels, refer to the calibration standard levels for guidance.)

- 13.4.1. Four aliquots of the QC check sample are analyzed using the same procedures used to analyze samples, including sample preparation.

- 13.4.2. Calculate the average recovery and standard deviation of the recovery for each analyte of interest. The %RSD should be  $\leq 15\%$  for each analyte, and the % recovery should be within 80-120%.

- 13.4.3. If any analyte does not meet the acceptance criteria, check the acceptance limits in the reference methods (Table 6 of method 8240B, paragraph 8.3.5 of method 8260A). If the recovery or precision is outside the limits in the reference methods, the test must be repeated. Only those analytes that did not meet criteria in the first test need to be evaluated. Repeated failure for any analyte indicates the need for the laboratory to evaluate the analytical procedure and take corrective action.

#### 13.4.4. Training Qualification

- 13.4.4.1. The group/team leader has the responsibility to ensure that this procedure is performed by an analyst who has been properly trained in its use and has the required experience.

---

## 14. POLLUTION PREVENTION

- 14.1. This method does not contain any specific modifications that serve to minimize or prevent pollution.

## 15. WASTE MANAGEMENT

- 15.1. Waste generated in the procedure must be segregated and disposed according to the facility hazardous waste procedures. The Health and Safety Director should be contacted if additional information is required.

## 16. REFERENCES

- 16.1. SW846, *Test Methods for Evaluating Solid Waste*, Third Edition, Gas Chromatography/Mass Spectrometry for Volatile Organics, Method 8260B, Update III, December 1996
- 16.2. SW846, *Test Methods for Evaluating Solid Waste*, Third Edition, Gas Chromatography/Mass Spectrometry for Volatile Organics, Method 8260A, Update II, September 1994.

## 17. MISCELLANEOUS

- 17.1. Modifications from the reference method
- 17.1.1. Ion 119 is used as the quantitation ion for chlorobenzene-d5 for 25 mL purge tests.
- 17.1.2. A retention time window of 0.2 minutes is used for all components, since some data systems do not have the capability of using the relative retention time units specified in the reference method.
- 17.1.3. The quantitation and qualifier ions for some compounds have been changed from those recommended in SW-846 in order to improve the reliability of qualitative identification.
- 17.1.4. Method 8260A recommends that the purge vessel is run through an additional purge cycle after 25 mL sample analysis to remove carryover. Instead, purge vessels are oven baked between analyses or disposable vessels are used one time only.
- 17.1.5. SW-846 recommends that a curve be used for any analytes with %RSD of the response factors > 15%. However, some industry standard data systems and forms generation software cannot report this data with the necessary information for data validation. In

addition most software available does not allow weighting of the curve. Unweighted curves may exhibit serious errors in quantitation at the low end, resulting in possible false positives or false negatives. Therefore, this SOP allows use of average response factors if the average %RSD for all compounds is  $\leq 15\%$ .

17.2. Modifications from previous revision

This SOP has been substantially revised to reflect the changes included in Update III to SW-846. Directions for method 524.2 and method 624 have also been added.

17.3. Facility specific SOPs

Each facility shall attach a list of facility-specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none.

17.4. Flow diagrams

17.4.1. Initial Demonstration and MDL

17.5. The following are protocols that must be followed when analyzing OhioVAP samples:

- Sections 9.5 and 9.6: n-Hexane must be spiked and reported for both the LCS and MS/MSD.
- Sections 10.4.6: All analytes must have a %RSDs  $\leq 15\%$ . Corrective action must be completed for any compounds failing the  $<15\%$  requirement.
- Section 11.1 and 17.1.5 (Method deviations) are not to be performed.
- Section 11.9.2: For OhioVAP projects, the laboratory will reanalyze any sample where the internal standard fails and there is no evidence of matrix interference.

17.6. The following are protocols that must be followed when analyzing BP Oil – Lima Refinery RFI work plan.

- Section 8.1 STL will continue to follow the 14 day holding time specified in the Corporate SOP.
- Delete for this project Section 8.3 At specific client request, unpreserved soil samples may be accepted.

- Delete for this project Section 8.8.1 At specific client request unpreserved soils packed into glass jars or brass tubes may be accepted and subsampled in the lab. This is the old procedure based on method 5030A. It is no longer included and is likely to generate results that are biased low, possibly by more than an order of magnitude.
- Modify Section 8.5.8 For the purpose of this project, the soil/methanol mixture may be stored for two days prior to analysis.
- Modify Section 8.6.12 For the purpose of this project, the soil/methanol mixture may be stored for two days prior to analysis.
- Modify (per discussion with Region V representative) to Section 10.4.6.2 Compounds with %RSD >15% are to be calibrated using an alternate calibration technique (e.g. linear or quadratic calibration curve). For poor responders, the alternate calibration technique requirements may not be met either. This sentence is added for those cases. If the correlation coefficient is < 0.990, then any hit for these compounds must be flagged as estimated.
- Modify Section 10.4.2 It is necessary to analyze the Appendix IX standard separately from the primary standard due to the presence of xylene solvent in the Appendix IX standard. Alternatively, STL will purchase the Appendix IX standard in a solvent other than xylene.
- Modify Section 10.4.9 For this project, this section will be modified to comply with the requirement of adding methanol to the calibration standards so that those standards contain the same amount of methanol as the diluted soil extracts.
- Modify Table 6
- For the project specific SOP, acetonitrile will be removed from table 6, page 49 and appended onto table 5, page 48. Acetonitrile will be calibrated as part of the STL primary standard, using a separate acetonitrile standard. This will ensure that the calibration curve for acetonitrile will be done free from any interference from allyl chloride.

17.7. The following are protocols that must be followed when analyzing South Carolina Projects only.

- **Delete** from Section 10.4.7 In a linear *or quadratic* calibration fit, the points at the lower end of the calibration curve have less weight in determining the curve generated than points at the high concentration end of the curve.
- **Delete** from Section 12.3.5 Concentrations in the sample may be determined from linear *or second order (quadratic)* curve fitted to the initial calibration points, or from the average response factor of the initial calibration points.
- **Delete** from Section 12.3.5.1

- Calculation of concentration using Quadratic fit

**Equation 13**

$$\text{Concentration } \mu\text{g} / \text{L} = A + Bx + Cx^2$$

$x$  is defined in equations 8, 9 and 10



$A$  is a constant defined by the intercept

$B$  is the slope of the curve

$C$  is the curvature

- **Change** Section 9.3 The compounds included in the surrogate spiking solutions are listed in Tables 8 and 9.

- 17.8. The following are protocols that must be followed to achieve the lower reporting limits required when analyzing Michigan projects.
- 17.8.1. Modify Section 8.5.4 and 8.6.8 (add 5 uL of 2500 ug/mL surrogate solution for a nominal 25 g sample).
- 17.8.2. Modify Section 8.5.5 and 8.6.9 (add 100 uL of 50 ug/mL spike solution for a nominal 25 g sample).
- 17.8.3. Modify Section 8.5.6 and 8.6.10 (add 100 uL of 50 ug/mL spike solution for a nominal 25g sample).
- 17.8.4. Michigan reporting limits for methanol preserved soils are achieved by injecting 100 uL of the methanol extract in a 5 mL purge. The instrument is calibrated using the recommended calibration levels in water of 1 ug/L, 2 ug/L (if a quadratic calibration is to be used), 5 ug/L, 10 ug/L, 20 ug/L and 40 ug/L. Some analytes are prepared at higher concentrations.
- 17.8.5. Samples for Michigan projects frequently require calibration for 2-Methylnaphthalene. Recommended calibration levels for this compounds are 2 ug/L, 10 ug/L, 20 ug/L, 40 ug/L and 80 ug/L.





8771290

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 45 of 67

Table 1

## STL Primary Standard and Reporting Limits

Compound	CAS Number	Reporting Limits <sup>1</sup>				
		5 mL Water µg/L	25 mL <sup>3</sup> water µg/L	Low soil µg/kg	8260B/ 5035 Soil ug/kg	8260A 5030A Med Level Soil µg/kg
Dichlorodifluoromethane	75-71-8	10	2	10	250	1200
Chloromethane	74-87-3	10	2	10	250	1200
Bromomethane	74-83-9	10	2	10	250	1200
Vinyl chloride	75-01-4	10	2	10	250	1200
Chloroethane	75-00-3	10	2	10	250	1200
Trichlorofluoromethane	75-69-4	10	2	10	250	1200
Acrolein	107-02-8	100	20	100	5000	12000
Acetone	67-64-1	20	10	20	1000	2500
Trichlorotrifluoroethane	76-13-1	5	1	5	250	620
Iodomethane	74-88-4	5	1	5	250	620
Carbon disulfide	75-15-0	5	1	5	250	620
Methylene chloride	75-09-2	5	1	5	250	620
tert-Butyl alcohol	75-65-0	200	50	200	10,000	25000
1,1-Dichloroethene	75-35-4	5	1	5	250	620
1,1-Dichloroethane	75-34-3	5	1	5	250	620
trans-1,2-Dichloroethene	156-60-5	5.0	1.0	5.0	250	310
Acrylonitrile	107-13-1	100	20	100	5000	12000
Methyl tert-butyl ether (MTBE)	1634-04-4	20	5	20	1000	2500
Hexane	110-54-3	5	1	5	250	620
cis-1,2-Dichloroethene	156-59-2	5.0	1.0	5.0	250	310
1,2-Dichloroethene (Total)	540-59-0	5	1	5	250	620
Tetrahydrofuran	109-99-9	20	5	20	1000	2500
Chloroform	67-66-3	5	1	5	250	620
1,2-Dichloroethane	107-06-2	5	1	5	250	620
Dibromomethane	74-95-3	5	1	5	250	620
2-Butanone	78-93-3	20	5	20	1000	2500
1,4-Dioxane	123-91-1	500	200	500	25000	62000
1,1,1-Trichloroethane	71-55-6	5	1	5	250	620
Carbon tetrachloride	56-23-5	5	1	5	250	620
Bromodichloromethane	75-27-4	5	1	5	250	620
1,2-Dichloropropane	78-87-5	5	1	5	250	620
cis-1,3-Dichloropropene	10061-01-5	5	1	5	250	620

**Table 1**  
**STL Primary Standard and Reporting Limits**

Compound	CAS Number	Reporting Limits <sup>1</sup>				
		5 mL Water µg/L	25 mL <sup>3</sup> water µg/L	Low soil µg/kg	8260B/ 5035 Soil ug/kg	8260A 5030A Med Level Soil µg/kg
Trichloroethene	79-01-6	5	1	5	250	620
Dibromochloromethane	124-48-1	5	1	5	250	620
1,2-Dibromochloroethane	106-93-4	5	1	5	250	620
1,2,3-Trichloropropane	96-18-4	5	1	5	250	620
1,1,2-Trichloroethane	79-00-5	5	1	5	250	620
Benzene	71-43-2	5	1	5	250	620
Ethylmethacrylate	97-63-2	5	1	5	250	620
trans-1,3-Dichloropropene	10061-02-6	5	1	5	250	620
Bromoform	75-25-2	5	1	5	250	620
4-Methyl-2-pentanone	108-10-1	20	5	20	1000	2500
2-Hexanone	591-78-6	20	5	20	1000	2500
Tetrachloroethene	127-18-4	5	1	5	250	620
Toluene	108-88-3	5	1	5	250	620
1,1,2,2-Tetrachloroethane	79-34-5	5	1	5	250	620
2-Chloroethyl vinyl ether	110-75-8	N/A <sup>2</sup>	N/A	50	1000	6200
Vinyl acetate	108-05-4	10	2	10	500	1200
Chlorobenzene	108-90-7	5	1	5	250	620
Ethylbenzene	100-41-4	5	1	5	250	620
Styrene	100-42-5	5	1	5	250	620
t-1,4-Dichloro-2-butene	110-57-6	5	1	5	250	620
m and p Xylenes		5.0	0.5	2.5	125	310
o-xylene	95-47-6	5.0	0.5	2.5	125	310
Total xylenes	1330-20-7	5	1	5	250	620
1,3-Dichlorobenzene	541-73-1	5	1	5	250	620
1,4-Dichlorobenzene	106-46-7	5	1	5	250	620
1,2-Dichlorobenzene	95-50-1	5	1	5	250	620
2,2-Dichloropropane	590-20-7	5	1	5	250	
Bromochloromethane	74-97-5	5	1	5	250	
1,1-Dichloropropene	563-58-6	5	1	5	250	
Bromodichloromethane	75-27-4	5	1	5	250	
1,2-Dichloropropane	78-87-5	5	1	5	250	
1,3-Dichloropropane	142-28-9	5	1	5	250	

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 47 of 67

Table 1  
 STL Primary Standard and Reporting Limits

Compound	CAS Number	Reporting Limits <sup>1</sup>				
		5 mL Water µg/L	25 mL <sup>3</sup> water µg/L	Low soil µg/kg	8260B/ 5035 Soil ug/kg	8260A 5030A Med Level Soil µg/kg
Isopropylbenzene	98-82-8	5	1	5	250	
Bromobenzene	108-86-1	5	1	5	250	
n-Propylbenzene	103-65-1	5	1	5	250	
2-Chlorotoluene	95-49-8	5	1	5	250	
4-Chlorotoluene	106-43-4	5	1	5	250	
1,3,5-Trimethylbenzene	108-67-8	5	1	5	250	
tert-Butylbenzene	98-06-6	5	1	5	250	
1,2,4-Trimethylbenzene	95-63-6	5	1	5	250	
sec-butylbenzene	135-98-8	5	1	5	250	
4-Isopropyltoluene	99-87-6	5	1	5	250	
n-Butylbenzene	104-51-8	5	1	5	250	
1,2,4-Trichlorobenzene	120-82-1	5	1	5	250	
Napthalene	91-20-3	5	1	5	250	
Hexachlorobutadiene	87-68-3	5	1	5	250	
1,2,3-Trichlorobenzene	87-61-6	5	1	5	250	
Acetonitrile	75-05-8	100	20	100	5000	

<sup>1</sup> Reporting limits listed for soil/sediment are based on wet weight. The reporting limits calculated by the laboratory for soil/sediment, calculated on dry weight basis, will be higher.

<sup>2</sup> 2-Chloroethyl vinyl ether cannot be reliably recovered from acid preserved samples

<sup>3</sup> Optionally, 5 mL purge volume if adequate sensitivity is obtained.

Table 2

STL Primary Standard Calibration Levels, 5 mL purge<sup>1</sup>

Compound	Calibration Level ug/L				
	Level 1	Level 2	Level 3	Level 4	Level 5
1,2-Dichloroethane-d4 (Surrogate)	5	20	50	100	200
Toluene-d8 (Surrogate)	5	20	50	100	200
4-Bromofluorobenzene (Surrogate)	5	20	50	100	200
Dichlorodifluoromethane	5	20	50	100	200
Chloromethane	5	20	50	100	200
Bromomethane	5	20	50	100	200
Vinyl chloride	5	20	50	100	200
Chloroethane	5	20	50	100	200
Trichlorofluoromethane	5	20	50	100	200
Acrolein	50	200	500	1000	2000
Acetone	5	20	50	100	200
Trichlorotrifluoroethane	5	20	50	100	200
Iodomethane	5	20	50	100	200
Carbon disulfide	5	20	50	100	200
Methylene chloride	5	20	50	100	200
tert-Butyl alcohol	100	400	1,000	2,000	4,000
1,1-Dichloroethene	5	20	50	100	200
1,1-Dichloroethane	5	20	50	100	200
trans-1,2-Dichloroethene	5	20	50	100	200
Acrylonitrile	50	200	500	1,000	2,000
Methyl tert-butyl ether (MTBE)	5	20	50	100	200
Hexane	5	20	50	100	200
cis-1,2-Dichloroethene	5	20	50	100	200
Tetrahydrofuran	5	20	50	100	200
Chloroform	5	20	50	100	200
1,2-Dichloroethane	5	20	50	100	200
Dibromomethane	5	20	50	100	200
2-Butanone	5	20	50	100	200
1,4-Dioxane	250	1000	2,500	5,000	10,000
1,1,1-Trichloroethane	5	20	50	100	200
Carbon tetrachloride	5	20	50	100	200
Bromodichloromethane	5	20	50	100	200
1,2-Dichloropropane	5	20	50	100	200
cis-1,3-Dichloropropene	5	20	50	100	200
Trichloroethene	5	20	50	100	200
Dibromochloromethane	5	20	50	100	200

8771294

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 49 of 67

Table 2

STL Primary Standard Calibration Levels, 5 mL purge<sup>1</sup>

Compound	Calibration Level ug/L				
	Level 1	Level 2	Level 3	Level 4	Level 5
1,2-Dibromoethane	5	20	50	100	200
1,2,3-Trichloropropane	5	20	50	100	200
Acetonitrile	50	200	500	1000	2000
1,1,2-Trichloroethane	5	20	50	100	200
Benzene	5	20	50	100	200
Ethylmethacrylate	5	20	50	100	200
trans-1,3-Dichloropropene	5	20	50	100	200
Bromoform	5	20	50	100	200
4-Methyl-2-pentanone	5	20	50	100	200
2-Hexanone	5	20	50	100	200
Tetrachloroethene	5	20	50	100	200
Toluene	5	20	50	100	200
1,1,2,2-Tetrachloroethane	5	20	50	100	200
2-Chloroethyl vinyl ether	10	40	100	200	400
Vinyl acetate	5	20	50	100	200
Chlorobenzene	5	20	50	100	200
Ethylbenzene	5	20	50	100	200
Styrene	5	20	50	100	200
t-1,4-Dichloro-2-butene	5	20	50	100	200
m and p Xylenes	10	40	100	200	400
o-xylene	5	20	50	100	200
1,3-Dichlorobenzene	5	20	50	100	200
1,4-Dichlorobenzene	5	20	50	100	200
1,2-Dichlorobenzene	5	20	50	100	200
2,2-Dichloropropane	5	20	50	100	200
Bromochloromethane	5	20	50	100	200
1,1-Dichloropropene	5	20	50	100	200
Bromodichloromethane	5	20	50	100	200
1,2-Dichloropropane	5	20	50	100	200
1,3-Dichloropropane	5	20	50	100	200
Isopropylbenzene	5	20	50	100	200
Bromobenzene	5	20	50	100	200
n-Propylbenzene	5	20	50	100	200
2-Chlorotoluene	5	20	50	100	200
4-Chlorotoluene	5	20	50	100	200
1,3,5-Trimethylbenzene	5	20	50	100	200
tert-Butylbenzene	5	20	50	100	200

Table 2

STL Primary Standard Calibration Levels, 5 mL purge<sup>1</sup>

Compound	Calibration Level ug/L				
	Level 1	Level 2	Level 3	Level 4	Level 5
1,2,4-Trimethylbenzene	5	20	50	100	200
sec-butylbenzene	5	20	50	100	200
4-Isopropyltoluene	5	20	50	100	200
n-Butylbenzene	5	20	50	100	200
1,2,4-Trichlorobenzene	5	20	50	100	200
Napthalene	5	20	50	100	200
Hexachlorobutadiene	5	20	50	100	200
1,2,3-Trichlorobenzene	5	20	50	100	200

<sup>1</sup> Levels for 25 mL purge are 5 times lower in all cases

Table 2A

8771296

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 51 of 67

STL Primary Standard Calibration Levels, Low Level<sup>1</sup>

Compound	Calibration Level ug/L				
	Level 1	Level 2	Level 3	Level 4	Level 5
Dibromofluoromethane (Surrogate)	1	5	10	20	40
1,2-Dichloroethane-d4 (Surrogate)	1	5	10	20	40
Toluene-d8 (Surrogate)	1	5	10	20	40
Bromofluorobenzene (Surrogate)	1	5	10	20	40
Dichlorodifluoromethane	1	5	10	20	40
Chloromethane	1	5	10	20	40
Vinyl Chloride	1	5	10	20	40
Bromomethane	1	5	10	20	40
Chloroethane	1	5	10	20	40
Trichlorofluoromethane	1	5	10	20	40
Acrolein	10	50	100	200	400
Acetone	2	10	20	40	80
1,1-Dichloroethene	1	5	10	20	40
Trichlorotrifluoroethane	1	5	10	20	40
Iodomethane	1	5	10	20	40
Carbon Disulfide	1	5	10	20	40
Methylene Chloride	1	5	10	20	40
Acetonitrile	10	50	100	200	400
Acrylonitrile	10	50	100	200	400
Methyl tert-butyl ether	1	5	10	20	40
trans-1,2-Dichloroethene	1	5	10	20	40
Hexane	1	5	10	20	40
Vinyl acetate	1	5	10	20	40
1,1-Dichloroethane	1	5	10	20	40
tert-Butyl Alcohol	20	100	200	400	800
2-Butanone	2	10	20	40	80
cis-1,2-dichloroethene	1	5	10	20	40
2,2-Dichloropropane	1	5	10	20	40
Bromochloromethane	1	5	10	20	40
Chloroform	1	5	10	20	40
Tetrahydrofuran	1	5	10	20	40
1,1,1-Trichloroethane	1	5	10	20	40
1,1-Dichloropropene	1	5	10	20	40
Carbon Tetrachloride	1	5	10	20	40
1,2-Dichloroethane	1	5	10	20	40
Benzene	1	5	10	20	40
Trichloroethene	1	5	10	20	40
1,2-Dichloropropane	1	5	10	20	40
1,4-Dioxane	50	250	500	1000	2000
Dibromomethane	1	5	10	20	40

**Table 2A**  
**STL Primary Standard Calibration Levels, Low Level<sup>1</sup>**

Compound	Calibration Level ug/L				
	Level 1	Level 2	Level 3	Level 4	Level 5
Bromodichloromethane	1	5	10	20	40
2-Chloroethyl vinyl ether	2	10	20	40	80
cis-1,3-Dichloropropene	1	5	10	20	40
4-Methyl-2-pentanone	2	10	20	40	80
Toluene	1	5	10	20	40
trans-1,3-Dichloropropene	1	5	10	20	40
Ethyl Methacrylate	1	5	10	20	40
1,1,2-Trichloroethane	1	5	10	20	40
1,3-Dichloropropane	1	5	10	20	40
Tetrachloroethene	1	5	10	20	40
2-Hexanone	2	10	20	40	80
Dibromochloromethane	1	5	10	20	40
1,2-Dibromoethane	1	5	10	20	40
Chlorobenzene	1	5	10	20	40
1,1,1,2-Tetrachloroethane	1	5	10	20	40
Ethylbenzene	1	5	10	20	40
m + p-Xylene	2	10	20	40	80
Xylene-o	1	5	10	20	40
Styrene	1	5	10	20	40
Bromoform	1	5	10	20	40
Isopropylbenzene	1	5	10	20	40
1,1,2,2-Tetrachloroethane	1	5	10	20	40
1,4-Dichloro-2-butene	1	5	10	20	40
1,2,3-Trichloropropane	1	5	10	20	40
Bromobenzene	1	5	10	20	40
n-Propylbenzene	1	5	10	20	40
2-Chlorotoluene	1	5	10	20	40
1,3,5-Trimethylbenzene	1	5	10	20	40
4-Chlorotoluene	1	5	10	20	40
tert-Butylbenzene	1	5	10	20	40
1,2,4-Trimethylbenzene	1	5	10	20	40
sec-Butylbenzene	1	5	10	20	40
4-Isopropyltoluene	1	5	10	20	40
1,3-Dichlorobenzene	1	5	10	20	40
1,4-Dichlorobenzene	1	5	10	20	40
n-Butylbenzene	1	5	10	20	40
1,2-Dichlorobenzene	1	5	10	20	40
1,2-Dibromo-3-chloropropane	1	5	10	20	40



8771298

Determination of Volatile analytes by GC/MS  
Analysis of Volatile Organics  
Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
Revision No. 2.3  
Revision Date: 05/23/01  
Page 53 of 67

Table 2A

STL Primary Standard Calibration Levels, Low Level<sup>1</sup>

Compound	Calibration Level ug/L				
	Level 1	Level 2	Level 3	Level 4	Level 5
1,2,4-Trichlorobenzene	1	5	10	20	40
Hexachlorobutadiene	1	5	10	20	40
Naphthalene	1	5	10	20	40
1,2,3-Trichlorobenzene	1	5	10	20	40
Cyclohexane	1	5	10	20	40
Methyl Acetate	2	10	20	40	80
Methylcyclohexane	1	5	10	20	40
1,3,5-Trichlorobenzene	1	5	10	20	40

<sup>1</sup> 25 mL purge samples analyzed at 5 mL purge on more sensitive equipment.

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 54 of 67

Table 3

## STL Appendix IX Standard and Reporting Limits, 5 mL purge

Compound	CAS Number	Reporting Limits			
		5 mL Water µg/L	25mL <sup>2</sup> water µg/L	Low Soil µg/kg	Medium Soil µg/mL
Allyl Chloride	107-05-1	10	2	10	500
Acetonitrile	75-05-8	100	20	100	5000
Dichlorofluoromethane		10	2	10	500
Isopropyl ether	108-20-3	10	2	10	500
Chloroprene	126-99-8	5	2	5	250
n-Butanol	71-36-3	200	50	200	10,000
Propionitrile	107-12-0	20	4	20	1000
Methacrylonitrile	126-98-7	5	2	5	250
Isobutanol	78-83-1	200	50	200	10,000
Methyl methacrylate	80-62-6	5	2	5	250
1,1,1,2-Tetrachloroethane	630-20-6	5	1	5	250
1,2-Dibromo-3-chloropropane	96-12-8	10	2	10	500
Ethyl ether	60-29-7	10	2	10	500
Ethyl Acetate	141-78-6	20	4	20	1,000
2-Nitropropane	79-46-9	10	4	10	500
Cyclohexanone	108-94-1	N/A <sup>1</sup>	N/A <sup>1</sup>	N/A <sup>1</sup>	N/A <sup>1</sup>
Isopropylbenzene	98-82-8	5	1	5	250
2-Methylnaphthalene (Michigan only)	91-57-6	NA	5	NA	330

<sup>1</sup> Cyclohexanone decomposes to 1,1-dimethoxycyclohexane in methanolic solution. Reporting limits cannot be accurately determined.

<sup>2</sup> Optionally, 5 mL purge volume if adequate sensitivity is obtained.

8771300

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 55 of 67

Table 4

Recommended/STL Appendix IX Standard Calibration Levels, µg/L

Compound	Level 1	Level 2	Level 3	Level 4	Level 5
Allyl Chloride	5	20	50	100	200
Acetonitrile	50	200	500	1,000	2,000
Dichlorofluoromethane	5	20	50	100	200
Isopropyl ether	5	20	50	100	200
Chloroprene	5	20	50	100	200
n-Butanol	100	400	1,000	2,000	4,000
Propionitrile	10	40	100	200	400
Methacrylonitrile	5	20	50	100	200
Isobutanol	100	400	1,000	2,000	4,000
Methyl methacrylate	5	20	50	100	200
1,1,1,2-Tetrachloroethane	5	20	50	100	200
1,2-Dibromo-3-chloropropane	10	40	100	200	400
Ethyl ether	5	20	50	100	200
Ethyl Acetate	10	40	100	200	400
2-Nitropropane	10	40	100	200	400
Cyclohexanone	50	200	500	1,000	2,000
2-Methylnaphthalene (Michigan only)	2	10	20	40	80

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 56 of 67

Table 5

## Reportable Analytes for STL Standard Tests, Primary Standard

Compound	CAS Number	STL Standard List	TCLP	TCL	Appendix IX	UTS
Dichlorodifluoromethane	75-71-8				X	X
Chloromethane	74-87-3	X		X	X	X
Bromomethane	74-83-9	X		X	X	X
Vinyl chloride	75-01-4	X	X	X	X	X
Chloroethane	75-00-3	X		X	X	X
Trichlorofluoromethane	75-69-4				X	X
Acrolein	107-02-8				X	X
Acetone	67-64-1	X		X	X	X
Trichlorotrifluoroethane	76-13-1					X
Ethanol	64-17-5					
Iodomethane	74-88-4				X	X
Carbon disulfide	75-15-0	X		X	X	X
Methylene chloride	75-09-2	X		X	X	X
tert-Butyl alcohol	75-65-0					
1,1-Dichloroethene	75-35-4	X	X	X	X	X
1,1-Dichloroethane	75-34-3	X		X	X	X
trans-1,2-Dichloroethene	156-60-5	X		X	X	X
Total dichloroethene		X		X	X	X
Acrylonitrile	107-13-1				X	X
Methyl tert-butyl ether (MTBE)	1634-04-4					
Hexane	110-54-3					
cis-1,2-Dichloroethene	156-59-2	X		X		
Tetrahydrofuran	109-99-9					
Chloroform	67-66-3	X	X	X	X	X
1,2-Dichloroethane	107-06-2	X	X	X	X	X
Dibromomethane	74-95-3				X	X
2-Butanone	78-93-3	X	X	X	X	X
1,4-Dioxane	123-91-1				X	X
1,1,1-Trichloroethane	71-55-6	X		X	X	X
Carbon tetrachloride	56-23-5	X	X	X	X	X
Bromodichloromethane	75-27-4	X		X	X	X
1,2-Dichloropropane	78-87-5	X		X	X	X
cis-1,3-Dichloropropene	10061-01-5	X		X	X	X
Trichloroethene	79-01-6	X	X	X	X	X

8771302

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 57 of 67

Table 5

## Reportable Analytes for STL Standard Tests, Primary Standard

Compound	CAS Number	STL Standard List	TCLP	TCL	Appendix IX	UTS
Dibromochloromethane	124-48-1	X		X	X	X
1,2-Dibromoethane	106-93-4				X	X
1,2,3-Trichloropropane	96-18-4				X	X
1,1,2-Trichloroethane	79-00-5	X		X	X	X
Benzene	71-43-2	X	X	X	X	X
Ethylmethacrylate	97-63-2				X	X
trans-1,3-Dichloropropene	10061-02-6	X		X	X	X
Bromoform	75-25-2	X		X	X	X
4-Methyl-2-pentanone	108-10-1	X		X	X	X
2-Hexanone	591-78-6	X		X	X	
Tetrachloroethene	127-18-4	X	X	X	X	X
Toluene	108-88-3	X		X	X	X
1,1,2,2-Tetrachloroethane	79-34-5	X		X	X	X
2-Chloroethyl vinyl ether	110-75-8					
Vinyl acetate	108-05-4				X	
Chlorobenzene	108-90-7	X	X	X	X	X
Ethylbenzene	100-41-4	X		X	X	X
Styrene	100-42-5	X		X	X	
t-1,4-Dichloro-2-butene	110-57-6				X	
m and p Xylenes		X		X	X	X
o-xylene	95-47-6	X		X	X	X
Total xylenes	1330-20-7	X		X	X	X
1,3-Dichlorobenzene	541-73-1					
1,4-Dichlorobenzene	106-46-7					
1,2-Dichlorobenzene	95-50-1					

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 58 of 67

Table 6

## Reportable Analytes for STL Standard Tests, Appendix IX standard

Compound	Number	STL Standard List	TCLP	TCL	Appendix IX	UTS
Allyl Chloride	107-05-1				X	
Acetonitrile	75-05-8				X	X
Dichlorofluoromethane	75-43-4					
Isopropyl ether	108-20-3					
Chloroprene	126-99-8				X	
n-Butanol	71-36-3					
Propionitrile	107-12-0				X	
Methacrylonitrile	126-98-7				X	X
Isobutanol	78-83-1				X	X
Methyl methacrylate	80-62-6				X	X
1,1,1,2-Tetrachloroethane	630-20-6				X	X
1,2-Dibromo-3-chloropropane	96-12-8				X	X
Ethyl ether	60-29-7					X
Ethyl Acetate	141-78-6					X
2-Nitrotoluene	79-46-9					
Cyclohexanone	108-94-1					X
Isopropylbenzene	98-82-8					

**Table 7**  
**Internal Standards**

	Standard Concentration µg/mL	Quantitation ion (5 mL purge)	Quantitation ion (25 mL purge)
Fluorobenzene	50	96	96
Chlorobenzene-d5	50	117	119
1,4-Dichlorobenzene-d4	50	152	152

Notes:

- 1) 5 µL of the internal standard is added to the sample. This results in a concentration of each internal in the sample of 50µg/L for a 5 mL purge or 10 µg/L for a 25 mL purge
- 2) Except for medium level soils, the surrogate and internal standards may be combined in one solution.

**Table 8**  
**Surrogate Standards**

Surrogate Compounds	Standard Concentration µg/mL
1,2-Dichloroethane-d <sub>4</sub>	50
Dibromofluoromethane	50
Toluene-d <sub>8</sub>	50
4-Bromofluorobenzene	50

Notes:

- 1) 5 µL of the surrogate standard is added to the sample. This results in a concentration of each surrogate in the sample of 50µg/L for a 5 mL purge or 10 µg/L for a 25 mL purge
- 2) Except for medium level soils, the surrogate and internal standards may be combined in one solution
- 3) Recovery limits for surrogates are generated from historical data and are maintained by the QA department

Table 9

## Matrix Spike / LCS Control Compounds

Compound	Standard Concentration $\mu\text{g/mL}$
1,1-Dichloroethene	50
Trichloroethene	50
Toluene	50
Benzene	50
Chlorobenzene	50
n-Hexane (Ohio VAP only)	50

## Notes:

- 1) 5  $\mu\text{L}$  of the standard is added to the LCS or matrix spiked sample. This results in a concentration of each spike analyte in the sample of 50  $\mu\text{g/L}$  for a 5 mL purge or 10  $\mu\text{g/L}$  for a 25 mL purge.
- 2) Recovery and precision limits for LCS and MS/MSD are generated from historical data and are maintained by the QA department.
- 3) Full analyte spikes may also be used at the laboratories option or at client request.

Table 10

## BFB Key Ion Abundance Criteria

Mass	Ion Abundance Criteria
50	15% to 40% of Mass 95
75	30% to 60% of Mass 95
95	Base Peak, 100% Relative Abundance
96	5% to 9% of Mass 95
173	Less Than 2% of Mass 174
174	Greater Than 50% of Mass 95
175	5% to 9% of Mass 174
176	Greater Than 95%, But Less Than 101% of Mass 174
177	5% to 9% of Mass 176



**Table 11**  
**SPCC Compounds and Minimum Response Factors**

Compound	8260B, 8260A Min. RF
Chloromethane	0.100
1,1-Dichloroethane	0.100
Bromoform	>0.100
1,1,2,2-Tetrachloroethane	0.300
Chlorobenzene	0.300

**Table 12**  
**CCC compounds**

Compound	Max. %RSD from Initial Calibration	Max. %D for continuing calibration
Vinyl Chloride	<30.0	<20.0
1,1-Dichloroethene	<30.0	<20.0
Chloroform	<30.0	<20.0
1,2-Dichloropropane	<30.0	<20.0
Toluene	<30.0	<20.0
Ethylbenzene	<30.0	<20.0
n-Hexane (Ohio VAP only)	<30.0	<20.0

**Table 13**  
**Characteristic ions**

Compound	Primary*	Secondary	Tertiary
1,2-Dichloroethane-d <sub>4</sub> (Surrogate)	65	102	
Dichlorodifluoromethane	85	87	50, 101, 103
Chloromethane	50	52	49
Vinyl chloride	62	64	61
Bromomethane	94	96	79
Chloroethane	64	66	49
Trichlorofluoromethane	101	103	66

Determination of Volatile analytes by GC/MS  
 Analysis of Volatile Organics  
 Based on Method 8260B, 8260A, and 624

SOP No. CORP-MS-0002NC  
 Revision No. 2.3  
 Revision Date: 05/23/01  
 Page 62 of 67

Table 13  
 Characteristic ions

Compound	Primary*	Secondary	Tertiary
1,1-Dichloroethene	96	61	98
Acrolein	56	55	58
Iodomethane	142	127	141
Carbon disulfide	76	78	
Trichlorotrifluoroethane	151	101	153
Ethanol	45	46	
Acetone	43	58	
Methylene chloride	84	49	51, 86
tert-Butyl alcohol	59	74	
trans-1,2-Dichloroethene	96	61	98
Acrylonitrile	53	52	51
Methyl tert butyl ether	73		
Hexane	57	43	
1,1-Dichloroethane	63	65	83
cis-1,2-Dichloroethene	96	61	98
2-Butanone	43	72**	
Tetrahydrofuran	42	71	
Chloroform	83	85	47
1,2-Dichloroethane	62	64	98
Dibromomethane	93	174	95, 172, 176
1,4-Dioxane	88	58	
Vinyl acetate	43	86	
1,1,1-Trichloroethane	97	99	117
Carbon tetrachloride	117	119	121
Benzene	78	52	77
Trichloroethene	130	95	97, 132
1,2-Dichloropropane	63	65	41
Bromodichloromethane	83	85	129
2-Chloroethyl vinyl ether	63	65	106
cis-1,3-Dichloropropene	75	77	39
trans-1,3-Dichloropropene	75	77	39
1,1,2-Trichloroethane	97	83	85, 99
Chlorodibromomethane	129	127	131
Bromoform	173	171	175, 252
1,2,3-Trichloropropane	75	110	77, 112, 97
Toluene-d <sub>8</sub> (Surrogate)	98	70	100
4-Bromofluorobenzene (Surrogate)	95	174	176
Toluene	91	92	65
4-Methyl-2-pentanone	43	58	57, 100

Table 13  
 Characteristic ions

Compound	Primary*	Secondary	Tertiary
Tetrachloroethene	164	166	131
Ethyl methacrylate	69	41	99, 86, 114
2-Hexanone	43	58	57, 100
Chlorobenzene	112	114	77
Ethylbenzene	106	91	
Xylenes	106	91	
Styrene	104	103	78, 51, 77
Dichlorobenzene (all isomers)	146	148	111
trans 1,4-Dichloro-2-butene	53	75	89, 77, 124
1,1,2,2-Tetrachloroethane	83	85	131, 133
Allyl Chloride	76	41	78
Acetonitrile	40	41	
Dichlorofluoromethane	67	69	
Isopropyl ether	87	59	45
Chloroprene	53	88	90
n-Butanol	56	41	42
Propionitrile	54	52	55
Methacrylonitrile	41	67	52
Isobutanol	41	43	74
Methyl methacrylate	41	69	100
1,1,1,2-Tetrachloroethane	131	133	119
1,2-Dibromo-3-chloropropane	157	155	75
Ethyl ether	59	74	
Ethyl Acetate	43	88	61
2-Nitropropane	41	43	46
Cyclohexanone	55	42	98
Isopropylbenzene	105	120	

\* The primary ion should be used for quantitation unless interferences are present, in which case a secondary ion may be used.

\*\* m/z 43 may be used for quantitation of 2-Butanone, but m/z 72 must be present for positive identification.

**1. REQUIREMENTS FOR EPA 624**

- 1.1. Method 624 is required for demonstration of compliance with NPDES wastewater discharge permits. This method can be applied only to aqueous matrices. The standard analyte list and reporting limits are listed in Table B-1.
- 1.2. The tune period for this method is defined as 24 hours.
- 1.3. The initial calibration curve for this method requires at least three points.
- 1.4. Sample concentrations are calculated using the average RRF from the initial calibration curve.
- 1.5. Each target analyte is assigned to the closest eluting internal standard.
- 1.6. Initial demonstration of Proficiency
  - 1.6.1. The spiking level for the four replicate initial demonstration of proficiency is 20 µg/L. The acceptance criteria are listed in Table B-2
- 1.7. Initial calibration curve requirements:
  - 1.7.1. Target compounds must have RSD ≤ 35%.
  - 1.7.2. If this requirement can not be met, a regression curve must be constructed for the non-compliant compounds. There is no correlation coefficient requirement for the regression curve.
- 1.8. Continuing calibration verification requirements:
  - 1.8.1. The continuing calibration standard is from a different source than the initial calibration standard. The acceptance criteria are listed in Table B-2.
- 1.9. Matrix Spike and LCS requirements
  - 1.9.1. The matrix spike and LCS are spiked at 20 µg/L. A matrix spike duplicate is not necessary for this method. The recovery limits for matrix spike and LCS recovery are listed in Table C-2.
- 1.10. Method clarifications, modifications and additions

- 1.10.1. Section 5.2.2 of the source method describes the trap packing materials as Tenax GC, Methyl silicone, silica gel and coconut charcoal. STL routinely employs the Supelco VOCARB 3000, which consists of Carbopack B and Carboxen 1000 and 1001.
- 1.10.2. Section 5.3.2 of the source method describes a packed analytical column. STL routinely employs capillary columns when performing this method.
- 1.10.3. The source method provides a suggested list of compounds for internal and surrogate standards. STL uses the following two compounds which are not on the table: Chlorobenzene- $d_5$  (internal standard) and 1,2-Difluorobenzene- $d_4$  (surrogate).

Table A-1.

Method 624 Analytes and Reporting Limits

Analytes	µg/L
Benzene	5
Bromodichloromethane	5
Bromoform	5
Bromomethane	5
Carbon tetrachloride	5
Chlorobenzene	5
Chloroethane	5
2-Chloroethyl vinyl ether	5
Chloroform	5
Chloromethane	5
Dibromochloromethane	5
1,2-Dichlorobenzene	5
1,3-Dichlorobenzene	5
1,4-Dichlorobenzene	5
1,1-Dichloroethane	5
1,2-Dichloroethane	5
1,1-Dichloroethene	5
trans-1,2-Dichloroethene	5
1,2-Dichloropropane	5
cis-1,3-Dichloropropene	5
trans-1,3-Dichloropropene	5
Ethylbenzene	5
Methylene chloride	5
1,1,2,2-Tetrachloroethane	5
Tetrachloroethene	5
Toluene	5
1,1,1-Trichloroethane	5
1,1,2-Trichloroethane	5
Trichloroethene	5
Trichlorofluoromethane	5
Vinyl chloride	5

Table A-2.

## Method 624 QC Acceptance Criteria

Analytes	Daily QC check acceptance criteria (20µg/L spike)	Mean recovery, 4 replicate initial demonstration acceptance criteria (20µg/L spike)	Standard deviation, 4 replicate initial demonstration acceptance criteria (20µg/L spike)	Matrix spike and LCS acceptance criteria (% recovery)
Benzene	12.8-27.2	15.2-26.0	6.9	37-151
Bromodichloromethane	13.1-26.9	10.1-28.0	6.4	35-155
Bromoform	14.2-25.8	11.4-31.1	5.4	45-169
Bromomethane	2.8-37.2	D-41 2	17.9	D-242
Carbon tetrachloride	14.6-25.4	17.2-23.5	5.2	70-140
Chlorobenzene	13.2-26.8	16.4-27.4	6.3	37-160
Chloroethane	7.6-32.4	8.4-40.4	11.4	14-230
2-Chloroethyl vinyl ether	D-44.8	D-50.4	25.9	D-305
Chloroform	13.5-26.5	13.7-24.2	6.1	51-138
Chloromethane	D-40.8	D-45.9	19.8	D-273
Dibromochloromethane	13.5-26.5	13.8-26.6	6.1	53-149
1,2-Dichlorobenzene	12.6-27.4	11.8-34.7	7.1	18-190
1,3-Dichlorobenzene	14.6-25.4	17.0-28.8	5.5	59-156
1,4-Dichlorobenzene	12.6-27.4	11.8-34.7	7.1	18-190
1,1-Dichloroethane	14.5-25.5	14.2-28.5	5.1	59-155
1,2-Dichloroethane	13.6-26.4	14.3-27.4	6.0	49-155
1,1-Dichloroethene	10.1-29.9	3.7-42.3	9.1	D-234
trans-1,2-Dichloroethene	13.9-26.1	13.6-28.5	5.7	54-156
1,2-Dichloropropane	6.8-33.2	3.8-36.2	13.8	D-210
cis-1,3-Dichloropropene	4.8-35.2	1.0-39.0	15.8	D-227
trans-1,3-Dichloropropene	10.0-30.0	7.6-32.4	10.4	17-183
Ethylbenzene	11.8-28.2	17.4-26.7	7.5	37-162
Methylene chloride	12.1-27.9	D-41.0	7.4	D-221
1,1,2,2-Tetrachloroethane	12.1-27.9	13.5-27.2	7.4	46-157
Tetrachloroethene	14.7-25.3	17.0-26.6	5.0	64-148
Toluene	14.9-25.1	16.6-26.7	4.8	47-150
1,1,1-Trichloroethane	15.0-25.0	13.7-30.1	4.6	52-162
1,1,2-Trichloroethane	14.2-25.8	14.3-27.1	5.5	52-150
Trichloroethene	13.3-26.7	18.6-27.6	6.6	71-157
Trichlorofluoromethane	9.6-30.4	8.9-31.5	10.0	17-181
Vinyl chloride	0.8-39.2	D-43.5	20.0	D-251

Controlled Copy

Copy No: \_\_\_\_\_

Implementation Date: 3/2/04

SOP No: CORP-MS-0001NC

Revision No: 2.10Revision Date: 02/24/04Page 1 of 50

## STL NORTH CANTON STANDARD OPERATING PROCEDURE

## TITLE: GC/MS ANALYSIS BASED ON METHODS 8270C

(SUPERSEDES: Revision 2.9, Dated 06/18/03)

Prepared by: Tom Hula 2/26/04

Date

Reviewed by: Mark Wlman 2/26/04  
Technology Specialist DateApproved by: Beth Lambert 2/27/04  
Quality Assurance Manager DateApproved by: Dorothy Dun 2/27/04  
Environmental Health and Safety DateApproved by: Paul M. Ehr 2/27/04  
Laboratory Director DateApproved by: Mark G. Bane 3/1/04  
Technical Director Date

## Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions.

THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY. ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.



**TABLE OF CONTENTS**

<b>1. SCOPE AND APPLICATION .....</b>	<b>4</b>
<b>2 SUMMARY OF METHOD.....</b>	<b>4</b>
<b>3 DEFINITIONS.....</b>	<b>5</b>
<b>4 INTERFERENCES .....</b>	<b>6</b>
<b>5 SAFETY PRECAUTIONS .....</b>	<b>6</b>
<b>6 EQUIPMENT AND SUPPLIES .....</b>	<b>8</b>
<b>7 REAGENTS AND STANDARDS.....</b>	<b>9</b>
<b>8 SAMPLE PRESERVATION AND STORAGE.....</b>	<b>9</b>
<b>9 QUALITY CONTROL .....</b>	<b>10</b>
<b>10 CALIBRATION AND STANDARDIZATION .....</b>	<b>14</b>
<b>11 PROCEDURE .....</b>	<b>18</b>
<b>12 DATA ANALYSIS AND CALCULATIONS.....</b>	<b>20</b>
<b>13 METHOD PERFORMANCE .....</b>	<b>27</b>
<b>14 POLLUTION PREVENTION.....</b>	<b>28</b>
<b>15 WASTE MANAGEMENT.....</b>	<b>28</b>
<b>16 REFERENCES .....</b>	<b>28</b>
<b>17 MISCELLANEOUS .....</b>	<b>29</b>

**LIST OF TABLES**

<b>TABLE 1</b>	STL North Canton Primary Standard and Standard Reporting Limits
<b>TABLE 2</b>	STL North Canton Appendix IX Standard and Standard Reporting Limits
<b>TABLE 2A</b>	STL North Canton Michigan Program
<b>TABLE 3</b>	Reportable Analytes for STL North Canton Standard Tests, Primary Standard
<b>TABLE 4</b>	Reportable analytes for STL North Canton Standard Tests, Appendix IX Standard
<b>TABLE 5</b>	Suggested Instrument Conditions
<b>TABLE 6</b>	DFTPP Ion Abundance Criteria
<b>TABLE 7</b>	Characteristic Ions, Primary Standard
<b>TABLE 8</b>	Characteristic Ions, Appendix IX Standard
<b>TABLE 9</b>	8270C LCS Control Compounds
<b>TABLE 9A</b>	8270 All-Analyte Spike Mix
<b>TABLE 10</b>	TCLP LCS Compounds
<b>TABLE 11</b>	8270C Surrogate Compounds
<b>TABLE 12</b>	Calibration Ranges
<b>TABLE 13</b>	Calibration Ranges, Appendix IX

## 1. SCOPE AND APPLICATION

- 1.1 This method is based upon SW846 8270C, and is applicable to the determination of the concentration of semivolatile organic compounds in extracts prepared from solid and aqueous matrices. Direct injection of a sample may be used in limited applications. Refer to Tables 1, 2, 3 and 4 for the list of compounds applicable for this method. Note that the compounds are listed in approximate retention time order. Additional compounds may be amenable to this method. If non-standard analytes are required, they must be validated by the procedures described in section 13 before sample analysis.
- 1.2 The following compounds may require special treatment when being determined by this method:
- Benzidine can be subject to oxidative losses during solvent concentration and exhibits poor chromatography. Neutral extraction should be performed if this compound is expected.
  - Hexachlorocyclopentadiene is subject to thermal decomposition in the inlet of the gas chromatograph, chemical reaction in acetone solution, and photochemical decomposition.
  - N-Nitrosodiphenylamine decomposes in the gas chromatographic inlet and cannot be distinguished from diphenylamine.
  - Pentachlorophenol, 2,4-dinitrophenol, 4-nitrophenol, 4,6-dinitro-2-methylphenol, 4-chloro-3-methylphenol, benzoic acid, 2-nitroaniline, 3-nitroaniline, 4-chloroaniline, and benzyl alcohol are subject to erratic chromatographic behavior, especially if the GC system is contaminated with high boiling material.
  - Hexachlorophene is not amenable to analysis by this method.
  - 3-Methylphenol cannot be separated from 4-methylphenol by the conditions specified in this method.
- 1.3 The standard reporting limit of this method for determining an individual compound is approximately 0.33 mg/kg (wet weight) for soil/sediment samples, 1 - 200 mg/kg for wastes (dependent on matrix and method of preparation), and 10 µg/L for groundwater samples. Some compounds have higher reporting limits. Refer to Tables 1 and 2 for specific SRLs. Reporting limits will be proportionately higher for sample extracts that require dilution.
- 1.4 The associated LIMS code is QL (8270C).

## 2 SUMMARY OF METHOD

- 2.1 Aqueous samples are extracted with methylene chloride using a separatory funnel, and/or a

continuous extractor. Solid samples are extracted with methylene chloride / acetone using sonication, soxhlet, accelerated soxhlet or pressurized fluid extraction. The extract is dried, concentrated to a final volume of 2 mL for waters and soils, and analyzed by GC/MS. Extraction procedures are detailed in SOP# CORP-OP-0001NC. Qualitative identification of the parameters in the extract is performed using the retention time and the relative abundance of characteristic ions. Quantitative analysis is performed using the internal standard technique with a single characteristic ion.

### 3 DEFINITIONS

- 3.1 CCC (Calibration Check Compounds) - A subset of target compounds used to evaluate the calibration stability of the GC/MS system. A maximum percent deviation of the CCC's is specified for calibration acceptance.
- 3.2 SPCC (System Performance Check Compounds) - Target compounds designated to monitor chromatographic performance, sensitivity, and compound instability or degradation on active sites. Minimum response factors are specified for acceptable performance.
- 3.3 Batch - The batch is a set of up to 20 samples of the same matrix processed using the same procedures and reagents within the same time period. The Quality Control batch must contain a matrix spike / spike duplicate (MS/MSD), a Laboratory Control Sample (LCS), and a method blank. Batches are defined at the sample preparation stage. Batches should be kept together through the whole analytical process to the extent possible, but it is not mandatory to analyze prepared extracts on the same instrument or in the same sequence. Refer to the STL North Canton QC Program document (QA-003) for further details of the batch definition.
- 3.4 Method Blank - An analytical control consisting of all reagents, internal standards and surrogate standards, that is carried through the entire analytical procedure. The method blank is used to define the level of laboratory background and reagent contamination.
- 3.5 LCS (Laboratory Control Sample) - A blank spiked with the parameters of interest that is carried through the entire analytical procedure. Analysis of this sample with acceptable recoveries of the spiked materials demonstrates that the laboratory techniques for this method are acceptable.
- 3.6 MS (Matrix Spike)- aliquot of a matrix (water or soil) fortified (spiked) with known quantities of specific compounds and subjected to the entire analytical procedure in order to indicate the appropriateness of the method for the matrix by measuring recovery.
- 3.7 MSD (Matrix Spike Duplicate)- a second aliquot of the same sample as the matrix spike (above) that is spiked in order to determine the precision of the method.

#### 4 INTERFERENCES

- 4.1 Method interferences may be caused by contaminants in solvents, reagents, glassware, and other processing apparatus that lead to discrete artifacts. All of these materials must be routinely demonstrated to be free from interferences under conditions of the analysis by running laboratory method blanks as described in the Quality Control section. Raw GC/MS data from all blanks, samples, and spikes must be evaluated for interferences. If an interference is detected it is necessary to determine if the source of interference is in the preparation and/or cleanup of the samples; then take corrective action to eliminate the problem.
- 4.2 The use of high purity reagents, solvents, and gases helps to minimize interference problems.
- 4.3 Matrix interferences may be caused by contaminants that are coextracted from the sample. The extent of matrix interferences will vary considerably from source to source, depending upon the nature of the sample.
- 4.4 Contamination by carryover can occur whenever high-level and low-level samples are sequentially analyzed. To reduce carryover, the sample syringe must be rinsed with solvent between samples. Whenever an unusually concentrated sample is encountered, it should be followed by the analysis of solvent to check for cross contamination.
- 4.5 Phthalate contamination is commonly observed in this analysis and its occurrence should be carefully evaluated as an indicator of a contamination problem in the sample preparation step of the analysis.

#### 5 SAFETY PRECAUTIONS

- 5.1 Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual and this document.
- 5.2 Eye protection that protects against splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents and reagents are being handled. Disposable gloves that have become contaminated will be removed and discarded; other gloves will be cleaned immediately.
- 5.3 Chemicals that have been classified as carcinogens, or potential carcinogens, under OSHA include: Benzo(a)anthracene, benzidine, 3,3'-dichlorobenzidine, benzo(a)pyrene, dibenzo(a,h)anthracene, and n-nitrosodimethylamine. Primary standards should be purchased in solution. If neat materials must be obtained, they shall be handled in a hood.

- 5.4 The following is a list of the materials used in this method, which have a serious or significant hazard rating. **NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table.** A complete list of materials used in the method can be found in the reagents and materials section. Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

Material (1)	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
Methylene Chloride	Carcinogen  Irritant	25 ppm-TWA  125 ppm-STEL	Causes irritation to respiratory tract. Has a strong narcotic effect with symptoms of mental confusion, light-headedness, fatigue, nausea, vomiting and headache. Causes irritation, redness and pain to the skin and eyes. Prolonged contact can cause burns. Liquid degrades the skin. May be absorbed through skin
Sulfuric Acid	Corrosive  Oxidizer  Dehydrator  Poison  Carcinogen	1 Mg/M3-TWA	Inhalation produces damaging effects on the mucous membranes and upper respiratory tract. Symptoms may include irritation of the nose and throat, and labored breathing. Symptoms of redness, pain, and severe burn can occur. Contact can cause blurred vision, redness, pain and severe tissue burns. Can cause blindness
1 – Always add acid to water to prevent violent reactions.			
2 – Exposure limit refers to the OSHA regulatory exposure limit.			

- 5.5 Exposure to chemicals must be maintained as low as reasonably achievable; therefore, unless they are known to be non-hazardous, all samples should be opened, transferred, and prepared in a fume hood, or under other means of mechanical ventilation. Solvent and waste containers should be kept closed unless transfers are being made.
- 5.6 The preparation of standards and reagents will be conducted in a fume hood with the sash closed as far as the operation will permit.
- 5.7 It is recommended that neat standards be purchased only as a last resort. The preparation of standards from neat materials and reagents {as well as glassware cleaning procedures that involved solvents such as methylene chloride} should be conducted in a fume hood with the sash closed as far as the operations will permit.
- 5.8 Standards in solution may be diluted in the open laboratory when syringes and the like are utilized.
- 5.9 All work must be stopped in the event of a known or potential compromise to the health and safety of a STL North Canton associate. The situation must be reported immediately to a laboratory supervisor.

## 6 EQUIPMENT AND SUPPLIES

- 6.1 Gas Chromatograph/Mass Spectrometer System: An analytical system complete with a temperature-programmable gas chromatograph suitable for split/splitless injection and all required accessories, including syringes, analytical columns, and gases. The capillary column should be directly coupled to the source.
- 6.2 Column: 20m x 0.18mm ID, 0.18 $\mu$ m film thickness silicon-coated fused-silica capillary column (J & W Scientific DB-5.625 or equivalent). Alternate columns are acceptable if they provide acceptable performance.
- 6.3 Mass Spectrometer: Capable of scanning from 35 to 500 AMU every one second or less, using 70 volts (nominal) electron energy in the electron impact ionization mode. The mass spectrometer must be capable of producing a mass spectrum for decafluorotriphenyl-phosphine (DFTPP) which meets all of the criteria in Table 6 when the GC/MS tuning standard is injected through the GC.
- 6.4 GC/MS Interface: Any GC-to-MS interface that gives acceptable calibration points and achieves acceptable tuning performance criteria may be used.
- 6.5 Data System: A computer system must be interfaced to the mass spectrometer. The

system must allow the continuous acquisition and storage on machine-readable media of all mass spectra obtained throughout the duration of the chromatographic program. The computer must have software that can search any GC/MS data file for ions of a specific mass and that can plot such ion abundances versus time or scan number. This type of plot is defined as the Extracted Ion Current Profile (EICP). Software must also be available that allows integrating the abundances in any EICP between specified time or scan-number limits. The most recent version of the EPA/NIH Mass Spectral Library is recommended.

6.6 Syringe: 5  $\mu$ L Hamilton Laboratory grade syringes or equivalent.

6.7 Carrier gas: Ultra high purity helium.

## 7 REAGENTS AND STANDARDS

7.1 A minimum five point calibration curve is prepared. If a quadratic regression is used, six points must be analyzed for the calibration curve. The low point should be at or below the reporting limit. Refer to Tables 12 and 13 for typical calibration levels for all analytes. Other calibration levels may be used, depending on instrument capability, but the low standard must support the reporting limit and the high standard defines the range of the calibration.

7.2 An Internal Standard solution is prepared by diluting a purchased standard. Compounds in the I.S. Mix are: acenaphthene-d10, chrysene-d12, 1,4-dichlorobenzene-d4, naphthalene-d8, perylene-d12, and phenanthrene-d10.

7.3 Surrogate Standard Spiking Solution: Prepare as indicated in the preparative methods. See appropriate preparation SOP. Surrogate compounds and levels are listed in Table 11.

7.4 GC/MS Tuning Standard: A methylene chloride solution containing decafluorotriphenylphosphine (DFTPP) is prepared. Pentachlorophenol, benzidine, and DDT, should also be included in the Tuning Standard. All components are at 25  $\mu$ g/mL.

7.5 The standards listed in 7.1 to 7.4 should be refrigerated at  $\leq 6^{\circ}\text{C}$  when not in use. Refrigeration at  $-10^{\circ}\text{C}$  to  $-20^{\circ}\text{C}$  may be used if it can be demonstrated that analytes do not fall out of solution at this temperature. The standards must be replaced at least once a year.

## 8 SAMPLE PRESERVATION AND STORAGE

8.1 Sample extracts are stored at  $4 \pm 2^{\circ}\text{C}$ . Samples and extracts should be stored in suitable glass containers with Teflon lined caps. (Extracts will normally be stored for 30 days after invoicing.)

8.3 Water samples are extracted within seven days of sampling and the extracts are analyzed



within forty days of extraction. Solids, sludges, and organic liquids are extracted within fourteen days of sampling and the extracts are analyzed within forty days of extraction.

## 9 QUALITY CONTROL

### 9.1 Initial Demonstration of Capability

- 9.1.1 For the standard analyte list, the initial demonstration and method detection limit (MDL) studies described in section 13 must be acceptable before analysis of samples may begin.
- 9.1.2 For non-standard analytes an MDL study should be performed and calibration curve generated before analyzing any samples, unless lesser requirements are previously agreed to with the client. In any event, the minimum initial demonstration required is analysis of an extracted standard at the reporting limit and a single point calibration.

### 9.2 Control Limits

In-house historical control limits must be determined for surrogates, matrix spikes, and laboratory control samples (LCS). These limits must be determined at least annually. The recovery limits are mean recovery  $\pm 3$  standard deviations for surrogates, MS and LCS. Precision limits for matrix spikes / matrix spike duplicates are mean relative percent difference  $\pm 3$  standard deviations.

- 9.2.1 These limits do not apply to dilutions (except for tests without a separate extraction), but surrogate and matrix spike recoveries will be reported.
- 9.2.2 All surrogate, LCS, and MS recoveries (except for dilutions) must be entered into QuantIMS (when available) or other database so that accurate historical control limits can be generated. For tests without a separate extraction, surrogates and matrix spikes will be reported for all dilutions.
- 9.2.3 Refer to the QC program document (QA-003) for further details of control limits.

### 9.3 Method Blank

A method blank is prepared and analyzed with each batch of samples. The method blank consists of reagent water for aqueous samples, and sodium sulfate for soil samples (Refer to SOP No. CORP-OP-0001NC for details). Surrogates are added and the method blank is carried through the entire analytical procedure. The method blank must not contain any analyte of interest at or above the reporting limit (except common laboratory contaminants, see below) or at or above 5% of the measured concentration of that analyte in the associated samples, whichever is higher.

- If the analyte is a common laboratory contaminant (phthalate esters), the data may be reported with qualifiers if the concentration of the analyte is less than five times the RL. Such action must be taken in consultation with the client.
- Reanalysis of any samples with reportable concentrations of analytes found in the method blank is required unless other actions are agreed with the client.
- If there is no target analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. Such action should be taken in consultation with the client.

9.3.1 The method blank must have acceptable surrogate recoveries. If surrogate recoveries are not acceptable, the data must be evaluated to determine if the method blank has served the purpose of demonstrating that the analysis is free of contamination. If surrogate recoveries are low and there are reportable analytes in the associated samples, re-extraction of the blank and affected samples will normally be required. Consultation with the client should take place.

9.3.2 If reanalysis of the batch is not possible due to limited sample volume or other constraints, the method blank is reported, all associated samples are flagged with a "B", and appropriate comments may be made in a narrative to provide further documentation.

9.3.3 Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective actions.

#### 9.4 Instrument Blank

9.4.1 Instruments must be evaluated for contamination during each 12 hour analytical run. This may be accomplished by analysis of a method blank. If a method blank is not available, an instrument blank must be analyzed. An instrument blank consists of methylene chloride with the internal standards added. It is evaluated in the same way as the method blank.

#### 9.5 Laboratory Control Sample (LCS)

9.5.1 A laboratory control sample (LCS) is prepared and analyzed with every batch of samples. All control analytes must be within established control limits. The LCS is spiked with the compounds listed in Tables 9 and 10 unless specified by a client or agency.

9.5.2 If any control analyte in the LCS is outside the laboratory established historical control limits, corrective action must occur. Corrective action may include re-extraction and

reanalysis of the batch.

- If the batch is not re-extracted and reanalyzed, the reasons for accepting the batch must be clearly presented in the project records and the report. (An example of acceptable reasons for not reanalyzing might be that the matrix spike and matrix spike duplicate are acceptable, and sample surrogate recoveries are good, demonstrating that the problem was confined to the LCS).
- If re-extraction and reanalysis of the batch is not possible due to limited sample volume or other constraints, the LCS is reported, all associated samples are flagged, and appropriate comments are made in a narrative to provide further documentation.

9.5.3 Ongoing monitoring of the LCS provides evidence that the laboratory is performing the method within accepted QC guidelines for accuracy and precision.

9.5.4 Additionally, if an all-analyte check sample is used, all non-controlling compounds must attain a recovery of 5% or greater if the compound is on the client's list.

#### 9.6 Matrix Spike/Matrix Spike Duplicate (MS/MSD)

A matrix spike/matrix spike duplicate (MS/MSD) is prepared and analyzed with every batch of samples. The MS/MSD is spiked with the same subset of analytes as the LCS (See Tables 9 and 10). Compare the percent recovery and relative percent difference (RPD) to that in the laboratory specific historically generated limits.

- If any individual recovery or RPD falls outside the acceptable range, corrective action must occur. The initial corrective action will be to check the recovery of that analyte in the Laboratory Control Sample (LCS). Generally, if the recovery of the analyte in the LCS is within limits, then the laboratory operation is in control and analysis may proceed. The reasons for accepting the batch must be documented.
- If the recovery for any component is outside QC limits for both the Matrix spike / spike duplicate and the LCS, the laboratory is out of control and corrective action must be taken. Corrective action will normally include re-preparation and reanalysis of the batch.
- If a MS/MSD is not possible due to limited sample, then a LCS duplicate should be analyzed. RPD of the LCS and LCSD are compared to the matrix spike limits.
- The matrix spike / duplicate must be analyzed at the same dilution as the unspiked sample, even if the matrix spike compounds will be diluted out.

## 9.7 Surrogates

9.7.1 Every sample, blank, and QC sample is spiked with surrogate standards. Surrogate spike recoveries must be evaluated by determining whether the concentration (measured as percent recovery) falls within the required recovery limits. The compounds routinely included in the surrogate spiking solution, along with recommended standard concentrations, are listed in Table 11.

9.7.2 If any surrogates are outside limits the following corrective actions must take place (except for dilutions):

- Check all calculations for error.
- Ensure that instrument performance is acceptable.
- Recalculate the data and/or reanalyze the extract if either of the above checks reveal a problem.

It is only necessary to reprepare / reanalyze a sample once to demonstrate that poor surrogate recovery is due to matrix effect, unless the analyst believes that the repeated out of control results are not due to matrix effect.

**Note:** If all associated QC meets criteria (blank, LCS/LCSD), up to one surrogate per fraction may be outside of acceptance criteria, as long as the recovery is greater than 10%.

**Note:** For Ohio VAP samples, all surrogates must be within acceptance criteria.

9.7.3 If the sample with surrogate recoveries outside the recovery limits was a sample used for an MS/MSD and the surrogate recoveries in the MS/MSD are also outside of the control limits, then the sample, the MS, and the MSD do not require reanalysis as this phenomenon would indicate a possible matrix problem.

9.7.4 If the sample is reanalyzed and the surrogate recoveries in the reanalysis are acceptable, then the problem was within the analyst's control and only the reanalyzed data should be reported. (Unless the reanalysis was outside holding times, in which case reporting both sets of results may be appropriate.)

9.7.5 If the reanalysis does confirm the original results, the original analysis is reported and the data flagged as estimated due to matrix effect.

## 9.8 Nonconformance and Corrective Action

- 9.8.1 Any deviations from QC procedures must be documented as a nonconformance, with applicable cause and corrective action approved by the facility QA Manager.

# 10 CALIBRATION AND STANDARDIZATION

## 10.1 Summary

- 10.1.1 The instrument is tuned for DFTPP, calibrated initially with a minimum five-point calibration curve, and verified each 12-hour shift with one or more continuing calibration standard(s). Recommended instrument conditions are listed in Table 5.

- 10.2 All standards and extracts are allowed to warm to room temperature before injecting.

## 10.3 Instrument Tuning

At the beginning of every twelve hour shift when analyses are to be performed, the GC/MS system must be checked to see if acceptable performance criteria (Table 6) is achieved for DFTPP (decafluorotriphenylphosphine).

- 10.3.1 Inject the GC/MS tuning standard (Section 7.4) into the GC/MS system. Obtain a background-corrected mass spectra of DFTPP and confirm that all the key  $m/z$  criteria in Table 6 are achieved. If all the criteria are not achieved, the analyst must retune the mass spectrometer and repeat the test until all criteria are achieved. The performance criteria must be achieved before any samples, blanks, or standards are analyzed.
- 10.3.2 The GC/MS tuning standard should also be used to evaluate the inertness of the chromatographic system. Benzidine and pentachlorophenol should not exhibit excessive tailing. If DDT is an analyte of interest, it must be included in the tuning standard, and its breakdown must be  $< 20\%$ . Refer to section 12 for the appropriate calculations.

## 10.4 Initial Calibration

- 10.4.1 Internal Standard Calibration Procedure: Internal standards are listed in Table 7. Use the base peak  $m/z$  as the primary  $m/z$  for quantitation of the standards. If interferences are noted, use one of the next two most intense masses for quantitation.
- 10.4.2 Compounds should be assigned to the IS with the closest retention time.
- 10.4.3 Prepare calibration standards at a minimum of five concentration levels for each parameter of interest. Six standards must be used for a quadratic least squares calibration. Quadratic fit may NOT be used for samples analyzed under South Carolina

Certification. It may also be useful to analyze six calibration levels and use the lower five for most analytes and the upper five for analytes that have poor response. Add the internal standard mixture to result in 2 ng on column. (For example, 5 uL of 80ppm IS mix is added to 100 uL of extract. This results in 4 ng, but only 0.5uL is injected, resulting in a final on column amount of 2 ng.) The concentration ranges of all analytes are listed in tables 12 and 13.

10.4.4 Analyze each calibration standard and tabulate the area of the primary characteristic m/z against concentration for each compound and internal standard. Calculate response factors (RF), average response factors, and the percent RSD of the response factors for each compound using the equations in section 12 and verify that the CCC and SPCC criteria in section 10.4.5 and 10.4.6 are met. **No sample analysis may be performed unless these criteria are met.**

10.4.5 System Performance Check Compounds (SPCCs): The minimum average RF for semivolatile SPCCs is 0.050. If the minimum response factors are not met, the system must be evaluated and corrective action must be taken before sample analysis begins. Some possible problems are standard mixture degradation, injection port inlet contamination, contamination at the front end of the analytical column, and active sites in the column or chromatographic system. This check must be met before analysis begins. SPCC Compounds:

N-nitroso-di-n-propylamine  
Hexachlorocyclopentadiene  
2,4-Dinitrophenol  
4-Nitrophenol

10.4.6 Calibration Check Compounds (CCCs): The %RSD of the response factors for each CCC in the initial calibration must be less than 30% for the initial calibration to be considered valid. This criterion must be met before sample analysis begins. Problems similar to those listed under SPCCs could affect this criterion.

10.4.6.1 If none of the CCCs are required analytes, project specific calibration specifications must be agreed with the client.

10.4.6.2 CCC Compounds:

Phenol  
Acenaphthene  
1,4-Dichlorobenzene  
N-nitrosodiphenylamine  
2-Nitrophenol

Pentachlorophenol  
2,4-Dichlorophenol  
Fluoranthene  
Hexachlorobutadiene  
Di-n-octylphthalate  
4-Chloro-3-methylphenol  
Benzo(a)pyrene  
2,4,6-Trichlorophenol

- 10.4.7 If the software in use is capable of routinely reporting curve coefficients for data validation purposes, and the necessary calibration reports can be generated, then the analyst should evaluate analytes with %RSD > 15% for calibration on a curve. If it appears that substantially better accuracy would be obtained using quantitation from a curve then the appropriate curve should be used for quantitation.
- 10.4.7.1 If an analyte in the initial calibration is > 15%, then calibration on a curve must be used. Linear or quadratic curve fits may be used. Linear curve fits only may be used for South Carolina Certification. The analyst should consider instrument maintenance to improve the linearity of response. Use of  $1/\text{Concentration}^2$  weighting is recommended to improve the accuracy of quantitation at the low end of the curve. If Relative Standard Error (RSE) is used to evaluate the curve it must be better than 15%. If the % RSD is >15%, the analyst may drop the low or high points in the ICAL, as long as a minimum of 5 points are maintained and the quantitation range is adjusted accordingly. If the % RSD is still >15%, a quadratic or linear curve may be used. The correlation coefficient (r) must be  $\geq 0.990$ . If the correlation coefficient is < 0.990, then any hits for these compounds must be flagged as estimated. If a curve is not linear for any compound that is found in a samples, the result must be flagged as estimated. Linear is defined as <15% RSD or a correlation coefficient of 0.990.
- 10.4.7.2 Note: Several components do not respond well by this method (poor linearity). These compounds are famphur, benzenethiol, kepone, and 2,4-toluenediamine. If these compounds are requested by a client and hits are found, alternate standards or methods will be needed for more accurate quantitation. Sensitivity as demonstrated by the low standard is sufficient to substantiate a non-detect.
- 10.4.8 If time remains in the 12 hour period initiated by the DFTPP injection before the initial calibration, samples may be analyzed. Otherwise, proceed to continuing calibration.

**10.4.9 Quantitation is performed using the calibration curve or average response factor from the initial curve, not the continuing calibration.**

**10.5 Continuing Calibration**

10.5.1 At the start of each 12-hour period, the GC/MS tuning standard must be analyzed. The injection of DFTPP must result in a mass spectrum for DFTPP which meets the criteria given in Table 6.

10.5.2 Following a successful DFTPP analysis the continuing calibration standard(s) are analyzed. The standards must contain all semivolatile analytes, including all required surrogates. A mid level calibration standard is used for the continuing calibration.

10.5.3 The following criteria must be met for the continuing calibration to be acceptable:

- The SPCC compounds must have a response factor of  $\geq 0.05$ .
- The percent difference or drift of the CCC compounds from the initial calibration must be  $\leq 20\%$ . (see section 12 for calculations) In addition, the percent difference or drift of all analytes must be  $\leq 50\%$ , with allowance for up to (4) compounds to be greater than 50%.
- The internal standard response must be within 50-200% of the response in the mid level of the initial calibration.
- The internal standard retention times must be within 30 seconds of the retention times in the mid-level of the initial calibration.
- NOTE: There is no internal standard criteria for samples. Criteria is only for continuing and initial calibrations.
- NOTE: Ohio VAP rules require that any sample with internal standard outliers be reanalyzed. The criteria for acceptance is between 50% and 200% of same internal standard in continuing calibration.

10.5.3.1 If none of the CCCs are required analytes, project specific calibration specifications must be agreed with the client.

10.5.4 Once the above criteria have been met, sample analysis may begin. Initial calibration average RFs (or the calibration curve) will be used for sample quantitation, not the continuing calibration RFs. Analysis may proceed until 12 hours from the injection of the DFTPP have passed. (A sample *injected* less than 12 hours after the DFTPP is acceptable.)



## 11 PROCEDURE

### 11.1 Sample Preparation

Samples are prepared following SOP CORP-OP-0001NC.

### 11.2 Sample Analysis Procedure

- 11.2.1 Calibrate the instrument as described in section 10. Depending on the target compounds required by the client, it may be necessary to use more than one calibration standard.
- 11.2.2 All samples must be analyzed using the same instrument conditions as the preceding continuing calibration standard.
- 11.2.3 Add internal standard to the extract to result in 2 ng injected on column. Mix thoroughly before injection into the instrument.
- 11.2.4 Inject the sample extract into the GC/MS system using the same injection technique as used for the standards.
- 11.2.5 The data system will determine the concentration of each analyte in the extract using calculations equivalent to those in section 12. Quantitation is based on the initial calibration, not the continuing calibration.
- 11.2.6 Identified compounds are reviewed for proper integration. Manual integrations are performed if necessary and are documented by the analyst or automatically by the data system.
- 11.2.7 Target compounds identified by the data system are evaluated using the criteria listed in section 12.1.
- 11.2.8 Library searches of peaks present in the chromatogram that are not target compounds (Tentatively Identified Compounds, TIC) may be performed if required by the client. They are evaluated using the criteria in section 12.3.

### 11.3 Dilutions

If the response for any compound exceeds the working range of the GC/MS system, a dilution of the extract is prepared and analyzed. An appropriate dilution should be in the

upper half of the calibration range. Samples may be screened to determine the appropriate dilution for the initial run. If the initial diluted run has no hits or hits below 20% of the calibration range and the matrix allows for analysis at a lesser dilution, the sample must be reanalyzed at a dilution targeted to bring the largest hit above 50% of the calibration range.

#### 11.3.1 Guidance for Dilutions Due to Matrix

If the sample is initially run at a dilution and the baseline rise is less than the height of the internal standards, or if individual non-target peaks are less than two times the height of the internal standards, the sample should be reanalyzed at a more concentrated dilution. This requirement is approximate and subject to analyst judgement. For example, samples containing organic acids may need to be analyzed at a higher dilution to avoid destroying the column.

#### 11.3.2 Reporting Dilutions

The most concentrated dilution with target compounds within the calibration range will be reported. Other dilutions will only be reported at client request.

- 11.4 Perform all qualitative and quantitative measurements. When the extracts are not being used for analyses, refrigerate them at  $4 \pm 2^{\circ}\text{C}$ , protected from light in screw cap vials equipped with unpierced Teflon lined septa.

#### 11.5 Retention time criteria for samples

If the retention time for any internal standard changes by more than 0.5 minutes from the last continuing calibration standard, the chromatographic system must be inspected for malfunctions and corrected. Reanalysis of samples analyzed while the system was malfunctioning is required.

- 11.5.1 If the retention time of any internal standard in any sample varies by more than 0.1 minute from the preceding continuing calibration standard, the data must be carefully evaluated to ensure that no analytes have shifted outside their retention time windows.

#### 11.6 Procedural Variations

- 11.6.1 One-time procedural variations are allowed only if deemed necessary in the professional judgment of supervision to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and approved by a Technical Specialist and QA Manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file. Any unauthorized deviations from this procedure must also be documented as a non-

conformance, with a cause and corrective action described.

## 11.7 Troubleshooting Guide

### 11.7.1 Daily Instrument Maintenance

In addition to the checks listed in the instrument maintenance schedule in the STL North Canton Quality Assurance Manual (LQM), current version, the following daily maintenance should be performed.

- Clip Column as necessary.
- Install new or cleaned injection port liner as necessary.
- Install new septum as necessary.
- Perform autotune.

### 11.7.2 Major Maintenance

A new initial calibration is necessary following major maintenance. Major maintenance includes changing the column, cleaning the source, and replacing the multiplier. Refer to the manufacturer's manual for specific guidance.

## 12 DATA ANALYSIS AND CALCULATIONS

### 12.1 Qualitative identification

An analyte is identified by retention time and by comparison of the sample mass spectrum with the mass spectrum of a standard of the suspected compound (standard reference spectrum). Mass spectra for standard reference may be obtained on the user's GC/MS by analysis of the calibration standards or from the NBS library. Two criteria must be satisfied to verify identification: (1) elution of sample component at the same GC retention time as the standard component; and (2) correspondence of the sample component and the standard component characteristic ions. (Note: Care must be taken to ensure that spectral distortion due to co-elution is evaluated.)

- The sample component retention time must compare to within  $\pm 0.2$  min. of the retention time of the standard component. For reference, the standard must be run within the same twelve hours as the sample.
- All ions present in the standard mass spectra at a relative intensity greater than 10% (most abundant ion in the spectrum equals 100%) should be present in the sample

spectrum.

- The characteristic ions of a compound must maximize in the same scan or within one scan of each other.
- The relative intensities of ions should agree to within  $\pm 30\%$  between the standard and sample spectra. (Example: For an ion with an abundance of 50% in the standard spectra, the corresponding sample abundance must be between 20% and 80%.)

12.1.1 If a compound cannot be verified by all the above criteria, but in the technical judgment of the analyst the identification is correct, the analyst shall report that identification and proceed with quantitation.

## 12.2 Mass chromatogram searches.

Certain compounds are unstable in the calibration standard and cannot be calibrated in the normal way. In particular, the compound hexachlorophene (CAS 70-30-4) falls into this category, and is required for Appendix IX analysis. For this analyte a mass chromatogram search is made.

### 12.2.1 Hexachlorophene

Display the mass chromatograms for mass 196 and mass 198 for the region of the chromatogram from at least 2 minutes before chrysene-d12 to at least 4 minutes after chrysene-d12. If peaks for both ions coincide then the analyst evaluates the spectrum for the presence of hexachlorophene. No quantitation is possible.

12.3 For samples containing components not associated with the calibration standards, a library search may be made for the purpose of tentative identification. The necessity to perform this type of identification will be determined by the type of analyses being conducted. Computer generated library search routines should not use normalization routines that would misrepresent the library or unknown spectra when compared to each other. Only after visual comparison of sample spectra with the nearest library searches shall the mass spectral interpretation specialist assign a tentative identification. Guidelines for making tentative identification are:

- Relative intensities of major ions in the reference spectrum (ions  $>10\%$  of the most abundant ion) should be present in the sample spectrum.
- The relative intensities of the major ions should agree within  $\pm 20\%$ . (Example: For an ion with an abundance of 50% in the standard spectrum, the corresponding sample ion abundance should be between 30% and 70%.)

- Molecular ions present in the reference spectrum should be present in the sample spectrum.
- Ions present in the sample spectrum, but not in the reference spectrum, should be reviewed for possible background contamination or presence of coeluting compounds.
- Ions present in the reference spectrum, but not in the sample spectrum, should be reviewed for possible subtraction from the sample spectrum because of background contamination or coeluting peaks. Data system library reduction programs can sometimes create these discrepancies.
- Automatic background subtraction can severely distort spectra from samples with unresolved hydrocarbons.

12.4 Anyone evaluating data is trained to know how to handle isomers with identical mass spectra and close elution times. These include:

Dichlorobenzenes  
Methylphenols  
Trichlorophenols  
Phenanthrene, anthracene  
Fluoranthene, pyrene  
Benzo(b) and (k)fluoranthene  
Chrysene, benzo(a)anthracene

Extra precautions concerning these compounds are to more closely scrutinize retention time vs. the calibration standard and also to check that all isomers have distinct retention times.

A second category of problem compounds would be the poor responders or compounds that chromatograph poorly. Included in this category would be:

Benzoic acid  
Chloroanilines  
Nitroanilines  
2,4-Dinitrophenol  
4-Nitrophenol  
Pentachlorophenol  
3,3'-Dichlorobenzidine  
Benzyl alcohol  
4,6-Dinitro-2-methylphenol

Manually checking the integrations would be appropriate for these compounds.

## 12.5 Calculations

## 12.5.1 Percent Relative Standard Deviation for Initial Calibration

$$\%RSD = \frac{SD}{RF} \times 100$$

$RF$  = Mean of RFs from initial calibration for a compound

$SD$  = Standard deviation of RFs from initial calibration for a compound,

$$= \sqrt{\frac{\sum_{i=1}^N (RF_i - \overline{RF})^2}{N - 1}}$$

$RF_i$  = RF for each of the calibration levels

$N$  = Number of RF values

## 12.5.2 Continuing calibration percent drift

$$\%Drift = \frac{C_{actual} - C_{found}}{C_{actual}} \times 100\%$$

$C_{actual}$  = Known concentration in standard

$C_{found}$  = Measured concentration using selected quantitation method

## 12.5.3 Concentration in the extract

The concentration of each identified analyte and surrogate in the extract is calculated from the linear or quadratic curve fitted to the initial calibration points, or from the average RF of the initial calibration. For South Carolina Certification the concentration of each identified analyte and surrogate in the extract is calculated from the linear curve only fitted to the initial calibration points, or from the average RF of the initial calibration.

## 12.5.3.1 Average response factor

If the average of all the %RSDs of the response factors in the initial calibration is  $\leq 15\%$ , the average response factor from the initial calibration may be used for quantitation.

$$C_{ex} = \frac{R_i C_s}{\overline{R_i RF}}$$

## 12.5.3.2 Linear fit (Use only Linear fit for South Carolina Certification)

$$C_{ex} = A + B \frac{(R_x C_{is})}{R_{is}}$$

$C_{ex}$  = Concentration in extract, µg/mL

$R_x$  = Response for analyte

$C_{is}$  = Concentration of internal standard

$A$  = Intercept

$B$  = Slope

## 12.5.3.3 Quadratic fit

$$C_{ex} = A + B \left( \frac{R_x C_{is}}{R_{is}} \right) + C \left( \frac{R_x C_{is}^2}{R_{is}} \right)$$

$C$  = Curvature

12.5.4 The concentration in the sample is then calculated.

#### 12.5.4.1 Aqueous Calculation

$$\text{Concentration, } \mu\text{g} / \text{L} = \frac{C_{ex} V_t}{V_o}$$

Where:

$V_t$  = Volume of total extract,  $\mu\text{L}$ , taking into account dilutions (i.e., a 1-to-10 dilution of a 1 mL extract will mean  $V_t = 10,000 \mu\text{L}$ . If half of the base/neutral extract and half of the acid extract are combined,  $V_t = 2,000$ .)

$V_o$  = Volume of water extracted (mL)

12.5.5 Sediment/Soil, Sludge (on a dry-weight basis) and Waste (normally on a wet-weight basis):

$$\text{Concentration, } \mu\text{g} / \text{kg} = \frac{C_{ex} V_t}{W_s D}$$

$W_s$  = Weight of sample extracted or diluted in grams

$D$  = (100 - % moisture in sample)/100, for a dry weight basis or 1 for a wet weight basis

12.6 MS/MSD percent recovery calculation.

$$\text{Matrix Spike Recovery} = \frac{S_{SR} - S_R}{S_A} \times 100\%$$

$S_{SR}$  = Spike sample result

$S_R$  = Sample result

$S_A$  = Spike added



## 12.7 Relative % Difference calculation for the MS/MSD

$$RPD = \frac{MS_R - MSD_R}{1/2(MS_R + MSD_R)} \times 100$$

$RPD$  = Relative percent difference

$MS_R$  = Matrix spike result

$MSD_R$  = Matrix spike duplicate result

## 12.8 Relative response factor calculation.

$$RF = \frac{A_x C_{is}}{A_{is} C_x}$$

$A_x$  = Area of the characteristic ion for the compound being measured

$A_{is}$  = Area of the characteristic ion for the specific internal standard

$C_x$  = Concentration of the compound being measured ( $\mu\text{g/L}$ )

$C_{is}$  = Concentration of the specific internal standard ( $\mu\text{g/L}$ )

## 12.9 Calculation of TICs: The calculation of TICs (tentatively identified compounds) is identical to the above calculations with the following exceptions:

$A_x$  = Area of the total ion chromatogram for the compound being measured

$A_{is}$  = Area of the total ion chromatogram for the nearest internal standard without interference

$$RF = 1$$

## 12.10 Percent DDT breakdown

$$\% \text{ DDT breakdown} = \frac{\text{DDEarea} + \text{DDDarea}}{\text{DDTarea} + \text{DDEarea} + \text{DDarea}}$$

The total ion current areas are used for this calculation

### 13 METHOD PERFORMANCE

#### 13.1 Method Detection Limit

Each laboratory must generate a valid method detection limit for each analyte of interest. The MDL must be below the reporting limit for each analyte. The procedure for determination of the method detection limit is given in 40 CFR Part 136, Appendix B, and further defined in QA Policy #: QA-005.

#### 13.2 Initial Demonstration

Each laboratory must make an initial demonstration of capability for each individual method. Demonstration of capability for both soil and water matrices is required. This requires analysis of QC check samples containing all of the standard analytes for the method. For some tests it may be necessary to use more than one QC check mix to cover all analytes of interest.

13.2.1 Four aliquots of the QC check sample are analyzed using the same procedures used to analyze samples, including sample preparation.

13.2.2 Calculate the average recovery and standard deviation of the recovery for each analyte of interest.

13.2.3 If any analyte does not meet the acceptance criteria the test must be repeated. Only those analytes that did not meet criteria in the first test need to be evaluated. Repeated failure for any analyte indicates the need for the laboratory to evaluate the analytical procedure and take corrective action.

#### 13.3 Non-standard analytes

For non-standard analytes, an MDL study must be performed and calibration curve generated before analyzing any samples, unless lesser requirements are previously agreed to with the client. In any event, the minimum initial demonstration required is analysis of an extracted standard at the reporting limit and a single point calibration.

#### 13.4 Training Qualification

The group/team leader has the responsibility to ensure that this procedure is performed by an analyst who has been properly trained in its use and has the required experience.

## 14 POLLUTION PREVENTION

- 14.1 This section is not applicable to this procedure.

## 15 WASTE MANAGEMENT

- 15.1 All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."
- 15.2 Laboratory personnel assigned to perform hazardous waste disposal procedures must have a working knowledge of the established procedures and practices of STL. They must have training on the hazardous waste disposal practices upon initial assignment to these tasks, followed by an annual refresher training.
- 15.3 Waste Streams Produced by the Method
- 15.3.1 **Vials containing sample extracts:** These vials are placed in the vial waste located in the GC/MS laboratory.

## 16 REFERENCES

- 16.1 References
- 16.1.1 SW846, Test Methods for Evaluating Solid Waste, Third Edition, Update II, October 1994, Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS): Capillary Column Technique, Method 8270C.
- 16.1.2 J. W. Eichelberger, L. E. Harris, and W. L. Budde, "Reference Compound to Calibrate Ion Abundance Measurement in Gas Chromatography/Mass Spectrometry," Analytical Chemistry, 47, 995 (1975)
- 16.1.3 Corporate Quality Management Plan (QMP), current version.
- 16.1.4 STL Laboratory Quality Manual (LQM), current version.
- 16.2 Associated SOPs and Policies, latest version
- 16.2.1 QA Policy, QA-003
- 16.2.2 Glassware Washing, NC-QA-0014

16.2.3 Statistical Evaluation of Data and Development of Control Charts, NC-QA-0018

16.2.4 Method Detection Limits and Instrument Detection Limits, NC-QA-0021

16.2.5 Navy/Army SOP, NC-QA-0016

**17 MISCELLANEOUS****17.1 Modifications from Reference Method**

17.1.1 A retention time window of 0.2 minutes is used for all components, since some data systems do not have the capability of using the relative retention time units specified in the reference method.

17.1.2 The quantitation and qualifier ions from coms compounds have been changed from those recommended in SW-846 in order to improve the reliability of qualitative identification.

**17.2 Tables****Table 1****STL North Canton Primary Standard and Standard Reporting Limits**

Analytes	CAS Number	Standard Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
Pyridine	110-86-1	20	660
N-nitrosodimethylamine	62-75-9	10	330
Aniline	62-53-3	10	330
Phenol	108-95-2	10	330
Bis(2-chloroethyl)ether	111-44-4	10	330
2-Chlorophenol	95-57-8	10	330
1,3-Dichlorobenzene	541-73-1	10	330
1,4-Dichlorobenzene	106-46-7	10	330
Benzyl alcohol	100-51-6	10	330
1,2-Dichlorobenzene	95-50-1	10	330
2-Methylphenol	95-48-7	10	330
2,2'-oxybis(1-chloropropane) <sup>1</sup>	108-60-1	10	330
4-Methylphenol	106-44-5	10	330
N-Nitroso-di-n-propylamine	621-64-7	10	330
Hexachloroethane	67-72-1	10	330
Nitrobenzene	98-95-3	10	330
Isophorone	78-59-1	10	330
2-Nitrophenol	88-75-5	10	330
2,4-Dimethylphenol	105-67-9	10	330
Benzoic acid	65-85-0	50	1600
Bis(2-chloroethoxy)methane	111-91-1	10	330

8771342

GC/MS ANALYSIS BASED ON METHOD 8270C

SOP No: CORP-MS-0001NC

Revision No: 2.10

Revision Date: 02/24/04

Page 30 of 50

Table 1

## STL North Canton Primary Standard and Standard Reporting Limits

Analytes	CAS Number	Standard Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
2,4-Dichlorophenol	120-83-2	10	330
1,2,4-Trichlorobenzene	120-82-1	10	330
Naphthalene	91-20-3	10	330
4-Chloroaniline	106-47-8	10	330
Hexachlorobutadiene	87-68-3	10	330
4-Chloro-3-methylphenol	59-50-7	10	330
2-Methylnaphthalene	91-57-6	10	330
Hexachlorocyclopentadiene	77-47-4	50	1600
2,4,6-Trichlorophenol	88-06-2	10	330
2,4,5-Trichlorophenol	95-95-4	10	330
2-Chloronaphthalene	91-58-7	10	330
2-Nitroaniline	88-74-4	50	1600
Dimethyl phthalate	131-11-3	10	330
Acenaphthylene	208-96-8	10	330
3-Nitroaniline	99-09-2	50	1600
Acenaphthene	83-32-9	10	330
2,4-Dinitrophenol	51-28-5	50	1600
4-Nitrophenol	100-02-7	50	1600
Dibenzofuran	132-64-9	10	330
2,4-Dinitrotoluene	121-14-2	10	330
2,6-Dinitrotoluene	606-20-2	10	330
Diethylphthalate	84-66-2	10	330
4-Chlorophenyl phenyl ether	7005-72-3	10	330
Fluorene	86-73-7	10	330
4-Nitroaniline	100-01-6	50	1600
4,6-Dinitro-2-methylphenol	534-52-1	50	1600
N-Nitrosodiphenylamine	86-30-6	10	330
Azobenzene	103-33-3	10	330
4-Bromophenyl phenyl ether	101-55-3	10	330
Hexachlorobenzene	118-74-1	10	330
Pentachlorophenol	87-86-5	50	1600
Phenanthrene	85-01-8	10	330
Anthracene	120-12-7	10	330
Carbazole	86-74-8	10	330
Di-n-butyl phthalate	84-74-2	10	330
Fluoranthene	206-44-0	10	330
Benzidine	92-87-5	100	3300
Pyrene	129-00-0	10	330
Butyl benzyl phthalate	85-68-7	10	330
3,3'-Dichlorobenzidine	91-94-1	50	1600
Benzo(a)anthracene	56-55-3	10	330
Bis(2-ethylhexyl)phthalate	117-81-7	10	330
Chrysene	218-01-9	10	330

Table 1

## STL North Canton Primary Standard and Standard Reporting Limits

Analytes	CAS Number	Standard Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
Di-n-octylphthalate	117-84-0	10	330
Benzo(b)fluoranthene	205-99-2	10	330
Benzo(k)fluoranthene	207-08-9	10	330
Benzo(a)pyrene	50-32-8	10	330
Indeno(1,2,3-cd)pyrene	193-39-5	10	330
Dibenz(a,h)anthracene	53-70-3	10	330
Benzo(g,h,i)perylene	191-24-2	10	330
Benzaldehyde	100-52-7	10	330
Caprolactam	105-60-2	10	330
1,1-Biphenyl	92-52-4	10	330
Atrazine	1912-24-9	10	330
Benzenethiol	108-98-5	10	330
Indene	95-13-6	10	330
Quinoline	91-22-5	10	330
1-Methyl Naphthalene	90-12-0	10	330

<sup>1</sup> 2,2'-oxybis(1-chloropropane) was formerly known as bis(2-chloroisopropyl)ether.

Table 2

STL North Canton Appendix IX<sup>1</sup> Standard Reporting Limits

Semivolatiles	CAS Number	Standard Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
2-Picoline	109-06-8	20	660
N-Nitrosomethylethylamine	10595-95-6	10	330
Methyl methanesulfonate	66-27-3	10	330
N-Nitrosodiethylamine	55-18-5	10	330
Ethyl methanesulfonate	62-50-0	10	330
Pentachloroethane	76-01-7	50	1600
Acetophenone	98-86-2	10	330
N-Nitrosopyrrolidine	930-55-2	10	330
N-Nitrosomorpholine	59-89-2	10	330
o-Toluidine	95-53-4	20	660
3-Methylphenol	108-39-4	10	330
N-Nitrosopiperidine	100-75-4	10	330
o,o,o-Triethyl-Phosphorothioate <sup>2</sup>	126-68-1	50	1600
a,a-Dimethyl-phenethylamine	122-09-8	50	1600
2,6-Dichlorophenol	87-65-0	10	330
Hexachloropropene	1888-71-7	100	3300
p-Phenylenediamine	106-50-3	100	3300

Table 2

STL North Canton Appendix IX<sup>1</sup> Standard Reporting Limits

Semivolatiles	CAS Number	Standard Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
n-Nitrosodi-n-butylamine	924-16-3	10	330
Safrole	94-59-7	20	660
1,2,4,5-Tetrachlorobenzene	95-94-3	10	330
Isosafrole	120-58-1	20	660
1,4-Dinitrobenzene	100-25-4	10	330
1,4-Naphthoquinone	130-15-4	50	1600
1,3-Dinitrobenzene	99-65-0	10	330
Pentachlorobenzene	608-93-5	10	330
1-Naphthylamine	134-32-7	10	330
2-Naphthylamine	91-59-8	10	330
2,3,4,6-Tetrachlorophenol	58-90-2	50	1600
5-Nitro-o-toluidine	99-55-8	20	660
Thionazin <sup>2</sup>	297-97-2	50	1600
1,3,5-Trinitrobenzene	99-35-4	50	1600
Sulfotepp <sup>2</sup>	3689-24-5	50	1600
Phorate <sup>2</sup>	298-02-2	50	1600
Phenacetin	62-44-2	20	660
Diallate <sup>2</sup>	2303-16-4	20	660
Dimethoate <sup>2</sup>	60-51-5	20	660
4-Aminobiphenyl	92-67-1	50	1600
Pentachloronitrobenzene	82-68-8	50	1600
Pronamide	23950-58-5	20	660
Disulfoton <sup>2</sup>	298-04-4	50	1600
2-secbutyl-4,6-dinitrophenol (Dinoseb <sup>2</sup> )	88-85-7	20	660
4-Nitroquinoline-1-oxide	56-57-5	100	3300
Methapyrilene	91-80-5	50	1600
Aramite	140-57-8	20	660
Famphur <sup>3</sup>	52-85-7	100	3300
p-(Dimethylamino)azobenzene	60-11-7	20	660
p-Chlorobenzilate	510-15-6	10	330
3,3'-Dimethylbenzidine	119-93-7	50	1600
2-Acetylaminofluorene	53-96-3	100	3300
Dibenz(a,j)acridine	224-42-0	20	660
7,12-Dimethylbenz(a)anthracene	57-97-6	20	660
3-Methylcholanthrene	56-49-5	20	660

<sup>1</sup> The Appendix IX standard contains additional analytes required for the Appendix IX list. The STL North Canton primary standard must also be analyzed to include all of the Appendix IX list

<sup>2</sup> May also be analyzed by method 8141, which can achieve lower reporting limits.

<sup>3</sup> It is highly recommended that Famphur is analyzed by method 8081. It is a poor responder by 8270C.

**Table 2A**  
**STL North Canton Michigan Program<sup>1</sup>**

Semivolatile	CAS Number	Michigan Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
Acenaphthene	83-32-9	5	330
Acenaphthylene	208-96-8	5	330
Acetophenone	98-86-2	5	330
Anthracene	120-12-7	5	330
Atrazine	1912-24-9	5	330
Benzaldehyde	100-52-7	10	330
Benzo(a)anthracene	56-55-3	1	330
Benzo(a)pyrene	50-32-8	2	330
Benzo(b)fluoranthene	205-99-2	2	330
Benzo(g,h,i)perylene	191-24-2	5	330
Benzo(k)fluoranthene	207-08-9	5	330
1,1'-Biphenyl	92-52-4	10	330
4-Bromophenylphenyl ether	101-55-3	5	330
Butylbenzylphthalate	85-68-7	5	330
di-n-Butylphthalate	84-74-2	5	330
Caprolactam	105-60-2	10	330
Carbazole	86-74-8	10	330
4-Chloroaniline	106-47-8	20	1700
bis(2-Chloroethoxy)methane	111-91-1	5	330
bis(2-Chloroethyl)ether	111-44-4	4	330
bis(2-Chloroisopropyl)ether	108-60-1	5	330
4-Chloro-3-Methylphenol	59-50-7	5	330
2-Chloronaphthalene	91-58-7	5	330
2-Chlorophenol	95-57-8	5	330
4-Chlorophenyl phenyl ether	7005-72-3	5	330
Chrysene	218-01-9	5	330
Dibenz(a,h)anthracene	53-70-3	2	330
Dibenzofuran	132-64-9	5	330
3,3'-Dichlorobenzidine	91-94-1	4	2000
2,4-Dichlorophenol	120-83-2	10	330
Diethylphthalate	84-66-2	5	330
2-4-Dimethylphenol	105-67-9	5	330
Dimethylphthalate	131-11-3	5	330
4,6-Dinitro-2-methylphenol	534-52-1	20	1700
2,4-Dinitrophenol	51-28-5	20	1700
2,4-Dinitrotoluene	121-14-2	5	330
2,6-Dinitrotoluene	606-20-2	5	330
bis(2-Ethylhexyl)phthalate	117-81-7	5	330
Fluoranthene	206-44-0	5	330
Fluorene	86-73-7	5	330



8771346

GC/MS ANALYSIS BASED ON METHOD 8270C

SOP No: CORP-MS-0001NC

Revision No: 2.10

Revision Date: 02/24/04

Page 34 of 50

Table 2A STL North Canton Michigan Program <sup>1</sup>			
Semivolatile	CAS Number	Michigan Reporting Limits	
		Aqueous µg/L	Low Soil/Sediment µg/kg
Hexachlorobenzene	118-74-1	5	330
Hexachlorobutadiene	87-68-3	5	330
Hexachlorocyclopentadiene	77-47-4	5	330
Hexachloroethane	67-72-1	5	330
Indeno(1,2,3-cd)pyrene	193-39-5	2	330
Isophorone	78-59-1	5	330
2-Methylnaphthalene	91-57-6	5	330
2-Methylphenol	95-48-7	5	330
4-Methylphenol	106-44-5	5	330
Naphthalene	91-20-3	5	330
2-Nitroaniline	88-74-4	20	1700
3-Nitroaniline	99-09-2	20	1700
4-Nitroaniline	100-01-6	20	1700
Nitrobenzene	95-95-3	4	330
2-Nitrophenol	88-75-5	5	330
4-Nitrophenol	100-02-7	20	1700
N-Nitroso-di-n-propylamine	621-64-7	5	330
N-Nitrosodiphenylamine (diphenylamine)	62-75-9	5	330
di-n-Octylphthalate	117-84-0	5	330
Pentachlorophenol	87-86-5	20	800
Phenanthrene	85-01-8	5	330
Phenol	108-95-2	5	330
Pyrene	129-00-0	5	330
2,4,5-Trichlorophenol	95-95-4	5	330
2,4,6-Trichlorophenol	88-06-2	4	330

<sup>1</sup> Reporting Limits are only for samples performed under the Michigan program.

Table 3

## Reportable Analytes for STL North Canton Standard Tests, Primary Standard

Analyte	CAS Number	TCLP	TCL	Appendix IX
Pyridine	110-86-1	X		X
N-nitrosodimethylamine	62-75-9			X
Aniline	62-53-3			X
Phenol	108-95-2		X	X
Bis(2-chloroethyl)ether	111-44-4		X	X
2-Chlorophenol	95-57-8		X	X
1,3-Dichlorobenzene	541-73-1		X	X
1,4-Dichlorobenzene	106-46-7	X	X	X
Benzyl alcohol	100-51-6			X
1,2-Dichlorobenzene	95-50-1		X	X
2-Methylphenol	95-48-7	X	X	X

Table 3

## Reportable Analytes for STL North Canton Standard Tests, Primary Standard

Analyte	CAS Number	TCLP	TCL	Appendix IX
2,2'-oxybis(1-chloropropane)	180-60-1		X	X
4-Methylphenol	106-44-5	X	X	X
N-Nitroso-di-n-propylamine	621-64-7		X	X
Hexachloroethane	67-72-1	X	X	X
Nitrobenzene	98-95-3	X	X	X
Isophorone	78-59-1		X	X
2-Nitrophenol	88-75-5		X	X
2,4-Dimethylphenol	105-67-9		X	X
Benzoic acid	65-85-0			
Bis(2-chloroethoxy)methane	111-91-1		X	X
2,4-Dichlorophenol	120-83-2		X	X
1,2,4-Trichlorobenzene	120-82-1		X	X
Naphthalene	91-20-3		X	X
4-Chloroaniline	106-47-8		X	X
Hexachlorobutadiene	87-68-3	X	X	X
4-Chloro-3-methylphenol	59-50-7		X	X
2-Methylnaphthalene	91-57-6		X	X
Hexachlorocyclopentadiene	77-47-4		X	X
2,4,6-Trichlorophenol	88-06-2	X	X	X
2,4,5-Trichlorophenol	95-95-4	X	X	X
2-Chloronaphthalene	91-58-7		X	X
2-Nitroaniline	88-74-4		X	X
Dimethyl phthalate	131-11-3		X	X
Acenaphthylene	208-96-8		X	X
3-Nitroaniline	99-09-2		X	X
Acenaphthene	83-32-9		X	X
2,4-Dinitrophenol	51-28-5		X	X
4-Nitrophenol	100-02-7		X	X
Dibenzofuran	132-64-9		X	X
2,4-Dinitrotoluene	121-14-2	X	X	X
2,6-Dinitrotoluene	606-20-2		X	X
Diethylphthalate	84-66-2		X	X
4-Chlorophenyl phenyl ether	7005-72-3		X	X
Fluorene	86-73-7		X	X
4-Nitroaniline	100-01-6		X	X
4,6-Dinitro-2-methylphenol	534-52-1		X	X
N-Nitrosodiphenylamine	86-30-6		X	X
Azobenzene <sup>1</sup>	103-33-3			
4-Bromophenyl phenyl ether	101-55-3		X	X
Hexachlorobenzene	118-74-1	X	X	X
Pentachlorophenol	87-86-5	X	X	X
Phenanthrene	85-01-8		X	X
Anthracene	120-12-7		X	X
Carbazole	86-74-8		X	
	84-74-2		X	X
Fluoranthene	206-44-0		X	X

Table 3

## Reportable Analytes for STL North Canton Standard Tests, Primary Standard

Analyte	CAS Number	TCLP	TCL	Appendix IX
Benzidine	92-87-5			
Pyrene	129-00-0		X	X
Butyl benzyl phthalate	85-68-7		X	X
3,3'-Dichlorobenzidine	91-94-1		X	X
Benzo(a)anthracene	56-55-3		X	X
Bis(2-ethylhexyl)phthalate	117-81-7		X	X
Chrysene	218-01-9		X	X
Di-n-octylphthalate	117-84-0		X	X
Benzo(b)fluoranthene	205-99-2		X	X
Benzo(k)fluoranthene	207-08-9		X	X
Benzo(a)pyrene	50-32-8		X	X
Indeno(1,2,3-cd)pyrene	193-39-5		X	X
Dibenz(a,h)anthracene	53-70-3		X	X
Benzo(g,h,i)perylene	191-24-2		X	X
Benzaldehyde	100-52-7		X	
Caprolactam	105-60-2		X	
1,1-Biphenyl	92-52-4		X	
Atrazine	1912-24-9		X	

<sup>1</sup> Azobenzene is formed by decomposition of 1,2-diphenylhydrazine. If 1,2-diphenylhydrazine is requested, it will be analyzed as azobenzene

Table 4

## Reportable analytes for STL North Canton Standard Tests, Appendix IX Standard

Semivolatiles	CAS Number	TCLP	TCL	Appendix IX
2-Picoline	109-06-8			X
N-Nitrosomethylethylamine	10595-95-6			X
Methyl methanesulfonate	66-27-3			X
N-Nitrosodiethylamine	55-18-5			X
Ethyl methanesulfonate	62-50-0			X
Pentachloroethane	76-01-7			X
Acetophenone	98-86-2		X	X
N-Nitrosopyrrolidine	930-55-2			X
N-Nitrosomorpholine	59-89-2			X
o-Toluidine	95-53-4			X
3-Methylphenol	108-39-4			X
N-Nitrosopiperidine	100-75-4			X
o,o,o-Triethyl-Phosphorothioate <sup>1</sup>	126-68-1			X
a,a-Dimethyl-phenethylamine	122-09-8			X
2,6-Dichlorophenol	87-65-0			X
Hexachloropropene	1888-71-7			X
p-Phenylenediamine	106-50-3			X
n-Nitrosodi-n-butylamine	924-16-3			X
Safrole	94-59-7			X
1,2,4,5-Tetrachlorobenzene	95-94-3			X

Table 4

## Reportable analytes for STL North Canton Standard Tests, Appendix IX Standard

Semivolatiles	CAS Number	TCLP	TCL	Appendix IX
Isosafrole	120-58-1			X
1,4-Dinitrobenzene	100-25-4			X
1,4-Naphthoquinone	130-15-4			X
1,3-Dinitrobenzene	99-65-0			X
Pentachlorobenzene	608-93-5			X
1-Naphthylamine	134-32-7			X
2-Naphthylamine	91-59-8			X
2,3,4,6-Tetrachlorophenol	58-90-2			X
5-Nitro-o-toluidine	99-55-8			X
Thionazin <sup>1</sup>	297-97-2			X
1,3,5-Trinitrobenzene	99-35-4			X
Sulfotepp <sup>1</sup>	3689-24-5			X
Phorate <sup>1</sup>	298-02-2			X
Phenacetin	62-44-2			X
Diallate	2303-16-4			X
Dimethoate <sup>1</sup>	60-51-5			X
4-Aminobiphenyl	92-67-1			X
Pentachloronitrobenzene	82-68-8			X
Pronamide	23950-58-5			X
Disulfoton <sup>1</sup>	298-04-4			X
2-secbutyl-4,6-dinitrophenol (Dinoseb) <sup>1</sup>	88-85-7			X
4-Nitroquinoline-1-oxide	56-57-5			X
Famphur <sup>2</sup>	52-85-7			X
Methapyrilene	91-80-5			X
Aramite	140-57-8			X
p-(Dimethylamino)azobenzene	60-11-7			X
p-Chlorobenzilate	510-15-6			X
3,3'-Dimethylbenzidine	119-93-7			X
2-Acetylaminofluorene	53-96-3			X
Dibenz(a,j)acridine <sup>3</sup>	224-42-0			X
7,12-Dimethylbenz(a)anthracene	57-97-6			X
3-Methylcholanthrene	56-49-5			X
Hexachlorophene <sup>4</sup>	70-30-4			X
Diphenylamine <sup>5</sup>	122-39-4			X

<sup>1</sup> May also be analyzed by method 8141, which can achieve lower reporting limits.

<sup>2</sup> May also be analyzed by method 8081, which can achieve lower reporting limits

<sup>3</sup> Skinner List Compound

<sup>4</sup> Hexachlorophene is a required analyte for Appendix IX. This compound is not stable, and therefore not included in the calibration standard. The characteristic ions for hexachlorophene are searched for in the chromatogram. (See section 12.2.1)

<sup>5</sup> Diphenylamine is a required compound for Appendix IX. N-nitrosodiphenylamine decomposes in the injection

port to form diphenylamine. Therefore these two compounds cannot be distinguished. Diphenylamine is not included in the calibration standard.

Table 5

## Suggested Instrumental Conditions

Mass Range	35-500 amu
Scan Time	≤1 second/scan
Initial Column Temperature/Hold Time	45°C for 1 minutes
Column Temperature Program	45- 100°C at 25°C/min for 0 min 100 - 280°C at 30°C/min for 0 min 280 - 100°C at 25°C/min for 2 min
Final Column Temperature/Hold Time	320°C (until at least one minute after benzo(g,h,i)perylene has eluted)
Injector Temperature	250 - 300°C
Transfer Line Temperature	250 - 300°C
Source Temperature	According to manufacturer's Specifications
Injector	Grob-type, split / splitless
Sample Volume	0.5 µl
Carrier Gas	Helium at 30 cm/sec

Table 6

## DFTPP Key Ions and Ion Abundance Criteria

Mass	Ion Abundance Criteria
51	30 – 60% of mass 198
68	<2% of mass 69
70	<2% of mass 69
127	40 – 60% of mass 198
197	<1% of mass 198
198	Base peak, 100% relative abundance
199	5 – 9% of mass 198
275	10 – 30% of mass 198
365	>1% of mass 198
441	Present, but less than mass 443
442	>40% of mass 198
443	17 – 23% of mass 442

Table 7

## Analytes in Approximate Retention Time Order and Characteristic Ions, Primary Standard

Analyte	Primary	Secondary	Tertiary
N-nitrosodimethylamine	74	42	
Pyridine	79	52	
<b>2-Fluorophenol (Surrogate Standard)</b>	112	64	63
<b>Phenol-d5 (Surrogate Standard)</b>	99	42	71
Benzaldehyde	77	105	106
Aniline	93	66	
Phenol	94	65	66
Bis(2-chloroethyl)ether	93	63	95
2-Chlorophenol	128	64	130
1,3-Dichlorobenzene	146	148	113
<b>1,4-Dichlorobenzene-d4 (Internal Standard)</b>	152	150	115
1,4-Dichlorobenzene	146	148	113
Benzyl Alcohol	108	79	77
1,2-Dichlorobenzene	146	148	113
2-Methylphenol	108	107	79
2,2'-oxybis(1-chloropropane) <sup>1</sup>	45	77	79
4-Methylphenol	108	107	79
N-Nitroso-di-n-propylamine	70	42	101,130
Hexachloroethane	117	201	199
<b>Nitrobenzene-d5 (Surrogate Standard)</b>	82	128	54
Nitrobenzene	77	123	65
Isophorone	82	95	138
2-Nitrophenol	139	65	109
2,4-Dimethylphenol	107	121	122
Benzoic Acid	122	105	77
Bis(2-chloroethoxy)methane	93	95	123
2,4-Dichlorophenol	162	164	98
1,2,4-Trichlorobenzene	180	182	145
<b>Naphthalene-d8 (Internal Standard)</b>	136	68	54
Naphthalene	128	129	127
4-Chloroaniline	127	129	65
Hexachlorobutadiene	225	223	227
Caprolactam	113	55	56
4-Chloro-3-methylphenol	107	144	142
2-Methylnaphthalene	142	141	115
Hexachlorocyclopentadiene	237	235	272
2,4,6-Trichlorophenol	196	198	200
2,4,5-Trichlorophenol	196	198	200
1,1'-Biphenyl	154	153	76
<b>2-Fluorobiphenyl (Surrogate Standard)</b>	172	171	170
2-Chloronaphthalene	162	164	127
2-Nitroaniline	65	92	138

8771352

## GC/MS ANALYSIS BASED ON METHOD 8270C

SOP No: CORP-MS-0001NC

Revision No: 2.10

Revision Date: 02/24/04

Page 40 of 50

Table 7

Analytes in Approximate Retention Time Order and Characteristic Ions, Primary Standard

Analyte	Primary	Secondary	Tertiary
Dimethylphthalate	163	194	164
Acenaphthylene	152	151	153
2,6-Dinitrotoluene	165	63	89
Acenaphthene-d10 (Internal Standard)	164	162	160
3-Nitroaniline	138	108	92
Acenaphthene	153	152	154
2,4-Dinitrophenol	184	63	154
Dibenzofuran	168	139	84
4-Nitrophenol	109	139	65
2,4-Dinitrotoluene	165	63	89
Diethylphthalate	149	177	150
Fluorene	166	165	167
4-Chlorophenylphenylether	204	206	141
4-Nitroaniline	138	92	108
4,6-Dinitro-2-methylphenol	198	182	77
N-Nitrosodiphenylamine	169	168	167
2,4,6-Tribromophenol (Surrogate Standard)	330	332	141
Azobenzene	77	182	105
4-Bromophenylphenylether	248	250	141
Hexachlorobenzene	284	142	249
Atrazine	200	173	215
Pentachlorophenol	266	264	268
Phenanthrene-d10 (Internal Standard)	188	94	80
Phenanthrene	178	179	176
Anthracene	178	179	176
Carbazole	167	166	139
Di-n-butylphthalate	149	150	104
Fluoranthene	202	101	100
Benzidine	184	92	185
Pyrene	202	101	100
Terphenyl-d14 (Surrogate Standard)	244	122	212
Butylbenzylphthalate	149	91	206
Benzo(a)Anthracene	228	229	226
Chrysene-d12 (Internal Standard)	240	120	236
3,3'-Dichlorobenzidine	252	254	126
Chrysene	228	226	229
Bis(2-ethylhexyl)phthalate	149	167	279
Di-n-octylphthalate	149	167	43
Benzo(b)fluoranthene	252	253	125
Benzo(k)fluoranthene	252	253	125
Benzo(a)pyrene	252	253	125
Perylene-d12 (Internal Standard)	264	260	265
Indeno(1,2,3-cd)pyrene	276	138	277
Dibenz(a,h)anthracene	278	139	279
Benzo(g,h,i)perylene	276	138	277

Table 8

## Additional Appendix IX Analytes in Approximate Retention Time Order and Characteristic Ions

Analyte	Primary	Secondary	Tertiary
2-Picoline	93	66	92
N-Nitrosomethylethylamine	88	42	43
Methyl methanesulfonate	80	79	65
N-Nitrosodiethylamine	102	44	57
Ethyl methanesulfonate	79	109	97
Pentachloroethane	117	119	167
Acetophenone	105	77	120
N-Nitrosopyrrolidine	100	41	42
N-Nitrosomorpholine	116	56	86
o-Toluidine	106	107	
3-Methylphenol	108	107	77
N-Nitrosopiperidine	114	42	55
o,o,o-Triethyl-Phosphorothioate	198	121	93
a,a-Dimethyl-phenethylamine	58	91	
2,6-Dichlorophenol	162	164	63
Hexachloropropene	213	215	211
p-Phenylenediamine	108	80	
n-Nitrosodi-n-butylamine	84	57	41
Safrole	162	104	77
1,2,4,5-Tetrachlorobenzene	216	214	218
Isosafrole 1	162	104	131
Isosafrole 2	162	104	131
1,4-Dinitrobenzene	168	75	122
1,4-Naphthoquinone	158	104	102
1,3-Dinitrobenzene	168	75	76
Pentachlorobenzene	250	248	252
1-Naphthylamine	143	115	
2-Naphthylamine	143	115	
2,3,4,6-Tetrachlorophenol	232	230	131
5-Nitro-o-toluidine	152	77	106
Thionazin	97	96	143
1,3,5-Trinitrobenzene	213	75	120
Sulfotepp	97	322	202
Phorate	75	97	121
Phenacetin	108	179	109
Diallate	86	234	
Dimethoate	87	93	125
4-Aminobiphenyl	169		
Pentachloronitrobenzene	237	142	214
Pronamide	173	175	255
Disulfoton	88	97	89
2-secbutyl-4,6-dinitrophenol (Dinoseb)	211	163	147
Methyl parathion	109	125	263
4-Nitroquinoline-1-oxide	190	128	160



8771354

GC/MS ANALYSIS BASED ON METHOD 8270C

SOP No: CORP-MS-0001NC

Revision No: 2.10

Revision Date: 02/24/04

Page 42 of 50

Table 8

## Additional Appendix IX Analytes in Approximate Retention Time Order and Characteristic Ions

Analyte	Primary	Secondary	Tertiary
Famphur	218	125	93
Methapyrilene	97	58	
Aramite 1	185	319	
Aramite 2	185	319	
p-(Dimethylamino)azobenzene	120	225	77
p-Chlorobenzilate	251	139	253
3,3'-Dimethylbenzidine	212	106	
2-Acetylaminofluorene	181	180	223
Dibenz(a,j)acridine	279	280	
7,12-Dimethylbenz(a)anthracene	256	241	120
3-Methylcholanthrene	268	252	253

Table 9

## 8270C LCS Control Compounds

LCS Compounds	Spiking Level, Conc. Added = 20 ug/L
1,2,4-Trichlorobenzene	20
Acenaphthene	20
2,4-Dinitrotoluene	20
Pyrene	20
N-Nitroso-di-n-propylamine	20
1,4-Dichlorobenzene	20
Pentachlorophenol	20
Phenol	20
2-Chlorophenol	20
4-Chloro-3-methylphenol	20
4-Nitrophenol	20

Table 9A 8270C All Analyte Spike Mix			
BNANPDES		Methanol	
	Acenaphthene		100
	Acenaphthylene		100
	Anthracene		100
	Benzo(a)anthracene		100
	Benzo(b)fluoranthene	Methanol	100
	Benzo(k)fluoranthene		100
	Benzo(a)pyrene		100
	Benzo(ghi)perylene		100
	Benzyl butyl phthalate		100
	Bis(2-chloroethyl)ether		100
	Bis(2-chloroethoxy)methane		100
	Bis(2-ethylhexyl)phthalate		100
	Bis(2-chloroisopropyl)ether		100
	4-Bromophenyl phenyl ether		100
	2-Chloronaphthalene		100
	4-Chlorophenyl phenyl ether		100
	Chrysene		100
	Dibenzo(a,h)anthracene		100
	Di-n-butylphthalate		100
	1,3-Dichlorobenzene		100
	1,2-Dichlorobenzene		100
	1,4-Dichlorobenzene		100

8771356

GC/MS ANALYSIS BASED ON METHOD 8270C

SOP No: CORP-MS-0001NC

Revision No: 2.10

Revision Date: 02/24/04

Page 44 of 50

Table 9A 8270C All Analyte Spike Mix			
BNANPDES	3,3'-Dichlorobenzidine	Methanol	100
	Diethyl phthalate		100
	Dimethyl phthalate		100
	2,4-Dinitrotoluene		100
	2,6-Dinitrotoluene		100
	Di-n-octylphthalate		100
	Fluoranthene		100
	Fluorene		100
	Hexachlorobenzene		100
	Hexachlorobutadiene		100
	Hexachloroethane		100
	Indeno(1,2,3-cd)pyrene		100
	Isophorone		100
	Naphthalene		100
	Nitrobenzene		100
	N-Nitrosodi-n-propylamine		100
	Phenanthrene		100
	Pyrene		100
	1,2,4-Trichlorobenzene		100
	4-Chloro-3-methylphenol		100
	2-Chlorophenol		100
	2,4-Dichlorophenol		100
	2,4-Dimethylphenol		100

Table 9A 8270C All Analyte Spike Mix			
BNANPDES		Methanol	
	2,4-Dinitrophenol		100
	2-Methyl-4,6-dinitrophenol		100
	2-Nitrophenol		100
	4-Nitrophenol		100
	Pentachlorophenol		100
	Phenol		100
	2,4,6-Trichlorophenol		100
	Acetophenone		100
	Atrazine		100
	Caprolactum		100
	Benzaldehyde		100
	1,1'-Biphenyl		100
	Safrole		100
	1,4-Dioxane		100
	Pronamide		100
	p-Chlorobenzilate		100
	Phenacetin		100
	Ethyl methanesulfonate		100
	2-Picoline		100
	Phorate		100
	Quinoline		100

Table 10

## TCLP LCS Compounds

LCS Compounds	Spiking Level, mg/L in extract
1,4-Dichlorobenzene	0.08
2,4-Dinitrotoluene	0.08
Hexachlorobenzene	0.08
Hexachlorobutadiene	0.08
Hexachloroethane	0.08
2-Methylphenol	0.08
3-Methylphenol	0.08
4-Methylphenol	0.08
Nitrobenzene	0.08
Pentachlorophenol	0.08
Pyridine	0.08
2,4,5-Trichlorophenol	0.08
2,4,6-Trichlorophenol	0.08

Recovery limits for the LCS and for matrix spikes are generated from historical data and are maintained by the QA department.

Table 11

## 8270C Surrogate Compounds

Surrogate Compounds	Spiking Level, Conc. Added = 20 ug/L / 30 ug/L
Nitrobenzene-d5	20
2-Fluorobiphenyl	20
Terphenyl-d14	20
1,2-Dichlorobenzene-d4 <sup>1</sup>	20
Phenol-d5	30
2-Fluorophenol	30
2,4,6-Tribromophenol	30
2-Chlorophenol-d4 <sup>1</sup>	30

Recovery limits for surrogates are generated from historical data and are maintained by the QA department.

**Table 12**  
**Calibration Ranges, µg/mL**

Analyte	Calibration Range
Pyridine	0.25-12.5 ug/mL
N-nitrosodimethylamine	0.25-12.5 ug/mL
Aniline	0.25-12.5 ug/mL
Phenol	0.25-12.5 ug/mL
Bis(2-chloroethyl)ether	0.25-12.5 ug/mL
2-Chlorophenol	0.25-12.5 ug/mL
1,3-Dichlorobenzene	0.25-12.5 ug/mL
1,4-Dichlorobenzene	0.25-12.5 ug/mL
Benzyl alcohol	0.25-12.5 ug/mL
1,2-Dichlorobenzene	0.25-12.5 ug/mL
2-Methylphenol	0.25-12.5 ug/mL
2,2'-oxybis(1-chloropropane) <sup>1</sup>	0.25-12.5 ug/mL
4-Methylphenol	0.25-12.5 ug/mL
N-Nitroso-di-n-propylamine	0.25-12.5 ug/mL
Hexachloroethane	0.25-12.5 ug/mL
Nitrobenzene	0.25-12.5 ug/mL
Isophorone	0.25-12.5 ug/mL
2-Nitrophenol	0.25-12.5 ug/mL
2,4-Dimethylphenol	0.25-12.5 ug/mL
Benzoic acid	0.25-12.5 ug/mL
Bis(2-chloroethoxy)methane	0.25-12.5 ug/mL
2,4-Dichlorophenol	0.25-12.5 ug/mL
1,2,4-Trichlorobenzene	0.25-12.5 ug/mL
Naphthalene	0.05-10 ug/mL
4-Chloroaniline	0.25-12.5 ug/mL
Hexachlorobutadiene	0.25-12.5 ug/mL
4-Chloro-3-methylphenol	0.25-12.5 ug/mL
2-Methylnaphthalene	0.05-10 ug/mL
Hexachlorocyclopentadiene	0.25-12.5 ug/mL
2,4,6-Trichlorophenol	0.25-12.5 ug/mL
2,4,5-Trichlorophenol	0.25-12.5 ug/mL
2-Chloronaphthalene	0.25-12.5 ug/mL
2-Nitroaniline	0.25-12.5 ug/mL
Dimethyl phthalate	0.25-12.5 ug/mL
Acenaphthylene	0.05-10 ug/mL
3-Nitroaniline	0.25-12.5 ug/mL
Acenaphthene	0.05-10 ug/mL
2,4-Dinitrophenol	0.25-12.5 ug/mL
4-Nitrophenol	0.25-12.5 ug/mL
Dibenzofuran	0.25-12.5 ug/mL
2,4-Dinitrotoluene	0.25-12.5 ug/mL
2,6-Dinitrotoluene	0.25-12.5 ug/mL
Diethylphthalate	0.25-12.5 ug/mL
4-Chlorophenyl phenyl ether	0.25-12.5 ug/mL
Fluorene	0.05-10 ug/mL

**Table 12**  
**Calibration Ranges, µg/mL**

Analyte	Calibration Range
4-Nitroaniline	0.25-12.5 ug/mL
4,6-Dinitro-2-methylphenol	0.25-12.5 ug/mL
N-Nitrosodiphenylamine	0.25-12.5 ug/mL
Azobenzene <sup>2</sup>	0.25-12.5 ug/mL
4-Bromophenyl phenyl ether	0.25-12.5 ug/mL
Hexachlorobenzene	0.25-12.5 ug/mL
Pentachlorophenol	0.25-12.5 ug/mL
Phenanthrene	0.05-10 ug/mL
Anthracene	0.05-10 ug/mL
Carbazole	0.05-10 ug/mL
Di-n-butyl phthalate	0.25-12.5 ug/mL
Fluoranthene	0.05-10 ug/mL
Benzidine	0.25-12.5 ug/mL
Pyrene	0.05-10 ug/mL
Butyl benzyl phthalate	0.25-12.5 ug/mL
3,3'-Dichlorobenzidine	0.25-12.5 ug/mL
Benzo(a)anthracene	0.05-10 ug/mL
Bis(2-ethylhexyl)phthalate	0.25-12.5 ug/mL
Chrysene	0.05-10 ug/mL
Di-n-octylphthalate	0.25-12.5 ug/mL
Benzo(b)fluoranthene	0.05-10 ug/mL
Benzo(k)fluoranthene	0.05-10 ug/mL
Benzo(a)pyrene	0.05-10 ug/mL
Indeno(1,2,3-cd)pyrene	0.05-10 ug/mL
Dibenz(a,h)anthracene	0.05-10 ug/mL
Benzo(g,h,i)perylene	0.05-10 ug/mL
Benzaldehyde	0.25-12.5 ug/mL
Caprolactam	0.25-12.5 ug/mL
1,1'-Biphenyl	0.25-12.5 ug/mL
Atrazine	0.25-12.5 ug/mL

<sup>1</sup> 2,2'-oxybis(1-chloropropane) was formerly known as bis(2-chloroisopropyl)ether

<sup>2</sup> Azobenzene is formed by decomposition of 1,2-diphenylhydrazine. If 1,2-diphenylhydrazine is requested, it will be analyzed as azobenzene.

**Note:** Nine calibration standards are prepared varying in concentration from 0.05 ug/mL to 12.5 ug/mL. A minimum of 5 calibration concentrations will be used for initial calibration. The concentration range of each analyte is listed in the table.

Table 13

## Calibration Ranges, Appendix IX, µg/mL

Semivolatiles	Calibration Range
2-Picoline	0.25-12.5 ug/mL
N-Nitrosomethylethylamine	0.25-12.5 ug/mL
Methyl methanesulfonate	0.25-12.5 ug/mL
N-Nitrosodiethylamine	0.25-12.5 ug/mL
Ethyl methanesulfonate	0.25-12.5 ug/mL
Pentachloroethane	0.25-12.5 ug/mL
Acetophenone	0.25-12.5 ug/mL
N-Nitrosopyrrolidine	0.25-12.5 ug/mL
N-Nitrosomorpholine	0.25-12.5 ug/mL
o-Toluidine	0.25-12.5 ug/mL
3-Methylphenol	0.25-12.5 ug/mL
N-Nitrosopiperidine	0.25-12.5 ug/mL
o,o,o-Triethyl-Phosphorothioate	0.25-12.5 ug/mL
a,a-Dimethyl-phenethylamine	0.25-12.5 ug/mL
2,6-Dichlorophenol	0.25-12.5 ug/mL
Hexachloropropene	0.25-12.5 ug/mL
p-Phenylenediamine	0.25-12.5 ug/mL
n-Nitrosodi-n-butylamine	0.25-12.5 ug/mL
Safrole	0.25-12.5 ug/mL
1,2,4,5-Tetrachlorobenzene	0.25-12.5 ug/mL
Isosafrole 1 + 2	0.25-12.5 ug/mL
1,4-Dinitrobenzene	0.25-12.5 ug/mL
1,4-Naphthoquinone	0.25-12.5 ug/mL
1,3-Dinitrobenzene	0.25-12.5 ug/mL
Pentachlorobenzene	0.25-12.5 ug/mL
1-Naphthylamine	0.25-12.5 ug/mL
2-Naphthylamine	0.25-12.5 ug/mL
2,3,4,6-Tetrachlorophenol	0.25-12.5 ug/mL
5-Nitro-o-toluidine	0.25-12.5 ug/mL
Thionazin	0.25-12.5 ug/mL
1,3,5-Trinitrobenzene	0.25-12.5 ug/mL
Sulfotepp	0.25-12.5 ug/mL
Phorate	0.25-12.5 ug/mL
Phenacetin	0.25-12.5 ug/mL
Diallate 1 + 2	0.25-12.5 ug/mL
Dimethoate	0.25-12.5 ug/mL
4-Aminobiphenyl	0.25-12.5 ug/mL
Pentachloronitrobenzene	0.25-12.5 ug/mL
Pronamide	0.25-12.5 ug/mL
Disulfoton	0.25-12.5 ug/mL
2-secbutyl-4,6-dinitrophenol (Dinoseb)	0.25-12.5 ug/mL
Methyl parathion	0.25-12.5 ug/mL
4-Nitroquinoline-1-oxide	0.25-12.5 ug/mL



Table 13

Calibration Ranges, Appendix IX,  $\mu\text{g/mL}$ 

Semivolatiles	Calibration Range
Parathion	0.25-12.5 $\mu\text{g/mL}$
Isodrin	0.25-12.5 $\mu\text{g/mL}$
Kepone	0.25-12.5 $\mu\text{g/mL}$
Famphur	0.25-12.5 $\mu\text{g/mL}$
Methapyriline	0.25-12.5 $\mu\text{g/mL}$
Aramite 1 and 2	0.25-12.5 $\mu\text{g/mL}$
p-(Dimethylamino)azobenzene	0.25-12.5 $\mu\text{g/mL}$
p-Chlorobenzilate	0.25-12.5 $\mu\text{g/mL}$
3,3'-Dimethylbenzidine	0.25-12.5 $\mu\text{g/mL}$
2-Acetylaminofluorene	0.25-12.5 $\mu\text{g/mL}$
Dibenz (a,j)acridine	0.25-12.5 $\mu\text{g/mL}$
7,12-Dimethylbenz(a)anthracene	0.25-12.5 $\mu\text{g/mL}$
3-Methylcholanthrene	0.25-12.5 $\mu\text{g/mL}$

**Note:** Nine calibrations standards are prepared varying in concentration from 0.05  $\mu\text{g/mL}$  to 12.5  $\mu\text{g/mL}$ . A minimum of 5 calibration concentrations will be used for initial calibration. The concentration range of each analyte is listed in the table.

Controlled Copy

Copy No: \_\_\_\_\_

Implementation Date 10/17/03

SOP No: CORP-GC-0001NC

Revision No: 5.7

Revision Date: 10/01/03

Page 1 of 84

## STL NORTH CANTON STANDARD OPERATING PROCEDURE

**TITLE: GAS CHROMATOGRAPHIC ANALYSIS BASED ON METHOD 8000B,  
8021B, 8081A, 8082, 8151A, 8310, 8141A, 8015B, and 615**

(SUPERSEDES: Revision 5.6 Dated 05/25/01)

Prepared by:	<u>Amey Duggal</u>	<u>10/9/03</u>
		Date
Approved by:	<u>Marcel Mol</u>	<u>10/9/03</u>
	Technology Specialist	Date
Approved by:	<u>Beth Lambert</u>	<u>10/17/03</u>
	Quality Assurance Manager	Date
Approved by:	<u>Deborah Dunn</u>	<u>10/16/03</u>
	Environmental Health and Safety	Date
Approved by:	<u>Carol Ann Felt</u>	<u>10/16/03</u>
	Laboratory Director	Date
Approved by:	<u>Mark Bunn</u>	<u>10/17/03</u>
	Technical Director	Date

## Proprietary Information Statement:

This documentation has been prepared by Severn Trent Laboratories (STL) solely for STL's own use and the use of STL's customers in evaluating its qualifications and capabilities in connection with a particular project. The user of this document agrees by its acceptance to return it to STL upon request and not to reproduce, copy, lend, or otherwise disclose its contents, directly or indirectly, and not to use it for any other purpose other than that for which it was specifically provided. The user also agrees that where consultants or other outside parties are involved in the evaluation process, access to these documents shall not be given to said parties unless those parties also specifically agree to these conditions. THIS DOCUMENT CONTAINS VALUABLE CONFIDENTIAL AND PROPRIETARY INFORMATION. DISCLOSURE, USE OR REPRODUCTION OF THESE MATERIALS WITHOUT THE WRITTEN AUTHORIZATION OF SEVERN TRENT LABORATORIES IS STRICTLY PROHIBITED. THIS UNPUBLISHED WORK BY STL IS PROTECTED BY STATE AND FEDERAL LAW OF THE UNITED STATES. IF PUBLICATION OF THIS WORK SHOULD OCCUR THE FOLLOWING NOTICE SHALL APPLY: ©2000 SEVERN TRENT LABORATORIES, INC. ALL RIGHTS RESERVED.

**TABLE OF CONTENTS**

1.	Scope and Application	4
2.	Summary of Method	4
3.	Definitions	4
4.	Interferences	4
5.	Safety	4
6.	Equipment and Supplies	5
7.	Reagents and Standards	5
8.	Sample Preservation and Storage	5
9.	Quality Control	6
10.	Calibration and Standardization	9
11.	Procedure	14
12.	Data Analysis and Calculations	16
13.	Method Performance	19
14.	Pollution Prevention	20
15.	Waste Management	20
16.	References	20
17.	Miscellaneous	20

**List of Appendices:**

Appendix A	Analysis of Volatile Organics based on Method 8021B
Appendix B	Analysis of Organochlorine Pesticides based on Method 8081A
Appendix C	Analysis of Organochlorine Pesticides and PCBs based on Method 8082
Appendix D	Analysis of Phenoxy Acid Herbicides based on Method 8151A
Appendix E	Analysis of Organophosphorus Pesticides Based on Method 8141A
Appendix F	Total Petroleum Hydrocarbons Based on 8015B
Appendix G	Analysis of Non-Halogenated Organic Compounds by Method 8015B, Direct Injection
Appendix H	Analysis of Phillips 66 Compounds

**GAS CHROMATOGRAPHIC ANALYSIS BASED ON  
METHOD 8000B, SW-846****SOP No. CORP-GC-0001NC**Revision No. 5.7Revision Date: 10/01/03Page 3 of 86

---

**List of Tables**

Table A1	Standard analyte lists for 8021B
Table A2-A4	Recommended conditions for method 8021B
Table A5	Surrogate and Internal standard concentrations for aqueous and low level soil samples, method 8021B
Table A6	Concentrations for LCS and MS/MSD compounds, low level soil and aqueous, method 8021B.
Table B1	Standard analyte and reporting limits, method 8081A
Table B2	Recommended conditions, method 8081A
Table B3	Calibration levels, method 8081A
Table B4	Column degradation evaluation mix, method 8081A
Table B5	LCS/matrix spike and surrogate levels, method 8081A
Table B6	LCS/matrix spike and surrogate levels for TCLP, method 8081A
Table B7	Suggested analytical sequence, method 8081A
Table B8	Performance limits, method 8081A
Table C1	Standard analyte list, method 8082
Table C2	Instrumental conditions, method 8082
Table C3	Calibration standards, method 8082
Table C4	LCS/Matrix spike and surrogate levels, method 8082
Table C5	Suggested Analytical Sequence, Method 8082
Table D1	Standard Analyte List, method 8151A
Table D2	Instrumental conditions, method 8151A
Table D3	LCS/Matrix spike and surrogate levels, method 8151B
Table D4	Performance limits, method 8151B
Table E1	Analyte and Reporting Limit List, Method 8141A
Table E2	Instrument Conditions, Method 8141A
Table E3	Initial Calibration, Method 8141A
Table E4	LCS/Matrix spike and surrogate compounds, Method 8141A
Table F1	Suggested GC Temperature Program for TPH Analysis
Table F2	Reporting Limits for TPH Analysis
Table G1	Non-Halogenated Organic Compounds Reporting Limits
Table G2	Non-Halogenated Organic Compounds Working Standards
Table H1	Phillips 66 Compounds Reporting Limits

**1. SCOPE AND APPLICATION**

This SOP describes procedures for analysis of organic analytes by Gas Chromatography (GC). The procedures are based on SW-846 methodology and are applicable for measurements made to comply with the Resource Conservation and Recovery Act (RCRA). Individual analytes and methods are described in the appendices. Appendix G describes procedures for the analysis of petroleum hydrocarbons by SW-846 8015B methodology. Appendix H includes criteria for the analysis of non-halogenated organic compounds by Method 8015B, Direct Injection. Appendix I describes the analysis of Phillips 66 analytes by Method 8015B.

**2. SUMMARY OF METHOD**

In general, semivolatile analytes in aqueous samples are prepared for analysis using continuous or separatory funnel liquid / liquid extraction or solid phase extraction (SOP # CORP-OP-0001NCNC). Solid samples are prepared using sonication, Soxhlet or pressurized fluid extraction (SOP # CORP-OP-0001NC). Volatile analytes are prepared for analysis using purge and trap methodology (Appendix A).

After the initial preparation step, the sample is introduced to the GC and concentrations of target analytes are measured by the detector response within a defined retention time window, relative to the response to standard concentrations. Internal or external standardization procedures are used as specified in the method appendices.

**3. DEFINITIONS**

Definitions of terms used in this SOP may be found in the glossary of the STL North Canton Laboratory Quality Manual (LQM), current version.

**4. INTERFERENCES**

Contamination by carryover can occur when a low concentration sample is analyzed after a high concentration sample. In addition, some purge and trap autosamplers are susceptible to port specific contamination. Co-elution of target analytes with non-targets can occur, resulting in false positives or biased high results. In particular, this is a problem with non-selective detectors such as the Flame Ionization Detector (FID). See the appendices for interferences specific to individual tests and suggested corrective actions.

**5. SAFETY**

- 5.1 Employees must abide by the policies and procedures in the Corporate Safety Manual, Radiation Safety Manual and this document.
- 5.2 Eye protection that prevents splash, laboratory coat, and appropriate gloves must be worn while samples, standards, solvents, and reagents are being handled. Disposable gloves that have become contaminated will be removed and discarded; other gloves will be cleaned immediately. Refer to the STL North Canton Chemical Hygiene plan for a complete description of personal protection equipment. Latex, Nitrile and vinyl gloves all provide adequate protection against the methanol used in this method.
- 5.3 The following is a list of the materials used in this method, which have a serious or significant hazard rating. **NOTE: This list does not include all materials used in the method. The table contains a summary of the primary hazards listed in the MSDS for each of the materials listed in the table.** A complete list of materials used in the method can be found in the reagents and materials section.

Employees must review the information in the MSDS for each material before using it for the first time or when there are major changes to the MSDS.

Material	Hazards	Exposure Limit (2)	Signs and symptoms of exposure
Methanol	Flammable Poison Irritant	200 ppm-TWA	A slight irritant to the mucous membranes. Toxic effects exerted upon nervous system, particularly the optic nerve. Symptoms of overexposure may include headache, drowsiness and dizziness. Methyl alcohol is a defatting agent and may cause skin to become dry and cracked. Skin absorption can occur; symptoms may parallel inhalation exposure. Irritant to the eyes.
Acetone	Flammable	1000 ppm-TWA	Inhalation of vapors irritates the respiratory tract. May cause coughing, dizziness, dullness, and headache.
Hexane	Flammable Irritant	500 ppm-TWA	Inhalation of vapors irritates the respiratory tract. Overexposure may cause lightheadedness, nausea, headache, and blurred vision. Vapors may cause irritation to the skin and eyes.
Methylene Chloride	Carcinogen Irritant	25 ppm-TWA 125 ppm-STEL	Causes irritation to respiratory tract. Has a strong narcotic effect with symptoms of mental confusion, light-headedness, fatigue, nausea, vomiting and headache. Causes irritation, redness and pain to the skin and eyes. Prolonged contact can cause burns. Liquid decreases the skin. May be absorbed through skin.
Sulfuric Acid	Corrosive Oxidizer Dehydrator Poison Carcinogen	1 Mg/M3-TWA	Inhalation produces damaging effects on the mucous membranes and upper respiratory tract. Symptoms may include irritation of the nose and throat, and labored breathing. Symptoms of redness, pain, and severe burn can occur. Contact can cause blurred vision, redness, pain and severe tissue burns. Can cause blindness.

8771368

**GAS CHROMATOGRAPHIC ANALYSIS BASED ON  
METHOD 8000B, SW-846**

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03Page 6 of 86

MTBE	Flammable Irritant	None established	Excessive exposure may cause irritation to the nose, throat, lungs, and respiratory tract. Central nervous system effects may include headache, dizziness, loss of balance and coordination.
THF	Flammable Irritant Possible Carcinogen	200 ppm TWA	Causes skin and eye irritation. Irritating to mucous membranes and upper respiratory tract.
Isooctane	Flammable Irritant	None Established	May cause eye, respiratory tract and skin irritation.
Hydrochloric Acid	Corrosive Poison	5 ppm-Ceiling	Inhalation of vapors can cause coughing, choking, inflammation of the nose, throat, and upper respiratory tract, and in severe cases, pulmonary edema, circulatory failure, and death. Can cause redness, pain, and severe skin burns. Vapors are irritating and may cause damage to the eyes. Contact may cause severe burns and permanent eye damage.
Sodium Bisulfate	Corrosive	None Listed	Contact may cause skin/eye burns. Inhalation can cause irritation of the respiratory tract with burning pain in the nose and throat, coughing, wheezing and shortness of breath. Causes chemical burns to the respiratory tract. May cause fatal spasms, inflammation or pulmonary/respiratory edema.
1 – Always add acid to water to prevent violent reactions.			
limit refers to the OSHA regulatory exposure limit.			

- 5.4. Opened containers of neat standards will be handled in a fume hood.
- 5.5. Sample extracts and standards, which are in a flammable solvent, shall be stored in an explosion-proof refrigerator.
- 5.6. When using hydrogen gas as a carrier, all precautions listed in the CSM shall be observed.
- 5.7. Standard preparation and dilution shall be performed inside an operating fume hood.
- 5.8. The gas chromatograph contains zones that have elevated temperatures. The analyst needs to be aware of the locations of those zones, and must cool them to room temperature prior to working on them.
- 5.9. There are areas of high voltage in both the gas chromatograph and the mass spectrometer. Depending on the type of work involved, either turn the power to the instrument off, or disconnect it from its source of power.

**6. EQUIPMENT AND SUPPLIES**

- 6.1. An analytical system complete with a gas chromatograph is required. A data system capable of measuring peak area and/or height is required. Recommended equipment and supplies for individual methods are listed in each method appendix.

**7. REAGENTS AND STANDARDS****7.1. Stock Standards**

Stock standards are purchased as certified solutions or prepared from pure solutions. Stock standards for method 8021B are stored at -10 to -20°C. Other stock standard solutions are stored as recommended by the manufacturer. All stock standards must be protected from light. Stock standard solutions should be brought to room temperature before using.

Semivolatile stock standard solutions must be replaced after one year. Stock standards of gases must be replaced at least every week, unless the acceptability of the standard is demonstrated (Less than 20% drift from the initial calibration is an acceptable demonstration). Other volatile stock standards must be replaced every 6 months or sooner if comparison with check standards prepared from an independent source indicates a problem.

- 7.1.1. Expiration times for all standards are measured from the time the standard is prepared or from the time that the standard ampoule is opened, if the standard is supplied in a sealed ampoule. If a vendor-supplied standard has an earlier expiration date then that date is used.

**7.2. Calibration Standards****7.2.1. Volatile Calibration Standards**

The procedure for preparation of volatile standards is given in Appendix A.

**7.2.2. Semivolatile Calibration Standards**

Semivolatile calibration standards are prepared as dilutions of the stock standards. Surrogates and internal standards are used as specified in the method appendices. Semivolatile calibration solutions must be refrigerated at  $\leq 6^{\circ}\text{C}$  and protected from light. The standards must be replaced at least every six months or sooner if comparison with check standards indicates a problem.

- 7.3. Gases for carrier and make-up: Hydrogen, Helium, Nitrogen, Argon/Methane.

**7.4. Quality control (QC) Standards**

QC standards (matrix spiking and LCS standards) are prepared and stored in the same way as calibration standards. They must be made from a stock independent from the calibration standards.

**8. SAMPLE PRESERVATION AND STORAGE**

Semivolatile extracts must be refrigerated at  $\leq 6^{\circ}\text{C}$  and analyzed within 40 days of the end of the extraction. Volatile sample storage conditions and holding times are given in Appendix A.

**9. QUALITY CONTROL****9.1. Initial Demonstration of Capability**

- 9.1.1. For the standard analyte list, the initial demonstration and method detection limit (MDL) studies described in Section 13 must be acceptable before analysis of samples may begin.



- 9.1.2. For non-standard analytes, a MDL study must be performed and calibration curve generated before analyzing any samples, unless lesser requirements are previously agreed to with the client. In any event the minimum initial demonstration required is analysis of an extracted standard at the reporting limit and a single point calibration.

9.2 Batch Definition

Batches are defined at the sample preparation stage. Batches should be kept together through the whole analytical process as far as possible, but it is not mandatory to analyze prepared extracts on the same instrument or in the same sequence. Refer to the STL North Canton QC Program document (QA-003) for further details of the batch definition.

9.2.1. Quality Control Batch

The batch is a set of up to 20 samples of the same matrix processed using the same procedures and reagents within the same time period. The Quality Control batch must contain a matrix spike / spike duplicate (MS/MSD), a Laboratory Control Sample (LCS), and a method blank. Laboratory generated QC samples (Blank, LCS, MS/MSD) do not count towards the maximum 20 samples in a batch. Field QC samples are included in the batch count. In some cases, at client request, the MS/MSD may be replaced with a matrix spike and sample duplicate. If insufficient sample is available for an MS/MSD a LCSD may be substituted.

9.3. Control Limits

In-house historical control limits must be determined for surrogates, matrix spikes, and laboratory control samples (LCS). These limits must be determined at least annually. The recovery limits are mean recovery  $\pm 3$  standard deviations, unless that limit is tighter than the calibration criteria, in which case limits may be widened. Refer to policy QA-003 for more details.

- 9.3.1. These limits do not apply to dilutions (except for tests without a separate extraction), but surrogate and matrix spike recoveries will be reported unless the dilution is more than 5X.

- 9.3.2. All surrogate, LCS, and MS recoveries (except for dilutions) must be entered into QuantIMS (when available) or other database so that accurate historical control limits can be generated. For tests without a separate extraction, surrogates and matrix spikes will be reported for all dilutions.

- 9.3.3. Refer to the QC Program document (QA-003) for further details of control limits.

9.4. Surrogates

All methods must use surrogates to the extent possible. Surrogate recoveries in samples and QC samples must be assessed to ensure that recoveries are within established limits. If any surrogates are outside limits, the following corrective actions must take place (except for dilutions):

- Check all calculations for error.
- Ensure that instrument performance is acceptable.
- Recalculate the data and/or reanalyze the extract if either of the above checks reveal a problem.
- Reprepare and reanalyze the sample or flag the data as "Estimated Concentration" if neither of the above resolves the problem. Repreparation is not necessary if there is obvious chromatographic interference.
- The decision to reanalyze or flag the data should be made in consultation with the client. It is only necessary to reprepare / reanalyze a sample once to demonstrate that poor surrogate recovery is due to matrix effect, unless the analyst believes that the repeated out of control results are not due to matrix effect.

**GAS CHROMATOGRAPHIC ANALYSIS BASED ON  
METHOD 8000B, SW-846**

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03Page 9 of 86

---

- 9.4.1. If dual column analysis is used the choice of which result to report is made in the same way as for samples (Section 12.1.2) unless one column is out of control, in which case the in-control result is reported.
- 9.4.2. If the surrogates are out of control for the sample, matrix spike, and matrix spike duplicate, then matrix effect has been demonstrated for that sample and reparation is not necessary. If the sample is out of control and the MS and/or MSD is in control, then reparation or flagging of the data is required.
- 9.4.3. Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective actions.

**9.5. Method Blanks**

For each batch of samples, analyze a method blank. The method blank consists of reagent water for aqueous semivolatile samples, and sodium sulfate for semivolatile soils tests (Refer to SOP No. CORP-OP-0001NCNC for details). For low level volatiles, the method blank consists of reagent water. For medium level volatiles, the method blank consists of methanol as described in Appendix A. Surrogates are added and the method blank is carried through the entire analytical procedure. The method blank must not contain any analyte of interest at or above the reporting limit (except common laboratory contaminants, see below) or at or above 5% of the measured concentration of that analyte in the associated samples, whichever is higher.

If the analyte is a common laboratory contaminant (methylene chloride, acetone, 2-butanone, phthalate esters) the data may be reported with qualifiers if the concentration of the analyte is less than five times the reporting limit. Such action must be taken in consultation with the client.

Re-extraction and reanalysis of samples associated with an unacceptable method blank is required when reportable concentrations are determined in the samples.

If there is no target analyte greater than the RL in the samples associated with an unacceptable method blank, the data may be reported with qualifiers. Such action should be taken in consultation with the client.

- 9.5.1. Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective actions.

**9.6. Instrument Blanks**

- 9.6.1. An instrument blank must be analyzed during any 12-hour period of analysis that does not contain a method blank.
- 9.6.2. An instrument blank consists of the appropriate solvent with internal standards added. If internal standards are not used the surrogates should be added.
- 9.6.3. Control criteria are the same as for the method blank, except that only reanalysis of affected samples would be required, not re-extraction.

**9.7. Laboratory Control Samples (LCS)**

For each batch of samples, analyze a LCS. The LCS contains a representative subset of the analytes of interest, and must contain the same analytes as the matrix spike. The LCS may also contain the full set of analytes with a subset of control analytes. If any control analyte or surrogate is outside established control limits, the system is out of control and corrective action must occur. Corrective action will normally be reparation and reanalysis of the batch; however, if the matrix spike and matrix spike duplicate are within limits; the batch may be acceptable

- 9.7.1. Refer to the STL North Canton QC Program document (QA-003) for further details of the corrective action.
- 9.7.2. If dual column analysis is used the choice of which result to report is made in the same way as for samples (Section 12.1.2) unless one column is out of control, in which case the in control result is reported.
- 9.7.3. LCS compound lists are included in the appendices.

9.8. Matrix Spikes

For each QC batch, analyze a matrix spike and matrix spike duplicate. Spiking compounds and levels are given in the appendices. Compare the percent recovery and relative percent difference (RPD) to those in the laboratory specific historically generated limits.

- If any individual recovery or RPD falls outside the acceptable range, corrective action must occur. The initial corrective action will be to check the recovery of that analyte in the Laboratory Control Sample (LCS). Generally, if the recovery of the analyte in the LCS is within limits, then the laboratory operation is in control and analysis may proceed.
- If the recovery for any component is outside QC limits for both the Matrix spike / spike duplicate and the LCS, the laboratory is out of control and corrective action must be taken. Corrective action will normally include reparation and reanalysis of the batch.
- If a MS/MSD is not possible due to limited sample, then a LCS duplicate should be analyzed.
- The matrix spike / duplicate must be analyzed at the same dilution as the unspiked sample, unless the matrix spike components would then be above the calibration range.

- 9.8.1. If dual column analysis is used the choice of which result to report is made in the same way as for samples (Section 12.1.2) unless one column is out of control, in which case the in control result is reported

9.9. Quality Assurance Summaries

Certain clients may require specific project or program QC that may supersede these method requirements. Quality Assurance Summaries should be developed to address these requirements.

9.10. STL North Canton QC Program

Further details of QC and corrective action guidelines are presented in the STL QC Program document (QA-003) Refer to this document if in doubt regarding corrective actions.

## 10. CALIBRATION AND STANDARDIZATION

Internal or external calibration may be used. Internal calibration is recommended unless the sample matrix is likely to interfere with the quantitation of the internal standard. In either event prepare standards containing each analyte of interest at a minimum of five concentration levels. The low level standard should be at or below the reporting limit. The other standards define the working range of the detector. Recommended calibration levels are given in the appendices.

- 10.1. A new calibration curve must be generated after major changes to the system or when the continuing calibration criteria cannot be met. Major changes include new columns, changing PID lamps or FID jets or replacing the ECD detector. A new calibration is not required after clipping the column, replacing the septum or syringe, or other minor maintenance.

- 10.2. With the exception of 10.3 below, it is NOT acceptable to remove points from a calibration curve for the purpose of meeting criteria, unless the points are the highest or lowest on the curve AND the reporting limit and/or linear range is adjusted accordingly. In any event, at least 5 points must be included in the calibration curve. Quadratic (second order) calibrations require at least six points. Third order calibrations require at least seven points.
- 10.3. A level may be removed from the calibration if the reason can be clearly documented, for example a broken vial or no purge run. A minimum of five levels must remain in the calibration. The documentation must be retained with the initial calibration. Alternatively, if the analyst believes that a point on the curve is inaccurate, the point may be reanalyzed and the reanalysis used for the calibration. All initial calibration points must be analyzed without any changes to instrument conditions, and all points must be analyzed within 24 hours.
- 10.4. External standard calibration
- Quantitation by the external standard method assumes a proportional relationship between the calibration run and the analyte in the sample. To use this approach, introduce each calibration standard into the GC using the technique that will be used for samples. The ratio of the peak height or area response to the mass or concentration injected may be used to prepare a calibration curve

$$\text{Calibration Factor (CF)} = \frac{\text{Area or Height of Peak}}{\text{Mass Injected (ng)}}$$

Some data systems may use the inverse of this formula. This is acceptable so long as the same formula is used for standards and samples. It is also possible to use the concentration of the standard rather than the mass injected. (This would require changes in the equations used to calculate the sample concentrations). Use of peak area or height must be consistent. However, if matrix interferences would make quantitation using peak area inaccurate for a particular sample, then peak height may be used as a substitute.

- 10.5. Internal standard calibration
- 10.5.1. The internal standard approach assumes that variations in instrument sensitivity, amount injected etc. can be corrected by determining the ratio of the response of the analyte to the response of an internal standard that has been added to the extract. To use this approach, select one or more internal standard(s) that are similar in analytical behavior to the compounds of interest. Recommended internal standards are given in the appendices. The analyst must demonstrate that the measurement of the internal standard is not affected by method or matrix interferences. If the sample matrix interferes with quantitation of the internal standard, then the external standard approach must be used instead. In this event use the response factors from the previous continuing calibration to quantitate the analytes in the sample with the interference (applies only to the sample with the interference).
- 10.5.2. Introduce each calibration standard into the GC using the technique that will be used for samples. Response factors (RF) for each compound are calculated as follows.

$$RF = \frac{A_s \times C_{is}}{A_{is} \times C_s}$$

Where:

$A_s$  = Response for the analyte to be measured

$A_{is}$  = Response for the internal standard

$C_{is}$  = Concentration of internal standard

$C_s$  = Concentration of the analyte to be determined in the standard

#### 10.6. Calibration curve fits

Average response factor, linear regression, or quadratic curves may be used to fit the data. Average response factor may be used if the average % RSD of the response factors or calibration factors of all the analytes in the calibration standard taken together is  $\leq 20\%$ . The average %RSD is calculated by summing the RSD value for each analyte and dividing by the total number of analytes.

10.6.1. In general, for environmental analysis, average response factors are the most appropriate calibration model. Linear or curved regression fits should only be used if the analyst has reason to believe that the average RF model does not fit the normal concentration/response behavior of the detector.

#### 10.6.2. Average response factor

The average response factor may be used if the average percent relative standard deviation (%RSD) of all the response factors taken together is  $\leq 20\%$ .

The equation for average response factor is.

$$\text{Average response factor} = \overline{RF} = \frac{\sum_{i=1}^n RF_i}{n}$$

Where:  $n$  = Number of calibration levels

$$\sum_{i=1}^n RF_i = \text{Sum of response factors for each calibration level}$$

## 10.6.3. Linear regression

The linear fit uses the following functions:

## 10.6.3.1. External Standard

$$y = ax + b$$

or

$$x = \frac{(y - b)}{a}$$

Where:  $y$  = Instrument response

$x$  = Concentration

$a$  = Slope

$b$  = Intercept

## 10.6.3.2. Internal Standard

$$C_s = \frac{\left[ \frac{A_s C_{is}}{A_{is}} - b \right]}{a}$$

Where:  $C_s$  = Concentration in the sample

$A_s$  = Area of target peak in the sample

$A_{is}$  = Area of internal standard in the sample

$C_{is}$  = Concentration of the internal standard

## 10.6.4. Quadratic curve

The quadratic curve uses the following functions:

## 10.6.4.1. External standard

$$y = ax + cx^2 + b$$

Where  $c$  is the curvature

## 10.6.4.2. Internal Standard

$$y = a \left( \frac{A_s \times C_s}{A_{is}} \right) + c \left( \frac{A_s \times C_s}{A_{is}} \right)^2 + b$$

## 10.7. Evaluation of calibration curves

- 10.7.1. The percent relative standard error (%RSE) from the calibration curve is used to evaluate the initial calibration. This provides a measure of how much error is associated with using the calibration curve for quantitation.
- 10.7.2. The least squares regression line is calculated and used to calculate the predicted concentration for each level. The percent relative standard error is calculated as follows:

$$\% RSE = 100\% \times \sqrt{\frac{\sum_{i=1}^N \left[ \frac{C_i - PC_i}{C_i} \right]^2}{(N - P)}}$$

Where

$N$  = Number of points in the curve

$P$  = Number of parameters in the curve (= 1 for average response factor, 2 for linear, 3 for quadratic)

$C_i$  = True concentration for level  $i$

$PC_i$  = Predicted concentration for level  $i$

Note that when average response factors are used, %RSE is equivalent to %RSD.

## 10.8. The following requirements must be met for any calibration to be used:

- Response must increase with increasing concentration.
- If a curve is used, the intercept of the curve at zero response must be less than  $\pm$  the reporting limit for the analyte.
- The average Relative Standard Error (RSD for average response factors) of the calibration points from the curve used must be  $\leq 20\%$ .
- Some data systems will not measure the %RSE from a linear or quadratic fit. For the linear case, the correlation coefficient may be used as an alternative to the %RSE, and must be greater than or equal to 0.990. For the quadratic case the Coefficient of Determination may be used, and must be greater or equal to 0.990.

**Note.** The Relative Standard Error (RSE) is superior to the Correlation Coefficient ( $r$ ) and Coefficient of Determination ( $r^2$ ) for testing the fit of a set of calibration points to a line. The lower points on a curve have little effect on  $r$ . As a result a curve may have a very good correlation coefficient ( $>0.995$ ), while also having  $> 100\%$  error at the low point

## 10.9. Weighting of data points

- 10.9.1. In linear and quadratic calibration fits, the points at the lower end of the calibration curve have less absolute variance than points at the high concentration end of the curve. This can cause severe errors in quantitation at the low end of the calibration. However, in environmental analysis, accuracy at the low end of the curve is very important. For this reason it is preferable to increase the weighting of the lower concentration points.  $1/\text{Concentration}^2$  weighting (often called  $1/X^2$  weighting) will improve accuracy at the low end of the curve and should be used if the data system has this capability.

- 10.10. Non-standard analytes are sometimes requested. For these analytes, it may be acceptable to analyze a single standard at the reporting limit with each continuing calibration rather than a five point initial calibration. This action must be with client approval. If the analyte is detected in any of the samples, a five point initial calibration must be generated and the sample(s) reanalyzed for quantitation.

## 10.11. Calibration Verification

## 10.11.1. 12 hour Calibration

The working calibration curve or RF must be verified by the analysis of a mid point calibration standard at the beginning, after every 12 hours, and at the end of the analysis sequence. The center of each retention time window is updated with each 12-hour calibration or calibration verification.

## 10.11.2. Calibration Verification

It may be appropriate to analyze a mid point standard more frequently than every 12 hours. If these calibration verification standards are analyzed, requirements are the same as the 12-hour calibration with the exception that retention times are not updated.

- 10.11.3. Any individual compounds with %D < 15% meet the calibration criteria. The calibration verification is also acceptable if the average of the %D for all the analytes is < 15%. This average is calculated by summing the entire absolute %D results in the calibration (including surrogates) and dividing by the number of analytes. Any analyte that is reportable as found must have a % difference of < 15% in the calibration verification or 12 hour calibration, on the column used for quantitation. Refer to section 12.1.2 for which result to report.

- 10.11.4. It is not necessary to run a calibration verification standard at the beginning of the sequence if samples are analyzed immediately after the completion of the initial calibration.

- 10.11.5. Samples quantitated by external standard methods must be bracketed by calibration verification standards that meet the criteria listed above. The bracketing standards on the column used for calibration must meet the same criteria as the opening standards. Bracketing is not necessary for internal standard methods.

- 10.11.6. If the analyst notes that a CCV has failed and can document the reason for failure (e.g. no purge, broken vial, carryover from the previous sample etc.) then a second CCV may be analyzed without any adjustments to the instrument. If this CCV meets criteria then the preceding samples have been successfully bracketed. If adjustments to the instrument are performed before the repeat CCV then the proceeding samples have not been successfully bracketed but analysis may continue.

- 10.11.7. In general, it is not advisable to analyze repeat CCVs on unattended runs. If repeat CCVs are analyzed then the first will serve as the bracketing standard for the preceding samples and the last will serve as the CCV for the following samples.



10.11.8. If highly contaminated samples are expected it is acceptable to analyze blanks or primers at any point in the run.

10.11.9 % Difference calculation

% Difference for internal and external methods is calculated as follows

Internal Standard:

External standard:

$$\%D = \frac{RF_c - \overline{RF}}{\overline{RF}} \times 100$$

$$\%D = \frac{CF_c - \overline{CF}}{\overline{CF}} \times 100$$

Where  $RF_c$  and  $CF_c$  are the response and calibration factors  
from the continuing calibration

$\overline{RF}$  and  $\overline{CF}$  are the average response and calibration factors  
from the initial calibration

10.11.10 % Drift calculation

% Drift is used for comparing the continuing calibration to a linear or quadratic curve. The criteria for % drift are the same as for % difference

$$\% \text{ Drift} = \frac{\text{Calculated Conc} - \text{Theoretical Conc}}{\text{Theoretical Conc}} \times 100\%$$

10.11.11. Corrective Actions for Continuing Calibration

If the overall average %D of all analytes is greater than  $\pm 15\%$  corrective action must be taken. This may include clipping the column, changing the liner or other minor instrument adjustments, followed by reanalyzing the standard. If the overall average %D still varies by more than  $\pm 15\%$ , a new calibration curve must be prepared.

10.11.12. Corrective Action for Samples

For internal standard methods, any samples injected after a standard not meeting the calibration criteria must be re-injected.

For external standard methods, any samples injected after the last good continuing calibration standard must be re-injected.

If the average %D for all the analytes in the calibration is over 15%, but all of the analytes requested for a particular sample have  $\%D \leq 15\%$ , then the analysis is acceptable for that sample.

## 11. PROCEDURE

### 11.1. Extraction

Extraction procedures are referenced in the appendices.

### 11.2. Cleanup

Cleanup procedures are referenced in the appendices.

## 11.3. Gas Chromatography

Chromatographic conditions for individual methods are presented in the appendices.

## 11.4. Sample Introduction

In general, volatile analytes are introduced using purge and trap as described in Appendix A.

Semivolatile analytes are introduced by direct injection of the extract. Samples, standards, and QC must be introduced using the same procedure.

## 11.5. Analytical Sequence

An analytical sequence starts with an initial calibration or a daily calibration. Refer to the individual method appendices for method specific details of daily calibrations and analytical sequences

11.5.1. The daily calibration includes analysis of standards containing all single response analytes and updating the retention time windows.

11.5.2. If there is a break in the analytical sequence of greater than 12 hours, a new analytical sequence must be started with a daily calibration.

## 11.6. Retention Time Windows

11.6.1. Retention time windows must be determined for all analytes. Make an injection of all analytes of interest each day over a three-day period. Calculate the standard deviation of the three retention times for each analyte (relative retention times may also be used). For multi-response analytes (e.g., Aroclors) use the retention time of major peaks. Plus or minus three times the standard deviation of the retention times of each analyte defines the retention time window.

11.6.2. The center of the retention time window is the retention time from the last of the three standards. The centers of the windows are updated with the mid- point of the initial calibration and each 12-hour calibration. The widths of the windows will remain the same until new windows are generated following the installation of a new column.

11.6.3. If the retention time window as calculated above is less than  $\pm 0.05$  minutes, use  $\pm 0.05$  minutes as the retention time window. This allows for slight variations in retention times caused by sample matrix.

11.6.4. The laboratory must calculate new retention time windows each time a new column is installed. The new windows must be generated within one week of the installation of the new column. Until these standards have been run on the new column, the retention time windows from the old column may be used, updated with the retention times from the new initial calibration

## 11.6.5. Corrective Action for Retention Times

The retention times of all compounds in the 12 hour calibration or calibration verification standard must be within the retention time window. If this condition is not met, all samples analyzed after the last compliant standard must be reanalyzed unless the following conditions are met for any compound that elutes outside the retention time window:

The retention time of that compound in the standard must be within a retention time range equal to twice the original window.

No peak that would be reportable may be present on the sample chromatogram within an elution time range equal to three times the original retention time window.

## 11.7. Daily Retention Time Windows

The center of the retention time windows determined in Section 11.6 are adjusted to the retention time of each analyte as determined in the 12 hour calibration standards or continuing calibration verification standards. (See the method 8081A and 8082 appendices for exceptions for multi-response components.) The retention time windows must be updated at the beginning of each analytical sequence and with each 12-hour calibration or continuing calibration verification.

11.8. Percent Moisture

Analytical results may be reported as dry or wet weight, as required by the client. Percent moisture must be determined if results will be reported as dry weight. Refer to SOP CORP-OP-0001NC for determination of percent moisture.

11.9. Procedural Variations

Procedural variations are allowed only if deemed necessary in the professional judgment of the supervisor to accommodate variation in sample matrix, radioactivity, chemistry, sample size, or other parameters. Any variation in procedure shall be completely documented using a Nonconformance Memo and approved by a supervisor and QA/QC manager. If contractually required, the client shall be notified. The Nonconformance Memo shall be filed in the project file. The nonconformance is also addressed in the case narrative. Any unauthorized deviations from this procedure must also be documented as a nonconformance, with a cause and corrective action described.

## 12. DATA ANALYSIS AND CALCULATIONS

### 12.1 Qualitative Identification

- 12.1.1 Tentative identification occurs when a peak is found within the retention time window for an analyte, at a concentration above the reporting limit, or above the MDL if J flags are required. Normally confirmation is required on a second column, but if the detector is sufficiently specific or if the sample matrix is well enough defined, single column analysis may be adequate. In some cases GC/MS confirmation may be required. Client specific requirements may also define the need for second column confirmation and / or GC/MS confirmation. Refer to the appendices for test specific requirements for confirmation. Identification is confirmed if a peak is also present in the retention time window for that analyte on the confirmatory column, at a concentration greater than the reporting limit (MDL if J flag confirmation required).

#### 12.1.2. Dual column quantitation

For confirmed results, two approaches are available to the analyst.

A) The primary column approach

Or

B) The better result approach

Both are acceptable to avoid the reporting of erroneous or unconfirmed data.

##### 12.1.2.1. Primary column approach:

The result from the primary column is normally reported. The result from the secondary column is reported if any of the following three bulleted possibilities are true.

- There is obvious chromatographic interference on the primary column
- The result on the primary column is > 40% greater than the result on the secondary column
- Continuing or bracketing standard fails on the primary column but is acceptable on the secondary column. (If the primary column result is > 40% higher than the secondary and the primary column calibration fails, then the sample must be evaluated for reanalysis.)

## 12.1.2.2. Better result approach

The lower of the two results is normally reported. The lower result is considered better because the higher result is generally higher because of chromatographic interference. The higher result is reported if any of the following two bulleted possibilities are true.

- There is obvious chromatographic interference on the column with the lower result
- The continuing or bracketing calibration on the column with the lower result fails. (If the higher result is > 40% higher and the calibration on the column with the lower result fails, then the sample must be evaluated for reanalysis.)

- 12.1.3. If the Relative percent difference (RPD) between the response on the two columns is greater than 40%, or if the opinion of an experienced analyst is that the complexity of the matrix is resulting in false positives, the confirmation is suspect and the results are qualified. RPD is calculated using the following formula:

$$RPD = \frac{R_1 - R_2}{\frac{1}{2}(R_1 + R_2)}$$

Where R=Result

## 12.1.4. Multi-response Analytes

For multi-response analytes, the analyst should use the retention time window, but should rely primarily on pattern recognition. The pattern of peaks will normally serve as confirmation.

- 12.1.5. The experience of the analyst should weigh heavily in the interpretation of the chromatogram. For example, sample matrix or laboratory temperature fluctuation may result in variation of retention times.

## 12.2. Calibration Range

If concentrations of any analytes exceed the working range as defined by the calibration standards, then the sample must be diluted and reanalyzed. Dilutions should target the most concentrated analyte in the upper half (over 50% of the high level standard) of the calibration range. It may be necessary to dilute samples due to matrix.

## 12.3. Dilutions

Samples may be screened to determine the appropriate dilution for the initial run. If the initial diluted run has no hits or hits below 20% of the calibration range and the matrix allows for analysis at a lesser dilution, then the sample must be reanalyzed at a dilution targeted to bring the largest hit above 50% of the calibration range.

## 12.3.1. Guidance for Dilutions Due to Matrix

If the sample is initially run at a dilution and only minor matrix peaks are, then the sample should be reanalyzed at a more concentrated dilution. Analyst judgement is required to determine the most concentrated dilution that will not result in instrument contamination.

## 12.3.2. Reporting Dilutions

The most concentrated dilution with no target compounds above the calibration range will be reported. Other dilutions may be reported at client request, if the lower dilutions will not cause detector saturation, column overload, or carryover. Analyst judgement and client site history will factors in the reporting of dual dilutions.

## 12.4. Interferences

If peak detection is prevented by interferences, further cleanup should be attempted. If no further cleanup is reasonable, then elevation of reporting levels and/or lack of positive identification must be addressed in the case narrative.

## 12.5. Internal Standard Criteria for Continuing Calibration

If internal standard calibration is used, then the internal standard response in a continuing calibration standard must be within 50 to 150% of the response in the mid level of the initial calibration.

## 12.6. Calculations

Capabilities of individual data systems may require the use of different formulas than those presented here. When this is the case, the calculations used must be shown to be equivalent and must be documented in an appendix attached to this document.

## 12.6.1. External Standard Calculations

## 12.6.1.1. Aqueous samples

$$\text{Concentration (mg / L)} = \frac{(A_s \times V_i \times D_f)}{(CF \times V_t \times V_s)}$$

Where:

$A_s$  = Response for the analyte in the sample

$V_i$  = Volume of extract injected,  $\mu\text{L}$

$D_f$  = Dilution factor

$V_t$  = Volume of total extract,  $\mu\text{L}$

$V_s$  = Volume of sample extracted or purged, mL

$CF$  = Calibration factor, area or height/ng, Section 10.1

## 12.6.1.2. Non-aqueous Samples

$$\text{Concentration (mg / kg)} = \frac{(A_s \times V_i \times D_f)}{(CF \times V_t \times W \times D)}$$

Where:

$W$  = Weight of sample extracted or purged, g

$$D = \frac{100 - \% \text{Moisture}}{100} \quad (D = 1 \text{ if wet weight is required})$$

## 12.6.2. Internal Standard Calculations

## 12.6.2.1. Aqueous Samples

$$\text{Concentration (mg / L)} = \frac{(A_s \times C_b \times D_f)}{(A_{is} \times RF \times V_s)}$$

Where:

$C_{is}$  = Amount of internal standard added, ng

$A_{is}$  = Response of the internal standard

$RF$  = Response factor for analyte

12.6.2.2. Non-aqueous Samples

$$\text{Concentration (mg / kg)} = \frac{(A_s \times C_{is} \times D_f)}{(A_{is} \times RF \times W \times D)}$$

12.6.3. Surrogate Recovery

Concentrations of surrogate compounds are calculated using the same equations as for the target compounds. The response factor from the initial calibration is used. Surrogate recovery is calculated using the following equation:

$$\% \text{ Recovery} = \frac{\text{Concentration (or amount) found}}{\text{Concentration (or amount) spiked}} \times 100$$

13. METHOD PERFORMANCE

13.1. Method Detection Limit

Each laboratory must generate a valid method detection limit for each analyte of interest. The MDL must be below the reporting limit for each analyte. The procedure for determination of the method detection limit is given in 40 CFR Part 136, Appendix B, and further defined in QA Policy #: QA-005.

13.2. Initial Demonstration

Each laboratory must make a one time initial demonstration of capability for each individual method. Demonstration of capability for both soils and water matrices is required. This requires analysis of QC check samples containing all of the standard analytes for the method. For some tests it may be necessary to use more than one QC check mix to cover all analytes of interest.

13.2.1. Four aliquots of the QC check sample are analyzed using the same procedures used to analyze samples, including sample preparation. The concentration of the QC check sample should be equivalent to a mid-level calibration.

13.2.2. Calculate the average recovery and standard deviation of the recovery for each analyte of interest. Compare these results with the acceptance criteria given in each appendix.

13.2.3. If any analyte does not meet the acceptance criteria, the test must be repeated. Only those analytes that did not meet criteria in the first test need to be evaluated. Repeated failure for any analyte indicates the need for the laboratory to evaluate the analytical procedure and take corrective action.

## 13.3. Training Qualification

The group/team leader has the responsibility to ensure that an analyst who has been properly trained in its use and has the required experience performs this procedure.

## 14. POLLUTION PREVENTION

This method does not contain any specific modifications that serve to minimize or prevent pollution.

## 15. WASTE MANAGEMENT

15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."

## 15.2. Waste Streams Produced by the Method.

15.2.1. The following waste streams are produced when this method is carried out.

15.2.1.1 **Vials containing sample extracts:** These vials are placed in the vial waste located in the GC/MS laboratory.

15.2.1.2 **Tubes containing sample extracts for TPH, Pesticides, PCBs and Herbicides.** These capped tubes are placed in the PCB/flammable waste located the GC prep laboratory

15.2.1.3 **Samples, standards, and all extraction materials contaminated with high levels (>50ppm) of PCB's must be segregated into their own waste stream.** PCB wastes are collected in one of three waste streams, solid PCB, liquid PCB and PCB vial waste. PCB containing samples are located through a LIMS query and disposed of as PCB containing.

15.2.1.4. **Extracted solid samples contaminated with methylene chloride/acetone or acetone/hexane.** These materials are disposed of in the solid waste and debris in a red container located in the extractions lab.

15.2.1.5. **Discarded samples** These samples are collected in the solid debris drum

## 16. REFERENCES

Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW846, 3rd Edition, Final Update III, December 1996, and Section 8000B

## 17. MISCELLANEOUS

## 17.1 Modifications from Reference Method

17.1.1. Chapter 1 of SW-846 states that the method blank should not contain any analyte of interest at or above the Method Detection Limit. This SOP states that the Method Blank must not contain any analyte of interest at or above the reporting limit. Common lab contaminants are allowed to be up to 5 times the reporting limit in the blank following consultation with the client.

**GAS CHROMATOGRAPHIC ANALYSIS BASED ON  
METHOD 8000B, SW-846**

**SOP No. CORP-GC-0001NC**

Revision No. 5.7

Revision Date: 10/01/03

Page 23 of 86

---

17.2. Modifications from Previous Revision



The calibration criteria in section 10.11 have been rewritten to improve consistency with SW-846 and to improve clarity.

17.3. Facility Specific SOPs


Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none.



## 17.4. Flow Diagrams

17.4.1. Initial demonstration and MDL<sup>1</sup>

<sup>1</sup> This flow diagram is for guidance and cannot cover all eventualities. Consult the SOP text and a supervisor if in doubt.

17.4.2. Sample Analysis<sup>1</sup>

<sup>1</sup> This flow diagram is for guidance and cannot cover all eventualities. Consult the SOP text and a supervisor if in doubt.

## APPENDIX A

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03Page A 1 of A17

---

**1. SCOPE AND APPLICATION**

- 1.1. This method describes sample preparation and extraction for the analysis of volatile organics by a purge and trap procedure, following method 8021B. However, where required by a client QAPP this section may also be used to analyze aromatic volatiles by discontinued methods 8020A and 8010B. All requirements of the 8000B section of this SOP must be met except when superseded by this Appendix. Refer to Table A-1 for the individual analytes normally determined by these procedures.
- 1.2. Compounds within the scope of this method have boiling points below 200°C and are soluble or slightly soluble in water. Classes of compounds best suited to purge-and-trap analysis include low molecular weight halogenated hydrocarbons, aromatics, ketones, nitriles, acetates, acrylates, ethers, and sulfides.
- 1.3. Water samples and soils samples with low levels of contamination may be analyzed directly by purge-and-trap extraction and gas chromatography. Higher concentrations of these analytes in soil may be determined by the medium level methanol extraction procedure.
- 1.4. This method also describes the preparation of water-miscible liquids, non-water-miscible liquids, solids, wastes, and soils/sediments for analysis by the purge-and-trap procedure
- 1.5. The associated LIMS method code is QR

**2. SUMMARY OF METHOD**

- 2.1. An inert gas is bubbled through the sample at ambient temperature or at 40°C (40°C required for low-level soils), and the volatile components are transferred from the aqueous phase to the vapor phase. The vapor is swept through a sorbent column where the volatile components are adsorbed. After purging is completed, the sorbent column is heated and back-flushed with inert gas to desorb the components onto a gas chromatographic column. Analytes are detected using a photoionization Detector, an electrolytic conductivity detector or a combination of both
- 2.2. For soil samples, a portion of the sample is dispersed in methanol to dissolve the volatile organic constituents. A portion of the methanolic solution is combined with water. It is then analyzed by purge-and-trap GC following the normal water method. If very low detection limits are needed for soil samples then direct purge using sodium bisulfate preservation may be necessary.

**3. DEFINITIONS**

Refer to the STL North Canton Laboratory Quality Manual (LQM), current version, for definitions of terms used in this SOP.

**4. INTERFERENCES**

- 4.1. Refer to section 4 of the method 8000B part of this SOP for general information on chromatographic interferences.
- 4.2. Impurities in the purge gas, and from organic compounds out-gassing from the plumbing ahead of the trap, account for the majority of contamination problems. The analytical system must be demonstrated to be free from contamination under the conditions of the analysis by running laboratory reagent blanks. The use of non-TFE plastic tubing, non-TFE thread sealants, or flow controllers with rubber components in the purging device should be avoided.
- 4.3. Samples can be contaminated by diffusion of volatile organics (particularly methylene chloride and fluorocarbons) through the septum seal of the sample vial during shipment and storage. A trip blank

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A2 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

prepared from organic-free reagent water and carried through sampling and handling protocols serves as a check on such contamination.

- 4.4. Contamination by carryover can occur whenever high-concentration and low-concentration samples are analyzed sequentially. Whenever an unusually concentrated sample is analyzed, it should be followed by an analysis of organic-free reagent water to check for cross-contamination. The trap and other parts of the system are subject to contamination. Therefore, frequent bake-out and purging of the system may be required.
- 4.5. When utilizing an autosampler system, which has multiple ports for sample analysis, it is likely that only a single stage or port may be contaminated by a highly concentrated sample. If a port is suspect, a water blank should be analyzed to verify lack of contamination. If the water blank and subsequent blanks on that port show contamination consistent with the concentrated sample, further maintenance is required. This may include replacing or cleaning the multi-port valve, transfer lines, etc.
- 4.6. A holding blank is kept in the sample refrigerator. This is analyzed and replaced every 14 days. If the holding blank does not meet the method blank criteria, the source of contamination must be found and corrected. Evaluation of all samples analyzed in the 14-day period prior to the analysis of the contaminated holding blank is required.
- 4.7. Acidification of samples may result in hydrolysis of 2-chloroethyl-vinyl ether.

## 5. SAFETY

- 5.1. Refer to section 5 of the Method 8000B section of this SOP for general safety requirements.
- 5.2. Often, purge vessels on purge-and-trap instrumentation are pressurized by the time analysis is completed. Therefore, vent the pressure prior to removal of these vessels to prevent the contents from spraying out.
- 5.3. The toxicity or carcinogenicity of each chemical used in this procedure has not been fully defined. Additional health and safety information can be obtained from the MSDS files maintained in the laboratory. The following specific hazards are known:

Methanol -- Flammable and toxic
- 5.4. The following method analytes have been tentatively classified as known or suspected human or mammalian carcinogens: Benzene, Carbon Tetrachloride, 1,4-Dichlorobenzene, 1,2-Dichloroethane, Hexachlorobutadiene, 1,1,2,2-Tetrachloroethane, 1,1,2-Trichloroethane, Chloroform, 1,2-Dibromochloroethane, Tetrachloroethene, Trichloroethene, Vinyl Chloride. Pure standard materials and stock standard solutions of these compounds should be handled in a hood.
- 5.5. Methanol shall not be used in a CaptAir hood.
- 5.6. GC VOA instruments use an ultraviolet (UV) light source, which must be shielded from view. There should also be a warning label/sticker on each instrument that identifies it as a UV light source.
- 5.7. Sodium bisulfate creates Sulfuric Acid when mixed with water.

## 6. EQUIPMENT AND SUPPLIES

- 6.1. Microsyringes -- 10µL, 25µL, 100µL, 250µL, 500µL, and 1000µL. These should be equipped with a 20 gauge (0.006" ID) needle. These will be used to measure and dispense methanolic solutions and aqueous samples.

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A3 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

- 6.2. Gas tight syringes -- 5 mL and 25 mL. Used for measuring sample volumes.
- 6.3. Purge and Trap Apparatus -- A device capable of extracting volatile compounds, trapping on a sorbent trap, and introducing onto a gas chromatograph.
- 6.4. Purge and Trap Autosampler -- In order to maintain high sample throughput, an autosampler is highly recommended.
- 6.5. Trap -- The trap used is dependent on the class of compound to be analyzed. Refer to Table A-2 for suggested traps for specific tests.
- 6.6. Purge Vessels -- These are dependent on the purge and trap unit/autosampler used. Both disposable culture tubes (needle sparge units) and specially designed vessels with fritted bottoms may be used. Follow the manufacturer's suggestions for configuration.
- 6.7. Columns - Refer to Table A-2 for details of columns.
- 6.8. Volumetric flasks, Class A. 5 mL to 250 mL
- 6.9. pH paper
- 6.10. Balance capable of weighing to 0.01g for samples.

## 7. REAGENTS AND SUPPLIES

- 7.1. Refer to the method 8000B section of this SOP for general requirements for reagents and supplies.
- 7.2. Organic Free Water

Organic free water is defined as water in which an interferent is not observed at the reporting limit of the compounds of interest. Suggested methods for generating organic free water include.

  - Filtration through a carbon bed
  - Continuously sparging water with helium or nitrogen.
  - Use of commercial water purification systems.

Other methods may be used, so long as the requirement that the water show no interference is met. The procedure used should be documented in a lab specific attachment.
- 7.3. Sodium Bisulfate
- 7.4. Methanol -- Purge and Trap Grade
- 7.5. Standards

Refer to tables A-5 and A-6 for details of surrogate, matrix spiking and internal standards. Calibration standard levels are not specified, since they may depend on the sensitivity and linear range of specific detectors. However, the low level standard must be equivalent to the reporting limits specified in Table A-1.

  - 7.5.1. Volatile standards are prepared by injecting a measured volume of the stock standard into a syringe containing the appropriate volume of organic free water. The calibration standard is then loaded into the purge device.

**ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**Revision No. 5.7Revision Date: 10/01/03

Page A4 of A17

**8. SAMPLE COLLECTION, PRESERVATION, AND STORAGE**

- 8.1. Holding times for all volatile analysis are 14 days from sample collection.
- 8.2. Water samples are normally preserved at pH < 2 with 1:1 hydrochloric acid. If residual chlorine is present, 2 drops of 10% sodium thiosulfate are added.
- 8.3. Solid samples are field preserved with sodium bisulfate solution for low level analysis, or with methanol for medium level analysis. Soil samples can also be taken using the EnCore™ sampler and preserved in the lab within 48 hours of sampling. At specific client request, unpreserved soil samples may be accepted.
- 8.4. There are several methods of sampling soil. The recommended method, which provides the minimum of field difficulties, is to take an EnCore sample. (The 5 g or 25 g sampler can be used, depending on client preference). Following shipment back to the lab the soil is preserved in methanol. This is the medium level procedure. If very low detection limits are needed (< 50 µg/kg for most analytes) then it will be necessary to use two additional 5 g EnCore samplers or to use field preservation.
- 8.5. Sample collection for medium level analysis using EnCore samplers.
  - 8.5.1. Ship one 5 g (or 25 g) EnCore sampler per field sample position.
  - 8.5.2. An additional bottle must be shipped for percent moisture determination.
  - 8.5.3. When the samples are returned to the lab, extrude the (nominal) 5g (or 25 g) sample into a tared VOA vial containing 5 mL methanol (25 mL methanol for the 25 g sampler). Obtain the weight of the soil added to the vial and note on the label.
  - 8.5.4. Add the correct amount of surrogate spiking mixture. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.)
  - 8.5.5. Add the correct amount of matrix spiking solution to the matrix spike and matrix spike duplicate samples. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.) The addition of spike introduces a slight error, (0.4%) which can be neglected, into the calculations.
  - 8.5.6. Prepare an LCS for each batch by adding the correct amount of matrix spiking solution to clean methanol. (100 µL of spike to 25 mL methanol or 20 µL spike to 5 mL methanol).
  - 8.5.7. Shake the samples for two minutes to distribute the methanol throughout the soil.
  - 8.5.8. Allow to settle, then remove a portion of methanol and store in a clean Teflon capped vial at 4±2°C until analysis.
- 8.6. Sample collection for medium level analysis using field methanol preservation
  - 8.6.1. Prepare a VOA vial by adding 5 mL purge and trap grade methanol. (If a 25 g sample is to be used, add 25 mL methanol to the VOA vial).
  - 8.6.2. Seal the bottle and attach a label.
  - 8.6.3. Weigh the bottle to the nearest 0.01g and note the weight on the label.
  - 8.6.4. Ship with appropriate sampling instructions

**ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**

- 
- 8.6.5. Each sample will require an additional bottle with no preservative for percent moisture determination.
- 8.6.6. At client request, the methanol addition and weighing may also be performed in the field.
- 8.6.7. When the samples are returned to the lab, obtain the weight of the soil added to the vial and note on the label.
- 8.6.8. Add the correct amount of surrogate spiking mixture. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.)
- 8.6.9. Add the correct amount of matrix spiking solution to the matrix spike and matrix spike duplicate samples. (Add 100 µL of 250 µg/mL solution for a nominal 25 g sample, 20µL for a nominal 5 g sample.) The addition of spike introduces a slight error, (0.4%) which can be neglected, into the calculations.
- 8.6.10. Prepare an LCS for each batch by adding the correct amount of matrix spiking solution to clean methanol. (100 µL of spike to 25 mL methanol or 20 µL spike to 5 mL methanol).
- 8.6.11. Shake the samples for two minutes to distribute the methanol throughout the soil.
- 8.6.12. Allow to settle, then remove a portion of methanol and store in a clean Teflon capped vial at 4±2°C until analysis
- 8.7 Low level procedure
- 8.7.1. If low detection limits are required (typically < 50 µg/kg) sodium bisulfate preservation must be used. However, it is also necessary to take a sample for the medium level (methanol preserved) procedure, in case the concentration of analytes in the soil is above the calibration range of the low-level procedure.
- 8.7.2. A purge and trap autosampler capable of sampling from a sealed vial is required for analysis of samples collected using this method. (Varian Archon or O.I 4552).
- 8.7.3. The soil sample is taken using a 5g EnCore sampling device and returned to the lab. It is recommended that two EnCore samplers be used for each field sample position, to allow for any reruns than may be necessary. A separate sample for % moisture determination is also necessary.
- 8.7.4. Prepare VOA vials by adding a magnetic stir bar, approximately 1 g of sodium bisulfate and 5 mL of reagent water.
- 8.7.5. Seal the vial and attach a label. The label must not cover the neck of the vial or the autosampler will malfunction.
- 8.7.6. Weigh the vial to the nearest 0.01g and note the weight on the label.
- 8.7.7. Extrude the soil sample from the EnCore sampler into the prepared VOA vial. Reweigh the vial to obtain the weight of soil and note on the label.
- 8.7.8. **Note:** Soils containing carbonates may effervesce when added to the sodium bisulfate solution. If this is the case at a specific site, add 5 mL of water instead, and freeze at ≥ -10°C until analysis.
- 8.7.9. Alternatively the sodium bisulfate preservation may be performed in the field. Ship at least two vials per sample. The field samplers must determine the weight of soil sampled. Each sample will

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A6 of A17



ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

require an additional bottle with no preservative for percent moisture determination, and an additional bottle preserved with methanol for the medium level procedure.

- 8.8. Aqueous samples are stored in glass containers with Teflon lined septa at 4°C +/- 2°C, with minimum headspace.
- 8.9. Medium level solid extracts are aliquoted into 2 - 5 mL glass vials with Teflon lined caps and stored at 4°C +/- 2°C. The extracts are stored with minimum headspace.
- 8.10. The maximum holding time is 14 days from sampling until the sample is analyzed. (Samples that are found to be unpreserved still have a 14 day holding time. However they should be analyzed as soon as possible. The lack of preservation should be addressed in the case narrative). Maximum holding time for the EnCore sampler (before the sample is added to methanol or sodium bisulfate) is 48 hours.
- 8.11. A holding blank is stored with the samples. This is analyzed and replaced if any of the trip blanks show any contamination. Otherwise it is replaced every 14 days


**Note:** Freezing is not allowed for Ohio VAP solids.



ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**EnCore procedure when low level is not required (field steps  
in gray)****EnCore procedure when low level is required**

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

## Field methanol extraction procedure (field steps in gray)



**ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B**

---

**9. QUALITY CONTROL**

- 9.1. Refer to the method 8000B section of this SOP, section 9, for general quality control procedures, including batch definition, requirements for method blanks, LCS, matrix spikes, surrogates, and control limits.

**10. CALIBRATION AND STANDARDIZATION**

- 10.1. Refer to the method 8000B section of this SOP, section 10, for general calibration procedures.

10.2. Gas Chromatograph Operating Conditions

Various column configurations are possible. If dual column confirmation is necessary, the sample may be split using a Y splitter at the injector end to direct the sample to two columns and two detectors. For simultaneous determination of aromatic and halogenated volatiles, a single column is used and the PID and ELCD detectors are connected in series

- 10.2.1. Refer to Table A-2, A-3 and A-4 for GC operating conditions.

10.3. Initial Calibration

- 10.3.1. Refer to Section 10 of the 8000B section of this SOP for details of initial calibration criteria.

- 10.3.2. Low level soil samples must be purged at 40°C; therefore the calibration curve must also be purged at 40°C. In addition, the low level soil calibration solutions should contain approximately the same amount of sodium bisulfate as the samples.

- 10.3.3. The low level calibration must be at the reporting limit or below. The remaining standards encompass the working range of the detector.

- 10.3.4. Calibrate the instrument using the same volume that will be used during sample analysis.

10.4. Calibration Verification

- 10.4.1. A mid level calibration standard is used for the calibration verification. The gases have 20 % D criteria rather than the 15% used for other analytes.

- 10.4.2. A calibration verification run is performed after every 10 samples for this method.

- 10.4.3. Bracketing of samples with calibration verification runs is only necessary for external standard analysis.

**11. PROCEDURE**

- 11.1. Refer to the method 8000B section of this SOP for general procedural requirements.

11.2. Analytical Sequence

The analytical sequence starts with an initial calibration of at least five points, or a 12 hour calibration that meets % difference criteria from an existing initial calibration.

11.3. Confirmation

The PID and ELCD detectors are sufficiently selective that second column confirmation is not always necessary. Requirements for second column confirmation should be decided in consultation with the

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A 10 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

client. If the PID and ELCD are used in series confirmatory information for many analytes can be gained by comparing the relative response from the two detectors.

## 11.4. Aqueous Sample Analysis (Purge and Trap units using sparge vessels)

- 11.4.1. Depending on the sensitivity of the instrument and capabilities of the purge and trap device, 5, 10, 20, or 25 mL sample volumes may be analyzed. A 5 mL sample volume is recommended.
- 11.4.2. Rinse a 5 mL (or 25 mL for larger sample volumes) gas-tight syringe with organic free water. Fill the syringe with the sample to be analyzed, and compress to volume.
- 11.4.3. Check and document the pH of the sample remaining in the VOA vial after loading the syringe.
- 11.4.4. This procedure invalidates the contents of the VOA vial for further analysis, unless an aliquot is transferred to a smaller VOA vial with no headspace (e.g., 20 mL) at the same time the analysis aliquot is removed.
- 11.4.5. Spike with the appropriate volume of surrogate/internal standard solution and spike solution (if required) through the barrel of the syringe. The method blank is spiked with surrogates only, the LCS and matrix spikes with the surrogate and matrix spiking solutions. Refer to Tables A-5 and A-6 for volumes and concentrations of spiking solutions.
- 11.4.6. Load onto the purge and trap device and start the run.
- 11.4.7. If the initial analysis of a sample or a dilution of the sample has a concentration of analytes that exceeds the initial calibration range, the sample must be reanalyzed at a higher dilution. When a sample has a high response for a compound, an organic free water blank should follow analysis. It is recognized that during automated unattended analysis, this may not occur. If any potential carryover hits are present in samples following highly contaminated samples, the sample must be reanalyzed to determine if any of these hits are a result of carryover or are actually present in the sample.
- 11.4.8. Dilutions may be made in gas tight syringes unless the volume of sample used is less than 5  $\mu$ L, in which case dilution in volumetric flasks will be necessary.
  - 11.4.8.1. Spike with the same volume of surrogate/internal standard solution as used for undiluted samples prior to loading onto the purge and trap device.
  - 11.4.8.2. For Matrix spike / matrix spike duplicates where the sample requires dilution, the sample is spiked after the dilution is performed.

## 11.5. Aqueous and Soil Sample Analysis (Purge and Trap units that sample directly from the VOA vial)

- 11.5.1. Units, which sample from the VOA vial, should be equipped with a module, which automatically adds surrogate and internal standard solution to the sample prior to purging the sample.
- 11.5.2. If the autosampler uses automatic IS/SS injection, no further preparation of the VOA vial is needed. Otherwise the internal and surrogate standards must be added to the vial. *Note:* Aqueous samples with high amounts of sediment present in the vial may not be suitable for analysis on this instrumentation, or they may need to be analyzed as soils.
- 11.5.3. Sample remaining in the vial after sampling with one of these mechanisms is no longer valid for further analysis. A fresh VOA vial must be used for further sample analysis.
- 11.5.4. Check the pH of the sample remaining in the VOA vial after analysis is completed.

# ANALYSIS OF VOLATILE ORGANICS BASED ON METHOD 8021B

## 11.6. Low-Level Solids Analysis using discrete autosamplers

**Note:** This technique may seriously underestimate analyte concentration and must not be used except at specific client request for the purpose of comparability with previous data. It is no longer part of SW-846.

*This method is based on purging a heated sediment/soil sample mixed with reagent water containing the surrogate and, if applicable, internal and matrix spiking standards. Analyze all reagent blanks and standards under the same conditions as the samples (e.g., heated). The calibration curve is also heated during analysis. Purge temperature is 40°C.*

11.6.1. Do not discard any supernatant liquids. Mix the contents of the container with a narrow metal spatula.

11.6.2. Weigh out 5 g (or other appropriate aliquot) of sample into a disposable culture tube or other purge vessel. Record the weight to the nearest 0.1 g. If method sensitivity is demonstrated, a smaller aliquot may be used. Do not use aliquots less than 1.0 g. If the sample is contaminated with analytes such that a purge amount less than 1.0 g is appropriate, use the medium level method described in section 11.7.

11.6.3. Connect the purge vessel to the purge and trap device.

11.6.4. Rinse a 5 mL gas-tight syringe with organic free water, and fill. Compress to 5 mL. Add surrogate/internal standard (and matrix spike solutions if required) (See Tables A-5, A-6, A-7 and A-8). Add directly to the sample from 11.6.2.

11.6.5. The above steps should be performed rapidly and without interruption to avoid loss of volatile organics.

11.6.6. Add the heater jacket or other heating device and start the purge and trap unit.

11.6.7. Soil samples that have low IS recovery when analyzed (<50%) should be reanalyzed once to confirm matrix effect. If external standard calibration is used, samples with surrogate recovery below the control limit should be reanalyzed once to confirm matrix effect.

## 11.7. Methanol Extract Soils

11.7.1. Rinse a gas-tight syringe with organic free water. Fill the syringe with the same volume of organic free water as used in the calibrations. Add no more than 2% (v/v) (100 µL for a 5 mL purge) methanolic extract (from Section 8.5 or 8.6) to the syringe. Add internal standard (if used). Load the sample onto the purge and trap device and analyze as for aqueous samples. If less than 5 µL of methanolic extract is to be added to the water, dilute the methanolic extract such that a volume greater than 5 µL will be added to the water in the syringe.

## 12. DATA ANALYSIS AND CALCULATIONS

Refer to section 12 of the 8000B section of this SOP.

## 13. METHOD PERFORMANCE

13.1. Performance limits for the four replicate initial demonstration of capability required under Section 13.1 of the 8000B section of this SOP.

## 14. POLLUTION PREVENTION

This method does not contain any specific modifications that serve to minimize or prevent pollution.

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7Revision Date: 10/01/03

Page A12 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

## 15. WASTE MANAGEMENT

- 15.1. All waste will be disposed of in accordance with Federal, State and Local regulations. Where reasonably feasible, technological changes have been implemented to minimize the potential for pollution of the environment. Employees will abide by this method and the policies in section 13 of the Corporate Safety Manual for "Waste Management and Pollution Prevention."
- 15.2. Waste streams produced by the method.
  - 15.2.1. The following waste streams are produced when this method is carried out.
    - 15.2.1.1. **Acidic material from the auto-sampler.** Waste stream must be collected and neutralized before discharge to a sewer system if the pH is less than 4.
    - 15.2.1.2. **Methanol waste from rinses and standards.** Methanol waste is discarded as a flammable liquid.
    - 15.2.1.3. **All samples including purged and extracted soils and waters:** Samples are collected in boxes and removed from the lab to storage. The waste coordinator handles crushing the vials and proper disposal.

## 16. REFERENCES

- 16.1. Test Methods for Evaluating Solid Waste, Physical/Chemical Methods, SW846, 3rd Edition, Final Update III, December 1996, Sections 5000, 5030B, 5035 and 8021B

## 17. MISCELLANEOUS

- 17.1. Modifications from Reference Method
- 17.2. Modifications from previous revision
  - 17.2.1. No revisions were made to this appendix.
- 17.3. Facility Specific SOPs

Each facility shall attach a list of facility specific SOPs or approved attachments (if applicable) which are required to implement this SOP or which are used in conjunction with this SOP. If no facility specific SOPs or amendments are to be attached, a statement must be attached specifying that there are none.

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7

Revision Date: 10/01/03

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

Page A13 of A17

## 17.4. Tables

Table A-1					
Standard Analyte List					
Test	Compound	CAS number	Reporting Limit, µg/L or µg/kg		
			Aqueous	Low Soil	Medium Soil
Halogenated volatiles by 8021B	Bromodichloromethane	75-27-4	1.0	1.0	50
	Bromoform	75-25-2	1.0	1.0	50
	Bromomethane	74-83-9	1.0	1.0	50
	Carbon Tetrachloride	56-23-5	1.0	1.0	50
	Chlorobenzene	108-90-7	1.0	1.0	50
	Chloroethane	70-00-3	1.0	1.0	50
	2-Chloroethyl vinyl ether	110-75-8	5.0	5.0	250
	Chloroform	67-66-3	1.0	1.0	50
	Chloromethane	74-87-3	1.0	1.0	50
	Dibromochloromethane	124-48-1	1.0	1.0	50
	1,2-Dichlorobenzene	95-50-1	1.0	1.0	50
	1,3-Dichlorobenzene	541-73-1	1.0	1.0	50
	1,4-Dichlorobenzene	106-46-7	1.0	1.0	50
	Dichlorodifluoromethane	75-71-8	1.0	1.0	50
	1,1-Dichloroethane	75-34-3	1.0	1.0	50
	1,2-Dichloroethane	107-06-2	1.0	1.0	50
	1,1-Dichloroethene	75-45-4	1.0	1.0	50
	cis-1,2 Dichloroethene	156-59-4	1.0	1.0	50
	trans-1,2-Dichloroethene	156-60-5	1.0	1.0	50
	Dichloromethane(DCM)	75-09-2	5.0	5.0	250
	1,2-Dichloropropane	78-87-5	1.0	1.0	50
	cis-1,3-Dichloropropene	10061-01-5	1.0	1.0	50
	trans-1,3-Dichloropropene	10061-02-6	1.0	1.0	50
	1,1,2,2-Tetrachloroethane	79-34-5	1.0	1.0	50
	Tetrachloroethene	127-18-4	1.0	1.0	50
	1,1,1-Trichloroethane	71-55-6	1.0	1.0	50
	1,1,2-Trichloroethane	79-00-5	1.0	1.0	50
	Trichloroethene	79-01-6	1.0	1.0	50
	Trichlorofluoromethane	75-69-4	1.0	1.0	50
	Vinyl Chloride	75-01-4	1.0	1.0	50
Additional halogenated volatiles	Benzyl Chloride	100-44-7	5.0	5.0	250
	Bromobenzene	108-86-1	1.0	1.0	50
	Dibromomethane	74-95-3	1.0	1.0	50
	1,1,1,2-Tetrachloroethane	630-20-6	1.0	1.0	50
	Freon 113	76-13-1	1.0	1.0	50
	1,2,3-Trichloropropane	96-18-4	1.0	1.0	50
BTEX by 8021B	Benzene	71-43-2	1.0	1.0	50
	Ethyl Benzene	100-41-4	1.0	1.0	50
	Toluene	108-88-3	1.0	1.0	50
	Xylenes (total)	1330-20-7	1.0	1.0	50
Aromatic	Benzene	71-43-2	1.0	1.0	50

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 57

Revision Date: 10/01/03

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

Page A14 of A17

Table A-1					
Standard Analyte List					
Test	Compound	CAS number	Reporting Limit, µg/L or µg/kg		
			Aqueous	Low Soil	Medium Soil
volatiles by 8021B					
	Chlorobenzene	108-90-7	1.0	1.0	50
	1,2-Dichlorobenzene	75-34-3	1.0	1.0	50
	1,3-Dichlorobenzene	107-06-2	1.0	1.0	50
	1,4-Dichlorobenzene	75-45-4	1.0	1.0	50
	Ethyl Benzene	100-41-4	1.0	1.0	50
	Toluene	108-88-3	1.0	1.0	50
	Xylenes (total)	1330-20-7	1.0	1.0	50
Additional aromatic and unsaturated volatiles	1,2,4 Trimethylbenzene	95-63-6	1.0	1.0	50
	1,3,5 Trimethylbenzene	108-67-8	1.0	1.0	50
	Acetone	67-64-1	10	10	500
	MEK (2-butanone)	78-93-3	5.0	5.0	250
	MIBK (4-methyl-2-pentanone)	108-10-1	5.0	5.0	250
	Naphthalene	91-20-3	2.0	2.0	250
	Styrene	100-42-5	1.0	1.0	50
	Methyl tert-butyl ether (MTBE)	1634-04-4	1.0	1.0	50
Combined halogenated and aromatic volatiles by 8021B	Benzene	71-43-2	1.0	1.0	50
	Bromobenzene	108-86-1	1.0	1.0	50
	Bromochloromethane	74-97-5	1.0	1.0	50
	Bromodichloromethane	75-27-4	1.0	1.0	50
	Bromoform	75-25-2	1.0	1.0	50
	Bromomethane	74-83-9	1.0	1.0	50
	n-butylbenzene	104-51-8	1.0	1.0	50
	sec-Butylbenzene	135-98-8	1.0	1.0	50
	tert-Butylbenzene	98-06-6	1.0	1.0	50
	Carbon Tetrachloride	56-23-5	1.0	1.0	50
	Chlorobenzene	108-90-7	1.0	1.0	50
	Chlorodibromomethane	124-48-1	1.0	1.0	50
	Chloroethane	75-00-3	1.0	1.0	50
	Chloroform	67-66-3	1.0	1.0	50
	Chloromethane	74-87-3	1.0	1.0	50
	2-Chlorotoluene	95-49-8	1.0	1.0	50
	4-Chlorotoluene	106-43-4	1.0	1.0	50
	1,2-Dibromo-3-Chloropropane(DBCP)	96-12-8	1.0	1.0	50
	1,2-Dibromoethane(EDB)	106-93-4	1.0	1.0	50
	Dibromomethane	74-95-3	1.0	1.0	50
	1,2-Dichlorobenzene	95-50-1	1.0	1.0	50
	1,3-Dichlorobenzene	541-73-1	1.0	1.0	50
	1,4-Dichlorobenzene	106-46-7	1.0	1.0	50
	Dichlorodifluoromethane	75-71-8	1.0	1.0	50
	1,1-Dichloroethane	75-34-3	1.0	1.0	50
	1,2-Dichloroethane	107-06-2	1.0	1.0	50
	1,1-Dichloroethene	75-35-4	1.0	1.0	50



8771402

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7

## ANALYSIS OF VOLATILE ORGANICS BASED ON

Revision Date: 10/01/03

## METHOD 8021B

Page A15 of A17

Table A-1					
Standard Analyte List					
Test	Compound	CAS number	Reporting Limit, µg/L or µg/kg		
			Aqueous	Low Soil	Medium Soil
	cis-1,2-Dichloroethene	156-59-4	1.0	1.0	50
	trans-1,2-Dichloroethene	156-60-5	1.0	1.0	50
	1,2-Dichloropropane	78-87-5	1.0	1.0	50
	1,3-Dichloropropane	142-28-9	1.0	1.0	50
	2,2-Dichloropropane	590-20-7	1.0	1.0	50
	1,1-Dichloropropene	563-58-6	1.0	1.0	50
	cis-1,3-Dichloropropene	10061-01-5	1.0	1.0	50
	trans-1,3-Dichloropropene	10061-02-6	1.0	1.0	50
	Ethylbenzene	100-41-4	1.0	1.0	50
	Hexachlorobutadiene	87-68-3	1.0	1.0	50
	Isopropylbenzene	98-82-8	1.0	1.0	50
	p-Isopropyltoluene	99-87-6	1.0	1.0	50
	Methylene Chloride	75-09-2	5.0	5.0	250
	Naphthalene	91-20-3	2.0	2.0	250
	n-Propylbenzene	10306501	1.0	1.0	50
	Styrene	100-42-5	1.0	1.0	50
	1,1,1,2-Tetrachloroethane	630-20-6	1.0	1.0	50
	1,1,2,2-Tetrachloroethane	79-34-5	1.0	1.0	50
	Tetrachloroethene	127-18-4	1.0	1.0	50
	Toluene	108-88-3	1.0	1.0	50
	1,2,3-Trichlorobenzene	87-61-6	1.0	1.0	50
	1,2,4-Trichlorobenzene	120-82-1	1.0	1.0	50
	1,1,1-Trichloroethane	71-55-6	1.0	1.0	50
	1,1,2-Trichloroethane	79-00-5	1.0	1.0	50
	Trichloroethene	79-01-6	1.0	1.0	50
	Trichlorofluoromethane	75-69-4	1.0	1.0	50
	1,2,3-Trichloropropane	96-18-4	1.0	1.0	50
	1,2,4-Trimethylbenzene	95-63-6	1.0	1.0	50
	1,3,5-Trimethylbenzene	108-67-8	1.0	1.0	50
	Vinyl Chloride	75-01-4	1.0	1.0	50
	Xylenes (total)	1330-20-7	1.0	1.0	50

Table A-2	
Recommended Conditions for Aromatic Volatiles	
Parameter	Recommended Conditions
Temperature program	50°C, 1min, 10°C/min to 200°C, 1min
Column 1	Rtx-502.2 or DB-502.2 60m x 0.53mm 3.0µm
Column 2	Rtx-1 or DB-1 60m x 0.53mm 3.0 µm
Carrier gas	Helium or hydrogen
Purge Flow / time	40 mL/min, 11 minutes
Desorb Temp / time	180°C, 2 minutes (220°C for Vocab 3000)
Bake Time / temp	200°C, 12 minutes (230°C for Vocab 3000)
Transfer line / valve temp	115°C

## APPENDIX A

SOP No. CORP-GC-0001NC

Revision No. 5.7

Revision Date: 10/01/03

Page A16 of A17

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

Table A-3 Recommended Conditions for Method Halogenated Volatiles	
Parameter	Recommended Conditions
Temperature program	35°C, 12 min, then 4°C/min to 200°C, hold for 5 min
Column 1	DB-VRX or RTX-502.2 105m x 0.53 mm id df = 3.0um
Column 2	DB-1 or RTX-1 105m x 0.53 mm ID df = 3.0um
Column 3	Rtx - Volatiles 120m x 0.53mm ID df=2.0um
Carrier gas	Helium
Purge Flow / time	40 mL/min, 11 minutes
Desorb Temp / time	180°C, 2 minutes (220°C for Vocab 3000)
Bake Time / temp	200°C, 12 minutes (230°C for Vocab 3000)
Transfer line / valve temp	115°C

Table A-4 Recommended Conditions for Method Combined Aromatic and Halogenated Volatiles	
Parameter	Recommended Conditions
Temperature program	35°C, 12 min, then 4°C/min to 200°C, hold for 5 min
Column 1	DB-VRX or RTX-502.2 105m x 0.53 mm id df = 3.0um
Column 2	DB-1 or RTX-1 105m x 0.53 mm ID df = 3.0um
Column 3	Rtx - Volatiles 120m x 0.53mm ID df=2.0um
Carrier gas	Helium
Purge Flow / time	40 mL/min, 11 minutes
Desorb Temp / time	180°C, 2 minutes (220°C for Vocab 3000)
Bake Time / temp	200°C, 12 minutes (230°C for Vocab 3000)
Transfer line / valve temp	115°C

## APPENDIX A

ANALYSIS OF VOLATILE ORGANICS BASED ON  
METHOD 8021B

SOP No. CORP-GC-0001NC

Revision No. 57

Revision Date: 10/01/03

Page A17 of A17

Table A-5 Surrogate and Internal Standard Concentrations				
Standard	Components	Working Solution µg/mL	Spike amount µL (for 5 mL purge)	Final concentration µg/L (µg/kg)
Aromatic volatiles IS/SS	4-Chlorotoluene (SS)	20	5	20
	1-Chloro-4-fluorobenzene (IS)	40		40
Halogenated volatiles IS/SS	4-chlorotoluene (SS)	20	5	20
	1-Chloro-4-fluorobenzene (IS)	40		40
Combined Aromatic and halogenated volatiles IS/SS	Fluorobenzene (SS)	20	5	20
	1,4-Dichlorobutane (SS)	20	5	20
	1-Chloro-4-fluorobenzene (IS)	40	10	40

It may be necessary to select different surrogates in order to minimize sample interferences. 1-chloro-4-fluorobenzene and 4-chlorotoluene are fairly well resolved from analytes listed in this SOP. However 4-chlorotoluene may sometimes be requested as a target analyte. Other surrogates that may be considered, and issues associated with their use are:

Bromochloromethane:

Elutes very close to chloroform and cis-1, 2-dichloroethene on the 502.2 column.  
May be a target analyte.

1,2-Bromochloroethane:

1-Chloro-2-fluorobenzene.

Elutes close to ethylbenzene on DB-1 or Rtx-1 and close to m,p-xylene on 502.2

a,a,a-Trifluorotoluene:

Good for aromatic volatiles, coelutes or very close to trichloroethene

Bromofluorobenzene:

Close to 1,1,2,2-trichloroethane and 1,2,3-trichloropropane on the 502.2 column.

Good on DB-1 or Rtx-1.

2-Bromo-1-chloropropane.

May coelute with 1,1,2-trichloroethane

Table A-6 Concentrations for LCS and MS/MSD compounds				
Standard	Components	Working Solution µg/mL	Spike amount µL (5 mL purge)	Final concentration µg/L (µg/kg)
Aromatic	Benzene	20	5	20
	Toluene	20		20
	Chlorobenzene	20		20
Halogenated	Chlorobenzene	20	5	20
	1,1-Dichloroethene	20		20
	Trichloroethene	20		20
Combination aromatic / halogenated	Benzene	20	5	20
	Toluene	20		20
	Chlorobenzene	20		20
	1,1-Dichloroethene	20		20
	Trichloroethene	20		20

**FINAL PAGE**

**PART II**

**ADMINISTRATIVE RECORD**

**PART II**

**FINAL PAGE**